# PRIME SURGICAL INTERNATIONAL LAPROSCOPIC INSTRUMENTS



# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### **Table of Contents**

	INTRODUCTION	4
	History	4
	Surgical Quality Management	7
.596	INSTRUMENTS	11
151	Scissors, Graspers, Dissectors	12-27
28	Modular System	12-23
8	Reposable Scissors	27
	Monopolar Electrodes	28
	Suturing	29
	Bipolar Instruments	30-31
	Specialty-Locking Graspers	32
	Specialty Instrument Sets - Advanced Procedure	33-36
	CLIP APPLIERS	37-40
de	SLS Single-Fire Ligation System	38
1	Challenger Ti-P Multi-Fire Ligation System	39
-	Vascular Clip Applier/Remover	40
1	ENDOSCOPES AND ACCESSORIES	41-43
47.17	HD Endoscopes	42
111.00	Unitrac® Retraction & Holding System	43-44
	TROCARS	45-51
43th	Reusable Trocar Systems	46-50
-89	Flexible Trocars	51
•		
-25-17	VIDEO SYSTEMS	53-55
<b>经</b>	Light Sources	54
H-A	Light Cables	55
	MISCELLANEOUS	57-62
The same	IVIISCELLAINEUUS	57-02
0		
1	GYNECOLOGICAL	_
		63-72
	PROCEDURE	
	INDEX	73-78

### History

PRIME SURGICAL INTERNATIONAL GENERAL SURGICAL INSTRUMETNS prime.instruments@hotmail.com

PRIME SURGICAL INTERNATIONAL GENERAL SURGICAL INSTRUMETNS prime.instruments@hotmail.com

PRIME SURGICAL INTERNATIONAL GYNECOLOGIC SURGICAL INSTRUMETNS prime.instruments@hotmail.com

PRIME SURGICAL INTERNATIONAL GYNECOLOGIC SURGICAL INSTRUMETNS prime.instruments@hotmail.com

PRIME SURGICAL INTERNATIONAL GYNECOLOGIC SURGICAL INSTRUMETNS prime.instruments@hotmail.com

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com




#### PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com



#### QUICK-SNAP™ Laparoscopic Instruments

line of instruments for minimally invasive surgery combines precision craftsmanship with modern materials and design.

Optional Handles for Customizing Instruments - Interchangeable for 5 mm and 10 mm instruments











PSI-L I -0001(Non-Ratchet)

PSI-L I -0002 (Ratchet)


PSI-L I -0003 (Non-Ratchet)


PSI-L I -0004 (Ratchet)


**PSI-L I -0005** (In-Line)

			COMPONENTS			
	COMPLETE	JAW INSERT	INSULATED	STANDARD	OPTIONAL	
	INSTRUMENT		OUTER TUBES	HANDLE	HANDLES	
1/1	PSI-L I -0006					
	Hook Scissors, 5 mm, 31 cm	PSI-L I -0008	PSI-L I -0009			
	PSI-L I -0007					
	Hook Scissors, 5 mm, 42 cm	PSI-L I -0010	PSI-L I -0011			
1/1	PSI-L I -0012	PSI-L I -00945	PSI-L I -00946			
	Micro-Dissecting Scissors, 5 mm, 31 cm					
	PSI-LI-00965 Order by component	PSI-L I -00947	PSI-L I -00948			
	Micro-Dissecting Scissors, 5 mm, 42 cm					
1/1	PSI-L I -0013	PSI-L I -00949				
	Metzenbaum Scissors, serrated,	PS1-L 1 -00949	PSI-L I -00950			
	Tungsten Carbide inserts, 5 mm, 31 cm					
6	PSI-L I -0014	PSI-L I -00951	PSI-L I -00952			
	Metzenbaum Scissors, serrated,					
	Tungsten Carbide inserts, 5 mm, 42 cm					
1/1	- 100 000 000	_				
	PSI-L I -0015	PSI-L I -00953	PSI-L I -00954			
0	Metzenbaum Scissors, insulated to tip, Tungsten Carbide Inserts, 5 mm, 31 cm					
	rangeen caronac macra, o min, or em			_	_	

www.primeinstrumentintl.com, prime.instruments@hotmail.com +92 300 2621673 , +92 313 5710887.








			COMPONE			
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES	
1/1	PST-LI-0971  CROCE/OLMI Dissecting and Suture Scissors with Knot Pusher, 5 mm, 31 cm	PSI-LI-0972	PSI-LI-0973			
1/1	PSI-LI-0974  Durogrip TC Needle Holder, straight, 5mm, 31cm	PSI-LI-0975	PSI-LI-0976		•	
1/1	PSI-LI-0977  Durogrip TC Needle Holder, curved, 5mm, 31cm	PSI-LI-0978	PSI-LI-097s		**	
1/1	PSI-LI-0980 Atraumatic Fixation Forceps, serrated, 5 mm, 31 cm	PSI-LI-0981	PSI-LI-	0982	**	
	Atraumatic Fixation Forceps, serrated, 5 mm, 42 cm	PSI-LI-0984	PSI-LI-	0985	**	
1/1	PSI-LI-0986 Tubal Fixation Forceps, 5 mm, 31 cm	PSI-LI-0987	PSI-L	1-0988		
	PSI-LI-0989  Tubal Fixation Forceps, 5 mm, 42 cm	PSI-LI-0990	PSI-I	.I-0991		
1/1	Oval Grasper, 5 mm, 31 cm	PSI-LI-0993	PSI	-LI-0994		
3.3	PSI-LI-0995 Oval Grasper, 5 mm, 42 cm	PSI-LI-0996	PSI-	LI-0997		

### www.primeinstrumentintl.com

QUICK-SNAP™ design









Axial design

PSI-L I -00 73 (Non-Ratchet)

PSI-L I -00 74 (Ratchet)

PSI-LI-00 75 (Non-Ratchet)


PSI-LI-00 76 (Ratchet)

PSI-L I -00 77 (In-Line)

		COMPONENTS			
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	PSI-LI-00 78  Curved Grasper, 5 mm, 31 cm	PSI-L I -00 79	PSI-L I -00 80		
(IIIIIIIIII	Curved Grasper, 5 mm, 42 cm	PSI-L I -00 82	PSI-L I -00 83		
1/1	PSI-LI-00 84 Fixation Forceps, fenestrated, 5 mm, 31 cm	PSI-L I -00 85	PSI-L I -00 86		
	PSI-LI-00 87 PSI-LI-0088  Fixation Forceps, fenestrated, serrated, 5 mm, 31 cm	PSI-LI-0998	PSI-LI-0999		
1/1	Atraumatic Dissecting Forceps, fine tip, 5 mm, 31 cm	PSI-L I -0090	PSI-L I -0091		
1/1	Maryland Dissecting Forceps, 5 mm, 31 cm	PSI-L I -0093	PSI-L I -0094		
The state of the s	Maryland Dissecting Forceps, 5 mm, 42 cm	PSI-L I -0096	PSI-L I -0097		
1/1	Maryland Dissecting Forceps, insulated to tip, 5 mm, 31 cm	PSI-L I -0099	PSI-L I -00100		::
1/1	PSI-L I -00101  Crile Dissecting Forceps, Curved, 5 mm, 31 cm	PSI-L I -00102	PSI-L I -00103		
	PSI-L I -00104  Crile Dissecting Forceps, Curved, 5 mm, 37 cm	PSI-L I -00105	PSI-L I -00106		•

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com









PSI-LI-00107 (Non-Ratchet)

PSI-LI-00108 (Ratchet)

PSI-L I -00109 (Non-Ratchet)

PSI-LI-00110 (Ratchet)

PSI-LI-00111 (In-Line)


			COMPONE	ENTS	
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	Overholt 60° Dissecting Forceps, 10 mm, 31 cm	PSI-L I -00113	PSI-L I -00114	÷	::
1/1	PSI-LI-00115 Right Angle, 90° Dissecting Forceps, 10 mm, 31 cm	PSI-L I -00116	PSI-L I -00117		::
<i>y</i> <sub>4</sub>	Mixter Dissector, Right Angle, 90° 5 mm, 31 cm	PSI-L I -00119	PSI-L I -00120		•
1/1	Right Angle, 90° Dissecting Forceps, 5 mm, 31 cm	PSI-L I -00122	PSI-L I -00123	•	
1/1	PSI-LI-00124  Dolphin Nose Dissecting Forceps, 5 mm, 31 cm	PSI-L I -00125	PSI-L I -00126		::
1/1	PSI-LI-00127  Duck Bill Dissecting Forceps, 5 mm, 31 cm	PSI-L I -00128	PSI-L I -00129		

## www.primeinstrumentintl.com



1

QUICK-SNAP™ design







PSI-L I -00130 (Non-Ratchet)

PSI-L I -00131 (Ratchet)


PSI-L I -00132 (Non-Ratchet)

PSI-L I -00133 (Ratchet)

**PSI-LI-01000** (In-Line)


		COMPONENTS			
	COMPLETE INSTRUMENT	JAW INSERT		STANDARD HANDLE	OPTIONAL HANDLES
1/1	PSI-L I -00134	PSI-L I -00135	PSI-L I -00136		
500	Duck Bill Blunt Grasper,				
\$-0.00m	5 mm, 37 cm				
1/1	PSI-L I -00137	PSI-L I -00138	PSI-L I -00139		
	Bullet Nose Grasper, 5 mm, 31 cm				
William Co.	PSI-L I -00140	PSI-L I -00141	PSI-L I -00142		
	Bullet Nose Grasper, 5 mm, 42 cm				
1/1	PSI-L I -00143	PSI-L I -00144	PSI-L I -00145		
	Aggressive Grasping Forceps,				
1110	5 mm, 31 cm				
1/1	PSI-L I -00146	PSI-L I -00147	PSI-L I -00148		.,1
	Alligator Grasping Forceps,				
00/00/00/00	5 mm, 31 cm				
1/1	PSI-L I -00149	PSI-L I -00150	PSI-L I -00151		
	Bowel Grasping Forceps, 20 mm Jaw,				
	5 mm, 31 cm				
	PSI-L I -00152				
	Bowel Grasping Forceps, 20 mm Jaw,	PSI-L I -00153	PSI-L I -00154		
	5 mm, 42 cm				
1/1		PSI-L I -00156	PSI-L I -00157		
., .	PSI-L I -00155	75555			
	Dorsey Intestinal Forceps, 40 mm Jaw, 5 mm, 31 cm				
	5 mm, 31 cm	PSI-L I -00159	PSI-L I -00160		
dinking and the same and the same		F51-L 1 -00159			
	Dorsey Intestinal Forceps, 40 mm Jaw, 5 mm, 42 cm				
	o min, T2 Cili				

## PRIME SURGICAL INTERNATIONAL













PSI-LI-00161 (Non-Ratchet)

PSI-LI-00162 (Ratchet)

PSI-L I -00163 (Non-Ratchet)


PSI-L I -00164 (Ratchet)

PSI-LI-010001 (In-Line)

			COMPONE	ENTS	
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	Universal Grasper, Fenestrated 5 mm, 31 cm	PSI-L I -00166	PSI-L I -00167	÷	::
1/2	PSI-L I -00168  Glassman Clamping Forceps, 60 mm Jaw, 10 mm, 31 cm	PSI-L I -00169	PSI-L I -00170	÷	::
1/2	PSI-LI-00171  DeBakey Atraumatic Forceps, 50 mm Jaw, 10 mm, 37 cm	PSI-L I -00172	PSI-L I -00173		::
1/1	DeBakey Atraumatic Forceps, 50 mm Jaw, 5 mm, 42 cm	PSI-L I -00175	PSI-L I -00176	÷	::
1/1	PSI-LI-00177  Babcock Grasping Forceps,  Single Action, 5 mm, 31 cm	PSI-L I -00178	PSI-L I -00179		::
1/1	PSI-L I -00180  Babcock Grasping Forceps  Double Action, 5 mm, 31 cm	PSI-L I -00181	PSI-L I -00182 PSI-L I -00185		
	Babcock Grasping Forceps, Double Action, 5 mm, 42 cm	PSI-L I -00184	, 54 24 -40493		•

#### www.primeinstrumentintl.com











PSI-L I -00186 (Non-Ratchet)

PSI-L I -00187 (Ratchet)

PSI-LI-00188 (Non-Ratchet)

PSI-LI-00189 (Ratchet)

**PSI-LI-010002** (In-Line)

			COMPONE	ENTS	
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	PSI-LI-00190 Babcock Grasping Forceps, Single Action, 10 mm, 31 cm	PSI-L I -00191	PSI-L I -00192		
10	Babcock Grasping Forceps, Single Action, 10 mm, 42 cm	PSI-L I -00194	PSI-L I -00195		
1/2	Babcock Grasping Forceps, Double Action, 10 mm, 37 cm	PSI-L I -00197	PSI-L I -00198		•
1/1	PSI-L I -00199  Claw Grasper, 2 x 3 teeth, 10 mm, 31 cm	PSI-L I -00200	PSI-L I -00201		::
1/1	Allis Atraumatic Grasping Forceps, 5 mm, 31 cm	PSI-L I -00203	PSI-L I -00204		
1/1	PSI-L I -00205 Allis Grasping Forceps, 5 mm, 31 cm	PSI-L I -00206	PSI-L I -00207		**
1/1	Universal Grasping Forceps, Fenestrated 5 mm, 31 cm	PSI-L I -00209 PS	I-L I -00210	•	
	Universal Grasping Forceps, Fenestrated 5 mm, 42 cm	PSI-L I -00212	PSI-L I -00213		

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com












Instruments may be ordered COMPLETE or assembled from COMPONENTS ordered separately

			COMPONE	ENTS	
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	Grasping Forceps, 2 x 4 teeth, serrated, 5 mm, 31 cm	PSI-L I -00219	PSI-L I -00220		•
Sterion S	Grasping Forceps, 2 x 4 teeth, serrated, 5 mm, 42 cm	731 27 3322			
1/2	Grasping Forceps, fine-toothed, non-serrated, 5 mm, 31 cm	PSI-L I -00225	PSI-L I -00226		•
1/1	PSI-LI-00227  Heywood-Smith Tissue Forceps, fenestrated, 10 mm, 31 cm	PSI-L I -00228	PSI-L I -00229	÷	::
1/1	Atraumatic Grasping Forceps, fenestrated, 5 mm, 31 cm	PSI-L I -00231	PSI-L I -00232		
1/1	Fundus Grasper, 5 mm, 31 cm,	PSI-L I -00234	PSI-L I -00235	•	
1/1	Fundus Grasper, 5 mm, 42 cm,  PSI-LI-00239  Short Fenestrated Grasper, 5 mm, 31 cm	PSI-L I -00240	PSI-L I -00241		-
	Short Fenestrated Grasper, 5 mm, 42 cm	PSI-L I -00243	PSI-L I -00244		

### www.primeinstrumentintl.com











PSI-L I -00245 (Non-Ratchet)

PSI-L I -00246 (Ratchet)

PSI-L I -00247 (Non-Ratchet)

PSI-L I -00248 (Ratchet)

**PSI-LI-010004** (In-Line)


		COMPONENTS			
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	PSI-L I -00249 Tenaculum Forceps 5 mm, 31 cm	PSI-L I -00250	PSI-L I -00251		
Comme	PSI-LI-00252 Tenaculum Forceps 10 mm, 31 cm	PSI-L I -00253	PSI-L I -00254		
1/1	PSI-L I -00255	PSI-L I -00256	PSI-L I -00257		
	Large Biopsy Spoon Forceps, 10 mm, 31 cm				•
1/1	PST-L I -00258 Biopsy Forceps, Cup, fenestrated, 5 mm, 31 cm	PSI-L I -00259	PSI-L I -00260		
	PSI-LI-00261 Biopsy Forceps, Cup, fenestrated, 5 mm, 42 cm	PSI-L I -00262	PSI-L I -00263		
1/1	PSI-LI-00264 Biopsy Forceps, Box, 5 mm, 31 cm	PSI-L I -00265	PSI-L I -00266		
	PSI-LI-00267 Biopsy Forceps, Box, 5 mm, 42 cm	PSI-L I -00268	PSI-L I -00269		
1/1	PSI-L I -00270 Biopsy Forceps, with teeth, 5 mm, 31 cm	PSI-L I -00271	PSI-L I -00272		
	PSI-LI-00273 Biopsy Forceps, with teeth, 5 mm, 42 cm	PSI-L I -00274	PSI-L I -00275		
1/1	PSI-LI-00276 Biopsy Forceps, Hook, 5 mm, 42 cm	PSI-LI-010005	PSI-LI-010006	•	

#### PRIME SURGICAL INTERNATIONAL

#### www.primeinstrumentintl.com











Instruments may be ordered COMPLETE or assembled from COMPONENTS ordered separately

	ORDER BY COMPONENT	JAW INSERT	COMPONENTS INSULATED OUTER TUBES	OPTIONAL HANDLES
1/1	Atraumatic Parenchyma Fixation Forceps, 5 mm, 31 cm	PSI-L I -00281	PSI-L I -00282	••
1/1	Dissecting Forceps, delicate, 5 mm, 31 cm	PSI-L I -00283	PSI-L I -00284	••
1/1	Pike Mouth Aggressive Grasper, 10 mm, 31 cm	PSI-L I -00285	PSI-L I -00286	••
1/1	Large Metzenbaum Scissors, TC, serrated 10 mm, 31 cm	PSI-LI-010008	PSI-LI-010009	

#### www.primeinstrumentintl.com













PSI-LI-00287 (Non-Ratchet)

PSI-LI-00288 (Ratchet)

PSI-LI-00289 (Non-Ratchet)

PSI-L I -00290 (Ratchet)

PSI-L I -00291 (In-Line)

Instruments may be ordered COMPLETE or assembled from COMPONENTS ordered separately


	ORDER BY COMPONENT	JAW INSERT	COMPONENTS INSULATED OUTER TUBES	OPTIONAL HANDLES
1/2	Endo Lymph Node Fixation Forceps, 10 mm, 31 cm	PSI-L I -00292	PSI-L I -00293	••
1/2	Pennington Tissue Grasping Forceps, 10 mm, 37 cm	PSI-L I -00294	PSI-L I -00295	••
1/2	Duval Atraumatic Fixation Forceps, 10 mm, 37 cm	PSI-L I -00296	PSI-L I -00297	••
1/2	DeBakey Atraumatic Forceps, 40 mm Jaw, 5 mm, 31 cm	PSI-LI-010010	PSI-LI-010011	

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com


## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

Laparoscopic Instruments
Instructions for quick and easy assembly and disassembly
For use with PM950R and PM951R, insulated metal handles

#### **ASSEMBLY**




Slide the inner metal tube onto the jaw insert with its four spring tongues towards the working end until the tongues snap into place.



Slide the outer insulation tube as far as it will go over the inner metal tube.



Hold the shaft firmly by its tip and pull back the handle rotation knob with fingers and maintain hold. Fully insert the shaft and release rotation knob.




Test that the instrument is secured in the handle by pulling on the shaft. If shaft pulls out of handle, then repeat step 3 and step 4.



Connect handle with jaw insert ball joint by depressing the black buttons on both sides of handle, continue to depress and lift rear ring of handle under ball joint of jaw insert until captured by handle. Release buttons.




Test instrument by opening and closing handle to verify that the jaw has been inserted properly.




#### **DISASSEMBLY**



Depress black buttons on both sides of handle, continue to depress, pull rear ring down to release ball of jaw insert. Release buttons.



Pull back rotation knob and pull out shaft with other hand. Release rotation knob.



Hold the jaw insert at the tip and remove the outer insulation tube.



Hold the jaw insert at the tip and remove the inner metal tube.


### www.primeinstrumentintl.com

## Laparoscopic Instruments Instructions for quick and easy assembly and disassembly For use with QUICK-SNAP™ PO958R and PO959R handles

#### **ASSEMBLY**




Slide the inner metal tube onto the jaw insert with its four spring tongues towards the working end until the tongues snap into place.



Slide the outer insulation tube as far as it will go over the inner metal tube.



If the red ring is not visible, pull the rotation knob back. If you are using the ratchet style handle, deactivate ratchet by flipping down the red toggle switch.




Hold the shaft firmly by its tip and the handle behind the rotation knob. Push the shaft into the handle until it locks in place and the red ring disappears.

Note: Do not hold onto the movable handle component.



Test the assembly of the instrument by opening and closing the handle to ensure the jaw is opening and closing.




#### **DISASSEMBLY**




If ratchet style handle, deactivate ratchet by flipping down red toggle switch.



Pull rotation knob back until it locks into place. A red ring will become visible. Once visible remove the shaft by pulling the instrument shaft out of the handle.



Hold the jaw insert at the tip and remove the outer insulation tube.



Hold the jaw insert at the tip and remove the inner metal tube.

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

#### **MODULAR INSTRUMENT SPARE PARTS**

#### **Insulated Outer Tubes**



**PSI-LI-010012** mm, 31 cm **PSI-LI-00298** mm, 37 cm **PSI-LI-010016** mm, 42 cm



PSI-LI-010013 10 mm, 31 cm PSI-LI-010017 10 mm, 42 cm



PSI-LI-010014 5/10 mm, 31 cm PSI-LI-010018 5/10 mm, 37 cm



PSI-LI-010015 10 mm, 37 cm

PSI-L I -00301	Cautery Pin for PO handles, 5 pkg.
PSI-LI-010019	Cautery Pin for PM handles, 5 pkg.
PSI-L I -00302	Rotation Knobs for PM instruments, grey, 5 pkg.
PSI-LI-010020	Disposable Monopolar Cable, 12 ft.
PSI-L I -00303	Reusable Monopolar Cable, ERBE T 4 mm diam., 12 ft.
PSI-LI-010021	Reusable Monopolar Cable, VALLEYLAB, CONMED, BIRCHER, BOVIE, ERBE ICC, 4 mm diam., 10 ft.
PSI-L I -00304	Cleaning Brush, 3.7 mm diam. for Inner Tube and Outer Tube, 5 mm diam., 12 pkg.
PSI-LI-010022	Cleaning Brush, 10 mm diam., 6 pkg.
PSI-L I -00305	Inner Metal Tube, 31 cm
PSI-LI-010023	Inner Metal Tube, 37 cm
PSI-L I -00306	Inner Metal Tube, 42 cm

### Reposable Scissors

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

#### Reposable Scissors



#### Reusable Handle and Disposable Jaw Insert

- Sharp scissors tip EVERY TIME with disposable jaw insert
  - Cost effective alternative compared to 100% disposable scissors
  - Three different scissors options
  - Convenient handle choices

	ITEM NO.	DESCRIPTION
1/1	PSI-L I -00307	Disposable Slightly Curved Metzenbaum Scissors Jaw Insert, 5 mm, 31 cm, 10 per box
1/1	PSI-L I -00308	Disposable Hook Scissors Jaw Insert, 5 mm, 31 cm, 10 per box
1/1	PSI-L I -00309	Disposable Mini-Metz Scissors, 5 mm, 31 cm, 10 per box
1/1	PSI-L I -00310 PSI-L I -00311	Disposable Metzenbaum Scissors Jaw Insert, 5 mm, 31 cm, 10 per box Disposable Metzenbaum Scissors Jaw Insert, 5 mm, 42 cm, 10 per box
10	PSI-L I -00312	Quick-Snap™ Handle
	PSI-L I -00313	Insulated Metal Handle

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

#### **Monopolar Electrodes**

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

#### Modular Monopolar Electrode System

QUALITY - German craftsmanship and design

**DURABILITY** - PEEK\* insulation is resistant to cracks, peeling, shrinkage and malformations

MODULARITY - Exchange of tips and replacements as needed

ENHANCED SAFETY - Ceramic insulation to the working tip provides controlled energy delivery

#### SYSTEM COMPONENTS



Monopolar Handle - Complete instrument includes both Outer Tube and Inner Tube.

PSI-LI-00314 5 mm diameter, 33 cm working length PSI-LI-00315 5 mm diameter, 42 cm working length



#### **Outer Tube**

**PSI-LI-00316** 5 mm, 33 cm **PSI-LI-00317** 5 mm, 42 cm



#### Inner Tube and Handle

PSI-LI-00318 5 mm, 33 cm PSI-LI-00319 5 mm, 42 cm

#### MONOPOLAR ELECTRODE TIPS

P	SI-L I -00320	Ceramic Electrode Tip, J-Hook
	PSI-L I -00321	Ceramic Electrode Tip, L-Hook
-	PSI-L I -00322	Button Electrode Tip
	PSI-L I -00323	Spatula Electrode Tip

#### **REUSABLE MONOPOLAR CABLES**

	PSI-L I -00324	ERBE T Reusable Monopolar Cable, 4 mm diameter, 12 ft.
0-	PSI-L I -00325	VALLEYLAB, CONMED, BIRCHER, BOVIE, ERBE ICC Reusable Monopolar Cable,
		4 mm diameter, 10 ft.

#### **SPARE PARTS**

PSI-LI-010024	Contact Pin for HF Connection	
PSI-LI-010025	Wire brush for cleaning the working tips	

## **Suturing**



INTRACORPOREAL		31 CM	37 CM	42 CM
1/1	Durogrip TC Needle Holder, straight, 5 mm	PSI-L I -00326	PSI-LI-01010026	PSI-L I -00327
1/1	Durogrip TC Needle Holder, curved left, 5 mm	PSI-L I -00328	PSI-LI-01010027	PSI-L I -00329
MEASURE	Durogrip TC Needle Holder, curved right, 5 mm	PSI-L I -00330	PSI-LI-01010028	
1/1	"The Viper" Durogrip TC Suture Assistant /Needle Holder, $2 \times 3$ teeth, $5 \text{ mm}$	PSI-L I -00331		
1/1 WAFSPIL A	Super Righting Needle Holder, 5 mm <sup>SI-L I -00965</sup>	PSI-L I -00333		
1/1	Durogrip TC Needle Holder, Fenestrated, 5 mm	PSI-L I -00335		
1/1	Durogrip TC Needle Holder, Crimped/Fenestrated, 5 mm	PSI-L I -00337		
EXTRACORPOREAL				
1/1	CROCE/OLMI Dissecting and Suture Scissors with Knot Pusher, 5 mm	PSI-L I -00339		
1/1		PSI-L I -00341		Jaw Insert
	*Knot Pushing/Knot Tightening Forceps, 5 mm	PSI-L I -00342 PSI-L I -00343		Outer Tube Handle
	Knot Pusher, 5 mm, 33 cm	PSI-L I -00344		
51	Knot Pusher, 5 mm, 33 cm	PSI-LI-00332		

<sup>\*</sup> Instrument sold by component

## **Bipolar Instruments**

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### The enhanced safety of Ceramic Insulation in Reusable Bipolar Laparoscopic Modular Instruments




	COMPLETE INSTRUMENT	JAW INSERT	OUTER TUBES	STANDARD HANDLE
1/1	PSI-LI-00345 Bipolar Metzenbaum Scissors, 5 mm, 31 cm	: PSI-L I -00346	PSI-L I -00347	PSI-L I -00348
7	PSI-LI -00349 Bipolar Metzenbaum Scissors, 5 mm, 42 cm	PSI-L I -00350	PSI-L I -00351	PSI-L I -00352
1/1	PSI-LI-00353 Bipolar Mini Metzenbaum Scissors, 5 mm, 31 cm	PSI-L I -00354	PSI-L I -00355	PSI-L I -00356
1/1	PSI-L I -00357 Bipolar Maryland Dissector, 5 mm, 31 cm	PSI-L I -00358	PSI-L I -00359	PSI-L I -00360
	PSI-LI-00361 Bipolar Maryland Dissector, 5 mm, 42 cm	PSI-L I -00362	PSI-L I -00363	PSI-L I -00364
1/1	PSI-LI-00365 Bipolar Maryland Dissector, fenestrated 5 mm, 31 cm	PSI-L I -00366	PSI-L I -00367	PSI-L I -00368
	PSI-L I -00369  Bipolar Maryland Dissector, fenestrated 5 mm, 42 cm	PSI-L I -00370	PSI-L I -00371	PSI-L I -00372
1/1	PSI-LI-00373 Bipolar Micro Forceps, 5 mm, 31 cm	PSI-L I -00374	PSI-L I -00955	PSI-L I -00956
1/1	PSI-LI -00957 Bipolar Macro Forceps, 5 mm, 31 cm	PSI-L I -00958	PSI-L I -00959	PSI-L I -00960
	PSI-LI-00961  Bipolar Macro Forceps, 5 mm, 42 cm	PSI-L I -00962	PSI-L I -00963	PSI-L I -00964

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

#### Disposable Bipolar Instruments



П	EM NO.	DISPOSABLE BIPOLAR SCISSORS
There	PSI-L I -00375	Metzenbaum Scissors, 5 mm, 31 cm, Six Single Sterile Scissors per Box
COLLAP.	PSI-L I -00376	Maryland Grasping Forceps, 5 mm, 31 cm, Six Single Sterile Graspers per Box
R CABLES & ACCE	SSORIES	
	PSI-L I -00377	Bipolar Cable for Aesculap, Wolf, Martin, Berchtold, 12 ft.
	PSI-L I -00378	4 m length, suitable for units like: Aesculap, BOWA, EMC, Lamidey, Söring, Valleylab, 28.6 mm Pin Distance
	PSI-L I -00379	4 m length, suitable for units like: Codman CMC II, CMC III, 38 mm Pin Distance
	PSI-L I -00380	Bipolar Cable for Valleylab, Conmed, Bovie, 12 ft. (Disposable, Single-use - 10 units/box)
Cable for Erbe		
Cable for Codman		
Wire Brush for cleaning	jaws	
	Cable for Erbe	PSI-LI-00376  R CABLES & ACCESSORIES  PSI-LI-00377  PSI-LI-00379  PSI-LI-00380  Cable for Erbe



	ITEM NO.	CLASSIC BIPOLAR TUBAL FORCEPS
-	PSI-L I -00384 PSI-L I -00385	32 cm, Syringe Handle Forceps 45 cm, Syringe Handle Forceps
	PSI-L I -00386 PSI-L I -00387	32 cm, Scissors Handle Forceps 45 cm, Scissors Handle Forceps
- 100	PSI-L I -00388	32 cm, Jaw Insert Only
	PSI-L I -00389	45 cm, Jaw Insert Only
	PSI-L I -00390	Reusable Cable for Wolf Generator
	PSI-L I -00391	Reusable Cable for Valley Lab Generator

### **Specialty-Locking Graspers**

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

#### Slide-Lock Graspers

- Single- and double-action models
- Atraumatic and retraction jaw patterns
- Flush ports for optimal cleaning
- Excellent tissue grasping properties

ITEM NO.	DESCRIPTION	HANDLE	SHAFT
PSI-L I -00392	Sovereign Atraumatic Grasper, single action, 5 mm, 42 cm	Slide Lock	Non-Rotatable
PSI-L I -00393	Sovereign Atraumatic Grasper, double action, 5 mm, 42 cm	Slide Lock	Rotatable
PSI-L I -00394	Sovereign Atraumatic Grasper, single action, 5 mm, 36 cm	Slide Lock	Non-Rotatable
PSI-L I -00395	Sovereign Atraumatic Grasper, double action, 5 mm, 36 cm	Slide Lock	Rotatable
PSI-L I -00396	Sovereign Retraction Grasper, single action, 5 mm, 36 cm	Slide Lock	Non-Rotatable
PSI-L I -00397	Sovereign Retraction Grasper, double action, 5 mm, 36 cm	Slide Lock	Rotatable
PSI-L I -00398	Prestige Atraumatic Grasper, single action, 5 mm, 36 cm	Slide Lock	Non-Rotatable
PSI-L I -00399	Prestige Atraumatic Grasper, double action, 5 mm, 36 cm	Slide Lock	Rotatable
PSI-L I -00400	Prestige Retraction Grasper, single action, 5 mm, 36 cm	Slide Lock	Non-Rotatable
PSI-LI-00334	Prestige Retraction Grasper, double action, 5 mm, 36 cm	Slide Lock	Rotatable



Sovereign® Graspers



Prestige<sup>™</sup> Graspers



Sovereign PSI-LI-00401 Prestige



Sovereign PSI-L I -00402 Prestige




Sovereign PSI-L I -00403 Prestige



Sovereign PSI-L I -00404 Prestige

#### **Specialty Instruments Sets - Advanced Procedure**

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com



#### Atraumatic Endo Vessel Clips (for temporary occlusion)

		CLOSING FORCE	JAW LENGTH	JAW SHAPE
	PSI-L I -00408	350 g	25 mm	short/straight
	PSI-L I -00409	450 g	45 mm	long/straight
	PSI-L I -00410	350 g	25 mm	short/curved
	PSI-L I -00411	450 g	45 mm	long/curved


#### **KIDNEY CLAMP**



### **Specialty Instruments Sets - Advanced Procedure**

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

Optional Handles for Customizing Instruments - Interchangeable for 5 mm and 10 mm instruments



		COMPONENTS				
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES	
	Potts-DeMartell Scissors, 10 mm, 31 cm  Potts-DeMartell Scissors, 10 mm, 31 cm  Potts-Li-00423 Protection Sleeve for safe introduction through the 12.5 mm trocar	PSI-L I -00421	PSI-L I -00422	•	•	
1/2	PSI-LI-00424  Satinsky Atraumatic Clamp, 38 mm jaw, 5 mm, 31 cm  PSI-LI-00427  Satinsky Atraumatic Clamp, 51 mm jaw, 5 mm, 31 cm	PSI-L I -00425	PSI-L I -004 PSI-L I -004			
1/2	Order by component DeBakey Atraumatic Clamp, straight, 50 mm Jaw, 10 mm, 22 cm	PSI-L I -00430	PSI-L I -00	0431		
1/2	PSI-LI-00432  DeBakey Atraumatic Clamp, curved, 5 mm, 31 cm	PSI-L I -00433	PSI-L I -00	1434		

<sup>\*</sup> For use with 12.5 mm flexible trocars

## PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

Optional Handles for Customizing Instruments - Interchangeable for 5 mm and 10 mm instruments



Instruments may be ordered COMPLETE or assembled from COMPONENTS ordered separately

PSI-L I -00440 (Ratchet)

PSI-L I -00441

441 (Ratchet) with lock

		COMPONENTS			
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/2	PSI-L I -00442	PSI-L I -00443	PSI-L I -00444		
	DeBakey Atraumatic Clamp, S-shaped, 5 mm, 31 cm				
1/2	PSI-LI-00445 Crafoord Endo-vascular Clamp, short jaw,	PSI-L I -00446	PSI-L I -00	447	
	5 mm, 31 cm PSI-LI-00448 Crefood Folo vectoriar Clamp long inv	PSI-L I -00449	PSI-L I -004	450	-
	Crafoord Endo-vascular Clamp, long jaw, 5 mm, 31 cm				
1/1	PSI-L I -00451	PSI-L I -00452	PSI-L I -00-	453	
The same of the sa	Overholt Dissecting Clamp, 90°, 10 mm, 31 cm				
1/1	PSI-L I -00454	PSI-L I -00455	PSI-L I -00456	5	
	Crile Dissecting Forceps, Curved 5 mm, 31 cm				
1	PSI-L I -00457	PSI-L I -00458	PSI-L I -00459		
	Crile Dissecting Forceps, Curved 5 mm, 37 cm				

<sup>\*</sup> For use with 12.5 mm flexible trocars

PSI-L I -00439 (In-Line)

# Specialty Instruments Sets - Advanced Procedure

# PRIME SURGICAL INTERNATIONAL


### www.primeinstrumentintl.com


### **FLEXIBLE TROCARS**



### Complete trocar -100 mm length with valve and obturator


- Shaft sufficiently flexible to permit insertion of curved instruments and sufficiently firm to ensure access to the peritoneal orthoracic cavity
- For the use of curved instruments or special mouth shapes
- Trocar shaft in two different diameters and two different lengths







# **SLS Single-Fire Ligation System**



### A Reliable and Cost-Efficient System that Delivers Superior Performance

- Clip cartridge is designed for easy loading and precise clip pickup
- Rotation knob allows for one finger, 360° rotation of applier shaft
- Modular clip applier can be disassembled for easy cleaning

ORDERING INFORMATION					
	Small	Small-Medium	Medium	Medium-Large	Large
Titanium Ligation Clips	Available in boxes of 30	Available in boxes of 30	Available in boxes of 30	Available in boxes of 20	Available in boxes of 20
	cartridges of 6 clips = 180 clips	cartridges of 6 clips = 180 clips	cartridges of 6 clips = 180 clips	cartridges of 6 clips = 120 clips	cartridges of 6 clips = 120 clips
Laparoscopic	PSI-L I -00476	PSI-LI-00340	PSI-L I -00477	PSI-LI-010026	PSI-L I -00478
Clip Appliers	Diameter	Diameter	Diameter	Diameter	Diameter
13" (33 cm)	5.5 mm	5.5 mm	6.5 mm	10.0 mm	12.5 mm
-8	Length	Length	Length	Length	Length
71	30 cm	30 cm	30 cm	33 cm (shown left)	33 cm
/ (			PSI-L I -00479		
			Diameter		
			10.0 mm Length 33 cm		
			33 CIII		

# Challenger Ti-P Multi-Fire Ligation System

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### **MULTI-FIRE**

- Fast, precise and reliable clip application complimented by an ergonomic delivery
- Shafts available in 5 mm and 10 mm diameters
- Superior quality ligation clip available in S/M and M/L sizes
- Reusable clip applier, disposable clip cartridges
- Easy and simple handling

### **CHALLENGER TI-P PNEUMATIC CLIP APPLIERS**



### PSI-LI-010027

5 mm diam., 31 cm length - 12 S/M clips per cartridge

### ERERRESE ST

### PSI-LI-010028

10 mm diam., 37 cm length
 8 M/L clips per cartridge



CLIP CARTRIDGES & CO2 CYLINDERS

S/M Clips

M/L Clips



### PSI-LI-010029

- 12 sterile clip cartridges
  - 1 CO<sub>2</sub> cylinder per cartridge
  - For use with 5 mm clip applier



### PSI-LI-010030

- 12 sterile clip cartridges
  - 1 CO<sub>2</sub> cylinder per cartridge
  - For use with 10 mm clip applier





	Complete Shaft	PSI-L I -00480			PSI-L I -00481		PSI-L I -00482	
		Shaft Leng 310 mm	th Diameter 5 mm		Shaft Length 370 mm	Diameter 10 mm	Shaft Length 260 mm	Diameter 10 mm
	Handle	PSI-L I -00483			PSI-L I -00484		PSI-L I -00485	
	Inner Tube	PSI-L I -0048	36		PSI-L I -00487		PSI-L I -00488	
		Shaft Leng 310 mm	th		Shaft Length 370 mm		Shaft Length 260 mm	
ts	Push Rod	PSI-L I -0048	9		PSI-L I -00490		PSI-L I -00491	
re Parts		Shaft Leng 310 mm	th		Shaft Length 370 mm		Shaft Length 260 mm	
Spare	Locking Ring	PSI-L I -00492			PSI-L I -00493		PSI-L I -00494	
•								
	CO <sub>2</sub> Cylinder Only	PSI-L I -00495	5		PSI-L I -00496		PSI-L I -00497	
	Pack of 3 Single Sterile Packed			PRIN	ME SURGIC	AL INTER	NATIONAI	_

# Vascular Clip Applier/Remover

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com



### **CLIP APPLIER/REMOVER**



PSI-L I -00498

Articulating Applicator for temporary endoscopic vessel clips, 12.5 mm, 34 cm length



SI-L I -00499

Fixed Applicator for temporary endoscopic vessel clips, 12.5 mm, 34 cm length

PSI-LI-010031

Rotation Knob, 5/package

### Atraumatic Endo Vessel Clips (for temporary occlusion)

ITEM NO.	CLOSING FORCE	JAW LENGTH	JAW SHAPE
PSI-L I -00500	350 g	25 mm	short/straight
PSI-L I -00501	450 g	45 mm	long/straight
PSI-L I -00502	350 g	25 mm	short/curved
PSI-L I -00503	450 g	45 mm	long/curved

### INTESTINAL BULLDOG CLIP



PSI-L I -00504

Intestinal bulldog clip,

Jaw length: 70 mm, Closing force: 3.94 mm



# **HD Endoscopes**

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com



### **LAPAROSCOPES**

ITEM NO.	DESCRIPTION	DIAMETER	VIEWING DIRECTION	LENGTH	Steam Autoclave	STERRAD®
PSI-L I -00505	Needlescopic, Semi-rigid	2 mm	0°	30 cm	•	•
PSI-L I -00506	Needlescopic, Semi-rigid	2.5 mm	0°	22 cm	•	•
PSI-L I -00507	Full HD Endoscope	4 mm	0°	180 mm	•	•
PSI-L I -00508	Full HD Endoscope	4 mm	30°	180 mm	•	•
PSI-L I -00509	Full HD Endoscope	4 mm	0°	300 mm	•	•
PSI-L I -00510	Full HD Endoscope	5 mm	30°	300 mm	•	•
PSI-L I -00511	Full Screen Technology	5 mm	0°	30 cm	•	•
PSI-L I -00512	Full Screen Technology	5 mm	30°	30 cm	•	•
PSI-L I -00513	Full Screen Technology	10 mm	0°	33 cm	•	•
PSI-L I -00514	Full Screen Technology	10 mm	30°	33 cm	•	•
PSI-L I -00515	Full Screen Technology	10 mm	45°	33 cm	•	•
PSI-L I -00516	Operative Laparoscope,	10 mm	0°	27 cm	•	•
	with working channel for 5 mm instruments					

### STORAGE SYSTEMS

PSI-L I -00517	Single Scope Tray, 17", for 5 mm/10 mm diameter scopes
PSI-L I -00518	Double Scope Tray, 17", for 5 mm/10 mm diameter scopes
PSI-L I -00519	Double Scope Tray, 17", for 2 mm/3 mm diameter scopes
PSI-L I -00520	Half-Size Storage Basket for 4 mm Perneczky angled scopes



# Unitrac® Retraction & Holding System

### PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

The Unitrac is a pneumatically-powered retraction and scope holding system that allows you to experience a more efficient procedure and obstruction-free surgical site. Precise adjustments are possible with the push of a single button, thus freeing the surgical assistant to carry out higher priority tasks.

### **Benefits Features** Single-piece construction utilizes compressed air Easy setup and disassembly source in the OR Extensive range of retractor blade, hook, Applications in Minimally Invasive and Open Surgical procedures camera and endoscope holders Reposition with push of a single button Increases efficiency during procedures Maintains full holding force in event of OR Protects patient and OR personnel safety compressed air failure

### UNITRAC ARM CONFIGURATIONS Unitrac Arm



### **OR Table Clamp**



OR Table Clamp w/ Flexible Joint



Compressed Air Hose



Y-connector for attaching 2 Unitrac arms to same compressed air source



# Unitrac® Retraction & Holding System


# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

The Unitrac is available in several configurations to ensure you have exactly what you need to achieve absolute rigid retraction and stable positioning. All configurations include the Unitrac pneumatic arm, a standard or flexible joint OR table clamp, and a compressed air hose or Y-connector for 2 Unitrac arms to be in use during the procedure.

Use the below ordering guide to select the Unitrac configuration that is right for you.

	Item No.			
	PSI-L I -00521	PSI-L I -00522	PSI-L I -00523	PSI-L I -00524
Unitrac Arm	•	•	•	•
OR Table Clamp	•	•		
OR Table Clamp w/ Flexible Joint			•	•
Compressed Air Hose	•		•	
Y-connector for attaching 2 Unitrac arms to same compressed air source		•		•

# Item No. PSI-L1-00525 Adjustable holder for endoscopes with 3.0 - 7.5 mm diam., with PEEK insert RT055P Replacement adjustable PEEK insert, 3.0 - 7.5 mm diam. PSI-L1-00526 PSI-L1-00527 Holder for endoscopes with 10 mm diam., with PEEK insert RT086P Replacement PEEK insert, 10 mm diam.



# **FROCARS**

# Reusable Trocar Systems

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### 5.5 mm TROCARS WORKING **COMPLETE CANNULA PYRAMIDAL TIP CONICAL TIP BLUNT TIP TROCAR ONLY OBTURATOR LENGTH OBTURATOR OBTURATOR** 60 mm PSI-LI-010032 PSI-LI-010033 PSI-LI-010034 PSI-L I -00529 PSI-L I -00531 PSI-L I -00532 PSI-L I -00533 110 mm PSI-L I -00530 Smooth Cannula PSI-LI-010037 PSI-LI-010036 PSI-LI-010035 150 mm 60 mm PSI-LI-010040 PSI-LI-010038 PSI-LI-010039 PSI-L I -00534 PSI-L I -00535 PSI-L I -00536 PSI-L I -00537 110 mm Threaded Cannula 150 mm 60 mm PSI-LI-010042 PSI-LI-010041 PSI-LI-010043 PSI-L I -00539 PSI-L I -00540 PSI-L I -00541 110 mm PSI-L I -00538 Radiolucent Threaded Cannula PSI-LI-010045 PSI-LI-010044 150 mm





(order EJ700251 for insufflation ring)



PSI-L I -00543 Manual Opening Trocar Body

### Sealing Caps for Reducing Sleeve 20/pkg

5.5 mm/10 mm 5.5 mm/12.5 mm 10 mm/12.5 mm
PSI-LI-010046 PSI-LI-010047 PSI-LI-010048

# ROCARS


# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### 10 mm TROCARS WORKING **PYRAMIDAL TIP CANNULA CONICAL TIP BLUNT TIP COMPLETE** LENGTH **TROCAR ONLY OBTURATOR OBTURATOR OBTURATOR** 60 mm 110 mm PSI-L I -00545 PSI-L I -00546 PSI-L I -00548 PSI-L I -00549 PSI-L I -00547 Smooth Cannula 150 mm 60 mm 110 mm PSI-L I -00550 PSI-L I -00551 PSI-L I -00552 PSI-L I -00553 Threaded Cannula 150 mm 60 mm PSI-L I -00554 PSI-L I -00555 110 mm PSI-L I -00556 PSI-L I -00557





Radiolucent Threaded Cannula



PSI-L I -00559 Manual Opening Trocar Body



**Reducing Converter** 

ITEM NO. TROCAR DIAMETER CONVERTS TO:

10 mm & 12.5 mm 5 mm/10 mm



PSI-LI-00560 10 mm PSI-LI-00561 180 mm PSI-LI-00562 5.5 mm

150 mm

Sealing Caps for Reducing Sleeve 20/pkg

5.5 mm/10 mm 5.5 mm/12.5 mm 10 mm/12.5 mm

# **Reusable Trocar Systems**

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com





Standard Trocar Body, without insufflation ring (order EJ750251 for insufflation ring)



### Reducing Sleeve

ITEM NO.	DIAMETER	LENGTH	GUIDE
PSI-L I -00575	13 mm	180 mm	5.5 mm
PSI-L I -00576	13 mm	180 mm	10 mm



### **Reducing Converter**

ITEM NO.	TROCAR DIAMETER	CONVERTS TO:
PSI-LI-010049	13 mm	5 mm, 10 mm, 13 mm

### Sealing Caps for Reducing Sleeve 20/pkg

5.5 mm/10 mm	5.5 mm/12.5 mm	10 mm/12.5 mm	5.5 mm/13 mm	10 mm/13 mm
PSI-LI-010050	PSI-LI-010051	PSI-LI-010052	PSI-LI-010053	PSI-LI-010054

# **TROCARS**

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### HASSON/MINI LAPAROTOMY STYLE TROCARS



COMPLETE TROCAR W/OBTURATOR CANNULA **ONLY** 

**BLUNT TIP OBTURATOR** 



100 mm

PSI-L I -00577

PSI-L I -00578

PSI-L I -00579









PSI-L I -00581 Stability Cone

### **TROCAR SPARE PARTS**





NO.	O.D. CANNULATED NEEDLE	O.D. CANNULA NEEDLE	LENGTH
PSI-L I -00584	1.6 mm	2.1 mm	100 mm
PSI-L I -00585 PSI-L I -00586	1.6 mm	2.1 mm	120 mm
	1.6 mm	2.1 mm	150 mm

# Reusable Trocar Systems

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### **Trocar Spare Parts (continued)**



PSI-LI-00588 Silicone Flap Valve for 5.5 mm or 7 mm trocars, 20/pkg.
PSI-LI-00588 Silicone Flap Valve for 10 mm or 12.5 mm trocars, 20/pkg.
PSI-LI-00589 Silicone Flap Valve for 13 mm trocars, 20/pkg.



PSI-LI-00590 Sealing Cap for Insufflation Connection for modular Trocars, 20/pkg.



PSI-LI-00591 Trocar Insufflation Ring with Luer Lock, 10 mm/12.5 mm, one piece
PSI-LI-00592 Trocar Insufflation Ring with Luer Lock, 5 mm/7 mm, one piece



PSI-LI-00593 Trocar Ring without Luer Lock Connector, 10 mm/12.5 mm, one piece
PSI-LI-00594 Trocar Ring without Luer Lock Connector, 5 mm/7 mm, one piece



PSI-LI-00595 Plastic Ring for 5 mm trocars, 20/pkg.
PSI-LI-00596 Plastic Ring for 10 mm trocars, 20/pkg.
PSI-LI-00597 Plastic Ring for 13 mm trocars, 20/pkg.



### Metal Valve for Insufflation Connection, one piece



PSI-LI-00599 Sealing Cap for Modular Trocars, 5.5 mm, 20/pkg.
PSI-LI-00599 Sealing Cap for Modular Trocars, 7 mm, 20/pkg.
PSI-LI-00601 Sealing Cap for Modular Trocars, 10 mm, 20/pkg.
PSI-LI-00601 Sealing Cap for Modular Trocars, 12.5 mm, 20/pkg.
PSI-LI-00602 Sealing Cap for Modular Trocars, 13 mm, 20/pkg.



PSI-L I -00603 Reducing Cap for Modular Trocars, 5/pkg.



PSI-L I -00604 Reducing Converter; 5 mm, 10 mm, 12.5 mm

### Flexible Trocars


# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com


### **FLEXIBLE TROCARS**



### Complete trocar -100 mm length with valve and obturator

- Shaft sufficiently flexible to permit insertion of curved instruments and sufficiently firm to ensure access to the peritoneal orthoracic cavity
- For the use of curved instruments or special mouth shapes
- Trocar shaft in two different diameters and two different lengths





PRIME SURGICAL INTERNATIONAL GYNECOLOGIC SURGICAL INSTRUMETNS prime.instruments@hotmail.com

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com



# **Light Sources**

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### **AXeL 180 XENON LIGHTSOURCE**





TE	CLIN	HCAL	DA-	ГΛ
ır	C HI	JICAI	UA	I A

Lamp	Xenon
Lamp Power	180 W
Lamp Voltage	15 V/12.0A

Light Intensity Continuously Adjustable

Color Temperature Approx. 6000 K

Dimensions (w x h x d) 12" x 5" x 12" (305 x 125 x 305 mm)

Weight 17 lbs., 17 oz. (6.5 kg)

### ITEM NO. DESCRIPTION

PSI-LI-010055 AXeL Xenon Lightsource (180 W)

Consisting of: PSI-L I -00618

PSI-L I -00619 Xenon Lamp, 15V, 180 W
PSI-L I -00619 Standard Adapter, Aesculap

### SYSTEM ACCESSORIES

PSI-L I -00620	Wolf LIGHT SOURCE Adapter
PSI-L I -00621	Olympus Lightsource Adapter
PSI-L I -00622	Circon ACMI Lightsource Adapter
PSI-L I -00623	Light Guide Cable, diameter: 5 mm, length: 250 cm (8.2 ft.) with additional bending protection
PSI-LI-010056	Xenon Spare Lamp

**VIDEO SYSTE** 

# **Light Cables**

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com


### **AESCULAP LIGHT CABLES and ADAPTERS**



### **Light Cable Ordering Information**

- Follow these steps to use the table below to identify the appropriate light cable:
  - Step 1 Select your light source from the horizontal column.
  - Step 2 Select your scope from the vertical column.
  - Step 3 The box where they meet represents the appropriate light cable.

(Example: Using an PRIME light source and an PRIME scope - PSI-LI-00624




• The universal light cable and all adapters can also be ordered separately

PSI-L I -00644	Universal light cable with no end fittings, Light fiber bundle: 5 mm, Length: 10 ft.
PSI-L I -00645	Aesculap / Storz adapter for light cable (LS end)
PSI-L I -00646	Olympus adapter for light cable (LS end)
PSI-L I -00647	Circon / ACMI adapter for light cable (LS end)
PSI-L I -00648	Stryker adapter for light cable (LS end)
PSI-L I -00649	Wolf/Dyonics adapter for light cable (LS end)
PSI-L I -00650	Aesculap / Storz / Olympus adapter for light cable (S end)
PSI-L I -00651	Circon / ACMI / Stryker adapter for light cable (S end)
PSI-L I -00652	Wolf/Dyonics adapter for light cable (S end)
PSI-L I -00653	Light cable tester



PRIME SURGICAL INTERNATIONAL GYNECOLOGIC SURGICAL INSTRUMETNS prime.instruments@hotmail.com



### Miscellaneous PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### **DISSECTING SPATULA** PSI-L I -00654 Dissecting Spatula, straight, 5 mm, 33 cm length PSI-LI-010057 Dissecting Spatula, curved, 5 mm, 33 cm length (shown) SCALPEL HANDLE for KNIFE BLADES Handle Features a retractable blade mechanism PSI-L I -00655 Scalpel Handle, 5 mm, 33 cm length KNIFE BLADES (order separately 10 per package) PSI-L I -00656 PSI-LI-010058 PSI-L I -00657 PSI-LI-010059 REUSABLE SUCTION IRRIGATION SYSTEMS PSI-L I -00658 Suction Tube With Luer-Lock connection and finger valve. Can be used with Y-connector, SK330C, 5 mm diam., 33 cm length. PSI-L I -00659 Suction/Irrigation Tube With Luer-Lock connection. Can be used with Y-connector, SK330C, 5 mm diam., 33 cm length. PSI-L I -00660 For suction and irrigation in connection with SR550R, SR551R, SR552R, SR553R, GF320R and GF321R Y-Connector PSI-L I -00661 Puncture Needle With Luer-Lock connection and finger valve, 5 mm diam., 33 cm length PSI-L I -00662 Injection Cannula, Sharp With Luer-Lock connection, 5 mm diam., 33 cm length

With Luer-Lock connection, 5 mm diam., 45 cm length

Injection Cannula, Sharp

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### LIGATURE INSTRUMENTS



PSI-L I -00664

Loop Applier 13 cm length Outside diam. - 5 mm Inside diam. - 4.4 mm

PSI-L I -00665

Loop Applier Sealing Cap

20/pkg.



PSI-L I -00666

Fascia Port-Closure Instrument 2.7 mm diam., 15.2 cm length, shaft only PL163R

### ADDITIONAL ITEMS



PSI-L I -00667

Myoma Fixation Screw

5 mm diam., 31 cm length



PSI-L I -00668

**Cholangiography Clamp** with Catheter Guide

Gall Bladder Extractor

(Maximum catheter diameter-5 French), ratcheted, 5 mm diam., 33 cm length

PSI-L I -00669

Sealing Cap for PL100R 20/pkg.



PSI-L I -00670

Three-blade design, 55 mm blade length.

Use with 10 mm trocars.

### Miscellaneous

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### **ADDITIONAL ITEMS**

3\_\_\_\_

PSI-L I -00671

**Exploratory Probe** 5 mm diam., 40 cm length

PSI-L I -00672

Exploratory Probe 5 mm diam., 50 cm length



PSI-L I -00673

Anvil Grasper, 10 mm, 37 cm

### **LUBRICANTS**



PSI-L I -00674

Aesculap Instrument Oil For box locks, ratchets, etc. Physiologically inert.



PSI-L I -00675

Aesculap Oil Spray
For maintenance of instruments prior to sterilization
Physiologically inert.

### **ENDOSCOPY TRAYS and BASKETS\***



PSI-L I -00676

Single Scope Tray, 10"
Two fit securely into
a Quarter Size container



PSI-L I -00677 PSI-L I -00678

Double Scope Tray, 17", with or without silicone supports Fits into a full size basket  $17\frac{1}{2}$ " long x  $3\frac{3}{4}$ " wide x  $1\frac{3}{6}$ " deep

# PRIME SURGICAL INTERNATIONAL

### www.primeinstrumentintl.com

### **ENDOSCOPY TRAYS and BASKETS\***



PSI-L I -00679

Two Scope Basket Includes: JF222R basket, JG691 (4) instrument holding clamps

### PSI-L I -00683

Three Scope Basket Includes: JF222R basket, JG691 (6) instrument holding clamps

### PSI-L I -00684

Four Scope Basket Includes: JF222R basket, JG691 (8) instrument holding clamps



PSI-L I -00680

Endorack™ with Silicone Cushioning Pad and Basket Holds 14 endoscopic instruments



PSI-L I -00685

Endorack, Full-Size Storage System Endorack with silicone mesh mat and basket; holds 18 endoscopic instruments



PSI-LI-00681 - Solid Bottom PSI-LI-00682 - Perforated Bottom **Full Size Container** 

## Miscellaneous

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### **EXTRA LONG CONTAINER SYSTEM\***



PSI-L I -00686

Perforated Bottom with Retention Plates




PSI-L I -00687

Lid with Retention Plates



PSI-L I -00688

Endorack™ with Silicone Cushioning Pad and Basket Holds 14 endoscopic instruments



PSI-L I -00689

Endorack, Extra-Long Storage System Endorack with silicone mesh mat and basket; holds 18 endoscopic instruments

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com Laparoscopic Instruments

For Gynecologic Procedures



PRIME SURGICAL INTERNATIONAL offers a wide variety of laparoscopic GYN instruments that improve surgical performance and enhance patient care during minimally invasive surgery.

also provides many innovative solutions for laparoscopy in addition to standard instrument patterns demanded by all surgeons. The Viper Needle Holder and Bulldog Vascular Clips are two of the innovative products Aesculap offers GYN surgeons today.

**PRIME'S** "Viper" Needle Holder is a truly unique instrument that was specifically designed with GYN surgeons in mind. The 2 x 3 teeth jaw design assists surgeons in firmly holding tough tissue, primarily the vaginal cuff, during suturing. The Viper not only provides great grasping and holding capabilities but it allows surgeons to securely grasp needles without having to interchange instruments throughout the suturing process.

**PRIME'S** Vascular Bulldog Clips are another specialty instrument that GYN surgeons utilize in hysterectomy and myomectomy procedures. The Bulldog Clips are primarily used to clamp uterine arteries to minimize back bleeding when dissecting the uterus away from the uterine wall or during the removal of fibroids.

Whether you are looking for a grasper or specialty instrument be assured that **PRIME'S** main focus is on providing high quality instruments and superior service to all customers.



# **Grasping and Dissecting**

### PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com













PSI-L I -00789 (Non-Ratchet)

PSI-L I -00790 (Ratchet)

PSI-LI-00791 (Non-Ratchet)

PSI-LI-00792 (Ratchet)

PSI-L I -00793 (In-Line)

Instruments may be ordered COMPLETE or assembled from COMPONENTS ordered separately

			COMPONE	NTS	
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	PSI-LI-00794 Tubal Fixation Forceps, 5 mm, 31 cm PSI-LI-00797	PSI-L I -00795	PSI-L I -00796		
	Tubal Fixation Forceps, 5 mm, 42 cm	PSI-L I -00798	PS1-L 1 -00/99		
1/1	PSI-LI-00800  Maryland Dissecting Forceps, 5 mm, 31 cm	PSI-L I -00801	PSI-L I -00802		
annumini like	PSI-LI-00803  Maryland Dissecting Forceps, 5 mm, 42 cm	PSI-L I -00804	PSI-L I -00805		**
1/1	PSI-LI-00806 Right Angle, 90° Dissecting Forceps 10 mm, 31 cm	PSI-L I -00807	PSI-L I -00808		•
1/1	PST-L I -00809  Dolphin Nose Dissecting Forceps 5 mm, 31 cm	PSI-L I -00810	PSI-L I -0081:		
1/1	Duck Bill Dissecting Forceps 5 mm, 31 cm	PSI-L I -00813	PSI-L I -004	<b>314</b>	

# **Grasping and Dissecting**

### PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com













PSI-L I -00815 (Non-Ratchet)

PSI-L I -00816 (Ratchet)

PSI-L I -00819 (In-Line)

Instruments may be ordered COMPLETE or assembled from COMPONENTS ordered separately

			COMPONE	NTS	
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	PSI-L I -00820  Bullet Nose Grasper, 5 mm, 31 cm  PSI-L I -00823  Bullet Nose Grasper, 5 mm, 42 cm	PSI-L I -00821 PSI-L I -00824	PSI-L I -00822 PSI-L I -00825		
1/1	Aggressive Grasping Forceps 5 mm, 31 cm	PSI-L I -00827	PSI-L I -00828		::
1/1	PSI-LI-00829 Claw Grasper, 2 x 3 teeth 10 mm, 31 cm	PSI-L I -00830	PSI-L I -00831	÷	::
1/1	PSI-LI-00832 Babcock Grasping Forceps Single Action, 5 mm, 31 cm	PSI-L I -00833	PSI-L I -00834	•	
1/1	PSI-LI-00835 Babcock Grasping Forceps Double Action, 5 mm, 31 cm	PSI-L I -00836	PSI-L I -008	37	
	PSI-LI-00838  Babcock Grasping Forceps  Double Action, 5 mm, 42 cm	PSI-L I -00839	PSI-L I -00:	840	**

# **Grasping and Dissecting**

### PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### SOVEREIGN® QUICK-SNAP™ design





Insulated metal design PSI-LI-00843 (Non-Ratchet)





PSI-LI-00841 (Non-Ratchet)

PSI-L I -00842 (Ratchet)

PSI-L I -00845 (|n-Line)

Instruments may be ordered COMPLETE or assembled from COMPONENTS ordered separately

			COMPONE	NTS	
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	STANDARD HANDLE	OPTIONAL HANDLES
1/1	PSI-L I -00846  Large Biopsy Spoon Forceps 10 mm, 31 cm	PSI-L I -00847	PSI-L I -00848		::
1/1	PSI-L I -00849 Biopsy Forceps, Box, 5 mm, 31 cm	PSI-L I -00850	PSI-L I -00851	•	
	Biopsy Forceps, Box, 5 mm, 42 cm	PSI-L I -00853	PSI-L I -00854	•	
1/1	PSI-LI-00855 Biopsy Forceps, with teeth, 5 mm, 31 cm	PSI-L I -00856	PSI-L I -00857	•	
1/1 AFSCULAP	Tenaculum Forceps, 5 mm, 31 cm	PSI-L I -00859	PSI-L I -00860	•	
	Tenaculum Forceps, 10 cm, 31 cm	PSI-L I -00862	PSI-L I -00863	•	-
1/1	Fundus Grasper, 5 mm, 31 cm	PSI-L I -00865	PSI-L I -00866	•	
1/1	Short Fenestrated Forceps, 5 mm, 31 cm	PSI-L I -00868	PSI-L I -00869	•	

# **Cutting**

### PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

Optional Handles for Customizing Instruments - Interchangeable for 5 mm and 10 mm instruments





PSI-L I -00870 (Non-Ratchet)

PSI-LI-00871 (Non-Ratchet)

Instruments may be ordered COMPLETE or assembled from COMPONENTS ordered separately

### SCISSORS, MONOPOLAR

	ITEM NO.	DESCRIPTION
1/1	PSI-L I -00872	Disposable Hook Scissors Jaw Insert, 5 mm, 31 cm, 10 per box for PO/PM handles
1/1	PSI-L I -00873	Disposable Mini-Metz Scissors Jaw Insert, 5 mm, 31 cm, 10 per box for PO/PM handles
1/1	PSI-L I -00874	Disposable Metzenbaum Scissors Jaw Insert, 5 mm, 31 cm, 10 per box for PO/PM handles Disposable Metzenbaum Scissors Jaw Insert, 5 mm, 42 cm, 10 per box for PO/PM handles

### MONOPOLAR ELECTRODES

### **SYSTEM COMPONENTS** PSI-L I -00876 Complete Monopolar Handle, 5 mm, 33 cm Complete Monopolar Handle, 5 mm, 42 cm PSI-L I -00877 MONOPOLAR ELECTRODE TIPS PSI-L I -00878 Ceramic Electrode Tip, J-Hook PSI-L I -00879 Ceramic Electrode Tip, L-Hook

# Hemostasis

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

**BIPOLAR INSTRUMENTS** 



	-		COMPONENTS	
	COMPLETE INSTRUMENT	JAW INSERT	INSULATED OUTER TUBES	HANDLE
1/1	PSI-LI-00880  Bipolar Metzenbaum Scissors, 5 mm, 31 cm	PSI-L I -00881	PSI-L I -00882	PSI-L I -00883
1/1	Bipolar Maryland Dissector, 5 mm, 31 cm	PSI-L I -00885	PSI-L I -00886	PSI-L I -00887
1/1	PSI-LI-00888 Bipolar Micro Forceps, 5 mm, 31 cm	PSI-L I -00889	PSI-L I -00890	PSI-L I -00891
1/1	PSI-LI-00892 Bipolar Macro Forceps, 5 mm, 31 cm	PSI-L I -00893	PSI-L I -00894	PSI-L I -00895
1/1	PSI-LI-00896 Bipolar Maryland Dissector, fenestrated 5 mm, 31 cm	PSI-L I -00897	PSI-L I -00898	PSI-L I -00899
	ITEM NO. DISPOSABL	E BIPOLAR		
1/1 AESCI	Metzenbaum Scissors, 5 mm, 31 cm Six single sterile scissors per box			
1/1 AESCU	Maryland Grasping Forceps, 5 mm, 31 cm Six single sterile scissors per box			

# **Suturing**

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com



### Durogrip TC Needle Holder, straight, 5 mm PSI-L I -00909 PSI-L I -00911 Durogrip TC Needle Holder, curved left, 5 mm PSI-L I -00910 PSI-L I -00912 Durogrip TC Needle Holder, curved right, 5 mm PSI-L I -00913 "The Viper" Durogrip TC Suture Assistant/Needleholder, PSI-L I -00914 2 x 3 teeth, 5 mm Durogrip TC Needle Holder, fenestrated, 5 mm PSI-L I -00915 Durogrip TC Needle Holder, crimped/fenestrated, 5 mm PSI-L I -00916 PSI-L I -00917 Super Righting Needle Holder, 5 mm

### **KNOT PUSHERS**


	31 cm	
1/1 Knot Pusher, 5 mm	PSI-L I -00918	
1/1 Knot Pusher, 5 mm	PSI-L I -00919	
Knot Pushing/Knot Tightening Forceps, 5 mm*	PSI-L I -00920 PSI-L I -00921 PSI-L I -00922	Jaw Insert Outer Tube Handle
CROCE/OLMI Dissecting and Suture Scissors with Knot Pusher, 5 mm	PSI-L I -00923	

<sup>\*</sup>Instrument sold by component

# **Specialty Instruments**

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com





### **BULLDOG CLAMPS FOR TEMPORARY OCCLUSION**

**CLIP APPLICATORS AND REMOVERS** 

# Articulating applicator for temporary endoscopic vessel clips Shaft length: 34 cm Diameter: 12.5 mm PSI-LI-00931 Fixed applicator for temporary endoscopic vessel clips Shaft length: 34 cm Diameter: 12.5 mm

# WAESCULAR AZIN

ATRAUMATIC ENDO VESSEL CLIPS (for temporary occlusion)

Item No.	Closing force	Jaw length	Jaw shape
PSI-L I -00933	350 g	25 mm	short/straight
PSI-L I -00934	450 g	45 mm	long/straight
PSI-L I -00935	350 g	25 mm	short/curved
PSI-L I -00936	450 g	45 mm	long/curved

# **Accessories**

# PRIME SURGICAL INTERNATIONAL www.primeinstrumentintl.com

### **CABLES & ACCESSORIES**

INSTRUMENT	
	PSI-LI-00937 Bipolar Cable for Aesculap, Wolf, Martin, Berchtold, 12 ft.
	PSI-L I -00938 Bipolar Cable for Valleylab, Conmed, Bovie, 12 ft.
	PSI-L I -00939  Bipolar Cable for Valleylab, Conmed, Bovie, 12 ft. (Disposable, Single-use, 10 units/box)

PSI-L I -00940	Reusable Monopolar Cable, Valleylab, Conmed, Bircher, Bovie, Erbe ICC, 4 mm diam., 10 ft.		
PSI-LI-00941 Cleaning Brush for 5 mm instruments, 12 pkg			
PSI-L I -00942	Cleaning Brush for 10 mm, 6 pkg		
PSI-L I -00943	Wire Brush for cleaning jaws		
PSI-L I -00944	Sterilit Instrument Oil		
PSI-LI-0930	Sterilit Instrument Oil Spray		