

PRIME Surgical International (PSI)

COMMON EXTERNAL FIXATOR

PRIME SURGICAL INTERNATIONAL

COMMON EXTERNAL FIXATOR

PRIME Surgical International (PSI)

COMMON EXTERNAL FIXATOR

PRIME SURGICAL INTERNATIONAL

PRIME Surgical International (PSI)

COMMON EXTERNAL FIXATOR

Common External Fixator System				
Product	REF	Picture	Specification	Accessories
Straight Type External Fixator	9100-0101		L	6*130mm Pin 4 PCS 6mm Hex Wrench 6mm Screwdriver
	9100-0102		M	6*120mm Pin 4 PCS 6mm Hex Wrench 6mm Screwdriver
	9100-0103		S	5*110mm Pin 4 PCS 6mm Hex Wrench 5mm Screwdriver
	9100-0104		XS	4*90mm Pin 4 PCS 4mm Hex Wrench 4mm Screwdriver
T-shaped Type External Fixator	9100-0201		L	6*130mm Pin 2 PCS HB6.0*150mm Pin 2 PCS 6mm Hex Wrench 6mm Screwdriver
	9100-0202		M	6*130mm Pin 2 PCS HB6.0*150mm Pin 2 PCS 6mm Hex Wrench 6mm Screwdriver
	9100-0203		S	5*110mm Pin 2 PCS HB5.0*120mm Pin 2 PCS 5mm Hex Wrench 5mm Screwdriver
	9100-0204		XS	4*90mm Pin 2 PCS 4mm Hex Wrench 4mm Screwdriver
Ankgle Joint External Fixator	9100-0301		L	6*130mm Pin 4 PCS 6mm Hex Wrench 6mm Screwdriver
	9100-0302		M	5*130mm Pin 4 PCS 6mm Hex Wrench 5mm Screwdriver
Straight Type with Middle Pins External Fixator	9100-0401		/	6*130mm Pin 6 PCS 6mm Hex Wrench 6mm Screwdriver
T-shaped Type with Middle Pins External Fixator	9100-0501		/	6*130mm Pin 4 PCS HB6.0*150mm Pin 2 PCS 6mm Hex Wrench 6mm Screwdriver
3D Type External Fixator	9100-0601		L	6*130mm Pin 4 PCS 6mm Hex Wrench 10-14mm Open Spanner
	9100-0602		M	5*120mm Pin 4 PCS 6mm Hex Wrench 10-14mm Open Spanner
	9100-0603		S	5*110mm Pin 4 PCS 5mm Hex Wrench 10-14mm Open Spanner
	9100-0604		XS	4*90mm Pin 4 PCS 4mm Hex Wrench 10mm Open Spanner

PRIME Surgical International (PSI)

COMMON EXTERNAL FIXATOR

Three-section Straight Type External Fixator	9100-0701		/	6*130mm Pin 6 PCS 6mm Hex Wrench 6mm Screwdriver
Three-section T-shaped Type External Fixator	9100-0801		/	6*130mm Pin 4 PCS HB6.0*150mm Pin 2 PCS 6mm Hex Wrench 6mm Screwdriver
Joint Type External Fixator	9100-0901		/	6*130mm Pin 2 PCS HB6.0*180mm Pin 2 PCS 6mm Hex Wrench 6mm Screwdriver
Oblique Type External Fixator	9100-1001		/	6*150mm Pin 2 PCS HB6.0*180mm Pin 2 PCS 6mm Hex Wrench 6mm Screwdriver
Dynamic External Fixator	9100-1101		/	4*250 Pin 2 PCS HB5*250mm Pin 3 PCS 10mm Open Spanner 4&5mm Screwdriver
External Fixator-I	9100-1201		L	6*130mm Pin 4 PCS 12-16mm Open Spanner 4&5mm Hex Wrench 6mm Screwdriver
	9100-1202		M	5*120mm Pin 4 PCS 10-14mm Open Spanner 4&5mm Hex Wrench 5mm Screwdriver
	9100-1203		S	5*110mm Pin 4 PCS 10-14mm Open Spanner 4&5mm Hex Wrench 5mm Screwdriver
External Fixator-II	9100-1301		L	6*130mm Pin 4 PCS 10-16mm Open Spanner 4mm Hex Wrench 6mm Screwdriver
	9100-1302		M	5*120mm Pin 4 PCS 10-14mm Open Spanner 4mm Hex Wrench 5mm Screwdriver
	9100-1303		S	5*110mm Pin 4 PCS 10-14mm Open Spanner 4mm Hex Wrench 5mm Screwdriver

PRIME Surgical International (PSI)

COMMON EXTERNAL FIXATOR

Mini Fragment External Fixator	9100-1401		/	/
Mini Fragment External Fixator	9100-1402		/	/
Wrist Joint External Fixator	9100-1501		/	/
Semi-circular External Fixator	9100-1601		/	Rod-rod Coupling 8/8mm 6PCS Pin-rod Coupling 8/4-5mm 6PCS Semi-circular Curved Rod 8mm Large 2PCS Connecting Rod 8*300mm 1PC Compression/distraction Rod 8mm 2PCS Cortical Pin 4*120mm 2PCS K-wire 4*250mm 2PCS
Pelvic External Fixator	9100-1701		/	Rod-rod Coupling 8/8mm 6PCS Pin-rod Coupling 8/4-5mm 6PCS Semi-circular Curved Rod 8mm Large 2PCS Connecting Rod 8*150mm 2PCS Compression/distraction Rod 8mm 1PCS Cortical Pin 5*130mm 6PCS
Pertrochanteric External Fixator	9100-1801		/	Double Rod-pin Clamp 8/4-5mm 5PCS Connecting Rod 8*300mm 2PCS Cortical Pin 5*130mm 2PCS Cortical Pin 5*180mm 3PCS

PRIME Surgical International (PSI)

COMMON EXTERNAL FIXATOR

Femur & Tibial External Fixator	9100-1901		/	Double Rod-pin Clamp 8/4-5mm 4PCS Connecting Rod 8*250mm 2PCS Cortical Pin 5*130mm 4PCS
Humeral Lateral Condylar Neck External Fixator	9100-2001		/	Pin-rod Coupling 8/4-5mm 4PCS Connecting Rod 8*200mm 1PC Connecting Rod 8*150mm 1PC Cortical Pin 5*130mm 4PCS
Humeral External Fixator	9100-2101		/	Rod-rod Coupling 8/8mm 2PCS Pin-rod Coupling 8/4-5mm 6PCS Connecting Rod 8*250mm 2PCS Connecting Rod 8*200mm 1PCS Cortical Pin 5*130mm 6PCS
Femoral Condylar External Fixator	9100-2201		/	Rod-rod Coupling 8/8mm 3PCS Double-pin Coupling 8/4-5mm 2PCS Pin-rod Coupling 8/4-5mm 5PCS Semi-circular Curved Rod 8mm Large 1PC Connecting Rod 8*300mm 2PCS Curved Rod 8mm 1PC Cortical Pin 5*130mm 7PCS
Femoral External Fixator	9100-2301		/	Rod-rod Coupling 8/8mm 3PCS Double-pin Coupling 8/4-5mm 4PCS Pin-rod Coupling 8/4-5mm 3PCS Semi-circular Curved Rod 8mm Large 2PCS Connecting Rod 8*250mm 2PCS Cortical Pin 5*130mm 4PCS K-wire 4*250mm 1PC

PRIME Surgical International (PSI)

COMMON EXTERNAL FIXATOR

Ankle External Fixator	9100-2401		/	Rod-rod Coupling 8/8mm 6PCS Pin-rod Coupling 8/4-5mm 8PCS Compression/distraction Rod 8mm 2PCS Connecting Rod 8*350mm 2PCS Connecting Rod 8*250mm 2PCS Cortical Pin 4*120mm 4PCS K-wire 4*250mm 2PCS
Proximal Tibial External Fixator	9100-2501		/	Rod-rod Coupling 8/8mm 3PCS Pin-rod Coupling 8/4-5mm 8PCS Connecting Rod 8*300mm 1PC Semi-circular Curved Rod 8mm Large 1PC Curved Rod 8mm 1PC Cortical Pin 4*120mm 3PCS K-wire 4*250mm 1PC
Humeral External Fixator	9100-2601		/	Rod-rod Coupling 8/8mm 3PCS Pin-rod Coupling 8/4-5mm 5PCS Connecting Rod 8*250mm 1PC Semi-circular Curved Rod 8mm Large 1PC Curved Rod 8mm 1PC Cortical Pin 4*120mm 3PCS K-wire 4*250mm 1PC
Knee Joint External Fixator	9100-2701		/	Rod-rod Coupling 8/8mm 8PCS Pin-rod Coupling 8/4-5mm 14PCS Semi-circular Curved Rod 8mm Large 2PCS Curved Rod 8mm 2PC Spanning Joint Connecting Rod 8mm 2PCS Cortical Pin 5*130mm 14PCS