

PRIME Surgical International (PSI)

SPINAL SYSTEM INSTRUMENT SET

PRIME SURGICAL INTERNATIONAL

SPINAL SYSTEM INSTRUMENT SET


PRIME Surgical International (PSI)

SPINAL SYSTEM INSTRUMENT SET

PRIME SURGICAL INTERNATIONAL



PRIME Surgical International (PSI)

SPINAL SYSTEM INSTRUMENT SET

Spinal Instrument Set						
Name	Picture		NO.	REF	Description	Qty.
2200-01 Spinal Pedicle Screw System Instrument Set			1	2200-0101	Trail Rod	1
			2	2200-0102	Screwdriver for Polyaxial Screw	1
			3	2200-0103	Screwdriver for Monoaxial Screw	1
			4	2200-0104	In-situ Bending Iron L	1
			5	2200-0105	In-situ Bending Iron R	1
			6	2200-0106	Crosslink Nut Holder Hex	1
			7	2200-0107	Screw Nut Holder Hex	1
			8	2200-0108	Screw Cutter for Long Arm Screw	1
			9	2200-0109	Feeler for Screw Channel Bent	1
			10	2200-0110	Tap $\phi 4.0$	1
			11	2200-0111	Tap $\phi 5.5$	1
			12	2200-0112	Screwdriver Hex 3.5mm for Crosslink Nut	1
			13	2200-0113	Tap $\phi 6.0$	1
			14	2200-0114	Tap $\phi 7.0$	1
			15	2200-0115	Hex Screwdriver for Srew Nut	1
			16	2200-0116	Feeler for Screw Channel Straight	1
			17	2200-0117	Rod Twist	1
			18	2200-0118	Parallel Compressor	1
			19	2200-0119	Parallel Spreader	1
			20	2200-0120	Rod Pusherial	1
			21	2200-0121	Insert Device for Fixation Pin	1
			22	2200-0122	Counter Torque for Screw Cutter	1
			23	2200-0123	Counter Torque	1
			24	2200-0124	Rod Bender	1
			25	2200-0125	Straight Handle Quick Coupling	1
			26	2200-0126	Rod Pushing Forcecp	1
			27	2200-0127	Rod Holding Forcecp	1
			28	2200-0128	T-handle Quick Coupling	1
			29	2200-0129	AWL	1
			30	2200-0130	Pedicle Probe Bent	1
			31	2200-0131	Pedicle Probe Straight	1
			32	2200-0132	Almunium Box	1
12 N.m Torque Wrench		/	/	2200-0133	/	1
Reduction Forcecp		/	/	2200-0134	/	1
T-handle Rechat Screwdriver		/	/	2200-0135	/	1




PRIME Surgical International (PSI)

SPINAL SYSTEM INSTRUMENT SET

2200-02 Anterior Cervical Plate System Instrument Set 	1	2200-0201	Plate and Screw Box	1
	2	2200-0202	AWL	1
	3	2200-0203	Tap	1
	4	2200-0204	Drill Bit	1
	5	2200-0205	Drill Bit with Straight Handle	1
	6	2200-0206	Hex Screwdriver	1
	7	2200-0207	Hex Screwdriver	1
	8	2200-0208	Nut Screwdriver	1
	9	2200-0209	Curette	1
	10	2200-0210	Brace Screwdriver	1
	11	2200-0211	Distractor	1
	12	2200-0212	Guider	1
	13	2200-0213	Location Screw Screwdriver	1
	14	2200-0214	Plate Holding Forcep	1
	15	2200-0215	Screwdriver Sleeve	1
	16	2200-0216	Screwdriver Sleeve	1
	17	2200-0217	Cervical Hook	1
	18	2200-0218	Plate Bender	1
	19	2200-0219	Location Screw	4
	20	2200-0220	Almunium Box	1
2200-03 Occipito-Cervical Fusion System Instrument Set 	1	2200-0301	Guider	1
	2	2200-0302	Screwdriver Hexagonal 2.5mm	1
	3	2200-0303	Rod Pushing Forcep	1
	4	2200-0304	Straight Handle Quick Coupling	1
	5	2200-0305	Insert Device for Fixation Pin	1
	6	2200-0306	Counter Torque	1
	7	2200-0307	Screwdriver for Polyaxial Screw	1
	8	2200-0308	Rod Bender	1
	9	2200-0309	AWL	1
	10	2200-0310	Protection Sleeve	1
	11	2200-0311	AWL	1
	12	2200-0312	Depth Gauge 0-40mm	1
	13	2200-0313	Drill Bit 2.7mm	1
	14	2200-0314	Drill Bit 2.4mm	1
	15	2200-0315	Trail Rod	1
	16	2200-0316	Screwdriver Hexagonal 2.5mm	1
	17	2200-0317	Tap ϕ 4.0	1
	18	2200-0318	Tap ϕ 3.5	1
	19	2200-0319	Nut Holder	1
	20	2200-0320	Feeler for Screw Channel Bent	1
	21	2200-0321	Feeler for Screw Channel Straight	1
	22	2200-0322	Distractor	1
	23	2200-0323	In-situ Bending Iron L	1
	24	2200-0324	In-situ Bending Iron R	1
	25	2200-0325	T-handle Quick Coupling	1
	26	2200-0326	Hook Holding Forcep	1
	27	2200-0327	Rod Holding Forcep	1
	28	2200-0328	Rod Twist	1
	29	2200-0329	Parallel Compression Forcep	1
	30	2200-0330	Rod Cutter	1
	31	2200-0331	Distraction Forcep	1
	32	2200-0332	Almunium Box	1


PRIME Surgical International (PSI)

SPINAL SYSTEM INSTRUMENT SET

2200-4 Instrument Set for Anterior Thoracic Plate System and Anterior Thoracolumbar Plate System 	1	2200-0401	Curved Guide Pin	1
	2	2200-0402	Straight Guide Pin	1
	3	2200-0403	Drill Bit	1
	4	2200-0404	Obturator	1
	5	2200-0405	Nut Driver SW8	1
	6	2200-0406	Nut Driver SW8	1
	7	2200-0407	Drill Guider	1
	8	2200-0408	Detector	1
	9	2200-0409	Hole Enlarger	1
	10	2200-0410	Tap 6.5mm	1
	11	2200-0411	Screw Holding Driver	1
	12	2200-0412	Tap 6mm	1
	13	2200-0413	Screwdriver	1
	14	2200-0414	Universal Screwdriver	1
	15	2200-0415	Screwdriver	1
	16	2200-0416	Calipers	1
	17	2200-0417	Compression Clamp 3mm	1
	18	2200-0418	Distractor Clamp	1
	19	2200-0419	Plate Holding Clamp	1
	20	2200-0420	Holding Bolt SW8	1
	21	2200-0421	Combination Wrench	1
	22	2200-0422	Wrench SW11	1
2200-5 Instrument Set for Titanium Mesh Cage 	1	2200-0501	Cage Stand	1
	2	2200-0502	Pressure 6mm	1
	3	2200-0503	Pressure 18mm	1
	4	2200-0504	Pusher Curved	1
	5	2200-0505	Osteotribe	1
	6	2200-0506	Pressure 12mm	1
	7	2200-0507	Pusher Straight	1
	8	2200-0508	Implant Measure 10-12mm	1
	9	2200-0509	Implant Measure 16-18mm	1
	10	2200-0510	Implant Measure 22-25mm	1
	11	2200-0511	Cage Holding Forcep	1
	12	2200-0512	Cage Cutter	1
	13	2200-0513	Aluminium Box	1
2200-06 Peek Cervical Cage Instrument Set 	1	2200-0601	Bone Graft Inserter	1
	2	2200-0602	Trail Holder	1
	3	2200-0603	Peek Cage Holder	1
	4	2200-0604	Bone Graft Cage	1
	5	2200-0605	Cage Box	1
	6	2200-0606	Trail 4mm	1
	7	2200-0607	Trail 5mm	1
	8	2200-0608	Trail 6mm	1
	9	2200-0609	Trail 7mm	1
	10	2200-0610	Trail 8mm	1
	11	2200-0611	Trail 9mm	1
	12	2200-0612	Trail 10mm	1
	13	2200-0613	Aluminium Box	1


PRIME Surgical International (PSI)

SPINAL SYSTEM INSTRUMENT SET

2200-07 PLIF PEEK Lumbar Cage & Posterior Lumbar Interbody Fusion Cage Instrument Set 	1	2200-0701	Trail 20*8mm	1
	2	2200-0702	Trail 20*10mm	1
	3	2200-0703	Trail 20*12mm	1
	4	2200-0704	Trail 20*14mm	1
	5	2200-0705	Trail 22*8mm	1
	6	2200-0706	Trail 22*10mm	1
	7	2200-0707	Trail 22*12mm	1
	8	2200-0708	Trail 22*14mm	1
	9	2200-0709	Trail 26*8mm	1
	10	2200-0710	Trail 26*10mm	1
	11	2200-0711	Trail 26*12mm	1
	12	2200-0712	Reamer 7mm	1
	13	2200-0713	Reamer 8mm	1
	14	2200-0714	Reamer 9mm	1
	15	2200-0715	Bone Graft Insert	1
	16	2200-0716	Hammer	1
	17	2200-0717	Trail 26*14mm	1
	18	2200-0718	Trail 32*8mm	1
	19	2200-0719	Trail 32*10mm	1
	20	2200-0720	Trail 32*12mm	1
	21	2200-0721	Trail 32*14mm	1
	22	2200-0722	Distracter	1
	23	2200-0723	Bone Graft Cage	1
	24	2200-0724	Cage Holder	1
	25	2200-0725	Curette	1
	26	2200-0726	Nerve Dissector	1
	27	2200-0727	Reamer 8*6.5mm	1
	28	2200-0728	Reamer 8.5*7.5mm	1
	29	2200-0729	Reamer 9.5*8.5mm	1
	30	2200-0730	Reamer 10*10.5mm	1
	31	2200-0731	Reamer 10.5*11.5mm	1
	32	2200-0732	Aluminium Box	1
2200-08 TLIF PEEK Cage Instrument Set 	1	2200-0801	Distracter 7mm	1
	2	2200-0802	Distracter 9mm	1
	3	2200-0803	Pusher Straight	1
	4	2200-0804	Pusher Bent	1
	5	2200-0805	Distracter 11mm	1
	6	2200-0806	Distracter 13mm	1
	7	2200-0807	Trail 7*25mm	1
	8	2200-0808	Trail 9*25mm	1
	9	2200-0809	Trail 11*25mm	1
	10	2200-0810	Trail 13*25mm	1
	11	2200-0811	Trail 15*25mm	1
	12	2200-0812	Trail 15*25mm	1
	13	2200-0813	Hook	1
	14	2200-0814	T-handle Quick Coupling	1
	15	2200-0815	Ring Curette	1
	16	2200-0816	Osteotome	1
	17	2200-0817	Cage Holder	1
	18	2200-0818	Hammer	1
	19	2200-0819	Curette Left	1
	20	2200-0820	Curette Right	1
	21	2200-0821	Scraper Large	1
	22	2200-0822	Scraper Small	1
	23	2200-0823	Aluminium Box	1

PRIME Surgical International (PSI)

SPINAL SYSTEM INSTRUMENT SET

2200-09 Laminar Shaping Plate Instrument Set		1	2200-0901	Plate Box	1
		2	2200-0902	Hand Drill Bit 1.5*4mm	1
		3	2200-0903	Hand Drill Bit 1.5*6mm	1
		4	2200-0904	Hand Drill Bit 1.5*8mm	1
		5	2200-0905	Hand Drill Bit 1.5*10mm	1
		6	2200-0906	Hand Drill Bit 1.5*12mm	1
		7	2200-0907	Mold Trial 4mm	1
		8	2200-0908	Mold Trial 6mm	1
		9	2200-0909	Mold Trial 8mm	1
		10	2200-0910	Mold Trial 10mm	1
		11	2200-0911	Mold Trial 12mm	1
		12	2200-0912	Bending Cutting Forcep	1
		13	2200-0913	Bending Cutting Forcep	1
		14	2200-0914	Plate Holding Forcep	1
		15	2200-0915	Plate Holding Forcep	1
		16	2200-0916	Plate Holding Clamp	1
		17	2200-0917	AWL	1
		18	2200-0918	Screwdriver	1
		19	2200-0919	Screwdriver	1
		20	2200-0920	Almunium Box	1