zia- ur-Rahman

Catalog of Strelization Boxes



Zia ur rahman- Address: Office No. 25 , Queen Centre M.T Khan Road Karachi- Paksitan e.mail: zia@primeinstrumentintl.com; prime.instruments@hotmail.com; Cell: 0092- 300-2621673, 0092-313-5710887


PRIME Containers.....Designed and Constructed for the Customer

- Container lid and bottom are both of one-piece construction and have a uniform, anodized aluminum alloy finish
- Lightweight, yet durable enough to accommodate acceptable weight loads of surgical instrumentation
- . Handles stop at 90° to aid lifting; low profile latch allows simple and secure lid closing
- . Retention plate employs an easy "click and lock" closure over disposable, quick-to-change filter
- Dual silicone gaskets within retention plate and around perimeter ensure complete Container and filter seal
- Basket offering includes both solid wire mesh and perforated metal designs to allow for maximum ultrasonic cleaning;

handles can be positioned to assist with secure and comfortable grasp for removal under sterile conditions



Zia ur rahman- Address: Office No. 25 , Queen Centre M.T Khan Road Karachi- Paksitan e.mail: zia@primeinstrumentintl.com; prime.instruments@hotmail.com; Cell: 0092- 300-2621673, 0092-313-5710887



LIMITED WARRENTY

PRIME hereby warrants all PRIME Containers purchased directly from PRIME or an Authorized Dealer in accordance with the following terms and conditions:

The PRIME Container is guaranteed to be free from defects in material and workmanship when used for its intended purpose. Any PRIME Container that proves defective in workmanship or material will either be repaired or replaced, at PRIME discretion, without charge. Containers not purchased from an authorized source may not be authentic or meet PRIME quality standards.

Container Gaskets Warranty

The warranty period for PRIME container gaskets is three years from the date of purchase.

Periodic Maintenance:

Periodic Maintenance is suggested to ensure all parts are in good working order.

Repair or Replacement:

Under this Warranty, and at its discretion,

PRIME will repair or replace any parts of the product found to be defective in workmanship or material.

PRIME will determine, at its discretion, the cause and nature of any defect, the necessity and manner of repair or replacement, and all other matters pertaining to the conditions of the PRIME Container System. PRIME, Inc. Returned Goods Policy applies to all returned items (refer to details in current price list).

Exclusions of warranty: This Warranty shall not apply to any condition(s) or damage resulting from negligence, misuse, improper handling, improper cleaning, improper opening techniques, or unauthorized repair work including but not limited to: Incidents of abuse such as denting of the PRIME Container due to dropping or other instances of mechanically applied pressure. Damage determined to be related to caustic or abrasive

cleaning agents. . Items modified by the customer. . Items modified or customized by PRIME at the request of

the customer. . Damage from fire, flood and other occurrences not under the control of PRIME, Inc.

Zia ur rahman- Address: Office No. 25, Queen Centre M.T Khan Road Karachi- Paksitan e.mail: zia@primeinstrumentintl.com; prime.instruments@hotmail.com; Cell: 0092- 300-2621673, 0092-313-5710887

ADVANTAGES



EASY HANDLING

- · Effortless removal and visual examination of the sterile barrier no filter exchange necessary
- · Trouble-free cleaning due to optimal microstructure of the surface



PROCESS SECURITY

- · Security in the reprocessing cycle aided by color indicators and one-way seals
- Simplified logistics with colored logistics framelets and label displays for record labels on top surface and front sides

zia- ur-Rahman

					WEIGHT
	CODE	ORDER NO.	DESCRIPTION	COLOR	kg
	1/1	K1100V	Lid perforated-Bottom non perforated, L23"X W11"X H4.5'	silver/gold/red/blue/green/black	3.55
		K1100VB	Lid perforated-Bottom perforated L23"X W11"X H5.3"		4.55
			Lid perforated-Bottom non perforated, L23"X W11"X H5.5		3.90
		K1135VB Lid perforated-Bottom perforated, L23"X W11" X H6.4"	311Verrgora/rea/blackgreen/black	5.05	
			Lid perforated-Bottom non perforated, L23"X W11"XH6.1	silver/gold/red/blue/green/black	4.00
			Lid perforated-Bottom perforated, L23"X W11"X H6.8"	311Ve17gOld/Ted/black	5.25
			Lid perforated-Bottom non perforated, L23"X W11"X H8"	silver/gold/red/blue/green/black	4.45
		K1200VB	Lid perforated-Bottom perforated, L23"X W11"X H9"	311Ve17gOld/1ed/blde/gleerl/black	5.95
			Lid perforated-Bottom non perforated, L23"X W11" X H10"	silver/gold/red/blue/green/black	4.95
		KI1260VB	Lid perforated-Bottom perforated, L23"X W11"X H11"	311Ve17gOld/Ted/black	6.75
П	3/4	K2100V	Lid perforated-Bottom non perforated, L18"X W11"X H4"		2.79
_		K2100VB	Lid perforated-Bottom perforated, L18"X W11"X H5.3"	silver/gold/red/blue/green/black	3.24
l S		K2135V	Lid perforated-Bottom non perforated, L18"X W11"X H5.5"	silver/gold/red/blue/green/black	2.94
PAR I		K2135VB	Lid perforated-Bottom perforated, L18"X W11"X H6.2"		3.39
R			Lid perforated-Bottom non perforated, L18"X W11"X H6.1"	silver/acid/sed/blve/areen/block	3.04
BIO BARRIER MODEL		K2150VB	Lid perforated-Bottom perforated, 118"X W11"X H7"		3.49
DE	1/2	K3100V	Lid perforated-Bottom non perforated, L11.2"X W11"X H4		2.12
۳		K3100VB	Lid perforated-Bottom perforated,L11.2"X W11"X H5.3"	silver/gold/red/blue/green/black	2.57
			Lid perforated-Bottom non perforated L11.2"X W11"X H5.5		2.22
		K3135VB	Lid perforated-Bottom perforated, L11.2"X W11"X H6.2"		2.67
		K3150V	Lid perforated-Bottom non perforated, L11.2"X W11"X H6.1	111	2.32
		K3150V K3150VB	Lid perforated-Bottom perforated, L11.2"X W11"X H7"	silver/gold/red/blue/green/black	2.77
		K3200V	Lid perforated-Bottom perforated, L11.2 × W11 × 117 Lid perforated-Bottom non perforated, L11.2 × W11 × H18		2.52
		K3200VB	Lid perforated-Bottom perforated, L11.2"X W11"X H9"	silver/gold/red/blue/green/black	2.98
		K3260VB		" silver/gold/red/blue/green/black	
			Lid perforated-Bottom non perforated L11.2"X W11"X H10		2.98
		K3260VB	Lid perforated-Bottom perforated, L11.2"X W11"X H11"		3.42
	XL SIZE	K7120V	Lid perforated-Bottom non perforated, L23"X W15"X H4.9'	silver/gold/red/blue/green/black	4.23
		K7120VB	Lid perforated-Bottom perforated, L23"X W15"X H5.3"		4.65
		K7180V	Lid perforated-Bottom non perforated, L23"X W15"X H7.2"		5.03
		K7180VB	Lid perforated-Bottom perforated, L23"X W15"X H7.6"	silver/gold/red/blue/green/black	5.47
$\overline{}$		10000	Lid politica Dottolii politicated, E23 x 1113 x 117.0		U. T/

zia- ur-Rahman

	FLAT CONTAINERS						
	CODE	ORDER NO.	DESCRIPTION	COLOR	WEIGHT kg		
		K4055R	Lid perforated-Bottom non perforated, L11.2"X W11"X H2'	silver/gold/red/blue/green/black	11.25		
		K4055RB	Lid perforated-Bottom perforated, L11.2"X W11"X H2"		1.35		
		K4085R	Lid perforated-Bottom non perforated,L11.2"X W11" X H3'	silver/gold/red/blue/green/black	11.35		
		K4085RB	Lid perforated-Bottom perforated, L11.2"X W11"X H3"	311VC17gOld/1Cd/bldc/glcc/1/bldck	1.55		
			MINI CONTAINERS				
	CODE	ORDER NO.	DESCRIPTION	COLOR	WEIGHT kg		
		K5040R	Lid perforated-Bottom non perforated, L12"X W5"X H1"	silver/gold/red/blue/green/black	0.55		
		K5040RB	Lid perforated-Bottom perforated, L12"X W5"X H1"	silver/gold/red/blue/green/black	0.65		
	2000	K5070R	Lid perforated-Bottom non perforated, L12"X W5"X H3"	silver/gold/red/blue/green/black	0.70		
_{DD}		K5070RB	Lid perforated-Bottom perforated,L12"X W5"X H3"		0.75		
ASI			Lid perforated-Bottom non perforated, L12"X W5"X H4"		0.80 0.85		
S ≥		K5100RB	Lid perforated-Bottom perforated L12"X W5"X H4"		0.65		
BASIC MOEDL	DENTAL CONTAINERS						
무	CODE	ORDER NO.	DESCRIPTION	COLOR	WEIGHT kg		
		K6050R	Lid perforated-Bottom non perforated, L12.5"X W7" X H2"	silver/gold/red/blue/green/black	0.85		
		K6050RB	Lid perforated-Bottom perforated, L12.5"X W7"X H2"	silver/gold/red/blue/green/black	0.95		
		K6065R	Lid perforated-Bottom non perforated, L12.5"X W7"X H2.5	silver/gold/red/blue/green/black silver/gold/red/blue/green/black	0.90		
		K6065RB	Lid perforated-Bottom perforated, L12.5"X W7"X H2.5"		1.00		
			Lid perforated-Bottom non perforated, L12.5"X W7"X H5"		1.15		
		K6130RB	Lid perforated-Bottom perforated, L12.5"X W7"X H5"		1.20		
	XL SIZE						
	CODE	ORDER NO.	DESCRIPTION	COLOR	WEIGHT		
	139	K7120R	Lid perforated-Bottom non perforated L23"X W14"X H4.7"	silver/gold/red/blue/green/black	3.85		
		K7120RB	Lid perforated-Bottom perforated, L23"X W14"X H4.7"		4.15		
		K7180R	Lid perforated-Bottom non perforated, L23"X W14"X H7"		4.65		
		K7180RB	Lid perforated-Bottom perforated, L23"X W14"X H7"	silver/gold/red/blue/green/black	4.90		

zia- ur-Rahman



sterile container and accessories



					_
	CODE	ORDER NO.	DESCRIPTION	COLOR	WEIGHT
		K1100R K1100RB	Lid perforated-Bottom non perforated_L23"X W11"X H 4" Lid perforated-Bottom perforated.L23"X W11" X H4"	silver/gold/red/blue/green/black	3.05
		K1135R K1135RB	Lid perforated-Bottom non perforated, 123"X W11" X H5" Lid perforated-Bottom perforated, 123" X W11" X H5"	silver/gold/red/blue/green/black	3.49 3.65
		K1150R K1150RB	Lid perforated-Bottom non perforated, L23" X W11" X H6" Lid perforated-Bottom perforated, L23" X W11" X H6"	silver/gold/red/blue/green/black	3.50 3.75
			Lid perforated-Bottom non perforated, L23" X W11" X H7" Lid perforated-Bottom perforated, L23" X W11" X H7"	silver/gold/red/blue/green/black	3.95
			Lid perforated-Bottom non perforated, L23" X W11" X H10	"silver/gold/red/blue/green/black	4.45
BASIC	3/4	K2100R	Lid perforated-Bottom perforated, 123" X W11" X H10" Lid perforated-Bottom non perforated, L18" X W11" X H4"	silver/gold/red/blue/green/black	2.45
SIC M		K2100RB K2135R	Lid perforated-Bottom perforated, L18" X W11" X H4" Lid perforated-Bottom non perforated, L18" X W11" X H5"		2.60
MOEDL		K2135RB	Lid perforated-Bottom perforated, L18" X W 11" X H5"	silver/gold/red/blue/green/black	2.75
		K2150R K2150RB	Lid perforated-Bottom non perforated, L18" X W11" X H6" Lid perforated-Bottom perforated, L18" X W11" X H6"	silver/gold/red/blue/green/black	2.70 2.90
	1/2	K3100R K3100RB	Lid perforated-Bottom non perforated, L11" X W11" X H 4 Lid perforated-Bottom perforated, L11" X W11" X H4"	silver/gold/red/blue/green/black	1.85 1.95
			Lid perforated-Bottom non perforated, L11" X W11" X H5" Lid perforated-Bottom perforated, L11" X W11" X H5"	silver/gold/red/blue/green/black	1 95
			Lid perforated-Bottom non perforated, L11" X W11" X H6"	silver/gold/red/blue/green/black	2.05
		K3200R	Lid perforated-Bottom perforated, L11" X W11" X H6" Lid perforated-Bottom non perforated, L11" X W11" X H7" Lid perforated Pattern perforated, L11" X W11" X H7"	silver/gold/red/blue/green/black	2.25
		K3200RB K3260R	Lid perforated-Bottom perforated, L11" X W11" X H7" Lid perforated-Bottom non perforated, L11" X W11" X H10	"silver/gold/red/blue/green/black	2.50
		K3260RB	Lid perforated-Bottom perforated, L11" X W11" X H10"		2.80

The Bar Basic Ser	K6203	Coloured identification labels. Engraved to your specifications. Maximum 16 ltters
	K6451	Rack,for endoscopy Holds 10 endoscopic instruments Fits into a 1/1 size 20cm deep container
	K6452	Rack,for arthroscopy Holds 20 arthroscopic punches Fits into a 1/1 size 20cm deep container
	K 6453	Rack,for Kerrison Holds 10 Kerrison punches Fits into a 1/1, 3/4, 1/2 size 20cm deep container

Zia ur rahman- Address: Office No. 25 , Queen Centre M.T Khan Road Karachi- Paksitan e.mail: zia@primeinstrumentintl.com; prime.instruments@hotmail.com; Cell: 0092- 300-2621673, 0092-313-5710887

	CODE	ORDER NO.	DESCRIPTION
		K6411	Wire Basket 540X250X30mm for 1/1 Containers
		K6412	Wire Basket 540X250X50mm for 1/1 Containers
		K6413	Wire Basket 540X250X70mm for 1/1 Containers
		K6414	Wire Basket 540X250X100mm for 1/1 Containers
		K6421	Wire Basket 405X250X30mm for 3/4 Containers
WIRE		K6422	Wire Basket 405X250X50mm for 3/4 Containers
WIRE BASKET		K6423	Wire Basket 405X250X70mm for 3/4 Containers
Ē		K6424	Wire Basket 405X250X100mm for 3/4 Containers
		K6431	Wire Basket 255X250X30mm for 1/2 Containers, Flat Containers
		K6432	Wire Basket 255X250X50mm for 1/2 Containers, Flat Containers
		K6433	Wire Basket 255X250X70mm for 1/2 Containers
		K6434	Wire Basket 255X250X100mm for 1/2 Containers
		K6441	Wire Basket 274X172X43mm for Dental Containers
		K6442	Wire Basket 267X125X50mm for Mini Containers

	K6501	Silicone Gasket for mini container 300X140mm
	K6502	Silicone Gasket for dental container 310X190mm
	K6601	Filter retainer, for 1/1, 3/4, 1/2, XL. For disposable paperfilter
	K6602	Filter retainer, for dental model
	K6603	Filter retainer, for mini model

zia- ur-Rahman

	CODE	ORDER NO.	DESCRIPTION		
	0	K6001	Single use paper filter, 100Pcs/Package, for 1/1,1/2,3/4 Containers		
		K6011	Teflon filter number of sterilization cycle>1200, for 1/1,1/2,3/4 Containers		
		K6002	Single use paper filter, 100Pcs/Package, for Mini Containers		
		K6022	Teflon filter number of sterilization cycle>1200, for Mini Containers		
		K6003	Single use paper filter, 100Pcs/Package, for Dental Containers		
	*	K6033	Teflon filter number of sterilization cycle>1200, for Dental Containers		
ACCE		K 6101	Plastic Security Seal 100Pcs/Package		
ACCESSORIES	THE	K 6201	Paper label with indicator 100Pcs/Package, for 1/1,1/2,3/4 Containers		
	DISTANT TO STANT THE STANT	K6202	Paper label with indicator 100Pcs/Package, for Mini & Dental Containers		
		K6301	Silicone Mat 520X230mm for 1/1 Containers		
		K6302	Silicone Mat 380X230mm for 3/4 Containers		
		K6303	Silicone Mat 250X235mm for 1/2, Flat Containers		
		K6304	Silicone Mat 260X180mm for Dental Containers		
		K6305	Silicone Mat 275X125 for Mini Containers		
	1/1	K6306	Silicone Handle in grey/yellow/red/blue/green/black		