

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

PRIME SURGICAL INTERNATIONAL

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

PRIME SURGICAL INTERNATIONAL

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

Styker Hoffmann Type					
Name	Picture	No.	REF	Products	Qty.
7100-01 Hoffmann External Fixation Small Fragment Set		1	7100-0101	Rod-Rod Coupling (5mm Rod)	12
		2	7100-0102	Rod-Pin Coupling (3/4mm Pin, 5mm Rod)	12
		3	7100-0103	4-Hole Pin Clamp (3/4mm Pin)	3
		4	7100-0104	Pin Clamp (3/4mm Pin)	1
		5	7100-0105	Straight Post	2
		6	7100-0106	30 Degree Angled Post	4
		7	7100-0107	Apex Pin Self-drilling 3*80mm	4
		8	7100-0108	Apex Pin Self-drilling 3*100mm	4
		9	7100-0109	Apex Pin Self-drilling 4*80mm	4
		10	7100-0110	Apex Pin Self-drilling 4*100mm	4
		11	7100-0111	Apex Pin Self-drilling 4*120mm	4
		12	7100-0112	Drill Sleeve 3mm	1
		13	7100-0113	Drill Sleeve 4mm	1
		14	7100-0114	Carbon Fiber Connecting Rod 5*120mm	4
		15	7100-0115	Carbon Fiber Connecting Rod 5*150mm	4
		16	7100-0116	Carbon Fiber Connecting Rod 5*180mm	2
		17	7100-0117	Carbon Fiber Connecting Rod 5*200mm	4
		18	7100-0118	T-shaped Clamp Wrench 5#	1
		19	7100-0119	Reduction Wrench	1
		20	7100-0120	Elbow Joint 5mm	1
		21	7100-0121	Apex Pin Chuck 3mm	1
		22	7100-0122	Apex Pin Chuck 4mm	1
		23	7100-0123	T-shaped Pin Wrench 3mm	1
		24	7100-0124	T-shaped Pin Wrench 4mm	1
		25	7100-0125	Arc Drilling Stand	1
		26	7100-0126	Adjust Wheel	1
		27	7100-0127	Aluminum Box	1
Hoffmann External Fixation Small Fragment Set Individuals		/	7100-0102	Rod-Pin Coupling (3/4mm Pin, 5mm Rod)	1
		/	7100-0101	Rod-Rod Coupling (5mm Rod)	1
		/	7100-0104	Pin Clamp (3/4mm Pin)	1
		/	7100-0103	4-Hole Pin Clamp (3/4mm Pin)	1
		/	7100-0127	Rod-Pin Coupling (2mm Pin, 3mm Rod)	1
		/	7100-0128	Rod-Rod Coupling (3mm Rod)	1

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

Hoffmann External Fixation Small Fragment Set Individuals		/	7100-0129	Tube-Rod Coupling 15/5mm	1
		/	7100-0105	Straight Post 5mm	1
		/	7100-0106	30 Degree Angled Post 5mm	1
		/	7100-0130	Aluminum Connecting Rod 5*200mm	1
			7100-0131	Aluminum Connecting Rod 5*250mm	1
			7100-0132	Aluminum Connecting Rod 5*300mm	1
		/	7100-0114	Carbon Fiber Connecting Rod 5*120mm	1
			7100-0115	Carbon Fiber Connecting Rod 5*150mm	1
			7100-0116	Carbon Fiber Connecting Rod 5*180mm	1
			7100-0117	Carbon Fiber Connecting Rod 5*200mm	1
		/	7100-0601	Wrist Joint 5mm	1
		/	7100-0120	Elbow Joint 5mm	1
		/	7100-1101	Compression Distraction Tube 15mm	1
		/	7100-0119	Reduction Wrench	1
		/	7100-0133	Wrench	1
		/	7100-0125	Arc Drilling Stand	1
		/	7100-0118	T-shaped Clamp Wrench 5#	1
		/	7100-0123	T-shaped Pin Wrench 3mm	1
		/	7100-0124	T-shaped Pin Wrench 4mm	1
		/	7100-0112	Drill Sleeve 3mm	1
		/	7100-0113	Drill Sleeve 4mm	1
		/	7100-0121	Apex Pin Chuck 3mm	1
		/	7100-0122	Apex Pin Chuck 4mm	1
		/	7100-0126	Adjust Wheel	1


PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

Ulna & Radius External Fixator 7100-03		/	7100-0301	Rod-Rod Coupling (5mm Rod)	4
		/	7100-0302	4-Hole Pin Clamp (3/4mm Pin)	2
		/	7100-0303	Post 5mm/30 Degree	4
		/	7100-0304	Aluminum Connecting Rod 5*200mm	2
		/	7100-0305	Apex Pin Self-drilling 3*120mm	4
Forearm External Fixator 7100-04		/	7100-0401	Rod-Rod Coupling (5mm Rod)	2
		/	7100-0402	4-Hole Pin Clamp (3/4mm Pin)	2
		/	7100-0403	Post 0 Degree 5mm	2
		/	7100-0404	Aluminum Connecting Rod 5*200mm	1
		/	7100-0405	Apex Pin Self-drilling 3*120mm	4
Ulna & Radius External Fixator 7100-05		/	7100-0501	Rod-Rod Coupling (5mm Rod)	1
		/	7100-0502	Rod-Pin Coupling (3/4mm Pin, 5mm Rod)	2
		/	7100-0503	Pin Clamp (3/4mm Pin)	1
		/	7100-0504	Aluminum Connecting Rod 5*200mm	1
		/	7100-0505	Apex Pin Self-drilling 3*120mm	4
Wrist Joint External Fixator 7100-06		/	7100-0601	Wrist Joint 5mm	1
		/	7100-0602	Rod-Pin Coupling (3/4mm Pin, 5mm Rod)	4
		/	7100-0603	Apex Pin Self-drilling 3*120mm	4
Elbow Joint External Fixator 7100-07		/	7100-0701	Elbow Joint 5mm	1
		/	7100-0702	Rod-Pin Coupling (3/4mm Pin, 5mm Rod)	4
		/	7100-0703	Apex Pin Self-drilling 4*120mm	4
Finger & Palm External Fixator 7100-08		/	7100-0801	Rod-Pin Coupling (3/2mm)	2
		/	7100-0802	Aluminum Connecting Rod 3*80mm	1
		/	7100-0803	Apex Pin Self-drilling 2*70mm	4
Humeral Shaft External Fixator 7100-09		/	7100-0901	Rod-Rod Coupling (5mm Rod)	2
		/	7100-0902	4-Hole Pin Clamp (3/4mm Pin)	2
		/	7100-0903	Post 5mm/30 Degree	2
		/	7100-0904	Aluminum Connecting Rod 5*200mm	1
		/	7100-0905	Apex Pin Self-drilling 4*120mm	4

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

Humeral Shaft External Fixator 7100-10		/	7100-1001	Rod-Rod Coupling (5mm Rod)	4
		/	7100-1002	Rod-Pin Coupling (3/4mm Pin, 5mm Rod)	2
		/	7100-1003	4-Hole Pin Clamp (3/4mm Pin)	1
		/	7100-1004	Aluminum Connecting Rod 5*250mm	2
		/	7100-1005	Aluminum Connecting Rod 5*200mm	1
		/	7100-1006	Post 5mm/30 Degree	2
		/	7100-1007	Apex Pin Self-drilling 4*120mm	4
Humeral Shaft External Fixator 7100-11		/	7100-1101	Compression Distraction Tube 15mm	1
		/	7100-1102	Tube-Rod Coupling 15/5mm	2
		/	7100-1103	4-Hole Pin Clamp (3/4mm Pin)	2
		/	7100-1104	Post 5mm/30 Degree	2
		/	7100-1105	Apex Pin Self-drilling 4*120mm	4
7100-02 Hoffmann External Fixation Large Fragment Set		1	7100-0201	Rod-Pin Coupling (5/6mm Pin, 8mm Rod)	12
		2	7100-0202	Rod-Rod Coupling (8mm Rod)	3
		3	7100-0203	5-Hole Pin Clamp (5/6mm Pin)	4
		4	7100-0204	30 Degree Angled Post 8mm	2
		5	7100-0205	Straight Post 8mm	12
		6	7100-0206	Apex Pin Self-drilling 5*120mm	4
		7	7100-0207	Apex Pin Self-drilling 5*150mm	2
		8	7100-0208	Apex Pin Self-drilling 5*180mm	1
		9	7100-0209	Drill Sleeve	8
		10	7100-0210	Carbon Fiber Connecting Rod 8*200mm	4
		11	7100-0211	Carbon Fiber Connecting Rod 8*250mm	2
		12	7100-0212	Carbon Fiber Connecting Rod 8*300mm	2
		13	7100-0213	Semi Circular Rod 8*180mm	2
		14	7100-0214	Semi Circular Rod 8*200mm	4
		15	7100-0215	Reduction Wrench	4
		16	7100-0216	Wrench with Extension Bar	4
		17	7100-0217	Adjust Wheel	2
		18	7100-0218	Apex Pin Chuck	1
		19	7100-0219	Arc Drilling Stand	1
		20	7100-0220	T-shaped Clamp Wrench 7#	1
		21	7100-0221	T-shaped Pin Wrench 5mm	1
		22	7100-0222	Aluminum Box	1
Hoffmann External Fixation Large Fragment Set Individuals		/	7100-0201	Rod-Pin Coupling (5/6mm Pin, 8mm Rod)	1
		/	7100-0202	Rod-Rod Coupling (8mm Rod)	1

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

Hoffmann External Fixation Large Fragment Set Individuals		/	7100-0203	5-Hole Pin Clamp (5/6mm Pin)	1
		/	7100-2001	Tube-Rod Coupling 20/8mm	2
		/	7100-0205	Straight Post 8mm	1
		/	7100-0204	30 Degree Angled Post 8mm	1
		/	7100-0234	Aluminum Connecting Rod 8*200mm	1
			7100-0235	Aluminum Connecting Rod 8*250mm	1
			7100-0236	Aluminum Connecting Rod 8*300mm	1
		/	7100-0210	Carbon Fiber Connecting Rod 8*200mm	1
			7100-0211	Carbon Fiber Connecting Rod 8*250mm	1
			7100-0212	Carbon Fiber Connecting Rod 8*300mm	1
		/	7100-0233	Semi Circular Rod 8*160mm	1
			7100-0213	Semi Circular Rod 8*180mm	1
			7100-0214	Semi Circular Rod 8*200mm	1
		/	7100-2002	Compression Distraction Tube 20mm	1
		/	7100-2028	Compression Distraction Single Rod 8mm	1
		/	7100-2029	Compression Distraction Double Rod 8mm	1
		/	7100-1903	Ankle Joint 8mm	1
		/	7100-1803	Knee Joint 8mm L/R	1
		/	7100-0216	Wrench with Extension Bar	1
		/	7100-0215	Reduction Wrench	1
		/	7100-0230	Wrench	1

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

Hoffmann External Fixation Large Fragment Set Individuals		/	7100-0219	Arc Drilling Stand	1		
				/	7100-0209	Drill Sleeve 5mm	1
	/	7100-0231	Drill Sleeve 6mm	1			
		/	7100-0218	Apex Pin Chuck 5mm	1		
		/	7100-0232	Apex Pin Chuck 6mm	1		
		/	7100-0217	Adjust Wheel	1		

Proximal Tibia External Fixator 7100-12		/	7100-1201	Rod-Rod Coupling (8mm Rod)	4
		/	7100-1202	Rod-Pin Coupling (5/6mm Pin, 8mm Rod)	3
		/	7100-1203	5-Hole Pin Clamp (5/6mm Pin)	1
		/	7100-1204	30 Degree Angled Post 8mm	2
		/	7100-1205	Semi Circular Rod 8*180mm	1
		/	7100-1206	Aluminum Connecting Rod 8*250mm	2
		/	7100-1207	Apex Pin Self-drilling 5*150mm	5

Femoral Shaft External Fixator 7100-13		/	7100-1301	Rod-Rod Coupling (8mm Rod)	2
		/	7100-1302	Rod-Pin Coupling (5/6mm Pin, 8mm Rod)	4
		/	7100-1303	Carbon Fiber Connecting Rod 8*200mm	1
		/	7100-1304	Carbon Fiber Connecting Rod 8*250mm	2
		/	7100-1305	Apex Pin Self-drilling 5*150mm	5

Tibia & Percone External Fixator 7100-14		/	7100-1401	Rod-Rod Coupling (8mm Rod)	6
		/	7100-1402	5-Hole Pin Clamp (5/6mm Pin)	2
		/	7100-1403	30 Degree Angled Post 8mm	4
		/	7100-1404	Aluminum Connecting Rod 8*250mm	2
		/	7100-1405	Aluminum Connecting Rod 8*300mm	1
		/	7100-1406	Apex Pin Self-drilling 5*150mm	4

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

Proximal Femur External Fixator 7100-15		/	7100-1501	Rod-Rod Coupling (8mm Rod)	4
		/	7100-1502	5-Hole Pin Clamp (5/6mm Pin)	2
		/	7100-1503	30 Degree Angled Post 8mm	4
		/	7100-1504	Aluminum Connecting Rod 8*300mm	2
		/	7100-1505	Apex Pin Self-drilling 5*150mm	2
		/	7100-1506	Apex Pin Self-drilling 5*180mm	2
Pelvic External Fixator 7100-16		/	7100-1601	Rod-Rod Coupling (8mm Rod)	6
		/	7100-1602	5-Hole Pin Clamp (5/6mm Pin)	2
		/	7100-1603	30 Degree Angled Post 8mm	2
		/	7100-1604	Aluminum Connecting Rod 8*300mm	2
		/	7100-1605	Aluminum Connecting Rod 8*200mm	2
		/	7100-1606	Apex Pin Self-drilling 5*150mm	4
Ankle Joint External Fixator 7100-17		/	7100-1701	Rod-Rod Coupling (8mm Rod)	3
		/	7100-1702	Rod-Pin Coupling (5/6mm Pin, 8mm Rod)	4
		/	7100-1703	Aluminum Connecting Rod 8*250mm	2
		/	7100-1704	Aluminum Connecting Rod 8*300mm	1
		/	7100-1705	Apex Pin Self-drilling 5*150mm	4
Knee Joint External Fixator 7100-18		/	7100-1801	Rod-Rod Coupling (8mm Rod)	2
		/	7100-1802	Rod-Pin Coupling (5/6mm Pin, 8mm Rod)	4
		/	7100-1803	Knee Joint 8mm L/R	1
		/	7100-1804	5-Hole Pin Clamp (5/6mm Pin)	1
		/	7100-1805	30 Degree Angled Post 8mm	1
		/	7100-1806	Aluminum Connecting Rod 8*150mm	1
		/	7100-1807	Apex Pin Self-drilling 5*150mm	6
Ankle Joint External Fixator 7100-19		/	7100-1901	Rod-Rod Coupling (8mm Rod)	1
		/	7100-1902	Rod-Pin Coupling (5/6mm Pin, 8mm Rod)	4
		/	7100-1903	Ankle Joint 8mm	1
		/	7100-1906	Apex Pin Self-drilling 5*150mm	4
		/	7100-1907	Apex Pin Self-drilling 3*120mm	1
Femoral External Fixator 7100-20		/	7100-2001	Tube-Rod Coupling 20/8mm	2
		/	7100-2002	Compression Distraction Tube 20mm	1
		/	7100-2003	5-Hole Pin Clamp (5/6mm Pin)	2
		/	7100-2004	30 Degree Angled Post 8mm	2
		/	7100-2005	Apex Pin Self-drilling 6*150mm	4

PRIME Surgical International (PSI)

STYKER HOFFMANN TYPE EXTERNAL FIXTOR

Femoral External Fixator 7100-21		/	7100-2101	Tube-Rod Coupling 20/8mm	2
		/	7100-2102	Compression Distraction Tube 20mm	1
		/	7100-2103	5-Hole Pin Clamp (5/6mm Pin)	2
		/	7100-2104	30 Degree Angled Post 8mm	2
		/	7100-2105	Apex Pin Self-drilling 6*150mm	4

PRIME SURGICAL