

## PRIME GYNECOLOGY & OBSTETRICS CATALOG - A



PRIME SURGICAL  
GYNECOLOGY  
INTERNATIONAL



PRIME SURGICAL  
INTERNATIONAL



PRIME SURGICAL  
INTERNATIONAL

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

---

PRIME SURGICAL INTERNATIONAL



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

The advertisement features a collection of various surgical instruments, including scalpels, forceps, and scissors, arranged on a dark background. The text is overlaid on the right side of the image.

**Zia ur Rahman**  
0092-300-2621673  
0092-313-5710887

**PRIME**  
**Surgical International**  
*Deals in All Kinds of Quality Instruments*

Ph: +92 213 5643935  
E.mail id: [prime.instruments@hotmail.com](mailto:prime.instruments@hotmail.com)  
[zia@primeinstrumentintl.com](mailto:zia@primeinstrumentintl.com) [www.primeinstrumentintl.com](http://www.primeinstrumentintl.com)  
Office # 25, Ground Floor, Queen Center M.T. Khan Road, Karachi.

Office address: QUEEN CENTER, OFFICE No.206, 2ND Floor, MT. KHAN Road , Karachi. Page 3 of 469

[Prime.instruments@hotmail.com](mailto:Prime.instruments@hotmail.com): [zia@primeinstrumentintl.com](mailto:zia@primeinstrumentintl.com) 0300-2621673

[www.primeinstrumentintl.com](http://www.primeinstrumentintl.com)

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

---

PRIME SURGICAL INTERNATIONAL

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

---

PRIME SURGICAL INTERNATIONAL

## OB/GYN

### Laser Instruments

#### PILLING® LASER SPECULA



Inner blades are matte finished to help reduce reflection of the laser beam. External blades are polished stainless steel for ease of insertion. One smoke evacuation tube is on the upper blade for rapid removal of steam and smoke from the ablated tissue. Available in three sizes.

| CATALOG # | DESCRIPTION                                                                                          |
|-----------|------------------------------------------------------------------------------------------------------|
| 508610    | Small, blades 3-1/4" (8.5 cm) long x 7/8" (2.2 cm) wide, suction tube 4-5/8" (12 cm) long (unbent)   |
| 508612    | Medium, blades 4-3/8" (11 cm) long x 1-3/8" (3.5 cm) wide, suction tube 4-5/8" (12 cm) long (unbent) |
| 508614    | Large, blades 5" (13 cm) long x 1-1/2" (3.7 cm) wide, suction tube 5-1/2" (14 cm) long (unbent)      |

#### SCHROEDER TENACULUM



Ebonized for use with laser surgery.

| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 508694    | 10" (25.5 cm) |

#### THOMS-ALLIS VULSELLUM FORCEPS




Ebonized for use with laser surgery.

| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 508696    | Straight, 6 x 7 Teeth, 8" (20.5 cm) |

## OB/GYN

### Laser Instruments

#### LASER MIRROR



Highly polished stainless steel front surface mirror.  
Suction to draw smoke and steam.

| CATALOG # | DESCRIPTION                                  |
|-----------|----------------------------------------------|
| 509209    | 6 mm diameter, working length 8-1/4" (21 cm) |


#### MILAN-LONDON VACUUM TUBE



Ebonized for use with laser surgery. Clips onto a standard metal Graves speculum for evacuation of smoke and plume.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 508602    | 4-7/8" (12.5 cm) |

#### OSSOFF MICRO SUCTION TUBE EBONIZED



Ebonized for use with laser surgery.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 509204    | 7 Fr., working length 7-7/8" (20 cm) |

## OB/GYN

### Retractors

#### SIMS VAGINAL RETRACTORS




212241

Single-ended.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 212241    | 1-1/8" x 3-5/8" (3 x 9 cm) |
| 212242    | 1-1/2" x 3-1/2" (4 x 9 cm) |

### Laser Instruments

#### BELLINA LASER RULERS




508600

508601

For determining depth and length of tissue for ablation purposes.

| CATALOG # | DESCRIPTION                                                                    |
|-----------|--------------------------------------------------------------------------------|
| 508600    | Lateral ruler, graduated in millimeters up to 15 mm                            |
| 508601    | Straight ruler, graduated in millimeters with a ridge indicating the 5 mm mark |

#### CERVICAL LASER HOOKS



508606

Ebonized for use with laser surgery.

| CATALOG # | DESCRIPTION                                                               |
|-----------|---------------------------------------------------------------------------|
| 508605    | Small, 2 mm hook, 6-1/2" (16.5 cm) working length, overall 9-1/2" (24 cm) |
| 508606    | Large, 4 mm hook, 7" (18 cm) working length, overall 9-7/8" (25 cm)       |



## OB/GYN

### Retractors

#### JACKSON VAGINAL RETRACTORS



212264

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 212263    | 1-9/16" x 3-1/4" (4 x 8.5 cm) |
| 212264    | 1-9/16" x 3-19/32" (4 x 9 cm) |
| 212265    | 1-19/32" x 4" (4 x 10 cm)     |

#### PRYOR-PEAN RETRACTOR



212295

| CATALOG # | DESCRIPTION                                                 |
|-----------|-------------------------------------------------------------|
| 212295    | 1" x 4-1/2" (2.5 x 11.5 cm), blade length 10-3/8" (26.5 cm) |

#### SIMS VAGINAL RETRACTORS



212228

Double-ended.

| CATALOG # | DESCRIPTION                                                  |
|-----------|--------------------------------------------------------------|
| 212225    | 7/8" x 2-5/8" (2 x 6.5 cm) and 1" x 3-1/8" (2.5 x 8 cm)      |
| 212226    | 1" x 2-1/2" (2.5 x 6.5 cm) and 1-1/4" x 3" (3 x 7.5 cm)      |
| 212227    | 1" x 3-1/4" (2.5 x 8.5 cm) and 1-3/16" x 3-1/2" (3 x 9 cm)   |
| 212228    | 1-1/4" x 3" (3 x 7.5 cm) and 1-3/8" x 3-1/4" (3.5 x 8.5 cm)  |
| 212229    | 1-3/16" x 3-1/4" (3 x 8.5 cm) and 1-7/16" x 4" (3.5 x 10 cm) |



## OB/GYN

### Retractors

#### DE LEE UNIVERSAL RETRACTOR



| CATALOG # | DESCRIPTION                                      |
|-----------|--------------------------------------------------|
| 212350    | 3" x 2-1/8" (7.5 x 5.5 cm) blade, 9-1/2" (24 cm) |

#### EASTMAN VAGINAL RETRACTORS



| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 212271    | 1-3/8" x 3-1/4" (3.5 x 8.5 cm) |
| 212275    | 1-3/8" x 3-7/8" (3.5 x 10 cm)  |

#### HEANEY-SIMON VAGINAL RETRACTOR



| CATALOG # | DESCRIPTION                                             |
|-----------|---------------------------------------------------------|
| 212305    | 1-1/8" x 4-3/4" (2.5 x 12 cm), blade length 11" (28 cm) |

## OB/GYN

### Specula

#### AUVARD WEIGHTED VAGINAL SPECULA



| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 212195    | 3 lbs., 3" (7.5 cm)         |
| 212196    | 2-1/2 lbs., 3-3/8" (8.5 cm) |
| 212197    | 3 lbs., 4" (10 cm)          |

#### GARRIGUE WEIGHTED VAGINAL SPECULUM



| CATALOG # | DESCRIPTION                                                                                    |
|-----------|------------------------------------------------------------------------------------------------|
| 212185    | Cut-out handle allows drainage. Blade 1-5/8" x 3-1/8" (4 x 8 cm), weight 2 lbs., chrome plated |

#### HENROTIN WEIGHTED VAGINAL SPECULUM



| CATALOG # | DESCRIPTION                                                     |
|-----------|-----------------------------------------------------------------|
| 212187    | Blade 1-1/2" x 3-1/2" (3.8 x 9 cm), weight 3 lbs, chrome plated |

## OB/GYN

### Specula

#### PEDERSON VAGINAL SPECULA



212011

Narrow, straight side blades.

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 212010    | Small, 3-1/4" x 1/2" (8.5 x 1.5 cm)     |
| 212011    | Medium, 3-1/4" x 1" (8.5 x 2.5 cm)      |
| 212012    | Large, 4-3/4" x 1" (12 x 2.5 cm) blades |

#### SLENDER PEDERSON SPECULUM



KM64288

| CATALOG # | DESCRIPTION                                                                                         |
|-----------|-----------------------------------------------------------------------------------------------------|
| KM64288   | Narrow blades, 4-11/32" (110 mm) long x 11/16" (17.5 mm) wide with a 2-7/8" (73 mm) maximum opening |

#### KOGAN ENDOCERVICAL SPECULUM



KM64830

| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| KM64830   | With gauge and set screw, 9-1/2" (24 cm) |

## OB/GYN

### Specula

#### GRAVES VAGINAL SPECULA



212001

| CATALOG # | DESCRIPTION                                                               |
|-----------|---------------------------------------------------------------------------|
| 212000    | Small blade 3-1/4" (8.5 cm) long x 7/8" (2.2 cm) wide                     |
| 212001    | Medium blade 4-3/8" (11 cm) long x 1-3/8" (3.5 cm) wide                   |
| 212002    | Large blade 5" (13 cm) long x 1-1/2" (3.7 cm) wide                        |
| 212005    | Medium, right open-sided blade 4-3/8" (11 cm) long x 1-3/8" (3.5 cm) wide |



212005

## OB/GYN

### Tenaculums and Vulsellums

#### M.G.H. VULSELLUM



| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 213175    | 3 x 3 teeth, 6-1/8" (15.5 cm) |

#### SCHROEDER VULSELLUM



| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 213180    | Straight, 9-7/8" (25 cm) |

#### TEALE UTERINE VULSELLUM



| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 754370    | Curved sideways, 3 x 4 teeth, 9-1/4" (23.5 cm) |

## OB/GYN

### Tenaculums and Vulsellums

#### HENROTIN VULSELLUM



| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 213150    | 8-1/4" (21 cm) |

#### JACOB'S VULSELLUM FORCEPS



| CATALOG # | DESCRIPTION                                            |
|-----------|--------------------------------------------------------|
| 213155    | Straight, 2 x 2 teeth, 8-1/4" (21 cm)                  |
| 213157    | Curved sideways, serrated, 2 x 2 teeth, 8-1/4" (21 cm) |



#### KELLY VULSELLUM



| CATALOG # | DESCRIPTION                     |
|-----------|---------------------------------|
| 213160    | Curved sideways, 9-1/2" (24 cm) |



## OB/GYN

### Tenaculums and Vulsellums

#### STAUDE TENACULUM




| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 213090    | 5-3/8" (13.5 cm) |

#### STAUDE-MOORE TENACULUM



| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 213095    | 6-3/8" (16 cm) |

#### WYLIE UTERINE TENACULUM



| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 754160    | 1 x 2 teeth, 9-1/2" (24 cm) |

## OB/GYN

### Tenaculums and Vulsellums

#### DUPLAY TENACULUM



| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 213085    | 9-7/8" (25 cm) |

#### HULKA TENACULUM AND PROBE



| CATALOG # | DESCRIPTION                                                         |
|-----------|---------------------------------------------------------------------|
| 754125    | Single tooth, probe 4-3/4" (12 cm) w/ 3.6 mm tip, 11-1/8" (28.5 cm) |

#### JARCHO TENACULUM



| CATALOG # | DESCRIPTION  |
|-----------|--------------|
| 213065    | 8" (20.5 cm) |

#### SCHROEDER TENACULUM



| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 213080    | 9-7/8" (25 cm) |

## OB/GYN

### Tenaculums and Vulsellums

#### ADAIR TENACULUM



| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 213120    | 6-7/8" (17.5 cm) |

#### BARRETT TENACULUM



| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 213110    | 7-1/4" (18.5 cm) |

#### BRAUN TENACULUM




| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 213100    | 9-5/8" (24.5 cm) |

## OB/GYN

### Curettes

#### RANDALL SUCTION BIOPSY CURETTE



| CATALOG # | DESCRIPTION                                  |
|-----------|----------------------------------------------|
| KM64808   | Sharp, curved, 4 mm diameter, 9-3/4" (25 cm) |

#### SIMS UTERINE CURETTES




Sharp blades. Malleable shafts.

| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 212749    | Size 0, 10-7/8" (27.5 cm) |
| 212750    | Size 1, 11-1/4" (28.5 cm) |
| 212751    | Size 2, 11-1/4" (28.5 cm) |
| 212752    | Size 3, 11-3/8" (29 cm)   |
| 212753    | Size 4, 11-1/2" (29 cm)   |
| 212754    | Size 5, 11-1/2" (29 cm)   |
| 212755    | Size 6, 11-5/8" (29.5 cm) |

## OB/GYN

### Curettes


#### KEVORKIAN-YOUNGE BIOPSY CURETTE



Pilling® Pattern.

| CATALOG # | DESCRIPTION     |
|-----------|-----------------|
| 212910    | 12-1/8" (31 cm) |

#### KEVORKIAN-YOUNGE CURETTES



Sharp tip.

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| KM64796   | Without basket, 12" (30.5 cm) |
| KM64802   | With basket, 12" (30.5 cm)    |

## OB/GYN

### Probes, Sounds and Dilators

#### PRATT UTERINE DILATORS



212647

Pilling® Pattern, satin finish. Double-ended, 11-5/8" (29.5 cm). Stainless steel.

| CATALOG # | DESCRIPTION                                           |
|-----------|-------------------------------------------------------|
| 212647    | Set of 8 in canvas case, sizes 13-15 Fr. to 41-43 Fr. |

#### PRATT UTERINE DILATORS




752394

Weck® Pattern, bright finish. Double-ended, 11-5/8" (29.5 cm).

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 752390    | Size 13-15 Fr. |
| 752391    | Size 17-19 Fr. |
| 752392    | Size 21-23 Fr. |
| 752393    | Size 25-27 Fr. |
| 752394    | Size 29-31 Fr. |
| 752395    | Size 33-35 Fr. |
| 752396    | Size 37-39 Fr. |
| 752397    | Size 41-43 Fr. |

## Curettes

#### HEANEY BIOPSY CURETTE



212905

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 212905    | 9-3/4" (25 cm) |

## OB/GYN

### Probes, Sounds and Dilators

#### HEGAR UTERINE DILATORS



212673

Pilling® Pattern, satin finish. Double-ended, 7-7/8" (20 cm).

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 212673    | Set of 8, sizes 3 mm to 18 mm |
| 212675    | Size 3-4 mm                   |
| 212676    | Size 5-6 mm                   |
| 212677    | Size 7-8 mm                   |
| 212678    | Size 9-10 mm                  |
| 212679    | Size 11-12 mm                 |
| 212680    | Size 13-14 mm                 |
| 212681    | Size 15-16 mm                 |
| 212682    | Size 17-18 mm                 |



212678

#### HEGAR UTERINE DILATORS



752354

Weck® Pattern, bright finish. Double-ended, hollow, 7-7/8" (20 cm).

| CATALOG # | DESCRIPTION          |
|-----------|----------------------|
| 752338    | Size 1/2 mm          |
| 752347    | Size 19/20 mm        |
| 752348    | Size 21/22 mm        |
| 752349    | Size 23/24 mm        |
| 752351    | Size 25/26 mm        |
| 752354    | Complete set in case |



## OB/GYN

### Probes, Sounds and Dilators

#### HANK UTERINE DILATORS



Pilling® Pattern, satin finish. Double-ended. Length 10-3/4" (27.5 cm). All sizes available individually.

| CATALOG # | DESCRIPTION                                                                     |
|-----------|---------------------------------------------------------------------------------|
| 212659    | Hank Uterine Dilators, set of 6, sizes 11-12 Fr. to 21-22 Fr. (3.66 to 7.33 mm) |
| 212661    | 9-10 Fr. (3 x 3.3 mm)                                                           |
| 212662    | 11-12 Fr. (3.6 x 4 mm)                                                          |
| 212663    | 13-14 Fr. (4.3 x 4.7 mm)                                                        |
| 212664    | 15-16 Fr. (5 x 5.3 mm)                                                          |
| 212665    | 17-18 Fr. (5.7 x 6 mm)                                                          |
| 212666    | 19-20 Fr. (6.3 x 6.7 mm)                                                        |
| 212667    | 21-22 Fr. (7 x 7.3 mm)                                                          |
| 212668    | 23-24 Fr. (7.7 x 8 mm)                                                          |
| 212669    | 25-26 Fr. (8.3 x 8.7 mm)                                                        |
| 212670    | 27-28 Fr. (9 x 9.3 mm)                                                          |
| 212671    | 29-30 Fr. (9.3 x 10 mm)                                                         |

#### HANK UTERINE DILATORS



Weck® Pattern, bright finish. Double-ended. End lengths 2-1/2" (6.5 cm). Length 10-3/8" (26.5 cm).

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 752380    | 9/10 American (4.5 x 5 mm)   |
| 752381    | 11/12 American (5.5 x 6 mm)  |
| 752382    | 13/14 American (6.5 x 7 mm)  |
| 752383    | 15/16 American (7.5 x 8 mm)  |
| 752384    | 17/18 American (8.5 x 9 mm)  |
| 752385    | 19/20 American (9.5 x 10 mm) |

## OB/GYN

### Probes, Sounds and Dilators

#### EMMETT UTERINE PROBE



212600

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 212600    | Sterling Silver, 11" (28 cm) |

#### SIMPSON UTERINE SOUND



212605

| CATALOG # | DESCRIPTION                                                            |
|-----------|------------------------------------------------------------------------|
| 212605    | 4 mm, graduated in inches, silver plated, malleable, 12-3/4" (32.5 cm) |

#### SIMS UTERINE SOUNDS



212612

| CATALOG # | DESCRIPTION                                                      |
|-----------|------------------------------------------------------------------|
| 212610    | 4 mm, non-graduated, silver plated, malleable, 12-3/4" (32.5 cm) |
| 212612    | 4 mm, graduated in cm, malleable, 12-3/4" (32.5 cm)              |

## OB/GYN

### Scissors

#### BUSCH UMBILICAL SCISSORS



| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| KM67294   | Curved sideways, 6-1/2" (16.5 cm) |

#### UMBILICAL SCISSORS



| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 214090    | 4-1/8" (10.5 cm)         |
| KM67292   | Straight, 4-3/8" (11 cm) |

#### BRAUN EPISIOTOMY SCISSORS



| CATALOG # | DESCRIPTION                                      |
|-----------|--------------------------------------------------|
| 214081    | Wide lower blade, angled to side, 8-3/4" (22 cm) |

## OB/GYN

### Delivery Forceps

#### DE LEE OBSTETRICAL FORCEPS



| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 214595    | 12" (30.5 cm) |

#### MCLEAN-TUCKER-LUIKART OBSTETRICAL FORCEPS



| CATALOG # | DESCRIPTION       |
|-----------|-------------------|
| 214605    | 15-1/2" (39.5 cm) |

#### SIMPSON OBSTETRICAL FORCEPS



| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 214550    | 14" (35.5 cm) |
| 214555    | 12" (30.5 cm) |

#### SIMPSON-LUIKART OBSTETRICAL FORCEPS



| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| 214570    | Fenestrated, smooth, 14" (35.5 cm) |

## OB/GYN

### Biopsy Punches

#### WITTNER BIOPSY FORCEPS



752526

8 x 6 mm toothed basket jaw, 8-7/8" (22.5 cm) stem.

| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 752526    | Straight, 11-1/2" (29 cm) |

### Circumcision Clamps

#### GOMCO® CIRCUMCISION CLAMPS




742110

| CATALOG # | DESCRIPTION        |
|-----------|--------------------|
| 742100    | Extra small, 11 mm |
| 742102    | Newborn, 13 mm     |
| 742104    | Infant, 15 mm      |
| 742106    | Child, 16 mm       |
| 742108    | Youth, 21 mm       |
| 742110    | Adult, 26 mm       |

## OB/GYN

### Biopsy Punches

#### THOMAS-GAYLOR CERVICAL BIOPSY FORCEPS



| CATALOG # | DESCRIPTION                     |
|-----------|---------------------------------|
| 213011    | Curved shanks, 9-1/4" (23.5 cm) |

213011

#### THOMAS-GAYLOR BIOPSY PUNCH



| CATALOG # | DESCRIPTION                     |
|-----------|---------------------------------|
| KM65084   | Interfitted jaw, 9-1/2" (24 cm) |

KM65084

#### TISCHLER BIOPSY PUNCHES



| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| KM65188   | Straight, 3 x 5 mm, 8" (20.5 cm)           |
| KM65186   | Straight, 3 x 7 mm, 9" (23 cm)             |
| KM65204   | Shaft only, 3 x 7 mm bite                  |
| KM65205   | Pistol grip handle only for Tischler punch |

KM65186

## OB/GYN

### Biopsy Punches

#### KEVORKIAN-YOUNGE CERVICAL BIOPSY FORCEPS



| CATALOG # | DESCRIPTION       |
|-----------|-------------------|
| 213036    | 10-1/8" (25.5 cm) |

#### SCHUBERT BIOPSY PUNCHES



| CATALOG # | DESCRIPTION                                      |
|-----------|--------------------------------------------------|
| 766520    | Angled to side, 5 x 8 mm punch, 9-5/8" (24.5 cm) |
| 766522    | Straight, 5 x 8 mm punch, 11" (28 cm)            |

#### SCHUBERT BIOPSY PUNCHES




| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| KM65115   | Straight, 8-1/2" (21.5 cm) |
| KM65118   | Angled, 10-1/2" (26.5 cm)  |



## OB/GYN

### Forceps and Clamps

#### HEANEY-BALLANTINE HYSTERECTOMY CLAMPS



Longitudinally serrated jaws. Single groove.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 214180    | Straight, 8-1/4" (21 cm) |
| 214182    | Curved, 8-1/8" (20.5 cm) |

#### WERTHEIM-CULLEN PEDICLE CLAMP



| CATALOG # | DESCRIPTION                                          |
|-----------|------------------------------------------------------|
| 214230    | Angled, longitudinally serrated jaws, 8-5/8" (22 cm) |

### Biopsy Punches

#### EPPENDORFER BIOPSY PUNCH




| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| KM65162   | Straight, 4 x 8 mm bite, 9" (23 cm) |

## OB/GYN

### Forceps and Clamps

#### DEBAKEY® HYSTERECTOMY CLAMP



| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 214215    | Curved, 9-3/8" (24 cm) |

#### HEANEY HYSTERECTOMY FORCEPS



Double serrated jaws with cross-groove.

| CATALOG # | DESCRIPTION                                       |
|-----------|---------------------------------------------------|
| 214170    | Curved, light jaws, single groove, 8-1/4" (21 cm) |
| 214172    | Curved, heavy jaws, single groove, 8-1/4" (21 cm) |
| 214175    | Curved, heavy jaws, double groove, 8-1/4" (21 cm) |



## OB/GYN

### Forceps and Clamps


#### PHANEUF UTERINE ARTERY FORCEPS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                                   |
|-----------|---------------------------------------------------------------|
| 214161    | Short, heavy jaws, 1 x 2 distal teeth, angled, 8-1/4" (21 cm) |

#### PHANEUF UTERINE ARTERY FORCEPS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                               |
|-----------|-----------------------------------------------------------|
| 762500    | Straight, serrated, 1 x 2 teeth, 4 mm tip, 8-1/4" (21 cm) |

#### BLAND CERVICAL TRACTION FORCEPS



| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 213140    | Straight, 10 mm tip, 9-1/4" (23.5 cm) |

## OB/GYN

### Forceps and Clamps


#### POLK PLACENTA FORCEPS



Fenestrated jaw with cross serrations.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 770730    | Straight, 9-5/8" (24.5 cm) |

#### SOMER UTERINE ELEVATING FORCEPS



| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 213230    | 9" (23 cm)  |

#### WALLACE CAESAREAN FORCEPS




| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 770405    | Serrated, 25 mm jaw, 8" (20.5 cm) |

## OB/GYN

### Forceps and Clamps

#### BOZEMAN UTERINE DRESSING FORCEPS




| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 212500    | Double-curved, one large ring, 10-1/4" (26 cm) |
| 212505    | Standard model, curved, 10-1/4" (26 cm)        |

#### DE LEE CERVICAL GRASPING FORCEPS



| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 213050    | 10" (25.5 cm) |

#### KELLY PLACENTA CLAMP



| CATALOG # | DESCRIPTION     |
|-----------|-----------------|
| 213055    | 12-5/8" (32 cm) |

## RETRACTORS, SELF-RETAINING

Abdominal/Balfour-Type Retractors

### WEXLER DEEP SPREADER



485262

For Wexler Retractors.

| CATALOG # | DESCRIPTION                                                 |
|-----------|-------------------------------------------------------------|
| 485262    | 4" x 1-7/8" blades (10 x 5 cm), 5-1/4" (13.5 cm), pair of 2 |

### WEXLER EXPANDABLE BLADE




485225

| CATALOG # | DESCRIPTION                                       |
|-----------|---------------------------------------------------|
| 485225    | 2-3/4" (7 cm) wide blade, expands to 5" (12.5 cm) |

## RETRACTORS, SELF-RETAINING

Abdominal/Balfour-Type Retractors

### WEXLER STANDARD OCTAGON FRAME



485205

Lateral catches for two blades.

| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 485205    | 11-5/16" x 9-3/16" (28.5 x 23.5 cm) |

### UNIVERSAL JOINT FOR WEXLER RETRACTORS



485210

| CATALOG # | DESCRIPTION                                                    |
|-----------|----------------------------------------------------------------|
| 485210    | Allows rotation on ball with locking arm for Wexler Retractors |

### WEXLER DEAVER BLADE



485220

For Wexler Retractors.

| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 485220    | 2" wide x 4-5/16" deep (5 x 11 cm), 9" (23 cm) |



## RETRACTORS, SELF-RETAINING

### Abdominal/Balfour-Type Retractors

#### BALFOUR FOURTH BLADE ATTACHMENTS



165170

Locks securely to arm of Balfour retractors, holding fourth blade. (Blade not included.)

| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 165170    | For right arm |
| 165171    | For left arm  |

#### O'SULLIVAN-O'CONNOR ABDOMINAL RETRACTORS



485450

| CATALOG # | DESCRIPTION                                                                                                                                                                                                                                                                                                    |
|-----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 165300    | Adjustable shape ring<br>Spread 7-1/2" (20 cm)<br>Two lateral blades, 1" wide x 3" deep (2.5 x 7.5 cm)<br>Set complete with three detachable blades:<br>2 small, 1-3/4" wide x 3" deep (4.5 x 7.5 cm), and<br>1 large, 2-3/4" wide x 4-1/8" deep (7.0 x 10.5 cm)<br>Positive action ratchet regulates diameter |
| 485450    | Blades on frame are 73 x 27 mm<br>Three detachable blades:<br>2 small, 1-3/4" wide x 3-1/8" deep (4.5 x 8 cm), and<br>1 large, 2-7/8" wide x 4" deep (7 x 10 cm)<br>Spread 12-1/2" (32 cm)                                                                                                                     |

#### WEXLER BANTAM RETRACTOR



486100

For pediatric and minor surgery use. Set comes complete with frame, 1 center blade and 2 side blades.

| CATALOG # | DESCRIPTION                                                                                                                   |
|-----------|-------------------------------------------------------------------------------------------------------------------------------|
| 486100    | Center blade is 1-1/4" deep x 3/4" wide (3 x 2 cm), side blade is 1-1/8" deep x 7/8" wide (3 x 2.5 cm), 3-7/8" (10 cm) spread |

## RETRACTORS, SELF-RETAINING

### Abdominal/Balfour-Type Retractors

#### BALFOUR CENTER BLADES



These blades are made for use with the Balfour Retractors 165100 through 165133.

| CATALOG # | DESCRIPTION                                              |
|-----------|----------------------------------------------------------|
| 165141    | Standard, 3-1/2" wide x 1-1/4" deep (8.8 cm x 3.2 cm)    |
| 165142    | Standard, 2-3/4" wide x 1-3/4" deep (7 cm x 4.4 cm)      |
| 165143    | Supra-Pubic, 3-1/2" wide x 2-5/8" deep (8.8 cm x 6.8 cm) |
| 165145    | Narrow, 1-5/8" wide x 3-1/2" deep (4.2 cm x 8.8 cm)      |
| 165146    | Wide, 2-1/2" wide x 2-5/8" deep (6.3 cm x 6.8 cm)        |
| 165147    | Wide, 2-1/2" wide x 3-1/2" deep (6.3 cm x 8.8 cm)        |
| 165149    | Horgan malleable 2" wide x 7-3/4" long (5 x 19.5 cm)     |
| 165150    | Horgan malleable 2-1/2" wide x 8" long (6.5 x 20.5 cm)   |

#### YU HOLTGREVE MALLEABLE BLADE




| CATALOG # | DESCRIPTION                                                |
|-----------|------------------------------------------------------------|
| 165160    | 2-1/2" wide x 5-7/8" deep x 9-3/8" long (6.5 x 15 x 24 cm) |

## RETRACTORS, SELF-RETAINING

Abdominal/Balfour-Type Retractors

### TRANSPLANT RECOVERY BALFOUR RETRACTORS



Transplant retractors for the recovery of the abdominal organs. Spiked lateral blades for maximum hold.

| CATALOG # | DESCRIPTION                                                                                                |
|-----------|------------------------------------------------------------------------------------------------------------|
| 165107    | Adult, 2 x 2 sharp prongs, blades 2-3/4" (7 cm), overall length of arms 7-7/8" (20 cm), 17" (43 cm) spread |

## RETRACTORS, SELF-RETAINING

### Abdominal/Balfour-Type Retractors

#### BALFOUR RETRACTORS PERMANENT SIDE BLADE MODEL



165105

| CATALOG # | DESCRIPTION                                                                                                                  |
|-----------|------------------------------------------------------------------------------------------------------------------------------|
| 165100    | Side blades 2-1/2" (6.5 cm) deep. Center blade 2-3/4" (7 cm) wide x 1-3/4" (4.5 cm) deep, standard size 8" (20.5 cm) opening |
| 165105    | Side blades 3-1/2" (9 cm) deep. Center blade 3-1/2" (9 cm) wide x 1-3/4" (4.5 cm) deep, standard size 8" (20.5 cm) opening   |

#### IMPROVED BALFOUR RETRACTORS INTERCHANGEABLE SIDE BLADE MODEL



165130

Lateral blades may be easily snapped on and off. Set includes 2 pairs of detachable wire lateral blades 2-1/2" (6.5 cm) and 3-1/2" (9 cm) deep.

| CATALOG # | DESCRIPTION                                                             |
|-----------|-------------------------------------------------------------------------|
| 165130    | Center blade 2-1/2" (6.5 cm) wide, standard size, 8" (20.5 cm) opening  |
| 165133    | Center blade 3-1/2" (8.8 cm) wide, extra wide 10-1/2" (26.5 cm) opening |
| 165135    | Extra solid lateral blades 2-1/2" (6.5 cm) deep, per pair               |
| 165137    | Extra wire lateral blades, 2-1/2" (6.5 cm) deep, per pair               |
| 165138    | Extra wire lateral blades, 3-1/2" (9 cm) deep, per pair                 |



165133

## RETRACTORS, SELF-RETAINING

### Abdominal/Balfour-Type Retractors


#### BALFOUR RETRACTORS



165200

Trigger controlled, positive automatic ratchet retaining design. Free, smooth easy movement of all parts. Large selection of blade designs and sizes. Stainless steel construction.

| CATALOG # | DESCRIPTION                                                                                                         |
|-----------|---------------------------------------------------------------------------------------------------------------------|
| 165200    | Wire lateral blades 2" deep (5 cm), center blade 3-5/8" wide x 1-3/4" deep (9 x 4.5 cm), 7" (18 cm) spread          |
| 165203    | Wire lateral blades 4-3/8" deep (11 cm), center blade 4-3/16" deep x 2-1/2" wide (10.5 x 6.5 cm), 7" (18 cm) spread |
| 165240    | Pediatric. Miniature model, arms 3-1/2" long (9 cm), lateral blades 3/4" (2 cm) deep, 4" (10 cm) spread             |



165240



## RETRACTORS, HAND-HELD

### Specialty Retractors

#### ALLISON LUNG RETRACTORS



Wire, spatula-shaped blade.

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 341678    | Child size, 10-1/4" (26 cm)   |
| 341679    | Adult Size, 12-3/8" (31.5 cm) |



#### DAVIDSON SCAPULA RETRACTORS



| CATALOG # | DESCRIPTION                                            |
|-----------|--------------------------------------------------------|
| 341650    | Small size, blade 2" (5 cm) wide, 7-3/4" (19.5 cm)     |
| 341651    | Large size, blade 3-1/2" (9 cm) wide, 7-3/4" (19.5 cm) |

#### DEVINE-MILLARD-AUFRECHT NASAL RETRACTOR FIBEROPTIC LIGHTING



| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 525150    | 7-3/4" (19.5 cm) |

#### FACE LIFT RETRACTOR FIBEROPTIC LIGHTING ON #2 DEAYER RETRACTOR



| CATALOG # | DESCRIPTION     |
|-----------|-----------------|
| 525215    | 10-5/8" (27 cm) |

## RETRACTORS, HAND-HELD

### Ribbon-Type Malleable Retractors

#### DAVIS BRAIN SPATULAS



380206

| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 380200    | 1/4" (6.5 mm) wide blade  |
| 380202    | 3/8" (9.5 mm) wide blade  |
| 380204    | 1/2" (12.5 mm) wide blade |
| 380206    | 5/8" (16 mm) wide blade   |
| 380208    | 3/4" (19 mm) wide blade   |
| 380210    | 1" (25.5 mm) wide blade   |

#### HAIGHT RETRACTORS



341694

Malleable "T"-shaped blades. For abdominal and intra-thoracic surgery.

| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 341690    | 3" wide x 11-1/2" long (7.5 x 29 cm)      |
| 341692    | 3-1/2" wide x 11-3/4" long (9 x 30 cm)    |
| 341694    | 4" wide x 12" long (10 x 30.5 cm)         |
| 341696    | 4-3/4" wide x 12-3/4" long (12 x 32.5 cm) |
| 341698    | 5-1/2" wide x 13-1/2" long (14 x 34.5 cm) |

#### OCHSNER RIBBON RETRACTORS



164170

Malleable.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 164164    | 1" x 13" (2.5 x 33 cm)     |
| 164166    | 1-1/4" x 13" (3 x 33 cm)   |
| 164168    | 1-1/2" x 13" (4 x 33 cm)   |
| 164170    | 2" x 13" (5 x 33 cm)       |
| 164172    | 2-1/2" x 13" (6.5 x 33 cm) |
| 164174    | 3" x 13" (7.5 x 33 cm)     |
| 442288    | 7/8" x 7" (2 x 18 cm)      |

#### RIBBON RETRACTORS



KM51766

Malleable.

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| KM51752   | 1/2" (12.5 mm) wide blade, 6" (15 cm)   |
| KM51755   | 1" (25.5 mm) wide blade, 7-1/2" (19 cm) |
| KM51766   | 3/4" (19 mm) wide blade, 13" (33 cm)    |



## RETRACTORS, HAND-HELD

### Abdominal-Type Retractors

#### RICHARDSON-EASTMAN RETRACTORS



| CATALOG # | DESCRIPTION                                                                                                                                                                  |
|-----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 480195    | Set of 2:<br>Small 9-7/8" (25 cm) blade A 30 x 29 mm,<br>blade B 24 x 22 mm<br>Large 10-5/8" (27 cm) blade A 63 x 42 mm,<br>blade B 39 x 37 mm, Weck® Pattern, bright finish |
| 480196    | Small only:<br>9-7/8" (25 cm) blade A 30 x 29 mm, blade B<br>24 x 22 mm, Weck® Pattern, bright finish                                                                        |
| KM51214   | Large only:<br>10-5/8" (27 cm) blade A 65 x 43 mm, blade B<br>39 x 38 mm, Satin finish                                                                                       |

#### WYLIE RENAL VEIN RETRACTORS



Designed for the retraction of the left renal vein during procedures on the abdominal aorta or renal arteries. Gently curved distal tip for atraumatic control. Stainless Steel.

| CATALOG # | DESCRIPTION                                                          |
|-----------|----------------------------------------------------------------------|
| 341795    | Short, 1" wide x 3-3/4" deep (2.5 x 9.5 cm) blade, 11-3/4" (30 cm)   |
| 341797    | Long, 1" wide x 6-7/8" deep (2.5 x 17.5 cm) blade, 12-3/4" (32.5 cm) |

#### WYLIE SPLANCHNIC RETRACTOR



| CATALOG # | DESCRIPTION                                                                                                                                                                                                                              |
|-----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 341780    | Rigid, acutely curved blade is 6-3/4" deep (17 cm)<br>Distal 1/2" is angled inward and deeply cross<br>serrated, which provides positive, non-slip traction<br>Blade is 2-1/2" wide (6.3 cm)<br>Large, fluted handle, overall length 13" |

## RETRACTORS, HAND-HELD

### Abdominal-Type Retractors

#### RICHARDSON RETRACTORS



Blade face is relatively flat with lateral edges curved upwards. Distal edge is straight across and turned downward.

| CATALOG # | DESCRIPTION                                               |
|-----------|-----------------------------------------------------------|
| 164542    | 3/4" x 1-1/4" blade (20 x 32 mm), length 9-3/4" (25 cm)   |
| 164543    | 1" x 1-1/4" blade (25 x 32 mm), length 9-3/4" (25 cm)     |
| 164544    | 1-1/2" x 1-1/2" blade (38 x 38 mm), length 9-3/4" (25 cm) |

#### RICHARDSON RETRACTORS



Stainless steel, hollow handle.

| CATALOG # | DESCRIPTION                                           |
|-----------|-------------------------------------------------------|
| 481616    | 1" x 3/4" blade (25 x 20 mm), 9-1/2" (24 .5 cm)       |
| 481617    | 1-1/4" x 1-1/4" blade (32 x 32 mm), 9-1/2" (24 .5 cm) |
| 481618    | 1-1/2" x 1-1/2" blade (38 x 38 mm), 9-1/2" (24 .5 cm) |
| 481629    | 2" x 3/4" blade (48 x 20 mm), 9-1/2" (24 .5 cm)       |

#### RICHARDSON-EASTMAN RETRACTORS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                                                                                                                       |
|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| 164550    | Set of 2:<br>Small 10" (25.5 cm), blade A 33 x 30 mm,<br>blade B 26 x 23 mm<br>Large 10-3/4" (27.5 cm), blade A 65 x 43 mm,<br>blade B 42 x 38 mm |

## RETRACTORS, HAND-HELD

### Abdominal-Type Retractors

#### HARRINGTON SPLANCHNIC RETRACTORS



Cross serrated retracting surface.

| CATALOG # | DESCRIPTION                                                  |
|-----------|--------------------------------------------------------------|
| 341770    | 1-1/2" wide x 4-3/4" deep blade (4 x 12 cm), 12-1/2" (32 cm) |
| 341773    | 2-1/2" wide x 5" deep blade (6.5 x 12.5 cm), 12" (30.5 cm)   |

#### KELLY RETRACTORS



Flat blade with slight upward turn on sides and downward turn on slightly reduced width distal edge.

| CATALOG # | DESCRIPTION                                       |
|-----------|---------------------------------------------------|
| KM51222   | 3" x 2-1/2" blade (76 x 64 mm), 9-3/4" (25 cm)    |
| 164525    | 1-1/2" x 2" blade (38 x 50 mm), 9" (23 cm)        |
| 164526    | 2" x 2-3/8" blade (50 x 61 mm), 9-1/2" (24 cm)    |
| 164527    | 2-1/2" x 3" blade (64 x 76 mm), 10" (25.5 cm)     |
| 164528    | 3" x 3-1/2" blade (76 x 90 mm), 10-1/2" (26.5 cm) |

#### MAYO ABDOMINAL RETRACTOR



| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 482350    | 40 x 70 mm blade, 11" (28 cm) |

#### RICHARDSON APPENDECTOMY RETRACTOR



| CATALOG # | DESCRIPTION                                          |
|-----------|------------------------------------------------------|
| 164540    | 3/4" wide x 2" deep blade (2 x 5 cm), 9-1/2" (24 cm) |

## RETRACTORS, HAND-HELD

### Abdominal-Type Retractors

#### DOYEN RETRACTORS



Concave blade.

| CATALOG # | DESCRIPTION                                        |
|-----------|----------------------------------------------------|
| 164520    | 1-3/4" x 2-3/4" blade (4.5 x 7 cm), 8-3/4" (22 cm) |
| 164521    | 2-3/8" x 3-3/8" blade (6 x 8.5 cm), 9" (23 cm)     |
| 164522    | 2-1/4" x 3-7/8" blade (5.5 x 10 cm), 9" (23 cm)    |

#### DEBAKEY®-COOLEY RETRACTOR



Extra-deep retractor of greater than usual width. Blade is carefully proportioned to provide maximum working depth and to prevent tendency toward slipping upward and out. Satin finish to prevent glare from overhead lights.

| CATALOG # | DESCRIPTION                                                                      |
|-----------|----------------------------------------------------------------------------------|
| 341740    | 7-1/4" (18.5 cm) deep blade with 2-1/4" (5.5 cm) Harrington end, 12-1/4" (31 cm) |

#### HARRINGTON EXTRA DEEP RETRACTORS



Extra-deep with fluted handle, satin finish.

| CATALOG # | DESCRIPTION                                                   |
|-----------|---------------------------------------------------------------|
| 341790    | 2-1/4" wide x 8" deep blade (5.5 x 20.5 cm), 11-3/4" (30 cm)  |
| 341791    | 4-1/2" wide x 7" deep blade (11.5 x 18 cm), 13-1/8" (33.5 cm) |

#### HARRINGTON RETRACTORS



Stainless steel, hollow handle.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 482370    | 135 x 44 mm blade, 12-3/4" (32.5 cm) |
| 482371    | 135 x 65 mm blade, 13" (33 cm)       |

## RETRACTORS, HAND-HELD

### Abdominal-Type Retractors

#### W.M. LEMMON EXTRA-DEEP DEAVER RETRACTORS



Distal blade section extended to provide working depth of 9-1/2" (24 cm). Blade length 15" (38 cm). Large fluted handles.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 341730    | 1" (2.5 cm) wide blade |
| 341732    | 2" (5 cm) wide blade   |

#### WILLAUER-DEAVER RETRACTORS



Deaver blades mounted on comfortable handle.

| CATALOG # | DESCRIPTION                                                        |
|-----------|--------------------------------------------------------------------|
| 341720    | Narrow, 1" wide x 6-1/4" deep (2.5 x 16 cm) blade, 14-1/2" (37 cm) |
| 341722    | Wide, 2" wide x 6-1/4" deep (5 x 16 cm) blade, 14-1/2" (37 cm)     |

#### DEEP DEAVER RETRACTORS



Stainless steel, hollow handle.

| CATALOG # | DESCRIPTION                                           |
|-----------|-------------------------------------------------------|
| 481825    | 6" deep x 1-1/2" (15 x 4 cm) blade, 13-3/4" (35 cm)   |
| 481826    | 6" deep x 2" (15 x 5 cm) blade, 13-3/4" (35 cm)       |
| 481827    | 6" deep x 2-1/2" (15 x 6.5 cm) blade, 13-3/4" (35 cm) |



## RETRACTORS, HAND-HELD

### Abdominal-Type Retractors

#### DEAVER RETRACTORS

Standard design. Rigid stainless steel.

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 164110    | 3/4" (2 cm) wide, 6-7/8" (17.5 cm)      |
| 164112    | 7/8" (2 cm) wide, 7-7/8" (20 cm)        |
| 164115    | 1" (2.5 cm) wide, 8-7/8" (22.5 cm)      |
| 164118    | #1, 1" (2.5 cm) wide, 8-1/8" (20.5 cm)  |
| 164120    | #2, 1" (2.5 cm) wide, 9-7/8" (25 cm)    |
| 164122    | #3, 1" (2.5 cm) wide, 12" (30.5 cm)     |
| 164125    | #4, 1" (2.5 cm) wide, 14-3/8" (36.5 cm) |
| 164130    | #5, 1-1/2" (4 cm) wide, 11-7/8" (30 cm) |
| 164132    | #6, 2" (5 cm) wide, 12" (30.5 cm)       |
| 164134    | #7, 3" (7.5 cm) wide, 12" (30.5 cm)     |



164110



164112



164115




164122



164118



164125



164120



164130



164134



## RETRACTORS, HAND-HELD

### Double-Ended Retractors

#### MILLER-SENN RETRACTORS



9.8 mm prongs, 21 x 6.5 mm blade, 6-1/2" (16.5 cm).

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 480125    | 3 sharp prongs |
| 480126    | 3 blunt prongs |

#### U.S. ARMY (ARMY-NAVY) RETRACTORS



Nested set of 2.

| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 164715    | 5/8" (1.5 cm) wide blade, 8-1/2" (21.5 cm) |

### Abdominal-Type Retractors

#### DEAVER RETRACTORS



Hollow grip handle.

| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| KM51230   | 1" wide x 10" long blade (2.5 x 25.5 cm)       |
| KM51232   | 1" wide x 12" long blade (2.5 x 30.5 cm)       |
| KM51234   | 1" wide x 13" long blade (2.5 x 33 cm)         |
| KM51240   | 2" wide x 12" long blade (5 x 30.5 mm)         |
| KM51236   | 2-1/2" wide x 14-1/2" long blade (6.5 x 37 cm) |
| KM51242   | 3" wide x 12" long blade (7.5 x 30.5 cm)       |
| KM51238   | 3-1/2" wide x 12" long blade (9 x 30.5 cm)     |

## RETRACTORS, HAND-HELD

### Double-Ended Retractors

#### RAGNELL RETRACTOR



KM52094

KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                                                                    |
|-----------|------------------------------------------------------------------------------------------------|
| KM52094   | Blade A 4 wide x 8 mm deep,<br>Blade B 6 wide x 15 mm deep, Overall length<br>5-3/4" (14.5 cm) |

#### ROUX RETRACTORS



480140



480145

7-1/4" (18.5 cm).

| CATALOG # | DESCRIPTION                                                                                                 |
|-----------|-------------------------------------------------------------------------------------------------------------|
| 480140    | Small:<br>Blade A 1-1/8" wide x 1" deep (27 x 25 mm),<br>Blade B 1-3/8" wide x 1-1/2" deep (34 x 38 mm)     |
| 480145    | Large:<br>Blade A 1-1/4" wide x 1-3/8" deep (31 x 36 mm),<br>Blade B 1-1/4" wide x 1-7/8" deep (31 x 47 mm) |

#### SENN RETRACTORS



164750

Pilling® Pattern, satin finish. 1/4" x 3/4" (6 x 19 mm) flat blade, 6-1/4" (16 cm).

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 164750    | 3 sharp prongs |
| 164752    | 3 blunt prongs |

#### SENN-KANAVAL DELICATE RETRACTOR



054590

| CATALOG # | DESCRIPTION                                                   |
|-----------|---------------------------------------------------------------|
| 054590    | 3 sharp 7 x 7 mm prongs, 21 x 6 mm blade,<br>6-1/2" (16.5 cm) |

## RETRACTORS, HAND-HELD

### Double-Ended Retractors

#### PARKER RETRACTORS



480170

| CATALOG # | DESCRIPTION                                                                                                                                                        |
|-----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 480170    | Small set of 2 retractors.<br>Short retractor: Blade A 7/8" wide x 5/8" deep (21 x 16 mm), Blade B 1" wide x 5/8" deep (25 x 16 mm), Overall length 5-1/2" (14 cm) |
| 480172    | Large set of 2 retractors.<br>Short retractor: Blade A 1" wide x 1" deep (25 x 24 mm), Blade B 1" wide x 1" deep (26 x 24 mm), Overall length 6-3/4" (17 cm)       |

#### PARKER-MOTT RETRACTORS



164698

Set of 2.

| CATALOG # | DESCRIPTION                                                         |
|-----------|---------------------------------------------------------------------|
| 164698    | 1" x 2-1/2" (25 x 65 mm) blades,<br>Overall length 5-3/4" (14.5 cm) |

#### RAGNELL RETRACTOR



054600

| CATALOG # | DESCRIPTION                                                                                          |
|-----------|------------------------------------------------------------------------------------------------------|
| 054600    | Blade A 3 mm wide x 8 mm deep,<br>Blade B 5 mm wide x 17 mm deep,<br>Overall length 5-3/4" (14.5 cm) |

## RETRACTORS, HAND-HELD

### Double-Ended Retractors

#### GOELET RETRACTOR



480150

| CATALOG # | DESCRIPTION                                                                                                                     |
|-----------|---------------------------------------------------------------------------------------------------------------------------------|
| 480150    | Blade A 1" wide x 1-1/8" deep (26 x 30 mm),<br>Blade B 1-1/4" wide x 1-7/8" deep (31 x 47 mm),<br>Overall length 7-1/2" (19 cm) |

#### INFANT RETRACTOR



480105

| CATALOG # | DESCRIPTION                                       |
|-----------|---------------------------------------------------|
| 480105    | Blades 3/8" x 1/2" (10 x 13 mm), 5-5/8" (14.5 cm) |

#### LUKENS THYMUS RETRACTOR



514412

| CATALOG # | DESCRIPTION                                                                |
|-----------|----------------------------------------------------------------------------|
| 514412    | Blade 3/8" wide x 1/4" deep (9 x 6 mm),<br>Overall length 6-7/8" (17.5 cm) |

#### MAYO-COLLINS RETRACTORS



164710


Nested set of 2.

| CATALOG # | DESCRIPTION                                                                                                                 |
|-----------|-----------------------------------------------------------------------------------------------------------------------------|
| 164710    | Blade A 3/4" wide x 1-1/8" deep (20 x 30 mm),<br>Blade B 3/4" wide x 1-1/4" deep (20 x 35 mm),<br>Overall length 6" (15 cm) |

## RETRACTORS, HAND-HELD

### Double-Ended Retractors

#### CRILE RETRACTOR



| CATALOG # | DESCRIPTION                                                                                                                           |
|-----------|---------------------------------------------------------------------------------------------------------------------------------------|
| 480100    | Large blade 1/4" wide x 3/4" deep (7.5 x 20 mm),<br>Small blade 1/4" wide x 5/8" deep (5 x 17 mm),<br>Overall length 4-1/8" (10.5 cm) |

#### DAVIS RETRACTORS




| CATALOG # | DESCRIPTION                                                                          |
|-----------|--------------------------------------------------------------------------------------|
| 442245    | Blade A 3/8" wide (10 mm), Blade B 1/4" wide (6 mm), Overall length 6-1/8" (15.5 cm) |

#### DAVIS RETRACTORS



| CATALOG # | DESCRIPTION                                                                                                                       |
|-----------|-----------------------------------------------------------------------------------------------------------------------------------|
| KM52076   | Blade A 1/4" wide x 5/8" deep blade (6 x 17 mm),<br>Blade B 3/8" wide x 3/4" deep blade (9 x 20 mm),<br>Overall length 6" (15 cm) |

#### FARABEUF RETRACTORS



Nested set of 2.

| CATALOG # | DESCRIPTION                                                                                                                   |
|-----------|-------------------------------------------------------------------------------------------------------------------------------|
| 164700    | Blade A 5/8" wide x 1-1/8" deep (15 x 28 mm),<br>Blade B 5/8" wide x 1-1/4" deep (15 x 30 mm),<br>Overall length 6" (15 cm)   |
| 164702    | Blade A 3/8" wide x 1" deep (10 x 25 mm),<br>Blade B 1/2" wide x 1-1/4" deep (12 x 30 mm),<br>Overall length 4-1/2" (11.5 cm) |

## RETRACTORS, HAND-HELD


### Double-Ended Retractors

#### KLEINERT-RAGNELL RETRACTOR



| CATALOG # | DESCRIPTION                                                                                                                     |
|-----------|---------------------------------------------------------------------------------------------------------------------------------|
| 442260    | Blade A 1/4" wide x 3/4" deep (5.6 x 20 mm),<br>Blade B 1/8" wide x 3/8" deep (3.7 x 10 mm),<br>Overall length 6-1/8" (15.5 cm) |

#### BREWSTER PHRENIC RETRACTORS



Nested set of 2.

| CATALOG # | DESCRIPTION                                                                                                                                               |
|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| 341635    | Large blade 1-3/8" x 1-1/2" wide (20 x 40 mm) x<br>1-3/8" deep (35 mm),<br>Small blade 1" wide x 3/8" deep (25 x 10 mm),<br>Overall length 5-7/8" (15 cm) |

#### CHILDREN'S HOSPITAL RETRACTOR



| CATALOG # | DESCRIPTION                                                                                                                  |
|-----------|------------------------------------------------------------------------------------------------------------------------------|
| 164754    | Blade A 3/8" wide x 3/4" deep (10 x 19 mm),<br>Blade B 1/4" wide x 3/4" deep (6 x 19 mm),<br>Overall length 6-1/8" (15.5 cm) |



## RETRACTORS, HAND-HELD

### Single-Ended Retractors

#### SAUERBRUCH RETRACTORS



Long, slender blades with slightly curved tips.  
Hollow handle.

| CATALOG # | DESCRIPTION                                        |
|-----------|----------------------------------------------------|
| 164635    | 5/8" x 1-3/4" (1.5 x 4.5 cm) blade, 9" (23 cm)     |
| 164636    | 7/8" x 2-1/2" (2 x 6.5 cm) blade, 9" (23 cm)       |
| 164637    | 3/4" x 2-7/8" (2 x 7.5 cm) blade, 8-7/8" (22.5 cm) |

#### SEMB RETRACTOR



For retracting the muscles during resection of the posterior medial ends of the ribs and during freeing of the transverse processes and separation of the vertebrocostal joints. Stainless steel blade.

| CATALOG # | DESCRIPTION                                       |
|-----------|---------------------------------------------------|
| 341640    | Small, blade 3/4" wide x 1-3/4" deep (2 x 4.5 cm) |

### Double-Ended Retractors

#### "S" RETRACTOR



| CATALOG # | DESCRIPTION                                                                                                                   |
|-----------|-------------------------------------------------------------------------------------------------------------------------------|
| 054620    | Blade A 1/4" wide x 3/8" deep (6 x 10 mm),<br>Blade B 1/2" wide x 5/8" deep (13.5 x 15 mm)<br>Overall length 5-1/4" (13.5 cm) |

## RETRACTORS, HAND-HELD

### Single-Ended Retractors

#### MEYERDING FINGER RETRACTORS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 164790    | 7.5 x 17 mm blade, 6-7/8" (17.5 cm) |

#### MEYERDING FINGER RETRACTORS



KMedic® Pattern, satin finish, 7" (18 cm).

| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| KM52116   | #1, 7 mm wide x 18 mm deep blade                |
| KM52118   | #2, 6 mm wide x 17 mm deep blade, with teeth    |
| KM52122   | #3, 4 mm wide x 8 mm deep blade, with teeth     |
| KM52124   | #4, 4 mm wide x 10 mm deep blade                |
| KM52128   | #5, 4 mm wide x 7 mm deep blade                 |
| KM52130   | #6, 8 mm wide x 7 mm deep blade, 3 sharp prongs |



## RETRACTORS, HAND-HELD


### Single-Ended Retractors

#### LOVE NERVE RETRACTORS



| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| KM51530   | Straight, 7 mm wide blade, 8-1/2" (21.5 cm)    |
| KM51532   | Angled 45 degrees, 7 mm wide blade, 7" (18 cm) |
| KM51534   | Angled 90 degrees                              |

#### LOVE NERVE RETRACTORS



| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 481250    | 4 x 6 mm blade, straight, 8" (20.5 cm) |

#### MCBURNEY RETRACTOR



| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| 481530    | Fenestrated, 20 x 16 mm blade, 7-5/8" (19.5 cm) |

## RETRACTORS, HAND-HELD

### Single-Ended Retractors

#### LAHEY RETRACTOR



| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 481610    | 30 x 5 mm blade, 7-3/4" (19.5 cm) |

#### LANGENBECK RETRACTORS



Distal blade tip, slightly curved. Flat handle.  
Length 8-1/2" (21.5 cm).

| CATALOG # | DESCRIPTION                      |
|-----------|----------------------------------|
| 164630    | 3/8" x 1-5/8" (1 x 4 cm) blade   |
| 164631    | 1/2" x 1-5/8" (1.3 x 4 cm) blade |
| 164632    | 5/8" x 1-5/8" (1.5 x 4 cm) blade |

#### LITTLER RETRACTOR




| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 481500    | Curved, 20 x 12 mm blade, 7-1/2" (19 cm) |

## RETRACTORS, HAND-HELD

### Single-Ended Retractors

#### CUSHING NERVE AND VEIN RETRACTOR



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| 164865    | 12 mm wide blade, 8-1/2" (21.5 cm) |

#### CUSHING NERVE AND VEIN RETRACTORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 481230    | 16 x 15 mm blade, 8-5/8" (22 cm)      |
| 481231    | 14.5 x 15 mm blade, 12-3/4" (32.5 cm) |

#### DESMARRES LID RETRACTORS



| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 421610    | Size 1, 12 mm wide blade, 5-3/8" (13.5 cm) |
| 421620    | Size 2, 13 mm wide blade, 5-3/8" (13.5 cm) |
| 421630    | Size 3, 15 mm wide blade, 5-3/8" (13.5 cm) |
| 421640    | Size 4, 17 mm wide blade, 5-3/8" (13.5 cm) |

#### GREEN RETRACTOR



| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| 481520    | Fenestrated, 21 x 24 mm blade, 8-7/8" (22.5 cm) |

## RETRACTORS, HAND-HELD

### Single-Ended Retractors

#### BREWSTER RETRACTORS



Curved, convex blades.

| CATALOG # | DESCRIPTION                                            |
|-----------|--------------------------------------------------------|
| 164650    | 3/4" x 1-3/4" (2 x 4.5 cm) blade, 8-1/2" (21.5 cm)     |
| 164652    | 1" x 2" (2.5 x 5 cm) blade, 9" (23 cm)                 |
| 164656    | 1-3/4" x 2-1/2" (4.5 x 6.5 cm) blade, 9-1/4" (23.5 cm) |
| 164658    | 1-7/8" x 2-7/8" (5 x 7.5 cm) blade, 9-1/2" (24 cm)     |

#### COTTLE ALAR RETRACTOR



| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 442210    | 1/2" (1.5 cm) wide blade, 5-1/2" (14 cm) |

#### CUSHING DECOMPRESSION RETRACTOR



| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| 481220    | 21 x 10 mm blade, 8-1/8" (20.5 cm) |



## RETRACTORS, HAND-HELD

Hook and Rake Retractors

### ISRAEL RETRACTOR



| CATALOG # | DESCRIPTION                                   |
|-----------|-----------------------------------------------|
| 481390    | Blunt, 4 prongs, 41 x 43 mm, 7-3/8" (18.5 cm) |

### OLLIER RAKE RETRACTOR



| CATALOG # | DESCRIPTION                                               |
|-----------|-----------------------------------------------------------|
| KM52282   | 4 blunt prongs, 5.5 cm wide x 4 cm deep blade, 9" (23 cm) |

### FREEMAN FACE LIFT RETRACTOR



| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 442295    | 1-5/8" (4 cm) wide blade, 7-1/8" (18 cm) |

## RETRACTORS, HAND-HELD

### Hook and Rake Retractors

#### VOLKMANN RAKE RETRACTORS



164403



164414

Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 164402    | Sharp, 2 prongs, 8-1/2" (21.5 cm) |
| 164403    | Sharp, 3 prongs, 8-1/2" (21.5 cm) |
| 164404    | Sharp, 4 prongs, 8-1/2" (21.5 cm) |
| 164406    | Sharp, 6 prongs, 8-1/2" (21.5 cm) |
| 164412    | Blunt, 2 prongs, 8-1/2" (21.5 cm) |
| 164413    | Blunt, 3 prongs, 8-1/2" (21.5 cm) |
| 164414    | Blunt, 4 prongs, 8-1/2" (21.5 cm) |
| 164416    | Blunt, 6 prongs, 8-1/2" (21.5 cm) |

#### VOLKMANN RETRACTORS



KM52264

KMedic® Pattern, satin finish. Sharp.

| CATALOG # | DESCRIPTION                                                            |
|-----------|------------------------------------------------------------------------|
| KM52266   | 6 prongs, 1-1/8" wide x 1/2" deep blade (3 x 1.5 cm), 8-1/2" (21.5 cm) |
| KM52254   | 1 prong, 1/2" deep blade (15 mm), 8-1/2" (21.5 cm)                     |
| KM52264   | 4 prongs, 3/4" wide x 1/2" deep blade (18 x 13 mm), 8-1/2" (21.5 cm)   |

#### MUELLER RAKE RETRACTOR



KM52338

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| KM52338   | 8 prongs, 1-1/2" (38 mm) wide |

#### DEEP RAKE RETRACTORS



KM52294

| CATALOG # | DESCRIPTION                                              |
|-----------|----------------------------------------------------------|
| KM52292   | 4 blunt prongs, 2-1/4" wide x 1-1/2" deep (5.5 x 4 cm)   |
| KM52294   | 4 blunt prongs, 2-1/4" wide x 2" deep (5.5 x 5 cm)       |
| KM52296   | 4 blunt prongs, 2-1/4" wide x 2-1/2" deep (5.5 x 6.5 cm) |
| KM52298   | 4 blunt prongs, 2-1/4" wide x 3" deep (5.5 x 7.5 cm)     |

## RETRACTORS, HAND-HELD

Hook and Rake Retractors

### BLAIR (ROLLET) RETRACTOR



KM54368

| CATALOG # | DESCRIPTION                                          |
|-----------|------------------------------------------------------|
| KM54366   | Rigid, 4 sharp prongs, 6 mm wide blade, 5" (12.5 cm) |
| KM54368   | Rigid, 4 blunt prongs, 5 mm wide blade, 5" (12.5 cm) |
| KM54372   | Flexible, 4 sharp prongs, 5-1/2" (14 cm)             |

### MURPHY RAKE RETRACTORS



164444



164454

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 164443    | Sharp, 3 prongs, 7-3/4" (19.5 cm) |
| 164444    | Sharp, 4 prongs, 8" (20.5 cm)     |
| 164446    | Sharp, 6 prongs, 7-7/8" (20 cm)   |
| 164453    | Blunt, 3 prongs, 8" (20.5 cm)     |
| 164454    | Blunt, 4 prongs, 8" (20.5 cm)     |
| 164456    | Blunt, 6 prongs, 7-3/4" (19.5 cm) |

### VOLKMANN POCKET RETRACTORS



164374

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 164374    | Sharp, 4 prongs, 4-1/2" (11.5 cm) |
| 164376    | Sharp, 6 prongs, 4-1/2" (11.5 cm) |
| 164384    | Blunt, 4 prongs, 4-1/2" (11.5 cm) |

## RETRACTORS, HAND-HELD

Hook and Rake Retractors

### JACKSON TRACHEA RETRACTOR



Weck® Pattern, bright finish.


| CATALOG # | DESCRIPTION          |
|-----------|----------------------|
| 070116    | 2 prongs, 7" (18 cm) |

### DINGMAN ZYGOMA HOOK



| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 026275    | Sharp, 13 mm hook, 7" (18 cm) |

### BERNAY FINGER RAKE TRACHEAL RETRACTORS



| CATALOG # | DESCRIPTION                     |
|-----------|---------------------------------|
| 514430    | 1 prong, sharp, 6-5/8" (17 cm)  |
| 514431    | 2 prongs, sharp, 6-5/8" (17 cm) |
| 514433    | 3 prongs, sharp, 6-5/8" (17 cm) |
| 514435    | 1 prong, blunt, 6-5/8" (17 cm)  |
| 514436    | 2 prongs, blunt, 6-5/8" (17 cm) |
| 514437    | 3 prongs, blunt, 6-5/8" (17 cm) |
| 514438    | 4 prongs, blunt, 6-5/8" (17 cm) |

### LAHEY FINGER RETRACTORS



| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| KM52082   | 1 prong, 4-3/4" (12 cm)      |
| KM52086   | 2 sharp prongs, 5" (12.5 cm) |
| KM52087   | 3 sharp prongs, 5" (12.5 cm) |

## RETRACTORS, HAND-HELD

### Hook and Rake Retractors

#### TRACHEA HOOK



| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 054165    | Sharp, 8.7 x 1.4 mm hook, 7-1/2" (19 cm)   |
| KM52098   | Hupp Trachea Hook, sharp, 6-1/2" (16.5 cm) |

#### SOLVAY SPAY HOOK



| CATALOG # | DESCRIPTION  |
|-----------|--------------|
| 801886    | 8" (20.5 cm) |

#### JACKSON TRACHEAL RETRACTOR, DOUBLE-ENDED




Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                                                   |
|-----------|-------------------------------------------------------------------------------|
| 514414    | Double-pronged, 1" wide x 1/4" deep (2.5 cm x 0.5 cm) blade, 6-7/8" (17.5 cm) |

## RETRACTORS, HAND-HELD

### Hook and Rake Retractors

#### DANDY NERVE HOOK



| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 185714    | Straight, 5.1 mm hook, 9" (23 cm) |

#### SISSON NERVE HOOK



| CATALOG # | DESCRIPTION                                          |
|-----------|------------------------------------------------------|
| 060120    | Blunt, 90 degree angled, 7 mm hook, 6-1/8" (15.5 cm) |

#### JACKSON TRACHEAL TENACULUM



| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 514413    | Hook depth 1/4" (0.5 cm), 5-1/4" (13.5 cm) |

#### NEWS TRACHEA RETRACTOR




| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 070114    | Sharp, 10 mm hook, 5-7/8" (15 cm) |



## RETRACTORS, HAND-HELD

### Skin Hooks

#### ROLLET RETRACTOR



481180

| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| 481180    | Sharp, 4 prongs, 2.8 x 5.5 mm, 5-1/4" (13.5 cm) |

#### SCHNITMAN SHARP SKIN HOOK WITH EYE



524230

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 524230    | 6-1/2" (16.5 cm) |

#### TYRELL HOOKS



KM54342

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| KM54338   | 1 sharp prong, 5" (12.5 cm) |
| KM54342   | 1 blunt prong, 5" (12.5 cm) |

#### TYRELL HOOKS



054015

Delicate, 1.3 mm hook, 4-3/4" (12 cm).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 054015    | Sharp       |
| 054016    | Blunt       |

#### WIENER SKIN HOOKS



KM54346

Single hook, 5" (12.5 cm).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| KM54346   | Sharp       |
| KM54348   | Blunt       |

## RETRACTORS, HAND-HELD

### Skin Hooks

#### KLEINERT-KUTZ HOOK



KM52060

KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| KM52060   | 3 mm hook, 6-1/4" (16 cm) |
| KM52064   | 5 mm hook, 6-1/4" (16 cm) |
| KM52068   | 6 mm hook, 6-1/4" (16 cm) |

#### KMEDIC® SKIN HOOK



KM52055

| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| KM52055   | 1 sharp prong, 2 mm hook, 4-3/4" (12 cm) |


#### O'CONNER SKIN HOOK



164786

| CATALOG # | DESCRIPTION             |
|-----------|-------------------------|
| 164786    | Sharp, 5-3/8" (13.5 cm) |

#### PECK-JOSEPH SKIN HOOK



442110

| CATALOG # | DESCRIPTION                                      |
|-----------|--------------------------------------------------|
| 442110    | 2 prongs, 10 mm between prongs, 5-5/8" (14.5 cm) |

#### RIGID RAKE RETRACTOR SKIN HOOKS



KM52172

6-1/2" (16.5 cm).

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| KM52156   | Sharp, 1 prong                    |
| KM52158   | Sharp, 2 prongs, 6-1/2" (16.5 cm) |
| KM52162   | Sharp, 3 prongs, 6-1/2" (16.5 cm) |
| KM52164   | Blunt, 1 prong, 6-1/2" (16.5 cm)  |
| KM52166   | Blunt, 2 prongs, 6-1/2" (16.5 cm) |
| KM52172   | Blunt, 3 prongs, 6-1/2" (16.5 cm) |

## RETRACTORS, HAND-HELD

### Skin Hooks

#### JOHNSON SKIN HOOK



| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| 054125    | Sharp, 5 mm hook, 6-1/2" (16.5 cm) |


#### JOSEPH SKIN HOOKS



Sharp, 4 mm hook. Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| 054100    | 1 prong, 6-1/8" (15.5 cm)                       |
| 054102    | 2 prongs, 2 mm between prongs, 6-1/4" (16 cm)   |
| 054105    | 2 prongs, 5 mm between prongs, 6-1/4" (16 cm)   |
| 054107    | 2 prongs, 6 mm between prongs, 6-1/4" (16 cm)   |
| 054110    | 2 prongs, 9 mm between prongs, 6-1/8" (15.5 cm) |

#### JOSEPH SKIN HOOKS



Pilling® Pattern, bright finish. 6-1/4" (16 cm).

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 442100    | 1 prong          |
| 442102    | 2 prongs, 3 mm   |
| 442104    | 2 prongs, 5.5 mm |
| 442106    | 2 prongs, 7 mm   |
| 442108    | 2 prongs, 10 mm  |

#### KLEINERT-KUTZ SKIN HOOKS




Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| 054140    | Small, 3 mm hook, 6-1/4" (16 cm)   |
| 054142    | Medium, 5 mm hook, 6-1/4" (16 cm)  |
| 054144    | Large, 7 mm hook, 6-1/8" (15.5 cm) |

## RETRACTORS, HAND-HELD

### Skin Hooks

#### FREER RETRACTOR



| CATALOG # | DESCRIPTION                     |
|-----------|---------------------------------|
| 481175    | Sharp, 2 prongs, 5-7/8" (15 cm) |

#### FREER SKIN HOOK




| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| KM52040   | Sharp, 2 prongs, 2.5 mm wide, 6" (15 cm) |

#### GILLIES CONVERSE SKIN HOOKS



| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 054128    | Small, sharp, 3 mm hook, 7-3/8" (18.5 cm) |
| 054130    | Large, sharp, 4 mm hook, 7-1/8" (18 cm)   |

#### GUTHRIE SKIN HOOK



| CATALOG # | DESCRIPTION                                          |
|-----------|------------------------------------------------------|
| 054126    | Sharp, 2 prongs, 2 mm between prongs, 5-1/8" (13 cm) |

## RETRACTORS, HAND-HELD

### Skin Hooks

#### CONVERSE SKIN HOOKS

| <br>442147 | CATALOG # | DESCRIPTION             |
|---------------------------------------------------------------------------------------------|-----------|-------------------------|
|                                                                                             | 442147    | Small, 7-1/4" (18.5 cm) |
|                                                                                             | 442145    | Large, 7" (18 cm)       |

#### COTTLE SKIN HOOKS

| <br>054135<br><br>054120 | CATALOG # | DESCRIPTION                                           |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|-------------------------------------------------------|
|                                                                                                                                                                                            | 054120    | Sharp, 2 prongs, 12 mm between prongs, 5-7/8" (15 cm) |
|                                                                                                                                                                                            | 054135    | Sharp, 3.5 mm hook, 4-7/8" (12.5 cm)                  |

#### FLEXIBLE RAKE RETRACTOR SKIN HOOKS

| <br>KM52182 | CATALOG # | DESCRIPTION                     |
|------------------------------------------------------------------------------------------------|-----------|---------------------------------|
|                                                                                                | KM52174   | Sharp, 1 prong, 6-1/4" (16 cm)  |
|                                                                                                | KM52176   | Sharp, 2 prongs, 6-1/4" (16 cm) |
|                                                                                                | KM52182   | Sharp, 3 prongs, 6-1/4" (16 cm) |
|                                                                                                | KM52184   | Blunt, 1 prong, 6-1/4" (16 cm)  |
|                                                                                                | KM52188   | Blunt, 2 prongs, 6-1/4" (16 cm) |
|                                                                                                | KM52192   | Blunt, 3 prongs, 6-1/4" (16 cm) |

#### FOMON RETRACTOR

| <br><br>054500 | CATALOG # | DESCRIPTION                                              |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|----------------------------------------------------------|
|                                                                                                                                                                                      | 054500    | Ball tip, 2 prongs, 10 mm between prongs, 6-3/8" (16 cm) |

## NEEDLE HOLDERS

Suture Tags

### GREGORY TECHNIQUE® SUTURE "STAY" CLAMPS



Smooth jaw hemostats for tagging sutures. May be used with or without rubber sheathing.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 351830    | Straight, 3-7/8" (10 cm), TC |
| 351831    | Curved, 3-7/8" (10 cm), TC   |



## NEEDLE HOLDERS

Wire Twisters

### BERRY STERNAL TECHNIQUE® NEEDLE HOLDER



A large, heavy-duty needle holder with 6 mm wide jaws designed to securely hold the large needles used to close sternotomies. LIFE-LOK® Box-Lock.

| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 343030    | 28 x 6 mm jaw, 7-3/4" (19.5 cm), TC |

### WIRE FORCEPS



| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 351348    | 5-3/8" (13.5 cm) |

## NEEDLE HOLDERS

Wire Twisters

### BAUMGARTEN WIRE TWISTER



| CATALOG # | DESCRIPTION                                         |
|-----------|-----------------------------------------------------|
| 142521    | Heavy square, 17.5 x 5 mm jaw, 8-1/8" (20.5 cm), TC |

### RUBIO TECHNIQUE® WIRE TWISTER



| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 351349    | 12 x 4.5 mm jaw, 5" (12.5 cm), TC |

## NEEDLE HOLDERS

Castro-Type Micro Needle Holders

### BARRAQUER DELICATE NEEDLE HOLDERS



With lock, round handle.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| KM44133   | Straight, smooth, 5-1/2" (14 cm)       |
| KM44129   | Curved, smooth, 5-1/2" (14 cm)         |
| KM44135   | Straight, smooth, 7" (18 cm)           |
| KM44131   | Curved, smooth, 7" (18 cm)             |
| KM44134   | Straight, smooth, 5-1/2" (14 cm), TC   |
| KM44130   | Curved, smooth, 5-1/2" (14 cm), TC     |
| KM44123   | Straight, serrated, 5-1/2" (14 cm), TC |
| KM44121   | Curved, serrated, 5-1/2" (14 cm), TC   |
| KM44136   | Straight, smooth, 7" (18 cm), TC       |
| KM44132   | Curved, smooth, 7" (18 cm), TC         |
| KM44124   | Straight, serrated, 7" (18 cm), TC     |
| KM44122   | Curved, serrated, 7" (18 cm), TC       |

### BARRAQUER DELICATE NEEDLE HOLDERS



Without lock, round handle.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| KM44151   | Straight, smooth, 5-1/2" (14 cm)       |
| KM44137   | Curved, smooth, 5-1/2" (14 cm)         |
| KM44153   | Straight, smooth, 7" (18 cm)           |
| KM44139   | Curved, smooth, 7" (18 cm)             |
| KM44152   | Straight, smooth, 5-1/2" (14 cm), TC   |
| KM44138   | Curved, smooth, 5-1/2" (14 cm), TC     |
| KM44127   | Straight, serrated, 5-1/2" (14 cm), TC |
| KM44125   | Curved, serrated, 5-1/2" (14 cm), TC   |
| KM44154   | Straight, smooth, 7" (18 cm), TC       |
| KM44150   | Curved, smooth, 7" (18 cm), TC         |
| KM44128   | Straight, serrated, 7" (18 cm), TC     |
| KM44126   | Curved, serrated, 7" (18 cm), TC       |

## NEEDLE HOLDERS

### Castro-Type Micro Needle Holders

#### JACOBSON MICRO NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                                                                   |
|-----------|-------------------------------------------------------------------------------|
| 358001    | Straight, round handle, with lock, 11 x 0.8 mm Diamo-grip jaw, 9" (23 cm)     |
| 358002    | Straight, round handle, with lock, 1.8 mm Diamo-grip jaw, 9" (23 cm)          |
| 358003    | Straight, round handle, with lock, 16 x 0.9 mm Diamo-grip jaw, 8-1/4" (21 cm) |

#### JACOBSON MICRO-SURGICAL NEEDLE HOLDERS



Castroviejo-type lock can be removed if not required.  
Stainless steel jaws.

| CATALOG # | DESCRIPTION                                                                  |
|-----------|------------------------------------------------------------------------------|
| 357601    | Straight, flat handle, with lock, 13 x 1.3 mm smooth jaw, 7" (18 cm)         |
| 357603    | Straight, flat handle, with lock, 13.5 x 1.2 mm smooth jaw, 8-1/2" (21.5 cm) |

#### BARRAQUER MICRO NEEDLE HOLDERS



Without lock, round handle.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| KM44175   | Straight, smooth, 5-1/2" (14 cm)       |
| KM44171   | Curved, smooth, 5-1/2" (14 cm)         |
| KM44177   | Straight, smooth, 7" (18 cm)           |
| KM44173   | Curved, smooth, 7" (18 cm)             |
| KM44176   | Straight, smooth, 5-1/2" (14 cm), TC   |
| KM44172   | Curved, smooth, 5-1/2" (14 cm), TC     |
| KM44161   | Straight, serrated, 5-1/2" (14 cm), TC |
| KM44159   | Curved, serrated, 5-1/2" (14 cm), TC   |
| KM44178   | Straight, smooth, 7" (18 cm), TC       |
| KM44174   | Curved, smooth, 7" (18 cm), TC         |
| KM44162   | Straight, serrated, 7" (18 cm), TC     |
| KM44160   | Curved, serrated, 7" (18 cm), TC       |

## NEEDLE HOLDERS

### Castro-Type Micro Needle Holders

#### MICRO SURGERY TECHNIQUE® NEEDLE HOLDERS



Flat handle, smooth jaws.

| CATALOG # | DESCRIPTION                                 |
|-----------|---------------------------------------------|
| 354966    | Straight, 13 x 1.4 mm jaw, 7" (18 cm), TC   |
| 354967    | Straight, 13 x 1.4 mm jaw, 8" (20.5 cm), TC |

#### SPENCER CORONARY TECHNIQUE® NEEDLE HOLDER



| CATALOG # | DESCRIPTION                                                   |
|-----------|---------------------------------------------------------------|
| 354969    | Straight, smooth, 12 x 1.2 mm jaw, 0.9 mm tip, 7" (18 cm), TC |

#### GIANNINI TECHNIQUE® NEEDLE HOLDER



| CATALOG # | DESCRIPTION                                      |
|-----------|--------------------------------------------------|
| 354971    | Straight, 17.5 x 1.8 mm jaw, 10-1/4" (26 cm), TC |

#### JACOBSON MICRO NEEDLE HOLDER



| CATALOG # | DESCRIPTION                                                                   |
|-----------|-------------------------------------------------------------------------------|
| 357603DG  | Straight with flat handle, lock, 18 x 1.3 mm Diamo-Grip jaw, 8-1/2" (21.5 cm) |



## NEEDLE HOLDERS

### Castro-Type Micro Needle Holders

#### CASTROVIEJO MICRO NEEDLE HOLDERS



357676DG

Round handle with Diama-Grip jaws.

| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 357656DG  | Straight, with lock, 1.2 mm, 7" (18 cm)   |
| 357676DG  | Straight, with lock, 1.2 mm, 8" (20.5 cm) |
| 357677DG  | Curved without lock, 1 mm, 8" (20.5 cm)   |

#### CASTROVIEJO NEEDLE HOLDERS



354964DG

Straight, with flat handle and lock. Diama-Grip jaws.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 354964DG  | 1.8 mm, 7-1/4" (18.5 cm) |
| 354965DG  | 8" (20.5 cm)             |

#### CASTROVIEJO TECHNIQUE® NEEDLE HOLDERS



354964

| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| 354964    | Straight, 13 x 1.4 mm jaw, 7-1/4" (18.5 cm), TC |
| 354965    | Straight, 13 x 1.4 mm jaw, 8" (20.5 cm), TC     |



## NEEDLE HOLDERS

### Castro-Type Micro Needle Holders

#### CASTROVIEJO MICRO NEEDLE HOLDERS



KM44219

Without lock.

| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| KM44218   | Curved, serrated, 5-1/2" (14 cm)   |
| KM44219   | Straight, serrated, 5-1/2" (14 cm) |
| KM44220   | Curved, smooth, 7" (18 cm)         |
| KM44221   | Curved, serrated, 7" (18 cm)       |
| KM44222   | Straight, smooth, 7" (18 cm)       |
| KM44223   | Straight, serrated, 7" (18 cm)     |

#### CASTROVIEJO MICRO NEEDLE HOLDERS TC



KM44210

Without lock.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| KM44210   | Curved, serrated, 5-1/2" (14 cm), TC   |
| KM44211   | Curved, serrated, 7" (18 cm), TC       |
| KM44212   | Straight, serrated, 5-1/2" (14 cm), TC |
| KM44213   | Straight, serrated, 7" (18 cm), TC     |
| KM44216   | Curved, smooth, 7" (18 cm), TC         |
| KM44217   | Straight, smooth, 7" (18 cm), TC       |

#### CASTROVIEJO MICRO NEEDLE HOLDERS



355559DG


Fine, straight, with round handle and lock. 1 mm  
Diam-Grip jaws.

| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 355558DG  | 9" (23 cm)    |
| 355559DG  | 10" (25.5 cm) |
| 355560DG  | 12" (30.5 cm) |

## NEEDLE HOLDERS

### Castro-Type Micro Needle Holders

#### CASTROVIEJO NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                                   |
|-----------|-----------------------------------------------|
| 152818    | Straight, 11 x 1.4 mm jaw, 6" (15 cm), TC     |
| 152819    | Curved, 12 x 1.4 mm jaw, 5-3/4" (14.5 cm), TC |

#### MICRO SURGERY NEEDLE HOLDERS



Flat handle, smooth jaws.

| CATALOG # | DESCRIPTION                                 |
|-----------|---------------------------------------------|
| 152825    | Straight, 12.5 x 1.4 mm jaw, 6" (15 cm), TC |
| 152827    | Curved, 13 x 1.4 mm jaw, 6" (15 cm), TC     |

#### CASTROVIEJO MICRO NEEDLE HOLDERS



With lock.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| KM44201   | Straight, serrated, 5-1/2" (14 cm)     |
| KM44200   | Curved, serrated, 5-1/2" (14 cm)       |
| KM44204   | Straight, smooth, 7" (18 cm)           |
| KM44202   | Curved, smooth, 7" (18 cm)             |
| KM44205   | Straight, serrated, 7" (18 cm)         |
| KM44203   | Curved, serrated, 7" (18 cm)           |
| KM44208   | Straight, serrated, 5-1/2" (14 cm), TC |
| KM44215   | Straight, smooth, 7" (18 cm), TC       |
| KM44214   | Curved, smooth, 7" (18 cm), TC         |
| KM44207   | Curved, serrated, 7" (18 cm), TC       |

## NEEDLE HOLDERS

### Castro-Type Micro Needle Holders

#### MICRO VASCULAR ROUND HANDLE NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                                                 |
|-----------|-------------------------------------------------------------|
| 357654    | Straight, without lock, 1.1 mm smooth jaw, 7-3/8" (18.5 cm) |
| 357655    | Curved, without lock, smooth, 7-3/8" (18.5 cm)              |
| 357656    | Straight, with lock, 1.2 mm smooth jaw, 7-1/2" (19 cm)      |
| 357657    | Curved, with lock, 0.6 mm smooth jaw, 7-3/8" (18.5 cm)      |
| 357667    | Straight, with lock, 7-3/8" (18.5 cm), TC                   |
| 357672    | Straight, with lock, 0.6 mm smooth jaw, 8-1/8" (20.5 cm)    |
| 357674    | Straight, with lock, 0.6 mm smooth jaw, 9" (23 cm)          |
| 357675    | Curved, with lock, 0.6 mm smooth jaw, 9" (23 cm)            |
| 357676    | Straight, with lock, 8" (20.5 cm), TC                       |
| 357677    | Curved, with lock, 1 mm jaw, 8" (20.5 cm), TC               |

#### BARRAQUER MICRO NEEDLE HOLDERS



With lock, round handle.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| KM44167   | Straight, smooth, 5-1/2" (14 cm)       |
| KM44163   | Curved, smooth, 5-1/2" (14 cm)         |
| KM44169   | Straight, smooth, 7" (18 cm)           |
| KM44165   | Curved, smooth, 7" (18 cm)             |
| KM44168   | Straight, smooth, 5-1/2" (14 cm), TC   |
| KM44164   | Curved, smooth, 5-1/2" (14 cm), TC     |
| KM44157   | Straight, serrated, 5-1/2" (14 cm), TC |
| KM44155   | Curved, serrated, 5-1/2" (14 cm), TC   |
| KM44170   | Straight, smooth, 7" (18 cm), TC       |
| KM44166   | Curved, smooth, 7" (18 cm), TC         |
| KM44158   | Straight, serrated, 7" (18 cm), TC     |
| KM44156   | Curved, serrated, 7" (18 cm), TC       |

## NEEDLE HOLDERS

CVT-Type Needle Holders

### W.M. LEMMON TECHNIQUE® NEEDLE HOLDERS



354916

| CATALOG # | DESCRIPTION                                       |
|-----------|---------------------------------------------------|
| 354915    | Delicate, 15.5 x 2.3 mm jaw, 5-3/4" (14.5 cm), TC |
| 354916    | 14.5 x 1.5 mm jaw, 6-3/4" (17 cm), TC             |
| 354917    | 15.5 x 2.3 mm jaw, 7-1/2" (19 cm), TC             |

### WANGENSTEEN NEEDLE HOLDER



152850

| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 152850    | 29.4 x 2.3 mm jaw, 10-3/4" (27.5 cm), TC |

## NEEDLE HOLDERS

CVT-Type Needle Holders

### PILLING® INTRACARDIAC TECHNIQUE® NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 354985    | 17 x 1.8 mm jaw, 7" (18 cm), TC       |
| 354986    | 11.5 x 1.6 mm jaw, 8" (20.5 cm), TC   |
| 354987    | 20.5 x 2 mm jaw, 9-1/4" (23.5 cm), TC |
| 354988    | 20.5 x 2 mm jaw, 10-1/4" (26 cm), TC  |

### PILLING® THIN VASCULAR AND CORONARY TECHNIQUE® NEEDLE HOLDERS

Delicately proportioned, slender needle holders.



| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 354961    | 16 x 1.7 mm jaw, 7" (18 cm), TC        |
| 354960    | 16 x 1.8 mm jaw, 6" (15 cm), TC        |
| 354962    | 20 x 1.8 mm jaw, 8-1/4" (21 cm), TC    |
| 354958    | 20.5 x 1.8 mm jaw, 9" (23 cm), TC      |
| 354959    | 20.5 x 2.2 mm jaw, 10-1/4" (26 cm), TC |



## NEEDLE HOLDERS

CVT-Type Needle Holders

### JAMISON TECHNIQUE® NEEDLE HOLDER



Long, slender needle holder. Jaws are lightly serrated.

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 354941    | 19.8 x 1.8 mm jaw, 9" (23 cm), TC |

### JULIAN NEEDLE HOLDER



| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 510250    | Straight, 21 x 2 mm jaw, 8" (20.5 cm), TC |

### PILLING® EXTRA-LIGHT TECHNIQUE® NEEDLE HOLDERS



Thinner shanks and box-locks than Pilling® Thin Needle Holders for reduced weight and more "feel."

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 354981    | 17 x 1.8 mm jaw, 6" (15 cm), TC   |
| 354982    | 17.5 x 1.6 mm jaw, 7" (18 cm), TC |



## NEEDLE HOLDERS

### CVT-Type Needle Holders

#### GERBODE TECHNIQUE® NEEDLE HOLDER



Medium-light construction. A strong full-sized needle holder with jaws tapered to 2 mm width for holding fine needles.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 343000    | 16.5 x 2 mm jaw, 10-1/2" (26.5 cm), TC |

#### HOWELL CORONARY TECHNIQUE® NEEDLE HOLDER



| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 354920    | 14 x 1.6 mm jaw, 7-1/4" (18.5 cm), TC |

#### HOWELL VALVE TECHNIQUE® NEEDLE HOLDER



| CATALOG # | DESCRIPTION                     |
|-----------|---------------------------------|
| 354928    | 19 x 2 mm jaw, 8" (20.5 cm), TC |

## NEEDLE HOLDERS

CVT-Type Needle Holders

### DEBAKEY® TECHNIQUE® NEEDLE HOLDERS



Carefully designed with light, short Tungsten Carbide jaws. Constructed in weight proportional to total length. May be purchased individually or as matched set of four.

| CATALOG # | DESCRIPTION                                             |
|-----------|---------------------------------------------------------|
| 354900    | No. 1. Length, 18.2 x 2 mm jaw, 7" (18 cm), TC          |
| 354902    | No. 2. Length, 21.5 x 1.8 mm jaw, 9-1/4" (23.5 cm), TC  |
| 354904    | No. 3. Length, 21.5 x 2.3 mm jaw, 10-1/2" (26.5 cm), TC |
| 354906    | No. 4. Length, 19 x 2.6 mm jaw, 11-3/4" (30 cm), TC     |
| 354908    | Matched set of four, TC                                 |

### FINOCHIETTO NEEDLE HOLDER



| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 152842    | 25.5 x 2.5 mm jaw, 10-1/2" (26.5 cm), TC |

## NEEDLE HOLDERS


### General Needle Holders

#### STRATTE NEEDLE HOLDERS



### CVT-Type Needle Holders

#### BAILEY CARDIAC TECHNIQUE® NEEDLE HOLDER



## NEEDLE HOLDERS

### General Needle Holders

#### SAROT NEEDLE HOLDERS



Weck® Pattern, bright finish. Parallel shanks with offset finger rings.

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 510467    | 16.5 x 1.8 mm jaw, 7" (18 cm), TC |
| 510469    | 17 x 1.8 mm jaw, 9" (23 cm), TC   |

#### SAROT NEEDLE HOLDER



| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 152796    | 19.5 x 1.8 mm jaw, 10-1/4" (26 cm), TC |

## NEEDLE HOLDERS

### General Needle Holders

#### HEANEY NEEDLE HOLDERS



#### HEANEY NEEDLE HOLDERS




#### MCALLISTER BULLDOG® NEEDLE HOLDER



## NEEDLE HOLDERS

### General Needle Holders

#### MASSON NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 511110    | 15.5 x 2.4 mm jaw, 7" (18 cm), TC    |
| 511112    | 16 x 3.2 mm jaw, 11-3/4" (30 cm), TC |

#### HEANEY BULLDOG® NEEDLE HOLDER



| CATALOG # | DESCRIPTION                                         |
|-----------|-----------------------------------------------------|
| 511251    | 19.3 x 3.2 mm curved Bulldog® Jaw, 8" (20.5 cm), TC |

#### HEANEY NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                                                      |
|-----------|------------------------------------------------------------------|
| 511250    | 20.5 x 2.6 mm curved jaw with cross serrations, 8" (20.5 cm), TC |



## NEEDLE HOLDERS

### General Needle Holders

#### MAYO HEGAR NEEDLE HOLDERS



Serrated, stainless steel jaw.

| CATALOG # | DESCRIPTION       |
|-----------|-------------------|
| KM41270   | 5-1/2" (14 cm)    |
| KM41272   | 6" (15 cm)        |
| KM41278   | 7" (18 cm)        |
| KM41284   | 8" (20.5 cm)      |
| KM41286   | 10-1/2" (26.5 cm) |
| KM41300   | 5" (12.5 cm), TC  |

#### MAYO HEGAR BROAD BULLDOG® NEEDLE HOLDERS



Weck® Pattern, bright finish, Bulldog® Jaw, broad.

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 510151    | 12 x 4 mm jaw, 7" (18 cm), TC     |
| 510161    | 17 x 3.5 mm jaw, 8" (20.5 cm), TC |

#### MASSON BULLDOG® NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                                   |
|-----------|-----------------------------------------------|
| 511101    | Heavy, 18.5 x 3 mm jaw, 10-1/4" (26 cm), TC   |
| 511102    | Light, 18.5 x 3.2 mm jaw, 10-1/4" (26 cm), TC |

## NEEDLE HOLDERS

### General Needle Holders


#### MAYO HEGAR NEEDLE HOLDERS



Weck® Pattern, bright finish, narrow.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 510100    | 14.6 x 2.6 mm jaw, 6" (15 cm), TC      |
| 510110    | 16.2 x 2.8 mm jaw, 7" (18 cm), TC      |
| 510120    | 15.7 x 3 mm jaw, 8" (20.5 cm), TC      |
| 510130    | 16.4 x 3.2 mm jaw, 9" (23 cm), TC      |
| 510140    | 16.3 x 3.5 mm jaw, 11-3/4" (30 cm), TC |

#### MAYO HEGAR NEEDLE HOLDERS



Weck® Pattern, bright finish, Bulldog® Jaw, narrow.

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 510101    | 15.5 x 3 mm jaw, 6" (15 cm), TC   |
| 510111    | 16 x 2.8 mm jaw, 7" (18 cm), TC   |
| 510121    | 19 x 2.7 mm jaw, 8" (20.5 cm), TC |
| 510131    | 20 x 3 mm jaw, 9" (23 cm), TC     |

## NEEDLE HOLDERS

### General Needle Holders

#### HEGAR BAUMGARTNER NEEDLE HOLDERS



Serrated stainless steel jaw.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| KM41162   | Serrated, 5" (12.5 cm)     |
| KM41164   | Serrated, 5-3/4" (14.5 cm) |

#### MAYO HEGAR NEEDLE HOLDERS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                      |
|-----------|--------------------------------------------------|
| 152720    | 14.5 x 3 mm jaw, 6-1/2" (16.5 cm), TC            |
| 152402    | 23 x 2.8 mm jaw, 7-1/4" (18.5 cm), TC            |
| 152722    | 16.5 x 3.3 mm jaw, 7-1/2" (19 cm), TC            |
| 152724    | 17 x 2.8 mm jaw, 8" (20.5 cm), TC                |
| 152406    | 17 x 2.8 mm jaw, 9-1/4" (23.5 cm), TC            |
| 152772    | Left-handed, 13 x 2.4 mm jaw, 6-1/4" (16 cm), TC |

## NEEDLE HOLDERS

### General Needle Holders

#### CRILE MURRAY BULLDOG® NEEDLE HOLDERS



510241

Weck® Pattern, bright finish. Long, slender Bulldog® Jaw.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 510241    | 15.4 x 1.8 mm jaw, 6-1/4" (16 cm), TC |
| 510251    | 15.5 x 2.4 mm jaw, 7" (18 cm), TC     |

#### CRILE WOOD BULLDOG® NEEDLE HOLDERS



510211

Weck® Pattern, bright finish. Bulldog® Jaw.

| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 510201    | Light, 6-1/4" (16 cm), TC |
| 510211    | Heavy, 6-1/4" (16 cm), TC |

## NEEDLE HOLDERS

### General Needle Holders

#### CRILE WOOD NEEDLE HOLDERS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 152771    | Left-handed, 7" (18 cm), TC  |
| 152780    | 6" (15 cm), TC               |
| 152781    | Smooth, 6" (15 cm), TC       |
| 152782    | 7" (18 cm), TC               |
| 152784    | 8" (20.5 cm), TC             |
| 152786    | 9" (23 cm), TC               |
| 152790    | Curved, 7" (18 cm), TC       |
| 152792    | Curved, 9-1/4" (23.5 cm), TC |

#### CRILE WOOD NEEDLE HOLDERS



Weck® Pattern, bright finish. Delicate, tapered tip.

| CATALOG # | DESCRIPTION         |
|-----------|---------------------|
| 510400    | 6" (15 cm), TC      |
| 510405    | 7" (18 cm), TC      |
| 510410    | 8" (20.5 cm), TC    |
| 510415    | 9" (23 cm), TC      |
| 510420    | 10-1/4" (26 cm), TC |
| 510430    | 12-1/4" (31 cm), TC |

## NEEDLE HOLDERS

### General Needle Holders


#### RYDER BULLDOG® NEEDLE HOLDERS



Weck® Pattern, bright finish. French eye with fine serrations and delicate, tapered, Bulldog® Jaw.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 511115    | 9.5 x 2 mm jaw, 5" (12.5 cm), TC     |
| 511130    | 12 x 2 mm jaw, 6" (15 cm), TC        |
| 511132    | 13 x 1.7 mm jaw, 7" (18 cm), TC      |
| 511134    | 13 x 2 mm jaw, 7-3/4" (19.5 cm), TC  |
| 511138    | 12.5 x 1.8 mm jaw, 10" (25.5 cm), TC |

#### RYDER NEEDLE HOLDERS



Weck® Pattern, bright finish. French eye with fine serrations and delicate, tapered jaw.

| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 511131    | 12.5 x 2 mm jaw, 6-1/4" (16 cm), TC |
| 511133    | 13 x 2.2 mm jaw, 7" (18 cm), TC     |
| 511135    | 12.5 x 1.7 mm jaw, 8" (20.5 cm), TC |
| 511139    | 12 x 2 mm jaw, 10" (25.5 cm), TC    |

#### RYDER EXTRA-FINE TECHNIQUE® NEEDLE HOLDERS



Pilling® Pattern, satin finish. Sculptured 1 mm wide jaws designed for the finest vascular, plastic and nerve suture needles. Lightweight throughout (average weight 1 oz).

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 354990    | 11.5 x 1.3 mm jaw, 5-1/4" (13.5 cm), TC |
| 354992    | 12.5 x 1.2 mm jaw, 6" (15 cm), TC       |
| 354994    | 14.3 x 1.2 mm jaw, 7" (18 cm), TC       |
| 354996    | 14.5 x 1.2 mm jaw, 8" (20.5 cm), TC     |
| 354998    | 15 x 1.3 mm jaw, 9" (23 cm), TC         |



## NEEDLE HOLDERS

### General Needle Holders

#### OLSEN-HEGAR NEEDLE HOLDERS



Combination needle holder and suture scissors.  
KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| KM41424   | Delicate, 4-3/4" (12 cm)       |
| KM41428   | 5" (12.5 cm)                   |
| KM41440   | Delicate, 4-1/2" (11.5 cm), TC |
| KM41452   | 6-1/2" (16.5 cm), TC           |

#### RYDER FRENCH EYE NEEDLE HOLDERS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 152760    | 20 x 2 mm jaw, 6" (15 cm), TC         |
| 152762    | 20.5 x 2 mm jaw, 7-1/4" (18.5 cm), TC |
| 152764    | 20 x 2 mm jaw, 8" (20.5 cm), TC       |

## NEEDLE HOLDERS

### General Needle Holders

#### WEBSTER NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| KM41068   | Smooth, stainless steel jaw, 5" (12.5 cm) |

#### OLSEN-HEGAR NEEDLE HOLDER AND SCISSORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                      |
|-----------|--------------------------------------------------|
| 510440    | 11 x 3.2 mm jaw, 5-3/4" (14.5 cm), TC            |
| 510441    | 15 x 3.1 mm jaw, 6" (15 cm), stainless steel jaw |
| 510445    | 11.5 x 3.2 mm jaw, 7" (18 cm), TC                |

## NEEDLE HOLDERS

### General Needle Holders


#### WEBSTER BULLDOG® NEEDLE HOLDER



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                                    |
|-----------|----------------------------------------------------------------|
| 510321    | Delicate, smooth, 13.5 x 2 mm Bulldog® jaw, 4-3/4" (12 cm), TC |

#### WEBSTER NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                                          |
|-----------|------------------------------------------------------|
| 152448    | 16 x 2.2 mm smooth stainless steel jaw, 5" (12.5 cm) |

## NEEDLE HOLDERS

### General Needle Holders

#### WEBSTER NEEDLE HOLDERS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 152808    | 10.2 x 1.8 mm smooth jaw, 4-1/2" (11.5 cm), TC |

#### WEBSTER NEEDLE HOLDERS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                           |
|-----------|-------------------------------------------------------|
| 510325    | Delicate, smooth, 13 x 2.2 mm jaw, 4-3/4" (12 cm), TC |

## NEEDLE HOLDERS

### General Needle Holders

#### NEIVERT NEEDLE HOLDER



Off-set ring handle.

| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 152805    | 16 x 2 mm jaw, 5-1/4" (13.5 cm), TC |

#### PAR NEEDLE HOLDERS



| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 510335    | 13 x 2 mm jaw, 5" (12.5 cm), TC     |
| 510395    | 11 x 1.8 mm jaw, 5-1/2" (14 cm), TC |

## NEEDLE HOLDERS

### General Needle Holders

#### HALSEY NEEDLE HOLDERS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| 152802    | Smooth, 11.5 x 1.8 mm jaw, 5-1/4" (13.5 cm), TC |
| 152800    | Serrated, 11 x 2.2 mm jaw, 5-1/4" (13.5 cm), TC |

#### JOHNSON BULLDOG® NEEDLE HOLDER



| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 510311    | 9.7 x 2.5 mm, Bulldog® Jaw, 4-3/4" (12 cm), TC |

#### JOHNSON NEEDLE HOLDER



| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 510312    | 12 x 2 mm jaw, 4-3/4" (12 cm), TC |



## NEEDLE HOLDERS

### General Needle Holders

#### GILLIES NEEDLE HOLDER SCISSORS



Large offset thumb ring.

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 152496    | 2.2 mm jaw, 6-1/4" (16 cm), TC |

#### HALSEY NEEDLE HOLDERS



KMedic® Pattern, satin finish. Stainless steel jaw.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| KM41064   | Smooth, 2.5 mm jaw, 5" (12.5 cm)       |
| KM41066   | Serrated, 2.2 mm jaw, 5-1/4" (13.5 cm) |

## NEEDLE HOLDERS

### General Needle Holders

#### DERF NEEDLE HOLDERS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 152798    | 7.5 x 1.7 mm jaw, 4-1/2" (11.5 cm), TC |

#### DERF NEEDLE HOLDER STANDARD – STAINLESS STEEL JAW



| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 443120    | 9.5 x 2 mm jaw, 5" (12.5 cm) |

## NEEDLE HOLDERS

### General Needle Holders

#### BROWN NEEDLE HOLDER



| CATALOG # | DESCRIPTION                                                                                      |
|-----------|--------------------------------------------------------------------------------------------------|
| 055200    | Regular, convex jaw with fine oblique serrations, 15 x 2.7 mm, 5" (12.5 cm), stainless steel jaw |

#### CONVERSE NEEDLE HOLDER



| CATALOG # | DESCRIPTION                                                     |
|-----------|-----------------------------------------------------------------|
| 443100    | Large rings, 14 x 2.3 mm jaw, 5" (12.5 cm), stainless steel jaw |

#### DERF NEEDLE HOLDERS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 055135    | 10.2 x 2.4 mm jaw, 4-3/4" (12 cm), TC |

## NEEDLE HOLDERS

### General Needle Holders

#### BAUMGARTNER NEEDLE HOLDER



| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 152710    | 11.7 x 2.2 mm jaw, 5-1/4" (13.5 cm), TC |

#### BAUMGARTNER BULLDOG® NEEDLE HOLDER



| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 510181    | 1.6 mm jaw, 5" (12.5 cm), TC |

## SCISSORS

Super-Cut and Ultimate-Cut Scissors

### ULTIMATE-CUT STEVENS TENOTOMY SCISSORS



| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 463920    | Straight, 4-1/2" (11.5 cm) |
| 463921    | Curved, 4-1/2" (11.5 cm)   |

### IRIS SCISSORS



Sharp/Sharp.

| CATALOG # | DESCRIPTION                                        |
|-----------|----------------------------------------------------|
| 463930    | Straight with one serrated blade, 4-1/2" (11.5 cm) |
| 463932    | Curved with one serrated blade, 4-1/2" (11.5 cm)   |

PR

## SCISSORS

Super-Cut and Ultimate-Cut Scissors

### ULTIMATE-CUT MAYO SCISSORS



| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 463900    | Straight, 6-3/4" (17 cm) |
| 463901    | Curved, 6-3/4" (17 cm)   |
| 463903    | Curved, 9-1/8" (23 cm)   |



### ULTIMATE-CUT METZENBAUM SCISSORS



| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 463910    | Straight, 5-3/4" (14.5 cm) |
| 463911    | Curved, 5-3/4" (14.5 cm)   |
| 463912    | Straight, 6-7/8" (17.5 cm) |
| 463913    | Curved, 7" (18 cm)         |
| 463915    | Curved, 9" (23 cm)         |

### ULTIMATE-CUT POTTS SCISSORS



| CATALOG # | DESCRIPTION                                                  |
|-----------|--------------------------------------------------------------|
| 463941    | Angled 45 degrees, with one serrated blade, 7-1/4" (18.5 cm) |



## SCISSORS

Super-Cut and Ultimate-Cut Scissors

### SUPER-CUT STEVENS SCISSORS



Blunt points. Delicate short narrow tips for dissecting fine tissue.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| KM35484   | Straight, 4-1/2" (11.5 cm) |
| KM35490   | Curved, 4-1/2" (11.5 cm)   |
| KM35494   | Curved, 5-1/4" (13.5 cm)   |

### SUPER-CUT IRIS SCISSORS



Sharp points.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| KM35450   | Straight, 3-1/2" (9 cm)    |
| KM35452   | Curved, 3-1/2" (9 cm)      |
| KM35454   | Straight, 4-1/2" (11.5 cm) |
| KM35460   | Curved, 4-1/2" (11.5 cm)   |

## SCISSORS

Super-Cut and Ultimate-Cut Scissors

### SUPER-CUT (BABY) METZENBAUM SCISSORS



Blunt points.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| KM35364   | Straight, 4-1/2" (11.5 cm) |
| KM35370   | Curved, 4-1/2" (11.5 cm)   |

### SUPER-CUT METZENBAUM SCISSORS



Blunt points.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| KM35380   | Straight, 6" (15 cm)         |
| KM35382   | Curved, 5-5/8" (14.5 cm)     |
| KM35392   | Straight, 7" (18 cm)         |
| KM35394   | Curved, 7" (18 cm)           |
| KM35396   | Curved, delicate, 7" (18 cm) |
| KM35406   | Curved, 8" (20.5 cm)         |
| KM35418   | Curved, 9" (23 cm)           |
| KM35430   | Curved, 11-1/4" (28.5 cm)    |

### SUPER-CUT JAMISON SCISSORS



| CATALOG # | DESCRIPTION                                        |
|-----------|----------------------------------------------------|
| KM35496   | Curved blunt points, very delicate, 6-1/4" (16 cm) |

## SCISSORS

### Super-Cut and Ultimate-Cut Scissors

#### SUPER-CUT OPERATING SCISSORS



Sharp/Blunt.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| KM35220   | Straight, 5-1/2" (14 cm) |
| KM35250   | Curved, 5-3/4" (14.5 cm) |

#### SUPER-CUT MAYO SCISSORS



Beveled blades, blunt points.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| KM35280   | Straight, 5-3/4" (14.5 cm) |
| KM35286   | Curved, 5-3/4" (14.5 cm)   |
| KM35292   | Straight, 6-3/4" (17 cm)   |
| KM35298   | Curved, 6-3/4" (17 cm)     |
| KM35316   | Straight, 9-1/4" (23.5 cm) |
| KM35322   | Curved, 9" (23 cm)         |
| KM35340   | Curved, 11" (28 cm)        |

#### ULTIMATE-CUT COOLEY SCISSORS



| CATALOG # | DESCRIPTION                                       |
|-----------|---------------------------------------------------|
| 463950    | Curved, with one serrated blade, 7-3/4" (19.5 cm) |

## SCISSORS

### Super-Cut and Ultimate-Cut Scissors

Our KMedic® Super-Cut Scissors and Pilling® Ultimate-Cut Scissors are manufactured with the finest surgical stainless steel, the latest heat treating methods and master craftsmanship.


Both style scissors are designed with one razor-sharp blade edge. The Super-Cut scissors' opposite blade has a standard smooth edge, while the Ultimate-Cut scissors' opposite blade incorporates a micro-serrated edge.

#### SUPER-CUT LISTER BANDAGE SCISSORS



| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| KM35210   | 5-3/4" (14.5 cm) |
| KM35212   | 7-1/4" (18.5 cm) |


#### SUPER-CUT HI-LEVEL BANDAGE SCISSORS



The delicate offset shanks and the sharp Super-Cut edge provide better leverage. The fine knob slides easily under finger bandages.

| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| KM35202   | Straight, angled on side, 5-1/2" (14 cm) |

#### SUPER-CUT LISTER BANDAGE SCISSORS WITH LARGE FINGER RING



Razor-sharp scissors ideal for cutting heavy material including: moleskin, drapes and stockinette. Not recommended for fiberglass or plaster splints. Oversized finger ring helps relieve pressure on the third metacarpal.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| KM35214   | Straight, 8" (20.5 cm) |

## SCISSORS

Delicate and Iris-Type Scissors

### WUTZLER SCISSORS



28 mm dissecting blades.

| CATALOG # | DESCRIPTION                                  |
|-----------|----------------------------------------------|
| 462540    | Straight, Blunt/Blunt blades, 5" (12.5 cm)   |
| 462541    | Curved, Blunt/Blunt blades, 5" (12.5 cm)     |
| 462550    | Straight, Sharp/Sharp blades, 4-3/4" (12 cm) |

### GOLDMAN-FOX SCISSORS

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| KM34648   | Curved, serrated, 5" (12.5 cm) |

## SCISSORS

Delicate and Iris-Type Scissors

### STEVENS TENOTOMY



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 462680    | Straight, 4-1/2" (11.5 cm) |
| 462681    | Curved, 4-1/2" (11.5 cm)   |
| 056110    | Curved, 5-1/8" (13 cm)     |

### STRABISMUS SCISSORS



Pilling® Pattern, satin finish. Rounded blades and tips.

| CATALOG # | DESCRIPTION                                            |
|-----------|--------------------------------------------------------|
| 144319    | Straight, 4-1/2" (11.5 cm), TC                         |
| 144321    | Curved, 4-1/2" (11.5 cm), TC                           |
| 144330    | Straight, round shank, 4-3/4" (12 cm), stainless steel |
| 144332    | Curved, round shank, 4-3/4" (12 cm), stainless steel   |

### STRABISMUS SCISSORS



Weck® Pattern, bright finish. Rounded 30 mm dissecting blades.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 462690    | Straight, 4-1/2" (11.5 cm) |
| 462691    | Curved, 4-1/2" (11.5 cm)   |



## SCISSORS

Delicate and Iris-Type Scissors

### STEVENS TENOTOMY



KMedic® Pattern, satin finish, Sharp/Sharp points.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| KM53174   | Straight, 4-3/8" (11 cm) |
| KM53176   | Curved, 4-3/8" (11 cm)   |

### STEVENS TENOTOMY



Pilling® Pattern, satin finish, slender style.

| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 144350    | Straight, 4-1/2" (11.5 cm)               |
| 144352    | Slightly curved blades, 4-1/2" (11.5 cm) |

### STEVENS TENOTOMY, DELICATE



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION          |
|-----------|----------------------|
| 442515    | Curved, 5" (12.5 cm) |

## SCISSORS

Delicate and Iris-Type Scissors

### RAGNELL SCISSORS, DELICATE



| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 442576    | Half-curved, 4-3/4" (12 cm) |

### RAGNELL SCISSORS



Curved, flat blunt tips.

| CATALOG # | DESCRIPTION  |
|-----------|--------------|
| KM34560   | 5" (12.5 cm) |
| KM34562   | 7" (18 cm)   |

### SHEA VEIN GRAFT SCISSORS




| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 442525    | Flat tips, 4-3/4" (12 cm) |

## SCISSORS

Delicate and Iris-Type Scissors

### PLASTIC SCISSORS



### PLASTIC SCISSORS



SCISSORS

Delicate and Iris-Type Scissors

LA GRANGE SCISSORS



KM34621

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| KM34621   | Curved, 4-1/2" (11.5 cm) |

PAR SCISSORS



464810

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 464810    | Curved, 4-1/2" (11.5 cm), TC |

## SCISSORS

Delicate and Iris-Type Scissors

### KELLY FISTULA SCISSORS




Smooth, Sharp/Sharp.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| KM34602   | Straight, 6-1/4" (16 cm) |
| KM34607   | Curved, 6-1/4" (16 cm)   |
| KM34612   | Angled, 6-1/4" (16 cm)   |




### DUFFIELD SCISSORS



Sharp/Sharp.

| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 461430    | Straight, 28 mm blades, 6-1/2" (16.5 cm) |
| 461431    | Curved, 28 mm blades, 6-1/2" (16.5 cm)   |

### KNAPP IRIS SCISSORS



One sharp and one blunt point, 4-1/2" (11.5 cm).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 144338    | Straight    |
| 144340    | Curved      |

## SCISSORS

Delicate and Iris-Type Scissors

### IRIS SCISSORS, TC



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 462632    | Straight, 4-1/2" (11.5 cm), TC |
| 462633    | Curved, 4-1/2" (11.5 cm), TC   |

### JOSEPH SCISSORS



| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 442510    | Straight, 5-3/4" (14.5 cm) |
| 442512    | Curved, 5-3/4" (14.5 cm)   |

### KAYE SCISSORS



| CATALOG # | DESCRIPTION                                         |
|-----------|-----------------------------------------------------|
| 464220    | Curved, 32 mm serrated blades, 4-1/2" (11.5 cm), TC |



## SCISSORS

Delicate and Iris-Type Scissors

### IRIS RIBBON SCISSORS



Slightly curved, Sharp/Sharp.

| CATALOG # | DESCRIPTION     |
|-----------|-----------------|
| KM53265   | 3-3/4" (9.5 cm) |
| KM53267   | 4" (10 cm)      |

### IRIS SCISSORS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                 |
|-----------|---------------------------------------------|
| 144300    | Straight, 4-1/2" (11.5 cm), stainless steel |
| 144310    | Curved, 4-1/2" (11.5 cm), stainless steel   |
| 144315    | Straight, 4-1/2" (11.5 cm), TC              |
| 144317    | Curved, 4-1/2" (11.5 cm), TC                |

### IRIS SCISSORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 462600    | Curved, 25 mm blades, 4-1/2" (11.5 cm)   |
| 462605    | Straight, 25 mm blades, 4-1/2" (11.5 cm) |
| 462640    | Straight, 30 mm blades, 5" (12.5 cm)     |
| 462641    | Curved, 28 mm blades, 5" (12.5 cm)       |

## SCISSORS

Delicate and Iris-Type Scissors

### GRADLE SCISSORS



Slightly curved, Sharp/Sharp.

| CATALOG # | DESCRIPTION     |
|-----------|-----------------|
| KM53316   | 3-3/4" (9.5 cm) |
| KM53318   | 4" (10 cm)      |

### IRIS MINIATURE SCISSORS




| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 055100    | Straight, 15 mm blades, 3-1/2" (9 cm) |

## SCISSORS


Delicate and Iris-Type Scissors

### IRIS SCISSORS WITH BARRAQUER HANDLE



| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 055502    | Curved, 18 mm blades, 5-3/4" (14.5 cm) |

### BLEPHAROPLASTY SCISSORS



| CATALOG # | DESCRIPTION                                                                |
|-----------|----------------------------------------------------------------------------|
| KM35510   | Curved flat tips, one serrated and one carbide blade, 4-1/2" (11.5 cm), TC |

PRIME

## SCISSORS

Delicate and Iris-Type Scissors

### DIGITATE®-HEIFETZ SCISSORS



066102

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 066102    | Curved, 5-3/4" (14.5 cm)   |
| 066100    | Straight, 5-3/4" (14.5 cm) |
| 066104    | Straight, 4-3/4" (12 cm)   |



### MCCLURE IRIS SCISSORS



144295

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 144295    | Angled, 4-3/4" (12 cm) |

## SCISSORS

### Metzenbaum-Type Scissors

#### LINCOLN-METZENBAUM SCISSORS



Rounded, narrow dissecting tips.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 460883    | Curved, 5-3/8" (13.5 cm) |
| 460873    | Curved, 7" (18 cm)       |


#### LINCOLN SCISSORS



Small, rounded dissecting blades.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 461310    | Straight, 6-1/2" (16.5 cm) |
| 461311    | Curved, 6-1/2" (16.5 cm)   |
| 462201    | Curved, 11" (28 cm)        |

#### STRULLY PROBED TIP SCISSORS



| CATALOG # | DESCRIPTION          |
|-----------|----------------------|
| 462530    | Curved, 8" (20.5 cm) |

## SCISSORS

### Metzenbaum-Type Scissors

#### JAMISON-METZENBAUM TENOTOMY SCISSORS



The standard Metzenbaum-style curved blades have specially narrowed tenotomy-type distal tips.

| CATALOG # | DESCRIPTION  |
|-----------|--------------|
| 342220    | 6" (15 cm)   |
| 342221    | 7" (18 cm)   |
| 342222    | 8" (20.5 cm) |
| 342223    | 9" (23 cm)   |

#### POTTS TENOTOMY SCISSORS



| CATALOG # | DESCRIPTION        |
|-----------|--------------------|
| 640280    | Curved, 6" (15 cm) |



## SCISSORS

### Metzenbaum-Type Scissors

#### NELSON LOBECTOMY SCISSORS



Pilling® Pattern, satin finish, medium weight with rounded Metzenbaum-type curved blades.

| CATALOG # | DESCRIPTION       |
|-----------|-------------------|
| 342242    | 10-1/4" (26 cm)   |
| 342243    | 11-1/4" (28.5 cm) |
| 342245    | 14-1/4" (36 cm)   |

#### NELSON SCISSORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 462148    | Straight, 8" (20.5 cm) |
| 462149    | Curved, 8" (20.5 cm)   |
| 462150    | Straight, 9" (23 cm)   |
| 462151    | Curved, 9" (23 cm)     |

#### NELSON SCISSORS, TC




| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 464800    | Curved, 11-1/8" (28.5 cm), TC |
| 464801    | Curved, 11-1/4" (28.5 cm), TC |

## SCISSORS

### Metzenbaum-Type Scissors

#### DEMARTEL SCISSORS



Very slim, light construction.

| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 352120    | Straight, sharp blades, 6-3/4" (17 cm)         |
| 352106    | Angled on flat, blunt blades, 6-3/4" (17 cm)   |
| 352109    | Curved, blunt blades, 6-3/4" (17 cm)           |
| 352129    | Curved, sharp blades, 6-3/4" (17 cm)           |
| 352100    | Straight, blunt blades, 7" (18 cm)             |
| 352103P   | Angled on flat, blunt tips, 7" (18 cm)         |
| 352123    | Angled on flat, sharp blades, 7" (18 cm)       |
| 352126    | Angled to side, sharp blades, 6-7/8" (17.5 cm) |

#### FERGUSON LONG DISSECTING SCISSORS



With Metzenbaum-type, curved blades.

| CATALOG # | DESCRIPTION     |
|-----------|-----------------|
| 342275    | 10-1/4" (26 cm) |

## SCISSORS

### Metzenbaum-Type Scissors

#### DEBAKEY® SCISSORS



Sturdy Metzenbaum-type scissors. Black rings for quick identification.

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 342160    | Curved, 5-3/4" (14.5 cm), TC  |
| 342162    | Curved, 7" (18 cm), TC        |
| 342166    | Curved, 9-1/8" (23 cm), TC    |
| 342167    | Curved, 10-1/4" (26 cm), TC   |
| 342168    | Curved, 11-1/4" (28.5 cm), TC |
| 342169    | Curved, 12" (30.5 cm), TC     |

#### DEBAKEY®-METZENBAUM DISSECTING SCISSORS



Curved blades. These scissors are proportionately matched in weight.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 342200    | Curved, 5-7/8" (15 cm)   |
| 342201    | Curved, 7" (18 cm)       |
| 342202    | Curved, 8-1/8" (20.5 cm) |
| 342203    | Curved, 9-1/8" (23 cm)   |
| 342204    | Curved, 10-1/4" (26 cm)  |
| 342206    | 12-1/4" (31 cm)          |

#### DEMARTEL ACUTELY ANGLED SCISSORS



Lightweight construction.

| CATALOG # | DESCRIPTION                                                            |
|-----------|------------------------------------------------------------------------|
| 352210    | Angled on flat 50 degrees, short, 16 mm blades, sharp tips, 7" (18 cm) |

## SCISSORS

### Metzenbaum-Type Scissors

#### METZENBAUM SCISSORS, TC



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 464700    | Straight, 5-3/4" (14.5 cm), TC |
| 464705    | Curved, 5-3/4" (14.5 cm), TC   |
| 464710    | Straight, 7" (18 cm), TC       |
| 464715    | Curved, 7" (18 cm), TC         |
| 464720    | Straight, 9" (23 cm), TC       |
| 464725    | Curved, 9" (23 cm), TC         |

#### METZENBAUM SCISSORS, DELICATE, TC



Delicate. Blunt/Blunt.

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 464701    | Straight, 5-3/4" (14.5 cm), TC |
| 464706    | Curved, 5-3/4" (14.5 cm), TC   |
| 464711    | Straight, 7" (18 cm), TC       |
| 464716    | Curved, 7" (18 cm), TC         |
| 464721    | Straight, 9" (23 cm), TC       |
| 464726    | Curved, 9" (23 cm), TC         |

## SCISSORS

### Metzenbaum-Type Scissors

#### METZENBAUM SCISSORS



KMedic® Pattern, satin finish, 6" and 8".

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| KM34458   | Straight, 6" (15 cm)   |
| KM34462   | Curved, 6" (15 cm)     |
| KM34472   | Straight, 8" (20.5 cm) |
| KM34474   | Curved, 8" (20.5 cm)   |

#### METZENBAUM SCISSORS, TC



Sharp/Sharp.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 464707    | Curved, 5-3/4" (14.5 cm), TC |
| 464712    | Straight, 7" (18 cm), TC     |
| 464718    | Curved, 7" (18 cm), TC       |
| 464722    | Straight, 9" (23 cm), TC     |
| 464727    | Curved, 9" (23 cm), TC       |

## SCISSORS

Metzenbaum-Type Scissors

### METZENBAUM SCISSORS



Weck® Pattern, bright finish.


| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 460820    | Straight, 5" (12.5 cm)      |
| 460821    | Curved, 5" (12.5 cm)        |
| 460840    | Straight, 5-3/8" (13.5 cm)  |
| 460841    | Curved, 5-3/8" (13.5 cm)    |
| 460880    | Straight, 5-3/8" (13.5 cm)  |
| 460881    | Curved, 5-3/8" (13.5 cm)    |
| 460870    | Straight, 7" (18 cm)        |
| 460871    | Curved, 7" (18 cm)          |
| 460884    | Straight, 9" (23 cm)        |
| 460885    | Curved, 9" (23 cm)          |
| 462105    | Straight, 10-3/8" (26.5 cm) |
| 462110    | Curved, 10-1/4" (26 cm)     |
| 460889    | Curved, 11" (28 cm)         |



## SCISSORS

### Metzenbaum-Type Scissors

#### METZENBAUM SCISSORS, DELICATE



| CATALOG # | DESCRIPTION                                              |
|-----------|----------------------------------------------------------|
| 460830    | Curved, 29 mm delicate blades with cut-out, 5" (12.5 cm) |

#### METZENBAUM SCISSORS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 141460    | Straight, 7" (18 cm)   |
| 141462    | Curved, 7" (18 cm)     |
| 141467    | Curved, 8-1/4" (21 cm) |
| 141469    | Curved, 9" (23 cm)     |

## SCISSORS

### Metzenbaum-Type Scissors

#### BABY METZENBAUM SCISSORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 056310    | Curved, 5-3/4" (14.5 cm) |

#### LAHEY (BABY METZ) SCISSORS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 141450    | Straight, 5-3/4" (14.5 cm) |
| 141452    | Curved, 5-3/4" (14.5 cm)   |

#### METZENBAUM SCISSORS, DELICATE



Pilling® Pattern, satin finish. Delicate curved.

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 141448    | 5-1/2" (14 cm) |
| 352235    | 7" (18 cm)     |
| 352237    | 8" (20.5 cm)   |

## SCISSORS

### Mayo-Type Scissors

#### THOREK SCISSORS



Medium weight construction, round surface, acutely curved blades.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 342270    | 7-1/4" (18.5 cm) |
| 342272    | 10" (25.5 cm)    |

#### MARTIN CARTILAGE SCISSORS



| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 790320    | Curved, serrated blades, 8" (20.5 cm) |

#### DEAN SCISSORS

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| KM34642   | Angled, serrated, 7" (18 cm) |

#### GORNEY RHYTIDECTOMY SCISSORS



| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| 442560    | Straight, serrated, 6-3/4" (17 cm) |

## SCISSORS

### Mayo-Type Scissors

#### SISTRUNK FINE SCISSORS



Delicate, 5-1/2" (14 cm).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| KM34348   | Straight    |
| KM34354   | Curved      |

#### SISTRUNK HEAVY SCISSORS



| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| KM34356   | Straight, 5-1/2" (14 cm) |

#### JORGENSEN SCISSORS



Heavy, rigid body with acutely curved blades that taper to a rounded narrow tip.

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 342273    | 8-3/4" (22 cm) |

## SCISSORS

Mayo-Type Scissors

### DOYEN SCISSORS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 141395    | Straight, 7-1/4" (18.5 cm) |

### DOYEN FERGUSON SCISSORS



| CATALOG # | DESCRIPTION                                  |
|-----------|----------------------------------------------|
| 460720    | Angled on flat, 59 mm blades, 7-1/8" (18 cm) |

### SIMS SCISSORS



8" (20.5 cm).

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| KM67164   | Straight, Sharp/Sharp blades |
| KM67166   | Straight, Sharp/Blunt blades |
| KM67168   | Straight, Blunt/Blunt blades |
| KM67170   | Curved, Sharp/Sharp blades   |
| KM67172   | Curved, Sharp/Blunt blades   |
| KM67174   | Curved, Blunt/Blunt blades   |

## SCISSORS

Mayo-Type Scissors

### CHURCH SCISSORS



Lightly constructed, long scissors with a beveled surface. Mayo style blades. Blunt, dissecting tips.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 352250    | Curved, 6-3/4" (17 cm) |
| 352240    | Straight, 7" (18 cm)   |

### COOLEY-MAYO CURVED SCISSORS



| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 352090    | 7-1/2" (19 cm) |

### DOYEN ABDOMINAL SCISSORS



KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION          |
|-----------|----------------------|
| KM34334   | Straight, 7" (18 cm) |
| KM34336   | Curved, 7" (18 cm)   |



## SCISSORS

### Mayo-Type Scissors

#### MAYO NOBLE SCISSORS



6-1/2" (16.5 cm).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 460600    | Straight    |
| 460601    | Curved      |

#### MAYO-HARRINGTON SCISSORS



11" (28 cm).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 462120    | Straight    |
| 462121    | Curved      |

#### MAYO-STILLE SCISSORS



| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| KM34306   | Straight, 5-1/2" (14 cm) |
| KM34308   | Curved, 5-1/2" (14 cm)   |
| KM34310   | Straight, 6-3/4" (17 cm) |
| KM34312   | Curved, 6-3/4" (17 cm)   |

## SCISSORS

Mayo-Type Scissors

### MAYO DISSECTING SCISSORS, TC



Pilling® Pattern, satin finish. Beveled blade.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 142220    | Curved, 5-1/2" (14 cm), TC |
| 142235    | Curved, 9" (23 cm), TC     |
| 142237    | Straight, 9" (23 cm), TC   |

### MAYO DISSECTING SCISSORS, TC



Weck® Pattern, bright finish. Round blade.

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 464605    | Curved, 5-3/4" (14.5 cm), TC   |
| 464615    | Curved, 6-1/2" (16.5 cm), TC   |
| 464600    | Straight, 5-3/4" (14.5 cm), TC |
| 464610    | Straight, 6-1/2" (16.5 cm), TC |

## SCISSORS

Mayo-Type Scissors

### MAYO ROUND SURFACE BLADE SCISSORS



| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 141332    | Straight, 6-3/4" (17 cm) |

### MAYO DISSECTING SCISSORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 460410    | Curved, 43 mm blades, 5-1/2" (14 cm)       |
| 460400    | Straight, 43.5 mm blades, 5-5/8" (14.5 cm) |
| 460420    | Straight, 51.5 mm blades, 6-3/4" (17 cm)   |
| 460430    | Curved, 49.5 mm blades, 6-3/4" (17 cm)     |
| 460620    | Straight, 64 mm blades, 8-3/4" (22 cm)     |
| 460630    | Curved, 64 mm blades, 8-3/4" (22 cm)       |

### MAYO DISSECTING SCISSORS



Weck® Pattern, bright finish, lightweight, 6-1/2" (16.5 cm).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 460460    | Straight    |
| 460461    | Curved      |

## SCISSORS

### Operating Scissors

#### OPERATING SCISSORS, TC



5-1/2" (14 cm), TC.

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 464500    | Sharp/Blunt |
| 464502    | Sharp/Sharp |



#### OPERATING SCISSORS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 141146    | Sharp/Blunt, straight, 5-1/2" (14 cm)   |
| 141106    | Sharp/Sharp, straight, 5-3/4" (14.5 cm) |

### Mayo-Type Scissors

#### MAYO BEVELED BLADE SCISSORS



| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 141310    | Straight, 5-3/4" (14.5 cm) |
| 141315    | Curved, 5-3/4" (14.5 cm)   |
| 141312    | Straight, 6-3/4" (17 cm)   |
| 141317    | Curved, 6-3/4" (17 cm)     |

## SCISSORS

### Operating Scissors

#### OPERATING SCISSORS, B/B



Blunt/Blunt, stainless steel.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 460165    | Straight, 5" (12.5 cm)     |
| 460219    | Curved, 5-1/2" (14 cm)     |
| 460170    | Straight, 5-3/4" (14.5 cm) |
| 460179    | Straight, 6" (15 cm)       |
| 460177    | Straight, 7" (18 cm)       |

#### OPERATING SCISSORS, S/B



Sharp/Blunt, stainless steel.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 460100    | Straight, 4-3/4" (12 cm) |
| 460105    | Straight, 5" (12.5 cm)   |
| 460110    | Straight, 5-1/2" (14 cm) |
| 460220    | Curved, 5-1/2" (14 cm)   |
| 460175    | Straight, 6" (15 cm)     |
| 460115    | Straight, 6-3/4" (17 cm) |

#### OPERATING SCISSORS, S/S



Sharp/Sharp, stainless steel.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 460245    | Curved, 5-1/2" (14 cm)   |
| 460135    | Straight, 5" (12.5 cm)   |
| 460140    | Straight, 5-1/2" (14 cm) |

## SCISSORS

### Wire Cutting Scissors

#### WIRE CUTTING SCISSORS



KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                             |
|-----------|---------------------------------------------------------|
| KM35200   | Angular carbide, one serrated blade, 4-3/4" (12 cm), TC |



#### VERNON WIRE SCISSORS



Straight, serrated blades.

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 790310    | 32 mm blades, 5-3/4" (14.5 cm) |
| 790315    | 42 mm blades, 7-1/2" (19 cm)   |

### Operating Scissors

#### DELICATE OPERATING SCISSORS



Sharp/Blunt.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 460180    | Straight, 5" (12.5 cm) |



## SCISSORS

### Wire Cutting Scissors

#### HEATH SUTURE AND WIRE CUTTING SCISSORS



| CATALOG # | DESCRIPTION                                         |
|-----------|-----------------------------------------------------|
| 520120    | Curved sideways, serrated, 23 mm blades, 6" (15 cm) |

#### ROGER WIRE CUTTING SCISSORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                                    |
|-----------|----------------------------------------------------------------|
| 520130    | One serrated blade, 15 mm blades angled sideways, 5" (12.5 cm) |

#### WIRE SCISSORS



Pilling® Pattern, satin finish. Angular with notch for wire.

| CATALOG # | DESCRIPTION                                        |
|-----------|----------------------------------------------------|
| 142515    | Stainless steel blades, one serrated, 5" (12.5 cm) |
| 142518    | TC blades, 5" (12.5 cm)                            |

## SCISSORS

### Stitch and Suture Removal Scissors

#### SUTURE REMOVAL SCISSORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                   |
|-----------|-----------------------------------------------|
| 056120    | Curved, pointed, 22 mm blades, 4-3/4" (12 cm) |



### Wire Cutting Scissors

#### BANTAM WIRE CUTTING SCISSORS



Serrated 18 mm blades.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 520150    | Straight, 4-3/8" (11 cm) |
| 520152    | Curved, 4" (10 cm)       |

## SCISSORS

Stitch and Suture Removal Scissors

### SPENCER STITCH SCISSORS



KMedic® Pattern, satin finish. Delicate, straight.

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| KM35068   | 3-1/2" (9 cm)  |
| KM35092   | 4-3/4" (12 cm) |

### SPENCER STITCH SCISSORS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 520135    | Slightly curved, 20 mm blades, 3-1/2" (9 cm)   |
| 520140    | Straight, serrated, 30 mm blades, 5" (12.5 cm) |

## SCISSORS

### Stitch and Suture Removal Scissors

#### SHORTBENT STITCH SCISSORS



| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| KM35080   | Curved on side, delicate, 3-1/2" (9 cm) |



#### SPENCER SCISSORS



Pilling® Pattern, satin finish. Very delicate for fine sutures.

| CATALOG # | DESCRIPTION  |
|-----------|--------------|
| 142312    | 5" (12.5 cm) |

## SCISSORS

Stitch and Suture Removal Scissors

### LITTLER STITCH SCISSORS



| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| KM35078   | Curved, 4-5/8" (11.5 cm) |

### NORTHBENT STITCH SCISSORS




| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| KM35104   | Curved on side, 5-1/8" (13 cm) |

PK

## SCISSORS

### Stitch and Suture Removal Scissors

#### HILL STITCH SCISSORS



KM35084

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| KM35084   | Curved, 4-3/8" (11 cm) |

#### LITTAUER STANDARD SUTURE SCISSORS



142325

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 142325    | 5-1/2" (14 cm) |

#### LITTAUER STITCH SCISSORS



KM35120

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| KM35120   | Straight, 5-1/2" (14 cm) |



## SCISSORS

Bandage Scissors and Shears

### MAKI™ SCISSORS



Uniquely designed teeth to grip bottle caps, needles, tubing.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 463820    | 5-3/4" (14.5 cm) |

### UNIVERSAL SHEARS



Stainless steel, comfortable polymer handles, serrated.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 142452    | 7-1/2" (19 cm)              |
| KM11202   | Red handles, 7-1/2" (19 cm) |

## SCISSORS

Bandage Scissors and Shears

### LISTER BANDAGE SCISSORS



One large finger ring.

| CATALOG # | DESCRIPTION  |
|-----------|--------------|
| 463840    | 8" (20.5 cm) |



### MOLESKIN SCISSORS



One blade serrated.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| KM69820   | Straight, 7" (18 cm)   |
| KM69824   | Curved, 7-1/2" (19 cm) |

## SCISSORS

Bandage Scissors and Shears

### LISTER BANDAGE SCISSORS



KM31680

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| KM31680   | 4-1/2" (11.5 cm) |

### LISTER BANDAGE SCISSORS



142415

Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 142415    | 5-3/4" (14.5 cm) |
| 142420    | 7-1/4" (18.5 cm) |
| 142422    | 8" (20.5 cm)     |

### LISTER BANDAGE SCISSORS



463810

Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 463810    | 5-3/4" (14.5 cm) |
| 463830    | 7-1/4" (18.5 cm) |

## SCISSORS

### Bandage Scissors and Shears

#### HI-LEVEL BANDAGE SCISSORS



The delicate offset shanks provide better leverage. The fine knob slides easily under finger and toe bandages.

| CATALOG # | DESCRIPTION                                               |
|-----------|-----------------------------------------------------------|
| KM31792   | Shanks angled on side, one serrated blade, 5-1/2" (14 cm) |

#### KMEDIC® UTILITY SHEARS



Effective for cloth or cotton wadding. One smooth and one serrated blade.

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| KM31836   | 8-1/4" (21 cm) |

#### KNOWLES BANDAGE SCISSORS



For finger and toe bandages.

| CATALOG # | DESCRIPTION                                  |
|-----------|----------------------------------------------|
| KM31804   | Straight, one serrated blade, 5-1/2" (14 cm) |

## SCISSORS

Bandage Scissors and Shears

### ESMARCH BANDAGE SCISSORS



One large finger ring.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 142425    | 7-3/4" (19.5 cm) |

### ESMARCH PLASTER SHEARS



| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 463885    | 7-7/8" (20 cm) |

### HERCULES PLASTER SCISSORS



For lightweight plaster or aluminum casts.

| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| KM31834   | Angled, one blade serrated, 7-1/2" (19 cm) |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps

#### STONE INTESTINAL CLAMPS



Handle may be used at any angle. Narrow longitudinally serrated blades. Set consists of 2 clamps with tip locks and 2 applying forceps.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 185402    | Set. With 4-3/8" (11 cm) length clamps |
| 185404    | Clamp with tip lock, 3-3/8" (8.5 cm)   |

#### STRELINGER RIGHT ANGLE CLAMP




| CATALOG # | DESCRIPTION                                                                        |
|-----------|------------------------------------------------------------------------------------|
| 353872    | Large size, 3-1/4" (8.5 cm) DeBakey®-type jaws angled 90 degrees, 9-5/8" (24.5 cm) |



## FORCEPS AND CLAMPS


### Gastrointestinal Forceps and Clamps

#### PENNINGTON FORCEPS



| CATALOG # | DESCRIPTION                                            |
|-----------|--------------------------------------------------------|
| 475250    | Serrated, triangular, 14 x 11 mm jaw, 6-1/8" (15.5 cm) |
| 616370    | Serrated, triangular, 16 x 12.5 mm jaw, 8-1/4" (21 cm) |

#### PILLING® GUT AND COLON CLAMPS




90 degree angled, 1 x 2 DeBakey®-type jaws, curved shafts.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 183400    | 2-3/4" (7 cm) jaws, 9-3/8" (24 cm)   |
| 183410    | 3-5/8" (9 cm) jaws, 9-5/8" (24.5 cm) |

## FORCEPS AND CLAMPS


### Gastrointestinal Forceps and Clamps

#### MAYO (GUYON) VESSEL CLAMP



| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 172560    | 9-1/4" (23.5 cm) |

#### MAYO-ROBSON INTESTINAL CLAMPS




Longitudinal serrations.

| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 185220    | Straight, 10-1/4" (26 cm) |
| 185222    | Curved, 10" (25.5 cm)     |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps


#### KOCHER INTESITAL CLAMPS



Soft blades, longitudinal serrations.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 185210    | Straight, 10" (25.5 cm)     |
| 185212    | Curved, 9-7/8" (25 cm)      |
| 185214    | Straight, 10-3/4" (27.5 cm) |

#### LELAND JONES CLAMPS



Delicate DeBakey® Jaws, 2.4 mm wide. LIFE-LOK® Box-Lock. Versatile clamps which have been found useful in small bowel procedures.

| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 353490    | Straight, 8-1/8" (20.5 cm)          |
| 353491    | Angled 15 degrees, 8-1/8" (20.5 cm) |
| 353492    | Angled 30 degrees, 8" (20.5 cm)     |

#### MIXTER FORCEPS




Long model. Very useful for deep dissection and as a bronchus clamp.

| CATALOG # | DESCRIPTION                                                            |
|-----------|------------------------------------------------------------------------|
| 342606    | Strongly curved, serrated (3 cm) tip length, 10-3/4" (27.5 cm)         |
| 342605    | Strongly curved, DeBakey® Jaw, 1-3/8" (3.5 cm) tip length, 11" (28 cm) |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps


#### HAYES ANTERIOR RESECTION CLAMPS



70 degree-angled, 1 x 2 DeBakey® Jaws, angular shafts.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 183975    | 2-3/8" (6 cm) jaws, 9-7/8" (25 cm)    |
| 183977    | 3-1/8" (8 cm) jaws, 10-1/2" (26.5 cm) |

#### KLINTMALM™ LIVER TRANSPLANT CLAMPS



The KLINTMALM™ Liver Transplant Clamp is designed to securely occlude the suprahepatic inferior vena cava of the recipient during anastomosis of the donor liver.

- Unique "teething" design of a single row against a double row of teeth helps prevent slippage and scissoring of the jaw.
- An even-parallel closure of the jaw helps prevent blood seepage and helps minimize the risk of the vena cava slipping from the clamp.
- Jaw curvature is designed to meet the anatomical requirements of the upper abdomen and to reduce trauma to the liver.


A standard vascular clamp on the cuff of the vena cava and diaphragm may not always properly close due to the amount of tissue in the clamp. The use of a standard vascular clamp may result in clamp insecurity and slippage leading to the possible retraction of the vena cava into the diaphragm, and visual obstruction caused by pooling blood.

| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 184001    | Small 2-1/2" (6.5 cm) jaw, 7-1/8" (18 cm)  |
| 184002    | Medium 3-1/8" (8 cm) jaw, 8-3/8" (21.5 cm) |
| 184003    | Large 3-3/8" (8.5 cm) jaw, 9-7/8" (25 cm)  |

## FORCEPS AND CLAMPS


### Gastrointestinal Forceps and Clamps

#### DENNIS ANASTOMOSIS CLAMP




| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 185030    | 54 mm jaws, 9-1/4" (23.5 cm) |

#### DESJARDIN GALLSTONE GRASPING FORCEPS



| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 172625    | 9-1/4" (23.5 cm) |

#### DOYEN INTESTINAL CLAMPS



Soft blades, diagonal serrations.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 185202    | Straight, 9-1/8" (23 cm) |
| 185204    | Curved, 9" (23 cm)       |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps

#### DEBAKEY® THORACIC DISSECTING FORCEPS



DeBakey®-type jaws. The jaws are narrowed and thinned to a minimum size and double-curved. This provides a lightweight, sure-holding forceps for dissecting and passing tapes and ligatures around vessels.

| CATALOG # | DESCRIPTION                                                            |
|-----------|------------------------------------------------------------------------|
| 342560    | C curved, 3/8" (1 cm) tip length, 8-1/8" (20.5 cm)                     |
| 342561    | D curved, 1/2" (1.5 cm) tip length, 2.8 mm tip width, 8-1/4" (21 cm)   |
| 342562    | C curved, 5/8" (1.5 cm) tip length, 2.8 mm tip width, 9-5/8" (24.5 cm) |
| 342563    | D curved, 7/8" (2 cm) tip length, 9-1/2" (24 cm)                       |
| 342565    | D curved, 3/4" (2 cm) tip length, 10-1/4" (26 cm)                      |



342560



342561



342562



342563



342565



## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps

#### BLAKE GALLSTONE GRASPING FORCEPS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 172616    | Curved, serrated, 8" (20.5 cm) |

#### CYSTIC DUCT CATHETER CLAMPS




For open cholangiography.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 172570    | Accepts 4 Fr. catheter, 7-1/2" (19 cm) |
| 172571    | Accepts 6 Fr. catheter, 7-1/2" (19 cm) |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps

#### DEBAKEY® BAINBRIDGE CLAMPS




With narrow DeBakey®-type jaws.

| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 353560    | Straight, 40 x 2.4 mm jaws, 7" (18 cm)     |
| 353562    | Straight, 80 x 2.8 mm jaws, 7-7/8" (20 cm) |
| 353565    | Curved, 60 x 2.4 mm jaws, 7" (18 cm)       |



#### BEST COLON CLAMP




Right angle design for open anastomosis. Strong crushing jaws.

| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 185331    | 2-5/8" (6.5 cm) jaw size, 10-5/8" (27 cm) |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps


#### BABCOCK TISSUE FORCEPS



Tungsten carbide jaw tips.

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 500410    | 3/8" (9.5 mm) tip, 6-1/2" (16.5 cm), TC |
| 500412    | 7/16" (11 mm) tip, 7-7/8" (20 cm), TC   |

#### BAINBRIDGE FORCEPS



Pilling® Pattern, satin finish. Slender blades longitudinally serrated.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 171117    | Straight, 7-1/4" (18.5 cm) |
| 171120    | Curved, 6-3/4" (17 cm)     |
| 171122    | Curved, 7-1/4" (18.5 cm)   |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps

#### WILLAUER-ALLIS TISSUE FORCEPS



Extra-long, springy, Allis-type forceps with 5 x 6 teeth. Dull, non-glare finish stainless steel.

| CATALOG # | DESCRIPTION       |
|-----------|-------------------|
| 342720    | 10" (25.5 cm)     |
| 342721    | 11-1/4" (28.5 cm) |

#### BABCOCK TISSUE FORCEPS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 183010    | 6-1/4" (16 cm)           |
| 183020    | 7-7/8" (20 cm)           |
| 342725    | Serrated, 9-1/2" (24 cm) |
| 342726    | Serrated, 9-1/2" (24 cm) |

#### BABCOCK TISSUE FORCEPS




Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 500399    | 5.8 mm tip, 5-3/4" (14.5 cm) |
| 500400    | 8.3 mm tip, 6-3/8" (16 cm)   |
| 500402    | 10 mm tip, 7-7/8" (20 cm)    |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps

#### JUDD-ALLIS INTESTINAL FORCEPS



| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 182900    | 3 x 4 teeth, 6" (15 cm)       |
| 182905    | 3 x 4 teeth, 7-1/4" (18.5 cm) |

#### LOCKWOOD-ALLIS INTESTINAL TISSUE FORCEPS



| CATALOG # | DESCRIPTION                                              |
|-----------|----------------------------------------------------------|
| 182920    | Light model, flattened jaws, 5 x 6 teeth, 6-1/4" (16 cm) |

#### THOMS-ALLIS TISSUE AND TRACTION FORCEPS



| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 182910    | 6 x 7 teeth, 8" (20.5 cm) |

## FORCEPS AND CLAMPS

### Gastrointestinal Forceps and Clamps

#### ALLIS DELICATE TISSUE FORCEPS



Weck® Pattern, bright finish, 6-1/8" (15.5 cm).

| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 541100    | 2 x 3 teeth , 3.8 mm tip                   |
| 541120    | Straight, delicate, 3 x 4 teeth , 4 mm tip |

#### ALLIS TISSUE FORCEPS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 182795    | 3 x 4 teeth, 6" (15 cm)       |
| 182850    | 4 x 5 teeth, 6" (15 cm)       |
| 182860    | 5 x 6 teeth, 6" (15 cm)       |
| 182870    | 5 x 6 teeth, 7-1/2" (19 cm)   |
| 182875    | 5 x 6 teeth, 9" (23 cm)       |
| 182880    | 4 x 5 teeth, 5-3/4" (14.5 cm) |

#### ALLIS TISSUE FORCEPS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 541094    | 3 x 4 teeth, 5" (12.5 cm)     |
| 500310    | 3 x 4 teeth, 6-1/8" (15.5 cm) |
| 500311    | 4 x 5 teeth, 6-1/8" (15.5 cm) |
| 500312    | 5 x 6 teeth, 6-1/8" (15.5 cm) |
| 500316    | 5 x 6 teeth, 7-1/2" (19 cm)   |



## FORCEPS AND CLAMPS

### Right Angle Clamps

#### WILLAUER FORCEPS



Lightweight full curved, longitudinally serrated jaws.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 342580    | 2" (5 cm) tip length, 9-1/2" (24 cm) |


### Gastrointestinal Forceps and Clamps

#### ADAIR-ALLIS INTESTINAL TISSUE FORCEPS



| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 182915    | 9 x 10 teeth, 6-1/4" (16 cm) |

#### ALLEN INTESTINAL CLAMPS



Longitudinal serrations, 1 x 2 teeth.

| CATALOG # | DESCRIPTION                                                    |
|-----------|----------------------------------------------------------------|
| 185005    | Jaws 3-3/4" (9.5 cm), 9-1/2" (24 cm)                           |
| 185006    | Jaws 2-3/8" (6 cm), 8" (20.5 cm)                               |
| 185010    | Pediatric size, delicate, jaws 2-1/8" (5.5 cm), 6-1/4" (16 cm) |

## FORCEPS AND CLAMPS

### Right Angle Clamps

#### STOREY HEMOSTATIC FORCEPS



33 mm jaw with cross serrations.

| CATALOG # | DESCRIPTION                     |
|-----------|---------------------------------|
| 496215    | Slightly curved, 8-3/4" (22 cm) |
| 496216    | Strongly curved, 8-5/8" (22 cm) |
| 496217    | Right angle, 8-3/8" (21.5 cm)   |

#### WESTPHAL DISSECTING FORCEPS



Longitudinal serrations, narrow tip.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 496112    | Angled, 7-7/8" (20 cm) |

## FORCEPS AND CLAMPS

### Right Angle Clamps


#### RUMEL THORACIC FORCEPS



Furnished with DeBakey®-type jaw or standard serrated jaws. In addition to being one of the more popular thoracic artery forceps, they are especially suited for passing tapes around the cava or pulmonary artery in open heart surgery.

| CATALOG # | DESCRIPTION                                                                        |
|-----------|------------------------------------------------------------------------------------|
| 342501    | A curved, DeBakey® jaw, 2-3/4" (7 cm) tip length, 9-5/8" (24.5 cm)                 |
| 342504    | D curved, DeBakey® jaw, 1-1/8" (3 cm) tip length, 2.8 mm tip width, 9-3/4" (25 cm) |
| 342511    | A curved, serrated, 1/2" (1.5 cm) tip length, 2.4 mm tip width, 9-3/8" (24 cm)     |
| 342513    | C curved, serrated, 3/4" (2 cm) tip length, 2 mm tip width, 9-5/8" (24.5 cm)       |
| 342514    | D curved, serrated, 1" (2.5 cm) tip length, 2.4 mm tip width, 9" (23 cm)           |
| 342533    | Child size, fully curved, standard serrated, 2 mm tip width, 7-7/8" (20 cm)        |

#### SAWTELL HEMOSTATIC FORCEPS



| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 342583    | Full curved, heavy tips, 7-1/4" (18.5 cm) |

## FORCEPS AND CLAMPS

### Right Angle Clamps

#### MIXTER "PEDI" GALL DUCT FORCEPS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 172495    | Delicate, 5-1/8" (13 cm) |



#### RIENHOFF SWAN-NECK CLAMP



A double-curved forceps with delicately rounded tips.

| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 342625    | Serrated, 10-5/8" (27 cm) |

#### ROOS SCALENECTOMY FORCEPS



Designed to give a secure hold on the scalene muscle during scaleneotomy.

| CATALOG # | DESCRIPTION                                                                                   |
|-----------|-----------------------------------------------------------------------------------------------|
| 342575    | Right angle, longitudinal serrations, 1/2" (1.5 cm) tip length, 1.5 mm tip width, 11" (28 cm) |

## FORCEPS AND CLAMPS

### Right Angle Clamps

#### MIXTER FORCEPS



| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| KM36566   | Fully curved, serrated, 9" (23 cm) |

#### MIXTER GALL DUCT FORCEPS



Acutely curved jaws. Full serrations run entire length of jaw.

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 172500    | 6-1/4" (16 cm) |
| 172503    | 7-1/8" (18 cm) |
| 172506    | 8-3/4" (22 cm) |

#### PEDIATRIC MIXTER DISSECTING FORCEPS




Weck® Pattern, bright finish, fully curved.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 477090    | 24 mm jaw, 2 mm tip, 5-1/4" (13.5 cm) |

## FORCEPS AND CLAMPS

### Right Angle Clamps


#### MIXTER THORACIC FORCEPS



Regular model. Available with DeBakey®-type jaw or standard serrated jaw, 3/4" (2 cm) long.


| CATALOG # | DESCRIPTION                                                                           |
|-----------|---------------------------------------------------------------------------------------|
| 342603    | Half curved, serrated, 5/8" (1.5 cm) tip length, 2.3 mm tip width, 9-1/8" (23 cm)     |
| 342600    | Curved, DeBakey®-type jaw, 1" (2.5 cm) tip length, 2.8 mm tip width, 9-1/4" (23.5 cm) |
| 342601    | Strongly curved, serrated, 3/4" (2 cm) tip length, 8-3/4" (22 cm)                     |

#### MIXTER CLAMPS



| CATALOG # | DESCRIPTION                                                             |
|-----------|-------------------------------------------------------------------------|
| 477108    | 40 degree-angled, fully serrated, 50 mm jaw, 3 mm tip, 7-3/4" (19.5 cm) |

#### MIXTER DISSECTING FORCEPS



Two-thirds serrations.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 477091    | 26 mm jaw, 2 mm tip, 7-1/4" (18.5 cm) |
| 477092    | 40 mm jaw, 2.2 mm tip, 9" (23 cm)     |
| 477093    | 55 mm jaw, 2.4 mm tip, 11" (28 cm)    |



## FORCEPS AND CLAMPS

### Right Angle Clamps

#### LOWER GALL DUCT FORCEPS



| CATALOG # | DESCRIPTION                                                       |
|-----------|-------------------------------------------------------------------|
| 172525    | Short jaws, fully curved, longitudinal serrations, 7-1/2" (19 cm) |

#### MEEKER HEMOSTATIC DISSECTING FORCEPS



| CATALOG # | DESCRIPTION                                                        |
|-----------|--------------------------------------------------------------------|
| 540150    | 90 degree-angled, serrated, 23 mm jaw, 1.8 mm tip, 6-3/8" (16 cm)  |
| 504190    | 90 degree-angled, serrated, 28 mm jaw, 1.8 mm tip, 10-5/8" (27 cm) |

#### MACQUIGG-MIXTER FORCEPS



An excellent general purpose forceps and bronchus clamp. Longitudinally serrated jaws, right angle tip.

| CATALOG # | DESCRIPTION                                                               |
|-----------|---------------------------------------------------------------------------|
| 342611    | Curved, longitudinal serrations, 1-5/8" (4 cm) tip length, 8-5/8" (22 cm) |

## FORCEPS AND CLAMPS

### Right Angle Clamps

#### LAHEY DISSECTING FORCEPS



| CATALOG # | DESCRIPTION                                                      |
|-----------|------------------------------------------------------------------|
| 496105    | Longitudinal serrations, 2.4 mm tip, 7-7/8" (20 cm)              |
| 496110    | Horizontal serrations, 44 mm jaw, 2.4 mm tip, 7-3/4" (19.5 cm)   |
| 496111    | Longitudinal serrations, 46 mm jaw, 2.5 mm tip, 7-3/4" (19.5 cm) |

#### LAHEY GALL DUCT FORCEPS, SLENDER BLADES



| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 172530    | Longitudinal serrations, 7-1/2" (19 cm)   |
| 172531    | Standard serrations, 7-1/2" (19 cm)       |
| 172532    | Longitudinal serrations, 8-7/8" (22.5 cm) |

## FORCEPS AND CLAMPS


### Right Angle Clamps

#### JACKSON TRACHEAL HEMOSTAT



| CATALOG # | DESCRIPTION                      |
|-----------|----------------------------------|
| 514480    | Acutely curved, 5-3/8" (13.5 cm) |

#### KANTROWITZ FORCEPS



True right angle forceps with cross serrated jaws. Useful for grasping and dissecting in thoracic, vascular and general surgery.

| CATALOG # | DESCRIPTION                                                       |
|-----------|-------------------------------------------------------------------|
| 342557    | Serrated, 3/4" (2 cm) tip length, 2 mm tip width, 7-7/8" (20 cm)  |
| 342558    | Serrated, 9-1/2" (24 cm)                                          |
| 342559    | Serrated, 3/4" (2 cm) tip length, 2 mm tip width, 10-5/8" (27 cm) |

PR

## FORCEPS AND CLAMPS

### Right Angle Clamps

#### GEMINI FORCEPS



A more delicate modification of the Mixer Forceps. Slender, acutely curved jaws with full serrations.


| CATALOG # | DESCRIPTION       |
|-----------|-------------------|
| 342596    | 7-1/4" (18.5 cm)  |
| 342597    | 9" (23 cm)        |
| 342598    | 11-1/8" (28.5 cm) |

#### GRAY CYSTIC DUCT FORCEPS



| CATALOG # | DESCRIPTION                                                                |
|-----------|----------------------------------------------------------------------------|
| 172550    | Curved jaws and shafts.<br>Set of two: 8-1/4" (21 cm) and 9-1/4" (23.5 cm) |

#### HARRINGTON CLAMPS




Right angle jaw. Dull, non-glare finish stainless steel.

| CATALOG # | DESCRIPTION                                                                              |
|-----------|------------------------------------------------------------------------------------------|
| 342632    | Heavy construction, longitudinal serrations, 1-1/2" (4 cm) tip length, 12-3/8" (31.5 cm) |
| 342633    | Light construction, longitudinal serrations, 1-1/2" (4 cm) tip length, 12-3/8" (31.5 cm) |

## FORCEPS AND CLAMPS

### Right Angle Clamps

#### DIETRICH RIGHT ANGLE CLAMPS



| CATALOG # | DESCRIPTION                                                   |
|-----------|---------------------------------------------------------------|
| 355798    | Serrated jaw, 5-3/8" (13.5 cm)                                |
| 355799    | Serrated jaw, 4-7/8" (12.5 cm)                                |
| 355800    | Serrated jaw 8 mm past angle, 7-1/4" (18.5 cm)                |
| 355801    | Longitudinally serrated jaw 8 mm past angle, 7-1/4" (18.5 cm) |

#### FINOCHIETTO THORACIC FORCEPS




A utility clamp for use as ligature carrier, hemostat or dissector. Particularly good in ligation during esophageal surgery. Blades have longitudinal serrations and holes for threading ligature.

| CATALOG # | DESCRIPTION                                                               |
|-----------|---------------------------------------------------------------------------|
| 342615    | Angled, longitudinal serrations, 1-1/8" (3 cm) tip length, 9-1/2" (24 cm) |

## FORCEPS AND CLAMPS

### Right Angle Clamps

#### CROUCH CLAMP



| CATALOG # | DESCRIPTION     |
|-----------|-----------------|
| 355055    | 10-1/4" (26 cm) |

#### DEBAKEY® RIGHT ANGLE FORCEPS



With full length longitudinal serrations and cross serrations at tip.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 342566    | 7-5/8" (19.5 cm) |
| 342567    | 9" (23 cm)       |
| 342568    | 10-1/4" (26 cm)  |
| 342569    | 11-3/8" (29 cm)  |



## FORCEPS AND CLAMPS

Extra Delicate Hemostatic Forceps

### PLASTI STATS® CRILE HEMOSTATS



Full serrations.

| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 181800    | Straight, 1-1/4" (32 mm) jaw, 5-3/4" (14.5 cm) |
| 181802    | Curved, 1-1/4" (32 mm) jaw, 5-1/2" (14 cm)     |

### PLASTI STATS® KOCHER FORCEPS



| CATALOG # | DESCRIPTION                                            |
|-----------|--------------------------------------------------------|
| 181832    | Straight, 1 x 2 teeth, full serrations, 5-1/2" (14 cm) |

### PLASTI STATS® MOSQUITO FORCEPS



Full serrations.

| CATALOG # | DESCRIPTION                                  |
|-----------|----------------------------------------------|
| 181810    | Straight, 3/4" (2 cm) jaw, 4" (10 cm)        |
| 181812    | Curved, 5/8" (1.5 cm) jaw, 4" (10 cm)        |
| 181818    | Straight, 1/2" (12.5 mm) jaw, 4-3/4" (12 cm) |
| 181820    | Curved, 1/2" (12.5 mm) jaw, 4-5/8" (11.5 cm) |
| 181814    | Straight, 3/4" (2 cm) jaw, 5" (12.5 cm)      |
| 181816    | Curved, 3/4" (2 cm) jaw, 5" (12.5 cm)        |
| 181806    | Straight, 1" (2.5 cm) jaw, 5-1/2" (14 cm)    |
| 181808    | Curved, 1" (2.5 cm) jaw, 5-1/2" (14 cm)      |

## FORCEPS AND CLAMPS

Extra Delicate Hemostatic Forceps

### MICRO MOSQUITO FORCEPS




KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                         |
|-----------|-----------------------------------------------------|
| KM36164   | Straight, 1 mm tip, 5" (12.5 cm)                    |
| KM36165   | Curved, serrated, extremely fine tips, 5" (12.5 cm) |



### PLASTI STATS® ADSON FORCEPS



| CATALOG # | DESCRIPTION                                          |
|-----------|------------------------------------------------------|
| 181824    | Curved, 1/2" (12.5 mm) jaw serration, 5-1/2" (14 cm) |

### PLASTI STATS® COLLIER FORCEPS



| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 181804    | Curved, 7-1/2" (19 cm) |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### VANDERBILT VESSEL FORCEPS



This hemostatic-type forceps is particularly useful in cardiovascular and thoracic surgery but has been used extensively in general surgery as well. The elongated handles afford access to a number of surgical sites. Jaws are of a standard cross-serration design.

| CATALOG # | DESCRIPTION                                                           |
|-----------|-----------------------------------------------------------------------|
| 342590    | Slightly curved, 2" (5 cm) tip length, 2 mm tip width, 8-5/8" (22 cm) |

#### SERREFINES CROSS-ACTION CLAMPS



Cross-action, full serrations.

| CATALOG # | DESCRIPTION               |
|-----------|---------------------------|
| 182200    | Straight, 1-1/2" (4 cm)   |
| 182205    | Straight, 2-1/4" (5.5 cm) |
| 182210    | Curved, 1-3/8" (3.5 cm)   |

### Extra-Delicate Hemostatic Forceps

#### JACOBSON HEMOSTATIC FORCEPS



Pilling® Pattern, satin finish. These delicate forceps feature tips as fine as 0.8 mm and can be used in a number of microsurgical procedures.

| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 443500    | Straight, 0.8 mm tip, 3-7/8" (10 cm)      |
| 443502    | Curved, 0.8 mm tip, 3-7/8" (10 cm)        |
| 443505    | Straight, 0.8 mm tip, 5-1/8" (13 cm)      |
| 443507    | Curved, 0.8 mm tip, 5-1/8" (13 cm)        |
| 443508    | Full curved, 0.8 mm tip, 4-7/8" (12.5 cm) |
| 443510    | Curved, 1.2 mm tip, 7-1/8" (18 cm)        |
| 443512    | Full curved, 1 mm tip, 7-1/8" (18 cm)     |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### SAROT INTRA-THORACIC ARTERY FORCEPS



Slender, slightly curved forceps with serrated jaws. Specially designed for pulmonary resections, cardiovascular and esophageal surgery.

| CATALOG # | DESCRIPTION                                                                    |
|-----------|--------------------------------------------------------------------------------|
| 342571    | Medium, serrated, 1/2" (1.5 cm) tip length, 1.6 mm tip width, 7-5/8" (19.5 cm) |
| 342572    | Large, serrated, 2-1/4" (5.5 cm) tip length, 2 mm tip width, 9" (23 cm)        |

#### SCHNIDT FORCEPS (TONSIL)



Pilling® Pattern, satin finish. Half serrations.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 182390    | Half curved, 7-1/2" (19 cm) |

#### SCHNIDT FORCEPS (TONSIL)



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                         |
|-----------|-----------------------------------------------------|
| 075116    | Slightly curved, 7-1/4" (18.5 cm)                   |
| 075120    | Lightly curved, serrated, 2 mm tip, 7-1/2" (19 cm)  |
| 075122    | Strongly curved, serrated, 2 mm tip, 7-1/2" (19 cm) |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### PROVIDENCE HOSPITAL HEMOSTATS



Pilling® Pattern, satin finish. Half serrations.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 182140    | Straight, 5-3/4" (14.5 cm) |
| 182150    | Curved, 5-5/8" (14.5 cm)   |

#### RANKIN FORCEPS



Pilling® Pattern, satin finish. Half serrations.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 182120    | Straight, 6-1/4" (16 cm) |
| 182130    | Curved, 6-1/4" (16 cm)   |

#### RANKIN-KELLY HEMOSTAT FORCEPS



KMedic® Pattern, satin finish. Three-quarter serrations.


| CATALOG # | DESCRIPTION                      |
|-----------|----------------------------------|
| KM36196   | Straight, 6-1/4" (16 cm)         |
| KM36198   | Curved, serrated, 6-1/4" (16 cm) |



## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### ROCHESTER-OCHSNER FORCEPS (KOCHER)



182620

Pilling® Pattern, straight. Full serrations, 1 x 2 teeth.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 182600    | 5-1/2" (14 cm)   |
| 182605    | 6-1/4" (16 cm)   |
| 182615    | 7-1/4" (18.5 cm) |
| 182620    | 8" (20.5 cm)     |

#### ROCHESTER-OCHSNER FORCEPS (KOCHER)



182670

Pilling® Pattern, curved. Full serrations, 1 x 2 teeth.

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 182655    | 6-1/4" (16 cm)                    |
| 182665    | 7-1/4" (18.5 cm)                  |
| 182670    | Curved, 1 x 2 teeth, 8" (20.5 cm) |

#### PROVIDENCE HOSPITAL FORCEPS



475200

Weck® Pattern, bright finish. Two-thirds serrations.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 475200    | Straight, 2 mm tip, 5-3/8" (13.5 cm) |
| 475201    | Curved, 2 mm tip, 5-3/8" (13.5 cm)   |



## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps


#### ROCHESTER-PEAN FORCEPS, 12"



KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| KM36360   | Hemostatic, straight, 4.6 mm tip, 12" (30.5 cm) |
| KM36362   | Curved, serrated, 12" (30.5 cm)                 |

#### OCHSNER FORCEPS (KOCHER)



Weck® Pattern, bright finish. Straight hemostatic forceps with full serrations, 1 x 2 teeth.

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 476200    | 3 mm tip, 5-5/8" (14.5 cm)              |
| 476205    | 4.2 mm tip, 6-3/8" (16 cm)              |
| 476206    | Light model, 2.5 mm tip, 6-3/8" (16 cm) |
| 476210    | 3.3 mm tip, 7-1/4" (18.5 cm)            |
| 476215    | 4.5 mm tip, 8" (20.5 cm)                |
| 476220    | 5.5 mm tip, 8-7/8" (22.5 cm)            |
| 476225    | 5 mm tip, 10-3/8" (26.5 cm)             |

#### OCHSNER FORCEPS (KOCHER)




Weck® Pattern, bright finish. Curved hemostatic forceps with full serrations, 1 x 2 teeth.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 476250    | 2.7 mm tip, 5-1/2" (14 cm)   |
| 476255    | 2.6 mm tip, 6-1/4" (16 cm)   |
| 476260    | 3.6 mm tip, 7-1/4" (18.5 cm) |
| 476265    | 3.4 mm tip, 7-7/8" (20 cm)   |
| 476270    | 4 mm tip, 8-3/4" (22 cm)     |
| 476275    | 4 mm tip, 10-1/4" (26 cm)    |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### ROCHESTER-PEAN FORCEPS



Pilling® Pattern, satin finish. Straight, full serrations.

| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| 181830    | Curved, 2 mm tip, 6-1/2" (16.5 cm) |



#### ROCHESTER-PEAN FORCEPS



Pilling® Pattern, satin finish. Straight, full serrations.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 182400    | 2.7 mm tip, 5-1/2" (14 cm)   |
| 182405    | 2.8 mm tip, 6-1/2" (16.5 cm) |
| 182415    | 3.0 mm tip, 7-1/4" (18.5 cm) |
| 182420    | 3.2 mm tip, 8" (20.5 cm)     |

#### ROCHESTER-PEAN FORCEPS




Pilling® Pattern, satin finish. Curved, full serrations.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 182440    | 3.2 mm tip, 5-1/2" (14 cm)   |
| 182445    | 2.8 mm tip, 6-1/4" (16 cm)   |
| 182455    | 3.0 mm tip, 7-1/4" (18.5 cm) |
| 182460    | 3.2 mm tip, 8" (20.5 cm)     |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### MAYO-PEAN FORCEPS



Weck® Pattern, bright finish. Straight hemostatic forceps with full serrations.

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 476105    | 2.9 mm tip, 6-5/8" (17 cm)    |
| 476106    | 2.4 mm tip, 6-5/8" (17 cm)    |
| 476110    | 3.3 mm tip, 7-1/4" (18.5 cm)  |
| 476115    | 3.7 mm tip, 8-1/4" (21 cm)    |
| 476120    | 3.7 mm tip, 8-3/4" (22 cm)    |
| 476125    | 3.8 mm tip, 10-3/8" (26.5 cm) |

#### MAYO-PEAN FORCEPS



Weck® Pattern, bright finish. Curved hemostatic forceps with full serrations.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 476150    | 2.5 mm tip, 5-5/8" (14.5 cm) |
| 476155    | 2.9 mm tip, 6-1/2" (16.5 cm) |
| 476156    | 2.4 mm tip, 6-1/2" (16.5 cm) |
| 476160    | 3.3 mm tip, 7-1/8" (18 cm)   |
| 476165    | 3.7 mm tip, 8-1/4" (21 cm)   |
| 476170    | 3.5 mm tip, 9-1/4" (23.5 cm) |
| 476175    | 3.8 mm tip, 10-1/4" (26 cm)  |

#### PEAN CLAMPS




Pilling® Pattern, satin finish. Horizontal serrations.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 185060    | Straight, 8-1/8" (20.5 cm)  |
| 185065    | Curved, 8-1/8" (20.5 cm)    |
| 185066    | Curved, 9-1/4" (23.5 cm)    |
| 185067    | Straight, 10-1/8" (25.5 cm) |


## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps


#### LAHEY TRACTION FORCEPS



#### LAWRENCE HEMOSTATIC FORCEPS



#### MASTIN MUSCLE FORCEPS



## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### JONES HEMOSTAT



Sturdy construction, fully serrated.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 182160    | Straight, 5" (12.5 cm) |

#### KELLY FORCEPS



Pilling® Pattern, satin finish. Half serrations.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 182100    | Straight, 5-3/4" (14.5 cm) |
| 182110    | Curved, 5-1/2" (14 cm)     |

#### KELLY FORCEPS



Weck® Pattern, bright finish. Half serrations.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 475210    | Straight, 5-1/2" (14 cm) |
| 475211    | Curved, 5-1/2" (14 cm)   |

## FORCEPS AND CLAMPS

Hemostatic Forceps and Clamps

### MOSQUITO FORCEPS SERRATED



| CATALOG # | DESCRIPTION                                 |
|-----------|---------------------------------------------|
| KM36078   | Curved, serrated, 1 x 2 teeth, 5" (12.5 cm) |

### SKILLMAN MOSQUITO FORCEPS



Weck® Pattern, bright finish. Full serrations.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 475130    | Straight, 1.6 mm tip, 5-1/4" (13.5 cm) |
| 475131    | Curved, 1.6 mm tip, 5-1/8" (13 cm)     |



## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### HARTMANN MOSQUITO FORCEPS



Weck® Pattern, bright finish. Delicate, serrated jaw.

| CATALOG # | DESCRIPTION                                 |
|-----------|---------------------------------------------|
| 475080    | Straight, 20 mm jaw, 1.3 mm tip, 4" (10 cm) |
| 475082    | Curved, 20 mm jaw, 1.3 mm tip, 4" (10 cm)   |

#### MOSQUITO FORCEPS



| CATALOG # | DESCRIPTION                                              |
|-----------|----------------------------------------------------------|
| 075118    | Curved, fully serrated, 1.7 x 2 mm tip, 7-1/4" (18.5 cm) |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### HALSTED MOSQUITO FORCEPS



Pilling® Pattern, satin finish, full serrations.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 182300    | Straight, 5-1/4" (13.5 cm) |
| 182310    | Curved, 5-1/8" (13 cm)     |
| 182350    | Straight, 4" (10 cm)       |
| 182355    | Curved, 4" (10 cm)         |

#### HALSTED MOSQUITO FORCEPS



Weck® Pattern, bright finish, serrated jaw.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 475084    | Straight, 1.2 mm tip, 4-3/4" (12 cm)   |
| 475086    | Curved, 1.2 mm tip, 4-5/8" (11.5 cm)   |
| 475100    | Straight, 1.6 mm tip, 5-1/4" (13.5 cm) |
| 475101    | Curved, 1.6 mm tip, 5-1/8" (13 cm)     |

#### HALSTED MOSQUITO FORCEPS WITH TEETH



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| 475102    | Straight, 1 x 2 teeth, 1.5 mm tip, 5" (12.5 cm) |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### CRILE FORCEPS




475311

Weck® Pattern, bright finish, delicate, full serrations.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 475310    | Straight, 1.8 mm tip, 5-1/2" (14 cm) |
| 475311    | Curved, 1.8 mm tip, 5-3/8" (13.5 cm) |
| 475312    | Straight, 1.7 mm tip, 6-3/8" (16 cm) |
| 475313    | Curved, 2 mm tip, 6-1/4" (16 cm)     |

#### DANDY HEMOSTATIC FORCEPS



181822

Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                     |
|-----------|---------------------------------|
| 181822    | Curved sideways, 5-1/2" (14 cm) |

#### DANDY KOLODNY FORCEPS



477170

Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                                       |
|-----------|-------------------------------------------------------------------|
| 477170    | Curved sideways, one-half horizontal serrations, 5-3/8" (13.5 cm) |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### COLLIER HEMOSTATIC FORCEPS



A Crile-type forceps with a lighter distal end for more delicate procedures.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 182085    | Straight, 5-1/2" (14 cm)   |
| 182086    | Curved, 5-1/2" (14 cm)     |
| 182090    | Straight, 6-1/4" (16 cm)   |
| 182095    | Curved, 6-1/4" (16 cm)     |
| 182097    | Straight, 7-3/4" (19.5 cm) |
| 182098    | Curved, 7-1/4" (18.5 cm)   |

#### CRILE HEMOSTATS



Pilling® Pattern, satin finish, standard tip, full serrations.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 182050    | Straight, 5-3/4" (14.5 cm) |
| 182060    | Curved, 5-1/2" (14 cm)     |
| 182070    | Straight, 6-1/4" (16 cm)   |
| 182080    | Curved, 6-1/4" (16 cm)     |

#### CRILE FORCEPS




Weck® Pattern, bright finish, full serrations.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 475308    | Straight, 5-5/8" (14.5 cm) |
| 475309    | Curved, 5-3/4" (14.5 cm)   |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### COLLER FORCEPS




| CATALOG # | DESCRIPTION                                                             |
|-----------|-------------------------------------------------------------------------|
| 475305    | Skillman serrations, straight, 31 x 1.7 mm tip, 5-1/2" (14 cm)          |
| 475306    | Skillman serrations, curved, 29 x 1.7 mm tip, 5-3/8" (13.5 cm)          |
| 475307    | Skillman serrations, strongly curved, 26 x 1.7 mm tip, 5-3/8" (13.5 cm) |
| 477109    | Straight, 38 x 1.8 mm tip, 1 x 3 catch, 6-5/8" (17 cm)                  |
| 477110    | Curved, 38 x 1.8 mm tip, 1 x 3 catch, 6-5/8" (17 cm)                    |
| 477112    | Curved, 42 x 1.7 mm tip, 1 x 3 catch, 7-1/4" (18.5 cm)                  |
| 477113    | Strongly curved, 37 x 1.6 mm tip, 1 x 3 catch, 7" (18 cm)               |
| 477116    | Curved, 44 x 2 mm tip, 4 x 4 catch, 9-3/8" (24 cm)                      |
| 477117    | Strongly curved, 40 x 2 mm tip, 4 x 4 catch, 9-1/4" (23.5 cm)           |
| 477118    | Angled, 32 x 2 mm tip, 4 x 4 catch, 8-7/8" (22.5 cm)                    |



## FORCEPS AND CLAMPS

Hemostatic Forceps and Clamps

### CARMALT FORCEPS



Cross-serrated tip with longitudinally serrated jaw.


| CATALOG # | DESCRIPTION                                     |
|-----------|-------------------------------------------------|
| 476300    | Straight, 43 mm jaw, 3 mm tip, 6-1/2" (16.5 cm) |
| 476310    | Straight, 58 mm jaw, 4 mm tip, 8-1/4" (21 cm)   |
| 476330    | Curved, 43 mm jaw, 3 mm tip, 6-1/2" (16.5 cm)   |
| 476340    | Curved, 58 mm jaw, 4 mm tip, 8-1/4" (21 cm)     |



## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps

#### LEJEUNE THORACIC FORCEPS



| CATALOG # | DESCRIPTION                                                                             |
|-----------|-----------------------------------------------------------------------------------------|
| 342594    | Slightly curved, serrated, 1-3/8" (3.5 cm) tip length, 1.8 mm tip width, 9-3/8" (24 cm) |


#### BENGOLEA FORCEPS



Lightweight, with narrow, criss-cross serrated jaws.

| CATALOG # | DESCRIPTION                                                |
|-----------|------------------------------------------------------------|
| 342587    | Slightly curved, serrated, 2 mm wide tip, 8-1/2" (21.5 cm) |

#### BRIDGE HEMOSTATIC FORCEPS



32 mm serrated jaw.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 504180    | Straight, 2 mm tip, 11-1/8" (28.5 cm) |
| 504181    | Curved, 2 mm tip, 11" (28 cm)         |
| 540111    | Curved, 1.6 mm tip, 6-7/8" (17.5 cm)  |

## FORCEPS AND CLAMPS

### Hemostatic Forceps and Clamps


#### ADSON HEMOSTATIC FORCEPS



30 mm jaw, 2 mm tip, two-thirds serrations.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 477190    | Straight, 7-1/2" (19 cm) |
| 477191    | Curved, 7-1/8" (18 cm)   |

#### PILLING® THORACIC ARTERY FORCEPS




Slim, slightly curved jaws have very light criss-cross serrations making this an excellent instrument for handling suture material as well as delicate tissue.

| CATALOG # | DESCRIPTION                                                      |
|-----------|------------------------------------------------------------------|
| 342585    | Serrated, 2" (5 cm) tip length, 2 mm tip width, 9-1/4" (23.5 cm) |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### RUSSIAN TISSUE FORCEPS



| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 181360    | 6" (15 cm)    |
| 181365    | 8" (20.5 cm)  |
| 181370    | 10" (25.5 cm) |

#### C.W. MAYO TISSUE FORCEPS



Russian-type fenestrated tips are 6.7 mm wide and deeply serrated.

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 342910    | 9" (23 cm)  |

#### DAVIS TISSUE FORCEPS



Russian-type jaws, 7.5 mm wide.

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 342915    | Straight, serrated, 10" (25.5 cm) |

#### SELMAN TISSUE FORCEPS



Small diameter circular tips have meshing teeth for secure, non-damaging grip.

| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 351840    | 6" (15 cm)    |
| 351841    | 8" (20.5 cm)  |
| 351842    | 10" (25.5 cm) |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### CORONARY ULTRA-FINE TISSUE FORCEPS



Tungsten carbide.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 351837    | 7" (18 cm), TC   |
| 351838    | 8" (20.5 cm), TC |
| 351839    | 9" (23 cm), TC   |

#### WHEELER TECHNIQUE® PLAQUE AND VESSEL FORCEPS



Very powerful grip at tip.

| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 351825    | DeBakey®-type jaw, 1.6 mm, 5-1/2" (14 cm) |
| 351827    | DeBakey®-type jaw, 1.6 mm, 7-1/8" (18 cm) |
| 351829    | DeBakey®-type jaw, 1.6 mm, 8-3/4" (22 cm) |

#### CRAFOORD TISSUE FORCEPS



Serrated tips are built up to provide an unusually secure grip.

| CATALOG # | DESCRIPTION            |
|-----------|------------------------|
| 351898    | Serrated, 8" (20.5 cm) |

#### POTTS-SMITH TISSUE FORCEPS



Fine, serrated, 1 x 2 distal teeth.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 351870    | 7-1/4" (18.5 cm) |
| 351872    | 8-1/4" (21 cm)   |
| 351874    | 9-1/2" (24 cm)   |
| 351876    | 10-1/4" (26 cm)  |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### NELSON TISSUE FORCEPS



| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 342950    | 6 x 7 sharp, fine teeth, 9" (23 cm) |

#### SINGLEY TISSUE FORCEPS



Serrated ring end 8.5 mm wide. Medium-heavy construction.

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 342933    | 9" (23 cm)  |

#### ADAMS TISSUE FORCEPS



| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 351835    | Delicate, 1 mm tip, 8" (20.5 cm), TC |

#### TUTTLE THORACIC TISSUE FORCEPS



Light construction. Ring tips 6.8 mm wide.

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 342935    | 9" (23 cm)  |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### IRIS TISSUE FORCEPS

Delicate.



181610



181614

| CATALOG # | DESCRIPTION                              |
|-----------|------------------------------------------|
| 181610    | Straight, 1 x 2 teeth, 4" (10 cm)        |
| 181614    | Full curve, 1 x 2 teeth, 3-3/4" (9.5 cm) |
| 181604    | Full curve, 3-3/4" (9.5 cm)              |

#### SEMKEN TISSUE FORCEPS



181272

These delicate forceps are midway in size between Iris forceps and delicate standard tissue forceps.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 181272    | 1 x 2 teeth, 5" (12.5 cm)   |
| KM39088   | 1 x 2 teeth, 6-1/4" (16 cm) |



## FORCEPS AND CLAMPS

### Tissue Forceps

#### GERALD DEBAKEY® TISSUE FORCEPS



351901



351910

Delicate, tapered DeBakey®-type serrations.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 351901    | 7" (18 cm)               |
| 351908    | 8" (20.5 cm)             |
| 351909    | 9-1/2" (24 cm)           |
| 351910    | Angled, 6-7/8" (17.5 cm) |

#### GRAEFE TISSUE FORCEPS




181250

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 181250    | 5 x 6 teeth, 4-3/8" (11 cm) |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### DEBAKEY® TITANIUM TISSUE FORCEPS



| CATALOG # | DESCRIPTION                                               |
|-----------|-----------------------------------------------------------|
| 359024    | Straight, flat handle, 1.5 mm jaw width, 6" (15 cm)       |
| 359025    | Straight, flat handle, 1.5 mm jaw width, 7-5/8" (19.5 cm) |
| 359021    | Straight, flat handle, 2.5 mm jaw width, 6" (15 cm)       |
| 359022    | Straight, flat handle, 2.5 mm jaw width, 7-5/8" (19.5 cm) |
| 359027    | Straight, flat handle, 3 mm jaw width, 6" (15 cm)         |
| 359028    | Straight, flat handle, 3 mm jaw width, 7-5/8" (19.5 cm)   |
| 359029    | Straight, flat handle, 3 mm jaw width, 9-1/2" (24 cm)     |
| 359030    | Angular, flat handle, 2.5 mm jaw width, 5-7/8" (15 cm)    |
| 359031    | Angular, flat handle, 2.5 mm jaw width, 7-1/2" (19 cm)    |
| 359032    | Angular, flat handle, 2.5 mm jaw width, 9-1/4" (23.5 cm)  |

#### FERRIS SMITH TISSUE FORCEPS



| CATALOG # | DESCRIPTION                                                     |
|-----------|-----------------------------------------------------------------|
| 468190    | 2 x 3 teeth with 4 mm serrated tying platform. 6-7/8" (17.5 cm) |

#### GERALD TISSUE FORCEPS



These are full-size tissue forceps with only the distal portion narrowed to provide light, delicate tips that maintain alignment.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 351905    | 1 x 2 teeth, 7" (18 cm)     |
| 351907    | 1 x 2 teeth, 9-3/4" (25 cm) |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### DEBAKEY® CARBIDE FORCEPS



351817

DeBaKey®-type serrations at distal tip. Tungsten carbide insert behind the jaw for secure needle pulling actions.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 351815    | Regular forceps, 6" (15 cm), TC        |
| 351817    | Regular forceps, 7-3/4" (19.5 cm), TC  |
| 351818    | Delicate forceps, 7-3/4" (19.5 cm), TC |
| 351819    | Regular forceps, 9-1/2" (24 cm), TC    |
| 351824    | Delicate forceps, 9-1/2" (24 cm), TC   |



351818

#### VASCULAR DEBAKEY® BAYONET TISSUE FORCEPS



351845

DeBaKey®-type serrations with a bayonet shape. Angled handles for difficult access and increased visibility.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 351845    | Fine, 1.3 mm tip, 7-3/4" (19.5 cm)     |
| 351843    | Standard, 2.0 mm tip, 7-3/4" (19.5 cm) |



351843

## FORCEPS AND CLAMPS

### Tissue Forceps


#### DEBAKEY® VASCULAR TISSUE FORCEPS



Standard. The distal 2 cm of these forceps is a narrow DeBakey®-type serration with a width of 2 mm. The body of the instrument is constructed especially to prevent feathering, and all pins and projections have been omitted to help eliminate suture catchers. Jaws are very well-fitted for handling tiny vessels as well as the smallest suture material.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 351800    | Straight, 2 mm tip, 6" (15 cm)       |
| 351802    | Straight, 2 mm tip, 7-3/4" (19.5 cm) |
| 351804    | Straight, 2 mm tip, 9-1/2" (24 cm)   |
| 351808    | Straight, 2 mm tip, 12" (30.5 cm)    |
| 351810    | Angled, 2 mm tip, 6" (15 cm)         |
| 351812    | Angled, 2 mm tip, 7-3/4" (19.5 cm)   |
| 351814    | Angled, 2 mm tip, 9-1/2" (24 cm)     |

#### DEBAKEY®-DIETHRICH DELICATE VASCULAR TISSUE FORCEPS



Designed for micro-vascular and coronary artery surgery. The distal portion is tapered to a 1 mm extra-fine DeBakey®-type serration.

| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| 351801    | Straight, delicate, 1 mm tip, 6" (15 cm)       |
| 351803    | Straight, delicate, 1 mm tip, 7-3/4" (19.5 cm) |
| 351806    | Straight, delicate, 1 mm tip, 9-1/2" (24 cm)   |

#### DEBAKEY® ULTRALIGHT FORCEPS



Approximately one-half the weight of standard DeBakey® Forceps, similar to that of titanium but at a much lower price. Ebonized stainless steel to reduce glare. Reduced tension minimizes fatigue and allows better control of delicate tissues.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 351820    | 1.5 mm tip, 6" (15 cm)       |
| 351821    | 1.5 mm tip, 7-3/4" (19.5 cm) |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### CUSHING BAYONET TISSUE FORCEPS



468240

Serrated handles, scraper end.

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 468240    | 1 x 2 teeth, 7-1/4" (18.5 cm) |

#### CUSHING GUTSCH TISSUE FORCEPS



468230

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 468230    | 1 x 2 teeth, serrated, 7" (18 cm) |

#### CUSHING TISSUE FORCEPS



468220

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 468220    | 1 x 2 teeth, serrated, 7" (18 cm) |

#### DEBAKEY® THORACIC TISSUE FORCEPS



342980



342986

The row of 1 x 2 DeBakey®-type serrations on the tips are designed to be large and strong enough to handle relatively large masses of tissue. Because of the groove construction and light spring tension, they may be used to hold and handle suture material as well. These forceps are of heavier design than the DeBakey® Vascular Tissue Forceps (351800 through 351814 on pg. 50). The distal portion is tapered to a 3 mm DeBakey®-type tip.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 342980    | Straight, 3 mm tip, 6-1/8" (15.5 cm) |
| 342981    | Straight, 3 mm tip, 7-3/4" (19.5 cm) |
| 342982    | Straight, 3 mm tip, 9-1/2" (24 cm)   |
| 342984    | Straight, 3 mm tip, 12" (30.5 cm)    |
| 342986    | Angled, 3 mm tip, 12" (30.5 cm)      |



## FORCEPS AND CLAMPS

### Tissue Forceps

#### BEASLEY-BABCOCK TISSUE FORCEPS



| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 057302    | 6 x 6 teeth, 5-1/4" (13.5 cm) |

#### BONNEY TISSUE FORCEPS



Strong, rugged construction. Serrated tying platform.  
Non-slip handles.

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 640524    | 1 x 2 teeth, 3.3 mm tip, 6-3/4" (17 cm) |
| 181327    | 2 x 3 teeth, 6-3/4" (17 cm)             |

#### BROWN TISSUE FORCEPS



Small 9 x 9 sets of teeth run longitudinally along  
narrow jaws.

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 181255    | 6" (15 cm)     |
| 342965    | 8" (20.5 cm)   |
| 342967    | 9-1/2" (24 cm) |



## FORCEPS AND CLAMPS

### Tissue Forceps

#### ADSON BROWN LASER TISSUE FORCEPS



Ebonized.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 447030    | Straight, 9 x 9 teeth, 4-3/4" (12 cm) |

#### EWALD TISSUE FORCEPS



Pilling® Pattern, satin finish. Delicate tips designed for careful handling of tissue. Handles are broad for comfortable, secure grasp.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 181240    | 1 x 2 teeth, 4-3/4" (12 cm) |

#### HUDSON TISSUE FORCEPS



| CATALOG # | DESCRIPTION                                  |
|-----------|----------------------------------------------|
| 057215    | 1 x 2 teeth, heavy, 1.5 mm tip, 5" (12.5 cm) |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### CALLISON-ADSON TISSUE FORCEPS



057220

Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                                 |
|-----------|-------------------------------------------------------------|
| 057220    | 1 x 2 teeth with tying platform, 1.3 mm tip, 4-3/4" (12 cm) |

#### ADSON TISSUE FORCEPS



KM39057

KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                   |
|-----------|-----------------------------------------------|
| KM39056   | Micro 1 x 2 teeth, 0.8 mm tip, 4-3/4" (12 cm) |
| KM39057   | Micro 1 x 2 teeth, 0.9 mm tip, 6" (15 cm)     |
| KM39064   | Delicate 1 x 2 teeth, 4-3/4" (12 cm)          |
| KM39050   | 1 x 2 teeth, 1.2 mm tip, 6" (15 cm)           |

#### ADSON-DEBAKEY® TISSUE FORCEPS



KM39054

KMedic® Pattern, satin finish. DeBakey® Teeth.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| KM39054   | 1.5 mm tip, 4-3/4" (12 cm) |

#### BROWN-ADSON TISSUE FORCEPS



181235

Pilling® Pattern, satin finish. Delicate tips designed for careful handling of tissue. Handles are broad for comfortable, secure grasp.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 181235    | 9 x 9 teeth, 4-3/4" (12 cm) |

#### BROWN-ADSON TISSUE FORCEPS



057250


Weck® Pattern, satin finish.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 057250    | 7 x 7 teeth, 2 mm tip, 4-3/4" (12 cm) |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### ADSON TISSUE FORCEPS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                           |
|-----------|---------------------------------------|
| 181223    | 1 x 2 teeth, 4-3/4" (12 cm)           |
| 181225    | 2 x 3 teeth, 5" (12.5 cm)             |
| 181228    | 1 x 2 teeth, serrated, 4-3/4" (12 cm) |

#### ADSON TISSUE FORCEPS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                                   |
|-----------|---------------------------------------------------------------|
| 057207    | 1 x 2 teeth, 1.3 mm tip, 4-3/4" (12 cm)                       |
| 057210    | Standard, 1 x 2 teeth, 1.4 mm tip, 4-3/4" (12 cm)             |
| 057212    | Delicate, 1 x 2 teeth, 1.2 mm tip, 5" (12.5 cm)               |
| 057214    | Delicate, serrated, 1 x 2 teeth, 1.3 mm tip, 4-7/8" (12.5 cm) |
| 057230    | 2 x 3 teeth, 2 mm tip, 5" (12.5 cm)                           |

## FORCEPS AND CLAMPS

### Tissue Forceps

#### TISSUE FORCEPS



KM38304

KMedic® Pattern, satin finish.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| KM38304   | 3 x 4 teeth, 5-1/2" (14 cm) |

#### THUMB TISSUE FORCEPS



467115

Weck® Pattern, bright finish. 1 x 2 teeth.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 467100    | 2.0 mm tip, 4-1/2" (11.5 cm) |
| 467105    | 2.2 mm tip, 5" (12.5 cm)     |
| 467110    | 2.2 mm tip, 5-5/8" (14.5 cm) |
| 467115    | 2.3 mm tip, 6-1/4" (16 cm)   |
| 467120    | 2.7 mm tip, 8" (20.5 cm)     |
| 467125    | 3.4 mm tip, 9-5/8" (24.5 cm) |

#### TISSUE FORCEPS



181310

Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                   |
|-----------|-----------------------------------------------|
| 181310    | 1.5 mm tip, 1 x 2 teeth, serrated, 6" (15 cm) |

## FORCEPS AND CLAMPS

### Dressing Forceps

#### BRIGHAM DRESSING FORCEPS



466170

| CATALOG # | DESCRIPTION                                       |
|-----------|---------------------------------------------------|
| 466170    | Serrated, 1.7 mm dissecting tip, 6-1/8" (15.5 cm) |

#### SEMKEN DRESSING FORCEPS



KM39086

| CATALOG # | DESCRIPTION                      |
|-----------|----------------------------------|
| 181270    | Serrated, 5" (12.5 cm)           |
| KM39086   | Serrated, 1.8 mm tip, 6" (15 cm) |

### Tissue Forceps

#### THUMB TISSUE FORCEPS



181131

Pilling® Pattern, satin finish. 1 x 2 teeth.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 181127    | 2 mm tip, 5" (12.5 cm)       |
| 181129    | 2.6 mm tip, 5-1/2" (14 cm)   |
| 181131    | 2.4 mm tip, 6" (15 cm)       |
| 181133    | 3.3 mm tip, 7-1/4" (18.5 cm) |
| 181139    | 3.4 mm tip, 11-3/4" (30 cm)  |

#### THUMB TISSUE FORCEPS



181166

Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                        |
|-----------|----------------------------------------------------|
| 181155    | 1 x 2 teeth, extra-fine 2.6 mm tip, 5-1/2" (14 cm) |
| 181166    | 2 x 3 teeth, 6-1/4" (16 cm)                        |

## FORCEPS AND CLAMPS

### Dressing Forceps

#### NEVINS FORCEPS



| CATALOG # | DESCRIPTION                                   |
|-----------|-----------------------------------------------|
| 466160    | Straight, 1.5 mm serrated tip, 6-3/4" (17 cm) |


#### POTTS-SMITH DRESSING FORCEPS



Fine, serrated tips.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 351850    | 7-1/4" (18.5 cm) |
| 351852    | 8-1/4" (21 cm)   |
| 351854    | 9-5/8" (24.5 cm) |
| 351856    | 10-1/4" (26 cm)  |
| 351858    | 11-7/8" (30 cm)  |

#### POTTS-SMITH TECHNIQUE® DRESSING FORCEPS



Narrow, serrated, TC tips.

| CATALOG # | DESCRIPTION         |
|-----------|---------------------|
| 351862    | 7" (18 cm), TC      |
| 351863    | 8" (20.5 cm), TC    |
| 351864    | 9-3/4" (25 cm), TC  |
| 351865    | 10-1/4" (26 cm), TC |



## FORCEPS AND CLAMPS

### Dressing Forceps

#### JANSEN BAYONET DRESSING FORCEPS



466150

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| KM55100   | Serrated, 2.1 mm tip, 6-1/4" (16 cm)   |
| 466150    | Serrated, 2.4 mm tip, 7-5/8" (19.5 cm) |

#### LUCAE BAYONET DRESSING FORCEPS



027100

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 027100    | Delicate, 1.5 mm tip, 5-1/2" (14 cm) |

#### LUCAE BAYONET DRESSING FORCEPS



031100

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 031100    | Serrated, 2 mm tip, 6-1/2" (16.5 cm) |

#### THUMB BAYONET DRESSING FORCEPS



395118

| CATALOG # | DESCRIPTION                      |
|-----------|----------------------------------|
| 395118    | Angled, smooth, 7-1/4" (18.5 cm) |

## FORCEPS AND CLAMPS

### Dressing Forceps

#### IRIS DRESSING FORCEPS



| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 181600    | Straight, 4" (10 cm)        |
| 181602    | Slightly curved, 4" (10 cm) |

#### DRESSING FORCEPS



| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| KM38082   | Curved, fine point serrated, 5" (12.5 cm) |

#### GRUENWALD BAYONET FORCEPS



| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 466151    | Serrated, 2.3 mm tip, 8-3/8" (21.5 cm) |

## FORCEPS AND CLAMPS

### Dressing Forceps

#### CUSHING BAYONET DRESSING FORCEPS



| CATALOG # | DESCRIPTION        |
|-----------|--------------------|
| 351861    | 7-1/2" (19 cm), TC |

#### CUSHING GUTSCH BAYONET DRESSING FORCEPS



Serrated handles, scraper end.

| CATALOG # | DESCRIPTION                            |
|-----------|----------------------------------------|
| 466220    | Serrated, 1.5 mm tip, 7-1/4" (18.5 cm) |

#### GERALD DRESSING FORCEPS



These are full-size dressing forceps with only the distal portion narrowed to provide light, delicate tips that maintain alignment.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 351900    | Serrated, 7" (18 cm)     |
| 351902    | Serrated, 9-3/4" (25 cm) |

#### GERALD DRESSING FORCEPS TC



| CATALOG # | DESCRIPTION                                          |
|-----------|------------------------------------------------------|
| 466225    | Straight, serrated, 1.1 mm tip, 6-7/8" (17.5 cm), TC |

## FORCEPS AND CLAMPS

### Dressing Forceps

#### ADSON LASER TISSUE FORCEPS



Ebonized.

| CATALOG # | DESCRIPTION                 |
|-----------|-----------------------------|
| 447025    | 1 x 2 teeth, 4-3/4" (12 cm) |

#### CUSHING DRESSING FORCEPS



Serrated handles, scraper end.

| CATALOG # | DESCRIPTION                         |
|-----------|-------------------------------------|
| 181410    | Straight, serrated tips, 7" (18 cm) |

#### CUSHING DRESSING FORCEPS TC



| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| 351860    | Straight, 6-7/8" (17.5 cm), TC |

#### CUSHING DRESSING FORCEPS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 466200    | Straight, serrated, 1.5 mm tip, 7" (18 cm) |

#### CUSHING GUTSCH DRESSING FORCEPS



Tungsten carbide.

| CATALOG # | DESCRIPTION                                |
|-----------|--------------------------------------------|
| 466211    | Serrated, 1.8 mm tip, 6-7/8" (17.5 cm), TC |
| 466210    | Serrated, 1.6 mm tip, 7" (18 cm), TC       |

## FORCEPS AND CLAMPS

### Dressing Forceps

#### ADSON DRESSING FORCEPS



Weck® Pattern, serrated tip.

| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| 057200    | Standard, 1.3 mm tip, 4-3/4" (12 cm) |
| 057202    | Delicate, 1.2 mm tip, 5" (12.5 cm)   |

#### ADSON DRESSING FORCEPS TC



Weck® Pattern, tungsten carbide.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 057205    | Standard, 4-3/4" (12 cm), TC |

#### ADSON DRESSING FORCEPS



KMedic® Pattern.

| CATALOG # | DESCRIPTION                                    |
|-----------|------------------------------------------------|
| KM39058   | Micro tip, serrated, 1 mm tip, 4-3/4" (12 cm)  |
| KM39060   | Delicate, serrated, 1.2 mm tip, 4-3/4" (12 cm) |
| KM39062   | Standard tip, serrated, 4-3/4" (12 cm)         |
| KM39059   | Micro serrated, 1 mm tip, 6" (15 cm)           |
| KM39052   | Serrated, 1.2 mm tip, 6" (15 cm)               |

#### EWALD DRESSING FORCEPS



Pilling® Pattern, satin finish. Delicate tips designed for careful handling of tissue. Handles are broad for comfortable, secure grasp.

| CATALOG # | DESCRIPTION              |
|-----------|--------------------------|
| 181237    | Serrated, 4-3/4" (12 cm) |

## FORCEPS AND CLAMPS

### Dressing Forceps

#### THUMB DRESSING FORCEPS



181106

Pilling® Pattern, satin finish. Standard, serrated tips.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 181102    | 2.7 mm tip, 5" (12.5 cm)   |
| 181104    | 3.4 mm tip, 5-1/2" (14 cm) |
| 181106    | 3.5 mm tip, 6" (15 cm)     |
| 181108    | 4 mm tip, 7-1/4" (18.5 cm) |

#### THUMB DRESSING FORCEPS



465115

Weck® Pattern, bright finish. Standard, serrated tips.

| CATALOG # | DESCRIPTION                  |
|-----------|------------------------------|
| 465100    | 2.6 mm tip, 4-5/8" (11.5 cm) |
| 465105    | 2.7 mm tip, 5" (12.5 cm)     |
| 465110    | 3 mm tip, 5-5/8" (14.5 cm)   |
| 465115    | 3.5 mm tip, 6-1/4" (16 cm)   |
| 465120    | 3.6 mm tip, 7-7/8" (20 cm)   |
| 465125    | 4 mm tip, 9-7/8" (25 cm)     |
| 465130    | 4.7 mm tip, 12" (30.5 cm)    |

#### THUMB FORCEPS



465160

Weck® Pattern, bright finish. Narrow, serrated tips.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 465150    | 2 mm tip, 4-5/8" (11.5 cm) |
| 465155    | 2 mm tip, 5" (12.5 cm)     |
| 465160    | 2 mm tip, 5-5/8" (14.5 cm) |
| 465165    | 2 mm tip, 6" (15 cm)       |
| 465170    | 2.1 mm tip, 7-7/8" (20 cm) |

#### ADSON DRESSING FORCEPS



181220

Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION    |
|-----------|----------------|
| 181220    | 4-3/4" (12 cm) |



## FORCEPS AND CLAMPS

### Splinter Forceps

#### VANGUARD SPLINTER LIBERATOR



KM39534

| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| KM39534   | Trocar in reversible handle, 5" (12.5 cm) |
| KM39536   | Trocar point only                         |

#### VIRTUS SPLINTER FORCEPS



KM39478

| CATALOG # | DESCRIPTION                    |
|-----------|--------------------------------|
| KM39478   | Straight, serrated, 6" (15 cm) |

#### WILLIAMS SPLINTER FORCEPS



KM39509

KM39511

KM39513

| CATALOG # | DESCRIPTION                                  |
|-----------|----------------------------------------------|
| KM39509   | Straight, serrated, 5-1/2" (14 cm)           |
| KM39511   | 45 degree angled, serrated, 5-1/4" (13.5 cm) |
| KM39513   | 90 degree angled, serrated, 5-1/2" (14 cm)   |

## FORCEPS AND CLAMPS

### Splinter Forceps

#### CARMALT SPLINTER FORCEPS



| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| KM39410   | Straight, serrated, 4-3/8" (11 cm) |
| KM39412   | Curved, serrated, 4-3/8" (11 cm)   |

#### FEILCHENFELD SPLINTER FORCEPS



| CATALOG # | DESCRIPTION                          |
|-----------|--------------------------------------|
| KM39386   | Straight, serrated, 3-1/4" (8.5 cm)  |
| KM39390   | Straight, serrated, 4-1/2" (11.5 cm) |

#### PLAIN SPLINTER FORCEPS



| CATALOG # | DESCRIPTION                                      |
|-----------|--------------------------------------------------|
| KM39364   | Straight, serrated, 3-1/2" (9 cm)                |
| KM39366   | Straight, serrated, 0.6 mm tip, 4-1/2" (11.5 cm) |

#### SPLINTER FORCEPS



| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 181556    | Straight, serrated, 3-1/8" (8 cm) |
| 181557    | Straight, serrated, 4" (10 cm)    |
| 181559    | Straight, serrated, 5" (12.5 cm)  |

## FORCEPS AND CLAMPS

### Towel Forceps

#### JONES CROSS-ACTION TOWEL FORCEPS



121640

| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 121640    | 3-1/2" (9 cm) |

#### PEERS TOWEL FORCEPS



121650

Double opposing row of non-penetrating teeth.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 121650    | 5-3/4" (14.5 cm) |

#### ROEDER TOWEL FORCEPS



121610

| CATALOG # | DESCRIPTION                       |
|-----------|-----------------------------------|
| 121610    | With ball stops, 5-1/4" (13.5 cm) |

## FORCEPS AND CLAMPS

### Towel Clamps

#### BACKHAUS TOWEL CLAMPS



454305

Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 454300    | 3-1/2" (9 cm)    |
| 454305    | 5-1/4" (13.5 cm) |

#### BACKHAUS TOWEL FORCEPS



121600

Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 121605    | 3-1/2" (9 cm)    |
| 121600    | 5-1/4" (13.5 cm) |

#### BACKHAUS TOWEL FORCEPS



447050

| CATALOG # | DESCRIPTION   |
|-----------|---------------|
| 447050    | 3-1/2" (9 cm) |

## FORCEPS AND CLAMPS

### Towel Clamps

#### EDNA TOWEL FORCEPS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                                         |
|-----------|-----------------------------------------------------|
| 121655    | Single row of non-penetrating teeth, 5-1/2" (14 cm) |
| 121656    | 4" (10 cm)                                          |

#### EDNA NON-PERFORATING TOWEL CLAMPS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 454320    | 3-3/4" (9.5 cm)  |
| 454322    | 5-3/8" (13.5 cm) |

## FORCEPS AND CLAMPS

### Tubing Clamps

#### VORSE TUBING OCCLUDING CLAMPS



Heavy construction with jaws designed for positive occlusion of Silastic®, Tygon® and similar thick-walled tubings.

| CATALOG # | DESCRIPTION                                        |
|-----------|----------------------------------------------------|
| 351373    | For tubing 1/4" I.D. (3/8" O.D.), 6-1/8" (15.5 cm) |
| 351374    | For tubing 3/8" I.D. (5/8" O.D.), 7-3/8" (18.5 cm) |
| 351375    | For tubing 1/2" I.D. (7/8" O.D.), 8-1/8" (20.5 cm) |

#### VORSE TUBING OCCLUDING CLAMPS WITH SAFETY GUARD



The familiar Vorse-type tube-occluding clamp supplied with a safety guard to prevent the box-lock from puncturing the tubing. Heavy construction with jaw designed for positive non-slip occlusion of Silastic®, Tygon® and similar thick-walled tubing.


| CATALOG # | DESCRIPTION                                        |
|-----------|----------------------------------------------------|
| 351376    | For tubing 1/4" I.D. (3/8" O.D.), 6-1/8" (15.5 cm) |
| 351377    | For tubing 3/8" I.D. (5/8" O.D.), 7-3/8" (18.5 cm) |
| 351378    | For tubing 1/2" I.D. (7/8" O.D.), 8-1/8" (20.5 cm) |



## FORCEPS AND CLAMPS

### Tubing Clamps

#### NICHOLS TUBING HOLDER




351370

Specially designed to hold and support the weight of the various diameters of tubing used by the pump-oxygenator, keeping them "in-line."

| CATALOG # | DESCRIPTION                                                                                              |
|-----------|----------------------------------------------------------------------------------------------------------|
| 351370    | 4-1/2" x 7-1/4" (11.4 x 18.5 cm) metal plate is perforated so it may be held in place with towel clamps. |

#### PRESBYTERIAN HOSPITAL TUBING CLAMP



351358


Durable construction with smooth jaws for clamping off tubing without damage.

| CATALOG # | DESCRIPTION      |
|-----------|------------------|
| 351358    | 6-1/2" (16.5 cm) |

## FORCEPS AND CLAMPS

### Tubing Clamps

#### GUNNAR ROLLER CLAMPS



| CATALOG # | DESCRIPTION                                                       |
|-----------|-------------------------------------------------------------------|
| 351390    | For "milking" pump tubing 1/2" (12.5 mm) diameter, 7-1/8" (18 cm) |

#### LORE TUBE HOLDING FORCEPS



Full serrations, serrated jaws grasp drape without puncturing. Fitting on shanks holds tubing securely and without drag when not in use. Suction tube may be stored in protective metal sleeve.

| CATALOG # | DESCRIPTION                        |
|-----------|------------------------------------|
| 536275    | Angled, 2.6 mm tip, 6-1/4" (16 cm) |

#### MORSE TOWEL CLIP TUBING HOLDERS



For holding extracorporeal tubing to drapes giving positive fixation. The serrated inner surface of the tube holding ring closes down on the tubing but does not collapse or reduce the lumen. Sizes listed are for outside diameter of extracorporeal tubing.

| CATALOG # | DESCRIPTION                |
|-----------|----------------------------|
| 351351    | For 3/16" and 1/4" tubing  |
| 351352    | For 5/16" and 3/8" tubing  |
| 351353    | For 7/16" and 1/2" tubing  |
| 351354    | For 9/16" and 5/8" tubing  |
| 351355    | For 11/16" and 3/4" tubing |

## FORCEPS AND CLAMPS

### Sponge Forceps

#### BALLENGER SPONGE FORCEPS



| CATALOG # | DESCRIPTION                               |
|-----------|-------------------------------------------|
| 121410    | Straight, smooth jaws, 7-1/4" (18.5 cm)   |
| 121412    | Straight, serrated jaws, 7-1/4" (18.5 cm) |
| 121422    | Curved, serrated jaws, 7" (18 cm)         |


#### FOERSTER SPONGE FORCEPS



Pilling® Pattern, satin finish.

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 121415    | Straight, smooth jaws, 9-3/4" (25 cm)   |
| 121417    | Straight, serrated jaws, 9-3/4" (25 cm) |
| 121427    | Curved, serrated jaws, 9-1/2" (24 cm)   |

#### FOERSTER SPONGE FORCEPS



Weck® Pattern, bright finish.

| CATALOG # | DESCRIPTION                                   |
|-----------|-----------------------------------------------|
| 454100    | Straight, cross serrated jaws, 9-1/2" (24 cm) |

## KNIFE HANDLES

Specialty

### GOULIAN DEPTH GAUGE AND GUARDS



050203

2-3/8" (6 cm). Use with Webster knife handle (050230), Wecprep blades (450205).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 050201    | 0.004"      |
| 050202    | 0.006"      |
| 050203    | 0.008"      |
| 050204    | 0.010"      |
| 050205    | 0.012"      |
| 050206    | 0.014"      |
| 050207    | 0.016"      |
| 050208    | 0.018"      |
| 050209    | 0.020"      |
| 050211    | 0.024"      |
| 050212    | 0.026"      |
| 050213    | 0.028"      |

### WECREP BLADES



450205

Single edge carbon steel blade. Use with 050200 and 050230.

| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 450205    | 2-1/4" (5.5 cm), carton of 50 |

### HOLINGER LARYNGEAL KNIFE HANDLE



506705

For use with surgical blades no. 10, 11, 12 and 15.

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 506705    | 11" (28 cm) |

## KNIFE HANDLES

General

### ROUND SHAFT BEAVER BLADE HANDLES



| CATALOG # | DESCRIPTION             |
|-----------|-------------------------|
| 456002    | Replacement collet      |
| 456005    | #3KM, 2-31/32" (7.5 cm) |
| 456006    | #3K, 4-1/8" (10.5 cm)   |
| 456007    | #3KL, 6" (15 cm)        |

Specialty

### WEBSTER SKIN GRAFT KNIFE HANDLE



Use with blade no. 450205, 7-1/4" (18.5 cm).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 050230    | Handle only |

### GOULIAN SKIN GRAFT KNIFE SET



Includes: 1 Webster knife handle (050230), 5 Wecprep blades (450205), 3 depth gauges and guards (050203, 050204, 050205).

| CATALOG # | DESCRIPTION |
|-----------|-------------|
| 050200    | Complete    |

## KNIFE HANDLES

General

### KNIFE HANDLES



352950



352951



352952



352953



352954



352957



352959

Correctly fit knife handle to standard surgical blades. The distal ends are tapered and beveled for use in blunt dissection. Each handle is contoured for firm grip, strength, manipulation and comfort. They are sold one handle each per package.

| CATALOG # | DESCRIPTION                                                                                         |
|-----------|-----------------------------------------------------------------------------------------------------|
| 352950    | Size 3, fits blade sizes 10, 11, 12 and 15                                                          |
| 352951    | Size 3L, fits blade sizes 10, 11, 12 and 15. Measures 8-3/8" (21.5 cm) not including blade.         |
| 352952    | Size 3LA, fits blade sizes 10, 11, 12 and 15                                                        |
| 352953    | Size 4, fits blade sizes 20, 21, 22, 23, 24 and 25                                                  |
| 352954    | Size 4L, fits blade sizes 20, 21, 22, 23, 24 and 25. Measures 8-3/8" (21.5 cm) not including blade. |
| 352957    | Size 7, fits blade sizes 10, 11, 12 and 15                                                          |
| 352959    | Size 9, fits blade sizes 10, 11, 12 and 15                                                          |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## STILLE skalpellhandtag | STILLE Scalpel Handle

The most common scalpel handles made with highest precision.

|              |      |             |             |
|--------------|------|-------------|-------------|
| 119-50010-3  | No 3 | Blade 10-15 | 12 cm 5"    |
| 119-50010-13 | No 3 | Blade 10-15 | 21 cm 8.25" |
| 119-50010-4  | No 4 | Blade 18-25 | 13 cm 5.25" |
| 119-50010-14 | No 4 | Blade 18-25 | 21 cm 8.25" |
| 119-50010-7  | No 7 | Blade 10-15 | 16 cm 6.25" |



## STILLE Barron skalpellhandtag | STILLE Barron Scalpel Handle

8-sided handle for precise handling and rotation.

|           |             |            |
|-----------|-------------|------------|
| 119-50012 | Blade 10-15 | 14 cm 5.5" |
|-----------|-------------|------------|

## STILLE skalpellhandtag | STILLE Scalpel Handle

Round handle guarantees precise finger tip control.

|           |             |               |
|-----------|-------------|---------------|
| 119-50015 | Blade 10-15 | 14.5 cm 5.75" |
|-----------|-------------|---------------|



## STILLE uterussond | STILLE Uterine Sound

Uterine sound with malleable silver shaft. Graduated in cm. Drop-shaped tip.

|           |              |
|-----------|--------------|
| 119-28452 | 27 cm 10.75" |
|-----------|--------------|



## Bergström-STILLE muskelbiopsi kanyl | Bergström-STILLE Muscle biopsy Cannula

Muscle biopsy cannula. Consists of a sharptipped, hollow outer needle with a small opening near the tip. A punch with a sharp edge fits tightly into the needle. The stylet is used for pushing out biopsy specimens from the punch.

|              |           |
|--------------|-----------|
| 119-29187-50 | Diam 5 mm |
| 119-29187-60 | Diam 6 mm |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



Oretorp-STILLE sond | Oretorp-STILLE Probe

Blunt ended probe for arthroscopy. Graduated shaft.



STILLE sond | STILLE Probe


Blunt ended probe for arthroscopy. Graduated shaft.



Övriga instrument | Various instruments | Sonstige  
Instrumente | Altri strumenti | Instrumentos varios |  
様々な器具

Ergonomic design ensures  
comfortable grip and exact  
handling of delicate tissue.

All instruments are manufactured in  
selected top quality surgical steel.



NEW  
Titanium

**Mikronålförare** | Micro needle holders | Mikronadelhalter |  
Micro porta-agujas | マイクロ持針器

## STILLE Titan mikronålförare | STILLE Titanium Micro Needle Holder

Round handle, 1.2 mm Diamond Jaws and **Streamline Box Lock**  
for sutures 5-0 and smaller. With or without ratchet.

|                |              |                 |          |               |
|----------------|--------------|-----------------|----------|---------------|
| 117-C115-18-TI | Diamond Jaws | Without Ratchet | Straight | 18 cm 7 in    |
| 117-C115-21-TI | Diamond Jaws | Without Ratchet | Straight | 21 cm 8.25 in |
| 117-C115-23-TI | Diamond Jaws | Without Ratchet | Straight | 23 cm 9 in    |
| 117-C125-18-TI | Diamond Jaws | With Ratchet    | Straight | 18 cm 7 in    |
| 117-C125-21-TI | Diamond Jaws | With Ratchet    | Straight | 21 cm 8.25 in |
| 117-C125-23-TI | Diamond Jaws | With Ratchet    | Straight | 23 cm 9 in    |
| 117-C126-18-TI | Diamond Jaws | With Ratchet    | Curved   | 18 cm 7 in    |
| 117-C126-21-TI | Diamond Jaws | With Ratchet    | Curved   | 21 cm 8.25 in |
| 117-C126-23-TI | Diamond Jaws | With Ratchet    | Curved   | 23 cm 9 in    |



PRIME SURGICAL

NEW  
Titanium

**Mikronålförare** | Micro needle holders | Mikron  
Micro porta-agujas | マイクロ持針器

## STILLE Titan mikronålförare | STILLE Titanium Micro Needle Holder

Round handle, fine 0.8 mm Diamond Jaws and **Regular Box Lock**  
for sutures 7-0 and smaller. With or without ratchet.

|                |              |                 |          |               |
|----------------|--------------|-----------------|----------|---------------|
| 117-C110-18-TI | Diamond Jaws | Without Ratchet | Straight | 18 cm 7 in    |
| 117-C110-21-TI | Diamond Jaws | Without Ratchet | Straight | 21 cm 8.25 in |
| 117-C110-23-TI | Diamond Jaws | Without Ratchet | Straight | 23 cm 9 in    |
| 117-C111-15-TI | Diamond Jaws | Without Ratchet | Curved   | 15 cm 6 in    |
| 117-C111-18-TI | Diamond Jaws | Without Ratchet | Curved   | 18 cm 7 in    |
| 117-C111-21-TI | Diamond Jaws | Without Ratchet | Curved   | 21 cm 8.25 in |
| 117-C120-18-TI | Diamond Jaws | With Ratchet    | Straight | 18 cm 7 in    |
| 117-C120-21-TI | Diamond Jaws | With Ratchet    | Straight | 21 cm 8.25 in |
| 117-C120-23-TI | Diamond Jaws | With Ratchet    | Straight | 23 cm 9 in    |
| 117-C121-18-TI | Diamond Jaws | With Ratchet    | Curved   | 18 cm 7 in    |
| 117-C121-21-TI | Diamond Jaws | With Ratchet    | Curved   | 21 cm 8.25 in |
| 117-C121-23-TI | Diamond Jaws | With Ratchet    | Curved   | 23 cm 9 in    |



PRIME

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE mikro nålförare | STILLE Micro Needle Holder

Round handle, 1.2 mm Diamond Jaws and **Streamline Box Lock** for sutures 5-0 and smaller. With ratchet.

|             |              |              |          |               |
|-------------|--------------|--------------|----------|---------------|
| 117-C125-15 | Diamond Jaws | With Ratchet | Straight | 15 cm 6 in    |
| 117-C125-18 | Diamond Jaws | With Ratchet | Straight | 18 cm 7 in    |
| 117-C125-21 | Diamond Jaws | With Ratchet | Straight | 21 cm 8.25 in |
| 117-C125-23 | Diamond Jaws | With Ratchet | Straight | 23 cm 9 in    |
| 117-C126-15 | Diamond Jaws | With Ratchet | Curved   | 15 cm 6 in    |
| 117-C126-18 | Diamond Jaws | With Ratchet | Curved   | 18 cm 7 in    |
| 117-C126-21 | Diamond Jaws | With Ratchet | Curved   | 21 cm 8.25 in |
| 117-C126-23 | Diamond Jaws | With Ratchet | Curved   | 23 cm 9 in    |



## STILLE mikro nålförare | STILLE Micro Needle Holder

Extra heavy model for handling of larger needles with TC-inserts. Intended for 2/0 suture and finer.

|             |            |              |             |
|-------------|------------|--------------|-------------|
| 117-C130    |            |              |             |
| 117-C140    |            |              |             |
| 117-C130-19 | TC-inserts | Without lock | 19 cm 7.5"  |
| 117-C130-23 | TC-inserts | Without lock | 23 cm 9"    |
| 117-C130-25 | TC-inserts | Without lock | 25 cm 9.75" |
| 117-C140-19 | TC-inserts | With lock    | 19 cm 7.5"  |
| 117-C140-23 | TC-inserts | With lock    | 23 cm 9"    |
| 117-C140-25 | TC-inserts | With lock    | 25 cm 9.75" |



## STILLE mikro nålförare | STILLE Micro Needle Holder

Round handle, fine 0.8 mm Diamond Jaws and **Regular Box Lock** for sutures 7-0 and smaller. Without ratchet.

|                                                                                              |                                                                                              |          |             |
|----------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|----------|-------------|
|  117-C110 |  117-C111 |          |             |
| 117-C110-15                                                                                  | Diamond Jaws                                                                                 | Straight | 15 cm 6"    |
| 117-C110-18                                                                                  | Diamond Jaws                                                                                 | Straight | 18 cm 7"    |
| 117-C110-21                                                                                  | Diamond Jaws                                                                                 | Straight | 21 cm 8.25" |
| 117-C110-23                                                                                  | Diamond Jaws                                                                                 | Straight | 23 cm 9"    |
| 117-C111-15                                                                                  | Diamond Jaws                                                                                 | Curved   | 15 cm 6"    |
| 117-C111-18                                                                                  | Diamond Jaws                                                                                 | Curved   | 18 cm 7"    |
| 117-C111-21                                                                                  | Diamond Jaws                                                                                 | Curved   | 21 cm 8.25" |
| 117-C111-23                                                                                  | Diamond Jaws                                                                                 | Curved   | 23 cm 9"    |



## STILLE mikro nålförare | STILLE Micro Needle Holder

Round handle, fine 0.8 mm Diamond Jaws and **Regular Box Lock** for sutures 7-0 and smaller. With ratchet.

| 117-C120    | 117-C121     |          |             |
|-------------|--------------|----------|-------------|
| 117-C120-15 | Diamond Jaws | Straight | 15 cm 6"    |
| 117-C120-18 | Diamond Jaws | Straight | 18 cm 7"    |
| 117-C120-21 | Diamond Jaws | Straight | 21 cm 8.25" |
| 117-C120-23 | Diamond Jaws | Straight | 23 cm 9"    |
| 117-C121-15 | Diamond Jaws | Curved   | 15 cm 6"    |
| 117-C121-18 | Diamond Jaws | Curved   | 18 cm 7"    |
| 117-C121-21 | Diamond Jaws | Curved   | 21 cm 8.25" |
| 117-C121-23 | Diamond Jaws | Curved   | 23 cm 9"    |



## STILLE mikro nålförare | STILLE Micro Needle Holder

Round handle, 1.2 mm Diamond Jaws and **Streamline Box Lock** for sutures 5-0 and smaller. Without ratchet.

|             |              |                 |          |               |
|-------------|--------------|-----------------|----------|---------------|
| 117-C115-15 | Diamond Jaws | Without Ratchet | Straight | 15 cm 6 in    |
| 117-C115-18 | Diamond Jaws | Without Ratchet | Straight | 18 cm 7 in    |
| 117-C115-21 | Diamond Jaws | Without Ratchet | Straight | 21 cm 8.25 in |
| 117-C115-23 | Diamond Jaws | Without Ratchet | Straight | 23 cm 9 in    |
| 117-C116-15 | Diamond Jaws | Without Ratchet | Curved   | 15 cm 6 in    |
| 117-C116-18 | Diamond Jaws | Without Ratchet | Curved   | 18 cm 7 in    |
| 117-C116-21 | Diamond Jaws | Without Ratchet | Curved   | 21 cm 8.25 in |
| 117-C116-23 | Diamond Jaws | Without Ratchet | Curved   | 23 cm 9 in    |





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

NEW  
Titanium

**Mikropincetter** | Micro forceps | Mikropinzetten | Micro pinze | Micro pinzas | マイクロ鉗子

## STILLE Titan mikro ringpincett | STILLE Titanium Micro Ring Forceps

Round handle, with fine ring tips, counter-balanced.

|                |              |                      |          |               |
|----------------|--------------|----------------------|----------|---------------|
| 117-B140-18-Ti | Diamond Jaws | 1.0x0.5 mm ring tips | Straight | 18 cm 7 in    |
| 117-B140-21-Ti | Diamond Jaws | 1.0x0.5 mm ring tips | Straight | 21 cm 8.25 in |
| 117-B140-23-Ti | Diamond Jaws | 1.0x0.5 mm ring tips | Straight | 23 cm 9 in    |



## STILLE Titan DeBaKey mikropincett | STILLE Titanium Micro DeBaKey Forceps

Round handle, with 1.2 mm fine atraumatic DeBaKey tips.

|                |                     |          |               |
|----------------|---------------------|----------|---------------|
| 117-B160-18-Ti | 1.2 mm DeBaKey tips | Straight | 18 cm 7 in    |
| 117-B160-21-Ti | 1.2 mm DeBaKey tips | Straight | 21 cm 8.25 in |



NEW  
Titanium

**Mikropincetter** | Micro forceps | Mikropinzetten | Micro pinze | Micro pinzas | マイクロ鉗子

## STILLE Titan mikropincett | STILLE Titanium Micro Forceps

Round handle, with tying platform, straight or curved 0.4 mm tips.

|                |              |          |               |
|----------------|--------------|----------|---------------|
| 117-B120-18-Ti | Diamond Jaws | Straight | 18 cm 7 in    |
| 117-B120-21-Ti | Diamond Jaws | Straight | 21 cm 8.25 in |
| 117-B120-23-Ti | Diamond Jaws | Straight | 23 cm 9 in    |
| 117-B120-25-Ti | Diamond Jaws | Straight | 25 cm 10 in   |
| 117-B121-18-Ti | Diamond Jaws | Curved   | 18 cm 7 in    |
| 117-B121-21-Ti | Diamond Jaws | Curved   | 21 cm 8.25 in |
| 117-B121-23-Ti | Diamond Jaws | Curved   | 23 cm 9 in    |
| 117-B121-25-Ti | Diamond Jaws | Curved   | 25 cm 10 in   |



## STILLE Titan mikropincett | STILLE Titanium Micro Forceps

Round handle, with tying platform, counter-balanced, straight or curved 0.4 mm tips.

|                |              |          |               |
|----------------|--------------|----------|---------------|
| 117-B130-18-Ti | Diamond Jaws | Straight | 18 cm 7 in    |
| 117-B130-21-Ti | Diamond Jaws | Straight | 21 cm 8.25 in |
| 117-B131-18-Ti | Diamond Jaws | Curved   | 18 cm 7 in    |
| 117-B131-21-Ti | Diamond Jaws | Curved   | 21 cm 8.25 in |





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE mikro ringpincett | STILLE Micro Ring Forceps

Counter-balanced micro forceps with fine Diamond Jaws ring tips.

117-B140

|             |            |              |          |             |
|-------------|------------|--------------|----------|-------------|
| 117-B140-18 | 1.0x0.5 mm | Diamond Jaws | Straight | 18 cm 7"    |
| 117-B140-21 | 1.0x0.5 mm | Diamond Jaws | Straight | 21 cm 8.25" |
| 117-B140-23 | 1.0x0.5 mm | Diamond Jaws | Straight | 23 cm 9 in  |



## STILLE mikro klopincett | STILLE Micro Tissue Forceps

Round handle, with 1x2 fine teeth on 0.4 mm tips.

|             |                |          |             |
|-------------|----------------|----------|-------------|
| 117-B150-15 | 1x2 fine teeth | Straight | 15 cm 6"    |
| 117-B150-18 | 1x2 fine teeth | Straight | 18 cm 7"    |
| 117-B150-21 | 1x2 fine teeth | Straight | 21 cm 8.25" |



## STILLE DeBaKey mikropincett | STILLE Micro DeBaKey Forceps

Round handle, with 1.2 mm fine atraumatic DeBaKey tips.

|             |        |              |          |               |
|-------------|--------|--------------|----------|---------------|
| 117-B160-18 | 1.2 mm | DeBaKey Tips | Straight | 18 cm 7 in    |
| 117-B160-21 | 1.2 mm | DeBaKey Tips | Straight | 21 cm 8.25 in |



## STILLE DeBaKey mikropincett | STILLE Micro DeBaKey Forceps

Round handle, with 1.2 mm fine atraumatic DeBaKey tips, counter-balanced.

|             |        |              |          |               |
|-------------|--------|--------------|----------|---------------|
| 117-B165-18 | 1.2 mm | DeBaKey Tips | Straight | 18 cm 7 in    |
| 117-B165-21 | 1.2 mm | DeBaKey Tips | Straight | 21 cm 8.25 in |



## STILLE mikropincett | STILLE Micro Forceps

Well balanced and ergonomic micro forceps of ultra high quality. With tying platform and straight or curved jaws. 0.4 mm tips.

117-B110 117-B111

|             |          |             |
|-------------|----------|-------------|
| 117-B110-15 | Straight | 15 cm 6"    |
| 117-B110-18 | Straight | 18 cm 7"    |
| 117-B110-21 | Straight | 21 cm 8.25" |
| 117-B110-23 | Straight | 23 cm 9 in  |
| 117-B111-15 | Curved   | 15 cm 6"    |
| 117-B111-18 | Curved   | 18 cm 7"    |
| 117-B111-21 | Curved   | 21 cm 8.25" |



## STILLE mikropincett | STILLE Micro Forceps

Same fine micro forceps, with Diamond Jaws, for increased tip durability and secure grip of fine needles. 0.4 mm tips.

117-B120 117-B121

|             |              |          |             |
|-------------|--------------|----------|-------------|
| 117-B120-15 | Diamond Jaws | Straight | 15 cm 6"    |
| 117-B120-18 | Diamond Jaws | Straight | 18 cm 7"    |
| 117-B120-21 | Diamond Jaws | Straight | 21 cm 8.25" |
| 117-B120-23 | Diamond Jaws | Straight | 23 cm 9 in  |
| 117-B121-15 | Diamond Jaws | Curved   | 15 cm 6"    |
| 117-B121-18 | Diamond Jaws | Curved   | 18 cm 7"    |
| 117-B121-21 | Diamond Jaws | Curved   | 21 cm 8.25" |



## STILLE mikropincett | STILLE Micro Forceps

Micro forceps with Diamond Jaws and Counter-balance. 0.4 mm tips.

117-B130 117-B131

|             |              |          |             |
|-------------|--------------|----------|-------------|
| 117-B130-18 | Diamond Jaws | Straight | 18 cm 7"    |
| 117-B130-21 | Diamond Jaws | Straight | 21 cm 8.25" |
| 117-B130-23 | Diamond Jaws | Straight | 23 cm 9 in  |
| 117-B131-18 | Diamond Jaws | Curved   | 18 cm 7"    |
| 117-B131-21 | Diamond Jaws | Curved   | 21 cm 8.25" |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



STILLE mikronålförare | STILLE Micro Needle Holder



STILLE Titan mikronålförare | STILLE Titanium Micro Needle Holder

## MIKRONÅLFÖRARE


STILLES mikronålförare är designade för att hålla även de finaste mikrosuturerna. Bara ett minimalt tryck behövs från kirurgens fingrar för ett säkert grepp. Alla modeller är försedda med STILLE Diamond Jaws eller på de kraftigare modellerna hårdmetallplattor. De erbjuds som raka eller böjda, med- eller utan låsning och i ett flertal olika längder. Leddesign Regular för 7-0 suturer och finare eller Streamline för 5-0 och finare (kraftiga modellerna för 2-0 till 4-0).

## MICRO NEEDLE HOLDERS

The STILLE Micro Needle Holders are designed to hold the finest sutures in a tight and secure grip. Only minimal pressure will need to be applied by surgeon's fingers. All models are made with STILLE Diamond Jaws and the heavier models with TC-inserts. They can be delivered as straight or curved, with or without ratchet and in several different lengths. Regular box lock for 7-0 sutures and finer, Streamline box lock for 5-0 sutures and finer. (heavy models for 2-0 to 4-0).

## STILLE dissektor | STILLE Dissector

Double-ended dissector in several different models




|              |                     |               |
|--------------|---------------------|---------------|
| 111-22440-18 | Fine with hole      | 18 cm 7"      |
| 111-22440-24 | Fine with hole      | 24 cm 9.5"    |
| 111-22442-19 | Coarse with hole    | 19.5 cm 7.75" |
| 111-22442-24 | Coarse with hole    | 24 cm 9.5"    |
| 111-22444    | Fine without hole   | 24.5 cm 9.75" |
| 111-22446    | Coarse without hole | 24.5 cm 9.5"  |



## Récamier STILLE uteruskyrett | Récamier STILLE Uterine curette

Malleable uterine curettes in different sizes. With shaft in non-hardened stainless steel, ensures long-time function without any surface coating falling off. Higher patient safety and longer life-time expectancy.




|              |       |       |           |
|--------------|-------|-------|-----------|
| 111-20340-5  | Sharp | 5 mm  | 30 cm 12" |
| 111-20340-7  | Sharp | 7 mm  | 30 cm 12" |
| 111-20340-8  | Sharp | 8 mm  | 30 cm 12" |
| 111-20340-10 | Sharp | 10 mm | 30 cm 12" |
| 111-20350-7  | Blunt | 7 mm  | 30 cm 12" |
| 111-20350-8  | Blunt | 8 mm  | 30 cm 12" |
| 111-20350-10 | Blunt | 10 mm | 30 cm 12" |
| 111-20350-12 | Blunt | 12 mm | 30 cm 12" |



PRIME SURGICAL

## Adson STILLE självhållande hake | Adson STILLE Self Retaining Retractor

Self retaining retractor with sharp prongs. Available in two different sizes.



|           |                  |       |             |
|-----------|------------------|-------|-------------|
| 108-26610 | 4x4 sharp prongs | 25 mm | 32 cm 12.5" |
| 108-26615 | 5x5 sharp prongs | 45 mm | 32 cm 12.5" |



## Weitlaner STILLE hake | Weitlaner STILLE Retractor

General purpose self-retaining retractor used in many procedures.

|           |                  |               |
|-----------|------------------|---------------|
| 108-26252 | 2x3 sharp prongs | 10.5 cm 4.25" |
| 108-26253 | 2x3 blunt prongs | 10.5 cm 4.25" |
| 108-26254 | 3x4 sharp prongs | 13.5 cm 5.25" |
| 108-26255 | 3x4 blunt prongs | 13.5 cm 5.25" |
| 108-26256 | 3x4 sharp prongs | 16.5 cm 6.5"  |
| 108-26257 | 3x4 blunt prongs | 16.5 cm 6.5"  |



## Alm STILLE hake | Alm STILLE Retractor


Self-retaining retractor for fine incisions. Available in two different sizes.

|           |                  |            |
|-----------|------------------|------------|
| 108-26261 | 4x4 sharp prongs | 7 cm 2.75" |
| 108-26262 | 4x4 sharp prongs | 10 cm 4"   |



## Carlens STILLE Tracheotomihake | Carlens STILLE Trachea Retractor

Double-ended trachea retractor.



108-24805

108-24805 Double-ended 22 cm 8.75"



PRIME SURGICAL

## Kocher STILLE bukspatel | Kocher STILLE Abdominal Spatula

STILLE's classical abdominal spatula, available in three different sizes.  
Also used in orthopaedics.

|              |       |            |
|--------------|-------|------------|
| 108-23780-25 | 25 mm | 22 cm 8.5" |
| 108-23780-35 | 35 mm | 22 cm 8.5" |
| 108-23780-45 | 45 mm | 22 cm 8.5" |



## Perman STILLE bukchake | Perman STILLE Abdominal Retractor

Self-retaining abdominal retractor. Can be custom designed in various sizes.

|           |             |                     |
|-----------|-------------|---------------------|
| 108-26865 | Depth 65 mm | Arm length 18 cm 7" |
|-----------|-------------|---------------------|





## Doyen STILLE bukhake | Doyen STILLE Abdominal Retractor

The famous STILLE-design abdominal retractor with a variety of blade sizes.  
The handle is rotated 90° for maximum comfort.

|              |           |             |
|--------------|-----------|-------------|
| 108-24930-20 | 90x20mm   | 22 cm 8.75" |
| 108-24930-30 | 90x30mm   | 22 cm 8.75" |
| 108-24935-35 | 110x35mm  | 22 cm 8.75" |
| 108-24935-40 | 110x40mm  | 22 cm 8.75" |
| 108-24935-55 | 110x55mm  | 22 cm 8.75" |
| 108-24940-55 | 120x55 mm | 24 cm 9.5"  |



PRIME

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE bröstshake | STILLE Breast Retractor


Useful double-sided retractor for breast augmentation. Serrated blades in two different sizes and shapes. Ergonomic handle for maximum comfort.

108-26038 Double-ended black coated handle 23 cm 9"



## Desmarres STILLE ögonlockshake | Desmarres STILLE Eye-lid Retractor

For lid retraction in blepharoplasty. Available in four different sizes, with or without insulation.



|              |                 |            |
|--------------|-----------------|------------|
| 108-26046-8  | 8 mm            | 14 cm 5.5" |
| 108-26046-10 | 10 mm           | 14 cm 5.5" |
| 108-26046-12 | 12 mm           | 14 cm 5.5" |
| 108-26046-14 | 14 mm           | 14 cm 5.5" |
| 108-26826-8  | Insulated 8 mm  | 14 cm 5.5" |
| 108-26826-10 | Insulated 10 mm | 14 cm 5.5" |
| 108-26826-12 | Insulated 12 mm | 14 cm 5.5" |
| 108-26826-14 | Insulated 14 mm | 14 cm 5.5" |



## Converse STILLE näshake | Converse STILLE Nasal Retractor

Nasal retractor with concave blade. Curved handle for better grip.

108-26015 11.5 cm 4.5"



## Deaver STILLE hake | Deaver STILLE Retractor

Flexible standard retractor for breast- and facial surgery.

|              |       |               |
|--------------|-------|---------------|
| 108-26037-19 | 19 mm | 17.5 cm 6.75" |
| 108-26037-22 | 22 mm | 20 cm 8"      |
| 108-26037-25 | 25 mm | 30 cm 12"     |



PRIME

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## Ragnell STILLE hake | Ragnell STILLE Retractor

Double-ended retractor for plastic and reconstructive surgery.



108-24808

108-24808 Double ended 9x3/15x5 mm 15 cm 6"



## Langenbeck STILLE hake | Langenbeck STILLE Retractor


General operating retractor with a variety of blade sizes.

|              |          |             |
|--------------|----------|-------------|
| 108-24812-7  | 20x7 mm  | 21 cm 8.25" |
| 108-24814-10 | 30x10 mm | 21 cm 8.25" |
| 108-24814-15 | 30x15 mm | 21 cm 8.25" |
| 108-24814-20 | 30x20 mm | 21 cm 8.25" |
| 108-24816-20 | 45x20 mm | 21 cm 8.25" |



## STILLE hake | STILLE Retractor

Retractor with one, two or three prongs, sharp or blunt.



108-24400-1

108-24400-2

108-24400-3

108-24401-2

108-24401-3

108-24402-2

108-24402-3

108-24403-2


108-24403-3

|             |            |               |          |
|-------------|------------|---------------|----------|
| 108-24400-1 | Sharp      | 1 prong       | 15 cm 7" |
| 108-24400-2 | Sharp      | 2 prongs 5 mm | 15 cm 7" |
| 108-24400-3 | Sharp      | 3 prongs 9 mm | 15 cm 7" |
| 108-24401-2 | Ergo Sharp | 2 prongs 5 mm | 15 cm 7" |
| 108-24401-3 | Ergo Sharp | 3 prongs 9 mm | 15 cm 7" |
| 108-24402-2 | Blunt      | 2 prongs 5 mm | 15 cm 7" |
| 108-24402-3 | Blunt      | 3 prongs 9 mm | 15 cm 7" |
| 108-24403-2 | Ergo Blunt | 2 prongs 5 mm | 15 cm 7" |
| 108-24403-3 | Ergo Blunt | 3 prongs 9 mm | 15 cm 7" |



## Volkmann STILLE hake | Volkmann STILLE Retractor

Retractor with one to four prongs. Sharp or blunt.



108-24680-1

108-24680-2

108-24680-3

108-24680-4

108-24682-2




108-24682-3

108-24682-4

|             |       |               |              |
|-------------|-------|---------------|--------------|
| 108-24680-1 | Sharp | 1 prong Large | 21.5 cm 8.5" |
| 108-24680-2 | Sharp | 2 prongs      | 22 cm 8.75"  |
| 108-24680-3 | Sharp | 3 prongs      | 22 cm 8.75"  |
| 108-24680-4 | Sharp | 4 prongs      | 22 cm 8.75"  |
| 108-24682-2 | Blunt | 2 prongs      | 22 cm 8.75"  |
| 108-24682-3 | Blunt | 3 prongs      | 22 cm 8.75"  |
| 108-24682-4 | Blunt | 4 prongs      | 22 cm 8.75"  |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS


(Surgical Instruments, Medical Equipments, Hospital Supplies)





## Kleinert-Kutz STILLE hudhake | Kleinert-Kutz STILLE Skin Hook

Skin hook with a sharp single prong in three different sizes. Popular in facial and nasal surgery.



108-25031-3      108-25031-5      108-25031-7

|             |        |      |             |
|-------------|--------|------|-------------|
| 108-25031-3 | Small  | 3 mm | 16 cm 6.25" |
| 108-25031-5 | Medium | 5 mm | 16 cm 6.25" |
| 108-25031-7 | Large  | 7 mm | 16 cm 6.25" |

## STILLE-AI Hussein lambâhake | STILLE-AI Hussein Flap Retractor


Flap retractor with four or five sharp hooks, for secure handling of larger flaps. Developed in cooperation with Ali Al-Husseiny, M.D.

|           |               |       |
|-----------|---------------|-------|
| 108-24116 | 4 sharp hooks | 10 cm |
| 108-24117 | 5 sharp hooks | 15 cm |



## Gillies STILLE hudhake | Gillies STILLE Skin Hook

Skin hook with small or large hook on long flat shaft.



108-24110      108-24112

|           |       |            |
|-----------|-------|------------|
| 108-24110 | Small | 19 cm 7.5" |
| 108-24112 | Large | 19 cm 7.5" |



## Joseph STILLE hudhake | Joseph STILLE Skin Hook

Skin hook with sharp prongs. Available in five different models. Popular in plastic surgery.




108-25011-1      108-25012-2      108-25012-5      108-25012-10

|              |             |                        |
|--------------|-------------|------------------------|
| 108-25011-1  | Single hook | 16 cm 6.25"            |
| 108-25012-2  | Double hook | 2 mm      16 cm 6.25"  |
| 108-25012-5  | Double hook | 5 mm      16 cm 6.25"  |
| 108-25012-7  | Double hook | 7 mm      16 cm 6.25"  |
| 108-25012-10 | Double hook | 10 mm      16 cm 6.25" |

PRIME SURGICAL

Hakar | Retractors | Wundhaken | Divaricatori |  
Separadores | 開創器、鉤、レトラクター



## STILLE uterusmanipulator | STILLE Uterine Manipulator

Designed to give better control of the uterus for improved visualization during laparoscopic procedures. Slightly rounded probe with a drop-shaped tip for minimized risk of tissue trauma.


107-13290 Probe 60 mm Curved 25.5 cm 10.5"





## Schroeder STILLE klotång | Schroeder STILLE Tenaculum Forceps

Uterine tenaculum forceps.



107-13280

107-13280 24 cm 9.5"



## Foerster STILLE ringtång | Foerster STILLE Sponge Holding Forceps

Sponge holding forceps with several different loops and lengths.



107-12485



107-12487



107-12486



107-12488

|           |          |                |            |
|-----------|----------|----------------|------------|
| 107-12489 | Serrated | Straight       | 18 cm 7"   |
| 107-12490 | Serrated | Curved         | 18 cm 7"   |
| 107-12485 | Serrated | Straight small | 24 cm 9.5" |
| 107-12486 | Serrated | Curved small   | 24 cm 9.5" |
| 107-12487 | Serrated | Straight large | 24 cm 9.5" |
| 107-12488 | Serrated | Curved large   | 24 cm 9.5" |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Maier STILLE tvättång | Maier STILLE Sponge Holding Forceps


Sponge holding forceps with serrated jaws. Available straight and curved.

|           |              |          |              |
|-----------|--------------|----------|--------------|
| 107-12478 | With ratchet | Straight | 26 cm 10.25" |
| 107-12479 | With ratchet | Curved   | 26 cm 10.25" |



## Allis STILLE fattningstång | Allis STILLE Intestinal Grasping Forceps

Atraumatic tissue-grasping forceps.



|              |               |            |
|--------------|---------------|------------|
| 107-13130-15 | Non-traumatic | 15.5 cm 6" |
| 107-13130-19 | Non-traumatic | 19 cm 7.5" |
| 107-13130-25 | Non-traumatic | 25 cm 10"  |



## STILLE-Norell spiraltång | STILLE-Norell IUD Forceps

IUD forceps with slightly curved, flattened jaws.



|           |                 |             |
|-----------|-----------------|-------------|
| 107-12491 | slightly curved | 27 cm 10.5" |
|-----------|-----------------|-------------|



## Allis STILLE fattningstång | Allis STILLE Intestinal Grasping Forceps

Tissue-grasping forceps for secure grip of tissue.



107-13120

|              |           |             |
|--------------|-----------|-------------|
| 107-13120-15 | 5x6 teeth | 15 cm 6"    |
| 107-13120-21 | 5x6 teeth | 21 cm 8.25" |



## Babcock STILLE fattningstång | Babcock STILLE Intestinal Grasping Forceps

Intestinal grasping forceps in two different lengths.




107-12370

|              |             |
|--------------|-------------|
| 107-12370-16 | 16 cm 6.25" |
| 107-12370-20 | 20 cm 8"    |



## Duval STILLE fattningstång | Duval STILLE Intestinal Grasping Forceps

Intestinal grasping forceps.



107-12563

|           |          |
|-----------|----------|
| 107-12563 | 18 cm 7" |
|-----------|----------|



PRIME

## Backhaus STILLE dukklämmare | Backhaus STILLE Towel Clamps

Towel clamps in two different versions. Sharp/sharp or ball/socket for paper towels.

|           |             |               |
|-----------|-------------|---------------|
| 107-10910 | Sharp/sharp | 11 cm 4.25"   |
| 107-10912 | Ball/socket | 13.5 cm 5.25" |



## Ericsson-STILLE profunda arteria femoralis klämmare | Ericsson-STILLE Profunda Femoral Artery Clamps

Profunda femoral artery clamp with atraumatic DeBakey jaws.



107-10898

|           |            |          |
|-----------|------------|----------|
| 107-10898 | Atraumatic | 18 cm 7" |
|-----------|------------|----------|



PRIME

## Ericsson-STILLE carotisklämmare | Ericsson-STILLE Carotid Clamp

Carotid clamp with angled atraumatic DeBakey jaws.



107-10880

|           |            |            |
|-----------|------------|------------|
| 107-10880 | Atraumatic | 17.5 cm 7" |
|-----------|------------|------------|



## STILLE DeBakey aorta aneurysmklämmare | STILLE DeBakey Aortic-Aneurysm Clamp

Aortic aneurysm clamp with atraumatic DeBakey jaws.



107-10846



107-10849

|           |            |           |              |
|-----------|------------|-----------|--------------|
| 107-10846 | Atraumatic | 8 cm jaws | 26 cm 10.25" |
| 107-10849 | Atraumatic | 7 cm jaws | 28 cm 11"    |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Caspar STILLE bittång | Caspar STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft. Available with or without coating. Fenestrated and serrated jaws in three different sizes.

|                |         |          |        |
|----------------|---------|----------|--------|
| 105-13955-18   | 2x12 mm | 18 cm 7" |        |
| 105-13955-18-C | 2x12 mm | 18 cm 7" | Coated |
| 105-13956-18   | 3x12 mm | 18 cm 7" |        |
| 105-13956-18-C | 3x12 mm | 18 cm 7" | Coated |
| 105-13957-18   | 4x14 mm | 18 cm 7" |        |
| 105-13957-18-C | 4x14 mm | 18 cm 7" | Coated |



## Spurling STILLE bittång | Spurling STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft. Available with or without coating. Jaw size 4x10 mm.

|                |         |          |        |
|----------------|---------|----------|--------|
| 105-13941-18   | 4x10 mm | 18 cm 7" |        |
| 105-13941-18-C | 4x10 mm | 18 cm 7" | Coated |



## Cushing STILLE bittång | Cushing STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft. Available with or without coating. Jaw size 2x10 mm.

|                |         |          |        |
|----------------|---------|----------|--------|
| 105-13921-18   | 2x10 mm | 18 cm 7" |        |
| 105-13921-18-C | 2x10 mm | 18 cm 7" | Coated |



## Love-Gruenwald STILLE bittång | Love-Gruenwald STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft. Available with or without coating. Jaw size 3x10 mm.

|                |         |          |        |
|----------------|---------|----------|--------|
| 105-13931-18   | 3x10 mm | 18 cm 7" |        |
| 105-13931-18-C | 3x10 mm | 18 cm 7" | Coated |





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE ståltrådsroterare | STILLE Wire Twisting Forceps

A good tool for twisting steel wires. With TC jaws.

|           |              |               |
|-----------|--------------|---------------|
| 104-15832 | TC - inserts | 15 cm 5"      |
| 104-15835 | TC - inserts | 18.5 cm 7.25" |



## STILLE ståltrådshållare | STILLE Wire Holding Forceps

For holding steel wires. Equipped with double-action joint and TC jaws. The jaws close parallel to secure a firm grip of the wire.

|           |        |          |
|-----------|--------|----------|
| 104-15854 | Fine   | 18 cm 7" |
| 104-15855 | Medium | 18 cm 7" |
| 104-15856 | Large  | 18 cm 7" |



## Rapp-STILLE knuthållartång | Rapp-STILLE Knot-holding Forceps

A 13 cm knot-holding forceps with smooth jaws. Ideal instrument to use as a tool for holding the suture knot in narrow areas.



104-15140

|           |          |
|-----------|----------|
| 104-15140 | 13 cm 5" |
|-----------|----------|



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Kaye STILLE nålförare | Kaye STILLE Needle holder

Delicate needle holder with short, fine TC jaws and sharp scissors. The scissor part enables suturing and cutting with one instrument. Intended for 5/0 suture and finer.



104-15777

|           |          |         |             |
|-----------|----------|---------|-------------|
| 104-15777 | Serrated | TC jaws | 12 cm 4.75" |
|-----------|----------|---------|-------------|



## Converse STILLE nålförare | Converse STILLE Needle holder

Short and delicate needle holder for fine work. With TC jaws. Intended for 4/0 suture and finer.

|           |         |               |
|-----------|---------|---------------|
| 104-15785 | TC jaws | 10.5 cm 4.25" |
|-----------|---------|---------------|



## Olsen-Hegar STILLE nålförare | Olsen-Hegar STILLE Needle holder

Popular needle holder with scissors for suturing and cutting with the same instrument. Intended for 0 suture and finer. See table "Needle holder in relation to suture size" for suitable suture size for each article code. 3+3 ratchets, 104-15590-171 has 2+2 ratchets.

|               |          |               |
|---------------|----------|---------------|
| 104-15590-121 | Scissors | 12 cm 4.75"   |
| 104-15590-131 | Scissors | 13.5 cm 5.25" |
| 104-15590-171 | Scissors | 17 cm 6.75"   |
| 104-15590-181 | Scissors | 18.5 cm 7.25" |



## Mathieu STILLE nålförare | Mathieu STILLE Needle holder

Needle holder with straight TC jaws. Plier type handle with ratchets that allows the surgeon to open and close the needle holder rapidly. Intended for 1 suture and finer.



104-15621

|               |         |             |
|---------------|---------|-------------|
| 104-15621-170 | TC jaws | 17 cm 6.75" |
| 104-15621-209 | TC jaws | 20 cm 8"    |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE nålförare | STILLE Needle holder

TC-sprayed delicate jaws. Intended for 4/0 suture and finer.



104-15585

|               |         |             |
|---------------|---------|-------------|
| 104-15585-161 | DD jaws | 16 cm 6.25" |
| 104-15585-181 | DD jaws | 18 cm 7.1"  |
| 104-15585-221 | DD jaws | 22 cm 8.75" |
| 104-15585-251 | DD jaws | 25 cm 10"   |



## STILLE nålförare | STILLE Needle holder

Delicate needle holder with TC jaws. Intended for 4/0 suture and finer.



104-15586

|           |         |          |
|-----------|---------|----------|
| 104-15586 | TC jaws | 18 cm 7" |
|-----------|---------|----------|



## New Orleans STILLE nålförare | New Orleans STILLE Needle holder

Short tapered TC jaws, one jaw straight serrated and the other cross-serrated. Intended for 3/0 suture and finer.



104-15571

|           |         |          |
|-----------|---------|----------|
| 104-15571 | TC jaws | 18 cm 7" |
|-----------|---------|----------|



## Ryder STILLE nålförare | Ryder STILLE Needle holder

Needle holder with very narrow, non-tapered, serrated TC jaws. Popular in vascular surgery. Intended for 3/0 suture and finer.



104-15573

|               |         |             |
|---------------|---------|-------------|
| 104-15573-161 | TC jaws | 16 cm 6.25" |
| 104-15573-181 | TC jaws | 18 cm 7"    |
| 104-15573-201 | TC jaws | 20 cm 8"    |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Crile-Wood STILLE gracil nålförare | Crile-Wood STILLE Delicate Needle holder

Delicate version of the popular Crile-Wood needle holder. With TC jaws. Intended for 3/0 suture and finer.



104-15566

|               |          |         |          |
|---------------|----------|---------|----------|
| 104-15566-151 | Delicate | TC jaws | 15 cm 6" |
| 104-15566-181 | Delicate | TC jaws | 18 cm 7" |
| 104-15566-201 | Delicate | TC jaws | 20 cm 8" |
| 104-15566-231 | Delicate | TC jaws | 23 cm 9" |



## STILLE-Frimand DUO-GRIP nålförare | STILLE-Frimand DUO-GRIP Needle holder

Designed in cooperation with Carl-Fredrik Frimand, M.D. Intended for precise use in the palmar grip without the need of grip change during suturing. This new and innovative needle holder reduces the suturing time and tissue manipulation as well as offering better ergonomics for the surgeon. Intended for 2/0 suture and finer.



104-15565


|               |         |          |
|---------------|---------|----------|
| 104-15565-151 | TC jaws | 15 cm 6" |
|---------------|---------|----------|



## Mayo-Hegar STILLE kraftig nålförare | Mayo-Hegar STILLE Heavy Needle holder

Heavier version of the widely used Mayo-Hegar needle holder. with TC jaws. Intended for 1 suture and finer.

|               |               |         |               |
|---------------|---------------|---------|---------------|
| 104-15562-151 | Heavy Pattern | TC jaws | 15 cm 6.75"   |
| 104-15562-181 | Heavy Pattern | TC jaws | 18.5 cm 7.25" |
| 104-15562-201 | Heavy Pattern | TC jaws | 20.5 cm 8"    |
| 104-15562-241 | Heavy Pattern | TC jaws | 24 cm 9.5"    |
| 104-15562-271 | Heavy Pattern | TC jaws | 27 cm 10.75"  |
| 104-15562-301 | Heavy Pattern | TC jaws | 30.5 cm 12"   |



## STILLE Crile-Wood nålförare | STILLE Crile-Wood Needle holder

General purpose needle holder with TC jaws. Intended for 2/0 suture and finer.



104-15565

|               |                    |
|---------------|--------------------|
| 104-15565-151 | TC jaws 15 cm 6 in |
| 104-15565-181 | TC jaws 18 cm 7 in |
| 104-15565-201 | TC jaws 20 cm 8 in |
| 104-15565-231 | TC jaws 23 cm 9 in |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## Halsey STILLE nålförare | Halsey STILLE Needle Holder

Needle holder with TC jaws used in many procedures. One model with ergonomically designed finger ring for thumb. For slightly larger needles than the Webster needle holder. Intended for 3/0 suture and finer.



104-15560

|             |                   |            |
|-------------|-------------------|------------|
| 104-15560   | TC jaws           | 12.5 cm 5" |
| 104-15560-2 | Ergonomic TC jaws | 12.5 cm 5" |

## Mayo-Hegar STILLE nålförare | Mayo-Hegar STILLE Needle holder

Commonly used general purpose needle holder with TC jaws. Intended for 0 suture and finer.



104-15561

|               |         |             |
|---------------|---------|-------------|
| 104-15561-141 | TC jaws | 14 cm 5.5"  |
| 104-15561-161 | TC jaws | 16 cm 6.25" |
| 104-15561-181 | TC jaws | 18 cm 7"    |
| 104-15561-201 | TC jaws | 20 cm 8"    |
| 104-15561-251 | TC jaws | 25 cm 10"   |



## Gillies STILLE nålförare | Gillies STILLE Needle Holder

Needle holder with scissors. Ergonomic design on shanks and finger rings. One rather large ring handle is set at an angle for fitting the top of the thumb. Without ratchets. With or without TC jaws. The scissor part enables suturing and cutting with one instrument. Intended for 2/0 suture and finer.



104-15550




104-15551

|           |                       |             |
|-----------|-----------------------|-------------|
| 104-15550 | With scissors         | 16 cm 6.25" |
| 104-15551 | With scissors TC jaws | 16 cm 6.25" |



## Webster/Kleinert STILLE nålförare | Webster/Kleinert STILLE Needle holder

Delicate needle holder for fine needles. With smooth jaws or TC jaws. The Webster-Kleinert (104-15557) has offset bows for a comfortable grip. Intended for 4/0 suture and finer.




104-15556



104-15559

|           |             |         |             |
|-----------|-------------|---------|-------------|
| 104-15556 | Delicate    | Smooth  | 12 cm 4.75" |
| 104-15557 | Offset bows | Smooth  | 12 cm 4.75" |
| 104-15559 | Serrated    | TC jaws | 12 cm 4.75" |



## Webster STILLE nålförare | Webster STILLE Needle holder

Delicate needle holder for fine needles. With smooth TC jaws. Intended for 3/0 suture and finer.



104-15558

|           |                |             |
|-----------|----------------|-------------|
| 104-15558 | Smooth TC jaws | 12 cm 4.75" |
|-----------|----------------|-------------|





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## Derf STILLE nålförare | Derf STILLE Needle holder

Needle holder with rather short TC jaws, popular in plastic surgery. Intended for 3/0 suture and finer.



104-15542

104-15542 TC jaws 12 cm 4.75"



## Webster STILLE nålförare | Webster STILLE Needle holder

Very delicate needle holder for very fine sutures. Smooth jaws or DD jaws. Intended for 5/0 suture and finer.



104-15545

104-15546

104-15545 Very delicate DD jaws 12 cm 4.75"

104-15546 Very delicate Smooth 12 cm 4.75"

## Kilner STILLE nålförare | Kilner STILLE Needle holder

Delicate needleholder with bayonet shaped shanks. TC jaws. Intended for 4/0 suture and finer.

104-15549-131 13 cm 5"




| Art.no        | Description                                              | Max suture size |
|---------------|----------------------------------------------------------|-----------------|
| 104-15542     | RAPP-STILLE KNOT-HOLDING FORCEPS 13 cm 5"                | 3/0             |
| 104-15545     | WEBSTER STILLE very delicate DD jaws 12 cm 4.75"         | 5/0             |
| 104-15546     | WEBSTER STILLE very delicate smooth 12 cm 4.75"          | 5/0             |
| 104-15549-131 | KILNER STILLE TC jaws 13 cm 5"                           | 4/0             |
| 104-15550     | GILLIES STILLE serrated 16 cm 6.25"                      | 2/0             |
| 104-15551     | GILLIES STILLE TC jaws 16 cm 6.25"                       | 2/0             |
| 104-15556     | WEBSTER STILLE delicate smooth 12 cm 4.75"               | 4/0             |
| 104-15557     | WEBSTER - KLEINERT STILLE smooth offset bows 12 cm 4.75" | 4/0             |
| 104-15558     | WEBSTER STILLE smooth TC jaws 12 cm 4.75"                | 3/0             |
| 104-15559     | WEBSTER STILLE serrated TC jaws 12 cm 4.75"              | 4/0             |
| 104-15560     | HALSEY STILLE TC jaws 12.5 cm 5"                         | 3/0             |
| 104-15560-2   | HALSEY STILLE ERGONOMIC TC jaws 12.5 cm 5"               | 3/0             |
| 104-15561-140 | MAYO-HEGAR STILLE TC jaws 14 cm 5.5"                     | 0               |
| 104-15561-161 | MAYO-HEGAR STILLE TC jaws 16 cm 6.25"                    | 0               |
| 104-15561-181 | MAYO-HEGAR STILLE TC jaws 18 cm 7"                       | 0               |
| 104-15561-201 | MAYO-HEGAR STILLE TC jaws 20 cm 8"                       | 0               |
| 104-15561-251 | MAYO-HEGAR STILLE TC jaws 25 cm 10"                      | 0               |
| 104-15562-161 | MAYO-HEGAR STILLE Heavy Pattern TC jaws 16 cm 6.25"      | 1               |
| 104-15562-181 | MAYO-HEGAR STILLE Heavy Pattern TC jaws 18 cm 7"         | 1               |
| 104-15562-201 | MAYO-HEGAR STILLE Heavy Pattern TC jaws 20 cm 8"         | 1               |
| 104-15562-241 | MAYO-HEGAR STILLE Heavy Pattern TC jaws 24 cm 9.5"       | 1               |
| 104-15562-271 | MAYO-HEGAR STILLE Heavy Pattern TC jaws 27 cm 10.5"      | 1               |
| 104-15562-301 | MAYO-HEGAR STILLE Heavy Pattern TC jaws 30 cm 12"        | 1               |
| 104-15565-151 | CRILE-WOOD STILLE TC jaws 15 cm 6"                       | 2/0             |
| 104-15565-181 | CRILE-WOOD STILLE TC jaws 18 cm 7"                       | 2/0             |
| 104-15565-201 | CRILE-WOOD STILLE TC jaws 20 cm 8"                       | 2/0             |
| 104-15565-231 | CRILE-WOOD STILLE TC jaws 23 cm 9"                       | 2/0             |

| Art.no        | Description                                            | Max suture size |
|---------------|--------------------------------------------------------|-----------------|
| 104-15566-151 | BABY CRILE-WOOD STILLE delicate TC 15 cm 6"            | 3/0             |
| 104-15566-181 | CRILE-WOOD STILLE delicate TC jaws 18 cm 7"            | 3/0             |
| 104-15566-201 | CRILE-WOOD STILLE delicate TC jaws 20 cm 8"            | 3/0             |
| 104-15566-231 | CRILE-WOOD STILLE delicate TC jaws 23 cm 9"            | 3/0             |
| 104-15571     | NEW ORLEANS STILLE TC jaws 18 cm 7"                    | 3/0             |
| 104-15573-161 | RYDER STILLE TC jaws 16 cm 6.25"                       | 3/0             |
| 104-15573-181 | RYDER STILLE TC jaws 18 cm 7"                          | 3/0             |
| 104-15573-201 | RYDER STILLE TC jaws 20 cm 8"                          | 3/0             |
| 104-15585-161 | STILLE DO jaws 16 cm 6.25"                             | 4/0             |
| 104-15585-181 | STILLE DO jaws 18.5 cm 7.5"                            | 4/0             |
| 104-15585-221 | STILLE DO jaws 22 cm 8.75"                             | 4/0             |
| 104-15585-251 | STILLE DO jaws 25 cm 10"                               | 4/0             |
| 104-15586     | STILLE TC jaws 18 cm 7"                                | 4/0             |
| 104-15590-121 | OLSEN-HEGAR STILLE TC jaws with scissors 12 cm 4.75"   | 3/0             |
| 104-15590-131 | OLSEN-HEGAR STILLE TC jaws with scissors 13.5 cm 5.25" | 2/0             |
| 104-15590-171 | OLSEN-HEGAR STILLE TC jaws with scissors 17 cm 6.75"   | 0               |
| 104-15590-181 | OLSEN-HEGAR STILLE TC jaws with scissors 18.5 cm 7.25" | 0               |
| 104-15621-179 | MATHIEU STILLE TC jaws 17 cm 6.75"                     | 1               |
| 104-15621-209 | MATHIEU STILLE TC jaws 20 cm 8"                        | 1               |
| 104-15665-151 | Frimand-STILLE DUO-GRIP Needle holder TC jaws 15 cm 6" | 2/0             |
| 104-15777     | KAYE STILLE with scissors serrated TC jaws 12 cm 4.75" | 5/0             |
| 104-15785     | CONVERSE STILLE TC jaws 10.5 cm 4.25"                  | 4/0             |

Suture size indicates the largest recommended for each needle holder.




## Nålförare | Needle holders | Nadelhalter | Portaghi | Porta-agujas | 持針器



### STILLE artär- och ligaturklämmare | STILLE Artery and Ligature Clamp

Four very delicate models of dissecting forceps with different jaw curvature.



|             |                 |               |
|-------------|-----------------|---------------|
| 103-10423-1 | Slightly curved | 20 cm 8"      |
| 103-10423-2 | Curved          | 20 cm 8"      |
| 103-10423-3 | Strongly curved | 19.5 cm 7.75" |
| 103-10423-4 | Angled 90°      | 19 cm 7.5"    |




### STILLE-Ehrén ligaturklämmare | STILLE-Ehrén Ligature Clamp

For fine and precise handling of ligatures in very narrow areas, such as paediatric surgery. The jaws close only at the tips to minimize the risk for accidental clamping of fragile vessels etc. Developed together with Henrik Ehrén, M.D. at Karolinska University Hospital. Matted finish.



## Crafoord STILLE peang | Crafoord STILLE Artery Forceps

Designed for thoracic and cardiovascular surgery. The fine tips make this a very useful instrument. Offers the surgeon Mosquito type tips, but with a length that allows working in a deep cavity.



|               |        |            |
|---------------|--------|------------|
| 103-10416     | Curved | 18 cm 7"   |
| 103-10430-241 | Curved | 24 cm 9.5" |



## Heiss STILLE peang | Heiss STILLE Artery Forceps

Dissecting and ligature clamps with half-serrated, jaws.




|               |                 |             |
|---------------|-----------------|-------------|
| 103-10421-211 | Straight        | 21 cm 8.25" |
| 103-10422-211 | Curved          | 21 cm 8.25" |
| 103-10420-201 | Strongly curved | 20 cm 8"    |



## Overholt STILLE peang | Overholt STILLE Dissecting Forceps

Strongly curved dissecting forceps.



|           |                 |            |
|-----------|-----------------|------------|
| 103-10425 | Strongly curved | 14 cm 5.5" |
|-----------|-----------------|------------|



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Dandy STILLE peang | Dandy STILLE Artery Forceps

Side curved hemostat without teeth. Often used in neuro surgery.



103-10310



## Schmidt STILLE peang | Schmidt STILLE Artery Forceps

A hemostat popular in general surgery.



103-10411-201



## Olivecrona STILLE peang | Olivecrona STILLE Artery Forceps

16 cm straight hemostat. Useful as a long hemostat with fine points.



103-10240-161




## Adson STILLE peang | Adson STILLE Artery Forceps

Popular model when a little longer and fairly delicate hemostat is needed.



103-10260-191




# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Rochester-Ochsner STILLE peang | Rochester-Ochsner STILLE Artery Forceps

Sturdy hemostat with 1x2 teeth. Used when a strong and long forceps is needed, both to occlude blood vessels and to lift or retract organs.



|               |           |          |            |
|---------------|-----------|----------|------------|
| 103-10215-201 | 1x2 teeth | Straight | 20 cm 8"   |
| 103-10215-241 | 1x2 teeth | Straight | 24 cm 9.5" |
| 103-10217-201 | 1x2 teeth | Curved   | 20 cm 8"   |
| 103-10217-241 | 1x2 teeth | Curved   | 24 cm 9.5" |



## Bengolea STILLE peang | Bengolea STILLE Artery Forceps

A long hemostat with straight or curved jaws without teeth.




|               |          |          |
|---------------|----------|----------|
| 103-10220-201 | Straight | 20 cm 8" |
| 103-10222-201 | Curved   | 20 cm 8" |



## Rochester-Ochsner STILLE peang | Rochester-Ochsner STILLE Artery Forceps

Sturdy hemostat with 1x2 teeth. Used when a strong forceps is needed, both to occlude blood vessels and to lift or retract organs.




|               |           |          |             |
|---------------|-----------|----------|-------------|
| 103-10205-161 | 1x2 teeth | Straight | 16 cm 6.25" |
| 103-10205-181 | 1x2 teeth | Straight | 18 cm 7"    |
| 103-10207-161 | 1x2 teeth | Curved   | 16 cm 6.25" |
| 103-10207-181 | 1x2 teeth | Curved   | 18 cm 7"    |



## Rochester-Péan STILLE peang | Rochester-Péan STILLE Artery Forceps

Sturdy hemostats. The 20 cm has short and the 24 cm has long jaws. Also referred to as hysterectomy forceps.



|               |          |            |
|---------------|----------|------------|
| 103-10210-201 | Straight | 20 cm 8"   |
| 103-10210-241 | Straight | 24 cm 9.5" |
| 103-10212-201 | Curved   | 20 cm 8"   |
| 103-10212-241 | Curved   | 24 cm 9.5" |
| 103-10212-301 | Curved   | 30 cm 12"  |




# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Leriche STILLE peang | Leriche STILLE Artery Forceps

Fairly tapered hemostat. For use in small operative fields.




|           |           |          |          |
|-----------|-----------|----------|----------|
| 103-10190 |           | Straight | 15 cm 6" |
| 103-10192 |           | Curved   | 15 cm 6" |
| 103-10195 | 1x2 teeth | Straight | 15 cm 6" |
| 103-10197 | 1x2 teeth | Curved   | 15 cm 6" |



## Rochester-Péan STILLE peang | Rochester-Péan STILLE Artery Forceps

Sturdy hemostat. Relatively blunt tips. Used in all types of surgery.




|               |  |          |             |
|---------------|--|----------|-------------|
| 103-10200-161 |  | Straight | 16 cm 6.25" |
| 103-10200-181 |  | Straight | 18 cm 7"    |
| 103-10202-161 |  | Curved   | 16 cm 6.25" |
| 103-10202-181 |  | Curved   | 18 cm 7"    |



## Crile STILLE peang | Crile STILLE Artery Forceps

General purpose hemostat without teeth. Used in relatively small operating fields, when more strength is needed than a Mosquito can offer.



|               |  |          |             |
|---------------|--|----------|-------------|
| 103-10180-141 |  | Straight | 14 cm 5.5"  |
| 103-10180-161 |  | Straight | 16 cm 6.25" |
| 103-10180-181 |  | Straight | 18 cm 7"    |
| 103-10182-141 |  | Curved   | 14 cm 5.5"  |
| 103-10182-161 |  | Curved   | 16 cm 6.25" |
| 103-10182-181 |  | Curved   | 18 cm 7"    |



## Crile STILLE peang | Crile STILLE Artery Forceps

General purpose hemostat with 1x2 teeth. Used in relatively small operating fields, when more strength is needed than a Mosquito can offer.



|               |           |          |             |
|---------------|-----------|----------|-------------|
| 103-10185     | 1x2 teeth | Straight | 14 cm 5.5"  |
| 103-10185-161 | 1x2 teeth | Straight | 16 cm 6.25" |
| 103-10187-141 | 1x2 teeth | Curved   | 14 cm 5.5"  |
| 103-10187-161 | 1x2 teeth | Curved   | 16 cm 6.25" |





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



Kelly-Rankin STILLE peang | Kelly-Rankin STILLE Artery Forceps

Longer version of the Kelly hemostat with jaws half-serrated.

|               |          |             |
|---------------|----------|-------------|
| 103-10160-161 | Straight | 16 cm 6.25" |
| 103-10162-161 | Curved   | 16 cm 6.25" |



Christophe-STILLE peang | Christophe-STILLE Artery Forceps

Curved hemostat without teeth.

|           |        |             |
|-----------|--------|-------------|
| 103-10172 | Curved | 16 cm 6.25" |
|-----------|--------|-------------|

Kocher STILLE peang | Kocher STILLE Artery Forceps

Hemostat with 1x2 teeth. Widely used when a short instrument is needed. Can also be used as tissue forceps to retract or lift an organ.



|               |           |          |            |
|---------------|-----------|----------|------------|
| 103-10115-141 | 1x2 teeth | Straight | 14 cm 5.5" |
| 103-10117-141 | 1x2 teeth | Curved   | 14 cm 5.5" |



Kelly STILLE peang | Kelly STILLE Artery Forceps

Fairly short hemostat with jaws, half-serrated.

|               |          |            |
|---------------|----------|------------|
| 103-10150-141 | Straight | 14 cm 5.5" |
| 103-10152-141 | Curved   | 14 cm 5.5" |





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



**Pean STILLE peang | Pean STILLE Artery Forceps**

General utility forceps used in many surgical disciplines.

103-10015-121 Straight 13.5 cm 5.25"



**Halsted-Mosquito STILLE peang | Halsted-Mosquito STILLE Artery Forceps**

Hemostat without teeth, used in small operative fields. For plastic surgery, neuro, hand or any delicate surgery.



103-10030-121 Straight 12 cm 4.75"

103-10032-121 Curved 12 cm 4.75"

103-10032-181 Curved 18 cm 7"



**Halsted-Mosquito STILLE gracil peang | Halsted-Mosquito STILLE delicate Artery Forceps**

Delicate hemostat without teeth for small operating field.



103-10026-121 Straight 12 cm 4.75"

103-10028-121 Curved 12 cm 4.75"



**Halsted-Mosquito STILLE peang | Halsted-Mosquito STILLE Artery Forceps**


Hemostat with 1x2 teeth, used in small operative fields. For plastic surgery, neuro, hand or any delicate surgery.



103-10035-121 1x2 teeth Straight 12 cm 4.75"


103-10037-121 1x2 teeth Curved 12 cm 4.75"

## Peanger | Artery forceps | Arterienklemmen | Pinze emostatiche | Pinzas hemostáticas | 止血鉗子



## Björk-STILLE knutpincett | Björk-STILLE Knot-Holding Forceps

Same as Björk-STILLE diathermy forceps, but without insulation.



102-9846-20

102-9846-20 20 cm 8 in



## STILLE De Bakley kärlpincett | STILLE De Bakley Vascular Forceps

Atraumatic forceps with De Bakley jaws for handling of blood vessels and other delicate tissue. Can be insulated on request.

| 102-9300 102-9302 | 102-9312 102-9314 | 102-9320 102-9322 |
|-------------------|-------------------|-------------------|
| 102-9304          | 102-9316 102-9318 | 102-9322          |

|          |                |        |              |
|----------|----------------|--------|--------------|
| 102-9340 | De Bakley type | 1 mm   | 15 cm 6"     |
| 102-9342 | De Bakley type | 1 mm   | 20 cm 8"     |
| 102-9354 | De Bakley type | 1.5 mm | 25 cm 10"    |
| 102-9300 | De Bakley type | 2 mm   | 15 cm 6"     |
| 102-9302 | De Bakley type | 2 mm   | 20 cm 8"     |
| 102-9304 | De Bakley type | 2 mm   | 25 cm 10"    |
| 102-9312 | De Bakley type | 2.7 mm | 15 cm 6"     |
| 102-9314 | De Bakley type | 2.7 mm | 20 cm 8"     |
| 102-9316 | De Bakley type | 2.7 mm | 25 cm 10"    |
| 102-9318 | De Bakley type | 2.7 mm | 30 cm 11.75" |
| 102-9320 | De Bakley type | 3.2 mm | 20 cm 8"     |
| 102-9322 | De Bakley type | 3.2 mm | 25 cm 10"    |



## STILLE De Bakley kärlpincett | STILLE De Bakley Vascular Forceps

Angled atraumatic forceps with De Bakley jaws for delicate tissue.



102-9308

102-9308 De Bakley type 2 mm 20 cm 8"

## STILLE De Bakley kärlpincett | STILLE De Bakley Vascular Forceps

Atraumatic forceps with De Bakley jaws for delicate tissue, non serrated handles. Stop pin for preventing from gripping the tissue too hard.



102-9330

102-9330 De Bakley type 2 mm 20 cm 8"

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE irispincett | STILLE Iris Forceps

Short delicate forceps with either 1x2 fine teeth or serrated jaws.



|          |           |          |
|----------|-----------|----------|
| 102-9460 | Serrated  | 10 cm 4" |
| 102-9470 | 1x2 Teeth | 10 cm 4" |



## Olivecrona-STILLE durapincett | Olivecrona-STILLE Dura Forceps

Dura forceps with 2x3 fine teeth. The nonpenetrating tips give the surgeon maximum grip without damaging the tissue.



|             |                |           |
|-------------|----------------|-----------|
| 102-9244-15 | 2x3 fine teeth | 15 cm 6"  |
| 102-9244-20 | 2x3 fine teeth | 20 cm 8"  |
| 102-9244-25 | 2x3 fine teeth | 25 cm 10" |



## Crafoord-STILLE anatomisk pincett | Crafoord-STILLE Dressing Forceps

Dressing forceps with rounded points. The perfect balance and delicate grip gives this instrument a unique feeling.



|             |              |          |
|-------------|--------------|----------|
| 102-9077-20 | Serrated     | 20 cm 8" |
| 102-9078-20 | Diamond jaws | 20 cm 8" |



## Waugh-Brophy STILLE anatomisk pincett | Waugh-Brophy STILLE Dressing Forceps

Dressing forceps with very fine points and two different jaw surfaces. The cross-serration gives excellent holding power.



|             |                |           |
|-------------|----------------|-----------|
| 102-9254    | Cross-serrated | 20 cm 8"  |
| 102-9255-20 | Diamond jaws   | 20 cm 8"  |
| 102-9255-25 | Diamond jaws   | 25 cm 10" |



NATIONAL

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## Adson STILLE anatomisk pincett | Adson STILLE Dressing Forceps

General-purpose dressing forceps with various types of jaw surfaces.



|             |                        |             |
|-------------|------------------------|-------------|
| 102-9020    | Serrated               | 12 cm 4.75" |
| 102-9020-15 | Serrated               | 15 cm 6"    |
| 102-9020-1  | Diamond jaws           | 12 cm 4.75" |
| 102-9024    | TC jaws cross-serrated | 12 cm 4.75" |

## STILLE anatomisk pincett | STILLE Dressing Forceps

Dressing forceps with square shaped and serrated points used in all surgical specialties.



|             |               |             |
|-------------|---------------|-------------|
| 102-9030-13 | Square points | 13 cm 5"    |
| 102-9030-15 | Square points | 15 cm 6"    |
| 102-9030-17 | Square points | 17 cm 6.75" |
| 102-9030-20 | Square points | 20 cm 8"    |
| 102-9030-25 | Square points | 25 cm 10"   |



## McIndoe STILLE anatomisk pincett | McIndoe STILLE Dressing Forceps


Dressing forceps with stabilizing pin and delicate jaws for grasping fine tissue.



|            |              |               |
|------------|--------------|---------------|
| 102-9025   | Serrated     | 15.5 cm 6.25" |
| 102-9025-1 | Diamond Jaws | 15.5 cm 6.25" |

## STILLE anatomisk pincett | STILLE Dressing Forceps

Dressing forceps with rounded points and serrated jaws.



|             |                |             |
|-------------|----------------|-------------|
| 102-9040-15 | Rounded points | 15 cm 6"    |
| 102-9040-17 | Rounded points | 17 cm 6.75" |
| 102-9040-20 | Rounded points | 20 cm 8"    |
| 102-9040-25 | Rounded points | 25 cm 10"   |
| 102-9040-30 | Rounded points | 30 cm 12"   |



## Adlercreutz STILLE pincett | Adlercreutz STILLE Tissue Forceps

Tissue forceps with 4x5 fine teeth. Very useful in neurosurgery to hold the dura, and in cardiovascular surgery to hold the pericardium.



|             |                |           |
|-------------|----------------|-----------|
| 102-9248-15 | 4x5 fine teeth | 15 cm 6"  |
| 102-9248-20 | 4x5 fine teeth | 20 cm 8"  |
| 102-9248-25 | 4x5 fine teeth | 25 cm 10" |

## Waugh-Brophy STILLE pincett | Waugh-Brophy STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth. Very fine tips, like Adson. Useful when long and slender instruments are needed.



|             |           |                |           |
|-------------|-----------|----------------|-----------|
| 102-9252-20 | 1x2 teeth | Cross-serrated | 20 cm 8"  |
| 102-9252-25 | 1x2 teeth | Cross-serrated | 25 cm 10" |
| 102-9252-30 | 1x2 teeth | Cross-serrated | 30 cm 12" |



## Russian STILLE pincett | Russian STILLE Tissue Forceps

Forceps for secure grip of tissue. Used when a better grip than with a dressing forceps is needed, and when damage to the tissue must be avoided. For hard and slippery tissue.



|             |          |
|-------------|----------|
| 102-9270-20 | 20 cm 8" |
|-------------|----------|



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE pincett | STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth and stabilizing pin. Delicate model.



102-9028 Serrated 1x2 fine teeth 15.5 cm 6.25"



## STILLE pincett | STILLE Tissue Forceps

Tissue forceps with 1x2 teeth.



102-9060

102-9060-15 1x2 teeth 15 cm 6"

102-9060-17 1x2 teeth 17 cm 6.75"

102-9060-20 1x2 teeth 20 cm 8"



## STILLE pincett | STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth.



102-9050-13

102-9050-13 1x2 fine teeth 13 cm 5"



## STILLE pincett | STILLE Tissue Forceps

Tissue forceps with 2x3 teeth. Used when more grasping power than 1x2 teeth is needed.



102-9062

102-9062-15 2x3 teeth 15 cm 6"

102-9062-17 2x3 teeth 17 cm 6.75"

102-9062-20 2x3 teeth 20 cm 8"



## Adson STILLE pincett | Adson STILLE Tissue Forceps

Tissue forceps with very delicate points and 1x2 fine teeth. Used in most specialties for delicate work.



102-9019



102-9022  
102-9022-15

102-9019 Very delicate 1x2 teeth 12 cm 4.75"

102-9022 Cross-serrated 1x2 teeth 12 cm 4.75"

102-9022-15 Cross-serrated 1x2 teeth 15 cm 6"



## Adson-Callison STILLE pincett | Adson-Callison STILLE Tissue Forceps

Tissue forceps with tying platform and 1x2 fine teeth.



102-9021



102-9021-1

102-9021 Tying platform 1x2 teeth 12 cm 4.75"

102-9021-1 Delicate tying platform 1x2 teeth 12 cm 4.75"



## Adson-Brown STILLE pincett | Adson-Brown STILLE Tissue Forceps

Tissue forceps with side-grasping teeth for secure, atraumatic grasping of soft tissue or cartilage.



102-9023-12

102-9023-12 Side-grasping teeth 12 cm 4.75"



## Gillies STILLE pincett | Gillies STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth and stabilizing pin. An excellent forceps for the surgeon who wants an Adson-size tip, but with greater length.




102-9027

102-9027 Serrated 1x2 fine teeth 15.5 cm 6.25"





## Pincetter | Forceps | Pinzetten | Pinze | Pinzas | 鑷子 (ピンセット)



Littauer STILLE sutursax |  
Littauer STILLE Suture scissors

A somewhat larger model than Spencer. Used in many procedures.

101-8597-14 Delicate pattern 14 cm 5.5"



Universal STILLE ståltrådsax |  
Universal STILLE Wire Cutting scissors

Wire cutting scissors with TC-inserts. For soft wire up to 1.2 mm,  
for hard wire up to 0.8 mm.

101-8815-12 TC 12 cm 4.75"



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## STILLE sutursax | STILLE Suture scissors

Suture removal scissors.



101-8584

101-8584 Curved 13.5 cm 5.25"



## Spencer STILLE sutursax | Spencer STILLE Suture scissors

Popular suture scissors for many procedures. Delicate model.

101-8586-9 9 cm 3.5"



## Lister STILLE TC bandagesax | Lister STILLE TC Bandage Scissors

Bandage scissors with TC-inserts for improved durability.

101-8595-18 TC-inserts 18 cm 7"

## Littler STILLE trådföringssax | Littler STILLE Suture carrying scissors

Suture carrying scissors. Rounded blades.



101-8164

101-8164 Curved 12.5 cm 5"



PRIME

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



**Rees STILLE SuperCut XE Face-lift sax |  
Rees STILLE SuperCut XE Face-lift scissors**

SuperCut XE face-lift scissors, double-beveled with semi-sharp outer blades for undermining skin flaps.



101-8456-171

101-8456-171 Curved 17 cm 6.75"



**Matarasso STILLE SuperCut XE Face-lift sax |  
Matarasso STILLE SuperCut XE Face-lift scissors**

SuperCut XE face-lift scissors with double-beveled semi-sharp outer blades.



101-8458

101-8459-20

101-8458-17 Curved 17 cm 6.75"

101-8458-20 Curved 20 cm 8"

101-8459-20 Straight 20 cm 8"



**Matarasso STILLE SuperCut XE ergonomisk  
Face-lift sax | Matarasso STILLE SuperCut XE  
Ergonomic Face-lift scissors**

Ergonomic version of the Matarasso SuperCut XE face-lift scissors with double-beveled semi-sharp outer blades.



101-8458


101-8458-172 Curved 17 cm 6.75"

101-8458-202 Curved 20 cm 8"



**Aston STILLE SuperCut XE Face-lift sax |  
Aston STILLE SuperCut XE Face-lift scissors**

Double-beveled SuperCut XE face-lift scissors with semi-sharp outer blades for undermining flaps.



101-8440-201

101-8441-201

101-8442-171

101-8440-201 Straight 20 cm 8 in

101-8441-201 Curved 20 cm 8 in

101-8442-171 Curved Serrated 17 cm 6.75 in



**Davis STILLE SuperCut XE Face-lift sax |  
Davis STILLE SuperCut XE Face-lift Scissors**

SuperCut XE face-lift scissors with tapered blades and sharp tips for precise cutting and dissection.



101-8452


101-8452 Curved 18 cm 7"

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Gorney STILLE SuperCut XE Face-lift sax | Gorney STILLE SuperCut XE Face-lift scissors

Super Cut double-beveled face-lift scissors with flat tips and semi-sharp outer blades. Two models with serrated blades to prevent tissue slippage.




|              |          |             |
|--------------|----------|-------------|
| 101-8430-171 | Straight | 17 cm 6.75" |
| 101-8431-171 | Curved   | 17 cm 6.75" |
| 101-8433-20  | Curved   | 20 cm 8"    |
| 101-8435-20  | Curved   | 20 cm 8"    |
| 101-8436-20  | Straight | 20 cm 8"    |



## Gorney STILLE SuperCut XE ergonomisk Face-lift sax | Gorney STILLE SuperCut XE Ergonomic Face-lift scissors

Ergonomic design to avoid fatigue in the hand. One model with serrated blades to prevent tissue slippage. Flat tips and semi-sharp outer blades.



|              |          |             |
|--------------|----------|-------------|
| 101-8430-172 | Straight | 17 cm 6.75" |
| 101-8431-172 | Curved   | 17 cm 6.75" |
| 101-8433-202 | Curved   | 20 cm 8"    |
| 101-8435-202 | Curved   | 20 cm 8"    |



## STILLE Cottle SuperCut XE nässax | STILLE Cottle SuperCut XE Nasal Scissors

SuperCut XE scissors with heavy blades for cutting dorsal cartilage. Angled shanks for better visibility. Serrated blades.



|          |                 |               |             |
|----------|-----------------|---------------|-------------|
| 101-8604 | Serrated blades | Angled shanks | 16 cm 6.25" |
|----------|-----------------|---------------|-------------|



## STILLE SuperCut XE nässax | STILLE SuperCut XE Nasal Scissors

SuperCut XE nasal scissors for excision of cartilage.



|          |          |               |             |
|----------|----------|---------------|-------------|
| 101-8411 | Straight | Angled shanks | 16 cm 6.25" |
|----------|----------|---------------|-------------|



## STILLE Fomon SuperCut XE nässax | STILLE Fomon SuperCut XE Nasal Scissors

Serrated SuperCut XE scissors for cutting dorsal cartilage. Angled shanks for better visibility.



|          |                 |               |            |
|----------|-----------------|---------------|------------|
| 101-8425 | Serrated blades | Angled shanks | 14 cm 5.5" |
|----------|-----------------|---------------|------------|




# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Peck-Joseph STILLE nässax | Peck-Joseph STILLE Nasal Scissors

Double-beveled scissors for rhinoplasty.




|          |       |          |          |
|----------|-------|----------|----------|
| 101-8120 | Sharp | Straight | 15 cm 6" |
| 101-8122 | Bunt  | Straight | 15 cm 6" |
| 101-8123 | Sharp | Curved   | 15 cm 6" |
| 101-8125 | Bunt  | Curved   | 15 cm 6" |



## Peck-Joseph STILLE SuperCut XE nässax | Peck-Joseph STILLE SuperCut XE Nasal Scissors

Double-beveled scissors for rhinoplasty.



|            |       |          |          |
|------------|-------|----------|----------|
| 101-8120-1 | Sharp | Straight | 15 cm 6" |
| 101-8123-1 | Sharp | Curved   | 15 cm 6" |



## Fomon STILLE nässax | Fomon STILLE Nasal Scissors

Scissors for rhinoplasty. Double-beveled blades.



|          |        |            |
|----------|--------|------------|
| 101-8025 | Curved | 12.5 cm 5" |
|----------|--------|------------|



## Joseph STILLE nässax | Joseph STILLE Nasal Scissors

Nasal scissors for the lower lateral cartilage. Useful in dissecting and cutting tough tissue. Beveled blades.



|          |          |          |
|----------|----------|----------|
| 101-8110 | Straight | 15 cm 6" |
| 101-8115 | Curved   | 15 cm 6" |



## Joseph STILLE SuperCut XE nässax | Joseph STILLE SuperCut XE Nasal Scissors

Super Cut nasal scissors for the lower lateral cartilage. Useful in dissecting and cutting tough tissue. Beveled blades.



|            |          |          |
|------------|----------|----------|
| 101-8110-1 | Straight | 15 cm 6" |
| 101-8115-1 | Curved   | 15 cm 6" |



PRIME

# PRIME Surgical-GYNECOLOGY & OBSTETRICS


(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE-Lund kärlsax |

### STILLE-Lund Vascular Scissors

Vascular scissors with one sharp and one probe-pointed blade for precise handling of fine blood vessels. Popular model in pathology procedures.


|            |                       |          |          |
|------------|-----------------------|----------|----------|
| 101-8490   | 1 probe-pointed blade | Straight | 15 cm 6" |
| 101-8490-2 | 1 probe-pointed blade | Ergo     | 15 cm 6" |



## STILLE SuperCut XE kärlsax |

### STILLE SuperCut XE Vascular Scissors

Super Cut scissors with angled blades for better reach of vessels etc. Used in vascular surgery.




|              |             |             |
|--------------|-------------|-------------|
| 101-8200-25  | Angled 25°  | 18 cm 7"    |
| 101-8200-45  | Angled 45°  | 18 cm 7"    |
| 101-8200-60  | Angled 60°  | 18 cm 7"    |
| 101-8200-90  | Angled 90°  | 17 cm 6.75" |
| 101-8200-125 | Angled 125° | 17 cm 6.75" |



## STILLE SuperCut XE kärlsax |

### STILLE SuperCut XE Vascular Scissors

Extra delicate version of the popular STILLE Vascular Scissor with shorter and thinner blades.



|              |             |             |
|--------------|-------------|-------------|
| 101-8201-25  | Angled 25°  | 18 cm 7"    |
| 101-8201-45  | Angled 45°  | 18 cm 7"    |
| 101-8201-60  | Angled 60°  | 18 cm 7"    |
| 101-8201-90  | Angled 90°  | 17 cm 6.75" |
| 101-8201-125 | Angled 125° | 17 cm 6.75" |





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE-Ekberg Stevens SuperCut ergonomisk tenotomisax | STILLE-Ekberg Stevens Ergonomic SuperCut Tenotomy Scissors

Ergonomic STILLE-Ekberg version. Available with blades curved up or down, suitable for both left- and right-handed surgeons.




|               |             |          |
|---------------|-------------|----------|
| 101-8167-153D | Curved down | 15 cm 6" |
| 101-8167-153U | Curved up   | 15 cm 6" |



## STILLE-Ekberg Stevens Diamond SuperCut ergonomisk tenotomisax | STILLE-Ekberg Stevens Diamond Ergonomic SuperCut Tenotomy Scissors

Diamond version for increased durability.




|                 |             |          |
|-----------------|-------------|----------|
| 101-8167-153D-C | Curved down | 15 cm 6" |
| 101-8167-153U-C | Curved up   | 15 cm 6" |



## Stevens-STILLE SuperCut XE ergonomisk tenotomisax | Stevens-STILLE SuperCut XE Ergonomic Tenotomy Scissors

Ergonomic design of the popular Stevens SuperCut XE tenotomy scissors.




|              |          |            |
|--------------|----------|------------|
| 101-8166-2   | Straight | 12.5 cm 5" |
| 101-8167-2   | Curved   | 12.5 cm 5" |
| 101-8167-152 | Curved   | 15 cm 6"   |
| 101-8167-182 | Curved   | 18 cm 7"   |



## Stevens-STILLE Diamond SuperCut ergonomisk tenotomisax | Stevens-STILLE Diamond Ergonomic SuperCut Tenotomy Scissors

Diamond version for increased durability.



|              |        |            |
|--------------|--------|------------|
| 101-8167-2-C | Curved | 12.5 cm 5" |
|--------------|--------|------------|





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Stevens-STILLE Diamond SuperCut tenotomisax | Stevens-STILLE Diamond SuperCut Tenotomy Scissors

Diamond version for increased durability.

|                                                                                   |                                                                                   |                                                                                   |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
|  |  |  |
| 101-8167-1-C                                                                      | 101-8167-151-C                                                                    | 101-8167-181-C                                                                    |
| 101-8167-1-C                                                                      | Curved                                                                            | 12.5 cm 5"                                                                        |
| 101-8167-151-C                                                                    | Curved                                                                            | 15 cm 6"                                                                          |
| 101-8167-181-C                                                                    | Curved                                                                            | 18 cm 7"                                                                          |



## Stevens-STILLE tenotomisax | Stevens-STILLE Tenotomy Scissors

Tenotomy scissors with rounded, very thin blades and tapered tips. Used in plastic surgery and hand surgery for cutting with extreme precision in small incisions.

|                                                                                   |                                                                                   |            |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------|
|  |  |            |
| 101-8166                                                                          | 101-8167                                                                          |            |
| 101-8166                                                                          | Straight                                                                          | 12.5 cm 5" |
| 101-8167                                                                          | Curved                                                                            | 12.5 cm 5" |



## Stevens-STILLE SuperCut tenotomisax | Stevens-STILLE SuperCut Tenotomy Scissors


Very popular SuperCut XE scissors for precise cutting and dissection of fine tissues. Rounded blades.

|                                                                                   |                                                                                    |              |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|--------------|
|  |  |              |
| 101-8166-1                                                                        | 101-8167-1                                                                         |              |
|  |  |              |
| 101-8167-151                                                                      | 101-8170-1                                                                         |              |
|  |                                                                                    |              |
| 101-8171-1                                                                        |                                                                                    |              |
| 101-8166-1                                                                        | Straight                                                                           | 12.5 cm 5"   |
| 101-8167-1                                                                        | Curved                                                                             | 12.5 cm 5"   |
| 101-8167-151                                                                      | Curved                                                                             | 15 cm 6"     |
| 101-8167-181                                                                      | Curved                                                                             | 18 cm 7"     |
| 101-8170-1                                                                        | Straight                                                                           | 9.5 cm 3.75" |
| 101-8171-1                                                                        | Curved                                                                             | 9.5 cm 3.75" |



## Ragnell STILLE operationssax | Ragnell STILLE Operating Scissors

Mammary plastic scissors intended for breast surgery. Rounded blades.

|                                                                                     |        |             |
|-------------------------------------------------------------------------------------|--------|-------------|
|  |        |             |
| 101-8460                                                                            |        |             |
| 101-8460                                                                            | Angled | 17 cm 6.75" |



## König-STILLE nageluppklippringsax | König-STILLE Nail-splitting Scissors

Nail-splitting scissors, for cutting up ingrown nails. Beveled blades.

|                                                                                     |          |          |
|-------------------------------------------------------------------------------------|----------|----------|
|  |          |          |
| 101-8580                                                                            |          |          |
| 101-8580                                                                            | Straight | 15 cm 6" |




# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Mayo-Harrington STILLE SuperCut XE operationsax | Mayo-Harrington STILLE SuperCut XE Operating Scissors

Operating Super Cut scissors useful for several surgical disciplines. Beveled blades.



|             |          |           |
|-------------|----------|-----------|
| 101-8097-20 | Straight | 20 cm 8"  |
| 101-8097-23 | Straight | 23 cm 9"  |
| 101-8099-20 | Curved   | 20 cm 8"  |
| 101-8099-23 | Curved   | 23 cm 9"  |
| 101-8099-30 | Curved   | 30 cm 12" |



## Boyd-STILLE tonsillsax | Boyd-STILLE Tonsil Scissors

Scissors with sharp tips and S-shaped shanks for greater visibility. Rounded blades.




101-8420 Curved S-shaped shanks 17 cm 6.75"



## STILLE Cooper operationsax | STILLE Cooper Operating Scissors

General purpose utility scissors, often used for sutures, gauze etc. Beveled blades.




|             |             |          |               |
|-------------|-------------|----------|---------------|
| 101-8042-15 | Sharp/blunt | Straight | 15 cm 6"      |
| 101-8043-15 | Blunt/blunt | Straight | 15 cm 6"      |
| 101-8047-13 | Sharp/blunt | Curved   | 13.5 cm 5.25" |
| 101-8047-15 | Sharp/blunt | Curved   | 15 cm 6"      |
| 101-8048-15 | Blunt/blunt | Curved   | 15 cm 6"      |
| 101-8048-18 | Blunt/blunt | Curved   | 18 cm 7"      |




## Mayo-Harrington STILLE operationsax | Mayo-Harrington STILLE Operating Scissors

Very useful operating scissors for deep surgery like gynecology and thoracic surgery. Beveled blades.



|             |          |           |
|-------------|----------|-----------|
| 101-8096-20 | Straight | 20 cm 8"  |
| 101-8096-23 | Straight | 23 cm 9"  |
| 101-8098-20 | Curved   | 20 cm 8"  |
| 101-8098-23 | Curved   | 23 cm 9"  |
| 101-8098-30 | Curved   | 30 cm 12" |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Ragnell STILLE SuperCut XE ergonomisk dissektionssax | Ragnell STILLE SuperCut XE Ergonomic Dissecting Scissors

Ergo Super Cut versions of Ragnell scissors for precise undermining of flaps and comfortable grip.



|              |        |            |
|--------------|--------|------------|
| 101-8380-122 | Curved | 12.5 cm 5" |
| 101-8380-152 | Curved | 15 cm 6"   |



## Thorek STILLE dissektionssax | Thorek STILLE Dissecting Scissors

Strongly curved dissecting scissors for good visibility. Rounded blades.



|          |                 |             |
|----------|-----------------|-------------|
| 101-8470 | Strongly curved | 21 cm 8.25" |
|----------|-----------------|-------------|



## Ragnell STILLE dissektionssax | Ragnell STILLE Dissecting Scissors

Dissecting scissors with flat, rounded tips. Often used in plastic surgery. Is also called Killner by some manufacturers. The short models are used for undermining in plastic surgery and the long models are used for further dissection.



|             |          |            |
|-------------|----------|------------|
| 101-8378    | Straight | 12.5 cm 5" |
| 101-8380-12 | Curved   | 12.5 cm 5" |
| 101-8380-15 | Curved   | 15 cm 6"   |
| 101-8380-18 | Curved   | 18 cm 7"   |
| 101-8380-20 | Curved   | 20 cm 8"   |
| 101-8380-23 | Curved   | 23 cm 9"   |



## Ragnell STILLE SuperCut XE dissektionssax | Ragnell STILLE SuperCut XE Dissecting Scissors

Super Cut version of the Ragnell scissors. For precise undermining of flaps.



|              |        |            |
|--------------|--------|------------|
| 101-8380-121 | Curved | 12.5 cm 5" |
| 101-8380-151 | Curved | 15 cm 6"   |
| 101-8380-181 | Curved | 18 cm 7"   |
| 101-8380-201 | Curved | 20 cm 8"   |
| 101-8380-231 | Curved | 23 cm 9"   |




# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Metzenbaum STILLE SuperCut XE dissektionssax | Metzenbaum STILLE SuperCut XE Dissecting Scissors

The popular Metzenbaum dissecting scissors with Super Cut edges. Rounded/beveled blades.



|             |          |              |
|-------------|----------|--------------|
| 101-8175-15 | Curved   | 15 cm 6"     |
| 101-8175-18 | Curved   | 18 cm 7"     |
| 101-8175-20 | Curved   | 20 cm 8"     |
| 101-8175-23 | Curved   | 23 cm 9"     |
| 101-8175-26 | Curved   | 26 cm 10.25" |
| 101-8175-30 | Curved   | 30 cm 12"    |
| 101-8175-33 | Curved   | 33 cm 13"    |
| 101-8176-15 | Straight | 15 cm 6"     |
| 101-8176-18 | Straight | 18 cm 7"     |
| 101-8176-20 | Straight | 20 cm 8"     |
| 101-8176-23 | Straight | 23 cm 9"     |
| 101-8176-26 | Straight | 26 cm 10.25" |



## Metzenbaum STILLE Diamond SuperCut dissektionssax | Metzenbaum STILLE Diamond SuperCut Dissecting Scissors

Diamond version of the popular Metzenbaum STILLE SuperCut scissors, for increased durability.



|               |        |              |
|---------------|--------|--------------|
| 101-8175-15-C | Curved | 15 cm 6"     |
| 101-8175-18-C | Curved | 18 cm 7"     |
| 101-8175-20-C | Curved | 20 cm 8"     |
| 101-8175-23-C | Curved | 23 cm 9"     |
| 101-8175-26-C | Curved | 26 cm 10.25" |
| 101-8175-30-C | Curved | 30 cm 12"    |



## Metzenbaum STILLE gracil Diamond SuperCut dissektionssax | Metzenbaum STILLE Delicate Diamond SuperCut Dissecting Scissors

With Diamond coating for increased durability.



|               |                 |          |
|---------------|-----------------|----------|
| 101-8173-15-C | Curved delicate | 15 cm 6" |
| 101-8173-18-C | Curved delicate | 18 cm 7" |



## Metzenbaum STILLE dissektionssax | Metzenbaum STILLE Dissecting Scissors

Dissecting scissors with rounded blades and tips. Used in all specialties for fine dissection. Metzenbaum is one of STILLE's most popular scissors model and is available in a variety of shapes and sizes.



|                                                    |        |              |
|----------------------------------------------------|--------|--------------|
| 101-8174-15, 101-8174-18, 101-8174-20, 101-8174-23 |        |              |
| 101-8174-15                                        | Curved | 15 cm 6"     |
| 101-8174-18                                        | Curved | 18 cm 7"     |
| 101-8174-20                                        | Curved | 20 cm 8"     |
| 101-8174-23                                        | Curved | 23 cm 9"     |
| 101-8174-26                                        | Curved | 26 cm 10.25" |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Metzenbaum STILLE gracil dissektionssax | Metzenbaum STILLE Delicate Dissecting Scissors

Short, delicate dissecting scissors with rounded blades. Used in delicate surgery requiring small instruments.



101-8172

101-8172-18 Curved delicate 18 cm 7"  
101-8172-23 Curved delicate 23 cm 9"



## Metzenbaum STILLE gracil SuperCut XE dissektionssax | Metzenbaum STILLE Delicate SuperCut XE Dissecting Scissors

Very delicate version with Super Cut edges and fine rounded tips.



101-8173-15



101-8173-18  
101-8173-20

101-8173-15 Curved delicate 15 cm 6"  
101-8173-18 Curved delicate 18 cm 7"  
101-8173-20 Curved delicate 20 cm 8"



## Baby Metzenbaum STILLE SuperCut XE dissektionssax | Baby Metzenbaum STILLE SuperCut XE Dissecting Scissors

Super Cut dissecting scissors for tough tissue. Used in delicate surgery requiring small instruments.



101-8161



101-8163

101-8161 Straight 12.5 cm 5"  
101-8163 Curved 12.5 cm 5"



## Baby Metzenbaum STILLE Diamond SuperCut dissektionssax | Baby Metzenbaum STILLE Diamond SuperCut Dissecting Scissors

With Diamond coating for increased durability.



101-8163-C

101-8163-C Curved 12.5 cm 5"



## Baby Metzenbaum STILLE SuperCut XE ergono- misk dissektionssax | Baby Metzenbaum STILLE SuperCut XE Ergonomic Dissecting Scissors

Ergonomic model for better control. Rounded blades.



101-8163-2

101-8163-2 Curved 12.5 cm 5"





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



**Sims STILLE dissektionssax |  
Sims STILLE Dissecting Scissors**

Mainly used in gynecological and other pelvic surgery. Beveled blades.




101-8138

101-8138 Curved 20 cm 8"



**Baby Metzenbaum STILLE dissektionssax |  
Baby Metzenbaum STILLE Dissecting Scissors**

Short, delicate dissecting scissors with rounded blades. Used in delicate surgery requiring small instruments.



101-8150


101-8152

101-8150 Straight 12.5 cm 5"  
101-8152 Curved 12.5 cm 5"



**STILLE-Mayo SuperCut XE dissektionssax |  
STILLE-Mayo SuperCut XE Dissecting Scissors**

Sturdy Super Cut scissors with rounded blades for dissection of tough tissue.



101-8085-15

101-8087-15


101-8085-17

101-8087-17

101-8085-19


101-8087-19

|               |                |          |                     |
|---------------|----------------|----------|---------------------|
| 101-8085-15   | Rounded blades | Straight | 15 cm 6"            |
| 101-8085-17   | Rounded blades | Straight | 17 cm 6.75"         |
| 101-8085-17-C | Rounded blades | Straight | 17 cm 6.75" Diamond |
| 101-8085-19   | Rounded blades | Straight | 19 cm 7.5"          |
| 101-8087-15   | Rounded blades | Curved   | 15 cm 6"            |
| 101-8087-17   | Rounded blades | Curved   | 17 cm 6.75"         |
| 101-8087-17-C | Rounded blades | Curved   | 17 cm 6.75" Diamond |
| 101-8087-19   | Rounded blades | Curved   | 19 cm 7.5"          |



**Mayo-Noble STILLE dissektionssax |  
Mayo-Noble STILLE Dissecting Scissors**

Short and strong dissecting scissors mainly for orthopaedic surgeons. Beveled blades.



101-8090

101-8092

101-8090 Straight 18 cm 7"  
101-8092 Curved 18 cm 7"

**Baby Mayo STILLE dissektionssax |  
Baby Mayo STILLE Dissecting Scissors**

Delicate version of the regular Mayo scissors. Mainly for plastic- and hand surgery. Beveled blades.



101-8100-15

101-8102-15


101-8100-15 Straight 15 cm 6"  
101-8102-15 Curved 15 cm 6"






# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## Mayo STILLE dissektionssax | Mayo STILLE Dissecting Scissors

The most used scissors, for all types of surgery. Used for routine dissection and cutting in tissue. Beveled blades.



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Kaye STILLE dissektionssax | Kaye STILLE Dissecting Scissors

Scissors with serrated blades for "holding" the tissue, making it easier to cut exactly as far as intended.



|          |                |        |            |
|----------|----------------|--------|------------|
| 101-8020 | Blunt serrated | Curved | 10 cm 4"   |
| 101-8021 | Blunt serrated | Curved | 12.5 cm 5" |
| 101-8022 | Sharp serrated | Curved | 12.5 cm 5" |



## Kaye STILLE SuperCut XE dissektionssax | Kaye STILLE SuperCut XE Dissecting Scissors

SuperCut XE version for improved cutting performance.



|            |                |        |            |
|------------|----------------|--------|------------|
| 101-8021-1 | Blunt serrated | Curved | 12.5 cm 5" |
|------------|----------------|--------|------------|



## Jabaley STILLE dissektionssax | Jabaley STILLE Dissecting Scissors

Sharp scissors suitable for hand surgery, where very strong scissors are used. Beveled blades.



|          |          |            |
|----------|----------|------------|
| 101-8060 | Straight | 12.5 cm 5" |
| 101-8062 | Curved   | 12.5 cm 5" |



## STILLE SuperCut XE ögonlockssax | STILLE SuperCut XE Blepharoplasty Scissors

Short delicate SuperCut XE scissors suitable for blepharoplasty.



|            |          |              |
|------------|----------|--------------|
| 101-8018-1 | Straight | 9.5 cm 3.75" |
| 101-8019-1 | Curved   | 9.5 cm 3.75" |



## STILLE SuperCut XE ergonomisk ögonlockssax | STILLE SuperCut XE Ergonomic Blepharoplasty Scissors

Blepharoplasty SuperCut XE scissors with ergonomically designed shanks for a more relaxed grip.



|            |        |              |
|------------|--------|--------------|
| 101-8019-2 | Curved | 9.5 cm 3.75" |
|------------|--------|--------------|



## Ragnell STILLE SuperCut XE ögonlockssax | Ragnell STILLE SuperCut XE Blepharoplasty Scissors

Short delicate SuperCut XE scissors with blunt, flat tips for blepharoplasty.



|            |        |              |
|------------|--------|--------------|
| 101-8031-1 | Curved | 9.5 cm 3.75" |
|------------|--------|--------------|



## Ragnell STILLE SuperCut XE ergonomisk ögonlockssax | Ragnell STILLE SuperCut XE Ergonomic Blepharoplasty Scissors

Blepharoplasty SuperCut XE scissors with ergonomically designed shanks for a more relaxed grip.



|            |        |              |
|------------|--------|--------------|
| 101-8031-2 | Curved | 9.5 cm 3.75" |
|------------|--------|--------------|





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STILLE ögonsax | STILLE Eye Scissors


Delicate eye scissors with sharp or blunt tips and large finger rings. These are also called iris scissors and are used for nerves, tendons and vessels. They are intended for facial surgery and any type of delicate surgery. Beveled blades.



|          |          |       |              |
|----------|----------|-------|--------------|
| 101-8001 | Straight | Sharp | 9.5 cm 3.75" |
| 101-8002 | Straight | Blunt | 9.5 cm 3.75" |
| 101-8003 | Curved   | Sharp | 9.5 cm 3.75" |
| 101-8004 | Curved   | Blunt | 9.5 cm 3.75" |

## STILLE ögonsax | STILLE Eye Scissors


Eye scissors for delicate work. Large finger rings for a comfortable grip. Beveled blades.



|          |          |       |          |
|----------|----------|-------|----------|
| 101-8010 | Straight | Sharp | 10 cm 4" |
| 101-8011 | Straight | Blunt | 10 cm 4" |
| 101-8013 | Curved   | Sharp | 10 cm 4" |
| 101-8015 | Curved   | Blunt | 10 cm 4" |

## STILLE SuperCut XE ögonsax | STILLE SuperCut XE Eye Scissors


SuperCut XE eye scissors with sharp or blunt tips and extremely sharp edges. Beveled blades.



|              |          |       |                  |
|--------------|----------|-------|------------------|
| 101-8010-1   | Straight | Sharp | 10 cm 4"         |
| 101-8011-1   | Straight | Blunt | 10 cm 4"         |
| 101-8015-1   | Curved   | Sharp | 10 cm 4"         |
| 101-8015-1-C | Curved   | Sharp | 10 cm 4" Diamond |
| 101-8016-1   | Curved   | Blunt | 10 cm 4"         |



## Saxar | Scissors | Scheren | Forbici | Tijeras | 剪刀 (はさみ)



## OB/GYN

Miscellaneous

### BETHEA SHEET HOLDERS



214770

| CATALOG # | DESCRIPTION                             |
|-----------|-----------------------------------------|
| 214770    | 1" (2.5 cm) wide, 5" (12.5 cm) diameter |

### IOWA TRUMPET NEEDLE GUIDE




884227

For paracervical anesthesia.


| CATALOG # | DESCRIPTION                   |
|-----------|-------------------------------|
| 884227    | Ball tip 7 mm, 5-1/2" (14 cm) |

## Obstetrics and Gynecology Instruments - Forceps and Scissors

| Catalog # | Description                                                                                         |
|-----------|-----------------------------------------------------------------------------------------------------|
| MH30-1700 | Heaney Hysterectomy Forceps, 8½" (213 mm), heavy pattern, single tooth, curved                      |
| MH30-1710 | Heaney Hysterectomy Forceps, 8½" (215 mm), heavy pattern, double tooth, curved                      |
| MH30-1716 | Heaney-Ballentine Hysterectomy Forceps, 8½" (216 mm), longitudinal serrations, single tooth, curved |
| MH30-1718 | Ballentine Hysterectomy Forceps, 8¾" (220 mm), longitudinal serrations, single tooth, straight      |
| MH30-2190 | Braun Episiotomy Scissors, 5½" (140 mm), angled to side, guarded lower blade                        |
| MH30-2195 | Braun Episiotomy Scissors, 8½" (216 mm), angled to side, guarded lower blade                        |



MH30-1700, MH30-1710



MH30-1716, MH30-1718



MH30-2190, MH30-2195



## Obstetrics and Gynecology Instruments - Biopsy Forceps

| Catalog #  | Description                                                                                                   |
|------------|---------------------------------------------------------------------------------------------------------------|
| MH30-1482  | Kevorkian-Young Uterine Biopsy Forceps, 9 1/2" (241 mm) shaft, 3.5 x 8 mm bite, teeth on lower jaw            |
| MH30-1485  | Eppendorfer Uterine Biopsy Forceps, 8 1/2" (216 mm) shaft                                                     |
| MH30-1405  | Thoms-Gaylor Uterine Biopsy Forceps, 8 1/2" (216 mm), angled shanks, curved jaws, 5 mm bite, interfitted cups |
| MH301442WL | Tischler Biopsy Forceps, 8" (203 mm), with lock                                                               |
| MH301443WL | Burke Baby Tischler Biopsy Forceps, 7 3/4" (197 mm), with lock                                                |



MH30-1482



MH30-1485



MH30-1405



## Obstetrics and Gynecology Instruments - Biopsy Curettes

| Catalog # | Description                                                                                          |
|-----------|------------------------------------------------------------------------------------------------------|
| MH30-1382 | Kevorkian-Younger Endocervical Biopsy Curette, 12 1/4" (310 mm), loop size 3 x 12 mm, with basket    |
| MH30-1383 | Kevorkian-Younger Endocervical Biopsy Curette, 12 1/4" (310 mm), loop size 3 x 12 mm, without basket |
| MH30-1385 | Heaney Uterine Biopsy Curette, 9" (229 mm), serrated loop                                            |



MH30-1382



MH30-1383



MH30-1385

## Obstetrics and Gynecology Instruments

| Catalog #   | Description                                                                             |
|-------------|-----------------------------------------------------------------------------------------|
| MH30-997    | Teale Uterine Vulsellum Forceps, 9" (230 mm), curved sideways, 3 x 4 teeth, 5 mm wide   |
| MH30-1000   | Schroeder Uterine Vulsellum Forceps, 10" (254 mm), straight                             |
| MH30-1025   | Jacobs Uterine Vulsellum Forceps, 8 1/4" (210 mm), straight, serrated jaws, 2 x 2 teeth |
| MH30-1205-0 | Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 0             |
| MH30-1205-1 | Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 1             |
| MH30-1205-2 | Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 2             |
| MH30-1205-3 | Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 3             |
| MH30-1205-4 | Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 4             |
| MH30-1205-5 | Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 5             |
| MH30-1205-6 | Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 6             |




MH30-997



MH30-1000



MH30-1025



0

1

2

3

4

5


6




MH30-1205-0 to MH30-1205-6

## Obstetrics and Gynecology Instruments - Tenaculum Forceps

| Catalog #   | Description                                                                                      |
|-------------|--------------------------------------------------------------------------------------------------|
| MH30-965    | Schroeder-Braun Uterine Tenaculum Forceps, 9 3/4" (248 mm), straight                             |
| MH30-965ATR | Braun Tenaculum Forceps, 10" (254 mm), straight/squared jaw with non-overlapping atraumatic tips |
| MH30-966    | Schroeder Uterine Tenaculum Forceps, 10" (254 mm), round jaw                                     |
| MH30-966ATR | Schroeder Tenaculum Forceps, 10" (254 mm), rounded jaw with non-overlapping atraumatic tips      |
| MH30-967    | Hulka Tenaculum Forceps, 11" (279 mm), with 3.5 mm probe                                         |
| MH30-985    | Barrett Uterine Tenaculum Forceps, 7 1/4" (184 mm), straight                                     |




MH30-965



MH30-966



MH30-967



MH30-985

## Obstetrics and Gynecology Instruments - Dilators

| Catalog #     | Description                                                                                                                                     |
|---------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| MH30-561      | Pratt Uterine Dilators, 11½" (292 mm), double end, set of 8 Dilators in Khaki Roll, sizes 13-15 Fr. (4.3 to 4.6 mm) to 41-43 Fr. (13.5-14.2 mm) |
| MH30-560-1315 | Pratt Uterine Dilator, 11½" (292 mm), double end, 13-15 Fr. (4.3-4.6 mm)                                                                        |
| MH30-560-1719 | Pratt Uterine Dilator, 11½" (292 mm), double end, 17-19 Fr. (5.6-6.3 mm)                                                                        |
| MH30-560-2123 | Pratt Uterine Dilator, 11½" (292 mm), double end, 21-23 Fr. (7-7.6 mm)                                                                          |
| MH30-560-2527 | Pratt Uterine Dilator, 11½" (292 mm), double end, 25-27 Fr. (8.3-9 mm)                                                                          |
| MH30-560-2931 | Pratt Uterine Dilator, 11½" (292 mm), double end, 29-31 Fr. (9.6-10.2 mm)                                                                       |
| MH30-560-3335 | Pratt Uterine Dilator, 11½" (292 mm), double end, 33-35 Fr. (10.9-11.6 mm)                                                                      |
| MH30-560-3739 | Pratt Uterine Dilator, 11½" (292 mm), double end, 37-39 Fr. (12.2-12.9 mm)                                                                      |
| MH30-560-4143 | Pratt Uterine Dilator, 11½" (292 mm), double end, 41-43 Fr. (13.5-14.2 mm)                                                                      |



MH30-561



MH30-560-1315 to MH30-560-4143

## Obstetrics and Gynecology Instruments - Dilators

| Catalog #     | Description                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------|
| MH30-550      | Hegar Uterine Dilators, 7 $\frac{3}{4}$ " (197 mm), double end, set of 8 Dilators in Khaki Roll, sizes 3-4 mm to 17-18 mm |
| MH30-555-34   | Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 3-4 mm                                                     |
| MH30-555-56   | Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 5-6 mm                                                     |
| MH30-555-78   | Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 7-8 mm                                                     |
| MH30-555-910  | Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 9-10 mm                                                    |
| MH30-555-1112 | Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 11-12 mm                                                   |
| MH30-555-1314 | Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 13-14 mm                                                   |
| MH30-555-1516 | Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 15-16 mm                                                   |
| MH30-555-1718 | Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 17-18 mm                                                   |



MH30-550



MH30-555-34 to MH30-555-1718

## Obstetrics and Gynecology Instruments - Dilators

| Catalog #     | Description                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| MH30-500      | Hank Uterine Dilators, 10½" (267 mm), double end, set of 6 Dilators in Khaki Roll, Sizes 11-12, (3.6-4 mm), 13-14 (4.3-4.6 mm), 15-16 (5-5.3 mm), 17-18 (5.6-6 mm), 19-20 (6.3-6.6 mm), and 21-22 (7-7.3 mm) French |
| MH30-505-910  | Hank Uterine Dilator, 10½" (267 mm), double end, 9-10 Fr. (3-3.3 mm)                                                                                                                                                |
| MH30-505-1112 | Hank Uterine Dilator, 10½" (267 mm), double end, 11-12 Fr. (3.6-4 mm)                                                                                                                                               |
| MH30-505-1314 | Hank Uterine Dilator, 10½" (267 mm), double end, 13-14 Fr. (4.3-4.6 mm)                                                                                                                                             |
| MH30-505-1516 | Hank Uterine Dilator, 10½" (267 mm), double end, 15-16 Fr. (5-5.3 mm)                                                                                                                                               |
| MH30-505-1718 | Hank Uterine Dilator, 10½" (267 mm), double end, 17-18 Fr. (5.6-6 mm)                                                                                                                                               |
| MH30-505-1920 | Hank Uterine Dilator, 10½" (267 mm), double end, 19-20 Fr. (6.3-6.6 mm)                                                                                                                                             |



MH30-500



MH30-505-910 to MH30-505-1920



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Obstetrics and Gynecology Instruments - Specula

| Catalog # | Description                                                                                   |
|-----------|-----------------------------------------------------------------------------------------------|
| MH30-10   | Graves Vaginal Speculum, original model, small size, ¾" (19 mm) x 3" (76 mm)                  |
| MH30-15   | Graves Vaginal Speculum, original model, medium size, 1⅜" (35 mm) x 4" (102 mm)               |
| MH30-20   | Graves Vaginal Speculum, original model, large size, 1½" (38 mm) x 4½" (114 mm)               |
| MH30-30   | Graves Open Sided Vaginal Speculum, medium size, 1⅜" (35 mm) x 4" (102 mm), wide angle blades |
| MH30-50   | Pederson Vaginal Speculum, narrow blades, small size, ½" (13 mm) x 3" (76 mm)                 |
| MH30-55   | Pederson Vaginal Speculum, narrow blades, medium size, ⅞" (22 mm) x 4" (102 mm)               |
| MH30-60   | Pederson Vaginal Speculum, narrow blades, large size, 1" (25 mm) x 4¾" (121 mm)               |
| MH30-62   | Pederson Vaginal Speculum, extra narrow blades, ⅝" (16 mm) x 4" (102 mm)                      |



MH30-10 to MH30-20



MH30-50 to MH30-62

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Suction Tubes

| Catalog # | Description                                                                    |
|-----------|--------------------------------------------------------------------------------|
| MH2-104SS | Yankauer Suction Tube, 11½" (292 mm), with removable tip and tubing connection |
| MH2-109   | Yankauer Pediatric Suction Tube, 8½" (216 mm), with removable tip, delicate    |
| MH2-110SS | Andrews-Pynchon Suction Tube 10" (254 mm), delicate                            |
| MH10-310  | Poole Suction Tube, 24 French, (8 mm), curved                                  |
| MH10-312  | Poole Suction Tube, 30 French, (10 mm), straight                               |



MH2-104SS



MH2-110SS



MH10-310



MH10-312

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## Neurological Instruments

| Catalog # | Description                                    |
|-----------|------------------------------------------------|
| MH1-226   | Wartenberg Neurological Pin Wheel, 7" (178 mm) |



MH1-226

## BRUSHES

OB/GYN

### Sets


#### 10-1426 INSTRUMENT CLEANING BRUSH SET

| PRODUCT NUMBER | DESCRIPTION                           | QUANTITIES |
|----------------|---------------------------------------|------------|
| 10-1660        | Cannula Cleaning Brush 12" 3mm, 3/pk  | 1 pk       |
| 10-1362        | Cannula Cleaning Brush 18" 5mm, 3/pk  | 1 pk       |
| 10-1374        | Cannula Cleaning Brush 24" 7mm, 3/pk  | 1 pk       |
| 10-1376        | Cannula Cleaning Brush 24" 10mm, 3/pk | 1 pk       |
| 10-1396        | Suction Tube Brush 6" 1.6mm, 3/pk     | 1 pk       |
| 10-1400        | Suction Tube Brush 8" 2mm, 3/pk       | 1 pk       |
| 10-1402        | Suction Tube Brush 8" 2.6mm, 3/pk     | 1 pk       |
| 10-1404        | Suction Tube Brush 8" 3.3mm, 3/pk     | 1 pk       |
| 10-1406        | Suction Tube Brush 8" 4mm, 3/pk       | 1 pk       |
| 10-1418        | Surface Cleaning Brush, 3/pk          | 1 pk       |
| 10-1444        | Double End Brush Nylon, 3/pk          | 2 pks      |
| 10-1650        | Brush Stainless Steel Bristle, 3/pk   | 1 pk       |
| 10-1651        | Brush Nylon Bristle, 3/pk             | 2 pks      |



#### 10-1428 ENDOSCOPIC BRUSH SET

| PRODUCT NUMBER | DESCRIPTION                            | QUANTITIES |
|----------------|----------------------------------------|------------|
| 10-1350        | Cannula Cleaning Brush 12" 2mm, 3/pk   | 1 pk       |
| 10-1660        | Cannula Cleaning Brush 12" 3mm, 3/pk   | 1 pk       |
| 10-1354        | Cannula Cleaning Brush 12" 5mm, 3/pk   | 1 pk       |
| 10-1362        | Cannula Cleaning Brush 18" 5mm, 3/pk   | 1 pk       |
| 10-1364        | Cannula Cleaning Brush 18" 10mm, 3/pk  | 1 pk       |
| 10-1370        | Cannula Cleaning Brush 24" 2.5mm, 3/pk | 1 pk       |
| 10-1372        | Cannula Cleaning Brush 24" 5mm, 3/pk   | 1 pk       |
| 10-1376        | Cannula Cleaning Brush 24" 10mm, 3/pk  | 1 pk       |
| 10-1398        | Double End Valve Brush Nylon, 3/pk     | 2 pks      |
| 10-1444        | Double End Brush Nylon, 3/pk           | 2 pks      |



## OB/GYN

## BRUSHES

### General Cleaning and Accessories

#### SURFACE CLEANING BRUSH

FOR USE ON SINKS, AUTOCLAVES, ULTRASONIC CLEANERS, BOWLS, BASINS, ETC.

| PRODUCT NUMBER | FEATURES                                                                                 | QUANTITIES |
|----------------|------------------------------------------------------------------------------------------|------------|
| 10-1418        | 3 ¼" L x 3" W Brush Head<br>4 ½" Plastic handle<br>2" White medical grade Nylon bristles | 3/pk       |



#### LARGE INSTRUMENT CLEANING BRUSH

FOR USE ON RETRACTORS, SPECULA, ETC.

| PRODUCT NUMBER | FEATURES                                                                             | QUANTITIES |
|----------------|--------------------------------------------------------------------------------------|------------|
| 10-1420        | 5" L x 1" W Brush Head<br>5" Plastic handle<br>1" White medical grade Nylon bristles | 3/pk       |



#### AUTOCLAVE CLEANING BRUSH

FOR USE ON AUTOCLAVE CHAMBERS AND DOORS, SINKS, CARTS, ETC.

| PRODUCT NUMBER | FEATURES                                                                                | QUANTITIES |
|----------------|-----------------------------------------------------------------------------------------|------------|
| 10-1422        | 3 ¼" L x 3" W Brush Head<br>16" Plastic handle<br>2" White medical grade Nylon bristles | 3/pk       |



#### INSTRUMENT SPOT REMOVER PAD


FOR USE ON O.R. GRADE SATIN FINISH INSTRUMENTS

| PRODUCT NUMBER | FEATURES | QUANTITIES |
|----------------|----------|------------|
| 10-1414        | 6" x 9"  | 1 Each     |
| 10-1415        | 6" x 9"  | 10/pk      |



#### BRUSH DRYING RACK & BRUSH HOLDER

| PRODUCT NUMBER | FEATURES     | QUANTITIES |
|----------------|--------------|------------|
| 10-1399        | Brush Rack   | 1 Each     |
| 10-1545        | Brush Holder | 1 Each     |





## BRUSHES

OB/GYN

### Instrument Cleaning


#### SUCTION TUBE CLEANING BRUSHES (BARON & FRAZIER)

| PRODUCT NUMBER | FEATURES | QUANTITIES |
|----------------|----------|------------|
| 10-1396        | 6" 1.6mm | 3/pk       |
| 10-1400        | 8" 2mm   | 3/pk       |
| 10-1401        | 8" 2.3mm | 3/pk       |
| 10-1402        | 8" 2.6mm | 3/pk       |
| 10-1403        | 8" 3mm   | 3/pk       |
| 10-1404        | 8" 3.3mm | 3/pk       |
| 10-1405        | 8" 3.6mm | 3/pk       |
| 10-1406        | 8" 4mm   | 3/pk       |



#### SURGICAL SCRUB BRUSHES

| PRODUCT NUMBER | FEATURES           | QUANTITIES |
|----------------|--------------------|------------|
| 06-2880        | With Channel Back  | 1 Each     |
| 06-2870        | With Channel Back  | 12/pk      |
| 06-2882        | With Nail Scrubber | 1 Each     |
| 06-2873        | With Nail Scrubber | 12/pk      |
| 06-2881        | Standard           | 1 Each     |
| 06-2872        | Standard           | 12/pk      |



#### CANNULA INSTRUMENT (PIPE) CLEANERS

| PRODUCT NUMBER | FEATURES                 | QUANTITIES |
|----------------|--------------------------|------------|
| 10-1470        | 12" Length 1/8" Diameter | 3/pk       |
| 10-1472        | 12" Length 1/4" Diameter | 3/pk       |



#### DOUBLE END VALVE BRUSHES

| PRODUCT NUMBER | FEATURES                                             | QUANTITIES |
|----------------|------------------------------------------------------|------------|
| 10-1398        | 5 3/4" Length<br>9 mm and 5mm Diameter<br>Nylon Tips | 3/pk       |





## OB/GYN

## BRUSHES

### Instrument Cleaning

#### TRACHEAL TUBE BRUSHES

| PRODUCT NUMBER | FEATURES                   | QUANTITIES |
|----------------|----------------------------|------------|
| 78-2821        | Small, 6" x 1/2" Diameter  | 6/pk       |
| 78-2822        | Medium, 7" x 5/8" Diameter | 6/pk       |



#### INSTRUMENT CLEANING BRUSHES (STAINLESS BRISTLES)

| PRODUCT NUMBER | FEATURES                 | QUANTITIES |
|----------------|--------------------------|------------|
| 10-1650        | Stainless Steel Bristles | 3/pk       |
| 10-1652        | Stainless Steel Bristles | 12/pk      |
| 10-1658        | Stainless Steel Bristles | 24/pk      |



#### INSTRUMENT CLEANING BRUSHES (NYLON BRISTLES)

| PRODUCT NUMBER | FEATURES                   | QUANTITIES |
|----------------|----------------------------|------------|
| 10-1651        | Nylon Bristles             | 3/pk       |
| 10-1657        | Nylon Bristles             | 12/pk      |
| 10-1659        | Nylon Bristles             | 24/pk      |
| 10-1465        | Nylon Bristles, Disposable | 12/pk      |
| 10-1468        | Nylon Bristles, Disposable | 24/pk      |



#### INSTRUMENT CLEANING BRUSHES (DOUBLE ENDED)

| PRODUCT NUMBER | FEATURES       | QUANTITIES |
|----------------|----------------|------------|
| 10-1444        | Nylon Bristles | 3/pk       |



#### INSTRUMENT CLEANING BRUSHES (DOUBLE TAPERED)

| PRODUCT NUMBER | FEATURES                 | QUANTITIES |
|----------------|--------------------------|------------|
| 10-1384        | 6", Double Tapered Brush | 3/pk       |
| 10-1388        | 8", Double Tapered Brush | 3/pk       |



## BRUSHES Endoscopy

OB/GYN

### CANNULA CLEANING BRUSHES

| PRODUCT NUMBER | FEATURES |                | QUANTITIES |
|----------------|----------|----------------|------------|
| 10-1350        | 12"      | 2mm Diameter   | 3/pk       |
| 10-1660        | 12"      | 3mm Diameter   | 3/pk       |
| 10-1661        | 12"      | 3mm Diameter   | 12/pk      |
| 10-1663        | 12"      | 3mm Diameter   | 24/pk      |
| 10-1352        | 12"      | 4mm Diameter   | 3/pk       |
| 10-1354        | 12"      | 5mm Diameter   | 3/pk       |
| 10-1356        | 12"      | 6mm Diameter   | 3/pk       |
| 10-1664        | 15"      | 9mm Diameter   | 3/pk       |
| 10-1665        | 15"      | 9mm Diameter   | 12/pk      |
| 10-1666        | 15"      | 9mm Diameter   | 24/pk      |
| 10-1668        | 15"      | 25mm Diameter  | 3/pk       |
| 10-1669        | 15"      | 25mm Diameter  | 12/pk      |
| 10-1672        | 15"      | 25mm Diameter  | 24/pk      |
| 10-1362        | 18"      | 5mm Diameter   | 3/pk       |
| 10-1364        | 18"      | 10mm Diameter  | 3/pk       |
| 10-1366        | 18"      | 15mm Diameter  | 3/pk       |
| 78-2824        | 19"      | 8mm Diameter   | 3/pk       |
| 10-1370        | 24"      | 2.5mm Diameter | 3/pk       |
| 10-1372        | 24"      | 5mm Diameter   | 3/pk       |
| 10-1374        | 24"      | 7mm Diameter   | 3/pk       |
| 10-1376        | 24"      | 10mm Diameter  | 3/pk       |



## OB/GYN

### SklarTip™ Instrument Identification

Instrument recognition is key to efficient processing in a busy clinical setting. SklarTip™ Instrument ID products offer a convenient color coded system for identifying instruments. Offered in various combinations including tape sheets, tape rolls and tags, products last through many autoclaving and gas sterilization processes. Latex Free.

#### ID SHEET TAPE

| PRODUCT NUMBER | FEATURES                      | QUANTITIES |
|----------------|-------------------------------|------------|
| 96-1845        | Blue with White Stripes       | 1 Each     |
| 96-1846        | Green with White Stripes      | 1 Each     |
| 96-1847        | Red with White Stripes        | 1 Each     |
| 96-1848        | Yellow with White Stripes     | 1 Each     |
| 96-1849        | Orange with White Stripes     | 1 Each     |
| 96-1850        | Brown with White Stripes      | 1 Each     |
| 96-1851        | Lavender with White Stripes   | 1 Each     |
| 96-1852        | Light Blue with White Stripes | 1 Each     |
| 96-1853        | Pink with White Stripes       | 1 Each     |



#### ID ROLL TAPE

| PRODUCT NUMBER | FEATURES           | QUANTITIES |
|----------------|--------------------|------------|
| 96-1860        | White 300" x 1/4"  | 1 Each     |
| 96-1862        | Green 300" x 1/4"  | 1 Each     |
| 96-1864        | Yellow 300" x 1/4" | 1 Each     |
| 96-1866        | Orange 300" x 1/4" | 1 Each     |
| 96-1868        | Black 300" x 1/4"  | 1 Each     |
| 96-1870        | Purple 300" x 1/4" | 1 Each     |
| 96-1872        | Blue 300" x 1/4"   | 1 Each     |
| 96-1874        | Red 300" x 1/4"    | 1 Each     |
| 96-1876        | Brown 300" x 1/4"  | 1 Each     |



#### INSTRUMENT ID TAGS

| PRODUCT NUMBER | FEATURES                  | QUANTITIES |
|----------------|---------------------------|------------|
| 96-1900        | ID Repair Tags, Red       | 100/pk     |
| 96-1901        | ID Repair Tags, Red       | 1000/pk    |
| 96-1903        | Sklar Tags, White         | 100/pk     |
| 96-1905        | ID Tags for Sharpen, Blue | 100/pk     |
| 96-1906        | ID Tags for Sharpen, Blue | 1000/pk    |





## SklarTip™ Instrument Protectors

Sklar Instrument Protectors provide a convenient, color-coded means of protecting surgical instrument tips and edges. Autoclavable, steam and gas permeable, and available in a wide variety of shapes and sizes, SklarTip™ Instrument Protectors provide exceptional protection during handling and sterilization. Latex Free.



### INSTRUMENT PROTECTORS - NON-VENTED

| PRODUCT NUMBER | FEATURES                      | QUANTITIES    |
|----------------|-------------------------------|---------------|
| 96-1510        | 19 x 1.5mm, 3/4" x 1/16"      | White 100/pk  |
| 96-1508        | 19 x 2.0mm, 3/4" x 5/64"      | Blue 100/pk   |
| 96-1506        | 19 x 3.0mm, 3/4" x 1/8"       | Green 100/pk  |
| 96-1504        | 25 x 3.0mm, 1" x 1/8"         | Red 100/pk    |
| 96-1502        | 25 x 5.0mm, 1" x 3/16"        | Yellow 100/pk |
| 96-1500        | 19 x 1.0mm, 3/4" x 3/64"      | Maroon 100/pk |
| 96-1516        | 25 x 1.5mm, 1" x 3/8" x 1/16" | Orange 100/pk |
| 96-1514        | 25 x 2.0mm, 1" x 5/8" x 5/64" | Brown 100/pk  |
| 96-1512        | 25 x 3.0mm, 1" x 1" x 1/8"    | Black 100/pk  |
| 96-1518        | Assorted                      | 100/pk        |

### INSTRUMENT PROTECTORS - VENTED


| PRODUCT NUMBER | FEATURES                      | QUANTITIES    |
|----------------|-------------------------------|---------------|
| 96-1410        | 19 x 1.5mm, 3/4" x 1/16"      | White 100/pk  |
| 96-1408        | 19 x 2.0mm, 3/4" x 5/64"      | Blue 100/pk   |
| 96-1406        | 19 x 3.0mm, 3/4" x 1/8"       | Green 100/pk  |
| 96-1404        | 25 x 3.0mm, 1" x 1/8"         | Red 100/pk    |
| 96-1402        | 25 x 5.0mm, 1" x 3/16"        | Yellow 100/pk |
| 96-1400        | 19 x 1.0mm, 3/4" x 3/64"      | Maroon 100/pk |
| 96-1416        | 25 x 1.5mm, 1" x 3/8" x 1/16" | Orange 100/pk |
| 96-1414        | 25 x 2.0mm, 1" x 5/8" x 5/64" | Brown 100/pk  |
| 96-1412        | 25 x 3.0mm, 1" x 1" x 1/8"    | Black 100/pk  |
| 96-1418        | Assorted                      | 100/pk        |

### OSTEOTOME GUARDS - VENTED

| PRODUCT NUMBER | FEATURES        | QUANTITIES |
|----------------|-----------------|------------|
| 96-1440        | 3/4" x 1/4"     | 50/pk      |
| 96-1442        | 1" x 3/8"       | 50/pk      |
| 96-1444        | 1" x 1/2"       | 50/pk      |
| 96-1446        | 1 1/2" x 1 1/2" | 50/pk      |

### TIP GRIPS

| PRODUCT NUMBER | FEATURES     | QUANTITIES   |
|----------------|--------------|--------------|
| 96-1550        | .062" Yellow | 3pr/pk, 5pks |
| 96-1552        | .062" Yellow | 5pr/pk, 5pks |



## OB/GYN

## PLASTICWARE

### Trays

#### PLASTIC TRAYS

| STERILE        |                    |            |
|----------------|--------------------|------------|
| PRODUCT NUMBER | FEATURES           | QUANTITIES |
| 96-7025        | Small Shallow Prep | 25/cs      |
| 96-7033        | Half Shallow       | 25/cs      |
| 96-7024        | Half Deep          | 25/cs      |
| 96-7175        | Full Deep          | 25/cs      |
| 96-7177        | 2-Compartment      | 25/cs      |
| 96-7178        | 3-Compartment      | 25/cs      |



#### PROCEDURE TRAYS\*

| NON-STERILE    |                                                    |               |
|----------------|----------------------------------------------------|---------------|
| PRODUCT NUMBER | FEATURES                                           | QUANTITIES    |
| 10-1960        | Mini, Non-Perforated, Blue, 10 1/8" x 6 5/8"       | 1 Each        |
| 10-1965        | Mini, Perforated, Blue, 10 1/8" x 6 5/8"           | 1, 6 or 12/cs |
| 10-1961        | Small, Non-Perforated, Blue, 10 7/16" x 13 9/16"   | 1 Each        |
| 10-1966        | Small, Perforated, Blue, 10 7/16" x 13 9/16"       | 1, 6 or 12/cs |
| 10-1962        | Medium, Non-Perforated, Blue, 11 3/4" x 16"        | 1 Each        |
| 10-1967        | Medium, Perforated, Blue, 11 3/4" x 16"            | 1, 6 or 12/cs |
| 10-1964        | Large, Heavy Duty, Non-Perforated, Blue, 15" x 20" | 1 Each        |



#### SOAKING TRAYS WITH REMOVABLE BASKET

| NON-STERILE            |                                          |            |
|------------------------|------------------------------------------|------------|
| PRODUCT NUMBER         | FEATURES                                 | QUANTITIES |
| 10-2555                | Small Instrument Tray, 8" x 3 1/2" x 2"  | 1 Each     |
| 10-2556                | Large Instrument Tray, 10" x 5" x 2 1/2" | 1 Each     |
| 10-2557 (not pictured) | Round Cup Size Tray, 3" x 3", 4 oz.      | 1 Each     |



\*Select products are available individually or in bulk packaging. To order large quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

## PLASTICWARE

### General

OB/GYN

#### EMESIS BASINS


| STERILE        |          |            |
|----------------|----------|------------|
| PRODUCT NUMBER | FEATURES | QUANTITIES |
| 96-7543        | 700cc    | 25/cs      |



#### OR BASINS

| STERILE        |                    |            |
|----------------|--------------------|------------|
| PRODUCT NUMBER | FEATURES           | QUANTITIES |
| 96-1679        | Round, 7 qt., Blue | 5/cs       |

| NON-STERILE    |                    |            |
|----------------|--------------------|------------|
| PRODUCT NUMBER | FEATURES           | QUANTITIES |
| 08-1011        | Round, 5 qt., Pink | 50/cs      |



\*Select products are available individually or in bulk packaging. To order large quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

#### GUIDEWIRE BOWLS

| STERILE        |             |            |
|----------------|-------------|------------|
| PRODUCT NUMBER | FEATURES    | QUANTITIES |
| 96-1630        | 8 1/4" x 3" | 10/cs      |



#### UTILITY BOWLS\*

| NON-STERILE    |                                      |            |
|----------------|--------------------------------------|------------|
| PRODUCT NUMBER | FEATURES                             | QUANTITIES |
| 10-1502        | Sponge/Utility Bowl, Aqua, 1 1/4 qt. | 1 or 12/cs |



| NON-STERILE    |                                          |            |
|----------------|------------------------------------------|------------|
| PRODUCT NUMBER | FEATURES                                 | QUANTITIES |
| 10-1504        | Solution Bowl, Reusable, Lt Green, 7 qt. | 1 or 12/cs |





## OB/GYN

## PLASTICWARE

### General

#### MEDICINE CUPS

STERILE

| PRODUCT NUMBER | FEATURES | QUANTITIES |
|----------------|----------|------------|
| 96-7027        | 2 oz.    | 25/cs      |



#### PITCHERS\*

STERILE

| PRODUCT NUMBER | FEATURES            | QUANTITIES |
|----------------|---------------------|------------|
| 96-7268        | Open Handle, 1200ml | 5/cs       |

NON-STERILE

| PRODUCT NUMBER | FEATURES            | QUANTITIES |
|----------------|---------------------|------------|
| 96-7270        | Open Handle, 1200ml | 35/cs      |



NON-STERILE

| PRODUCT NUMBER | FEATURES                            | QUANTITIES |
|----------------|-------------------------------------|------------|
| 10-1546        | Graduated Measure, 16 oz., Reusable | 1 or 6/cs  |



#### SPECIMEN TUBES

STERILE

| PRODUCT NUMBER | FEATURES      | QUANTITIES            |
|----------------|---------------|-----------------------|
| 96-7459        | 3 1/2" x 1/2" | 3/pk, 50pk/cs, 150/cs |



#### PLASTIC BOWLS

STERILE

| PRODUCT NUMBER | FEATURES | QUANTITIES |
|----------------|----------|------------|
| 96-7028        | 8 oz.    | 25/cs      |
| 96-7020        | 16 oz.   | 10/cs      |
| 96-7174        | 32 oz.   | 20/cs      |
| 96-7451        | 48 oz.   | 5/cs       |



\*Select products are available individually or in bulk packaging. To order large quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

## STEELWARE

## OB/GYN

### Sterilization Containers

## SKLARLITE™

### Sterilization Container System

The SklarLite™ Rigid Sterilization Container System provides instruments an organized, routine and secure approach to the sterilization of surgical instruments. Customize your SklarLite™ System by selecting colored perforated container lids, solid safety lids and ID labels. Available in your choice of six vibrant colors for easy identification of instrument sets by procedure and/or department. Choose from Silver, Yellow, Green, Blue, Black and Red.

#### BASKETS

SklarLite™ features stackable inner baskets to optimize internal space while silicone mats securely protect delicate and laparoscopic instruments during sterilization, transport and storage.

| PRODUCT NUMBER | MODEL     | FEATURES        |
|----------------|-----------|-----------------|
| 10-4784        | Half Size | 10" x 9½" x 1"  |
| 10-4786        | Half Size | 10" x 9½" x 2"  |
| 10-4788        | Half Size | 10" x 9½" x 3"  |
| 10-4790        | Half Size | 10" x 9½" x 4"  |
| 10-4478        | Mid Size  | 16" x 10" x 1"  |
| 10-4480        | Mid Size  | 16" x 10" x 2"  |
| 10-4482        | Mid Size  | 16" x 10" x 3"  |
| 10-4484        | Mid Size  | 16" x 10" x 4"  |
| 10-4284        | Full Size | 21¼" x 10" x 1" |
| 10-4286        | Full Size | 21¼" x 10" x 2" |
| 10-4288        | Full Size | 21¼" x 10" x 3" |
| 10-4290        | Full Size | 21¼" x 10" x 4" |

#### SILICONE MATS

| PRODUCT NUMBER | MODEL     | FEATURES  |
|----------------|-----------|-----------|
| 10-4970        | Half Size | 8½" x 9"  |
| 10-4975        | Mid Size  | 15" x 9"  |
| 10-4985        | Full Size | 20½" x 9" |



#### ACCESSORIES

| PRODUCT NUMBER | MODEL                   | FEATURES                                      |
|----------------|-------------------------|-----------------------------------------------|
| 10-4900        | Process Indicator Cards | Steam, 250/box                                |
| 10-4910        | Filters                 | Disposable, 100/pack                          |
| 10-4925        | Security Seals          | 100/pack                                      |
| 10-4930        | ID Labels               | Silver                                        |
| 10-4930ENG     | ID Labels               | Silver Custom Engraved, Maximum 16 Characters |



## OB/GYN

## STEELWARE

### Sterilization Containers

## SKLARLITE™

### Sterilization Container System

#### HALF SIZE

| NON-PERFORATED<br>BOTTOM | PERFORATED<br>BOTTOM | MODEL    | FEATURES         |
|--------------------------|----------------------|----------|------------------|
| 10-4500                  | 10-4512              | Standard | 11¼" x 11" x 4"  |
| 10-4548                  | 10-4560              | Standard | 11¼" x 11" x 5"  |
| 10-4596                  | 10-4608              | Standard | 11¼" x 11" x 6"  |
| 10-4644                  | 10-4656              | Standard | 11¼" x 11" x 8"  |
| 10-4692                  | 10-4704              | Standard | 11¼" x 11" x 10" |
| 10-4524                  | 10-4536              | Safe     | 11¼" x 11" x 4"  |
| 10-4572                  | 10-4584              | Safe     | 11¼" x 11" x 5"  |
| 10-4620                  | 10-4632              | Safe     | 11¼" x 11" x 6"  |
| 10-4668                  | 10-4680              | Safe     | 11¼" x 11" x 8"  |
| 10-4716                  | 10-4728              | Safe     | 11¼" x 11" x 10" |

#### MID SIZE

| NON-PERFORATED<br>BOTTOM | PERFORATED<br>BOTTOM | MODEL    | FEATURES        |
|--------------------------|----------------------|----------|-----------------|
| 10-4300                  | 10-4312              | Standard | 18¼" x 11" x 4" |
| 10-4348                  | 10-4360              | Standard | 18¼" x 11" x 5" |
| 10-4394                  | 10-4406              | Standard | 18¼" x 11" x 6" |
| 10-4324                  | 10-4336              | Safe     | 18¼" x 11" x 4" |
| 10-4370                  | 10-4382              | Safe     | 18¼" x 11" x 5" |
| 10-4418                  | 10-4430              | Safe     | 18¼" x 11" x 6" |

#### FULL SIZE

| NON-PERFORATED<br>BOTTOM | PERFORATED<br>BOTTOM | MODEL    | FEATURES        |
|--------------------------|----------------------|----------|-----------------|
| 10-4000                  | 10-4012              | Standard | 23" x 11" x 4"  |
| 10-4048                  | 10-4060              | Standard | 23" x 11" x 5"  |
| 10-4096                  | 10-4108              | Standard | 23" x 11" x 6"  |
| 10-4144                  | 10-4156              | Standard | 23" x 11" x 8"  |
| 10-4192                  | 10-4204              | Standard | 23" x 11" x 10" |
| 10-4024                  | 10-4036              | Safe     | 23" x 11" x 4"  |
| 10-4072                  | 10-4084              | Safe     | 23" x 11" x 5"  |
| 10-4120                  | 10-4132              | Safe     | 23" x 11" x 6"  |
| 10-4168                  | 10-4180              | Safe     | 23" x 11" x 8"  |
| 10-4216                  | 10-4228              | Safe     | 23" x 11" x 10" |

Containers and accessories for the SklarLite™ Sterilization Container System are available in a variety of sizes and colors. Contact your Sklar Account Manager for information on the complete product line.

## SKLARLITE™

### Sterilization Container System

The SklarLite™ Rigid Sterilization Container System is a reusable rigid container made from anodized aluminum alloy, designed to be strong and durable yet lightweight and easy to handle.

- Stackable Containers
- Stackable Baskets
- 6 Month Shelf Life
- Load Limit 25lbs.



#### TWO MODELS TO CHOOSE FROM:

##### STANDARD

Consists of a body and perforated lid and used with Sklar's disposable filters and other compatible accessories for sterilization. Choose from either perforated or non-perforated bottoms.

##### SAFE

Consists of identical body and lid, but features an additional solid cover that safeguards the perforation field and contents during transportation and storage.

## STEELWARE

OB/GYN

### Instrument Stringers

#### WEISENBACH INSTRUMENT HOLDER

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1775        | 6"          | 1 Each     |
| 10-1776        | 12"         | 1 Each     |



#### WORCESTER INSTRUMENT HOLDER

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1839        | 13" x 3"    | 1 Each     |



#### BUNT INSTRUMENT STRINGER

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1770        |             | 1 Each     |




#### SAUNDERS INSTRUMENT RACK

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1780        |             | 1 Each     |



#### MAYO INSTRUMENT HOLDER

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1772        |             | 1 Each     |



#### POSILOCK INSTRUMENT RACK

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1763        |             | 1 Each     |





## OB/GYN

## STEELWARE

### Instrument Stringers

#### WEINSTEIN INSTRUMENT STRINGERS

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1800        | 5½" x 3¼"   | 1 Each     |
| 10-1805        | 11¼" x 3¼"  | 1 Each     |



#### INSTRUMENT STRINGERS

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1340        | 2½" x 6"    | 1 Each     |
| 10-1342        | 2½" x 8"    | 1 Each     |
| 10-1344        | 2½" x 10"   | 1 Each     |
| 10-1346        | 2½" x 12"   | 1 Each     |
| 10-1347        | 2½" x 14"   | 1 Each     |
| 10-1348        | 2½" x 16"   | 1 Each     |
| 10-1349        | 2½" x 18"   | 1 Each     |



#### E-Z LOCK INSTRUMENT STRINGERS

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 10-1806        | 2½" x 6"    | 1 Each     |
| 10-1808        | 2½" x 8"    | 1 Each     |
| 10-1810        | 2½" x 10"   | 1 Each     |
| 10-1812        | 2½" x 12"   | 1 Each     |
| 10-1814        | 2½" x 14"   | 1 Each     |
| 10-1816        | 2½" x 16"   | 1 Each     |
| 10-1818        | 2½" x 18"   | 1 Each     |
| 10-1826        | 5" x 6"     | 1 Each     |
| 10-1828        | 5" x 8"     | 1 Each     |
| 10-1830        | 5" x 10"    | 1 Each     |
| 10-1832        | 5" x 12"    | 1 Each     |
| 10-1834        | 5" x 14"    | 1 Each     |
| 10-1836        | 5" x 16"    | 1 Each     |
| 10-1838        | 5" x 18"    | 1 Each     |



## STEELWARE Instrument Trays

OB/GYN

### LIFTING HANDLES

| PRODUCT NUMBER | DESCRIPTION             | QUANTITIES |
|----------------|-------------------------|------------|
| 10-1725        | For Sterilization Trays | 1 Pair     |



### FLASH STERILIZATION TRAY

| PRODUCT NUMBER | DESCRIPTION                                        | QUANTITIES |
|----------------|----------------------------------------------------|------------|
| 10-1711        | 10" x 6½" x 2"                                     | 1 Each     |
| 10-1712        | Single Handle for Flash Sterilization Tray 10-1711 | 1 Each     |




### DOUBLE DECKER MICRO TRAY

| PRODUCT NUMBER | DESCRIPTION     | QUANTITIES |
|----------------|-----------------|------------|
| 10-1717        | 10" x 15" x 2¼" | 1 Each     |



### MICRO TRAY WITH HINGED LID

| PRODUCT NUMBER | DESCRIPTION                                 | Quantities |
|----------------|---------------------------------------------|------------|
| 10-1706        | 8½" x 8" x 1½" with 2½" x 2½" Compartment   | 1 Each     |
| 10-1707        | 8½" x 12" x 1½"                             | 1 Each     |
| 10-1708        | 8½" x 14" x 1½"                             | 1 Each     |
| 10-1702        | 10½" x 10" x 1½"                            | 1 Each     |
| 10-1716        | 10½" x 20" x 2½" with 10½" x 5" Compartment | 1 Each     |



| PRODUCT NUMBER | DESCRIPTION      | QUANTITIES |
|----------------|------------------|------------|
| 10-1700        | 8½" x 8" x 1¼"   | 1 Each     |
| 10-1703        | 10½" x 15" x 1¼" | 1 Each     |



### MICRO TRAY WITH HINGED LID (Not pictured)

| PRODUCT NUMBER | DESCRIPTION      | Quantities |
|----------------|------------------|------------|
| 10-1701        | 10" x 6½" x 7⁄8" | 1 Each     |

### STERILIZATION BOX (Not pictured)

| PRODUCT NUMBER | DESCRIPTION     | QUANTITIES |
|----------------|-----------------|------------|
| 10-1705        | 6¾" x 3¾" x 1¾" | 1 Each     |



## OB/GYN

## STEELWARE

### Instrument Trays

#### STERILIZATION TRAY / BATH SETS

| PRODUCT NUMBER | DESCRIPTION         |         |                   | QUANTITIES |
|----------------|---------------------|---------|-------------------|------------|
| 10-1947        | Large Set Includes: |         |                   |            |
|                | Solid Tray          | 10-1948 | 20¾" x 12¾" x 6"  | 1          |
|                | Perforated Tray     | 10-1949 | 20¾" x 12¾" x 4"  | 1          |
|                | Cover               | 10-1535 | 20¾" x 12¾"       | 1          |
| 10-1679        | Small Set Includes: |         |                   |            |
|                | Solid Tray          | 10-1747 | 12¾" x 10½" x 4"  | 1          |
|                | Perforated Tray     | 10-1697 | 12¾" x 10½" x 2½" | 1          |
|                | Cover               | 10-1882 | 12¾" x 10½"       | 1          |



#### STERILIZATION TRAYS\*

| PRODUCT NUMBER | DESCRIPTION                             | COVERS  | QUANTITIES     |
|----------------|-----------------------------------------|---------|----------------|
| 10-1850        | 4 qt., 12¾" x 10½" x 2½"                | 10-1864 | 1 Each         |
| 10-1860        | 7 qt., 12¾" x 10½" x 4"                 | 10-1864 | 1 Each         |
| 10-1750        | 5 ¼ qt., 20¾" x 6¾" x 4"                | 10-1854 | 1, 6 or 12/box |
| 10-1752        | 9 qt., 20¾" x 6¾" x 6"                  | 10-1854 | 1, 6 or 12/box |
| 10-1953        | 8 ½ qt., 20¾" x 12¾" x 2½", Heavyweight | 10-1535 | 1, 6 or 12/box |
| 10-1954        | 14 ½ qt., 20¾" x 12¾" x 4", Heavyweight | 10-1535 | 1, 6 or 12/box |
| 10-1948        | 24 qt., 20¾" x 12¾" x 6"                | 10-1535 | 1, 6 or 12/box |



#### PERFORATED STERILIZATION TRAYS

| PRODUCT NUMBER | DESCRIPTION                  | QUANTITIES |
|----------------|------------------------------|------------|
| 10-1710        | 10" x 6½" x 2½"              | 1 Each     |
| 10-1713        | 10" x 10½" x 3½"             | 1 Each     |
| 10-1714        | 15" x 10½" x 3½"             | 1 Each     |
| 10-1720        | 20" x 10½" x 3½"             | 1 Each     |
| 10-1957        | 20¾" x 12¾" x 4" Mesh bottom | 1 Each     |



\*All products are packaged individually. Instrument trays and covers are sold separately unless otherwise noted. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.



## STEELWARE Instrument Trays

OB/GYN

### INSTRUMENT TRAYS, PERFORATED BOTTOM

| PRODUCT NUMBER | DESCRIPTION                                                 | QUANTITIES |
|----------------|-------------------------------------------------------------|------------|
| 10-1740        | 8 $\frac{1}{8}$ " x 5" x 2"                                 | 1 Each     |
| 10-1874        | 10" x 6 $\frac{1}{2}$ " x 2"                                | 1 Each     |
| 10-1892        | 12" x 7 $\frac{3}{4}$ " x 2"                                | 1 Each     |
| 10-1879        | 12" x 7 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ "                 | 1 Each     |
| 10-1758        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ " | 1 Each     |
| 10-1862        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 4"                | 1 Each     |
| 10-1894        | 16 $\frac{1}{2}$ " x 10" x 2 $\frac{1}{2}$ "                | 1 Each     |
| 10-1895        | 20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ " | 1 Each     |

### INSTRUMENT TRAYS, PERFORATED BOTTOM / SIDES\*

| PRODUCT NUMBER | DESCRIPTION                                                 | QUANTITIES     |
|----------------|-------------------------------------------------------------|----------------|
| 10-1757        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 1 $\frac{1}{4}$ " | 1 Each         |
| 10-1697        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ " | 1 Each         |
| 10-1943        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 4"                | 1 Each         |
| 10-1944        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 6"                | 1 Each         |
| 10-1941        | 20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 1 $\frac{1}{4}$ " | 1 Each         |
| 10-1759        | 20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ " | 1 Each         |
| 10-1949        | 20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 4"                | 1, 6 or 12/box |
| 10-1942        | 20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 6"                | 1 Each         |

### INSTRUMENT TRAYS, SOLID\*

| PRODUCT NUMBER | DESCRIPTION                                                             | COVERS  | QUANTITIES     |
|----------------|-------------------------------------------------------------------------|---------|----------------|
| 10-1755        | 6 $\frac{1}{8}$ " x 4 $\frac{1}{4}$ " x 2"                              | 10-1737 | 1, 6 or 12/box |
| 10-1741        | 6 $\frac{1}{8}$ " x 6 $\frac{1}{8}$ " x 4" Heavyweight                  | 10-1835 | 1 or 6/box     |
| 10-1764        | 10" x 6 $\frac{1}{2}$ " x 2"                                            | 10-1766 | 1 Each         |
| 10-1742        | 10 $\frac{3}{8}$ " x 6 $\frac{3}{8}$ " x 2 $\frac{1}{2}$ "              | 10-1542 | 1, 6 or 12/box |
| 10-1877        | 12" x 7 $\frac{3}{4}$ " x 2 $\frac{1}{4}$ "                             | 10-1878 | 1, 6 or 12/box |
| 10-1744        | 12 $\frac{3}{4}$ " x 6 $\frac{7}{8}$ " x 2 $\frac{1}{2}$ "              | 10-1724 | 1, 6 or 12/box |
| 10-1746        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "             | 10-1882 | 1, 6 or 12/box |
| 10-1950        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ " Heavyweight | 10-1882 | 1, 6 or 12/box |
| 10-1747        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 4"                            | 10-1882 | 1, 6 or 12/box |
| 10-1952        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 4" Heavyweight                | 10-1882 | 1, 6 or 12/box |
| 10-1870        | 12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 6"                            | 10-1882 | 1, 6 or 12/box |
| 10-1753        | 13 $\frac{7}{8}$ " x 12 $\frac{7}{8}$ " x 2 $\frac{1}{2}$ "             | 10-1754 | 1, 6 or 12/box |
| 10-1748        | 16 $\frac{1}{2}$ " x 10" x 2 $\frac{1}{2}$ "                            | 10-1756 | 1, 6 or 12/box |
| 10-1749        | 16 $\frac{1}{2}$ " x 10" x 4"                                           | 10-1756 | 1 Each         |
| 10-1953        | 20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ " Heavyweight | 10-1535 | 1, 6 or 12/box |
| 10-1954        | 20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 4" Heavyweight                | 10-1535 | 1, 6 or 12/box |
| 10-1948        | 20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 6"                            | 10-1535 | 1, 6 or 12/box |



\*All products are packaged individually. Instrument trays and covers are sold separately unless otherwise noted. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

## OB/GYN

## STEELWARE

### Instrument Trays

#### MESH CLEANING BASKET\*

| PRODUCT NUMBER | DESCRIPTION      | QUANTITIES  |
|----------------|------------------|-------------|
| 10-1729        | 4¼" x 7¼" x 2"   | 1 Each      |
| 10-1731        | 6" x 9" x 3"     | 1 Each      |
| 10-1732        | 8½" x 12½" x 3½" | 1 Each      |
| 10-1733        | 9½" x 16" x 4"   | 1 Each      |
| 10-1734        | 12" x 21" x 5"   | 1 Each      |
| 10-1736        | 14" x 20" x 5"   | 1 or 12/box |



#### PERFORATED CLEANING BASKET

| PRODUCT NUMBER | DESCRIPTION      | QUANTITIES |
|----------------|------------------|------------|
| 10-1560        | 12¾" x 10½" x 6" | 1 Each     |



#### INSTRUMENT DRYING TRAYS\*

| PRODUCT NUMBER | DESCRIPTION      | QUANTITIES     |
|----------------|------------------|----------------|
| 10-1876        | 10" x 6½" x ¾"   | 1, 6 or 12/box |
| 10-1883        | 13⅛" x 9⅜" x ¾"  | 1, 6 or 12/box |
| 10-1884        | 15½" x 10⅝" x ¾" | 1, 6 or 12/box |
| 10-1683        | 17" x 11⅝" x ¾"  | 1, 6 or 12/box |
| 10-1888        | 19⅞" x 12⅝" x ¾" | 1, 6 or 12/box |
| 10-1680        | 21¼" x 16¼" x ⅞" | 1, 6 or 12/box |



## STEELWARE Instrument Trays

OB/GYN

\*All products are packaged individually. Instrument trays and covers are sold separately unless otherwise noted. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

### MAYO TRAYS (SOLID)\*

| PRODUCT NUMBER | DESCRIPTION       | QUANTITIES     |
|----------------|-------------------|----------------|
| 10-1690        | 10" x 6½" x ¾"    | 1, 6 or 12/box |
| 10-1691        | 13⅞" x 9⅜" x ¾"   | 1, 6 or 12/box |
| 10-1692        | 15⅞" x 10⅞" x ¾"  | 1, 6 or 12/box |
| 10-1693        | 17" x 11⅞" x ¾"   | 1, 6 or 12/box |
| 10-1694        | 19 ⅞" x 12⅞" x ¾" | 1, 6 or 12/box |
| 10-1695        | 21¼" x 16¼" x ⅞"  | 1, 6 or 12/box |



### MAYO TRAYS (PERFORATED)\*

| PRODUCT NUMBER | DESCRIPTION      | QUANTITIES     |
|----------------|------------------|----------------|
| 10-1684        | 10" x 6½" x ¾"   | 1, 6 or 12/box |
| 10-1685        | 13¾" x 9¾" x ¾"  | 1, 6 or 12/box |
| 10-1686        | 15⅞" x 10⅞" x ¾" | 1, 6 or 12/box |
| 10-1687        | 17" x 11⅞" x ¾"  | 1, 6 or 12/box |
| 10-1688        | 19⅞" x 12⅞" x ¾" | 1, 6 or 12/box |
| 10-1689        | 21¼" x 16¼" x ⅞" | 1 or 6/box     |



### CATHETER INSTRUMENT TRAYS\*

| PRODUCT NUMBER | DESCRIPTION     | COVERS         | QUANTITIES     |
|----------------|-----------------|----------------|----------------|
| 10-1955        | 8¾" x 3¼" x 1½" | 10-1956 - Dome | 1, 6 or 12/box |
| 10-1738        | 8⅞" x 5" x 2"   | 10-1722 - Dome | 1, 6 or 12/box |
| 10-1760        | 8⅞" x 5" x 3"   | 10-1722 - Dome | 1, 6 or 12/box |
| 10-1886        | 17¼" x 4½" x 2" | 10-1887 - Flat | 1, 6 or 12/box |





## OB/GYN

## STEELWARE

### Surgical / Lab Products

#### OINTMENT JAR WITH LID\*

| PRODUCT NUMBER | DESCRIPTION            | QUANTITIES  |
|----------------|------------------------|-------------|
| 10-1525        | 9 oz., 3 1/8" x 2 3/4" | 1 or 12/box |



#### MINI TRAY WITH LID

| PRODUCT NUMBER | DESCRIPTION          | QUANTITIES |
|----------------|----------------------|------------|
| 10-1768        | 8 1/2" x 3" x 1 1/2" | 1 Each     |





**CONVENIENT DISPENSING CANISTER**

SEE PAGE 265 FOR ORDERING DETAILS

## SKLAR DISINFECTANT™ SURFACE WIPES

BACTERICIDAL • VIRUCIDAL • FUNGICIDAL • TUBERCULOCIDAL

**3 MINUTE KILL TIME!**

- READY TO USE, NON-TOXIC, NON-IRRITATING FORMULA
- SPECIALLY FORMULATED FOR THE CONTROL OF CROSS CONTAMINATION
- EFFECTIVE AGAINST C.DIFF, TB, HBV, HCV, VRE, HYDROPHILIC AND LIPOPHILIC VIRUSES INCLUDING HIV-1 (AIDS VIRUS), BACTERIA AND FUNGI ON INANIMATE OBJECTS
- SAFE FOR ALL AREAS OF THE FACILITY INCLUDING THE NICU, OPERATING ROOMS, ISOLATION ROOMS, PATIENT CARE AREAS, LABORATORIES AND DIALYSIS
- SATISFIES OSHA'S BLOODBORNE PATHOGEN STANDARDS
- DESTROYS ODOR-CAUSING BACTERIA
- 7" X 10" NON-WOVEN, TEAR RESISTANT TOWELETTE
- CLEAN SCENT
- EPA REGISTERED

## STEELWARE

OB/GYN

### Surgical / Lab Products

#### DRESSING JAR WITH LID\*

| PRODUCT NUMBER | DESCRIPTION         | QUANTITIES     |
|----------------|---------------------|----------------|
| 10-1523        | ½ qt., 4½" x 2¼"    | 1 or 12/box    |
| 10-1533        | 1 ¼ qt., 4" x 5"    | 1 or 12/box    |
| 10-1538        | 1 ¾ qt., 4½" x 6 ¼" | 1 or 12/box    |
| 10-1527        | 3 qt., 5½" x 6¾"    | 1 or 12/box    |
| 10-1543        | 4¾ qt., 6¾" x 7¾"   | 1 or 12/box    |
| 10-1544        | 6 qt., 7¼" x 8½"    | 1, 6 or 12/box |
| 10-1547        | 8 qt., 8" x 9½"     | 1, 6 or 12/box |



#### FORCEP/THERMOMETER JAR\*

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES     |
|----------------|-------------|----------------|
| 10-1550        | 1½" x 3"    | 1, 6 or 12/box |
| 10-1541        | 2¼" x 4½"   | 1, 6 or 12/box |
| 10-1545        | 2¼" x 7"    | 1, 6 or 12/box |



#### FUNNELS WITH HANGING HOOK

| PRODUCT NUMBER | DESCRIPTION        | QUANTITIES |
|----------------|--------------------|------------|
| 10-1157        | 2½ oz., 3" x 3¾"   | 1 Each     |
| 10-1158        | 6½ oz., 4½" x 4¾"  | 1 Each     |
| 10-1159        | 13 oz., 5½" x 5½"  | 1 Each     |
| 10-1160        | 22¼ oz., 6¼" x 6½" | 1 Each     |
| 10-1170        | 32 oz., 7" x 8½"   | 1 Each     |
| 10-1171        | 64 oz., 8½" x 9½"  | 1 Each     |



#### LAB / UTILITY SCOOP

| PRODUCT NUMBER | DESCRIPTION                             | QUANTITIES |
|----------------|-----------------------------------------|------------|
| 10-1192        | 5 oz., 5" x 2½"<br>(Overall length 9½") | 1 Each     |





## OB/GYN

## STEELWARE

### Surgical / Lab Products

#### EMESIS BASINS\*

| PRODUCT NUMBER | DESCRIPTION                                                        | QUANTITIES  |
|----------------|--------------------------------------------------------------------|-------------|
| 10-1497        | 10 oz., 6 $\frac{1}{8}$ " x 3 $\frac{1}{2}$ " x 1 $\frac{1}{2}$ "  | 1 or 12/box |
| 10-1698        | 12 oz., 8 $\frac{1}{8}$ " x 3 $\frac{3}{8}$ " x 1 $\frac{1}{8}$ "  | 1 or 12/box |
| 10-1699        | 26 oz., 9 $\frac{7}{8}$ " x 4 $\frac{1}{2}$ " x 2 $\frac{1}{8}$ "  | 1 or 12/box |
| 10-1487        | 48 oz., 12 $\frac{1}{8}$ " x 5 $\frac{3}{4}$ " x 2 $\frac{1}{4}$ " | 1 Each      |



#### HOSPITAL BUCKETS

| PRODUCT NUMBER | DESCRIPTION                                              | QUANTITIES |
|----------------|----------------------------------------------------------|------------|
| 10-1182        | 13 qt., 11 $\frac{1}{8}$ " x 9 $\frac{1}{4}$ " w/o lid   | 1 Each     |
| 10-1450        | 13 qt., 11 $\frac{1}{8}$ " x 9 $\frac{1}{4}$ " with lid  | 1 Each     |
| 10-1183        | 16 qt., 11 $\frac{3}{4}$ " x 10 $\frac{1}{8}$ " w/o lid  | 1 Each     |
| 10-1451        | 16 qt., 11 $\frac{3}{4}$ " x 10 $\frac{1}{8}$ " with lid | 1 Each     |
| 10-1184        | 20 qt., 13 $\frac{3}{4}$ " x 11 $\frac{1}{8}$ " w/o lid  | 1 Each     |



#### HOSPITAL KICK BUCKET WITHOUT LID

| PRODUCT NUMBER | DESCRIPTION                                    | QUANTITIES |
|----------------|------------------------------------------------|------------|
| 10-1551        | 13 qt., 11 $\frac{1}{8}$ " x 9 $\frac{1}{4}$ " | 1 Each     |



#### STORAGE CONTAINERS\*

| PRODUCT NUMBER | DESCRIPTION                                                             | QUANTITIES     |
|----------------|-------------------------------------------------------------------------|----------------|
| 10-1540        | $\frac{1}{2}$ gal., 5 $\frac{7}{8}$ " x 6 $\frac{5}{8}$ "<br>no handles | 1, 6 or 12/box |
| 10-1186        | 3 gal., 10 $\frac{3}{8}$ " x 9"                                         | 1 Each         |
| 10-1188        | 4 gal., 13" x 9"                                                        | 1 Each         |
| 10-1189        | 5 gal., 13" x 11"                                                       | 1 Each         |
| 10-1190        | 6 gal., 13" x 13"                                                       | 1 Each         |



\*All products are packaged individually. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## STEELWARE

OB/GYN

### Surgical / Lab Products

#### GRADUATED MEASURE\*

| PRODUCT NUMBER | DESCRIPTION       | QUANTITIES     |
|----------------|-------------------|----------------|
| 10-1552        | 16 oz., 3¾" x 4⅝" | 1, 6 or 12/box |
| 10-1510        | 32 oz., 4¾" x 5¼" | 1, 6 or 12/box |
| 10-1674        | 64 oz., 5½" x 6½" | 1, 6 or 12/box |



#### SPONGE BOWLS\*

| PRODUCT NUMBER | DESCRIPTION       | QUANTITIES  |
|----------------|-------------------|-------------|
| 10-1496        | 12 oz., 4¼" x 2"  | 1 or 12/box |
| 10-1498        | 22 oz., 5" x 2⅜"  | 1 or 12/box |
| 10-1484        | 29 oz., 5⅞" x 2⅞" | 1 or 12/box |
| 10-1485        | 36 oz., 6¼" x 2¾" | 1 or 12/box |
| 10-1483        | 48 oz., 7⅞" x 2¾" | 1 or 12/box |
| 10-1486        | 64 oz., 7⅞" x 3"  | 1 or 12/box |



#### MIXING/SOLUTION BOWLS\*

| PRODUCT NUMBER | DESCRIPTION        | QUANTITIES     |
|----------------|--------------------|----------------|
| 10-1488        | ⅞ qt., 6⅞" x 2½"   | 1 or 12/box    |
| 10-1491        | 1⅞ qt., 8⅞" x 3"   | 1 or 12/box    |
| 10-1492        | 3 qt., 9¾" x 4"    | 1 or 12/box    |
| 10-1493        | 5 qt., 11⅞" x 4½"  | 1 or 12/box    |
| 10-1494        | 6¾ qt., 12¾" x 5¼" | 1, 6 or 12/box |
| 10-1495        | 8 qt., 13⅞" x 5⅞"  | 1, 6 or 12/box |



#### SOLUTION BASINS\*

| PRODUCT NUMBER | DESCRIPTION       | QUANTITIES     |
|----------------|-------------------|----------------|
| 10-1489        | 4½ qt., 12" x 4¼" | 1 or 12/box    |
| 10-1505        | 7 qt., 13⅞" x 4⅞" | 1, 6 or 12/box |
| 10-1507        | 9½ qt., 14½" x 5" | 1, 6 or 12/box |



## OB/GYN

## STEELWARE

### Surgical / Lab Products

All Sklar Stainless Steelware is produced according to the strictest standards to ensure superior performance and durability.

\*All products are packaged individually. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

#### IODINE CUP\*

| PRODUCT NUMBER | DESCRIPTION                                   | QUANTITIES          |
|----------------|-----------------------------------------------|---------------------|
| 10-1670        | 6 oz., 3 $\frac{7}{16}$ " x 2"                | 1, 12, 24 or 48/box |
| 10-1671        | 14 oz., 4 $\frac{3}{8}$ " x 2 $\frac{5}{8}$ " | 1, 12, 24 or 48/box |



#### MEDICINE CUP\*

| PRODUCT NUMBER | DESCRIPTION                                        | QUANTITIES          |
|----------------|----------------------------------------------------|---------------------|
| 10-1675        | 2 oz., 2 $\frac{1}{8}$ " x 2"<br>Graduated Measure | 1, 12, 24 or 48/box |



#### TUMBLER\*

| PRODUCT NUMBER | DESCRIPTION                                  | QUANTITIES |
|----------------|----------------------------------------------|------------|
| 10-1536        | 7 oz., 2 $\frac{1}{8}$ " x 3 $\frac{1}{8}$ " | 1, 12/box  |



#### PITCHER / BEAKER\*

| PRODUCT NUMBER | DESCRIPTION                                                      | QUANTITIES     |
|----------------|------------------------------------------------------------------|----------------|
| 10-1172        | 125ml, 2 $\frac{1}{8}$ " x 2 $\frac{1}{8}$ " w/o Handle          | 1 Each         |
| 10-1173        | 250ml, 2 $\frac{3}{8}$ " x 3 $\frac{1}{4}$ " w/o Handle          | 1 Each         |
| 10-1174        | 600ml, 3 $\frac{3}{8}$ " x 4 $\frac{1}{8}$ " w/o Handle          | 1 Each         |
| 10-1176        | 2000ml, 4 $\frac{7}{8}$ " x 7 $\frac{1}{8}$ " w/o Handle         | 1 Each         |
| 10-1177        | 3000ml, 5 $\frac{1}{2}$ " x 8 $\frac{1}{4}$ " w/o Handle         | 1 Each         |
| 10-1178        | 4000ml, 6 $\frac{1}{8}$ " x 8 $\frac{7}{8}$ " w/o Handle         | 1 Each         |
| 10-1175        | 1200ml, 4 $\frac{1}{8}$ " x 6 $\frac{1}{8}$ " with closed Handle | 1 Each         |
| 10-1539        | 2000ml, 5 $\frac{1}{8}$ " x 6 $\frac{3}{8}$ " with closed Handle | 1, 6 or 12/box |
| 10-1180        | 3000ml, 6 $\frac{1}{8}$ " x 7" with closed Handle                | 1 Each         |
| 10-1181        | 4000ml, 6 $\frac{3}{8}$ " x 8 $\frac{1}{8}$ " with closed Handle | 1 Each         |



## OB/GYN

## BRUSHES

### Sets


#### 10-1424 AUTOCLAVE CLEANING BRUSH SET

| PRODUCT NUMBER | DESCRIPTION                                                   | QUANTITIES |
|----------------|---------------------------------------------------------------|------------|
| 10-1424HDL     | Telescopic Handle extends from 24" to 48"                     | 1 Each     |
| 10-1424PBR     | Large Cleaning Brush 7½" L x 5½" W with round corner bristles | 1 Each     |
| 10-1424BHD     | Surface Cleaning Brush, 4½" handle                            | 1 Each     |



#### 10-1432 MICRO-OPTHALMIC BRUSH SET

| PRODUCT NUMBER | DESCRIPTION                          | QUANTITIES |
|----------------|--------------------------------------|------------|
| 10-1350        | Cannula Cleaning Brush 12" 2mm, 3/pk | 1 pk       |
| 10-1660        | Cannula Cleaning Brush 12" 3mm, 3/pk | 1 pk       |
| 10-1396        | Suction Tube Brush 6" 1.6mm, 3/pk    | 1 pk       |
| 10-1400        | Suction Tube Brush 8" 2mm, 3/pk      | 1 pk       |
| 10-1402        | Suction Tube Brush 8" 2.6mm, 3/pk    | 1 pk       |
| 10-1404        | Suction Tube Brush 8" 3.3mm, 3/pk    | 1 pk       |
| 10-1406        | Suction Tube Brush 8" 4mm, 3/pk      | 1 pk       |
| 10-1444        | Double End Brush Nylon, 3/pk         | 1 pk       |



#### 10-1434 ALL PURPOSE BRUSH SET

| PRODUCT NUMBER | DESCRIPTION                          | QUANTITIES |
|----------------|--------------------------------------|------------|
| 10-1650        | Brush Stainless Steel Bristle 3/pk   | 1 pk       |
| 10-1651        | Brush Nylon Bristle 3/pk             | 1 pk       |
| 10-1420        | Large Instrument Cleaning Brush 3/pk | 1 pk       |
| 10-1418        | Surface Cleaning Brush 3/pk          | 1 pk       |





## MedGyn Needle Adapters - Extenders

An effective means of extending ordinary hypodermic needles when longer needles are required. The paracervical and uterosacral areas can easily be reached with an extended needle.

### Needle Adapters

- 030501 Needle Adapter 3" (76mm)
- 030502 Needle Adapter 4" (102mm)
- 030503 Needle Adapter 5" (127mm)
- 030504 Needle Adapter 6" (152mm)




## Miscellaneous Products

- 040004 Emesis Basin Medium 8" (203mm)
- 040005 Emesis Basin Large 10" (254mm)
- 040003 Emesis Basin Small 6" (152mm)
- 040008 Iodine Cup 6oz
- 040001 Medicine Cup (2 oz.)
- 040009 Instrument Cleaner 8oz Bottle
- 040000 Instrument Lubricant Spray - 8oz
- 040010 Instrument Stain Remover 3oz
- 040002 Instrument Steel Tray - 12¼" X 7⅝" X 2⅛"  
(Graham 3259)








**MAYO-HEGAR NEEDLE HOLDERS**


|        |        |
|--------|--------|
| 97-842 | 5 1/4" |
| 97-847 | 6"     |
| 97-853 | 7"     |
| 97-869 | 8"     |



**WEBSTER NEEDLE HOLDER**

|        |    |        |
|--------|----|--------|
| 97-809 | 5" | Smooth |
|--------|----|--------|

## SKLAR CERVIX BRUSH CERVICAL CELL SAMPLER



### SKLAR CERVIX BRUSH

#### NON-STERILE

| PRODUCT NUMBER | QUANTITIES |
|----------------|------------|
| 96-4140NS      | 100/bx     |

#### STERILE

| PRODUCT NUMBER | QUANTITIES |
|----------------|------------|
| 96-4140        | 100/bx     |

Sklar's Cervix Brush is made of polyethylene and has a total length of 20 centimeters. The top portion uses a soft flexible brush to obtain cell samples. The brush consists of 57 semi-circular plastic bristles of various lengths. The shape is such that the top edges of the brush follow the contours of the cervix. The longer middle bristles reach deep into the endocervical canal. The shorter bristles touch both the ectocervical area and transformation zone at the same time.





## COTTLE TENACULUM HOOK

| ITEM #        | DESCRIPTION |
|---------------|-------------|
| HNM-1-20-4016 | 6 1/4", ANG |



## EMMETT TENACULUM HOOK

| ITEM #        | DESCRIPTION |
|---------------|-------------|
| HNM-1-20-4201 | 8 3/4", #1  |
| HNM-1-20-4202 | #2          |
| HNM-1-20-4203 | 8 3/4", #3  |
| HNM-1-20-4204 | #4          |
| HNM-1-20-4205 | 8 3/4", #5  |

## ASPIRATION SET



### ASPIRATION SET W/LOCKING 60CC SYRINGE AND CANNULA

| ITEM #      | DESCRIPTION           |
|-------------|-----------------------|
| HBI-R57-980 | 3mm, FLEXIBLE CURETTE |
| HBI-R57-981 | 4mm, FLEXIBLE CURETTE |
| HBI-R57-982 | 5mm, FLEXIBLE CURETTE |
| HBI-R57-983 | 6mm, FLEXIBLE CURETTE |
| HBI-R57-985 | 7mm, FLEXIBLE CURETTE |

## GYN SUPPLIES

### JACKSON TRACHEAL TENACULUM

| ITEM #        | DESCRIPTION |
|---------------|-------------|
| HNM-1-20-4014 | 5 3/4"      |



UTERINE SOUND



UTERINE SOUND

| ITEM #     | DESCRIPTION                                 |
|------------|---------------------------------------------|
| HBI-BA-100 | 9 3/4", STR, SINGLE USE, STERILE, BOX OF 25 |






## SIMS UTERINE CURETTE

| ITEM #         | DESCRIPTION                                 |
|----------------|---------------------------------------------|
| HNM-1-52-8150  | 11", #0, MALLEABLE SHAFT, SHARP             |
| HNM-1-52-8151  | 11", #1, MALLEABLE SHAFT, SHARP             |
| HNM-1-52-8152  | 11", #2, MALLEABLE SHAFT, SHARP             |
| HNM-1-52-8153  | 11", #3, MALLEABLE SHAFT, SHARP             |
| HNM-1-52-8154  | 11", #4, MALLEABLE SHAFT, SHARP             |
| HNM-1-52-8155  | 11", #5, MALLEABLE SHAFT, SHARP             |
| HNM-1-52-8156  | 11", #6, MALLEABLE SHAFT, SHARP             |
| HNM-1-52-8150S | 11", #0, MALLEABLE SHAFT, SHARP, STND GRADE |
| HNM-1-52-8152S | 11", #2, MALLEABLE SHAFT, SHARP, STND GRADE |
| HNM-1-52-8156S | 11", #6, MALLEABLE SHAFT, SHARP, STND GRADE |

## SIZE CHART




## CURETTES

### KEVORKIAN YOUNGE BIOPSY CURETTE

| ITEM #        | DESCRIPTION         |
|---------------|---------------------|
| HNM-1-52-7010 | W/O BASKET, 11 3/4" |
| HNM-1-52-7015 | W/BASKET, 12 1/8"   |

### NOVAK ENDOMETRIUM BIOPSY SUCTION CURETTE


| ITEM #        | DESCRIPTION                                |
|---------------|--------------------------------------------|
| HNM-1-52-7048 | 9", DIA 2MM                                |
| HNM-1-52-7050 | 9 7/8", DIA 4MM, W/LUER LOCK HUB AND TEETH |



AUVARD VAGINAL SPECULUM

| ITEM #        | DESCRIPTION                   |
|---------------|-------------------------------|
| HNM-1-52-2786 | 2 3/4" x 1 1/2" BLADES, 2LB   |
| HNM-1-52-2787 | 2 3/4" x 1 1/2" BLADES, 2.5LB |
| HNM-1-52-2788 | 2 3/4" x 1 1/2" BLADES, 3LB   |
| HNM-1-52-2793 | 4" x 1 3/4" BLADES, 2.5LB     |

## ENDOCERVICAL SPECULUM



KOGAN ENDOCERVICAL

| ITEM #        | DESCRIPTION               |
|---------------|---------------------------|
| HNM-1-52-3955 | W/GAUGE AND SET SCREW, 9" |
| HNM-1-52-3950 | W/RACHET, 9"              |





## WEISMAN GRAVES VAGINAL SPECULUM

| ITEM #        | DESCRIPTION                                       |
|---------------|---------------------------------------------------|
| HNM-1-52-0769 | SMALL, 3 1/4" x 1/4" BLADES, RIGHT SIDE OPENING   |
| HNM-1-52-0400 | LEFT SIDE OPENING                                 |
| HNM-36-150759 | MEDIUM, RIGHT SIDE OPENING                        |
| HNM-36-150760 | MEDIUM, 1 1/4" x 4" BLADES, LEFT SIDE OPENING     |
| HNM-1-52-0405 | LARGE, LEFT SIDE OPENING                          |
| HNM-1-52-0770 | LARGE, 4 3/4" x 1 1/2" BLADES, RIGHT SIDE OPENING |
| HNM-36-150761 | EXTRA LARGE, RIGHT SIDE OPENING                   |
| HNM-36-150762 | EXTRA LARGE, LEFT SIDE OPENING                    |

## WEISMAN PEDERSON VAGINAL SPECULUM

| ITEM #        | DESCRIPTION                |
|---------------|----------------------------|
| HNM-36-150763 | MEDIUM, RIGHT SIDE OPENING |
| HNM-36-150764 | MEDIUM, LEFT SIDE OPENING  |
| HNM-36-150766 | LARGE, RIGHT SIDE OPENING  |
| HNM-36-150765 | LARGE, LEFT SIDE OPENING   |


## VAGINAL SPECULUM

### GRAVES VAGINAL SPECULUM


| ITEM #        | DESCRIPTION                                    |
|---------------|------------------------------------------------|
| HNM-1-52-0607 | SMALL, 3/4" x 3" BLADES                        |
| HNM-1-52-0610 | MEDIUM, 1 1/4" x 4" BLADES                     |
| HNM-1-52-0612 | LARGE, 1 3/8" x 4 1/2" BLADES                  |
| HNM-36-150755 | EXTRA LARGE, 1 1/2" x 6 3/4" BLADES            |
| HNM-36-150756 | IMPROVED, MEDIUM, 1 3/8" x 4" BLADES, WIDE ANG |

### PEDERSON VAGINAL SPECULUM

| ITEM #        | DESCRIPTION                                     |
|---------------|-------------------------------------------------|
| HNM-1-52-0507 | SMALL, 1/2" x 3" BLADES                         |
| HNM-1-52-0510 | MEDIUM, 7/8" x 4" BLADES                        |
| HNM-1-52-0512 | LARGE, 1" x 4 3/4" BLADES                       |
| HNM-36-151353 | VIRGIN SIZE, 5/8" x 2 3/4" BLADES, EXTRA NARROW |



GRAVES



PEDERSON




TISCHLER-KEVORKIAN BIOPSY PUNCH

| ITEM #        | DESCRIPTION                         |
|---------------|-------------------------------------|
| HNM-36-150754 | 8", 3mm x 9.5mm BITE, SPRING HANDLE |




KEVORKIAN-PACIFIC BIOPSY PUNCH

| ITEM #     | DESCRIPTION                       |
|------------|-----------------------------------|
| H-233-5365 | 8", 3mm x 8mm BITE, SPRING HANDLE |




## BIOPSY PUNCHES



TISCHLER-MORGAN BIOPSY PUNCH

| ITEM #        | DESCRIPTION                        |
|---------------|------------------------------------|
| HNM-36-150753 | 8", 3mm x 7mm BITE, SPRING HANDLE  |
| HNM-1-52-6556 | 10", 3mm x 7mm BITE, SPRING HANDLE |



BABY TISCHLER BIOPSY PUNCH

| ITEM #        | DESCRIPTION                        |
|---------------|------------------------------------|
| HNM-36-150750 | 8", 2mm x 4mm BITE, SPRING HANDLE  |
| HNM-36-150752 | 10", 2mm x 4mm BITE, SPRING HANDLE |



## BRAUN TENACULUM FORCEPS

| ITEM #        | DESCRIPTION |
|---------------|-------------|
| HNM-1-52-4426 | 10", STR    |



## LAUFE UTERINE POLYP FORCEPS

| ITEM #        | DESCRIPTION |
|---------------|-------------|
| HNM-1-30-5640 | 8"          |



## FORCEPS



### SCHROEDER TENACULUM FORCEPS


| ITEM #         | DESCRIPTION              |
|----------------|--------------------------|
| HNM-36-151325  | 10", TIPS OVERLAP        |
| HNM-36-151324  | 10", TIPS DO NOT OVERLAP |
| HNM-1-52-4425  | 10"                      |
| HNM-1-52-4425S | 10", STND GRADE          |





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## SMOKE EVACUATOR TUBING

| ITEM #  | DESCRIPTION                            |
|---------|----------------------------------------|
| H-786T  | 7/8" x 6" TUBE, NON STERILE, BOX OF 24 |
| H-786TS | 7/8" x 6" TUBE, STERILE, BOX OF 24     |



## REDUCER

| ITEM #        | DESCRIPTION                                           |
|---------------|-------------------------------------------------------|
| HNM-36-150722 | 7/8", TO PRE-FILTER OR CANISTER (FLUTED TUBING) 10/PK |



## CYTOLOGY BRUSH

| ITEM #       | DESCRIPTION |
|--------------|-------------|
| HBI-TEAR1011 | 100/PK      |

OB/GYN

DISPOSABLES



## SKLARSPEC™

### DISPOSABLE GRAVES VAGINAL SPECULA

SKLARSPEC plastic specula are clear, disposable, and comfortable. Plastic material allows 360 degree visibility and reduces thermal shock. Provides protection of the doctor, the nurse and the patient from cross contamination. Clean up is eliminated by single-use only. Plastic material avoids unpleasant temperature shock associated with metal. One hand operation with wide range adjustment.

### PLASTIC VAGINAL SPECULUMS, GRAVES STYLE

#### STERILE

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 96-2602        | Small       | 50/cs      |
| 96-2604        | Medium      | 50/cs      |
| 96-2606        | Large       | 50/cs      |
| 96-1857        | Small       | 100/cs     |
| 96-1867        | Medium      | 100/cs     |
| 96-1877        | Large       | 100/cs     |

#### NON-STERILE

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 96-0857        | Small       | 100/cs     |
| 96-0867        | Medium      | 100/cs     |
| 96-0877        | Large       | 100/cs     |

## DISPOSABLES

## OB/GYN



### PLASTIC VAGINAL SPECULUMS, LIGHT SOURCE STYLE

#### NON-STERILE

| PRODUCT NUMBER | DESCRIPTION   | QUANTITIES |
|----------------|---------------|------------|
| 96-0827        | Small, Blue   | 100/cs     |
| 96-0837        | Medium, Green | 100/cs     |
| 96-0847        | Large, White  | 100/cs     |



### PLASTIC VAGINAL SPECULUMS, CUSCO STYLE STERILE

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES |
|----------------|-------------|------------|
| 96-2642        | Small       | 50/cs      |
| 96-2644        | Medium      | 50/cs      |
| 96-2648        | Large       | 50/cs      |

#### NON-STERILE

| PRODUCT NUMBER | DESCRIPTION | QUANTITIES             |
|----------------|-------------|------------------------|
| 96-0557        | Small       | 10/pk, 10pk/cs, 100/cs |
| 96-0567        | Medium      | 10/pk, 10pk/cs, 100/cs |
| 96-0577        | Large       | 10/pk, 10pk/cs, 100/cs |

## DISPOSABLES



### MedGyn Single-Use Forceps

029701-D Tischler Biopsy Single-Use Forceps, 9 3/4" (10/Pack)



029702-D Baby Tischler Biopsy Single-Use Forceps, 9 3/4" (10/Pack)



029703-D Kevorkian-Younger Biopsy Single-Use Forceps 9 3/4" (10/Pack)



030112-D Alligator Single-Use Forceps 8" (10/Pack)



### MedGyn Rotating Biopsy Forceps

MedGyn's Rotating Biopsy Punch is our latest single-use device with a plastic handle and stainless steel shaft to sample cervical tissue.

030100-D Disposable Rotating Biopsy Forceps



### MedGyn Needle Extender

An effective means of extending ordinary hypodermic needles when longer needles are required.

022600 MedGyn Disposable Needle Extender





## MedGyn Aspiration Kits

MedGyn Aspiration Kit and Aspiration Kit IV are effective devices for simultaneous endometrial curettage and aspiration.

022514 MedGyn Aspiration Kit with 4mm Cannula-Curette (10/Pack)

022515 MedGyn Aspiration Kit with 5mm Cannula-Curette (10/Pack)

022516 MedGyn Aspiration Kit with 6mm Cannula-Curette (10/Pack)

022517 MedGyn Aspiration Kit with 7mm Cannula-Curette (10/Pack)

022511 MedGyn Aspiration Kit IV with 4mm, 5mm, 6mm and 7mm Cannula-Curette (10/Pack)



## MedGyn Aspiration Kit Double Valve

MedGyn Aspiration Kit Double Valve is used for uterine aspiration/uterine evacuation in obstetrics and gynecologic procedures. The device fits the MedGyn Flexible Karman Curettes and a wide variety of MedGyn Rigid and Straight Curettes.

022526 MedGyn Aspiration Kit Double Valve



## MedGyn s IUD Retriever

MedGyn s IUD Retriever is designed to assist in the removal of the intrauterine device from the uterus.

022366 Disposable IUD Retriever (25/Box)



## MedGyn Os Finders

Ideal for cervical dilation, uterine sounding or other intrauterine procedures. Made of semi-flexible plastic and packaged individually sterile.

022731 MedGyn Disposable Os Finder (50/Pack)

## DISPOSABLES



## MedGyn Collection Sets, Collection Jars and Bags

MedGyn Disposable Collection Sets include 6' tubing (either 3/8" or 1/2" I.D.) attached to a control handle with slip-ring for finger-tip vacuum control. The Collection Sets are individually packaged and sterile.

### Collection Sets & Bottles

- 022310 3/8" Collection Set with Handle (50/Pack)
- 022310NS 3/8" Collection Set w/Handle & Adapter (Non-Sterile)
- 022310-WO 3/8" Collection Set with Handle (Without Adapter, Sterile) (50/Pack)
- 022314 3/8" Collection Set x 8 Ft with Handle & Adaptor
- 022312 1/2" Collection Set with Handle & Adapter (25/Pack)
- 022312NS 1/2" Collection Set w/Handle & Adapter (Non-Sterile)
- 022331 3/8" Disp. Collection Jar - with Cap (10/Box)
- 022332 1/2" Disp. Collection Jar - with Cap (10/Box)
- 022323 Disp. Tops Large Non-Conductive (10/Box)
- 022324 Disp. Tops Small Non-Conductive (10/Box)

### Metal Swivel Handles

Stainless steel, autoclavable, with slip-ring control

- 022450 3/8" Handle for Curettes up to 14mm with small base
- 022451 1/2" Handle for 14mm, 15mm & 16mm Curettes
- 022340 Stainless Steel 3/8" Female Fitting
- 022341 Stainless Steel 3/8" Male Fitting
- 022343 Stainless Steel 1/2" Male Fitting
- 022352 Stainless Steel 3/8" Handle Extender
- 022353 Stainless Steel 1/2" Handle Extender

### Accessories


- 022316-3/8-6 Tubing 3/8" Precut to 6'
- 022311 Tubing 3/8", Custom Cut Per Foot
- 022316-1/2-6 Tubing 1/2" Precut to 6'
- 022313 Tubing 1/2", Custom Cut Per Foot
- 022317 Tubing 16" Connector
- 022318 Tubing 20" Connector
- 022317-ADAP Tube Fitting - Male 3/8" Plastic
- 022317-INT Tube - Interconnecting with Adapter (10/Box)
- 022319-ADAP Tube Fitting - Male 1/2" Plastic
- 022350 Plastic Adapter for 1/2" Base Curettes (Sterile)
- 022351 Plastic Adapter for 1/2" Base Curettes (Non-Sterile)
- O-RING-FAD 3/8" O Ring for Female Adapter (10/Pack)
- 022320 Disposable Filter Assembly (10/Box)
- ADAP-REDC 1/2" to 3/8" Adapter for Aspiration Handle (For Small Curette) (10Box)
- 022354 3/8" to 1/2" Adapter for Aspiration Handle (For Large Curette)

### Permeable Gauze Sacks

Gauze collection sack easily attaches to collection canister for simple and efficient separation of tissue from aspirated fluid.

- 022460 Gauze Sacks Tissue Trap (10/Box)
- O-RING3/8 O Ring for 3/8" Sack (10/Box)
- O-RING1/2 O Ring for 1/2" Sack (10/Pack)





## MedGyn Vacuum Cannulas-Curettes

MedGyn offers a complete line of high-quality, disposable vacuum cannulas-curettes in the industry. Available in Rigid Curved, Rigid Straight, Flexible and Flexible Karman options, our cannulas-curettes are individually packaged and sterilized.

### Rigid Curved Cannula-Curette

|                       |                                            |
|-----------------------|--------------------------------------------|
| 022106 6mm (50/Pack)  | 022113 13mm - 3/8 Base (25/Pack)           |
| 022107 7mm (50/Pack)  | 022114 14mm - 1/2 Base (25/Pack)           |
| 022108 8mm (50/Pack)  | 022144 14mm - 1/2 Base w/Adapter (25/Pack) |
| 022109 9mm (50/Pack)  | 022145 14mm - 3/8 Base (25/Pack)           |
| 022110 10mm (50/Pack) | 022115 15mm - 1/2 Base (25/Pack)           |
| 022111 11mm (50/Pack) | 022116 16mm - 1/2 Base (25/Pack)           |
| 022112 12mm (50/Pack) |                                            |

### Rigid Straight Cannula-Curette

|                       |                                            |
|-----------------------|--------------------------------------------|
| 022206 6mm (50/Pack)  | 022213 13mm - 3/8 Base (25/Pack)           |
| 022207 7mm (50/Pack)  | 022214 14mm - 1/2 Base (25/Pack)           |
| 022208 8mm (50/Pack)  | 022244 14mm - 1/2 Base w/Adapter (25/Pack) |
| 022209 9mm (50/Pack)  | 022245 14mm - 3/8 Base (25/Pack)           |
| 022210 10mm (50/Pack) | 022215 15mm - 1/2 Base (25/Pack)           |
| 022211 11mm (50/Pack) | 022216 16mm - 1/2 Base (25/Pack)           |
| 022212 12mm (50/Pack) |                                            |



### Flexible/Karman Cannula-Curette

|                        |
|------------------------|
| 022004K 4mm (10/Pack)  |
| 022005K 5mm (10/Pack)  |
| 022006K 6mm (10/Pack)  |
| 022007K 7mm (10/Pack)  |
| 022008K 8mm (10/Pack)  |
| 022009K 9mm (10/Pack)  |
| 022010K 10mm (10/Pack) |
| 022011K 11mm (10/Pack) |
| 022012K 12mm (10/Pack) |

### Flexible Cannula-Curette with Dual Port Design

|                       |
|-----------------------|
| 022004 4mm (50/Pack)  |
| 022005 5mm (50/Pack)  |
| 022006 6mm (50/Pack)  |
| 022007 7mm (50/Pack)  |
| 022008 8mm (50/Pack)  |
| 022009 9mm (50/Pack)  |
| 022010 10mm (50/Pack) |
| 022011 11mm (50/Pack) |
| 022012 12mm (50/Pack) |



## DISPOSABLES



### MedGyn Laminaria

Laminaria (type of hygroscopic dilator) is a 100% natural sterile cervical dilator made from a sea-grown plant containing no synthetic materials. Safely utilized for over a century, Laminaria provides gentle, gradual dilation. After harvesting, the plant is dehydrated. Once inserted into the cervical canal, Laminaria absorbs liquid, swelling to several times its initial dried diameter, gently dilating the cervix.

#### Standard, approximate length, 60mm

- 021002 2mm Extra Small/Extra Thin (25/Box)
- 021003 3mm Small (25/Box)
- 021004 4mm Medium (25/Box)
- 021005 5mm Large/Thick (25/Box)
- 021006 6mm Extra Large/Extra Thick (25/Box)
- 021008 8mm Jumbo (25/Box)
- 021010 10mm Extra Jumbo (25/Box)

#### Extra Long, approximate length, 70mm

- 021014 4mm Long 70mm (25/Box)
- 021015 5mm Long 70mm (25/Box)
- 021016 6mm Long 70mm (25/Box)
- 021018 8mm Long 70mm (25/Box)
- 021003L Laminaria #3 Long-String
- 021004L Laminaria #4 Long-String
- 021005L Laminaria #5 Long-String
- 021006L Laminaria #6 Long-String
- 021008L Laminaria #8 Long-String

#### With buttons to prevent over-insertion

- 021022 2mm Extra Small with Buttons (12/Box)
- 021023 3mm Small with Buttons (12/Box)
- 021024 4mm Medium with Buttons (12/Box)
- 021025 5mm Large with Buttons (12/Box)
- 021026 6mm Extra Large with Buttons (12/Box)

#### Laminaria Inserter

- 031201 MedGyn Laminaria Inserter
- 031202 McMahon Laminaria Inserter



## MedGyn Pregnancy Test Kits

MedGyn One-Step Pregnancy Test Kit is an easy-to-use clinical pregnancy test that delivers fast and accurate results. Our tests have built-in test controls and highly visible color.

022704 MedGyn Pregnancy HCG Cassette Urine Test (25/Box)

022705 MedGyn Pregnancy HCG Combo Urine/Serum Test (25/Box)

022706 MedGyn Pregnancy HCG Dipstick Urine Test (25/Box)



## MedGyn Amnio Hook

MedGyn Amnio Hook is designed to perforate the amniotic membrane with a rigid talon tip.

020301 MedGyn Disposable Amnio Hook (50/Pack)



## MedGyn Umbilical Cord

MedGyn Umbilical Cord Clamp has grooved clamping and lock to secure the umbilical cord prior to cutting.

020302 MedGyn Disposable Umbilical Cord Clamp (100/Pack)



## MedGyn Uterine Manipulator

This device is used for laparoscopic tubal occlusion, diagnostic laparoscopy, mini laparotomy, salpingoplasty and fertility exams. Pear-shaped balloon provides optimal hold and a perfect seal at the internal OS. Smaller cannula size allows for easier insertion and comfort.

022701 MedGyn Uterine Manipulator(12/Box)

## DISPOSABLES



### MedGyn LED Speculum

MedGyn LED Speculum has a cordless LED light to improve visualization of the exam area.

020227 MedGyn Graves/Pederson LED Large (80/Box)

020226 MedGyn Graves/Pederson w/ LED Medium (100/Box)

020225 MedGyn Graves/Pederson w/ LED Small (100/Box)

020225STRL MedGyn Sterile Graves/Pederson w/ LED Small (100/Box)

020226STRL MedGyn Sterile Graves/Pederson w/ LED Medium (100/Box)

020227STRL MedGyn Sterile Graves/Pederson w/ LED Large (80/Box)



Large



Medium



Small



### MedGyn Vaginal Speculum

MedGyn Disposable Speculum is available in both Graves and Pederson styles and easily locks in a variety of positions to reduce examination time.

020205 MedGyn Disposable Graves Large (100/Box)

020204 MedGyn Disposable Graves Medium (100/Box)

020201 MedGyn Disposable Pederson Small (100/Box)



Large



Medium



Small



### MedGyn Balloon Catheter

MedGyn Postpartum Balloon Catheter is a disposable, multiple lumen catheter attached to an inflatable balloon system designed to provide tamponade for controlling hemorrhage from the uterus and vagina.

022703 MedGyn Balloon Catheter





## MedGyn HSG Catheters

MedGyn HSG Catheter is used to inject contrast media or saline into the uterus and fallopian tubes for hysterosalpingography and sonohysterography procedures. It has a latex-free balloon catheter with a rigid insertion sheath.

022726 MedGyn HSG Catheters - 5Fr (10/Box)

022728 MedGyn HSG Catheters - 7Fr (10/Box)



## MedGyn HVP-100 Hand-Held Vacuum Delivery Pump and Vacuum Assisted Delivery Cups

MedGyn Hand-Held Vacuum Delivery Pump may be used to safely assist medical personnel during vaginal deliveries. MedGyn Vacuum Delivery Cups offer optimal adhesion to ensure proper vacuum levels.

020322 MedGyn Vacuum Delivery Cups 50mm with Release Valve

020323 MedGyn Vacuum Delivery Cups 60mm with Release Valve

020316 MedGyn Vacuum Delivery Cups 70mm with Release Valve

020320 Vacuum Delivery Cup 50mm

020321 Vacuum Delivery Cup 60mm

016160 MedGyn HVP-100 Hand-Held Vacuum Delivery Pump

## DISPOSABLES



### MedGyn Endocervical Curette

The MedGyn Endocervical Curette provides an efficient sampling of the cervical canal to test for abnormal, precancerous conditions, or cervical cancer.

The Kevorkian-Young endocervical curette is a sterile, single-use device with a basket and a square tip.



022732 MedGyn Disposable Endocervical Curette (25/Box)




### MedGyn IUI Catheters

MedGyn Insemination Catheter is a sterile, flexible catheter designed for use in intrauterine insemination. Tip of catheter is rounded and beveled for patient comfort.

022723 MedGyn Intrauterine Insemination Catheter - Regular/Straight (25/Box)

022724 MedGyn Intrauterine Insemination Catheter - Mini (25/Box)

022725 MedGyn Intrauterine Insemination Catheter - Curved (25/Box)



### MedGyn Ezinject Catheter

A simple and effective way to allow visibility of the uterine cavity.

Balloonless catheter allows for a more complete view

Long tip at distal end improves placement in the cervical canal

022797 MedGyn Ezinject Catheter (10/box)







## MedGyn Uterine Sounds

MedGyn Disposable Uterine Sounds are made of malleable plastic with centimeter graduations for accurate measurement of uterine cavity.

020250 MedGyn Disposable Uterine Sounds O.D. 3mm 10" Length with Marking (25/Box)



## MedGyn Cell Sweep®

The MedGyn Cell Sweep® makes pap smears easier than ever by simultaneously sampling both the endocervix and ectocervix, eliminating the need for both a cytobrush and spatula.

- Optimal cell harvesting and transfer
- Separation of cells in one-slide technique
- Minimizes false negatives, red blood cells, and blood artifacts

- 022362 MedGyn Cell Sweep® Slide Holder (100/pack)
- 022361 MedGyn Cell Sweep® Non-Sterile (100/pack)
- 022361K MedGyn Cell Sweep® Non-Sterile with Kit (Slide and Slide Holder) (100/pack)
- 022363 MedGyn Cell Sweep® Sterile (100/pack)
- 022363K MedGyn Cell Sweep® Sterile with Kit (Slide and Slide Holder) (100/pack)



## MedGyn Cytobrush and Pap Cell

- 022710 MedGyn Cytobrush (100/Box)
- 022364 MedGyn Pap Cell (50/Box)

## DISPOSABLES



### MedGyn Endosampler™

The unique design of the MedGyn Endosampler™ allows for superior accuracy and a more reliable sample. Available with one or two holes for the right amount of suction for each patient.

- 3mm curette and 10cc syringe
- Highly-visible graduated markings
- Semi-rigid curette
- Curved for improved manipulation
- Prevent back-flow with syringe lock-spring
- Sharp slot in curette tip
- Closed system prevents sample contamination



Endosampler™ EndosamplerII™

022720 MedGyn Endosampler™ (25/Box)

022722 MedGyn EndosamplerII™ (25/Box)



### MedGyn Pipette™

The MedGyn Pipette™, Pipette II™, and Pipette IV™ are unique, sterile, disposable cannulas for sampling the endometrium.

- Semi-rigid and tapered for easy insertion and safety
- Flexible enough collect tissue from all areas

022721 MedGyn Pipette 3.1mm O.D. (25/box)

022734 MedGyn Pipette II 3.1mm O.D. (25/box)

022730 MedGyn Pipette IV 3.1mm O.D. (25/box)



Pipette™



Pipette II™



Pipette IV™



### MedGyn Mucous Sampler

The MedGyn Mucous Sampler is a single-use, sterile, endocervical aspirator.

- 23 cm long, clear, flexible polypropylene tube
- 2.2mm diameter circular opening
- Distal end rounded for safety



022729 MedGyn Disposable Mucous Sampler (25/pack)



## MedGyn Monsel's Solution

Ferric Subsulfate Solution 20% is an effective hemostatic solution for use after cervical biopsies.

022610 Monsel's Solution - 8mL bottles  
(12/Box)

022612 Monsel's Solution - 500 mL bottle



## MedGyn's Disposable Derma Curette

022797-3 Disposable Derma- Curette 3.0mm- STERILE

022797-4 Disposable Derma- Curette 4.0mm- STERILE

022797-7 Disposable Derma- Curette 7.0mm- STERILE



## MedGyn's Disposable Derma (Keyes Biopsy) Punch

MedGyn's latest products include a disposable derma (keyes biopsy) punch made of quality stainless steel and a ribbed sturdy handle for the ejection of lodged specimen.

022727-1 Disposable Derma (Keyes Biopsy) Punch 1mm- STERILE

022727-2 Disposable Derma (Keyes Biopsy) Punch 2mm- STERILE

022727-3 Disposable Derma (Keyes Biopsy) Punch 3mm- STERILE

022727-4 Disposable Derma (Keyes Biopsy) Punch 4mm- STERILE

022727-5 Disposable Derma (Keyes Biopsy) Punch 5mm- STERILE



## MedGyn's Disposable Scalpel

MedGyn's latest products include disposable scalpels made of quality stainless steel blades with sturdy, non-slip plastic handles.

022798-11 Disposable Scalpel #11- STERILE

022798-15 Disposable Scalpel #15- STERILE



## MedGyn Word Bartholin's Catheter

For use following an incision and drainage of a Bartholin's cyst or abscess. Catheter kit includes balloon catheter, syringe for balloon inflation and scalpel.

022719 MedGyn Word Bartholin's Catheter (6/Box)



## Other

040065 Video Print Paper - 110mm x 20mm rolls  
(10/Box)

## FORCEPS



### MedGyn Placenta Forceps

#### Hern-Barret Placenta Forceps

Slightly curved placenta forceps with oval jaws.

031151 Hern-Barret Placenta Forceps 13" (330mm)  
16mm Serrated Jaws

031152 Hern-Barret Placenta Forceps 13" (330mm)  
16mm Smooth Jaws



#### Kelly Placenta Forceps

Slightly curved placenta forceps with serrated jaws.

031140 Kelly Placenta Forceps 12.5" (318mm)  
19mm Serrated Jaws



#### Hern-Kelly Placenta Forceps

Slightly curved placenta forceps with serrated jaws.

031145 Hern-Kelly 13" (330mm)  
19mm Serrated Jaws



### MedGyn Myoma Grasping Forceps

#### Myoma Grasping Forceps

Ideal for transcervical management of polyps and pedunculated myoma.

031210 Myoma Grasping Forceps 13" (330mm) 7mm Jaws

031211 Myoma Grasping Forceps 13" (330mm) 11mm  
Jaws

031212 Myoma Grasping Forceps 13" (330mm) 17mm  
Jaws



### MedGyn Tissue Forceps

#### Ochsner Tissue Forceps

031169 Ochsner Tissue Forceps





## Boston Forceps

Straight forceps with large oval jaws and coarse teeth.

031156 Boston Evacuation Forceps 13½" (343mm)  
27mm Jaws - High Fulcrum

031157 Boston Evacuation Forceps 13½" (343mm)  
27mm Jaws - Low Fulcrum



## Meltzer Forceps

031143 Meltzer Evacuation Forceps 13" (330mm)  
19mm Jaws



## McClintock Forceps

Slightly curved with fine-tooth jaws.


031122 McClintock Ovum Forceps 9½"  
(241mm) 14mm Jaws



## Fink Forceps

Slightly curved with jaws and concave tip.

031176 Fink Ovum Forceps 10½" (267mm)  
13mm Jaws



## Lisa Forceps

031173 Lisa Forceps 10" (254mm)

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## FORCEPS



### Stubblefield Forceps

Slightly curved, very long forceps with oval jaws.

031144 Stubblefield Forceps 15" (381mm) 14mm Jaws



### Agadi Evacuation Forceps

Smaller forceps available in both straight and curved versions with oval jaws.

031170 Agadi Evacuation Forceps 9¾" (248mm) Curved

031171 Agadi Evacuation Forceps 9½" (241mm) Straight



### Blumenthal Forceps

Slightly curved forceps with large oval jaws and coarse teeth.

031141 Blumenthal Evacuation Forceps 13" (330mm)  
20mm Jaws



### Saenger Forceps

Forceps with unique inner-diameter teeth within oval jaws.

031163 Saenger Ovum Forceps 11¼" (286mm) Curved

031164 Saenger Ovum Forceps 10¾" (273mm) Straight  
with Ratchet

031165 Saenger Ovum Forceps 10½" (267mm) Curved  
with Ratchet

031172 Saenger Ovum Forceps 13" (330mm) Straight  
with Ratchet




### Jacobson Forceps

Slightly curved forceps, with coarse teeth, oval jaw and enlarged ring handles. With concave tip.

031142 Jacobson Evacuation Forceps 13" (330mm)  
Large with 20mm Jaws





## Hern-Patterson Forceps

Straight forceps with coarse teeth, oval jaws and ratchet.

031139 Hern Patterson Evacuation Forceps 13" (330mm) 17mm Jaws




## Hern Van-Lyth Forceps

Curved forceps with fine teeth.

031166 Hern Van Lyth Evacuation Forceps 9½" (241mm) 9mm Jaws

031167 Hern Van Lyth Evacuation Forceps 10" (254mm) 11mm Jaws



## Agadi-Bierer Forceps

031153 Agadi-Bierer Forceps 13" (330mm) - 14mm without Ratchet

031154 Agadi-Bierer Forceps 13" (330mm) - 14mm with Ratchet



## Ramplly Ring Forceps

031007 Ramplly with Oval Jaws 10" (254mm) Straight



## West-Gluck Forceps

031168 West-Gluck Forceps 11" (279mm)

## FORCEPS



### MedGyn Ovum Forceps

#### Sopher Ovum Forceps

Slightly curved forceps with oval, serrated jaws.

- 031123 Sopher Ovum Forceps 11" (279mm) 12mm Jaws
- 031123R Sopher Ovum Forceps 11" (279mm) 12mm Jaws with Ratchet
- 031126 Sopher Ovum Forceps 11" (279mm) 14mm Jaws
- 031149 Sopher Ovum Forceps 11" (279mm) 14mm Jaws with Ratchet
- 031124 Sopher Ovum Forceps 13" (330mm) 12mm Jaws
- 031125 Sopher Ovum Forceps 13" (330mm) 12mm Jaws with H Ratchet
- 031127 Sopher Ovum Forceps 13" (330mm) 14mm Jaws
- 031128 Sopher Ovum Forceps 13" (330mm) 14mm Jaws with Ratchet



#### Bierer Ovum Forceps

Slightly curved forceps with oval, coarse, serrated jaws.

- 031131 Bierer Ovum Forceps 11" (279mm) 16mm Jaws
- 031132 Bierer Ovum Forceps 11" (279mm) 16mm Jaws with Ratchet
- 031135 Bierer Ovum Forceps 11" (279mm) 19mm Jaws
- 031136 Bierer Ovum Forceps 11" (279mm) 19mm Jaws with Ratchet
- 031133 Bierer Ovum Forceps 13" (330mm) 16mm Jaws
- 031134 Bierer Ovum Forceps 13" (330mm) 16mm Jaws with Ratchet
- 031137 Bierer Ovum Forceps 13" (330mm) 19mm Jaws
- 031138 Bierer Ovum Forceps 13" (330mm) 19mm Jaws with Ratchet



#### Sopher Original Forceps

Slightly curved forceps with oval jaw. Concave tip design to prevent grasping of fundus.

- 031129 Sopher Original Forceps 10½" (267mm) Small Jaws
- 031130 Sopher Original Forceps 10½" (267mm) Medium Jaws



#### Hern Evacuation Forceps

Slightly curved forceps with coarse teeth, oval jaws.

- 031146 Hern Evacuation Forceps 13" (330mm) Small 12mm Jaws
- 031147 Hern Evacuation Forceps 13" (330mm) Medium 15mm Jaws
- 031148 Hern Evacuation Forceps 13" (330mm) Large 19mm Jaws



- 030316 Spiral Endocervical Curette Blue Titanium
- 030301 Kevorkian-Young Malleable Curette w/o Basket 12" (305mm)
- 030302 Kevorkian-Young Malleable Curette with Basket 12" (305mm)
- 030303 Kevorkian-Young Rigid Curette w/o Basket 12" (305mm)
- 030304 Kevorkian-Young Rigid Curette with Basket 12" (305mm)
- 030303-H Kevorkian-Young Rigid Curette w/o Basket (Square Tip, Gold Handle, Heavy) 12" (305mm)
- 030304-H Kevorkian-Young Rigid Curette with Basket (Square Tip, Gold Handle, Heavy) 12" (305mm)
- 030315 Kevorkian-Young Rigid Curette w/o Basket (Square Tip, Titanium) 12" (305mm)
- 030311 Kevorkian-Young Rigid Curette with Basket (Square Tip, Titanium) 12" (305mm)
- 030306 Kevorkian-Young Rigid Curette with Basket (Round Tip, Gold Handle, Heavy) 12" (305mm)
- 030310 Kevorkian-Young Rigid Curette with Basket (Round Tip, Titanium) 12" (305mm)

## Townsend Curette

Round, tapered tip allows for entry through stenotic os.



- 030305 Townsend Curette 12" (305mm)

## Metal Uterine Sounds

MedGyn reusable, graduated metal sounds made of highly malleable metal.



- 030704 Simpson Uterine Sound 3mm Diameter, Graduated in inches, 13" (330mm)
- 030705 Simpson Uterine Sound 3mm Diameter, Graduated in cm, 13" (330mm)

## Simpson Uterine Sound



- 030700 Sims Uterine Sound 3mm Diameter, Graduated in inches, 13" (330mm)
- 030701 Sims Uterine Sound 3mm Diameter, Graduated in cm, 13" (330mm)
- 030702 Sims Uterine Sound 2mm Diameter, Graduated in cm, 13" (330mm)
- 030703 Sims Uterine Sound 3mm Plain Probe, Malleable, 13" (330mm)

## METAL CURETTES



### Sims Curette

Metal curette with sharp blade on malleable or rigid shanks 11" (279mm).

030640 Sims Malleable Curette 8 Piece Set Sizes #00 - #6

030640 12 Sims Malleable Curette 12 Piece Set Sizes #00-10

030641 Sims Malleable Curette #00

030642 Sims Malleable Curette #0

030643 Sims Malleable Curette #1

030644 Sims Malleable Curette #2

030645 Sims Malleable Curette #3

030646 Sims Malleable Curette #4

030647 Sims Malleable Curette #5

030648 Sims Malleable Curette #6

030649 Sims Malleable Curette #7

030637 Sims Malleable Curette #8

030638 Sims Malleable Curette #9

030639 Sims Malleable Curette #10

030650 Sims Rigid Curette 8 Piece Set Sizes #00 - #6

030651 Sims Rigid Curette #00

030652 Sims Rigid Curette #0

030653 Sims Rigid Curette #1

030654 Sims Rigid Curette #2

030655 Sims Rigid Curette #3

030656 Sims Rigid Curette #4

030657 Sims Rigid Curette #5

030658 Sims Rigid Curette #6

030659 Sims Rigid Curette #7

030660 Sims Rigid Curette #8

030661 Sims Rigid Curette #9

030662 Sims Rigid Curette #10

## MedGyn Endocervical Curettes



### Kevorkian-Younger Curette

Metal curette available with or without basket.



Kevorkian-Younger w/o Basket



Kevorkian-Younger with Basket



## MedGyn Endometrial Curettes



### Novak Curette

Popular metal curette available in a 3 piece set or in individual sizes.

030600 Novak Curette 3 Piece Set 2mm - 4mm 9¾" (248mm)

030602 Novak Curette 2mm 9¾" (248mm)

030603 Novak Curette 3mm 9¾" (248mm)

030604 Novak Curette 4mm 9¾" (248mm)



### Evans Curette

030630 Evans Curette 6 Piece Set Sizes #0 - #5 11" (279mm)

030631 Evans Curette #0-4mm 11" (279mm)

030632 Evans Curette #1-6mm 11" (279mm)

030633 Evans Curette #2-8mm 11" (279mm)

030634 Evans Curette #3-10mm 11" (279mm)

030635 Evans Curette #4-12mm 11" (279mm)

030636 Evans Curette #5-14mm 11" (279mm)



### Randall Curette

030610 Randall Curette 3 Piece Set #2 - #4 9¾" (248mm)

030612 Randall Curette #1 9¾" (248mm)

030613 Randall Curette #2 9¾" (248mm)

030614 Randall Curette #3 9¾" (248mm)



### Bumm Curette

Extra-large metal curettes on malleable shanks

030620 4 Piece Set - Sizes #1 - #4 12" (305mm)

030621 Bumm - Size #1 12" (305mm)

030622 Bumm - Size #2 12" (305mm)

030623 Bumm - Size #3 12" (305mm)

030624 Bumm - Size #4 12" (305mm)



### Banjo/Hunter Curette

030619 Banjo/Hunter Curette Size 14 - 30mm



### Heaney Curette

030680 Heaney Curette Small 9½" (241mm)

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## DILATORS



### McMahan Cervico Meter

Angled, used for measuring dilation of cervix.

031300 McMahan Cervico Meter



### Mini Dilator Set

6 piece set (5 dilators & handle).

030779 6 Piece Mini Dilator Set:

1mm, 1.5mm, 2mm, 2.5mm, 3mm



030779BT Mini Dilator Set Blue Titanium



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## DILATORS



### Hegar Dilator

Stainless steel single-end dilator 7½" (191mm).

030841-1 Hegar Dilator 1mm  
030841-1.5 Hegar Dilator 1.5mm  
030841-2 Hegar Dilator 2mm  
030841-2.5 Hegar Dilator 2.5mm  
030841-3 Hegar Dilator 3mm  
030841-3.5 Hegar Dilator 3.5mm  
030841-4 Hegar Dilator 4mm  
030841-4.5 Hegar Dilator 4.5mm  
030841-5 Hegar Dilator 5mm  
030841-5.5 Hegar Dilator 5.5mm  
030842-6 Hegar Dilator 6mm  
030842-6.5 Hegar Dilator 6.5mm  
030842-7 Hegar Dilator 7mm  
030842-7.5 Hegar Dilator 7.5mm  
030842-8 Hegar Dilator 8mm  
030842-8.5 Hegar Dilator 8.5mm  
030842-9 Hegar Dilator 9mm  
030842-9.5 Hegar Dilator 9.5mm



030842-10 Hegar Dilator 10mm  
030842-10.5 Hegar Dilator 10.5mm  
030843-11 Hegar Dilator 11mm  
030843-11.5 Hegar Dilator 11.5mm  
030843-12 Hegar Dilator 12mm  
030843-12.5 Hegar Dilator 12.5mm  
030843-13 Hegar Dilator 13mm  
030843-13.5 Hegar Dilator 13.5mm  
030843-14 Hegar Dilator 14mm  
030843-14.5 Hegar Dilator 14.5mm  
030843-15 Hegar Dilator 15mm  
030843-15.5 Hegar Dilator 15.5mm  
030844-16 Hegar Dilator 16mm  
030844-16.5 Hegar Dilator 16.5mm  
030844-17 Hegar Dilator 17mm  
030844-17.5 Hegar Dilator 17.5mm  
030844-18 Hegar Dilator 18mm  
030844-18.5 Hegar Dilator 18.5mm



### Pratt Dilator Plastic

Light-weight, reusable 5 piece plastic dilator set.

030890 Pratt Dilator Plastic 5 Piece Set 5-6mm-  
13-14mm



### Os Cervical Finder Set

Made of soft, malleable Teflon™. 3 piece set  
(Os Finder, Canal Finder & Fundus Finder)

030803 Os Cervical Finder 3 Piece Set

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## DILATORS



### Hegar Dilator

Stainless steel double-end dilator 7½" (191mm).

030830 Hegar 8 Piece Dilator Set 3-4mm - 17-18mm

030831 Hegar Dilator 1-2mm

030832 Hegar Dilator 3-4mm

Hegar Dilator 5-6mm

030834 Hegar Dilator 7-8mm

030835 Hegar Dilator 9-10mm

030836 Hegar Dilator 11-12mm

030837 Hegar Dilator 13-14mm

030838 Hegar Dilator 15-16mm

030839 Hegar Dilator 17-18mm



### Kleegman Dilator

Double-end dilator available in stainless steel and Teflon™. (Stainless steel pictured) (1.0mm/2.0mm)

030804 Kleegman Stainless Steel Dilator 11" (279mm)

030805 Kleegman Teflon™ Dilator 12" (305mm)

## MedGyn Single-End Dilators



### Hern Dilator

Teflon single-end dilator with handle 14¾" (375mm).

030807 Hern 21 Piece Dilator Set Sizes 44 Fr - 104 Fr

030808 Hern Dilator 44 Fr

030809 Hern Dilator 47 Fr

030810 Hern Dilator 50 Fr

030811 Hern Dilator 53 Fr

030812 Hern Dilator 56 Fr

030813 Hern Dilator 59 Fr

030814 Hern Dilator 62 Fr

030815 Hern Dilator 65 Fr

030816 Hern Dilator 68 Fr

030817 Hern Dilator 71 Fr

030818 Hern Dilator 74 Fr

030819 Hern Dilator 77 Fr

030820 Hern Dilator 80 Fr

030821 Hern Dilator 83 Fr

030822 Hern Dilator 86 Fr

030823 Hern Dilator 89 Fr

030824 Hern Dilator 92 Fr

030825 Hern Dilator 95 Fr

030826 Hern Dilator 98 Fr

030827 Hern Dilator 101 Fr

030828 Hern Dilator 104 Fr

030829-107 Hern Dilator

107 Fr

030829-110 Hern Dilator

110 Fr

030829-113 Hern Dilator

113 Fr



### Hawkin Ambler

Stainless steel single-end, dilator 7¼" (184mm).

030780 16 Piece Dilator Set Sizes 3-6mm - 18-21mm

030781 Hawkin Ambler Dilator 3-6mm

030782 Hawkin Ambler Dilator 4-7mm

030783 Hawkin Ambler Dilator 5-8mm

030784 Hawkin Ambler Dilator 6-9mm

030785 Hawkin Ambler Dilator 7-10mm

030786 Hawkin Ambler Dilator 8-11mm

030787 Hawkin Ambler Dilator 9-12mm

030788 Hawkin Ambler Dilator 10-13mm

030789 Hawkin Ambler Dilator 11-14mm

030790 Hawkin Ambler Dilator 12-15mm

030791 Hawkin Ambler Dilator 13-16mm

030792 Hawkin Ambler Dilator 14-17mm

030793 Hawkin Ambler Dilator 15-18mm

030794 Hawkin Ambler Dilator 16-19mm

030795 Hawkin Ambler Dilator 17-20mm

030796 Hawkin Ambler Dilator 18-21mm

## MedGyn Double-End Dilators



### Pratt Dilator

Stainless steel double-end dilator 11½" (292mm).

030850 Pratt 8 Piece Dilator Set 13-15 Fr - 41-43 Fr

030850-ECO Pratt 8 Piece Dilator Set 13-15 Fr -41-43 Fr

030859 Pratt Dilator 9-11Fr

030851 Pratt Dilator 13-15 Fr

030852 Pratt Dilator 17-19 Fr

030853 Pratt Dilator 21-23 Fr

030854 Pratt Dilator 25-27 Fr

030855 Pratt Dilator 29-31 Fr

030856 Pratt Dilator 33-35 Fr

030857 Pratt Dilator 37-39 Fr

030858 Pratt Dilator 41-43 Fr



### Hank Dilator

Designed with guide rings for ease-of-use 11" (279mm)

030801 Hank 6 Piece Dilator Set 11-12 Fr - 21-22 Fr

030801-7 Hank 7 Piece Dilator Set 9-10 Fr - 21-22 Fr

030801-7/8 Hank Dilator 7-8 Fr

030801-9/10 Hank Dilator 9-10 Fr

030801-11/12 Hank Dilator 11-12 Fr

030801-13/14 Hank Dilator 13-14 Fr

030801-15/16 Hank Dilator 15-16 Fr

030801-17/18 Hank Dilator 17-18 Fr

030801-19/20 Hank Dilator 19-20 Fr

030801-21/22 Hank Dilator 21-22 Fr



### Pratt Extra Large Dilator

Pratt style in larger sizes 11½" (292mm).

030870 Pratt Extra Large 9 Piece Dilator Set 45-47 Fr - 77-79 Fr

030870-ECO Pratt Extra Large 9 Piece Dilator Set 45-47 Fr - 77-79 Fr

030871 Pratt Dilator Extra Large 45-47 Fr

030872 Pratt Dilator Extra Large 49-51 Fr

030873 Pratt Dilator Extra Large 53-55 Fr

030874 Pratt Dilator Extra Large 57-59 Fr

030875 Pratt Dilator Extra Large 61-63 Fr

030876 Pratt Dilator Extra Large 65-67 Fr

030877 Pratt Dilator Extra Large 69-71 Fr

030878 Pratt Dilator Extra Large 73-75 Fr

030879 Pratt Dilator Extra Large 77-79 Fr

030880 Pratt Extra Large, Light Weight 9 Piece Dilator Set 45-47 Fr - 77-79 Fr

030881 Pratt Extra Large, Light Weight 45-47 Fr

030882 Pratt Extra Large, Light Weight 49-51 Fr

030883 Pratt Extra Large, Light Weight 53-55 Fr

030884 Pratt Extra Large, Light Weight 57-59 Fr

030885 Pratt Extra Large, Light Weight 61-63 Fr

030886 Pratt Extra Large, Light Weight 65-67 Fr

030887 Pratt Extra Large, Light Weight 69-71 Fr

030888 Pratt Extra Large, Light Weight 73-75 Fr

030889 Pratt Extra Large, Light Weight 77-79 Fr

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## TENACULUMS



### Atraumatic Tenaculum

- 031001 Atraumatic Straight 1/4" (6mm) Jaws, 9" (229mm)
- 031002 Atraumatic Straight 1/8" (3mm) Jaws, 9" (229mm)
- 031003 Atraumatic Angled 1/4" (6mm) Jaws, 9 1/2" (241mm)



### Shah Tenaculum

- 031016 Shah Tenaculum 9 1/2" (241mm)



### Iris Hook

- 030406 Iris Hook 10" (254mm)



### Emmett Tenaculum Hooks

- 031401 Emmett Curved Hook 9 1/2" (241mm)
- 031402 Emmett Half-Curved Hook 9 3/4" (248mm)
- 031403 Emmett Hook L-Shaped 9" (229mm)
- 031404 Emmett Hook Acute Angle 9" (229mm)
- 031405 Emmett Hook Double Angled 9" (229mm)







## Schroeder Double-Tooth Tenaculum

- 031011 Schroeder Double-Tooth - 10" Straight (254mm)
- 031012 Schroeder Double-Tooth - 10" Curved (254mm)
- 031013 Schroeder Double-Tooth - 10" Side Curve (254mm)



## Jacobs Double-Tooth Tenaculum

- 031197 Jacobs Double-Tooth 8½" Straight - 2x2 Teeth (215mm)
- 031198 Jacobs Double-Tooth 8½" Side Curve - 2x2 Teeth (215mm)



## Bierer Tenaculum

- 031010 Bierer Angled 5x6 Teeth, 9" (229mm)
- 031033 Bierer Angled 5x6 Teeth, Economy 9" (229mm)



## Teale Tenaculum

- 031005 Teale Curved, 3x4 Teeth, 5mm Wide, 10" (254mm)



## Kelly Tenaculum

- 031017 Kelly Tenaculum 4X5 Teeth, 9" (229mm)

## TENACULUMS



### MedGyn Tenaculums

#### Schroeder Braun Single-Tooth Tenaculum

031004 Schroeder Braun Single-Tooth - 9¾" (248mm)

031019 Schroeder Braun Tenaculum Light-Touch Straight - 9¾" (248mm)

031034 Schroeder Braun Tenaculum Economy - 9¾" (248 mm)

031018 Braun Tenaculum Light-Touch Round - 7¼" (184mm)



#### Skene Single-Tooth Tenaculum

031006 Skene Single-Tooth-Side-Curve - 10" (254mm)



#### Duplay Single-Tooth Tenaculum

031008 Duplay Single-Tooth Double-Curve - 11" (279mm)



#### Barrett Tenaculum

031009 Barrett Straight - 7¼" (184mm)



#### Jacobson Double-Tooth Tenaculum

031199 Jacobson Double-Tooth 10" Straight - 2x2 Teeth with Larger Rings (254mm)

031200 Jacobson Double-Tooth 10" Side Curve - 2x2 Teeth with Larger Rings (254mm)





**Van Doren Biopsy Forceps**

030155 Van Doren Biopsy Forceps with Ratchet 10" (254mm)

030166 Van Doren Biopsy Forceps without Ratchet 10" (254mm)



**Schubert Biopsy Forceps**

030157 Schubert Biopsy Forceps 10" (254mm) Curved



**Wittner Biopsy Punch**

030117 Wittner Biopsy Punch 8½" (216mm)



**Rotating Biopsy Punch Options**

030778 Vulvar Keyes Biopsy Punch Set (Handle + 6 Tips)

030104 Rotating Tischler Biopsy Punch

030145 Rotating Tischler Biopsy Punch - Angled Down

030146 Rotating Tischler Biopsy Punch - Angled Up

030105 Rotating Baby Tischler Biopsy Punch

030105-DOWN Rotating Baby Tischler Biopsy Punch - Angled Down

030105-UP Rotating Baby Tischler Biopsy Punch - Angled Up

030106 Rotating Kevorkian Biopsy Punch

030107 Rotating Burke Biopsy Punch

030111-BB Rotating Shaft - Baby Tischler Tip Only

030111-BRK Rotating Shaft - Burke Tip Only

030111-KEV Rotating Shaft - Kevorkian Tip Only

030111-TIS Rotating Shaft - Tischler Tip Only

030147 Rotating Shaft - Tischler Tip - Up

030148 Rotating Shaft - Tischler Tip - Down

030108 Rotating Biopsy Punch Handle Only

030100 Rotating Shaft (with any tip)



## BIOPSY PUNCHES



### Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma.

Bite Size: 4.2mm x 2.3mm

- 030102-10 Baby Tischler Biopsy Punch 10 (254mm)
- 029811 Baby Tischler Biopsy Punch 10 (254mm) Angled Up
- 029812 Baby Tischler Biopsy Punch 10 (254mm) Angled Down
- 029813 Baby Tischler Biopsy Punch 10 (254mm) Angled Left
- 029814 Baby Tischler Biopsy Punch 10 (254mm) Angled Right
- 030102 Baby Tischler Biopsy Punch 8½ (216mm)
- 030140 Baby Tischler Biopsy Punch 8½ (216mm) Angled Up
- 030141 Baby Tischler Biopsy Punch 8½ (216mm) Angled Down
- 030142 Baby Tischler Biopsy Punch 8½ (216mm) Angled Left
- 030143 Baby Tischler Biopsy Punch 8½ (216mm) Angled Right



### Kevorkian Biopsy Punch

(Kevorkian-Youngs Biopsy Forceps displayed)

- 030103-10 Kevorkian Biopsy Punch 10 (254mm)
- 030103 Kevorkian Biopsy Punch 8½ (216mm)
- 030169 Kevorkian-Youngs Biopsy Forceps without Spring 8" (203mm)



### Tischler-Kevorkian Cervical Biopsy Punch

- 030116 Tischler-Kevorkian Biopsy Punch 10 (254mm)
- 030115 Tischler-Kevorkian Biopsy Punch 8½ (216mm)



### Burke Biopsy Punch

Open-handle design for optimal control with small cutting jaws and sharp teeth.

- 030109 Burke Biopsy Punch 8½



### Eppendorfer Biopsy Forceps

- 030110 Eppendorfer Biopsy Forceps 8" (203mm)

## BIOPSY PUNCHES



### Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

- 030173 Baby Tischler Biopsy Punch Titanium 10 (254mm)  
Heavy-Grip Handle
- 030121 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- 030134 Baby Tischler Biopsy Punch Titanium 10 (254mm)  
- Angled Up
- 030135 Baby Tischler Biopsy Punch Titanium 10 (254mm)  
- Angled Down
- 030161 Baby Tischler Biopsy Punch Titanium 10 (254mm)  
- Angled Left
- 030160 Baby Tischler Biopsy Punch Titanium 10 (254mm)  
- Angled Right
- 030119 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- 030132 Baby Tischler Biopsy Punch Titanium 8½ (216mm)  
- Angled Up
- 030133 Baby Tischler Biopsy Punch Titanium 8½ (216mm)  
-Angled Down
- 030159 Baby Tischler Biopsy Punch Titanium 8½ (216mm)  
- Left
- 030158 Baby Tischler Biopsy Punch Titanium 8½ (216mm)  
- Right



## MedGyn Stainless Steel

### Tischler Biopsy Punch

Classic Tischler bite designed for drawing in maximum tissue near Os with minimal trauma.

Bite Size 7mm x 3mm

- 030101-10 Tischler Biopsy Punch 10 (254mm)
- 029801 Tischler Biopsy Punch 10 (254mm) Angled Up
- 029802 Tischler Biopsy Punch 10 (254mm) Angled Down
- 029803 Tischler Biopsy Punch 10 (254mm) Angled Left
- 029804 Tischler Biopsy Punch 10 (254mm) Angled Right
- 030101 Tischler Biopsy Punch 8½ (216mm)
- 030113 Tischler Biopsy Punch 8½ (216mm) Angled Up
- 030114 Tischler Biopsy Punch 8½ (216mm) Angled Down
- 030136 Tischler Biopsy Punch 8½ (216mm) Angled Left
- 030137 Tischler Biopsy Punch 8½ (216mm) Angled Right

## BIOPSY PUNCHES



### MedGyn Titanium

#### Tischler Biopsy Punch

Classic Tischler bite, with durable coating.

Bite Size: 7mm x 3mm

030172 Tischler Biopsy Punch Titanium 10 (254mm) Heavy-Grip Handle

030120 Tischler Biopsy Punch Titanium 10 (254mm)

030124 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Up

030125 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Down

030162 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Left

030163 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Right

030118 Tischler Biopsy Punch Titanium 8½ (216mm)

030122 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Up

030123 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Down

030138 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Right

030139 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Left



#### Baby Tischler Biopsy Punch (Forceps Style)

Classic Baby Tischler bite with durable coating. Bite Size:

4.2mm x 2.3mm

030170 Baby Tischler Biopsy Punch Titanium Gold Handle (Forceps Style) 10" (254mm)



#### Oval Tischler Biopsy Punch

030171 Oval Tischler Biopsy Punch Titanium Heavy-Grip Handle 10 (254mm)



#### Kevorkian Biopsy Punch

Popular design with strong, narrow jaw with durable coating. Bite Size: 8mm x 3mm.

030174 Kevorkian Biopsy Punch Titanium Heavy-Grip Handle 10" (254mm)

030127 Kevorkian Biopsy Punch Titanium 10 (254mm)

030126 Kevorkian Biopsy Punch Titanium 8½ (216mm)





## MedGyn Blue Titanium

Classic Tischler bite with durable coating.  
Bite Size: 7mm x 3mm

030176 Tischler Biopsy Punch Blue Titanium  
Heavy-Grip Handle 10 (254mm)

030181 Tischler Biopsy Punch Blue Titanium 10  
(254mm)



## Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

030177 Baby Tischler Biopsy Punch Blue Titanium  
Heavy-Grip Handle 10 (254mm)

030129 Baby Tischler Biopsy Punch Blue Titanium  
8½ (216mm)

030180 Baby Tischler Biopsy Punch Blue Titanium  
10 (254mm)



## Baby Tischler Biopsy Punch (Forceps Style)

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

030179 Baby Tischler Biopsy Punch Blue Titanium  
Gold Handle (Forceps Style) 10" (254mm)



## Kevorkian Biopsy Punch

Popular design with strong, narrow jaw with durable coating and a heavy-grip handle. Bite Size: 8mm x 3mm

030178 Kevorkian Biopsy Punch Blue Titanium  
Heavy-Grip Handle 10 (254mm)



## Oval Tischler Biopsy Punch

Oval Tischler bite with durable coating and a heavy-grip handle.

030175 Oval Tischler Biopsy Punch Blue Titanium  
Heavy-Grip Handle 10 (254mm)

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## SPECULUMS

### Weighted Speculum



- 030921 Auward 2.0lbs, 75° angle,  
Blade 1½" x 3¼" (38mm x 83mm)
- 030922 Auward 2.5lbs, 75° angle,  
Blade 1½" x 3¼" (38mm x 83mm)
- 030923 Auward 3.0lbs, 75° angle,  
Blade 1½" x 3¼" (38mm x 83mm)
- 030924 Auward 2.5lbs, 45° angle,  
Blade 1½" x 4" (35mm x 102mm)
- 030926 Auward 2.5lbs, 90° angle,  
Blade 1" x 5½" (25mm x 140mm)

### Collins Speculum

The unique hands-free design offers the advantages of an open-sided speculum while permitting multiple positioning for optimal exposure with blades that open laterally instead of vertically.



- 030946 1¼" x 3¾" (32mm x 95mm)
- 030947 1½" x 4" (35mm x 102mm)
- 030948 1½" x 4½" (41mm x 114mm)


### Lateral Vaginal Wall Retractor

Used for the retraction of vaginal wall and made to fit inside any standard speculum to expand past blades of speculum. Allows for greater exposure of cervix for safer, more accurate procedures.




- 030970 Lateral Vaginal Wall Retractor Blue Titanium
- 030966 Open ½" (13mm), 10" (254mm),  
Blade Size 1" x 2½" (25mm x 64mm)
- 030967 Open 1" (25mm), 10" (254mm),  
Blade Size 1" x 2½" (25mm x 64mm)
- 030968 Wide Open, 10" (254mm),  
Blade Size 1" x 2½" (25mm x 64mm)
- 030969 Closed, 10" (254mm),  
Blade Size 1" x 2½" (25mm x 64mm)

### Collins Speculum



### Endocervical Speculum



- 030205 Endocervical Speculum Blue Titanium
- 030915 Kogan with Screw 2-4mm
- 030916 Kogan with Screw 4-8mm
- 030201 Kogan with Gauge & Ratchet 2-4mm
- 030202 Kogan with Gauge & Ratchet 4-8mm
- 030203 Kogan with Gauge & Ratchet 8-12mm
- 030204 Townsend Endocervical Speculum 8-12mm



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## SPECULUMS

### Klopper Speculum

Wider opening (4cm) and short blade



030931 1 3/8" x 3 1/2" (35mm x 89mm)

### Grundman Speculum

Wide open design with screw adjustment/lever



030949 Small 1 1/8" x 3 1/8" (34mm x 85mm)

030950 Medium 1 3/8" x 3 3/4" (35mm x 95mm)

030951 Large 1 1/2" x 4 3/4" (38mm x 121mm)

### Cusco Speculum

Unique handle-design speculum with center screw or side screw options.



030960 Small w/Center Screw 1 1/8" x 3 1/2" (34mm x 89mm)

030961 Medium w/Center Screw 1 3/8" x 4" (35mm x 102mm)


030962 Large w/Center Screw 1 3/8" x 4 1/2" (35mm x 114mm)

030963 Small w/Side Screw 1 1/8" x 3 1/2" (34mm x 89mm)

030964 Medium w/Side Screw 1 3/8" x 4" (35mm x 102mm)

030965 Large w/Side Screw 1 3/8" x 4 1/2" (35mm x 114mm)

### Klopper, Grundman, Cusco, For-Way Speculum



### Four-Way Speculum

The design of this speculum combines the characteristics of both a standard speculum and a vaginal wall retractor.



038095 4-Way Lateral Expander Speculum - Medium  
7/8" x 4 1/4" (22mm x 108mm)

038096 4-Way Lateral Expander Speculum - Large  
7/8" x 5 1/4" (22mm x 133mm)

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## SPECULUMS

### Pederson Vu-Max Speculum

Largest view area (7cm vs. standard 3cm) with durable coating.



030983 Medium 4½" (114mm)

030984 Large 5" (127mm)

### Pederson Speculum

Narrow blade design



030941 XNarrow ⅝" x 4" (17mm x 102mm)

030913 XSmall ½" x 2¾" (13mm x 70mm)

030906 Small ½" x 3" (13mm x 76mm)


030907 Medium ⅞" x 4" (22mm x 102mm)

030908 Large 1" x 4¾" (25mm x 121mm)

030914 XLarge 1" x 6½" (25mm x 165mm)

030917 Extra Long Snowman Clinton 1" x 6" (25mm x 152mm)

### Pederson Speculum



### Pederson Vu-More Speculum

Larger view area (4cm vs. standard 3cm) with durable coating.



030981 Medium 4½" (114mm)

030982 Large 5" (127mm)

### Pederson - Open Side Speculum

Narrow blade design with open side



030972 45° Medium - Right ⅞" x 4" (22mm x 102mm)

030973 45° Large - Right 1" x 4¾" (25mm x 121mm)

030942 XNarrow - Right ⅝" x 4" (17mm x 102mm)

030943 XNarrow - Left ⅝" x 4" (17mm x 102mm)

030939 Medium - Right ⅞" x 4" (22mm x 102mm)

030939-LT Medium - Left ⅞" x 4" (22mm x 102mm)

030908R Large - Right 1" x 4¾" (25mm x 121mm)

030908-LT Large - Left 1" x 4¾" (25mm x 121mm)

030914-L-XL XLarge - Left 1" x 6½" (25mm x 165mm)

030914-R-XL XLarge - Right 1" x 6½" (25mm x 165mm)

(Surgical Instruments, Medical Equipments, Hospital Supplies)



www.primeinstrumentintl.com

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

## SPECULUMS

### Graves Vu-More Speculum

Wider blade design with larger view area (4cm vs. standard 3cm) with durable coating.



030985 Medium 4½" (114mm)  
030986 Large 5" (127mm)


### Graves Vu-More Speculum

Wider blade design with larger view area (4cm vs. standard 3cm).



030952 Small 3¾" (95mm)  
030953 Medium 4½" (114mm)  
030954 Large 5" (127mm)  
030955 XLong 6" (152mm)  
030956 XLong 6½" (161mm)  
030957 XLong 6¾" (171mm)  
030958 XLong 7" (178mm)

### Graves Speculum



### Graves Vu-Max Speculum

Wider blade design with largest view area (7cm vs. standard 3cm) with durable coating.



030987 Medium 4½" (114mm)  
030988 Large 5" (127mm)

### Graves Vu-Max Speculum

Wider blade design with largest view area (7cm vs. standard 3cm).



030975 Medium 4½" (114mm)  
030976 Large 5" (127mm)  
030977 XLong 6" (152mm)  
030978 XLong 6½" (161mm)  
030979 XLong 6¾" (171mm)  
030980 XLong 7" (178mm)

### Weisman Graves Speculum

Open-side speculum



030911-SM Small - Right ¾" x 3" (10mm x 76mm)  
030911 Medium - Right 1⅜" x 4" (35mm x 102mm)  
030911-LG Large - Right 1⅜" x 4¾" (35mm x 121mm)  
030911-XL XLarge - Right 1½" x 6½" (38mm x 165mm)  
030912 Medium - Left 1⅜" x 4" (35mm x 102mm)  
030912-LG Large - Left 1⅜" x 4¾" (35mm x 121mm)  
030912-XL XLarge - Left 1½" x 6½" (38mm x 165mm)






## Baby Allis Tissue Forceps

- 031111 Baby Allis Forceps 5½" (140mm)  
4x5 Teeth-Extra Delicate Jaws
- 031112 Baby Allis Forceps 5" (127mm)  
4x5 Teeth-Extra Delicate Jaws



## Adson Tissue Forceps

- 031196 Adson Tissue Forceps 4¾" (121mm)  
1x2 Teeth



## Tissue Forceps

- 031180 Tissue Forceps 8" (203mm) Small 1x2 Teeth
- 031181 Tissue Forceps 10" (254mm) Large 1x2 Teeth




## MedGyn Tissue Forceps

### Allis Tissue Forceps

- 031107 Allis Tissue Forceps 4" (104mm) 4x5 Teeth
- 031105 Allis Tissue Forceps 5" (127mm) 4x5 Teeth
- 031108 Allis Tissue Forceps 6" (152mm) 4x5 Teeth
- 031109 Allis Tissue Forceps 6" (152mm) 5x6 Teeth
- 031110 Allis Tissue Forceps 7½" (190mm) 5x6 Teeth
- 031103 Allis Tissue Forceps 9½" (241mm) 5x6 Teeth
- 031106 Allis Tissue Forceps 10" (254mm) 5x6 Teeth
- 031104 Allis Tissue Forceps 6" (152mm) 4x5 Teeth  
with Angular Jaws
- 031101 Allis Tissue Forceps 6" (152mm) 5x6 Teeth  
with Angular Jaws
- 031102 Allis Tissue Forceps 7¼" (184mm) 5x6 Teeth  
with Angular Jaws

## OB/GYN

## DISPOSABLES



### SKLAR CERVIX BRUSH CERVICAL CELL SAMPLER

Sklar's Cervix Brush is made of polyethylene and has a total length of 20 centimeters. The top portion uses a soft flexible brush to obtain cell samples. The brush consists of 57 semi-circular plastic bristles of various lengths. The shape is such that the top edges of the brush follow the contours of the cervix. The longer middle bristles reach deep into the endocervical canal. The shorter bristles touch both the ectocervical area and transformation zone at the same time.

### SKLAR CERVIX BRUSH

#### STERILE

| PRODUCT NUMBER | QUANTITIES |
|----------------|------------|
| 96-4140        | 100/bx     |

#### NON-STERILE

| PRODUCT NUMBER | QUANTITIES |
|----------------|------------|
| 96-4140NS      | 100/bx     |




# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Nabelschnurklemmen  
Umbilical cord clamps

Clamps pour le cordon ombilical  
Clamps umbilicales  
Pinze ombelicali



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Abortuszangen – Abortusküretten  
Abortion forceps – Ovum curettes

Pinces à faux germes – Curette d'avortement  
Pinzas para aborto – Curetas para aborto  
Pinze per placenta ed aborto – Curette per placenta ed aborto



SAENGER


270 mm – 10 3/4"

tk 35720-27


tk 35721-27



1:1



1:1



1:3

WALLICH

420 mm – 16 1/2"

tk 35735-02

stumpf  
blunt  
mousse  
roma  
smussa



1:3

PESTALOZZA

300 mm – 12"

tk 35732-30

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Abortuszangen – Abortusküretten  
Abortion forceps – Ovum curettes

Pinces à faux germes – Curette d'avortement  
Pinzas para aborto – Curetas para aborto  
Pinze per placenta ed aborto – Curette per placenta ed aborto



WINTER

|               | Fig. 1      | Fig. 2      | Fig. 3      |
|---------------|-------------|-------------|-------------|
| 280 mm – 11 " | tk 35701-28 | tk 35702-28 | tk 35703-28 |

gerade  
straight  
droit  
recta  
rette

|             |             |             |
|-------------|-------------|-------------|
| tk 35711-28 | tk 35712-28 | tk 35713-28 |
|-------------|-------------|-------------|

gebogen  
curved  
courbe  
curva  
curve



C UZZI

|               | 14.5 mm     | 18.0 mm     | 21.0 mm     |
|---------------|-------------|-------------|-------------|
| 300 mm – 12 " | tk 35730-01 | tk 35730-02 | tk 35730-03 |

stumpf, Rückseite gerieft  
blunt, back serrated  
mousse, dos strié  
romo, dorso estriado  
smusso, dorso zigrinato

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Geburtshilfe  
Obstetrics

Obstétrique  
Obstetricia  
Ostetricia



SMELLIE

250 mm – 10" **tk 35741-25**

Perforatorium  
Perforator  
Perforateur  
Perforador  
Perforatore



BRAUN

420 mm – 16½" **tk 35750-42**

Kranioklast, ein Maul gefenstert  
Cranioclast, one jaw fenestrated  
Cranioclasse, une partie fenestrée  
Craneolasto, una boca fenestrada  
Cranioclasse, un lato fenestrato



TARNIER

460 mm – 17¼" **tk 35748-46**

Basiotrib  
Basiotribe  
Basiotribe  
Basiotribo  
Basiotribo

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Geburtshilfe  
*Obstetrics*

Obstétrique  
*Obstetricia*  
Ostetricia



SELLHEIM 1:2

280 mm – 11 " **tk 35664-28**

Geburtslöffel  
*Obstetric lever*

Cuillère obstétricale  
*Cuchara obstétrica*  
Cucchaio obstetrico



BRAUN 1:2

325 mm – 12 3/4 " **tk 35745-30**

Dekapitationshaken  
*Decapitating hook*

Crochet de décapitation  
*Granchio de decapitación*  
Uncino di decapitazione



GREEN-  
ARMYTAGE 1:2

210 mm – 8 1/4 " **tk 35665-21**

Uterusklemmzange  
*Uterine elevating forceps*

Pince utérine  
*Pinza para el útero*  
Pinza uterina

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Geburtszangen  
Obstetrical forceps

Forceps obstétriques  
Forceps obstétricos  
Forcipi



ELLIOT

|               |             |
|---------------|-------------|
| 320 mm - 12½" | tk 35620-32 |
| 380 mm - 15"  | tk 35620-38 |



MCLEAN

|               |             |
|---------------|-------------|
| 360 mm - 14¼" | tk 35622-36 |
|---------------|-------------|




KIELLAND

|               |             |
|---------------|-------------|
| 410 mm - 16¼" | tk 35635-41 |
|---------------|-------------|



KIELLAND-  
LUIKART

|              |             |
|--------------|-------------|
| 400 mm - 16" | tk 35637-40 |
|--------------|-------------|



BILL

tk 35640-00

Zughebel, passend für alle Geburtszangen mit gehorntem Standard-Griff  
*Traction handle for all the forceps with horned standard handle*

Dispositif de traction pour tous les forceps avec manche normal à cornes  
*Dispositivo de tracción para todos los forceps con mango standard cornudo*  
Impugnatura di trazione adatto a tutti i forcipi con manico standard con corna



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Geburtszangen  
Obstetrical forceps

Forceps obstétriques  
Forceps obstétricos  
Forcipi



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Geburtszangen  
Obstetrical forceps

Forceps obstétriques  
Forceps obstétricos  
Forcipi



SIMPSON

1:3

300 mm - 12" tk 35600-30



SIMPSON-  
BRAUN

330 mm - 13" tk 35601-33  
360 mm - 14 1/4" tk 35601-36



NAEGELE

360 mm - 14 1/4" tk 35605-36  
400 mm - 16" tk 35605-40



PIPER

1:3

440 mm - 17 1/4" tk 35607-44

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scharfe Löffel  
Uterine spoons

Curettes tranchantes  
Cucharillas cortantes  
Cucchiai taglienti



| SIMON                                                                                       |  |  |  |
|---------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
|                                                                                             | Fig. 1                                                                            | Fig. 2                                                                            | Fig. 3                                                                             |
| 250 mm - 9 3/4"                                                                             | tk 24445-01                                                                       | tk 24445-02                                                                       | tk 24445-03                                                                        |
| scharf, fest<br>sharp, rigid<br>tranchante, rigide<br>cortante, rigida<br>taglienti, rigido |                                                                                   |                                                                                   |                                                                                    |
|                                                                                             |  |  |  |
|                                                                                             | Fig. 4                                                                            | Fig. 5                                                                            | Fig. 6                                                                             |
|                                                                                             | tk 24445-04                                                                       | tk 24445-05                                                                       | tk 24445-06                                                                        |




| SCHRÖDER     |  |  |  |
|--------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
|              | Fig. 1                                                                              | Fig. 2                                                                              | Fig. 3                                                                               |
| 300 mm - 12" | tk 24448-01                                                                         | tk 24448-02                                                                         | tk 24448-03                                                                          |
|              |  |  |  |
|              | Fig. 4                                                                              | Fig. 5                                                                              | Fig. 6                                                                               |
|              | tk 24448-04                                                                         | tk 24448-05                                                                         | tk 24448-06                                                                          |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterusküretten  
Uterine curettes

Curettes utérines  
Curetas uterinas  
Curette uterine




# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterusküretten  
Uterine curettes

Curettes utérines  
Curetas uterinas  
Curette uterine



| RECAMIER      |      |      | scharf, biegsam<br><i>sharp, malleable</i><br>tranchante, malléable<br><i>cortante, maleable</i><br>tagliente, flessibile | stumpf, biegsam<br><i>blunt, malleable</i><br>mousse, malléable<br><i>roma, maleable</i><br>smussa, flessibile | scharf, fest<br><i>sharp, rigid</i><br>tranchante, rigide<br><i>cortante, rigida</i><br>tagliente, rigida | stumpf, fest<br><i>blunt, rigid</i><br>mousse, rigid<br><i>roma, rigida</i><br>smussa, rigida |
|---------------|------|------|---------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|
| y=340 mm-13½" | Fig. | x=mm |                                                                                                                           |                                                                                                                |                                                                                                           |                                                                                               |
|               | 13   | 25   | tk 35004-13                                                                                                               | tk 35005-13                                                                                                    | tk 35006-13                                                                                               | tk 35007-13                                                                                   |
|               | 14   | 30   | tk 35004-14                                                                                                               | tk 35005-14                                                                                                    | tk 35006-14                                                                                               | tk 35007-14                                                                                   |
|               | 15   | 35   | tk 35004-15                                                                                                               | tk 35005-15                                                                                                    | tk 35006-15                                                                                               | tk 35007-15                                                                                   |
|               | 16   | 40   | tk 35004-16                                                                                                               | tk 35005-16                                                                                                    | tk 35006-16                                                                                               | tk 35007-16                                                                                   |




# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterusküretten  
Uterine curettes

Curettes utérines  
Curetas uterinas  
Curette uterine



scharf, biegsam  
sharp, malleable  
tranchante, malléable  
cortante, maleable  
tagliente, flessibile

stumpf, biegsam  
blunt, malleable  
mousse, malléable  
roma, maleable  
smussa, flessibile

scharf, fest  
sharp, rigid  
tranchante, rigide  
cortante, rigida  
tagliente, rigida

stumpf, fest  
blunt, rigid  
mousse, rigid  
roma, rigida  
smussa, rigida

RECAMIER

| y=260mm-10 3/4" | Fig. | x=mm |             |             |             |             |
|-----------------|------|------|-------------|-------------|-------------|-------------|
|                 | 00   | 5    | tk 35000-20 | tk 35001-20 | tk 35002-20 | tk 35003-20 |
|                 | 0    | 6    | tk 35000-00 | tk 35001-00 | tk 35002-00 | tk 35003-00 |
|                 | 1    | 7    | tk 35000-01 | tk 35001-01 | tk 35002-01 | tk 35003-01 |
|                 | 2    | 8    | tk 35000-02 | tk 35001-02 | tk 35002-02 | tk 35003-02 |
|                 | 3    | 9    | tk 35000-03 | tk 35001-03 | tk 35002-03 | tk 35003-03 |
|                 | 4    | 11   | tk 35000-04 | tk 35001-04 | tk 35002-04 | tk 35003-04 |
|                 | 5    | 12   | tk 35000-05 | tk 35001-05 | tk 35002-05 | tk 35003-05 |
|                 | 6    | 14   | tk 35000-06 | tk 35001-06 | tk 35002-06 | tk 35003-06 |
|                 | 7    | 15   | tk 35000-07 | tk 35001-07 | tk 35002-07 | tk 35003-07 |
|                 | 8    | 16   | tk 35000-08 | tk 35001-08 | tk 35002-08 | tk 35003-08 |
|                 | 9    | 19   | tk 35000-09 | tk 35001-09 | tk 35002-09 | tk 35003-09 |
|                 | 10   | 20   | tk 35000-10 | tk 35001-10 | tk 35002-10 | tk 35003-10 |
|                 | 11   | 21   | tk 35000-11 | tk 35001-11 | tk 35002-11 | tk 35003-11 |
|                 | 12   | 23   | tk 35000-12 | tk 35001-12 | tk 35002-12 | tk 35003-12 |

| y=310mm-12 1/4" | Fig. | mm |             |             |             |             |
|-----------------|------|----|-------------|-------------|-------------|-------------|
|                 | 00   | 5  | tk 35004-20 | tk 35005-20 | tk 35006-20 | tk 35007-20 |
|                 | 0    | 6  | tk 35004-00 | tk 35005-00 | tk 35006-00 | tk 35007-00 |
|                 | 1    | 7  | tk 35004-01 | tk 35005-01 | tk 35006-01 | tk 35007-01 |
|                 | 2    | 8  | tk 35004-02 | tk 35005-02 | tk 35006-02 | tk 35007-02 |
|                 | 3    | 9  | tk 35004-03 | tk 35005-03 | tk 35006-03 | tk 35007-03 |
|                 | 4    | 11 | tk 35004-04 | tk 35005-04 | tk 35006-04 | tk 35007-04 |
|                 | 5    | 12 | tk 35004-05 | tk 35005-05 | tk 35006-05 | tk 35007-05 |
|                 | 6    | 14 | tk 35004-06 | tk 35005-06 | tk 35006-06 | tk 35007-06 |
|                 | 7    | 15 | tk 35004-07 | tk 35005-07 | tk 35006-07 | tk 35007-07 |
|                 | 8    | 16 | tk 35004-08 | tk 35005-08 | tk 35006-08 | tk 35007-08 |
|                 | 9    | 19 | tk 35004-09 | tk 35005-09 | tk 35006-09 | tk 35007-09 |
|                 | 10   | 20 | tk 35004-10 | tk 35005-10 | tk 35006-10 | tk 35007-10 |
|                 | 11   | 21 | tk 35004-11 | tk 35005-11 | tk 35006-11 | tk 35007-11 |
|                 | 12   | 23 | tk 35004-12 | tk 35005-12 | tk 35006-12 | tk 35007-12 |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Salpingograph  
Salpingograph

Salpingographe  
Salpingógrafo  
Salpingografo



SCHULTZE

295 mm – 11 1/2" **tk 34880-00**

Salpingograph, komplett im Holzetui  
Salpingograph, complete in wooden case  
Salpingographe, complet en etui de bois  
Salpingógrafo, completa en caja de madera  
Salpingografo, completo in astuccio di legno



4 mm Ø

x = 230 mm – 9" **tk 34872-00**

Zange zum Entfernen von Spiralen  
Forceps for I.U.D. removings

Pince pour enlever des stérilets  
Pinza para remover la espiral de esterilidad  
Pinza per distaccare la spirale

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Absaugküretten  
*Aspirating curettes*

Curettes d'aspiration  
*Curetas de aspiración*  
Curette per aspirazione



|         |                            |
|---------|----------------------------|
|         | 200 mm – 11½" <sup>a</sup> |
| 3 mm Ø  | <b>tk 35530-03</b>         |
| 4 mm Ø  | <b>tk 35530-04</b>         |
| 6 mm Ø  | <b>tk 35530-06</b>         |
| 8 mm Ø  | <b>tk 35530-08</b>         |
| 10 mm Ø | <b>tk 35530-10</b>         |
| 12 mm Ø | <b>tk 35530-12</b>         |
| 14 mm Ø | <b>tk 35530-14</b>         |

Probe- und Absaugkürette  
*Specimen and aspirating curette*  
Curette à biopsies et d'aspiration  
*Cureta de biopsia y de aspiración*  
Curette per biopsia ed per aspirazione



|                           |                    |
|---------------------------|--------------------|
| 120 mm – 4¾" <sup>a</sup> | <b>tk 35530-50</b> |
|---------------------------|--------------------|

Kupplungsstück  
*Coupling piece*  
Pièce d'accouplement  
*Pieza de conexión*  
Pezzo di connessione

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Biopsieküretten  
Biopsy curettes


Curettes pour biopsie  
Curetas para biopsía  
Curette per biopsia



| RANDALL         |                    |
|-----------------|--------------------|
|                 | 4 mm Ø             |
| 240 mm - 9 1/2" | <b>tk 34800-24</b> |



| KEVORKIAN    |                                                        |                                                       |
|--------------|--------------------------------------------------------|-------------------------------------------------------|
| 300 mm - 12" | <b>tk 34810-30</b>                                     | <b>tk 34813-30</b>                                    |
|              | fest, ohne Korb<br><i>rigid, without basket</i>        | fest, mit Korb<br><i>rigid, with basket</i>           |
|              | rigide, sans panier<br><i>rigida, sin cesta</i>        | rigide, avec panier<br><i>rigida, con cesta</i>       |
|              | rigida, senza cestino                                  | rigida, con cestino                                   |
|              | <b>tk 34811-30</b>                                     | <b>tk 34812-30</b>                                    |
|              | biegsam, ohne Korb<br><i>malleable, without basket</i> | biegsam, mit Korb<br><i>malleable, with basket</i>    |
|              | malléable, sans panier<br><i>malleable, sin cesta</i>  | malléable, avec panier<br><i>malleable, con cesta</i> |
|              | flessibile, senza cestino                              | flessibile, con cestino                               |



| NOVAK       |                    |
|-------------|--------------------|
|             | 4 mm Ø             |
| 230 mm - 9" | <b>tk 34805-23</b> |

| MINI-NOVAK |                    |
|------------|--------------------|
|            | 235 mm - 9 1/4"    |
| 1 mm Ø     | <b>tk 34807-23</b> |
| 2 mm Ø     | <b>tk 34808-23</b> |
| 3 mm Ø     | <b>tk 34809-23</b> |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Probeexzisionszangen  
Cervical biopsy and specimen forceps

Pincès emporte-pièces à biopsie  
Pinzas saca-bocados para biopsia  
Pinze taglienti per biopsia



UNI-TISCHLER

3 x 7 mm

tk 34492-01

gerade, Aufsatz allein  
straight, Tip only

droit, mors seul  
recta, boca sola  
retta, punta sola



tk 34490-00

Uni-Griff allein  
Handle only

Poignee seule  
Mango solo  
Manico solo



UNI-TISCHLER

3 x 7 mm

tk 34492-02

aufwärts, Aufsatz allein  
up, Tip only

haut, mors seul  
arriba, boca sola  
in alto, punta sola



UNI-TISCHLER

3 x 7 mm

tk 34492-03

abwärts, Aufsatz allein  
down, Tip only

bas, mors seul  
abajo, boca sola  
in basso, punta sola



UNI-TOWNSEND

2.3 x 4 mm

tk 34492-04

gerade, Aufsatz allein  
straight, Tip only

droit, mors seul  
recta, boca sola  
retta, punta sola



UNI-TOWNSEND

2.3 x 4 mm

tk 34492-05

aufwärts, Aufsatz allein  
up, Tip only

haut, mors seul  
arriba, boca sola  
in alto, punta sola



UNI-TOWNSEND

2.3 x 4 mm

tk 34492-06

abwärts, Aufsatz allein  
down, Tip only

bas, mors seul  
abajo, boca sola  
in basso, punta sola



UNI-KEVORKIAN

3 x 9 mm

tk 34492-07

Aufsatz allein  
Tip only

Mors seul  
Boca sola  
Punta sola



tk 34491-00

Uni-Griff allein mit Ring  
Handle only

Poignee seule  
Mango solo  
Manico solo



UNI-OVAL

3 x 5 mm

tk 34492-08

Aufsatz allein  
Tip only

Mors seul  
Boca sola  
Punta sola



UNI-BURKE

3 x 5 mm

tk 34492-09

Aufsatz allein  
Tip only

Mors seul  
Boca sola  
Punta sola

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Probeexzisionszangen  
*Cervical biopsy and specimen forceps*

Pinces emporte-pièces à biopsie  
*Pinzas saca-bocados para biopsia*  
Pinze taglianti per biopsia



MINI-TISCHLER  
USA MODELL

|                     |                |
|---------------------|----------------|
| x = 200 mm – 8"     | tk 34492-20 US |
| x = 250 mm – 9 7/8" | tk 34492-25 US |

MINI-  
TISCHLER

|             |
|-------------|
| tk 34492-20 |
| tk 34492-25 |

Instrument wie abgebildet, jedoch ohne Feststellschraube.  
*Instrument as illustrated, but without locking screw.*


Instrument comme illustré, mais sans vis de fixation.  
*Instrumento como ilustrado, pero sin tornillo de fijación.*  
Pinza come illustrata però senza vite di fissaggio.

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Probeexzisionszangen  
Cervical biopsy and specimen forceps

Pinces emporte-pièces à biopsie  
Pinzas saca-bocados para biopsia  
Pinze taglianti per biopsia



TISCHLER-MORGAN  
USA MODELL

|                     |                |
|---------------------|----------------|
| x = 200 mm - 8"     | tk 34491-20 US |
| x = 250 mm - 9 3/4" | tk 34491-25 US |

TISCHLER-MORGAN

|             |
|-------------|
| tk 34491-20 |
| tk 34491-25 |

Instrument wie abgebildet, jedoch ohne Feststellschraube.

Instrument as illustrated, but without locking screw.

Instrument comme illustré, mais sans vis de fixation.

Instrumento como ilustrado, pero sin tornillo de fijación.

Pinza come illustrata però senza vite di fissaggio.



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Probeexzisionszangen  
Cervical biopsy and specimen forceps

Pinces emporte-pièces à biopsie  
Pinzas saca-bocados para biopsia  
Pinze taglienti per biopsia



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Probeexzisionszangen  
Cervical biopsy and specimen forceps


Pinces emporte-pièces à biopsie  
Pinzas saca-bocados para biopsia  
Pinze taglianti per biopsia



VAN DOREN  
240 mm - 9 1/2"  
tk 34470-25



SCHUMACHER  1:1  
Fig. 1  
240 mm - 9 1/2" tk 34473-01



SCHUMACHER  1:1  1:1  
Fig. 1 Fig. 2  
x = 240 mm - 9 1/2" tk 34475-01 tk 34475-02

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Probeexzisionszangen  
Cervical biopsy and specimen forceps

Pinces emporte-pièces à biopsie  
Pinzas saca-bocados para biopsia  
Pinze taglianti per biopsia



ALEXANDER


|                 |              |
|-----------------|--------------|
|                 | 280 mm - 11" |
| x = 200 mm - 8" | tk 34450-28  |



WITTNER

|                 |             |             |
|-----------------|-------------|-------------|
| x = 230 mm - 9" | tk 34460-23 | tk 34461-23 |
|-----------------|-------------|-------------|

|          |         |
|----------|---------|
| gerade   | gebogen |
| straight | curved  |
| droit    | courbe  |
| recta    | curva   |
| rette    | curve   |



SCHUMACHER

|                 |             |             |
|-----------------|-------------|-------------|
| x = 230 mm - 9" | tk 34470-00 | tk 34470-01 |
|-----------------|-------------|-------------|

|          |         |
|----------|---------|
| komplett | Aufsatz |
| complete | Tip     |
| complète | Embout  |
| completa | Tija    |
| completa | Punta   |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Probeexzisionszangen  
*Cervical biopsy and specimen forceps*

Pinces emporte-pièces à biopsie  
*Pinzas saca-bocados para biopsia*  
Pinze taglienti per biopsia



GELLHORN

230 mm - 9"

tk 34402-23



THOMS-GAYLOR

220 mm - 8 3/4"

tk 34420-22

240 mm - 9 1/2"

tk 34420-24

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Probeexzisionszangen  
Cervical biopsy and specimen forceps

Pinces emporte-pièces à biopsie  
Pinzas saca-bocados para biopsia  
Pinze taglienti per biopsia



SCHUBERT  
210 mm - 8½" tk 34400-21  
280 mm - 11" tk 34400-29



FAURE  
210 mm - 8½" tk 34410-24




SCHUBERT  
260 mm - 10¼" tk 34440-26

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Fistelhaken  
Fistula hooks

Crochets pour fistules  
Ganchos para fístulas  
Uncini per fistole






# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Myom-Instrumente  
*Myoma instruments*

Instruments pour myomes  
*Instrumentos para fibroma*  
Strumenti per miomi



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterusfazzangen  
Uterine seizing forceps

Pinces hystérolabes  
Histerólabos  
Pinze isterolabi



**DARTIGUES**

270 mm – 10<sup>3</sup>/<sub>4</sub>" **tk 34394-27**



**SOMER**

230 mm – 9" **tk 34393-23**



**COLLIN**

250 mm – 9<sup>3</sup>/<sub>4</sub>" **tk 34395-25**

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterus-Polypenzangen  
*Uterine polypus forceps*

Pinces à polypes utérins  
*Pinzas para pólipos uterinos*  
Pinze per polipi uterini



**KELLY**

320 mm - 12½" **tk 34380-32**



**NOTO**

270 mm - 10¾" **tk 34375-27**



**DOYEN**

200 mm - 8" **tk 34370-20**

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Hakenzangen  
*Tenaculum forceps*

Pinces pour abaissement du col de l'utérus  
*Pinzas para coger el cuello del útero*  
Pinze per abbassare il collo dell'utero



PRATT

260 mm - 10 1/4" tk 34364-26



CZERNY


205 mm - 8" tk 34366-21

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Hakenzangen  
*Tenaculum forceps*

Pinces pour l'abaissement du col de l'utérus  
*Pinzas para coger el cuello del útero*  
Pinze per abbassare il collo dell'utero



MUSEUX

|       | 240 mm - 9 1/2" | 240 mm - 9 1/2" |
|-------|-----------------|-----------------|
| 8 mm  | tk 34352-24     | tk 34357-24     |
| 6 mm  | tk 34353-24     | tk 34358-24     |
| 10 mm | tk 34351-24     | tk 34356-24     |



JACOBS

|                 |             |             |
|-----------------|-------------|-------------|
| 210 mm - 8 3/4" | tk 34360-21 | tk 34361-21 |
|-----------------|-------------|-------------|



TEALE

230 mm - 9"  
tk 34328-23

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Hakenzangen  
*Tenaculum forceps*

Pinces pour l'abaissement du col de l'utérus  
*Pinzas para coger el cuello del útero*  
Pinze per abbassare il collo dell'utero



**DUPLAY**

280 mm - 11 "

**tk 34320-28**



**IOWA**

260 mm - 10 1/4 "

**tk 34329-26**



**BEACHAM**

260 mm - 10 1/4 "

**tk 34322-26**

Amniotom

*Amniotome*

Amniotomo

*Amniotomo*

Amniotomo



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Hakenzangen  
*Tenaculum forceps*

Pinces pour l'abaissement du col de l'utérus  
*Pinzas para coger el cuello del útero*  
Pinze per abbassare il collo dell'utero



POZZI

255 mm - 10"

tk 34311-26



SCHROEDER

250 mm - 9 3/4"

tk 34310-25



SCHROEDER

240 mm - 9 1/2"



1:2



1:2

tk 34314-24

tk 34315-24

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Hakenzangen  
*Tenaculum forceps*

Pinces pour l'abaissement du col de l'utérus  
*Pinzas para coger el cuello del útero*  
Pinze per abbassare il collo dell'utero



STAUDE-  
MOORE

170mm - 6 <sup>3</sup>/<sub>4</sub>"

tk 34302-17



BARRETT

180mm - 7"

tk 34305-18



SKENE

240mm - 9 <sup>1</sup>/<sub>2</sub>"

tk 34301-24



BRAUN

250mm - 9 <sup>3</sup>/<sub>4</sub>"

tk 34303-25

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterus-Dilatatoren und -Depressoren  
*Uterine dilators and depressors*

Dilatateurs du col de l'utérus et dépresseurs utérins  
*Dilatadores uterinos y depresores uterinos*  
Dilatatori uterini e depressori uterini



## HANK

|            |                 |             |
|------------|-----------------|-------------|
| Fig. 19/10 | 4.5 / 5.0 mm Ø  | tk 34280-10 |
| Fig. 11/12 | 5.5 / 6.0 mm Ø  | tk 34280-12 |
| Fig. 13/14 | 6.5 / 7.0 mm Ø  | tk 34280-14 |
| Fig. 15/16 | 7.5 / 8.0 mm Ø  | tk 34280-16 |
| Fig. 17/18 | 8.5 / 9.0 mm Ø  | tk 34280-18 |
| Fig. 19/20 | 9.5 / 10.0 mm Ø | tk 34280-20 |

## HANK

|                            |   |              |
|----------------------------|---|--------------|
| 4.5 / 5 mm – 9.5 / 10 mm Ø | 6 | tk 34281-80  |
| Stück per Satz             |   | Satz, allein |
| Pieces per set             |   | Set, only    |
| Pièces per jeu             |   | Jeu, seul    |
| Piezas para juego          |   | Juego, solo  |
| Pezzi per serie            |   | Serie, solo  |



1:3

## HANK

|                            |   |                            |
|----------------------------|---|----------------------------|
| 4.5 / 5 mm – 9.5 / 10 mm Ø | 6 | tk 34281-82                |
| Stück per Satz             |   | Satz, in Stoffrolltasche   |
| Pieces per set             |   | Set, in canvas roll        |
| Pièces per jeu             |   | Jeu, en poche de textile   |
| Piezas para juego          |   | Juego, en bolsa de tela    |
| Pezzi per serie            |   | Serie, in astuccio di tela |



## BRAUN

280 mm – 11 "

tk 34295-26

1:3

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterus-Dilatatoren  
Uterine dilators

Dilatateurs du col de l'utérus  
Dilatadores uterinos  
Dilatatori uterini



## HEGAR

|              |             |
|--------------|-------------|
| 1 / 2 mm Ø   | tk 34262-02 |
| 3 / 4 mm Ø   | tk 34262-04 |
| 5 / 6 mm Ø   | tk 34262-06 |
| 7 / 8 mm Ø   | tk 34262-08 |
| 9 / 10 mm Ø  | tk 34262-10 |
| 11 / 12 mm Ø | tk 34262-12 |
| 13 / 14 mm Ø | tk 34262-14 |
| 15 / 16 mm Ø | tk 34262-16 |
| 17 / 18 mm Ø | tk 34262-18 |
| 19 / 20 mm Ø | tk 34262-20 |
| 21 / 22 mm Ø | tk 34262-22 |
| 23 / 24 mm Ø | tk 34262-24 |
| 25 / 26 mm Ø | tk 34262-26 |

## HEGAR

|                         |    |             |
|-------------------------|----|-------------|
| 3 / 4 mm - 17 / 18 mm Ø | 14 | tk 34265-05 |
| 1 / 2 mm - 25 / 26 mm Ø | 26 | tk 34265-06 |

Stück per Satz      Satz, allein  
Pieces per set      Set, only  
Pièces per jeu      Jeu, seul  
Piezas para juego      Juego, solo  
Pezzi per serie      Serie, solo



## HEGAR

|                         |    |             |             |
|-------------------------|----|-------------|-------------|
| 3 / 4 mm - 17 / 18 mm Ø | 8  | tk 34265-00 | tk 34265-08 |
| 1 / 2 mm - 25 / 26 mm Ø | 13 |             | tk 34265-13 |

Stück per Satz      Satz in Stoffrolltasche  
Pieces per set      Set in canvas roll  
Pièces per jeu      Jeu en poche de textile  
Piezas para juego      Juego en bolsa de tela  
Pezzi per serie      Serie in astuccio di tela


Satz, in Metalletui  
Set, in metal case  
Jeu, en boîte de métal  
Juego, en caja de metal  
Serie, in astuccio metallico

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterus-Dilatatoren  
*Uterine dilators*

Dilatateurs du col de l'utérus  
*Dilatadores uterinos*  
Dilatatori uterini



## HEGAR

|                |    |             |
|----------------|----|-------------|
| 4 mm - 17 mm Ø | 14 | tk 34261-14 |
| 1 mm - 26 mm Ø | 26 | tk 34261-26 |

|                   |                               |
|-------------------|-------------------------------|
| Stück per Satz    | Satz, in Metalletui           |
| Pieces per set    | Set, in metal case            |
| Pièces per jeu    | Jeu, en boîte de métal        |
| Piezas para juego | Juego, en caja de metal       |
| Pezzi per serie   | Serie, con astuccio metallico |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterus-Dilatatoren  
*Uterine dilators*

Dilateurs du col de l'utérus  
*Dilatadores uterinos*  
Dilatatori uterini



## HEGAR

|          |                    |
|----------|--------------------|
| 1.0 mm Ø | <b>tk 34250-10</b> |
| 1.5 mm Ø | <b>tk 34250-15</b> |
| 2.0 mm Ø | <b>tk 34250-20</b> |
| 2.5 mm Ø | <b>tk 34250-25</b> |
| 3.0 mm Ø | <b>tk 34250-30</b> |
| 3.5 mm Ø | <b>tk 34250-35</b> |
| 4.0 mm Ø | <b>tk 34250-40</b> |
| 4.5 mm Ø | <b>tk 34250-45</b> |
| 5.0 mm Ø | <b>tk 34250-50</b> |
| 5.5 mm Ø | <b>tk 34250-55</b> |
| 6.0 mm Ø | <b>tk 34250-60</b> |
| 6.5 mm Ø | <b>tk 34250-65</b> |
| 7.0 mm Ø | <b>tk 34250-70</b> |
| 7.5 mm Ø | <b>tk 34250-75</b> |
| 8.0 mm Ø | <b>tk 34250-80</b> |
| 8.5 mm Ø | <b>tk 34250-85</b> |
| 9.0 mm Ø | <b>tk 34250-90</b> |
| 9.5 mm Ø | <b>tk 34250-95</b> |

|           |                    |
|-----------|--------------------|
| 10.0 mm Ø | <b>tk 34251-00</b> |
| 10.5 mm Ø | <b>tk 34251-05</b> |
| 11.0 mm Ø | <b>tk 34251-10</b> |
| 11.5 mm Ø | <b>tk 34251-15</b> |
| 12.0 mm Ø | <b>tk 34251-20</b> |
| 12.5 mm Ø | <b>tk 34251-25</b> |
| 13.0 mm Ø | <b>tk 34251-30</b> |
| 13.5 mm Ø | <b>tk 34251-35</b> |
| 14.0 mm Ø | <b>tk 34251-40</b> |
| 14.5 mm Ø | <b>tk 34251-45</b> |
| 15.0 mm Ø | <b>tk 34251-50</b> |
| 15.5 mm Ø | <b>tk 34251-55</b> |
| 16.0 mm Ø | <b>tk 34251-60</b> |
| 16.5 mm Ø | <b>tk 34251-65</b> |
| 17.0 mm Ø | <b>tk 34251-70</b> |
| 17.5 mm Ø | <b>tk 34251-75</b> |
| 18.0 mm Ø | <b>tk 34251-80</b> |
| 18.5 mm Ø | <b>tk 34251-85</b> |
| 19.0 mm Ø | <b>tk 34251-90</b> |
| 19.5 mm Ø | <b>tk 34251-95</b> |

|           |                    |
|-----------|--------------------|
| 20.0 mm Ø | <b>tk 34252-00</b> |
| 21.0 mm Ø | <b>tk 34252-10</b> |
| 22.0 mm Ø | <b>tk 34252-20</b> |
| 23.0 mm Ø | <b>tk 34252-30</b> |
| 24.0 mm Ø | <b>tk 34252-40</b> |
| 25.0 mm Ø | <b>tk 34252-50</b> |
| 26.0 mm Ø | <b>tk 34252-60</b> |
| 27.0 mm Ø | <b>tk 34252-70</b> |
| 28.0 mm Ø | <b>tk 34252-80</b> |
| 29.0 mm Ø | <b>tk 34252-90</b> |
| 30.0 mm Ø | <b>tk 34253-00</b> |

## HEGAR

|                |    |                    |
|----------------|----|--------------------|
| 4 mm - 17 mm Ø | 14 | <b>tk 34260-14</b> |
| 1 mm - 26 mm Ø | 26 | <b>tk 34260-26</b> |

|                   |                                 |
|-------------------|---------------------------------|
| Stück per Satz    | Satz, ohne Metalletui           |
| Pieces per set    | Set, without metal case         |
| Pièces per jeu    | Jeu, sans boîte en métal        |
| Piezas para juego | Juego, sin caja de metal        |
| Pezzi per serie   | Serie, senza astuccio metallico |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Uterus-Sonden  
Uterine sounds

Hystéromètres du col de l'utérus  
Histerómetros uterinos  
Isterometri uterini



|                                                                                                                           |                                                                                    |                                                                                      |
|---------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
|                                         |  |  |
| VALLEIX 1:2                                                                                                               | MARTIN                                                                             | SIMS                                                                                 |
| 250 mm – 9¾" tk 34197-25                                                                                                  | 320 mm – 12½" tk 34198-32 tk 34199-32                                              | 320 mm – 12½" tk 34200-32 tk 34201-32                                                |
| graduirt, mit Schieber<br>graduated, with slider<br>gradué, avec curseur<br>graduado, con cursor<br>graduato, con cursore | fest<br>rigid<br>rigide<br>rigido<br>rigido                                        | biegsam<br>malleable<br>malléable<br>maleable<br>malleabile                          |
|                                                                                                                           |                                                                                    | fest<br>rigid<br>rigide<br>rigido<br>rigido                                          |
|                                                                                                                           |                                                                                    | biegsam<br>malleable<br>malléable<br>maleable<br>malleabile                          |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Hysterektomieklemmen  
Hysterectomy forceps

Pinces à hystérectomie  
Pinzas-clamps para histerectomia  
Pinze per isterectomia








# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Parametriumklemmen  
*Parametrium forceps*

Pinces parametrium  
*Pinzas parametrium*  
Pinze parametrium



|                  |             |             |             |             |             |
|------------------|-------------|-------------|-------------|-------------|-------------|
| 210 mm - 8 1/4"  | tk 34155-21 | tk 34156-21 | tk 34157-21 | tk 34158-21 | tk 34159-21 |
| 240 mm - 9 1/2"  | tk 34155-24 | tk 34156-24 | tk 34157-24 | tk 34158-24 | tk 34159-24 |
| 300 mm - 12"     | tk 34155-30 | tk 34156-30 | tk 34157-30 | tk 34158-30 | tk 34159-30 |
| 350 mm - 13 3/4" | tk 34155-35 | tk 34156-35 | tk 34157-35 | tk 34158-35 | tk 34159-35 |

S-förmig gebogen  
*S-curved*

coube en S  
*curva en S*  
curva a forma di S

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Endospekula  
Endospecula

Endospéculums  
Endospéculos  
Endospecoli



KOGAN

240 mm - 9½"

tk 34153-24

mit Skalensperre und Stellschraube  
with scale-ratchet and fixation screw

avec crémaillère à échelle et vis à fixation  
con cremallera graduada y tornillo fijación  
con cremagliera graduata e vite di fissaggio



KOGAN

240 mm - 9½"

tk 34154-24



KOGAN-MODIF.

280 mm - 11"

x = 4 mm y = 8 mm tk 34150-01

x = 2 mm y = 4 mm tk 34150-02

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scheidenspekula  
Vaginal specula

Spéculums vaginaux  
Espéculos vaginales  
Specoli vaginali



|               |               |               |                       |                       |
|---------------|---------------|---------------|-----------------------|-----------------------|
| Fig. 1        | Fig. 2        | Fig. 3        | Fig. 4                | Fig. 5                |
| 80 mm x 30 mm | 85 mm x 35 mm | 90 mm x 40 mm | 55 mm x 40 mm / 60 mm | 70 mm x 45 mm / 70 mm |
| tk 34090-01   | tk 34090-02   | tk 34090-03   | tk 34090-04           | tk 34090-05           |

SCHERBACK

tk 34090-00

Satz mit 400 g Gewicht (tk 34090-10)  
Set with weight 400 grs. (tk 34090-10)

Jeu avec poid 400 gr. (tk 34090-10)  
Juego con peso 400 gr. (tk 34090-10)  
Serie con peso di 400 gr. (tk 34090-10)



|       |             |
|-------|-------------|
| 400 g | tk 34090-10 |
| 800 g | tk 34090-12 |



tk 34090-14



DOYEN

|                |                 |                 |
|----------------|-----------------|-----------------|
|                | 240 mm - 9 1/2" | 240 mm - 9 1/2" |
| 55 mm x 35 mm  | tk 34100-01     | tk 34110-01     |
| 75 mm x 35 mm  | tk 34100-02     | tk 34110-02     |
| 115 mm x 35 mm | tk 34100-03     | tk 34110-03     |
| 55 mm x 45 mm  | tk 34102-01     | tk 34112-01     |
| 75 mm x 45 mm  | tk 34102-02     | tk 34112-02     |
| 115 mm x 45 mm | tk 34102-03     | tk 34112-03     |
| 160 mm x 45 mm | tk 34102-04     |                 |
| 55 mm x 60 mm  | tk 34104-01     | tk 34114-01     |
| 75 mm x 60 mm  | tk 34104-02     | tk 34114-02     |
| 115 mm x 60 mm | tk 34104-03     | tk 34114-03     |
| 160 mm x 60 mm | tk 34104-04     |                 |

leighthohl  
slightly concave

légèrement concave  
ligeramente concava  
legermente concavo

halbhohl  
medium concave

semi concave  
semi concava  
semi concavo

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scheidenspekula  
Vaginal specula

Spéculums vaginaux  
Espéculos vaginales  
Specoli vaginali



LANDAU

|                |                    |
|----------------|--------------------|
|                | 290 mm – 11½"      |
| 100 mm x 28 mm | <b>tk 34128-01</b> |
| 110 mm x 38 mm | <b>tk 34128-02</b> |

BREISKY

|                |                    |
|----------------|--------------------|
|                | 300 mm – 12"       |
| 130 mm x 20 mm | <b>tk 34130-20</b> |
| 130 mm x 25 mm | <b>tk 34130-25</b> |
| 130 mm x 30 mm | <b>tk 34130-30</b> |
| 130 mm x 35 mm | <b>tk 34130-35</b> |
| 130 mm x 40 mm | <b>tk 34130-40</b> |

|                |                    |
|----------------|--------------------|
|                | 300 mm – 12"       |
| 100 mm x 20 mm | <b>tk 34142-20</b> |
| 100 mm x 25 mm | <b>tk 34142-25</b> |
| 100 mm x 30 mm | <b>tk 34142-30</b> |
| 100 mm x 35 mm | <b>tk 34142-35</b> |
| 100 mm x 40 mm | <b>tk 34142-40</b> |

|                |                    |
|----------------|--------------------|
|                | 350 mm – 13¾"      |
| 160 mm x 40 mm | <b>tk 34140-40</b> |
| 180 mm x 40 mm | <b>tk 34140-41</b> |

Wiener Modell  
Vienna pattern

Modèle Vienne  
Modelo Vienna  
Modello Vienna



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scheidenspekula  
Vaginal specula

Spéculums vaginaux  
Espéculos vaginales  
Specoli vaginali



KALLMORGEN

|               |             |
|---------------|-------------|
| 70 mm x 40 mm | tk 34072-01 |
| 90 mm x 40 mm | tk 34072-02 |

Spekulum  
Speculum

Spéculum  
Espéculo  
Specolo



|               |             |
|---------------|-------------|
| 70 mm x 33 mm | tk 34073-01 |
| 90 mm x 33 mm | tk 34073-02 |

Halter  
Retractor

Valve  
Valva  
Valve



|             |
|-------------|
| tk 34074-01 |
| tk 34074-02 |

Satz  
Set

Jeu  
Juego  
Serie

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scheidenspekula  
Vaginal specula

Spéculums vaginaux  
Espéculos vaginales  
Specoli vaginali



JACKSON

|                |             |
|----------------|-------------|
| 75 mm x 35 mm  | tk 34070-10 |
| 75 mm x 40 mm  | tk 34070-15 |
| 90 mm x 35 mm  | tk 34070-20 |
| 90 mm x 40 mm  | tk 34070-25 |
| 100 mm x 35 mm | tk 34070-30 |
| 100 mm x 40 mm | tk 34070-35 |



AUVARD

|                 |
|-----------------|
| 220 mm - 8 3/4" |
| 80 mm x 38 mm   |
| tk 34060-00     |



SIMS

|                       |             |
|-----------------------|-------------|
| 60 mm x 25 mm / 30 mm | tk 34075-01 |
| 70 mm x 30 mm / 35 mm | tk 34075-02 |
| 80 mm x 35 mm / 40 mm | tk 34075-03 |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



KRISTELLER

|                |                    |
|----------------|--------------------|
| 105mm x 18 mm  | <b>tk 34055-01</b> |
| 105 mm x 21 mm | <b>tk 34055-02</b> |

Spekulum konisch, für Jungfrauen  
*Speculum tapered, for virgins*  
Spéculum conique, pour vierges  
*Espéculo conico, para virgenes*  
Specolo conico, per virgini



|                |                    |
|----------------|--------------------|
| 110 mm x 14 mm | <b>tk 34056-01</b> |
| 115 mm x 17 mm | <b>tk 34056-02</b> |

Halter  
*Retractor*  
Valve  
*Valva*  
Valve



|                    |
|--------------------|
| <b>tk 34057-01</b> |
| <b>tk 34057-02</b> |

Satz  
*Set*  
Jeu  
*Juego*  
Serie

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scheidenspekula  
*Vaginal specula*

Spéculums vaginaux  
*Espéculos vaginales*  
Specoli vaginali



## KRIESTELLER

|                |             |
|----------------|-------------|
| 70 mm x 27 mm  | tk 34050-01 |
| 80 mm x 30 mm  | tk 34050-02 |
| 90 mm x 36 mm  | tk 34050-03 |
| 110 mm x 27 mm | tk 34050-04 |
| 110 mm x 30 mm | tk 34050-05 |
| 110 mm x 36 mm | tk 34050-06 |

Spekulum  
*Speculum*

Spéculum  
*Espéculo*  
Specolo

|                |             |
|----------------|-------------|
| 75 mm x 23 mm  | tk 34051-01 |
| 85 mm x 26 mm  | tk 34051-02 |
| 100 mm x 32 mm | tk 34051-03 |
| 115 mm x 23 mm | tk 34051-04 |
| 115 mm x 26 mm | tk 34051-05 |
| 115 mm x 32 mm | tk 34051-06 |

Halter  
*Retractor*

Valve  
*Valva*  
Valve

|             |
|-------------|
| tk 34052-01 |
| tk 34052-02 |
| tk 34052-03 |
| tk 34052-04 |
| tk 34052-05 |
| tk 34052-06 |

Satz  
*Set*


Jeu  
*Juego*  
Serie

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scheidenspekula  
Vaginal specula

Spéculums vaginaux  
Espéculos vaginales  
Specoli vaginali



**COLLIN**

|                |                    |
|----------------|--------------------|
| 85 mm x 30 mm  | <b>tk 34018-01</b> |
| 100 mm x 35 mm | <b>tk 34018-02</b> |
| 110 mm x 40 mm | <b>tk 34018-03</b> |



**SEMM**

|                     |                    |
|---------------------|--------------------|
| 90 mm x 14 / 18 mm  | <b>tk 34020-00</b> |
| 100 mm x 17 / 20 mm | <b>tk 34020-01</b> |
| 100 mm x 25 / 30 mm | <b>tk 34020-02</b> |



**TRELAT**

|                |                    |
|----------------|--------------------|
| 85 mm x 30 mm  | <b>tk 34030-01</b> |
| 95 mm x 35 mm  | <b>tk 34030-02</b> |
| 115 mm x 43 mm | <b>tk 34030-03</b> |



**PEDERSON**

|                |                    |
|----------------|--------------------|
| 75 mm x 13 mm  | <b>tk 34032-01</b> |
| 100 mm x 22 mm | <b>tk 34032-02</b> |
| 120 mm x 25 mm | <b>tk 34032-03</b> |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scheidenspekula  
*Vaginal specula*

Spéculums vaginaux  
*Espéculos vaginales*  
Specoli vaginali



## CUSCO

|                |             |
|----------------|-------------|
| 80 mm x 22 mm  | tk 34012-01 |
| 100 mm x 25 mm | tk 34012-02 |
| 105 mm x 27 mm | tk 34012-03 |

## Schweizer Modell *Swiss pattern*

Modèle suisse  
*Modelo suizo*  
Modello svizzero



## GRAVE

|                |             |
|----------------|-------------|
| 75 mm x 20 mm  | tk 34015-01 |
| 95 mm x 35 mm  | tk 34015-02 |
| 115 mm x 35 mm | tk 34015-03 |



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Scheidenspekula  
Vaginal specula

Spéculums vaginaux  
Espéculos vaginales  
Specoli vaginali



Cusco

|                |             |
|----------------|-------------|
| 75 mm x 32 mm  | tk 34001-01 |
| 85 mm x 35 mm  | tk 34001-02 |
| 100 mm x 35 mm | tk 34001-03 |

|               |             |
|---------------|-------------|
| 75 mm x 17 mm | tk 34005-04 |
|---------------|-------------|

für Jungfrauen  
for virgins

pour vierges  
para virgenes  
per vergine



Cusco

|                |             |
|----------------|-------------|
| 75 mm x 32 mm  | tk 34007-01 |
| 85 mm x 35 mm  | tk 34007-02 |
| 100 mm x 37 mm | tk 34007-03 |

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

Beckenzirkel  
Pelvimeters

Pelvimètres  
Pelvimetros  
Pelvimetri



MARTIN

340 mm - 13 1/2" tk35800-34

cm / Zoll-Scala  
cm / inch-scale

échelle en cm / pouces  
escala cm / pulgadas  
scala cm / pollici



COLLIN

350 mm - 13 3/4" tk35805-35

cm-Scala  
cm-scale

échelle en cm  
escala en cm  
scala in cm

# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

---

PRIME SURGICAL INTERNATIONAL

## *electrosurgery*

### **Sklar Blue™ Electrosurgical Instruments**

Sklar Blue™ Electrosurgery Instruments are designed to provide an effective, non-conductive, patient-safe surface. Sklar's line of electrosurgical instruments features our exclusive Sklar Blue™ insulated coating that prevents transmission of RF current through the instrument. Manufactured of the highest quality stainless steel, these instruments are fully autoclavable. All Sklar Blue™ electrosurgical coatings are guaranteed for one year against manufacturing defects.

Speculums .....

Retractors .....

Endocervical Speculums .....

Forceps .....

Hooks .....

### **Electrosurgical Accessories**

Sklar's long held tradition of providing quality surgical instruments is carried through with the establishment of the Sklar line of electrosurgical pencils, grounding pads, active electrodes and accessories. Offering unsurpassed performance and design features, Sklar Electrosurgical products delivery quality you can feel, see and trust for any challenging surgical procedure.

Electrosurgical Pencils .....

Electrodes .....

Accessories .....

Bipolar Forceps .....

## ELECTROSURGERY

OB/GYN

### Bipolar Forceps

#### JANSEN BAYONET FORCEP

| PRODUCT NUMBER | FEATURES |
|----------------|----------|
| 91-6106        | 1 MM, 6" |



#### JEWELER #5 FORCEP

| PRODUCT NUMBER | FEATURES  |
|----------------|-----------|
| 91-6107        | .5 MM 4½" |



#### MALIS FORCEP

| PRODUCT NUMBER | FEATURES |
|----------------|----------|
| 91-6108        | 1 MM, 7" |



#### MCPHERSON FORCEP

| PRODUCT NUMBER | FEATURES          |
|----------------|-------------------|
| 91-6109        | .5 MM, 3½" Curved |



#### MCPHERSON FORCEP

| PRODUCT NUMBER | FEATURES            |
|----------------|---------------------|
| 91-6110        | .5 MM, 3½" Straight |



#### SEMKIN FORCEP

| PRODUCT NUMBER | FEATURES   |
|----------------|------------|
| 91-6111        | .5 MM, 5½" |



#### TENZEL FORCEP

| PRODUCT NUMBER | FEATURES  |
|----------------|-----------|
| 91-6112        | .5 MM, 4" |



www.primeinstrumentintl.com

## OB/GYN

## ELECTROSURGERY

### Bipolar Forceps

#### ADSON FORCECEP

| PRODUCT NUMBER | FEATURES  |
|----------------|-----------|
| 91-6100        | 1 MM, 4¾" |



#### COAPTATION FORCECEP

| PRODUCT NUMBER | FEATURES  |
|----------------|-----------|
| 91-6101        | .5 MM, 4" |



#### CUSHING FORCECEP

| PRODUCT NUMBER | FEATURES |
|----------------|----------|
| 91-6102        | 1 MM, 6" |



#### GERALD BAYONET FORCECEP

| PRODUCT NUMBER | FEATURES |
|----------------|----------|
| 91-6103        | 1 MM, 7" |



#### GERALD FORCECEP

| PRODUCT NUMBER | FEATURES |
|----------------|----------|
| 91-6104        | 1 MM, 7" |



#### HARDY BAYONET FORCECEP

| PRODUCT NUMBER | FEATURES  |
|----------------|-----------|
| 91-6105        | 1 MM, 8½" |



#### WARNING: STERILIZATION GUIDELINES

Follow manufacturer's recommendations and AAMI standards for sterilization.

STEAM: steam autoclave at 250°F / 121°C for 20 minutes.

FLASH: steam autoclave at 275°F / 134°C for 10 minutes.

GAS (ETO): follow manufacturer's guide for the unit which you are using.



## ELECTROSURGERY Accessories

OB/GYN

### CORDS AND CABLES

| PRODUCT NUMBER | FEATURES                                                                | QUANTITIES |
|----------------|-------------------------------------------------------------------------|------------|
| 96-9698        | Reusable Return Pad Cable For Use With Non-Corded Pads, 15' Non-Sterile | 10/Box     |

| PRODUCT NUMBER | FEATURES                                                    | QUANTITIES           |
|----------------|-------------------------------------------------------------|----------------------|
| 96-9700        | Single Use, Sterile, Bipolar Forceps Cord, Flying Lead, 12' | 10/Box               |
| 96-9701        | Single Use, Sterile, Bipolar Forceps Cord, Flying Lead, 12' | 10/bx, 5bx/cs, 50/cs |

| PRODUCT NUMBER | FEATURES                                                        | QUANTITIES           |
|----------------|-----------------------------------------------------------------|----------------------|
| 96-9702        | Single Use, Sterile, Bipolar Forceps Cord, Solid Connector, 12' | 10/Box               |
| 96-9703        | Single Use, Sterile, Bipolar Forceps Cord, Solid Connector, 12' | 10/bx, 5bx/cs, 50/cs |



### TIP CLEANERS-DISPOSABLE

| PRODUCT NUMBER | FEATURES    | QUANTITIES |
|----------------|-------------|------------|
| 96-9704        | Sterile     | 50/Box     |
| 96-9706        | Non-Sterile | 500/Box    |

42.5 MM X 42.5 MM  
(1.7" X 1.7")



### HOLSTERS-DISPOSABLE

| PRODUCT NUMBER | FEATURES    | QUANTITIES |
|----------------|-------------|------------|
| 96-9708        | Sterile     | 100/Box    |
| 96-9710        | Non-Sterile | 500/Box    |



## OB/GYN

## ELECTROSURGERY

### Grounding Pads | Single Use | Sterile

Sklar Premium Patient Return Pads are manufactured using an advanced design hydro-gel adhesive that, incorporated with a gel overlay system, provide a much safer return pad for both the patient and surgeon. Sklar offers Patient Return Pads in five different configurations.

#### RETURN ELECTRODES / GROUNDING PADS

| PRODUCT NUMBER | FEATURES                                 | QUANTITIES            |
|----------------|------------------------------------------|-----------------------|
| 96-9688        | Dual Element, 9' Cord Attached, Standard | 50/Box                |
| 96-9690        | Dual Element, Non-Corded, Standard       | 5/pk, 40pk/cs, 200/bx |

96-9688 and 96-9690 are dual element return pads designed to be used with electrosurgical units that are equipped with REM or Neutral Electrode safety systems that continually monitor impedance levels. The 96-9688 comes with a pre-attached 9 foot cord with standard plug. The 96-9690 requires the use of a reusable cable, 96-9698.



| PRODUCT NUMBER | FEATURES                                   | QUANTITIES            |
|----------------|--------------------------------------------|-----------------------|
| 96-9692        | Single Element, 9' Cord Attached, Standard | 50/Box                |
| 96-9694        | Single Element, Non-Corded, Standard       | 5/pk, 40pk/cs, 200/bx |

96-9692 and 96-9694 are traditional single element return pads manufactured using the same hydro-gel technology. Both are designed for use with all electrosurgical generators. The 96-9692 comes with a pre-attached 9 foot cord, with standard connector and the 96-9694 requires the use of a reusable cable, 96-9698.



| PRODUCT NUMBER | FEATURES                                           | QUANTITIES |
|----------------|----------------------------------------------------|------------|
| 96-9696        | Dual Element, 9' Cord Attached, Pediatric / Infant | 50/Box     |



## ELECTROSURGERY

OB/GYN

Electrodes | Single Use | Sterile

Sklar Premium Straight Ball Electrodes feature our Hex Locking system to assure stability during difficult procedures. They also have our Ribbed Molded insulation for easy removal and extended insulation for protection of the surgeon and patient.

### STRAIGHT BALL ELECTRODES

| PRODUCT NUMBER | FEATURES                            | QUANTITIES             |
|----------------|-------------------------------------|------------------------|
| 96-9676        | Stainless Steel 4mm Ball, 2"        | 10/Box                 |
| 96-9677        | Stainless Steel 4mm Ball, 2"        | 10/bx, 10bx/cs, 100/cs |
| 96-9678        | Stainless Steel 4mm Coated Ball, 2" | 10/Box                 |
| 96-9679        | Stainless Steel 4mm Coated Ball, 2" | 10/bx, 10bx/cs, 100/cs |
| 96-9680        | Stainless Steel 5mm Ball, 2"        | 10/Box                 |
| 96-9681        | Stainless Steel 5mm Ball, 2"        | 10/bx, 10bx/cs, 100/cs |
| 96-9682        | Stainless Steel 5mm Coated Ball, 2" | 10/Box                 |
| 96-9683        | Stainless Steel 5mm Coated Ball, 2" | 10/bx, 10bx/cs, 100/cs |



|         |                                |                        |
|---------|--------------------------------|------------------------|
| 96-9684 | Stainless Steel 5mm Ball, 5.3" | 10/Box                 |
| 96-9685 | Stainless Steel 5mm Ball, 5.3" | 10/bx, 10bx/cs, 100/cs |



|         |                                       |                        |
|---------|---------------------------------------|------------------------|
| 96-9686 | Stainless Steel 5mm Coated Ball, 5.3" | 10/Box                 |
| 96-9687 | Stainless Steel 5mm Coated Ball, 5.3" | 10/bx, 10bx/cs, 100/cs |



## OB/GYN

## ELECTROSURGERY

Electrodes | Single Use | Sterile

Sklar Premium Needle Electrodes feature a durable and ultra fine needle point for precise locations during the surgical procedure. All electrodes feature a ribbed molded insulation for easy removal and extended insulation to protect the surgeon and patient.

### STRAIGHT NEEDLE ELECTRODES

| PRODUCT NUMBER | FEATURES                                  | QUANTITIES             |
|----------------|-------------------------------------------|------------------------|
| 96-9664        | Solid Stainless Steel Needle, 2.8"        | 10/Box                 |
| 96-9665        | Solid Stainless Steel Needle, 2.8"        | 10/bx, 10bx/cs, 100/cs |
| 96-9666        | Solid Stainless Steel Coated Needle, 2.8" | 10/Box                 |
| 96-9667        | Solid Stainless Steel Coated Needle, 2.8" | 10/bx, 10bx/cs, 100/cs |



|         |                                  |                        |
|---------|----------------------------------|------------------------|
| 96-9668 | Solid Stainless Steel Needle, 6" | 10/Box                 |
| 96-9669 | Solid Stainless Steel Needle, 6" | 10/bx, 10bx/cs, 100/cs |



|         |                                         |                        |
|---------|-----------------------------------------|------------------------|
| 96-9670 | Solid Stainless Steel Coated Needle, 6" | 10/Box                 |
| 96-9671 | Solid Stainless Steel Coated Needle, 6" | 10/bx, 10bx/cs, 100/cs |



|         |                                   |                        |
|---------|-----------------------------------|------------------------|
| 96-9672 | Solid Stainless Steel Needle, 10" | 10/Box                 |
| 96-9673 | Solid Stainless Steel Needle, 10" | 10/bx, 10bx/cs, 100/cs |



|         |                                          |                        |
|---------|------------------------------------------|------------------------|
| 96-9674 | Solid Stainless Steel Coated Needle, 10" | 10/Box                 |
| 96-9675 | Solid Stainless Steel Coated Needle, 10" | 10/bx, 10bx/cs, 100/cs |





## ELECTROSURGERY

Electrodes | Single Use | Sterile

OB/GYN

### LLETZ BALL ELECTRODES

| PRODUCT NUMBER | FEATURES      | QUANTITIES |
|----------------|---------------|------------|
| 96-1040        | 3 MM / 0.3 CM | 5/Box      |

### LLETZ SQUARE ELECTRODES

| PRODUCT NUMBER | FEATURES                 | QUANTITIES |
|----------------|--------------------------|------------|
| 96-1099        | 10 x 8 MM / 1.0 x 0.8 CM | 5/Box      |

### LLETZ LOOP ELECTRODES

| PRODUCT NUMBER | FEATURES                  | QUANTITIES |
|----------------|---------------------------|------------|
| 96-1060        | 15 x 8 MM / 1.5 x 0.8 CM  | 5/Box      |
| 96-1075        | 20 x 8 MM / 2.0 x 0.8 CM  | 5/Box      |
| 96-1085        | 20 x 20 MM / 2.0 x 2.0 CM | 5/Box      |



## OB/GYN

## ELECTROSURGERY

Electrodes | Single Use | Sterile

Sklar Premium Loop Electrodes feature a strong and durable tungsten steel loop incorporated into a reinforced T-molded head for added stability during difficult and delicate surgical procedures. All electrodes are a standard five inches in length and feature an extended insulation covering to protect the surgeon and patient.

### LOOP ELECTRODES

| PRODUCT NUMBER                                                                       | FEATURES                         | QUANTITIES           |
|--------------------------------------------------------------------------------------|----------------------------------|----------------------|
| 96-9652                                                                              | Tungsten Loop, 10 MM x 5 MM, 5"  | 5/Box                |
| 96-9653                                                                              | Tungsten Loop, 10 MM x 5 MM, 5"  | 5/bx, 10bx/cs, 50/cs |
|    |                                  |                      |
| 96-9654                                                                              | Tungsten Loop, 10 MM x 10 MM, 5" | 5/Box                |
| 96-9655                                                                              | Tungsten Loop, 10 MM x 10 MM, 5" | 5/bx, 10bx/cs, 50/cs |
|    |                                  |                      |
| 96-9656                                                                              | Tungsten Loop, 15 MM x 10 MM, 5" | 5/Box                |
| 96-9657                                                                              | Tungsten Loop, 15 MM x 10 MM, 5" | 5/bx, 10bx/cs, 50/cs |
|   |                                  |                      |
| 96-9658                                                                              | Tungsten Loop, 15 MM x 15 MM, 5" | 5/Box                |
| 96-9659                                                                              | Tungsten Loop, 15 MM x 15 MM, 5" | 5/bx, 10bx/cs, 50/cs |
|  |                                  |                      |
| 96-9660                                                                              | Tungsten Loop, 20 MM x 10 MM, 5" | 5/Box                |
| 96-9661                                                                              | Tungsten Loop, 20 MM x 10 MM, 5" | 5/bx, 10bx/cs, 50/cs |
|  |                                  |                      |
| 96-9662                                                                              | Tungsten Loop, 20 MM x 15 MM, 5" | 5/Box                |
| 96-9663                                                                              | Tungsten Loop, 20 MM x 15 MM, 5" | 5/bx, 10bx/cs, 50/cs |
|  |                                  |                      |



## ELECTROSURGERY

OB/GYN

### Electrodes | Single Use | Sterile

Sklar Premium Blade Electrodes feature a ribbed molded insulator, which allows the surgeon or technician easy removal of the electrode during procedures. They are available with a hex lock feature that assures stability and prohibits rotation of the blade during difficult procedures. They are also available coated with an anti-stick material that allows easy clean up of accumulated debris.

#### BLADE ELECTRODES

| PRODUCT NUMBER | FEATURES                                      | QUANTITIES             |
|----------------|-----------------------------------------------|------------------------|
| 96-9632        | Flex Tip Blade, 2.4"                          | 10/Box                 |
| 96-9633        | Flex Tip Blade, 2.4"                          | 10/bx, 10bx/cs, 100/cs |
| 96-9634        | Flex Tip Blade, Hex Lock Hub, 2.4"            | 10/Box                 |
| 96-9635        | Flex Tip Blade, Hex Lock Hub, 2.4"            | 10/bx, 10bx/cs, 100/cs |
| 96-9636        | Solid Steel Blade, 2.75"                      | 10/Box                 |
| 96-9637        | Solid Steel Blade, 2.75"                      | 10/bx, 10bx/cs, 100/cs |
| 96-9638        | Solid Steel Coated Blade, 2.75"               | 10/Box                 |
| 96-9639        | Solid Steel Coated Blade 2.75"                | 10/bx, 10bx/cs, 100/cs |
| 96-9640        | Solid Steel Blade, Hex Lock Hub, 2.75"        | 10/Box                 |
| 96-9641        | Solid Steel Blade, Hex Lock Hub, 2.75"        | 10/bx, 10bx/cs, 100/cs |
| 96-9642        | Solid Steel Coated Blade, Hex Lock Hub, 2.75" | 10/Box                 |
| 96-9643        | Solid Steel Coated Blade, Hex Lock Hub, 2.75" | 10/bx, 10bx/cs, 100/cs |



| PRODUCT NUMBER | FEATURES              | QUANTITIES             |
|----------------|-----------------------|------------------------|
| 96-9644        | Solid Steel Blade, 6" | 10/Box                 |
| 96-9645        | Solid Steel Blade, 6" | 10/bx, 10bx/cs, 100/cs |



| PRODUCT NUMBER | FEATURES                     | QUANTITIES             |
|----------------|------------------------------|------------------------|
| 96-9646        | Solid Steel Coated Blade, 6" | 10/Box                 |
| 96-9647        | Solid Steel Coated Blade, 6" | 10/bx, 10bx/cs, 100/cs |



| PRODUCT NUMBER | FEATURES               | QUANTITIES             |
|----------------|------------------------|------------------------|
| 96-9648        | Solid Steel Blade, 10" | 10/Box                 |
| 96-9649        | Solid Steel Blade, 10" | 10/bx, 10bx/cs, 100/cs |



| PRODUCT NUMBER | FEATURES                      | QUANTITIES             |
|----------------|-------------------------------|------------------------|
| 96-9650        | Solid Steel Coated Blade, 10" | 10/Box                 |
| 96-9651        | Solid Steel Coated Blade, 10" | 10/bx, 10bx/cs, 100/cs |



## OB/GYN

## ELECTROSURGERY

### Electrosurgical Pencils | Single Use | Sterile

#### Sklar Premium

Electrosurgical Pencils are available in a push button, foot switch or rocker switch configuration. Sklar pencils are ergonomically designed to provide maximum balance and comfort for the surgeon. Each pencil features a removable 2.75" stainless steel blade electrode, with hex locking feature to provide stability during difficult surgical procedures.

Unsurpassed performance and design features, like sealed switching components, make the Sklar Premium Electrosurgical Pencil the right choice. Pencil tips are interchangeable with any of Sklar's electrodes and most other manufacturer's electrodes. Sklar Pencils come standard with a pre-attached 3 meter cord with a standard three prong plug. All Electrosurgical Pencils are provided in a single use, sterile pouch.

#### PUSH BUTTON SWITCH

| PRODUCT NUMBER | FEATURES                     | ACCESSORIES          | QUANTITIES            |
|----------------|------------------------------|----------------------|-----------------------|
| 96-9601        | 2.75" Blade Electrode        | –                    | 50/Box                |
| 96-9602        | 2.75" Blade Electrode        | –                    | 50/bx, 2bx/cs, 100/cs |
| 96-9603        | 2.75" Blade Electrode        | Holster              | 50/Box                |
| 96-9604        | 2.75" Blade Electrode        | Holster              | 50/bx, 2bx/cs, 100/cs |
| 96-9605        | 2.75" Blade Electrode        | Holster, Tip Cleaner | 50/Box                |
| 96-9606        | 2.75" Blade Electrode        | Holster, Tip Cleaner | 50/bx, 2bx/cs, 100/cs |
| 96-9607        | 2.75" Coated Blade Electrode | –                    | 50/Box                |
| 96-9608        | 2.75" Coated Blade Electrode | –                    | 50/bx, 2bx/cs, 100/cs |
| 96-9609        | 2.75" Coated Blade Electrode | Holster              | 50/Box                |
| 96-9610        | 2.75" Coated Blade Electrode | Holster              | 50/bx, 2bx/cs, 100/cs |



#### ROCKER SWITCH

| PRODUCT NUMBER | FEATURES                     | ACCESSORIES          | QUANTITIES            |
|----------------|------------------------------|----------------------|-----------------------|
| 96-9611        | 2.75" Blade Electrode        | –                    | 50/Box                |
| 96-9612        | 2.75" Blade Electrode        | –                    | 50/bx, 2bx/cs, 100/cs |
| 96-9613        | 2.75" Blade Electrode        | Holster              | 50/Box                |
| 96-9614        | 2.75" Blade Electrode        | Holster              | 50/bx, 2bx/cs, 100/cs |
| 96-9615        | 2.75" Blade Electrode        | Holster, Tip Cleaner | 50/Box                |
| 96-9616        | 2.75" Blade Electrode        | Holster, Tip Cleaner | 50/bx, 2bx/cs, 100/cs |
| 96-9617        | 2.75" Coated Blade Electrode | –                    | 50/Box                |
| 96-9618        | 2.75" Coated Blade Electrode | –                    | 50/bx, 2bx/cs, 100/cs |
| 96-9619        | 2.75" Coated Blade Electrode | Holster              | 50/Box                |
| 96-9620        | 2.75" Coated Blade Electrode | Holster              | 50/bx, 2bx/cs, 100/cs |



#### FOOT SWITCH


| PRODUCT NUMBER | FEATURES                     | ACCESSORIES          | QUANTITIES            |
|----------------|------------------------------|----------------------|-----------------------|
| 96-9621        | 2.75" Blade Electrode        | –                    | 50/Box                |
| 96-9622        | 2.75" Blade Electrode        | –                    | 50/bx, 2bx/cs, 100/cs |
| 96-9623        | 2.75" Blade Electrode        | Holster              | 50/Box                |
| 96-9624        | 2.75" Blade Electrode        | Holster              | 50/bx, 2bx/cs, 100/cs |
| 96-9625        | 2.75" Blade Electrode        | Holster, Tip Cleaner | 50/Box                |
| 96-9626        | 2.75" Blade Electrode        | Holster, Tip Cleaner | 50/bx, 2bx/cs, 100/cs |
| 96-9627        | 2.75" Coated Blade Electrode | –                    | 50/Box                |
| 96-9628        | 2.75" Coated Blade Electrode | –                    | 50/bx, 2bx/cs, 100/cs |
| 96-9629        | 2.75" Coated Blade Electrode | Holster              | 50/Box                |
| 96-9630        | 2.75" Coated Blade Electrode | Holster              | 50/bx, 2bx/cs, 100/cs |



## ELECTROSURGERY

OB/GYN

### Sklar Blue™ Hooks



#### IRIS HOOKS

##### QUICK RELEASE

|         |    |          |         |     |          |
|---------|----|----------|---------|-----|----------|
| 91-5244 | 9" | 1mm      | 91-5248 | 12" | 1mm      |
| 91-5243 | 9" | 3mm      | 91-5247 | 12" | 3mm      |
| 91-5242 | 9" | 4mm      | 91-5246 | 12" | 4mm      |
| 91-5241 | 9" | 6mm      | 91-5245 | 12" | 6mm      |
| 91-5236 | 9" | Set of 4 | 91-5238 | 12" | Set of 4 |

#### EMMETT HOOKS

|         |    |          |
|---------|----|----------|
| 91-5251 | 9" | Style 1  |
| 91-5252 | 9" | Style 2  |
| 91-5253 | 9" | Style 3  |
| 91-5254 | 9" | Style 4  |
| 91-5255 | 9" | Style 5  |
| 91-5250 | 9" | Set of 5 |

## OB/GYN

## ELECTROSURGERY

### Sklar Blue™ Forceps



**BOZEMAN FORCEP**

91-5310 10¼" S-Curve, One Large Ring



**BOZEMAN FORCEPS**

91-5307 10½" Straight  
91-5308 10¼" S-Curve



**ADSON FORCEPS**

91-5428 4¾" Serrated



**DRESSING FORCEPS**

91-5420 8" Serrated  
91-5422 10" Serrated  
91-5424 12" Serrated



**POTTS-SMITH FORCEPS**

91-5411 7" Serrated  
91-5412 10" Serrated  
91-5413 8¼" 1x2 Teeth  
91-5414 10" 1x2 Teeth



**TUTTLE/ SINGLEY FORCEPS**

91-5300 7" Tuttle  
91-5302 9" Singley



## ELECTROSURGERY

OB/GYN

### Sklar Blue™ Forceps



**SCHWEIZER FORCEP**

91-5404 9½" Straight



**CLEMETSON FORCEP**

WITHOUT RATCHET

91-5400 9½" Straight



**LAUE UTERINE POLYP FORCEP**

91-5408 7¾"



**SCHROEDER TENACULUM  
FORCEPS**

91-5450 10" Braun Flat Tips  
91-5452 10" Pozzi Round Tips



**BABCOCK FORCEPS**

91-5268 8"  
91-5270 9½"



**ALLIS TISSUE FORCEPS**

91-5278 7½" 5x6 Teeth  
91-5440 9½" 5x6 Teeth

## OB/GYN

## ELECTROSURGERY

### Sklar Blue™ Endocervical Speculums



**KOGAN ENDOCERVICAL SPECULUM**

WITH GAUGE AND SCREW

91-5230 Narrow Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**

WITH SPRING HANDLE


91-5228 Narrow Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**

WITH RATCHET HANDLE


91-5235 Narrow Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**

WITH GAUGE AND SCREW

91-5232 Wide Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**

WITH SPRING HANDLE

91-5226 Wide Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**

WITH RATCHET HANDLE

91-5234 Wide Wire Blades



## ELECTROSURGERY

OB/GYN

### Sklar Blue™ Retractors



**LATERAL VAGINAL RETRACTOR  
WITH RATCHET**

91-5190 Fits Inside Speculum



**SIMS VAGINAL RETRACTORS**

**SINGLE END**

91-5200 Medium  
91-5202 Medium with loop to hold and  
direct disposable Smoke  
Evacuation Tubing (not illustrated)



**SIMS VAGINAL RETRACTORS**

**DOUBLE END**

91-5210 Medium  
91-5212 Medium with loop to hold and  
direct disposable Smoke  
Evacuation Tubing (not illustrated)



**JACKSON VAGINAL RETRACTORS  
SINGLE END**

|         |           |        |
|---------|-----------|--------|
| 91-5215 | 3" x 1½"  | Small  |
| 91-5218 | 3½" x 1½" | Medium |
| 91-5220 | 4" x 1½"  | Large  |



**TAUBER VAGINAL SPATULA**

91-5176 3½"



**SECOND VAGINAL SPATULA**

91-5178 4¾"

## OB/GYN

## ELECTROSURGERY

### Sklar Blue™ Speculums



#### AUVARD VAGINAL SPECULUMS

##### MEDIUM

- 91-5120 2 3/4" x 1 1/2" 2 1/2 lbs. with permanent Smoke Evacuation Tube (not illustrated)
- 91-5122 2 3/4" x 1 1/2" 2 1/2 lbs. with loop to hold and direct disposable Smoke Evacuation Tubing



#### AUVARD VAGINAL SPECULUMS

##### MEDIUM LONG BLADE

- 91-5130 4" x 1 3/4" 2 1/2 lbs. without loop
- 91-5132 4" x 1 3/4" 2 1/2 lbs. with loop to hold and direct disposable Smoke Evacuation Tubing



#### STEINER-AUVARD VAGINAL SPECULUMS

##### MEDIUM EXTRA LONG BLADE

- 91-5140 5 1/2" x 1 1/4" 2 1/2 lbs. without loop
- 91-5142 5 1/2" x 1 1/4" 2 1/2 lbs. with loop to hold and direct disposable Smoke Evacuation Tubing



#### COLLIN VAGINAL SPECULUMS

- 91-5180 4 1/4" x 1 1/4" Medium, with permanent Smoke Evacuation Tube (not illustrated)
- 91-5185 4 1/4" x 1 1/4" Medium, with loops to hold and direct disposable Smoke Evacuation Tubing

## MONSEL'S & LUGOL'S

### IN-OFFICE GYNECOLOGICAL ELECTROSURGERY




MONSEL'S PASTE Item #96-1200  
LUGOL'S SOLUTION Item #96-1210  
8ml Ready to use bottles

## ELECTROSURGERY

OB/GYN

### Sklar Blue™ Speculums



**GRAVES VAGINAL SPECULUM  
WIDE VIEW PATTERN  
WITH LARGE PROXIMAL OPENING FOR  
INSTRUMENTATION**

91-5080 4" x 1 1/4" Medium, with loops  
to hold and direct  
disposable Smoke  
Evacuation Tubing



**GRAVES OB-VIEW VAGINAL  
SPECULUMS  
WITH PERMANENT SMOKE EVACUATION TUBE  
(NOT ILLUSTRATED)**

91-3740 4" x 1 1/2" Medium  
91-3741 4 3/4" x 1 1/2" Large



**FREEWAY GRAVES VAGINAL  
SPECULUMS  
OPEN SIDE, 90° ANGLE**

91-5160 4" x 1 1/4" Medium, with permanent  
Smoke Evacuation Tube  
(not illustrated)  
91-5161 4" x 1 1/4" Medium, without  
permanent Smoke  
Evacuation Tube



**FREEWAY GRAVES VAGINAL  
SPECULUM  
OPEN SIDE, 45° ANGLE**

91-5162 4" x 1 1/4" Medium, with  
permanent Smoke  
Evacuation Tube  
(not illustrated)



**CUSCO VAGINAL SPECULUMS  
WITH PERMANENT SMOKE EVACUATION TUBE**

91-5150 3" x 1 1/4" Small  
91-5152 3 1/2" x 1 1/2" Medium  
91-5154 3 3/4" x 1 3/8" Large



**DISPOSABLE SMOKE  
EVACUATION TUBING**

93-2160 10" Precut  
3/16" Inside Diameter  
1/4" Outside Diameter 50/cs

## OB/GYN

## ELECTROSURGERY

### Sklar Blue™ Speculums



#### GRAVES VAGINAL SPECULUMS WITH LOOPS FOR DISPOSABLE SMOKE EVACUATION TUBING

|         |                 |             |
|---------|-----------------|-------------|
| 91-5015 | 3" x 1/2"       | Extra Small |
| 91-5010 | 3" x 3/4"       | Small       |
| 91-5012 | 4" x 1 1/4"     | Medium      |
| 91-5014 | 4 1/2" x 1 3/8" | Large       |
| 91-5018 | 4 3/4" x 1 1/2" | Extra Large |



#### GRAVES VAGINAL SPECULUMS WITH PERMANENT SMOKE EVACUATION TUBE

|         |                 |             |
|---------|-----------------|-------------|
| 91-5020 | 3" x 3/4"       | Small       |
| 91-5022 | 4" x 1 1/4"     | Medium      |
| 91-5024 | 4 1/2" x 1 3/8" | Large       |
| 91-5028 | 4 3/4" x 1 1/2" | Extra Large |




#### GRAVES VAGINAL SPECULUMS WITHOUT TUBE OR LOOPS

|         |                 |             |
|---------|-----------------|-------------|
| 91-5066 | 3" x 1/2"       | Extra Small |
| 91-5060 | 3" x 3/4"       | Small       |
| 91-5062 | 4" x 1 1/4"     | Medium      |
| 91-5064 | 4 1/2" x 1 3/8" | Large       |
| 91-5068 | 4 3/4" x 1 1/2" | Extra Large |



#### PEDERSON VAGINAL SPECULUMS WITH LOOPS FOR DISPOSABLE SMOKE EVACUATION TUBING

|         |               |             |
|---------|---------------|-------------|
| 91-5038 | 2 3/4" x 1/2" | Extra Small |
| 91-5032 | 3" x 1/2"     | Small       |
| 91-5040 | 4" x 5/8"     | Slender     |
| 91-5034 | 4" x 7/8"     | Medium      |
| 91-5036 | 4 3/4" x 1"   | Large       |



#### PEDERSON VAGINAL SPECULUMS WITH PERMANENT SMOKE EVACUATION TUBE

|         |                 |             |
|---------|-----------------|-------------|
| 91-5042 | 3" x 1/2"       | Small       |
| 91-5050 | 4" x 5/8"       | Slender     |
| 91-5044 | 4" x 7/8"       | Medium      |
| 91-5046 | 4 3/4" x 1"     | Large       |
| 91-5047 | 4 3/4" x 1 5/8" | Extra Large |



#### PEDERSON VAGINAL SPECULUMS WITHOUT TUBE OR LOOPS

|         |               |             |
|---------|---------------|-------------|
| 91-5088 | 2 3/4" x 1/2" | Extra Small |
| 91-5082 | 3" x 1/2"     | Small       |
| 91-5092 | 4" x 5/8"     | Slender     |
| 91-5084 | 4" x 7/8"     | Medium      |
| 91-5086 | 4 3/4" x 1"   | Large       |



## ENDOCERVICAL SPECULA

Sklar Blue™ Electrosurgical

OB/GYN



**KOGAN ENDOCERVICAL SPECULUM**  
**WITH GAUGE AND SCREW**

91-5230 Narrow Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**  
**WITH SPRING HANDLE**

91-5228 Narrow Wire Blades




**KOGAN ENDOCERVICAL SPECULUM**  
**WITH RATCHET HANDLE**

91-5235 Narrow Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**  
**WITH GAUGE AND SCREW**

91-5232 Wide Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**  
**WITH SPRING HANDLE**

91-5226 Wide Wire Blades



**KOGAN ENDOCERVICAL SPECULUM**  
**WITH RATCHET HANDLE**

91-5234 Wide Wire Blades

## OB/GYN

### RETRACTORS

#### Sklar Blue™ Electrosurgical



##### LATERAL VAGINAL RETRACTOR

###### RATCHET HANDLE

91-5190 Blade fits inside Vaginal Speculum



##### SIMS VAGINAL RETRACTORS

###### SINGLE END

91-5200 Medium  
91-5202 Medium, with loop to hold and direct Disposable Smoke Evacuation Tubing



##### SIMS VAGINAL RETRACTORS

###### DOUBLE END

91-5210 Medium  
91-5212 Medium, with loop to hold and direct Disposable Smoke Evacuation Tubing



##### JACKSON VAGINAL RETRACTORS

###### SINGLE END

91-5215 3" x 1½" Small  
91-5218 3½" x 1½" Medium  
91-5220 4" x 1½" Large



##### TAUBER VAGINAL SPATULA

91-5176 10" Blade 3½" x 1⅝"



##### SECOND VAGINAL SPATULA



91-5178 10½" Blade 4¾" x 1⅝"



## VAGINAL RETRACTORS

### SIMS VAGINAL RETRACTOR

| ITEM #        | DESCRIPTION                                       |
|---------------|---------------------------------------------------|
| HNM-1-52-2020 | SMALL, HOLLOW GRIP HANDLE, 2 1/2" x 1" BLADES     |
| HNM-1-52-2021 | MEDIUM, HOLLOW GRIP HANDLE, 3" x 1 1/4" BLADES    |
| HNM-1-52-2022 | LARGE, HOLLOW GRIP HANDLE, 3 1/2" x 1 1/2" BLADES |



### LEEP COATED LATERAL WALL RETRACTOR

| ITEM #        | DESCRIPTION                                  |
|---------------|----------------------------------------------|
| HNM-36-150757 | CLOSE SHANKS 22mm x 57mm LONG BLADES         |
| HNM-36-150758 | CLOSE SHANKS, 22mm x 75mm LONG BLADES, LARGE |







## LEEP COATED GRAVES VAGINAL SPECULUM

| ITEM #        | DESCRIPTION                           |
|---------------|---------------------------------------|
| HNM-36-150770 | W/SMOKE EVACUATION TUBE, SMALL        |
| HNM-36-150768 | W/SMOKE EVACUATION TUBE, MEDIUM       |
| HNM-36-150769 | W/SMOKE EVACUATION TUBE, LARGE        |
| HNM-36-150773 | W/SMOKE EVACUATION TUBE GUIDE, SMALL  |
| HNM-36-150774 | W/SMOKE EVACUATION TUBE GUIDE, MEDIUM |
| HNM-36-150775 | W/SMOKE EVACUATION TUBE GUIDE, LARGE  |

## LEEP COATED PEDERSON VAGINAL SPECULUM

| ITEM #        | DESCRIPTION                                              |
|---------------|----------------------------------------------------------|
| HNM-36-151534 | W/SMOKE EVACUATION TUBE, SMALL                           |
| HNM-36-150776 | W/SMOKE EVACUATION TUBE, MEDIUM, 7/8" x 4" BLADES        |
| HNM-36-150767 | W/SMOKE EVACUATION TUBE, LARGE                           |
| HNM-36-151535 | W/SMOKE EVACUATION TUBE GUIDE, SMALL                     |
| HNM-36-151536 | W/SMOKE EVACUATION TUBE GUIDE, MEDIUM                    |
| HNM-36-150777 | W/SMOKE EVACUATION TUBE GUIDE, LARGE, 1" x 4 3/4" BLADES |





## DISPOSABLE SQUARE LOOP ELECTRODES

| ITEM #    | DESCRIPTION                    |
|-----------|--------------------------------|
| H-ES14    | 5mm x 5mm, STERILE, BOX OF 5   |
| H-ES16    | 10mm x 8mm, STERILE, BOX OF 5  |
| H-OB1010S | 10mm x 10mm, STERILE, BOX OF 5 |




## DISPOSABLE BALL ELECTRODES

| ITEM # | DESCRIPTION            |
|--------|------------------------|
| H-OB3B | 3mm, STERILE, BOX OF 5 |
| H-ES07 | 5mm, STERILE, BOX OF 5 |



## ELECTRODES



### DISPOSABLE LOOP ELECTRODES

| ITEM #    | DESCRIPTION                         |
|-----------|-------------------------------------|
| H-ES08    | 5mm x 5mm LOOP, STERILE, BOX OF 5   |
| H-ES46    | 10mm x 5mm, STERILE, BOX OF 5       |
| H-OB-1007 | 10mm x 7mm, STERILE, BOX OF 5       |
| H-ES52    | 10mm x 8mm, STERILE, BOX OF 5       |
| H-ES09    | 10mm x 10mm LOOP, STERILE, BOX OF 5 |
| H-OB-1507 | 15mm x 7mm, STERILE, BOX OF 5       |
| H-ES11    | 15mm x 10mm LOOP, STERILE, BOX OF 5 |
| H-ES53    | 15mm x 12mm, STERILE, BOX OF 5      |
| H-ES43    | 15mm x 15mm, STERILE, BOX OF 5      |
| H-ES12    | 20mm x 8mm LOOP, STERILE, BOX OF 5  |
| H-ES49    | 20mm x 10mm, STERILE, BOX OF 5      |
| H-ES42    | 20mm x 12mm, STERILE, BOX OF 5      |
| H-ES13    | 20mm x 15mm LOOP, STERILE, BOX OF 5 |
| H-ES31    | 20mm x 20mm LOOP, STERILE, BOX OF 5 |
| H-ES47    | 25mm x 10mm LOOP, STERILE, BOX OF 5 |

## OB/GYN

## FORCEPS

### Sklar Blue™ Electrosurgical



**ADSON FORCEP**  
SERRATED  
91-5428 4 3/4"



**DRESSING FORCEPS**  
SERRATED  
91-5420 8"  
91-5422 10"  
91-5424 12"



**POTTS-SMITH FORCEPS**  
91-5411 7" Serrated  
91-5412 10" Serrated  
91-5413 8 1/4" 1x2 Teeth  
91-5414 10" 1x2 Teeth



**TUTTLE/SINGLEY FORCEPS**  
91-5300 7" Tuttle  
91-5302 9" Singley

## SKLAR FACE MASKS

The Sklar Disposable Face Mask is offered in the most popular style to meet clinical requirements in the operating room, general medicine, or any other areas requiring facial protection. Protecting against infectious pathogens, Sklar masks are ≥98% BFE, non-woven, soft and breathable. Convenient ear loops and adjustable nose-bridge ensures a proper fit. Latex Free.

### NON-STERILE

| PRODUCT NUMBER | FEATURES                | QUANTITIES |
|----------------|-------------------------|------------|
| 96-5557        | Mask Earloop, Blue      | 50/bx      |
| 96-5557A       | Mask Earloop, Blue      | 500/cs     |
| 96-5557BNS     | Mask Earloop, Blue Bulk | 2000/cs    |





# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



## Waugh-Brophy STILLE isolerad diatermipincett | Waugh-Brophy STILLE insulated Diathermy Forceps

Insulated version of the popular Waugh-Brophy forceps. Used in monopolar diathermy with 1x2 teeth.



102-9852

|             |           |                |           |
|-------------|-----------|----------------|-----------|
| 102-9852-20 | 1x2 teeth | Cross-serrated | 20 cm 8"  |
| 102-9852-25 | 1x2 teeth | Cross-serrated | 25 cm 10" |
| 102-9852-30 | 1x2 teeth | Cross-serrated | 30 cm 12" |

## Björk-STILLE diatermipincetter | Björk-STILLE Diathermy Forceps

Insulated forceps for monopolar diathermy. Serrations allow grasping of delicate structures.




102-9845

|             |           |           |
|-------------|-----------|-----------|
| 102-9845-15 | Insulated | 15 cm 6"  |
| 102-9845-20 | Insulated | 20 cm 8"  |
| 102-9845-25 | Insulated | 25 cm 10" |

## Waugh-Brophy STILLE isolerad diatermipincett | Waugh-Brophy STILLE insulated Diathermy Forceps

Insulated version of the popular Waugh-Brophy forceps. Used in monopolar diathermy. With Diamond Jaws.




102-9855

|             |              |           |
|-------------|--------------|-----------|
| 102-9855-20 | Diamond jaws | 20 cm 8"  |
| 102-9855-25 | Diamond jaws | 25 cm 10" |






## ELECTRODES




### DISPOSABLE LOOP ELECTRODES

| ITEM #    | DESCRIPTION                         |
|-----------|-------------------------------------|
| H-ES08    | 5mm x 5mm LOOP, STERILE, BOX OF 5   |
| H-ES46    | 10mm x 5mm, STERILE, BOX OF 5       |
| H-OB-1007 | 10mm x 7mm, STERILE, BOX OF 5       |
| H-ES52    | 10mm x 8mm, STERILE, BOX OF 5       |
| H-ES09    | 10mm x 10mm LOOP, STERILE, BOX OF 5 |
| H-OB-1507 | 15mm x 7mm, STERILE, BOX OF 5       |
| H-ES11    | 15mm x 10mm LOOP, STERILE, BOX OF 5 |
| H-ES53    | 15mm x 12mm, STERILE, BOX OF 5      |
| H-ES43    | 15mm x 15mm, STERILE, BOX OF 5      |
| H-ES12    | 20mm x 8mm LOOP, STERILE, BOX OF 5  |
| H-ES49    | 20mm x 10mm, STERILE, BOX OF 5      |
| H-ES42    | 20mm x 12mm, STERILE, BOX OF 5      |
| H-ES13    | 20mm x 15mm LOOP, STERILE, BOX OF 5 |
| H-ES31    | 20mm x 20mm LOOP, STERILE, BOX OF 5 |
| H-ES47    | 25mm x 10mm LOOP, STERILE, BOX OF 5 |



## DISPOSABLE SQUARE LOOP ELECTRODES

| ITEM #    | DESCRIPTION                    |
|-----------|--------------------------------|
| H-ES14    | 5mm x 5mm, STERILE, BOX OF 5   |
| H-ES16    | 10mm x 8mm, STERILE, BOX OF 5  |
| H-OB1010S | 10mm x 10mm, STERILE, BOX OF 5 |



## DISPOSABLE BALL ELECTRODES

| ITEM # | DESCRIPTION            |
|--------|------------------------|
| H-OB3B | 3mm, STERILE, BOX OF 5 |
| H-ES07 | 5mm, STERILE, BOX OF 5 |





## Pederson Vu-More Speculum, (40mm) Opening

- 038100 Medium 4½" (114mm)
- 038101 Large 4¾" (124mm)
- 038102 Extra Large 5" (149mm)



## Pederson Vu-Max Speculum, (70mm) Opening

- 038110 Medium 4½" (114mm)
- 038111 Large 4¾" (124mm)
- 038112 Extra Large 5" (149mm)



## Weisman Open Side Speculum Closed Tube

- 038050 Right Side Open Medium 1⅝" x 4" (35mm x 102mm)
- 038050-ANG Right Side Open Medium Angled 1⅝" x 4" (35mm x 102mm)
- 038052 Right Side Open Large 1⅝" x 4¾" (35mm x 121mm)
- 038052-ANG Right Side Open Large Angled 1⅝" x 4¾" (35mm x 121mm)
- 038051 Left Side Open Medium 1⅝" x 4" (35mm x 102mm)
- 038051-ANG Left Side Open Medium Angled 1⅝" x 4" (35mm x 102mm)
- 038053 Left Side Open Large 1⅝" x 4¾" (35mm x 121mm)
- 038053-ANG Left Side Open Large Angled 1⅝" x 4¾" (35mm x 121mm)



## Collins Speculum Closed Tube

- 038047 Small 1¼" x 3¾" (32mm x 95mm)
- 038048 Medium 1⅝" x 4" (35mm x 102mm)
- 038049 Large 1⅝" x 4½" (41mm x 114mm)



## Four-Way Lateral Speculum

- 038093 Medium ⅞" x 4¼" (22mm x 108mm)
- 038094 Large ⅞" x 5¼" (22mm x 133mm)



## Lateral Vaginal Wall Retractor

- 038060 Closed, 10" (254mm), Blade Size 1" x 2½" (25mm x 64mm)
- 038061 Open ½" (13mm) 10" (254mm), Blade Size 1" x 2½" (25mm x 64mm)
- 038062 Open 1" (25mm) 10" (254mm) Blade Size 1" x 2½" (25mm x 64mm)
- 038063 Wide Open, 10" (254mm) Blade Size 1" x 2½" (25mm x 64mm)



## Klopfer Speculum

- 038054 1⅝" x 3½" (35mm x 89mm)

## ELECTROSURGERY

### MedGyn Lletz Instruments - Blue Coated

All MedGyn Lletz instruments have a durable coating to insulate against transmission of electric current, ensuring the highest level of protection for patients and physicians.



#### Graves Speculum Closed

- 038030 Small, ¾" x 3" (19mm x 76mm)
- 038031 Medium, 1⅜" x 4" (35mm x 102mm)
- 038032 Large, 1⅝" x 4¾" (35mm x 121mm)
- 038033 Extra Large, 1½" x 6½" (38mm x 165mm)
- 030900-LTZ Extra Long Snowman, Clinton 1½" x 6" (38mm x 152mm)



#### Graves Vu-More Speculum Closed

- 038055 Medium 1¼" x 4" (32mm x 102mm)
- 038056 Large 1⅜" x 4½" (35mm x 114mm)



#### Graves Vu-Max Speculum Closed Tube

- 038046 Large 1⅜" x 4½" (35mm x 114mm)



#### Graves Vu-More Speculum, (40mm) Opening

- 038115 Small 3⅝" (98mm)
- 038116 Medium 4½" (114mm)
- 038117 Large 4⅞" (124mm)
- 038118 Extra Large 5⅞" (149mm)
- 038119 Extra Large 6⅞" (161mm)
- 038120 Extra Large 6⅞" (169mm)
- 038121 Extra Large 7" (178mm)




#### Graves Vu-Max Speculum, (70mm) Opening

- 038125 Small 3⅝" (98mm)
- 038126 Medium 4½" (114mm)
- 038127 Large 4⅞" (124mm)
- 038128 Extra Large 5⅞" (149mm)
- 038129 Extra Large 6⅞" (161mm)
- 038130 Extra Large 6⅞" (169mm)
- 038131 Extra Large 7" (178mm)



## MedGyn Disposable Electrodes

Made of high-quality Tungsten wire, for superior shape and integrity, MedGyn offers a full range of shapes and sizes to suit most electrosurgical procedures.



- 028067 Disposable Loop Electrode 5mm x 5mm (5/Pack)
- 028069 Disposable Loop Electrode 10mm x 5mm (5/Pack)
- 028070 Disposable Loop Electrode 10mm x 7mm (5/Pack)
- 028071 Disposable Loop Electrode 10mm x 8mm (5/Pack)
- 028055 Disposable Loop Electrode 10mm x 10mm (5/Pack)
- 028056 Disposable Loop Electrode 10mm x 13mm (5/Pack)
- 028057 Disposable Loop Electrode 15mm x 5mm (5/Pack)
- 028087 Disposable Loop Electrode 15mm X 7mm (5/Pack)
- 028058 Disposable Loop Electrode 15mm x 8mm (5/Pack)
- 028068 Disposable Loop Electrode 15mm x 10mm (5/Pack)
- 028076 Disposable Loop Electrode 15mm x 13mm (5/Pack)
- 028059 Disposable Loop Electrode 15mm x 15mm (5/Pack)
- 028060 Disposable Loop Electrode 20mm x 8mm (5/Pack)
- 028063 Disposable Loop Electrode 20mm x 10mm (5/Pack)
- 028061 Disposable Loop Electrode 20mm x 12mm (5/Pack)
- 028062 Disposable Loop Electrode 20mm x 13mm (5/Pack)
- 028064 Disposable Loop Electrode 20mm x 15mm (5/Pack)
- 028065 Disposable Loop Electrode 20mm x 20mm (5/Pack)
- 028066 Disposable Loop Electrode 25mm x 10mm (5/Pack)

- 028041 Disposable Ball Electrode 3mm (5/Pack)
- 028042 Disposable Ball Electrode 5mm (5/Pack)
- 028081 Disposable Blade Electrode (5/Pack)
- 028053 Disposable Needle 5mm (5/Pack)


- 028075 Disposable Square Loop Electrode 5mm x 5mm (5/Pack)
- 028072 Disposable Square Loop Electrode 10mm x 4mm (5/Pack)
- 028073 Disposable Square Loop Electrode 10mm x 8mm (5/Pack)
- 028074 Disposable Square Loop Electrode 10mm x 10mm (5/Pack)



# PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)


## ELECTROSURGERY



### Lletz Kogan Endocervical Speculum

038057 Speculum 2-4mm, 11" (279mm)

038058 Speculum 4-8mm, 11" (279mm)



### Lletz Tenaculum


038070 Single-Tooth Tenaculum 9 3/4" (248mm)



### Lletz Tenaculum Hook

038080 Emmet Tenaculum 9" (229mm)



038081 Iris Tenaculum Hook 9" (229mm)



### Lletz Dressing Forceps

038082 Lletz Dressing Forceps Small 8" (203mm)


038092 Lletz Dressing Forceps Large 10" (254mm)



### Lletz Tissue Forceps

038090 Lletz Tissue Forceps Small 8" (203mm)

038091 Lletz Tissue Forceps Large 10" (254mm)



### Lletz Foerster Forceps

038071 Lletz Foerster Forceps - Straight 9 1/2" (241mm)

038072 Lletz Foerster Forceps - Curved 9 1/2" (241mm)

