

PRIME GYNECOLOGY & OBSTETRICS CATALOG - B



PRIME SURGICAL
GYNECOLOGY
INTERNATIONAL



PRIME SURGICAL
INTERNATIONAL



PRIME SURGICAL
INTERNATIONAL

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PRIME SURGICAL INTERNATIONAL

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



AL

PRIME SURGICAL

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



Zia ur Rahman
0092-300-2621673
0092-313-5710887

PRIME
Surgical International
Deals in All Kinds of Quality Instruments

Ph: +92 213 5643935
E.mail id: prime.instruments@hotmail.com
zia@primeinstrumentintl.com www.primeinstrumentintl.com
Office # 25, Ground Floor, Queen Center M.T. Khan Road, Karachi.

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PRIME SURGICAL INTERNATIONAL

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PRIME SURGICAL INTERNATIONAL

CARE AND CLEANING

OB/GYN

SKLAR INSTRUMENT CARE KIT

A Trio of Sklar's Care and Cleaning products to keep instruments in tiptop shape.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1625		
	Sklar Lube™ 8 oz. Spray Bottle	1 each
	Sklar Kleen™ Liquid 8 oz. Bottle	1 each
	Sklar Polish™ 8 oz. Jar	1 each
	Stainless Steel Bristle Cleaning Brush	1 each
	Nylon Bristle Cleaning Brush	1 each



TROUBLE SHOOTING

NO STEEL IS TRULY "STAINLESS"

It's a fact that so-called stainless steel is subject to both water-spotting and staining. However what many people identify as a "rust" problem usually turns out to be something else.

Most often the problem is actually a stain caused by a surface deposit. These deposits appear in a variety of colors depending upon the type of deposit involved. Once the cause is found and eliminated, most of the time the problem quickly diminishes or disappears completely.

As you probably know, stainless steel is ideally suited for the surgical suite because it is rust resistant, it can be honed to an extremely sharp edge or fine point,

and it can be hardened to maintain the delicate, yet precise, requirements of the surgeon. But it is this hardenable quality – the amount of carbon used during the manufacturing process – coupled with harsh chemicals and improper care that renders stainless steel susceptible to corrosion.

To protect your instrument investment, Sklar cleaning agents have been specifically formulated to provide cost effective instrument care solutions that help maintain the life and performance of your valuable stainless steel instruments.

OB/GYN

CARE AND CLEANING

SKLAR FRESH™

LIQUID SPRAY ODOR ABSORBER

SKLAR FRESH is an environmentally friendly blend of pure biodegradable organic essential oils and trace elements. It is a non-allergenic liquid spray that neutralizes odors at their source. Simply spray SKLAR FRESH directly in the air or on the source (i.e. bathrooms, refuse, soiled items from vomitus, incontinence odors, etc.)

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2842	22 oz. Spray Bottles	12/case

SKLAR POLISH™

INSTRUMENT STAIN REMOVER

Non-Abrasive

SKLAR POLISH is an effective stain remover for use on all your metal surgical instruments, including stainless steel, chrome and brass. It will brighten instrument trays, basins and racks and may also be used on porcelain and tile surfaces.



PRODUCT NUMBER	FEATURES	QUANTITIES
10-1626	8 oz. Jars	12/case
10-1629	8 oz. Jars	6/case
10-1927	8 oz. Jar	1 jar

ODOR ABSORBER / LIQUID SPRAY

SKLAR FRESH™



CARE AND CLEANING

OB/GYN

SKLAR STAT™

INSTANT FOAM HAND SANITIZER

KILLS 99% OF GERMS IN LESS THAN 15 SECONDS!

SKLAR-STAT Instant Foam Hand Sanitizer, based on the active ingredient Benzalkonium Chloride, is a unique formulation featuring exceptional skin feel, conditioning and moisturizing properties. The efficacy of this product has been confirmed to reduce S. Aureus 99.999% and S. Aureus MRSA 99.9% in as little as 15 seconds. SKLAR-STAT Instant Foaming Hand Sanitizer produces a fast drying, non-sticky foam that contains unique conditioning and moisturizing ingredients, leaving the skin with a soft, silky after-feel.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2850	50ml Pump Dispensers	24/case
10-2852	550ml Pump Dispensers	12/case

SKLAR STAT™

AUTOMATIC FOAM DISPENSER

Motion activated, battery powered and wall-mounted, Sklar's AUTOMATIC FOAM DISPENSER allows for easy access and reduced risk of cross contamination. Ideal for use in hospitals, medical centers, offices, schools and any high-traffic public area with multiple users.



PRODUCT NUMBER	FEATURES	QUANTITIES
10-2856	Wall Mount Dispenser	1 Each
10-2854	1000ml Refill Bottles	6/case

INSTANT FOAM HAND SANITIZER

SKLAR-STAT™

- Alcohol-Free
- Non-Sticky
- Moisturizes and Conditions Hands

1. Pump a small amount into palm of hand.
2. Rub thoroughly over both surfaces of both hands.
3. Rub hands together briskly until dry.



OB/GYN

CARE AND CLEANING

3 EASY STEPS FOR THE CARE AND CLEANING OF YOUR INSTRUMENTS

DECONTAMINATION

1

SKLAR SPRAY DISINFECTANT™ Disinfectant / Decontaminant Cleaner

SKLAR DISINFECTANT is a pre-cleaner and disinfectant for the safe decontamination of soiled surgical instruments, medical and surgical equipment, and inanimate surfaces including plastics, stainless steel, chrome, glass, Plexiglas and baked on painted surfaces. The use of SKLAR DISINFECTANT for 2 minutes kills infectious organisms including H1N1, HIV-1, HPV, VRE, Norovirus, Hepatitis B & C, as well as TB, MRSA, A-Baumannii, and C-Diff in its vegetative form.



CLEANING

2

SKLAR KLEEN™ DETERGENTS AND ENZYMATICS

Sklar has an array of agents for effectively cleaning not only your surgical instruments, but for other metals, glassware, plastic, and rubber. SKLAR KLEEN detergent is available as a powder or liquid. Our enzymatic products, specially formulated for processing instruments heavily laden with organic debris, are concentrated and available in a variety of sizes. Sklar also offers SPRAY-ZYME™ a new foam spray pre-cleaner.



Choose The Cleaning Product Right For You!

LUBRICATION

3

SKLAR LUBRICANT AND SPRAY LUBE™

SKLAR LUBE should be a part of your instrument processing regimen to guarantee extending the life and performance of this very important investment. Sklar offers several different types of lubricant to suit your needs. All designed specifically for surgical instruments, they each contain rust inhibitors to protect against corrosion, rusting, and staining.



YOUR FINAL STEP IS AUTOCLAVE STERILIZATION

Your Instruments Will Look And Work Like New!

CARE AND CLEANING

OB/GYN

SKLAR ENZYMATIC™ DETERGENT

100% Biodegradable

SKLAR ENZYMATIC Detergent is a low-foaming, fast acting detergent specifically designed for use when instruments and/or cloth materials become soiled with organic debris.

SKLAR ENZYMATIC Detergent may be used as a manual instrument cleaner as well as a pre-soak for blood-contaminated linens, clothing and instruments. It is also used in ultrasonic units and in washer/sterilizer/decontaminators. When used as directed, the enzymatic pH neutral formula is safe for use on rubber, plastic, stainless steel instruments, and surgical equipment as well as cloth.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2777	1 Gallon Bottles	4/case
10-2775	1 Gallon Bottles	2/case
10-2776	1 Gallon Bottle	1 bottle



OB/GYN

CARE AND CLEANING

SKLAR KLEEN™

REGULAR FORMULA FOR MANUAL AND ULTRASONIC CLEANING

SKLAR KLEEN is non-corrosive, non aggressive, phosphate free and biodegradable. We have specifically developed SKLAR KLEEN for manual and ultrasonic cleaning. It is moderately alkaline, affording optimum cleaning without etching glassware or attacking instrument metal.

LOW FOAM FORMULA FOR AUTOMATIC CLEANING

The SKLAR KLEEN Low Foam Liquid is concentrated in order to provide a very cost-effective solution. It is especially formulated to maintain low foam, ensuring effective and trouble-free action in washer / sterilizers and tunnel washers. Dispersing quickly, the SKLAR KLEEN Low Foam formula is also well suited for automatic injection systems.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2701	Low Foam 1 Gallon Bottles	4/case
10-1608	Powdered 5 lb. Containers	4/case
10-1607	Powdered 5 lb. Container	1 container
10-1630	Powdered 3.5 lb. Pitchers	6/case
10-2765	Powdered 3.5 lb. Pitcher	1 pitcher
10-1612	Liquid 8 oz. Bottles	6/case
10-1614	Liquid 8 oz. Bottle	1 bottle
10-1613	Liquid 1 Gallon Bottles	4/case
10-1627	Liquid 1 Gallon Bottle	1 bottle



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PRIME SURGICAL INTERNATIONAL

OB/GYN

CARE AND CLEANING

SKLAR INSTRU-GUARD™ ONE STEP

A "ONE STEP" INSTRUMENT CLEANER AND LUBRICANT

EXTENDS INSTRUMENT LIFE!

SKLAR INSTRU-GUARD ONE STEP contains a penetrating lubricant that protects and helps extend the life of surgical instruments. The detergents and emulsifiers quickly break up and remove organic soils from box locks and other hard to reach places.

SKLAR INSTRU-GUARD ONE STEP contains four rust inhibitors to control corrosion and eliminate scaling and spotting. It does not contain phosphates and can reduce the time and cost involved in processing your instruments by cleaning and lubricating in one easy step.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2800	1 Gallon Bottles	4/case
10-2801	1 Gallon Bottle	1 bottle



CARE AND CLEANING

OB/GYN

SKLAR SOAK™

INSTRUMENT PRE-SOAK & DECONTAMINANT

RUST RESISTANT!

Immersing soiled instruments in SKLAR SOAK effectively pre-cleans and decontaminates, without fear of rusting!

SKLAR SOAK is an effective pre-cleaner that breaks down bio-burden, lifting blood and other organic debris from the surface, while penetrating box locks and crevices. The rust inhibitor protects the delicate chromium oxide layer of the instrument, thus defending against rust and corrosion. Using SKLAR SOAK™ as a pre-cleaner is the first line of defense against rust, prolonging the life of the instrument and therefore protecting your investment!

Soaking for only two minutes thoroughly decontaminates instruments, reducing the risk of exposure to pathogens such as H1N1, HIV-1, HPV, VRE, Norovirus, Hep B & C as well as TB, MRSA, A. Baumannii, and C. Diff in its vegetative form. However, because SKLAR SOAK™ contains a rust inhibitor, instruments may soak for several hours without fear of rusting! Employees are protected from cross contamination, and at the same time can enjoy the freedom to reprocess instruments when it's convenient!

SKLAR SOAK is ideal for use on all stainless steel instruments:

- PH balanced, designed specifically for stainless steel instruments.
- Intermediate level disinfectant with broad spectrum of efficacy.
- Rust resistant.
- Non-caustic.
- Ready to use.
- No special handling or disposal.
- Fragrant, no offensive odors.
- Non-glutaraldehyde.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1603	1 Gallon Bottles	4/case
10-1609	1 Gallon Bottle	1 bottle



OB/GYN

CARE AND CLEANING

SKLAR SPRAY-ZYME™

MULTI-ENZYME BASED INSTRUMENT PRESOAK & CLEANER

FOAM SPRAY CLEANER

The newest product from Sklar for expert care of your surgical instruments. This multi-enzymatic spray creates a heavy foam which lasts for hours, keeping soiled instruments moist to avoid encrustation while also pre-cleaning prior to processing.

- Prevents Organic Encrustation - enzyme foam stops blood and bodily fluids from drying onto surgical instruments.
- No More Spills - messy spills between the operating room and central cleaning /processing department are eliminated.
- Neutral pH - safe and effective for stainless steel, tungsten carbide, aluminum and brass instruments.
- Bacteriostatic - microbiological control agents inhibit bacterial growth in the ready-to-use foam solution.
- Long Lasting Foam - saves time, foam lasts many hours after application.
- Multi-Enzymatic - contains a synergistic blend of lipase and protease enzymes.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2722	22 oz. Refillable Spray Bottles	12/case
10-2724	1 Gallon Bottles	4/case

SKLAR-ZYME™ POUCHES

CONVENIENT VALUE

Each pouch of SKLAR-ZYME produces one gallon of enzymatic detergent. This convenient, easy to use product guarantees very little waste and is ideal for the office or smaller surgical setting.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2774	½ oz. Pouches	32/box



CARE AND CLEANING

OB/GYN

SKLAR DISINFECTANT™ SURFACE WIPES

2 Minute Kill Time!

**SAME EFFECTIVE FORMULA AS
SKLAR DISINFECTANT!**

- Ready to use, non-toxic formula.
- Specially formulated for the control of cross contamination.
- Effective against a broad range of pathogens including H1N1, HIV-1, HPV, VRE, Norovirus, Hep B & C as well as TB, MRSA, A. Baumannii, and C. Diff in its vegetative form.
- Safe for all areas of the facility including the NICU, operating rooms, isolation rooms, patient care areas and laboratories.
- Satisfies OSHA's Bloodborne Pathogen Standards.
- Destroys odor-causing bacteria.
- 7" x 10" non-woven, tear resistant towelette.
- Clean scent.
- EPA registered.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1616	100 Towelettes / Canister	6/case
10-1617	100 Towelettes / Canister	12/case



OB/GYN

CARE AND CLEANING

SKLAR DISINFECTANT™

READY-TO-USE SPRAY
SURFACE DISINFECTANT

2 Minute Kill Time!

- Ready to use, non-toxic formula.
- Specially formulated for the control of cross contamination.
- Effective against a broad range of pathogens including H1N1, HIV-1, HPV, VRE, Norovirus, Hep B & C as well as TB, MRSA, A. Baumannii, and C. Diff in its vegetative form. Safe for all areas of the facility including the NICU, operating rooms, isolation rooms, patient care areas and laboratories.
- Satisfies OSHA's Bloodborne Pathogen Standards.
- Non-irritating formula.
- EPA registered.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1643	24 oz. Spray Bottles	6/case
10-1642	24 oz. Spray Bottle	1 bottle
10-1653	1 Gallon Bottles	4/case
10-1654	1 Gallon Bottle	1 bottle
10-1644	2.5 Gallon Bottles	2/case



CARE AND CLEANING

OB/GYN

SKLAR-CIDE™

**SPORICIDAL • BACTERICIDAL
VIRUCIDAL • FUNGICIDAL**

SKLAR-CIDE provides a safe and efficient means for the high-level disinfection or sterilization of medical, ophthalmic, dental and surgical devices including endoscopes and lensed instruments. It may be used as a standard soak solution or in automatic reprocessing units set at room temperature. Unlike other solutions, there is no need to increase the solution temperature to provide high-level disinfection in 25 minutes.

SKLAR-CIDE is a surfactant-free, activated glutaraldehyde liquid chemical sterilant and high-level disinfectant. Its active ingredient is glutaraldehyde in a nominal concentration of 3%. When the powdered activator is added, the solution turns green in color to indicate its readiness for use. In addition to the active ingredient, this ready-to-use solution contains buffers, wetting agents, and a rust and corrosion inhibitor and it has a pleasant fragrance.

SKLAR-CIDE produces high-level disinfection at room temperature in as little time as 25 minutes, and complete sterilization in 10 hours at room temperature.

- Reusable for 28 days instead of 14, safer and saves time and money.
- Surfactant free - Will not damage scope surfaces or cloud lenses.
- Can be used safely on all devices, including instruments and endoscopes. No need to stock several products.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1623	28 HLD 1 Gallon Bottles	4/case
10-1624	28 HLD 1 Gallon Bottle	1 bottle

GLUTARALDEHYDE STRIPS

PRODUCT NUMBER	FEATURES	QUANTITIES
96-4105	2.1% MEC 60 Strips / Bottle	2 bottles/case



OB/GYN

BRUSHES

Sets


10-1424 AUTOCLAVE CLEANING BRUSH SET

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1424HDL	Telescopic Handle extends from 24" to 48"	1 Each
10-1424PBR	Large Cleaning Brush 7½" L x 5½" W with round corner bristles	1 Each
10-1424BHD	Surface Cleaning Brush, 4½" handle	1 Each



10-1432 MICRO-OPTHALMIC BRUSH SET

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1350	Cannula Cleaning Brush 12" 2mm, 3/pk	1 pk
10-1660	Cannula Cleaning Brush 12" 3mm, 3/pk	1 pk
10-1396	Suction Tube Brush 6" 1.6mm, 3/pk	1 pk
10-1400	Suction Tube Brush 8" 2mm, 3/pk	1 pk
10-1402	Suction Tube Brush 8" 2.6mm, 3/pk	1 pk
10-1404	Suction Tube Brush 8" 3.3mm, 3/pk	1 pk
10-1406	Suction Tube Brush 8" 4mm, 3/pk	1 pk
10-1444	Double End Brush Nylon, 3/pk	1 pk



10-1434 ALL PURPOSE BRUSH SET

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1650	Brush Stainless Steel Bristle 3/pk	1 pk
10-1651	Brush Nylon Bristle 3/pk	1 pk
10-1420	Large Instrument Cleaning Brush 3/pk	1 pk
10-1418	Surface Cleaning Brush 3/pk	1 pk



BRUSHES

OB/GYN

Sets


10-1426 INSTRUMENT CLEANING BRUSH SET

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1660	Cannula Cleaning Brush 12" 3mm, 3/pk	1 pk
10-1362	Cannula Cleaning Brush 18" 5mm, 3/pk	1 pk
10-1374	Cannula Cleaning Brush 24" 7mm, 3/pk	1 pk
10-1376	Cannula Cleaning Brush 24" 10mm, 3/pk	1 pk
10-1396	Suction Tube Brush 6" 1.6mm, 3/pk	1 pk
10-1400	Suction Tube Brush 8" 2mm, 3/pk	1 pk
10-1402	Suction Tube Brush 8" 2.6mm, 3/pk	1 pk
10-1404	Suction Tube Brush 8" 3.3mm, 3/pk	1 pk
10-1406	Suction Tube Brush 8" 4mm, 3/pk	1 pk
10-1418	Surface Cleaning Brush, 3/pk	1 pk
10-1444	Double End Brush Nylon, 3/pk	2 pks
10-1650	Brush Stainless Steel Bristle, 3/pk	1 pk
10-1651	Brush Nylon Bristle, 3/pk	2 pks



10-1428 ENDOSCOPIC BRUSH SET

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1350	Cannula Cleaning Brush 12" 2mm, 3/pk	1 pk
10-1660	Cannula Cleaning Brush 12" 3mm, 3/pk	1 pk
10-1354	Cannula Cleaning Brush 12" 5mm, 3/pk	1 pk
10-1362	Cannula Cleaning Brush 18" 5mm, 3/pk	1 pk
10-1364	Cannula Cleaning Brush 18" 10mm, 3/pk	1 pk
10-1370	Cannula Cleaning Brush 24" 2.5mm, 3/pk	1 pk
10-1372	Cannula Cleaning Brush 24" 5mm, 3/pk	1 pk
10-1376	Cannula Cleaning Brush 24" 10mm, 3/pk	1 pk
10-1398	Double End Valve Brush Nylon, 3/pk	2 pks
10-1444	Double End Brush Nylon, 3/pk	2 pks



OB/GYN

BRUSHES

General Cleaning and Accessories

SURFACE CLEANING BRUSH

FOR USE ON SINKS, AUTOCLAVES, ULTRASONIC CLEANERS, BOWLS, BASINS, ETC.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1418	3 ¼" L x 3" W Brush Head 4 ½" Plastic handle 2" White medical grade Nylon bristles	3/pk



LARGE INSTRUMENT CLEANING BRUSH

FOR USE ON RETRACTORS, SPECULA, ETC.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1420	5" L x 1" W Brush Head 5" Plastic handle 1" White medical grade Nylon bristles	3/pk



AUTOCLAVE CLEANING BRUSH

FOR USE ON AUTOCLAVE CHAMBERS AND DOORS, SINKS, CARTS, ETC.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1422	3 ¼" L x 3" W Brush Head 16" Plastic handle 2" White medical grade Nylon bristles	3/pk



INSTRUMENT SPOT REMOVER PAD


FOR USE ON O.R. GRADE SATIN FINISH INSTRUMENTS

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1414	6" x 9"	1 Each
10-1415	6" x 9"	10/pk



BRUSH DRYING RACK & BRUSH HOLDER

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1399	Brush Rack	1 Each
10-1545	Brush Holder	1 Each



BRUSHES

OB/GYN

Instrument Cleaning

SUCTION TUBE CLEANING BRUSHES (BARON & FRAZIER)

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1396	6" 1.6mm	3/pk
10-1400	8" 2mm	3/pk
10-1401	8" 2.3mm	3/pk
10-1402	8" 2.6mm	3/pk
10-1403	8" 3mm	3/pk
10-1404	8" 3.3mm	3/pk
10-1405	8" 3.6mm	3/pk
10-1406	8" 4mm	3/pk



SURGICAL SCRUB BRUSHES

PRODUCT NUMBER	FEATURES	QUANTITIES
06-2880	With Channel Back	1 Each
06-2870	With Channel Back	12/pk
06-2882	With Nail Scrubber	1 Each
06-2873	With Nail Scrubber	12/pk
06-2881	Standard	1 Each
06-2872	Standard	12/pk



CANNULA INSTRUMENT (PIPE) CLEANERS

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1470	12" Length 1/8" Diameter	3/pk
10-1472	12" Length 1/4" Diameter	3/pk



DOUBLE END VALVE BRUSHES

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1398	5 3/4" Length 9 mm and 5mm Diameter Nylon Tips	3/pk



© 2018

OB/GYN

BRUSHES

Instrument Cleaning

TRACHEAL TUBE BRUSHES

PRODUCT NUMBER	FEATURES	QUANTITIES
78-2821	Small, 6" x 1/2" Diameter	6/pk
78-2822	Medium, 7" x 5/8" Diameter	6/pk



INSTRUMENT CLEANING BRUSHES (STAINLESS BRISTLES)

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1650	Stainless Steel Bristles	3/pk
10-1652	Stainless Steel Bristles	12/pk
10-1658	Stainless Steel Bristles	24/pk



INSTRUMENT CLEANING BRUSHES (NYLON BRISTLES)

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1651	Nylon Bristles	3/pk
10-1657	Nylon Bristles	12/pk
10-1659	Nylon Bristles	24/pk
10-1465	Nylon Bristles, Disposable	12/pk
10-1468	Nylon Bristles, Disposable	24/pk



INSTRUMENT CLEANING BRUSHES (DOUBLE ENDED)

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1444	Nylon Bristles	3/pk



INSTRUMENT CLEANING BRUSHES (DOUBLE TAPERED)

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1384	6", Double Tapered Brush	3/pk
10-1388	8", Double Tapered Brush	3/pk



BRUSHES

Endoscopy

OB/GYN

CANNULA CLEANING BRUSHES

PRODUCT NUMBER	FEATURES		QUANTITIES
10-1350	12"	2mm Diameter	3/pk
10-1660	12"	3mm Diameter	3/pk
10-1661	12"	3mm Diameter	12/pk
10-1663	12"	3mm Diameter	24/pk
10-1352	12"	4mm Diameter	3/pk
10-1354	12"	5mm Diameter	3/pk
10-1356	12"	6mm Diameter	3/pk
10-1664	15"	9mm Diameter	3/pk
10-1665	15"	9mm Diameter	12/pk
10-1666	15"	9mm Diameter	24/pk
10-1668	15"	25mm Diameter	3/pk
10-1669	15"	25mm Diameter	12/pk
10-1672	15"	25mm Diameter	24/pk
10-1362	18"	5mm Diameter	3/pk
10-1364	18"	10mm Diameter	3/pk
10-1366	18"	15mm Diameter	3/pk
78-2824	19"	8mm Diameter	3/pk
10-1370	24"	2.5mm Diameter	3/pk
10-1372	24"	5mm Diameter	3/pk
10-1374	24"	7mm Diameter	3/pk
10-1376	24"	10mm Diameter	3/pk



OB/GYN

SklarTip™ Instrument Identification

Instrument recognition is key to efficient processing in a busy clinical setting. SklarTip™ Instrument ID products offer a convenient color coded system for identifying instruments. Offered in various combinations including tape sheets, tape rolls and tags, products last through many autoclaving and gas sterilization processes. Latex Free.

ID SHEET TAPE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-1845	Blue with White Stripes	1 Each
96-1846	Green with White Stripes	1 Each
96-1847	Red with White Stripes	1 Each
96-1848	Yellow with White Stripes	1 Each
96-1849	Orange with White Stripes	1 Each
96-1850	Brown with White Stripes	1 Each
96-1851	Lavender with White Stripes	1 Each
96-1852	Light Blue with White Stripes	1 Each
96-1853	Pink with White Stripes	1 Each



ID ROLL TAPE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-1860	White 300" x 1/4"	1 Each
96-1862	Green 300" x 1/4"	1 Each
96-1864	Yellow 300" x 1/4"	1 Each
96-1866	Orange 300" x 1/4"	1 Each
96-1868	Black 300" x 1/4"	1 Each
96-1870	Purple 300" x 1/4"	1 Each
96-1872	Blue 300" x 1/4"	1 Each
96-1874	Red 300" x 1/4"	1 Each
96-1876	Brown 300" x 1/4"	1 Each



INSTRUMENT ID TAGS

PRODUCT NUMBER	FEATURES	QUANTITIES
96-1900	ID Repair Tags, Red	100/pk
96-1901	ID Repair Tags, Red	1000/pk
96-1903	Sklar Tags, White	100/pk
96-1905	ID Tags for Sharpen, Blue	100/pk
96-1906	ID Tags for Sharpen, Blue	1000/pk



SklarTip™ Instrument Protectors

Sklar Instrument Protectors provide a convenient, color-coded means of protecting surgical instrument tips and edges. Autoclavable, steam and gas permeable, and available in a wide variety of shapes and sizes, SklarTip™ Instrument Protectors provide exceptional protection during handling and sterilization. Latex Free.



INSTRUMENT PROTECTORS - NON-VENTED

PRODUCT NUMBER	FEATURES	QUANTITIES
96-1510	19 x 1.5mm, 3/4" x 1/16"	White 100/pk
96-1508	19 x 2.0mm, 3/4" x 5/64"	Blue 100/pk
96-1506	19 x 3.0mm, 3/4" x 1/8"	Green 100/pk
96-1504	25 x 3.0mm, 1" x 1/8"	Red 100/pk
96-1502	25 x 5.0mm, 1" x 3/16"	Yellow 100/pk
96-1500	19 x 1.0mm, 3/4" x 3/64"	Maroon 100/pk
96-1516	25 x 1.5mm, 1" x 3/8" x 1/16"	Orange 100/pk
96-1514	25 x 2.0mm, 1" x 5/8" x 5/64"	Brown 100/pk
96-1512	25 x 3.0mm, 1" x 1" x 1/8"	Black 100/pk
96-1518	Assorted	100/pk

INSTRUMENT PROTECTORS - VENTED


PRODUCT NUMBER	FEATURES	QUANTITIES
96-1410	19 x 1.5mm, 3/4" x 1/16"	White 100/pk
96-1408	19 x 2.0mm, 3/4" x 5/64"	Blue 100/pk
96-1406	19 x 3.0mm, 3/4" x 1/8"	Green 100/pk
96-1404	25 x 3.0mm, 1" x 1/8"	Red 100/pk
96-1402	25 x 5.0mm, 1" x 3/16"	Yellow 100/pk
96-1400	19 x 1.0mm, 3/4" x 3/64"	Maroon 100/pk
96-1416	25 x 1.5mm, 1" x 3/8" x 1/16"	Orange 100/pk
96-1414	25 x 2.0mm, 1" x 5/8" x 5/64"	Brown 100/pk
96-1412	25 x 3.0mm, 1" x 1" x 1/8"	Black 100/pk
96-1418	Assorted	100/pk

OSTEOTOME GUARDS - VENTED

PRODUCT NUMBER	FEATURES	QUANTITIES
96-1440	3/4" x 1/4"	50/pk
96-1442	1" x 3/8"	50/pk
96-1444	1" x 1/2"	50/pk
96-1446	1 1/2" x 1 1/2"	50/pk

TIP GRIPS

PRODUCT NUMBER	FEATURES	QUANTITIES
96-1550	.062" Yellow	3pr/pk, 5pks
96-1552	.062" Yellow	5pr/pk, 5pks



OB/GYN

PLASTICWARE

Trays

PLASTIC TRAYS

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-7025	Small Shallow Prep	25/cs
96-7033	Half Shallow	25/cs
96-7024	Half Deep	25/cs
96-7175	Full Deep	25/cs
96-7177	2-Compartment	25/cs
96-7178	3-Compartment	25/cs



PROCEDURE TRAYS*

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1960	Mini, Non-Perforated, Blue, 10 1/8" x 6 5/8"	1 Each
10-1965	Mini, Perforated, Blue, 10 1/8" x 6 5/8"	1, 6 or 12/cs
10-1961	Small, Non-Perforated, Blue, 10 7/16" x 13 9/16"	1 Each
10-1966	Small, Perforated, Blue, 10 7/16" x 13 9/16"	1, 6 or 12/cs
10-1962	Medium, Non-Perforated, Blue, 11 3/4" x 16"	1 Each
10-1967	Medium, Perforated, Blue, 11 3/4" x 16"	1, 6 or 12/cs
10-1964	Large, Heavy Duty, Non-Perforated, Blue, 15" x 20"	1 Each



SOAKING TRAYS WITH REMOVABLE BASKET

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2555	Small Instrument Tray, 8" x 3 1/2" x 2"	1 Each
10-2556	Large Instrument Tray, 10" x 5" x 2 1/2"	1 Each
10-2557 (not pictured)	Round Cup Size Tray, 3" x 3", 4 oz.	1 Each



*Select products are available individually or in bulk packaging. To order large quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

PLASTICWARE

General

OB/GYN

EMESIS BASINS

STERILE		
PRODUCT NUMBER	FEATURES	QUANTITIES
96-7543	700cc	25/cs



OR BASINS

STERILE		
PRODUCT NUMBER	FEATURES	QUANTITIES
96-1679	Round, 7 qt., Blue	5/cs

NON-STERILE		
PRODUCT NUMBER	FEATURES	QUANTITIES
08-1011	Round, 5 qt., Pink	50/cs



*Select products are available individually or in bulk packaging. To order large quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

GUIDEWIRE BOWLS

STERILE		
PRODUCT NUMBER	FEATURES	QUANTITIES
96-1630	8 1/4" x 3"	10/cs



UTILITY BOWLS*

NON-STERILE		
PRODUCT NUMBER	FEATURES	QUANTITIES
10-1502	Sponge/Utility Bowl, Aqua, 1 1/4 qt.	1 or 12/cs



NON-STERILE		
PRODUCT NUMBER	FEATURES	QUANTITIES
10-1504	Solution Bowl, Reusable, Lt Green, 7 qt.	1 or 12/cs



OB/GYN

PLASTICWARE

General

MEDICINE CUPS

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-7027	2 oz.	25/cs



PITCHERS*

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-7268	Open Handle, 1200ml	5/cs

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-7270	Open Handle, 1200ml	35/cs



NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1546	Graduated Measure, 16 oz., Reusable	1 or 6/cs



SPECIMEN TUBES

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-7459	3 1/2" x 1/2"	3/pk, 50pk/cs, 150/cs



PLASTIC BOWLS

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-7028	8 oz.	25/cs
96-7020	16 oz.	10/cs
96-7174	32 oz.	20/cs
96-7451	48 oz.	5/cs



*Select products are available individually or in bulk packaging. To order large quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

STEELWARE

OB/GYN

Sterilization Containers

SKLARLITE™

Sterilization Container System

The SklarLite™ Rigid Sterilization Container System provides instruments an organized, routine and secure approach to the sterilization of surgical instruments. Customize your SklarLite™ System by selecting colored perforated container lids, solid safety lids and ID labels. Available in your choice of six vibrant colors for easy identification of instrument sets by procedure and/or department. Choose from Silver, Yellow, Green, Blue, Black and Red.

BASKETS

SklarLite™ features stackable inner baskets to optimize internal space while silicone mats securely protect delicate and laparoscopic instruments during sterilization, transport and storage.

PRODUCT NUMBER	MODEL	FEATURES
10-4784	Half Size	10" x 9½" x 1"
10-4786	Half Size	10" x 9½" x 2"
10-4788	Half Size	10" x 9½" x 3"
10-4790	Half Size	10" x 9½" x 4"
10-4478	Mid Size	16" x 10" x 1"
10-4480	Mid Size	16" x 10" x 2"
10-4482	Mid Size	16" x 10" x 3"
10-4484	Mid Size	16" x 10" x 4"
10-4284	Full Size	21¼" x 10" x 1"
10-4286	Full Size	21¼" x 10" x 2"
10-4288	Full Size	21¼" x 10" x 3"
10-4290	Full Size	21¼" x 10" x 4"

SILICONE MATS

PRODUCT NUMBER	MODEL	FEATURES
10-4970	Half Size	8½" x 9"
10-4975	Mid Size	15" x 9"
10-4985	Full Size	20½" x 9"



ACCESSORIES

PRODUCT NUMBER	MODEL	FEATURES
10-4900	Process Indicator Cards	Steam, 250/box
10-4910	Filters	Disposable, 100/pack
10-4925	Security Seals	100/pack
10-4930	ID Labels	Silver
10-4930ENG	ID Labels	Silver Custom Engraved, Maximum 16 Characters



OB/GYN

STEELWARE

Sterilization Containers

SKLARLITE™

Sterilization Container System

HALF SIZE

NON-PERFORATED BOTTOM	PERFORATED BOTTOM	MODEL	FEATURES
10-4500	10-4512	Standard	11¼" x 11" x 4"
10-4548	10-4560	Standard	11¼" x 11" x 5"
10-4596	10-4608	Standard	11¼" x 11" x 6"
10-4644	10-4656	Standard	11¼" x 11" x 8"
10-4692	10-4704	Standard	11¼" x 11" x 10"
10-4524	10-4536	Safe	11¼" x 11" x 4"
10-4572	10-4584	Safe	11¼" x 11" x 5"
10-4620	10-4632	Safe	11¼" x 11" x 6"
10-4668	10-4680	Safe	11¼" x 11" x 8"
10-4716	10-4728	Safe	11¼" x 11" x 10"

MID SIZE

NON-PERFORATED BOTTOM	PERFORATED BOTTOM	MODEL	FEATURES
10-4300	10-4312	Standard	18¼" x 11" x 4"
10-4348	10-4360	Standard	18¼" x 11" x 5"
10-4394	10-4406	Standard	18¼" x 11" x 6"
10-4324	10-4336	Safe	18¼" x 11" x 4"
10-4370	10-4382	Safe	18¼" x 11" x 5"
10-4418	10-4430	Safe	18¼" x 11" x 6"

FULL SIZE

NON-PERFORATED BOTTOM	PERFORATED BOTTOM	MODEL	FEATURES
10-4000	10-4012	Standard	23" x 11" x 4"
10-4048	10-4060	Standard	23" x 11" x 5"
10-4096	10-4108	Standard	23" x 11" x 6"
10-4144	10-4156	Standard	23" x 11" x 8"
10-4192	10-4204	Standard	23" x 11" x 10"
10-4024	10-4036	Safe	23" x 11" x 4"
10-4072	10-4084	Safe	23" x 11" x 5"
10-4120	10-4132	Safe	23" x 11" x 6"
10-4168	10-4180	Safe	23" x 11" x 8"
10-4216	10-4228	Safe	23" x 11" x 10"

Containers and accessories for the SklarLite™ Sterilization Container System are available in a variety of sizes and colors. Contact your Sklar Account Manager for information on the complete product line.

SKLARLITE™

Sterilization Container System

The SklarLite™ Rigid Sterilization Container System is a reusable rigid container made from anodized aluminum alloy, designed to be strong and durable yet lightweight and easy to handle.

- Stackable Containers
- Stackable Baskets
- 6 Month Shelf Life
- Load Limit 25lbs.



TWO MODELS TO CHOOSE FROM:

STANDARD

Consists of a body and perforated lid and used with Sklar's disposable filters and other compatible accessories for sterilization. Choose from either perforated or non-perforated bottoms.

SAFE

Consists of identical body and lid, but features an additional solid cover that safeguards the perforation field and contents during transportation and storage.

STEELWARE

OB/GYN

Instrument Stringers

WEISENBACH INSTRUMENT HOLDER

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1775	6"	1 Each
10-1776	12"	1 Each



WORCESTER INSTRUMENT HOLDER

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1839	13" x 3"	1 Each



BUNT INSTRUMENT STRINGER

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1770		1 Each




SAUNDERS INSTRUMENT RACK

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1780		1 Each



MAYO INSTRUMENT HOLDER

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1772		1 Each



POSILOCK INSTRUMENT RACK

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1763		1 Each



OB/GYN

STEELWARE

Instrument Stringers

WEINSTEIN INSTRUMENT STRINGERS

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1800	5½" x 3¼"	1 Each
10-1805	11¼" x 3¼"	1 Each



INSTRUMENT STRINGERS

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1340	2½" x 6"	1 Each
10-1342	2½" x 8"	1 Each
10-1344	2½" x 10"	1 Each
10-1346	2½" x 12"	1 Each
10-1347	2½" x 14"	1 Each
10-1348	2½" x 16"	1 Each
10-1349	2½" x 18"	1 Each



E-Z LOCK INSTRUMENT STRINGERS

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1806	2½" x 6"	1 Each
10-1808	2½" x 8"	1 Each
10-1810	2½" x 10"	1 Each
10-1812	2½" x 12"	1 Each
10-1814	2½" x 14"	1 Each
10-1816	2½" x 16"	1 Each
10-1818	2½" x 18"	1 Each
10-1826	5" x 6"	1 Each
10-1828	5" x 8"	1 Each
10-1830	5" x 10"	1 Each
10-1832	5" x 12"	1 Each
10-1834	5" x 14"	1 Each
10-1836	5" x 16"	1 Each
10-1838	5" x 18"	1 Each



STEELWARE Instrument Trays

OB/GYN

LIFTING HANDLES

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1725	For Sterilization Trays	1 Pair



FLASH STERILIZATION TRAY

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1711	10" x 6½" x 2"	1 Each
10-1712	Single Handle for Flash Sterilization Tray 10-1711	1 Each




DOUBLE DECKER MICRO TRAY

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1717	10" x 15" x 2¼"	1 Each



MICRO TRAY WITH HINGED LID

PRODUCT NUMBER	DESCRIPTION	Quantities
10-1706	8½" x 8" x 1½" with 2½" x 2½" Compartment	1 Each
10-1707	8½" x 12" x 1½"	1 Each
10-1708	8½" x 14" x 1½"	1 Each
10-1702	10½" x 10" x 1½"	1 Each
10-1716	10½" x 20" x 2½" with 10½" x 5" Compartment	1 Each



PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1700	8½" x 8" x 1¼"	1 Each
10-1703	10½" x 15" x 1¼"	1 Each



MICRO TRAY WITH HINGED LID (Not pictured)

PRODUCT NUMBER	DESCRIPTION	Quantities
10-1701	10" x 6½" x 7⁄8"	1 Each

STERILIZATION BOX (Not pictured)

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1705	6¾" x 3¾" x 1¾"	1 Each

OB/GYN

STEELWARE

Instrument Trays

STERILIZATION TRAY / BATH SETS

PRODUCT NUMBER	DESCRIPTION			QUANTITIES
10-1947	Large Set Includes:			
	Solid Tray	10-1948	20¾" x 12¾" x 6"	1
	Perforated Tray	10-1949	20¾" x 12¾" x 4"	1
	Cover	10-1535	20¾" x 12¾"	1
10-1679	Small Set Includes:			
	Solid Tray	10-1747	12¾" x 10½" x 4"	1
	Perforated Tray	10-1697	12¾" x 10½" x 2½"	1
	Cover	10-1882	12¾" x 10½"	1



STERILIZATION TRAYS*

PRODUCT NUMBER	DESCRIPTION	COVERS	QUANTITIES
10-1850	4 qt., 12¾" x 10½" x 2½"	10-1864	1 Each
10-1860	7 qt., 12¾" x 10½" x 4"	10-1864	1 Each
10-1750	5 ¼ qt., 20¾" x 6¾" x 4"	10-1854	1, 6 or 12/box
10-1752	9 qt., 20¾" x 6¾" x 6"	10-1854	1, 6 or 12/box
10-1953	8 ½ qt., 20¾" x 12¾" x 2½", Heavyweight	10-1535	1, 6 or 12/box
10-1954	14 ½ qt., 20¾" x 12¾" x 4", Heavyweight	10-1535	1, 6 or 12/box
10-1948	24 qt., 20¾" x 12¾" x 6"	10-1535	1, 6 or 12/box



PERFORATED STERILIZATION TRAYS

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1710	10" x 6½" x 2½"	1 Each
10-1713	10" x 10½" x 3½"	1 Each
10-1714	15" x 10½" x 3½"	1 Each
10-1720	20" x 10½" x 3½"	1 Each
10-1957	20¾" x 12¾" x 4" Mesh bottom	1 Each



*All products are packaged individually. Instrument trays and covers are sold separately unless otherwise noted. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.



STEELWARE Instrument Trays

OB/GYN

INSTRUMENT TRAYS, PERFORATED BOTTOM

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1740	8 $\frac{3}{8}$ " x 5" x 2"	1 Each
10-1874	10" x 6 $\frac{1}{2}$ " x 2"	1 Each
10-1892	12" x 7 $\frac{3}{4}$ " x 2"	1 Each
10-1879	12" x 7 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ "	1 Each
10-1758	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	1 Each
10-1862	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 4"	1 Each
10-1894	16 $\frac{1}{2}$ " x 10" x 2 $\frac{1}{2}$ "	1 Each
10-1895	20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ "	1 Each

INSTRUMENT TRAYS, PERFORATED BOTTOM / SIDES*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1757	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 1 $\frac{1}{4}$ "	1 Each
10-1697	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	1 Each
10-1943	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 4"	1 Each
10-1944	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 6"	1 Each
10-1941	20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 1 $\frac{1}{4}$ "	1 Each
10-1759	20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ "	1 Each
10-1949	20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 4"	1, 6 or 12/box
10-1942	20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 6"	1 Each

INSTRUMENT TRAYS, SOLID*

PRODUCT NUMBER	DESCRIPTION	COVERS	QUANTITIES
10-1755	6 $\frac{7}{8}$ " x 4 $\frac{1}{4}$ " x 2"	10-1737	1, 6 or 12/box
10-1741	6 $\frac{7}{8}$ " x 6 $\frac{3}{8}$ " x 4" Heavyweight	10-1835	1 or 6/box
10-1764	10" x 6 $\frac{1}{2}$ " x 2"	10-1766	1 Each
10-1742	10 $\frac{3}{8}$ " x 6 $\frac{3}{8}$ " x 2 $\frac{1}{2}$ "	10-1542	1, 6 or 12/box
10-1877	12" x 7 $\frac{3}{4}$ " x 2 $\frac{1}{4}$ "	10-1878	1, 6 or 12/box
10-1744	12 $\frac{3}{4}$ " x 6 $\frac{7}{8}$ " x 2 $\frac{1}{2}$ "	10-1724	1, 6 or 12/box
10-1746	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	10-1882	1, 6 or 12/box
10-1950	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ " Heavyweight	10-1882	1, 6 or 12/box
10-1747	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 4"	10-1882	1, 6 or 12/box
10-1952	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 4" Heavyweight	10-1882	1, 6 or 12/box
10-1870	12 $\frac{3}{4}$ " x 10 $\frac{1}{2}$ " x 6"	10-1882	1, 6 or 12/box
10-1753	13 $\frac{7}{8}$ " x 12 $\frac{7}{8}$ " x 2 $\frac{1}{2}$ "	10-1754	1, 6 or 12/box
10-1748	16 $\frac{1}{2}$ " x 10" x 2 $\frac{1}{2}$ "	10-1756	1, 6 or 12/box
10-1749	16 $\frac{1}{2}$ " x 10" x 4"	10-1756	1 Each
10-1953	20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ " Heavyweight	10-1535	1, 6 or 12/box
10-1954	20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 4" Heavyweight	10-1535	1, 6 or 12/box
10-1948	20 $\frac{3}{4}$ " x 12 $\frac{3}{4}$ " x 6"	10-1535	1, 6 or 12/box



*All products are packaged individually. Instrument trays and covers are sold separately unless otherwise noted. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

OB/GYN

STEELWARE

Instrument Trays

MESH CLEANING BASKET*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1729	4¼" x 7¼" x 2"	1 Each
10-1731	6" x 9" x 3"	1 Each
10-1732	8½" x 12½" x 3½"	1 Each
10-1733	9½" x 16" x 4"	1 Each
10-1734	12" x 21" x 5"	1 Each
10-1736	14" x 20" x 5"	1 or 12/box



PERFORATED CLEANING BASKET

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1560	12¾" x 10⅝" x 6"	1 Each



INSTRUMENT DRYING TRAYS*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1876	10" x 6½" x ¾"	1, 6 or 12/box
10-1883	13⅛" x 9⅜" x ¾"	1, 6 or 12/box
10-1884	15⅝" x 10⅝" x ¾"	1, 6 or 12/box
10-1683	17" x 11⅝" x ¾"	1, 6 or 12/box
10-1888	19⅝" x 12⅝" x ¾"	1, 6 or 12/box
10-1680	21¼" x 16¼" x ⅞"	1, 6 or 12/box



STEELWARE Instrument Trays

OB/GYN

*All products are packaged individually. Instrument trays and covers are sold separately unless otherwise noted. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

MAYO TRAYS (SOLID)*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1690	10" x 6½" x ¾"	1, 6 or 12/box
10-1691	13⅞" x 9⅜" x ¾"	1, 6 or 12/box
10-1692	15⅞" x 10⅞" x ¾"	1, 6 or 12/box
10-1693	17" x 11⅞" x ¾"	1, 6 or 12/box
10-1694	19 ⅞" x 12⅞" x ¾"	1, 6 or 12/box
10-1695	21¼" x 16¼" x ⅞"	1, 6 or 12/box



MAYO TRAYS (PERFORATED)*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1684	10" x 6½" x ¾"	1, 6 or 12/box
10-1685	13¾" x 9¾" x ¾"	1, 6 or 12/box
10-1686	15⅞" x 10⅞" x ¾"	1, 6 or 12/box
10-1687	17" x 11⅞" x ¾"	1, 6 or 12/box
10-1688	19⅞" x 12⅞" x ¾"	1, 6 or 12/box
10-1689	21¼" x 16¼" x ⅞"	1 or 6/box



CATHETER INSTRUMENT TRAYS*

PRODUCT NUMBER	DESCRIPTION	COVERS	QUANTITIES
10-1955	8¾" x 3¼" x 1½"	10-1956 - Dome	1, 6 or 12/box
10-1738	8⅞" x 5" x 2"	10-1722 - Dome	1, 6 or 12/box
10-1760	8⅞" x 5" x 3"	10-1722 - Dome	1, 6 or 12/box
10-1886	17¼" x 4½" x 2"	10-1887 - Flat	1, 6 or 12/box



OB/GYN

STEELWARE

Surgical / Lab Products

OINTMENT JAR WITH LID*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1525	9 oz., 3 $\frac{1}{8}$ " x 2 $\frac{3}{4}$ "	1 or 12/box



MINI TRAY WITH LID

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1768	8 $\frac{1}{2}$ " x 3" x 1 $\frac{1}{2}$ "	1 Each





CONVENIENT
DISPENSING
CANISTER

SEE PAGE 265 FOR
ORDERING DETAILS

SKLAR DISINFECTANT™ SURFACE WIPES

BACTERICIDAL • VIRUCIDAL • FUNGICIDAL • TUBERCULOCIDAL

3 MINUTE KILL TIME!

- READY TO USE, NON-TOXIC, NON-IRRITATING FORMULA
- SPECIALLY FORMULATED FOR THE CONTROL OF CROSS CONTAMINATION
- EFFECTIVE AGAINST C.DIFF, TB, HBV, HCV, VRE, HYDROPHILIC AND LIPOPHILIC VIRUSES INCLUDING HIV-1 (AIDS VIRUS), BACTERIA AND FUNGI ON INANIMATE OBJECTS
- SAFE FOR ALL AREAS OF THE FACILITY INCLUDING THE NICU, OPERATING ROOMS, ISOLATION ROOMS, PATIENT CARE AREAS, LABORATORIES AND DIALYSIS
- SATISFIES OSHA'S BLOODBORNE PATHOGEN STANDARDS
- DESTROYS ODOR-CAUSING BACTERIA
- 7" X 10" NON-WOVEN, TEAR RESISTANT TOWELETTE
- CLEAN SCENT
- EPA REGISTERED

STEELWARE

OB/GYN

Surgical / Lab Products

DRESSING JAR WITH LID*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1523	½ qt., 4½" x 2¼"	1 or 12/box
10-1533	1 ¼ qt., 4" x 5"	1 or 12/box
10-1538	1 ¾ qt., 4½" x 6 ¼"	1 or 12/box
10-1527	3 qt., 5½" x 6¾"	1 or 12/box
10-1543	4¾ qt., 6¾" x 7¾"	1 or 12/box
10-1544	6 qt., 7¼" x 8½"	1, 6 or 12/box
10-1547	8 qt., 8" x 9½"	1, 6 or 12/box



FORCEP/THERMOMETER JAR*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1550	1½" x 3"	1, 6 or 12/box
10-1541	2¼" x 4½"	1, 6 or 12/box
10-1545	2¼" x 7"	1, 6 or 12/box



FUNNELS WITH HANGING HOOK

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1157	2½ oz., 3" x 3¾"	1 Each
10-1158	6½ oz., 4½" x 4¾"	1 Each
10-1159	13 oz., 5½" x 5¾"	1 Each
10-1160	22¼ oz., 6¼" x 6¾"	1 Each
10-1170	32 oz., 7" x 8¾"	1 Each
10-1171	64 oz., 8½" x 9½"	1 Each



LAB / UTILITY SCOOP

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1192	5 oz., 5" x 2½" (Overall length 9½")	1 Each



OB/GYN

STEELWARE

Surgical / Lab Products

EMESIS BASINS*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1497	10 oz., 6 $\frac{1}{8}$ " x 3 $\frac{1}{2}$ " x 1 $\frac{1}{2}$ "	1 or 12/box
10-1698	12 oz., 8 $\frac{1}{8}$ " x 3 $\frac{3}{8}$ " x 1 $\frac{1}{8}$ "	1 or 12/box
10-1699	26 oz., 9 $\frac{7}{8}$ " x 4 $\frac{1}{2}$ " x 2 $\frac{1}{8}$ "	1 or 12/box
10-1487	48 oz., 12 $\frac{1}{8}$ " x 5 $\frac{3}{4}$ " x 2 $\frac{1}{4}$ "	1 Each



HOSPITAL BUCKETS

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1182	13 qt., 11 $\frac{1}{8}$ " x 9 $\frac{1}{4}$ " w/o lid	1 Each
10-1450	13 qt., 11 $\frac{1}{8}$ " x 9 $\frac{1}{4}$ " with lid	1 Each
10-1183	16 qt., 11 $\frac{3}{4}$ " x 10 $\frac{1}{8}$ " w/o lid	1 Each
10-1451	16 qt., 11 $\frac{3}{4}$ " x 10 $\frac{1}{8}$ " with lid	1 Each
10-1184	20 qt., 13 $\frac{3}{4}$ " x 11 $\frac{13}{16}$ " w/o lid	1 Each



HOSPITAL KICK BUCKET WITHOUT LID

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1551	13 qt., 11 $\frac{1}{8}$ " x 9 $\frac{1}{4}$ "	1 Each



STORAGE CONTAINERS*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1540	$\frac{1}{2}$ gal., 5 $\frac{7}{8}$ " x 6 $\frac{5}{8}$ " no handles	1, 6 or 12/box
10-1186	3 gal., 10 $\frac{3}{8}$ " x 9"	1 Each
10-1188	4 gal., 13" x 9"	1 Each
10-1189	5 gal., 13" x 11"	1 Each
10-1190	6 gal., 13" x 13"	1 Each



*All products are packaged individually. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

STEELWARE

OB/GYN

Surgical / Lab Products

GRADUATED MEASURE*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1552	16 oz., 3 $\frac{3}{4}$ " x 4 $\frac{5}{8}$ "	1, 6 or 12/box
10-1510	32 oz., 4 $\frac{3}{4}$ " x 5 $\frac{1}{4}$ "	1, 6 or 12/box
10-1674	64 oz., 5 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ "	1, 6 or 12/box



SPONGE BOWLS*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1496	12 oz., 4 $\frac{1}{4}$ " x 2"	1 or 12/box
10-1498	22 oz., 5" x 2 $\frac{3}{8}$ "	1 or 12/box
10-1484	29 oz., 5 $\frac{9}{16}$ " x 2 $\frac{5}{16}$ "	1 or 12/box
10-1485	36 oz., 6 $\frac{1}{4}$ " x 2 $\frac{3}{4}$ "	1 or 12/box
10-1483	48 oz., 7 $\frac{1}{8}$ " x 2 $\frac{3}{4}$ "	1 or 12/box
10-1486	64 oz., 7 $\frac{1}{8}$ " x 3"	1 or 12/box



MIXING/SOLUTION BOWLS*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1488	$\frac{7}{8}$ qt., 6 $\frac{5}{8}$ " x 2 $\frac{1}{2}$ "	1 or 12/box
10-1491	1 $\frac{1}{8}$ qt., 8 $\frac{1}{8}$ " x 3"	1 or 12/box
10-1492	3 qt., 9 $\frac{3}{4}$ " x 4"	1 or 12/box
10-1493	5 qt., 11 $\frac{1}{8}$ " x 4 $\frac{1}{2}$ "	1 or 12/box
10-1494	6 $\frac{3}{4}$ qt., 12 $\frac{3}{4}$ " x 5 $\frac{1}{4}$ "	1, 6 or 12/box
10-1495	8 qt., 13 $\frac{1}{8}$ " x 5 $\frac{7}{8}$ "	1, 6 or 12/box



SOLUTION BASINS*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1489	4 $\frac{1}{2}$ qt., 12" x 4 $\frac{1}{4}$ "	1 or 12/box
10-1505	7 qt., 13 $\frac{3}{8}$ " x 4 $\frac{1}{8}$ "	1, 6 or 12/box
10-1507	9 $\frac{1}{2}$ qt., 14 $\frac{1}{2}$ " x 5"	1, 6 or 12/box



OB/GYN

STEELWARE

Surgical / Lab Products

All Sklar Stainless Steelware is produced according to the strictest standards to ensure superior performance and durability.

*All products are packaged individually. Select products are available in bulk. To order larger quantities, simply add the package quantity to the end of the product number when ordering, ie. 10-1670-12 for 12/pk or 10-1670-48 for 48/pk.

IODINE CUP*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1670	6 oz., 3 ⁷ / ₁₆ " x 2"	1, 12, 24 or 48/box
10-1671	14 oz., 4 ³ / ₁₆ " x 2 ⁵ / ₁₆ "	1, 12, 24 or 48/box



MEDICINE CUP*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1675	2 oz., 2 ¹ / ₁₆ " x 2" Graduated Measure	1, 12, 24 or 48/box



TUMBLER*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1536	7 oz., 2 ¹ / ₁₆ " x 3 ¹ / ₁₆ "	1, 12/box



PITCHER / BEAKER*

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1172	125ml, 2 ¹ / ₁₆ " x 2 ¹ / ₁₆ " w/o Handle	1 Each
10-1173	250ml, 2 ³ / ₁₆ " x 3 ¹ / ₁₆ " w/o Handle	1 Each
10-1174	600ml, 3 ³ / ₁₆ " x 4 ¹ / ₁₆ " w/o Handle	1 Each
10-1176	2000ml, 4 ⁷ / ₁₆ " x 7 ¹ / ₁₆ " w/o Handle	1 Each
10-1177	3000ml, 5 ¹ / ₂ " x 8 ¹ / ₄ " w/o Handle	1 Each
10-1178	4000ml, 6 ¹ / ₁₆ " x 8 ⁷ / ₁₆ " w/o Handle	1 Each
10-1175	1200ml, 4 ¹ / ₁₆ " x 6 ¹ / ₁₆ " with closed Handle	1 Each
10-1539	2000ml, 5 ¹ / ₁₆ " x 6 ³ / ₁₆ " with closed Handle	1, 6 or 12/box
10-1180	3000ml, 6 ¹ / ₁₆ " x 7" with closed Handle	1 Each
10-1181	4000ml, 6 ⁷ / ₁₆ " x 8 ⁷ / ₁₆ " with closed Handle	1 Each



BLADES AND SCALPELS

Specialty

OB/GYN

MINI BLADES AND SCALPELS

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
97-1822	#61	25/box
97-1721	#62	12/box
97-1821	#62	25/box
97-1820	#63	25/box
97-1826	#64	25/box
97-1825	#65	25/box
97-1824	#67	25/box
97-1823	#69	25/box



MINI PODIATRY HANDLE

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
97-1450	#3K	1 Each



SPECIALTY BLADES

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
97-1400	#318, Curved Edge	12/box
97-1405	#312, Curved Edge	12/box
97-1410	#313, Curved Edge	12/box
97-1415	#314, Curved Edge	12/box
97-1420	#316, Straight Edge	12/box
97-1425	#300, Straight Edge	12/box
97-1430	#361, Knife Edge	12/box



SELF-LOCKING HANDLE

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
97-1445	#6	1 Each



CARBON P-M BLADE

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3024	#60	50/box



MINI THUMB SCALPEL

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-7052	#11	1500/cs



STITCH CUTTER

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3030	Stainless Steel	100/box



OB/GYN

BLADES AND SCALPELS

SklarSafe™ Safety Scalpels



SKLARSAFE SAFETY SCALPES

- Protective, transparent sheath can be opened, closed, and locked with a single hand.
- Precision stainless steel blade made in England.
- Sturdy, "Sure-Grip" handle.
- Audible "click" when sheath advanced to locked position.
- OSHA Compliant.
- Individually sealed in sturdy film and Tyvek® pack.

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3094	#10	10/box
06-3095	#11	10/box
06-3096	#15	10/box

NON-STERILE

06-3130	#10	500/cs
06-3132	#11	500/cs
06-3134	#15	500/cs

BLADES AND SCALPELS

OB/GYN

General



BLADES, STAINLESS STEEL

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3100	#10	100/box
06-3101	#11	100/box
06-3102	#12	100/box
06-3103	#15	100/box
06-3108	#15C	100/box
06-3104	#20	100/box
06-3105	#21	100/box
06-3106	#22	100/box
06-3107	#23	100/box

BLADES, CARBON STEEL

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3010	#10	100/box
06-3011	#11	100/box
06-3012	#12	100/box
06-3013	#12D	100/box
06-3015	#15	100/box
06-3016	#15C	100/box
06-3020	#20	100/box
06-3021	#21	100/box
06-3022	#22	100/box
06-3023	#23	100/box

SCALPELS, DISPOSABLE

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3110	#10	10/box
06-3111	#11	10/box
06-3112	#12	10/box
06-3115	#15	10/box
06-3120	#20	10/box
06-3121	#21	10/box
06-3122	#22	10/box
06-3123	#23	10/box



OB/GYN

BLADES AND SCALPELS

Handles

SCALPEL HANDLES

PRODUCT NUMBER	FEATURES
----------------	----------

06-2910	#3 with metric scale
---------	----------------------



PRODUCT NUMBER	FEATURES
----------------	----------

06-2903	#3 fits Blade #10, 11, 12, 15
---------	-------------------------------

06-2915	#3L
---------	-----



PRODUCT NUMBER	FEATURES
----------------	----------

06-2904	#4 fits Blade #20, 21, 22, 23
---------	-------------------------------

06-2920	#4L
---------	-----



PRODUCT NUMBER	FEATURES
----------------	----------

06-2907	#7 fits Blade #10, 11, 12, 15
---------	-------------------------------



PRODUCT NUMBER	FEATURES
----------------	----------

06-2912	#8 fits Blade #60
---------	-------------------



PRODUCT NUMBER	FEATURES
----------------	----------

06-2925	Large Lab Type Handle #3
---------	--------------------------

06-2930	Large Lab Type Handle #4
---------	--------------------------



PRODUCT NUMBER	FEATURES
----------------	----------

06-2908	Round Scalpel Handle, Straight #3R, fits Blade #10, 11, 12, 15
---------	---



PRODUCT NUMBER	FEATURES
----------------	----------

06-3090	Blade Remover
---------	---------------



PROTECTIVE APPAREL

OB/GYN

SKLAR DISPOSABLE FACE MASKS

The Sklar Disposable Face Mask is offered in the most popular style to meet clinical requirements in the operating room, general medicine, or any other areas requiring facial protection. Protecting against infectious pathogens, Sklar masks are $\geq 98\%$ BFE, non-woven, soft and breathable. Convenient ear loops and adjustable nose-bridge ensures a proper fit. Latex Free.

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-5557	Mask Earloop, Blue	50/bx
96-5557A	Mask Earloop, Blue	500/cs
96-5557BNS	Mask Earloop, Blue, Bulk	2000/cs



SKLAR SURGICAL GOWNS

- Low linting, lightweight and breathable SMS fabric.
- Generous cut provides comfort and freedom of movement during various procedures.
- Superior liquid barrier protection provides confidence against strikethrough.
- Velcro-style neck closure easily adjusts for a perfect fit.
- Meets safety standards for Class I Flame Resistance.
- Meets AAMI Level I requirements.
- Disposable, single use.

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
96-5329	Non-Reinforced, Large	40/cs
96-5331	Non-Reinforced, X-Large	40/cs
96-5333	Non-Reinforced, XX-Large	40/cs

BULK, NON-STERILE

96-5003B	Non-Reinforced, Large	160/cs
96-5005B	Non-Reinforced, X-Large	160/cs
96-5006B	Non-Reinforced, XX-Large	120/cs



OB/GYN

STETHOSCOPES

General

BOWLES STETHOSCOPE

PRODUCT NUMBER

06-1000



PERRELL FETOSCOPE

PRODUCT NUMBER

06-1701

FEATURES

With Slip on Tubing Adaptor



PINARD FETOSCOPE

PRODUCT NUMBER

06-1101

FEATURES

Aluminum



STETHOSCOPES

OB/GYN

General

NURSES STETHOSCOPES

SINGLE SIDED

- Economic solution for situations that require a quality listening device.
- Latex Free.

PRODUCT NUMBER	FEATURES	
06-1004	30"	Blue
06-1635	30"	Black
06-1638	30"	Gray



TEACHING TRAINING STETHOSCOPE

SINGLE SIDED

- Specifically designed for multi-function teaching applications.
- Elongated tubing for comfortable two-person use.
- Latex Free.

PRODUCT NUMBER	FEATURES	
06-1668	40"	Black



DISPOSABLE STETHOSCOPES

SINGLE SIDED

- Designed for single use to prevent cross contamination in infectious situations.
- Plastic, ultra sensitive chest piece with remarkable acoustic performance.
- Latex Free.

PRODUCT NUMBER	FEATURES	
96-1201	30"	Blue
96-1202	30"	Yellow



OB/GYN

STETHOSCOPES

Pediatric and Infant

PEDIATRIC AND INFANT STETHOSCOPES

- Traditional bell/diaphragm chest pieces are specifically designed, sized and acoustically precise for children and infants.
- Solid Stainless Steel chest piece features a 1" Traditional Bell for Pediatric and ¾" Bell for Infant.
- Easy rotation between diaphragm or open bell side of stethoscope and sensitive to both high and low frequencies.
- Available as a dual-sided design, in a variety of attractive colors for use with children and infants.
- Excellent acoustic seal and comfortable fit with patented soft ear tips and a comfortably angled, anatomically correct headset.
- Latex Free.



PEDIATRIC STETHOSCOPES

DUAL HEAD

PRODUCT NUMBER	FEATURES
06-1640	30" Gray
06-2415	30" Pink
06-2417	30" Black
06-2418	30" Burgundy
06-2422	30" Blue



INFANT STETHOSCOPES

DUAL HEAD

PRODUCT NUMBER	FEATURES
06-1642	30" Gray
06-2426	30" Pink
06-2427	30" Light Blue
06-2428	30" Black

FETAL STETHOSCOPE

CHROME-PLATED INNER SPRING BINAURAL AND VINYL COATED HEADREST

PRODUCT NUMBER	FEATURES
06-1662	22" Y-Tubing Black



STETHOSCOPES

OB/GYN

Adult

ADULT STETHOSCOPES

- Unsurpassed quality and value.
- Available in single-sided or dual-sided design, in a variety of color choices, with very high acoustic response.
- Specifically designed, sized and acoustically precisioned for adults.
- Patient friendly non-chill rim and diaphragm, excellent acoustic seal and comfortable fit with patented soft ear tips and a comfortably angled, anatomically correct headset.
- Unlimited flexibility and portable convenience.
- Latex Free.

SINGLE SIDED

PRODUCT NUMBER	FEATURES
06-2760	30" Black



DUAL HEAD

Deluxe, stainless steel, two-sided traditional combination chest piece with high and low frequency capabilities.

PRODUCT NUMBER	FEATURES
06-1601	30" Gray
06-2298	30" Purple
06-2299	30" Hunter Green
06-2300	30" Black
06-2302	30" Blue



LIGHT ADULT STETHOSCOPES

DUAL HEAD

Lighter construction reduces fatigue while maintaining performance.

PRODUCT NUMBER	FEATURES
06-1630	30" Black
06-1632	30" Gray



OB/GYN

STETHOSCOPES

Cardiology

CARDIOLOGY STETHOSCOPES

- Superb craftsmanship, construction and materials to offer the ultimate performance and sound transmission.
- Available in single-sided or dual-sided design with very high acoustic response.
- Offered as a double-lumen tube design in a variety of color choices to meet individual preferences.
- Comfortably angled, anatomically correct headset with patented soft ear tips and excellent acoustic seal.
- Latex Free.


SINGLE SIDED

PRODUCT NUMBER	FEATURES	
06-2750	27"	Black
06-2800	28"	Black
06-2802	28"	Blue
06-2803	28"	Green

DUAL HEAD

Deluxe, two-sided traditional combination chest piece with patented turnable diaphragm to allow for high, medium and low frequency sounds to be heard.

PRODUCT NUMBER	FEATURES	
06-1670	28"	Black
06-1671	28"	Burgundy
06-1672	28"	Navy Blue
06-1673	28"	Green



STETHOSCOPES

OB/GYN

Sprague Rappaport

Designed to meet the exacting standards of health care professionals who demand precise acoustic response in various settings, Sklar offers a full range of Stethoscopes that feature superb craftsmanship and performance.

SPRAGUE RAPPAPORT STETHOSCOPES

- The most popular stethoscope style on the market, worldwide.
- Includes five easily interchangeable chest pieces: 3 Bells (Adult, Medium and Infant) and 2 Diaphragms (Small and Large) allowing for a full range of diagnostic capabilities. Also includes 3 different size ear tips.
- Superior high and low frequency sensitivity; designed to detect the entire range of heart and chest sounds. Clinically proven to deliver more faint heart sounds and murmurs than other stethoscopes.
- Chest piece is equipped with a full rotational dual-output acoustic valve mechanism that ensures smooth rotation and accurate acoustic transmission.
- All components are individually machined and precisely assembled to deliver the choice stethoscope of medical professionals.
- Comfortably angled, anatomically correct headset with patented soft ear tips and excellent acoustic seal.
- Latex Free.

SPRAGUE RAPPAPORT STETHOSCOPES

DUAL HEAD

PRODUCT NUMBER	FEATURES
----------------	----------

06-1625	28" Gray
---------	----------


06-1627	28" Black
---------	-----------

06-2304	28" Blue
---------	----------



OB/GYN

SUCTION TUBES



YANKAUER SUCTION TUBE

07-1801	11"	Chrome Plated
07-1805	11"	Stainless Steel
22-6585	8½"	Stainless Steel, Pediatric
07-1802	Tip	For use with 07-1801
07-1806	Tip	For use with 07-1805
22-6586	Tip	For use with 22-6585



POOLE SUCTION TUBE

22-6577	16 French	Curved	Pediatric
34-2023	23 French	Curved	
34-2230	30 French	Straight	Slip on End




BARON SUCTION TUBE

50-1953	3 French
50-1955	5 French
50-1957	7 French



ADSON SUCTION TUBE

50-1982	6"	12 French
50-1985	7½"	15 French



FRAZIER SUCTION TUBE

50-2006	6 French
50-2007	7 French
50-2008	8 French
50-2009	9 French
50-2010	10 French
50-2011	11 French
50-2012	12 French
50-2014	14 French



FRAZIER SUCTION TUBE 90 DEGREE

50-2206	7½"	6 French
50-2207	7½"	7 French
50-2208	7½"	8 French
50-2210	7½"	10 French
50-2212	7½"	12 French



BUIE SUCTION TUBE

80-2715	15 French	With Valve
80-2815	15 French	Without Valve

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

ECONO™ STERILE


OB/GYN

Wireforms




STITCH SCISSORS

96-8971 4½" Littauer 50/cs
96-8973 4½" Sharp/Blunt 50/cs
96-8972 5¼" Sharp/Blunt 50/cs
96-8975 5¼" Sharp/Sharp 50/cs



HALSTED FORCEPS

96-8964 5" Straight 50/cs
96-8965 5" Curved 50/cs




KELLY FORCEPS

96-8967 5½" Straight 50/cs
96-8966 5½" Curved 50/cs



SPONGE FORCEPS

96-8969 7" Serrated 50/cs
96-8970 9½" Serrated 50/cs



TOWEL FORCEP

96-8962 5¼" 50/cs



THUMB FORCEP

96-8963 4½" 1x2 Teeth 50/cs



TISSUE FORCEP

96-8976 5" 50/cs




NEEDLE HOLDERS

96-8968 5¼" Smooth 50/cs
96-8974 5¼" Serrated 50/cs

OB/GYN


ECONO™ STERILE

Needle Holders




BAUMGARTNER NEEDLE HOLDERS

96-2581 5½" Serrated 25/cs
96-2580 5½" Serrated 50/cs




MAYO-HEGAR NEEDLE HOLDERS

96-2599 5" Serrated 25/cs
96-2595 6" Serrated 25/cs
96-2594 6" Serrated 50/cs
96-2710 7" Serrated 50/cs
96-2712 8" Serrated 50/cs




CRILE-WOOD NEEDLE HOLDERS

96-2591 6" Serrated 25/cs
96-2590 6" Serrated 50/cs




WEBSTER NEEDLE HOLDERS

96-2587 5¼" Smooth 25/cs
96-2586 5¼" Smooth 50/cs
96-2659 5¼" Smooth 100/cs



DERF NEEDLE HOLDERS

96-2583 4¾" Serrated 25/cs
96-2582 4¾" Serrated 50/cs




HALSEY NEEDLE HOLDERS

96-2585 5" Smooth 25/cs
96-2584 5" Smooth 50/cs


ECONO™ STERILE Clamps

OB/GYN




BACKHAUS TOWEL CLAMPS

96-2632	3 1/2"	25/cs
96-2636	5 1/4"	25/cs
96-2630	3 1/2"	50/cs
96-2634	5 1/4"	50/cs




LORNA EDNA TOWEL CLAMPS

96-2529	3 1/2"	25/cs
96-2528	3 1/2"	50/cs




HOFF TOWEL CLAMPS

96-2565	5 1/4"	25/cs
96-2564	5 1/4"	50/cs



SPONGE FORCEPS

96-2533	9 1/2"	Straight, Serrated	25/cs
96-2532	9 1/2"	Straight, Serrated	50/cs



PRESBYTERIAN TUBE OCCLUDING FORCEPS

96-2408	7"	Smooth	25/cs
96-2406	7"	Smooth	50/cs


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN


ECONO™ STERILE

Forceps




ADSON BROWN FORCEPS

96-2559	4 3/4"	Side Grasping	25/cs
96-2558	4 3/4"	Side Grasping	50/cs




DRESSING FORCEPS

96-2569	5"	Serrated	25/cs
96-2921	8"	Serrated	25/cs
96-2688	4 1/2"	Serrated	50/cs
96-2568	5"	Serrated	50/cs
96-2660	5"	Serrated	100/cs



TISSUE FORCEPS

96-2579	5" 1x2 Teeth	25/cs
96-2578	5" 1x2 Teeth	50/cs
96-2661	5" 1x2 Teeth	100/cs



FINE POINT SPLINTER FORCEPS

96-2412	4 1/2"	25/cs
96-2410	4 1/2"	50/cs

SKLAR-STAT™

INSTANT FOAM HAND SANITIZER

- Alcohol-Free
- Non-Sticky
- Moisturizes and Conditions Hands

PRODUCT NUMBER

10-2850
50ml Pump Dispensers
24/case

10-2852
550ml Pump Dispensers
12/case


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

ECONO™ STERILE


OB/GYN

Forceps



HARTMANN FORCEPS

96-2556	3½" Straight	25/cs
96-2557	3½" Curved	25/cs
96-2554	3½" Straight	50/cs
96-2555	3½" Curved	50/cs




ROCHESTER PEAN FORCEPS

96-2545	6¼" Straight	25/cs
96-2547	6¼" Curved	25/cs
96-3006	7¼" Straight	25/cs
96-3008	7¼" Curved	25/cs
96-2548	8½" Curved	25/cs
96-2544	6¼" Straight	50/cs
96-2546	6¼" Curved	50/cs
96-2549	8½" Curved	50/cs
96-2664	8½" Curved	100/cs




ROCHESTER-OSHSNER FORCEPS

96-2566	6¼" Straight	25/cs
96-2567	6¼" Curved	25/cs
96-2707	6¼" Straight	50/cs
96-2708	6¼" Curved	50/cs




ALLIS TISSUE FORCEPS

96-2575	6"	4x5 Teeth	25/cs
96-2576	6"	4x5 Teeth	50/cs
96-2774	7¼"	6x7 Teeth	50/cs



BABCOCK FORCEPS

96-2706	6¼" Straight	25/cs
96-2705	6¼" Straight	50/cs



ADSON FORCEPS

96-2571	4¾"	Serrated	25/cs
96-2573	4¾"	1x2 Teeth	25/cs
96-2570	4¾"	Serrated	50/cs
96-2572	4¾"	1x2 Teeth	50/cs
96-2577	4¾"	Serrated Delicate	50/cs
96-2775	4¾"	1x2 Teeth Delicate	50/cs
96-2655	4¾"	Serrated	100/cs
96-2656	4¾"	1x2 Teeth	100/cs

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

ECONO™ STERILE

Scissors




UTILITY SCISSORS

96-2701	5½"	Black	25/cs
96-2703	7"	Black	25/cs
96-1220	7"	Yellow	25/cs
96-1221	7"	Orange	25/cs
96-1222	7"	Purple	25/cs


96-1223	7"	Burgundy	25/cs
96-1224	7"	Magenta	25/cs
96-1225	7"	Navy Blue	25/cs
96-1226	7"	Royal	25/cs
96-1227	7"	Teal Blue	25/cs

96-1228	7"	Lime Green	25/cs
96-1229	7"	Hunter Green	25/cs
96-2702	5½"	Black	50/cs
96-2704	7"	Black	50/cs




HALSTEAD MOSQUITO FORCEPS

96-2537	5"	Straight	25/cs
96-2539	5"	Curved	25/cs
96-2536	5"	Straight	50/cs
96-2538	5"	Curved	50/cs
96-2657	5"	Straight	100/cs
96-2658	5"	Curved	100/cs



CRILE FORCEPS

96-2551	5½"	Straight	25/cs
96-2553	5½"	Curved	25/cs
96-2550	5½"	Straight	50/cs
96-2552	5½"	Curved	50/cs




KELLY FORCEPS

96-2561	5½"	Straight	25/cs
96-2563	5½"	Curved	25/cs
96-2560	5½"	Straight	50/cs
96-2562	5½"	Curved	50/cs
96-2652	5½"	Straight	100/cs
96-2653	5½"	Curved	100/cs

ECONO™ STERILE


OB/GYN

Scissors




STRABISMUS SCISSORS

96-2542	4½" Straight	25/cs
96-2543	4½" Curved	25/cs
96-2540	4½" Straight	50/cs
96-2541	4½" Curved	50/cs




SPENCER STITCH SCISSORS

96-2598	4½"	25/cs
96-2672	3½"	50/cs
96-2596	4½"	50/cs




LITTAUER STITCH SCISSORS

96-2511	3½"	25/cs
96-2513	5½"	25/cs
96-2510	3½"	50/cs
96-2512	5½"	50/cs
96-2663	3½"	100/cs
96-2662	5½"	100/cs




WIRE CUTTING SCISSORS

96-2519	4¾" Angled	25/cs
96-2518	4¾" Angled	50/cs



LISTER BANDAGE SCISSORS

96-2503	5½"	25/cs
96-2501	4½"	50/cs
96-2502	5½"	50/cs
96-2654	5½"	100/cs



KNOWLES BANDAGE SCISSORS

96-2404	5½"	25/cs
96-2402	5½"	50/cs


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

ECONO™ STERILE

Scissors



OPERATING SCISSORS

STRAIGHT


96-2597	4½"	Sharp/Blunt	25/cs	96-2514	5½"	Blunt/Blunt	50/cs
96-2592	4½"	Sharp/Sharp	25/cs	96-2520	5½"	Sharp/Blunt	50/cs
96-2521	5½"	Sharp/Blunt	25/cs	96-2522	5½"	Sharp/Sharp	50/cs
96-2523	5½"	Sharp/Sharp	25/cs	96-2726	6½"	Sharp/Blunt	50/cs
96-2728	4½"	Sharp/Blunt	50/cs	96-2649	5½"	Sharp/Blunt	100/cs
96-2593	4½"	Sharp/Sharp	50/cs				



OPERATING SCISSORS


CURVED

96-2526	5½"	Sharp/Blunt	25/cs
96-2527	5½"	Sharp/Sharp	25/cs
96-2524	5½"	Sharp/Blunt	50/cs
96-2525	5½"	Sharp/Sharp	50/cs




MAYO SCISSORS

96-2531	5½"	Straight	25/cs
96-2530	5½"	Straight	50/cs



METZENBAUM SCISSORS

96-2517	5½"	Straight	25/cs
96-2720	5½"	Curved	25/cs
96-2682	7"	Curved	25/cs
96-2516	5½"	Straight	50/cs
96-2721	5½"	Curved	50/cs
96-2680	7"	Curved	50/cs



IRIS SCISSORS

96-2505	4½"	Straight	25/cs
96-2507	4½"	Curved	25/cs
96-2497	3½"	Straight	50/cs
96-2499	3½"	Curved	50/cs
96-2504	4½"	Straight	50/cs
96-2506	4½"	Curved	50/cs
96-2650	4½"	Straight	100/cs
96-2651	4½"	Curved	100/cs


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

ECONO™ STERILE

OB/GYN

Scissors



GRAVES SPECULUMS

96-2603	Small	25/cs	96-3000	Small	50/cs
96-2605	Medium	25/cs	96-3001	Medium	50/cs
96-2607	Large	25/cs	96-3002	Large	50/cs




PEDERSON SPECULUMS

96-2613	Small	25/cs	96-2612	Small	50/cs
96-2615	Medium	25/cs	96-2614	Medium	50/cs
96-2617	Large	25/cs	96-2616	Large	50/cs



BOZEMAN FORCEPS

96-2622	10¼" S-Curve	25/cs
96-2620	10¼" S-Curve	50/cs



SCHROEDER TENACULUM FORCEPS

96-2624	9½"	25/cs
96-2626	9½"	50/cs


OB/GYN

ECONO™ STERILE




Econo™ Sterile Instruments

- FACTORY IS ISO 9000 CERTIFIED AND CGMP COMPLIANT
- CHROMIUM CONTENT OF 11.5% TO 12.5% TO PREVENT RUSTING
- COMBINES HIGH QUALITY WITH LOW COST
- SAVES LABOR COSTS
- CONVENIENT & SAFE
- ELIMINATES COSTLY CLEANING, PROCESSING AND STERILIZATION
- MANUFACTURED IN PAKISTAN IN A FAMILY-OWNED FACTORY
- PACKAGED AND ETO STERILIZED IN OUR OWN US FACILITY
- STERILE ONE-TIME USE
- 101 PATTERNS



MAYO-HEGAR NEEDLE HOLDERS



97-842	5 1/4"
97-847	6"
97-853	7"
97-869	8"



WEBSTER NEEDLE HOLDER

97-809	5"	Smooth
--------	----	--------

SKLAR CERVIX BRUSH CERVICAL CELL SAMPLER



SKLAR CERVIX BRUSH

NON-STERILE

PRODUCT NUMBER	QUANTITIES
96-4140NS	100/bx

STERILE

PRODUCT NUMBER	QUANTITIES
96-4140	100/bx


Sklar's Cervix Brush is made of polyethylene and has a total length of 20 centimeters. The top portion uses a soft flexible brush to obtain cell samples. The brush consists of 57 semi-circular plastic bristles of various lengths. The shape is such that the top edges of the brush follow the contours of the cervix. The longer middle bristles reach deep into the endocervical canal. The shorter bristles touch both the ectocervical area and transformation zone at the same time.

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

MERIT™




HALSEY NEEDLE HOLDER

97-810 5" Smooth




OLSEN-HEGAR NEEDLE HOLDERS

97-872 5½"
97-877 6½"
97-881 7½"




CRILE-WOOD NEEDLE HOLDER

97-838 6"




ABBY NEEDLE HOLDER

97-894 5"



DERF NEEDLE HOLDER

97-826 4¾"



BAUMGARTNER NEEDLE HOLDER


97-832 5½"

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

MERIT™

OB/GYN




IRIS DRESSING/ TISSUE FORCEPS

DRESSING-SERRATED

98-123	4"	Straight
98-129	4"	Half Curved
98-135	4"	Full Curved


TISSUE-1X2 TEETH

98-142	4"	Straight
97-146	4"	Half Curved
97-151	4"	Full Curved




PLAIN SPLINTER FORCEPS

97-777	3½"
97-781	4½"




BACKHAUS TOWEL CLAMPS

97-606	3½"	1x2 Teeth
97-611	5¼"	1x2 Teeth
97-623	5¼"	Roeder With Ball Stops




VESSELS NON-PERFORATING TOWEL FORCEPS

97-645	4½"
97-655	5½"



LORNA (EDNA) TOWEL CLAMPS

97-629	3½"
97-635	5¼"



FORESTER SPONGE FORCEPS

SERRATED


97-571	7"	Straight
97-577	9½"	Straight
97-586	9½"	Curved

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)


OB/GYN

MERIT™




MAGILL CATHETER FORCEPS

97-285 Child
97-289 Adult




DRESSING FORCEPS

97-705 4½"
97-710 5"
97-727 5½"
97-731 6"
97-736 8"
97-742 10"
97-748 12"



TISSUE FORCEPS

97-753 4½" 1x2 Teeth
97-757 5½" 1x2 Teeth
97-762 5½" 1x2 Teeth
97-769 5½" 2x3 Teeth




ADSON DRESSING/ TISSUE FORCEPS

DRESSING-SERRATED

97-771 4¾"
97-772 4¾" Delicate

TISSUE-1X2 TEETH


97-775 4¾"
97-776 4¾" Delicate



ADSON-BROWN TISSUE FORCEP

STRAIGHT

97-778 4¾" Side Grasping



RUSSIAN TISSUE FORCEP


97-604 6"

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)


MERIT™

OB/GYN




CRILE FORCEPS

STRAIGHT		CURVED	
97-447	5½"	97-452	5½"
97-458	6¼"	97-463	6¼"




ROCHESTER-PEAN FORCEPS

STRAIGHT		CURVED	
97-470	5½"	97-471	5½"
97-473	6¼"	97-477	6¼"
97-481	7¼"	97-485	7¼"
97-492	8"	97-498	8"
97-483	9"	97-487	9"
97-494	10¼"	97-496	10¼"




ROCHESTER-OSHSNER FORCEPS

STRAIGHT		CURVED	
97-500	5½"	97-501	5½"
97-504	6¼"	97-512	6¼"
97-521	7¼"	97-525	7¼"
97-513	8"		
97-517	10"		




ROCHESTER-CARMALT FORCEPS

STRAIGHT		CURVED	
97-528	6¼"	97-534	6¼"
97-537	8"	97-542	8"



PRESBYTERIAN TUBE OCCLUDING FORCEP

97-592	6½"
--------	-----



ALLIS TISSUE FORCEPS

97-930	6"	4x5 Teeth
97-934	6"	5x6 Teeth
97-938	7½"	5x6 Teeth

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN


MERIT™



MULTI-CUT UTILITY SCISSORS


SERRATED

97-105	7½"	Black
97-107	7½"	Blue
97-108	7½"	Purple
97-109	7½"	Green



WIRE CUTTING SCISSORS

97-364 4¾"




LISTER BANDAGE SCISSORS

97-222	3½"	
97-226	4½"	
97-231	5½"	
97-237	5½"	With Clip
97-241	7¼"	




LISTER BAND SCISSORS

97-245 8" Large Ring



HALSTEAD MOSQUITO FORCEPS

STRAIGHT		CURVED	
97-426	5"	97-431	5"



KELLY FORCEPS


STRAIGHT		CURVED	
97-436	5½"	97-442	5½"

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

MERIT™


OB/GYN



OPERATING SCISSORS

STRAIGHT


97-263	4½"	Sharp/Blunt
97-272	5½"	Sharp/Blunt
97-281	6½"	Sharp/Blunt
97-266	4½"	Sharp/Sharp
97-275	5½"	Sharp/Sharp
97-284	6½"	Sharp/Sharp
97-269	4½"	Blunt/Blunt
97-278	5½"	Blunt/Blunt
97-288	6½"	Blunt/Blunt



OPERATING SCISSORS

CURVED

97-292	5½"	Sharp/Blunt
97-293	6½"	Sharp/Blunt
97-295	5½"	Sharp/Sharp
97-297	6½"	Sharp/Sharp
97-270	4½"	Blunt/Blunt
97-299	5½"	Blunt/Blunt
97-300	5½"	Sharp/Sharp




METZENBAUM SCISSORS

STRAIGHT

97-335	5¾"
97-339	7"

CURVED

97-336	5¾"	Sharp/Blunt
97-340	7"	Sharp/Blunt




MAYO DISSECTING SCISSORS

STRAIGHT

97-308	5½"
97-321	6¾"

CURVED

97-312	5½"
97-325	6¾"



IRIS SCISSORS

STRAIGHT

98-103	3½"
98-104	4½"

CURVED

98-107	3½"
98-109	4½"



SPENCER STITCH SCISSORS

97-346	3½"
97-352	4¾"

LITTAUR STITCH SCISSORS


97-358	5½"
--------	-----

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

MERIT™



HEGAR UTERINE DILATOR

SET OF 6
98-370




UMBILICAL SCISSOR

98-367 4"



BOZEMAN UTERINE DRESSING FORCEP

98-374 10¼"



SCHROEDER TENACULUM FORCEP

98-386 9½"

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

MERIT™

OB/GYN



GRAVES VAGINAL SPECULUMS

98-303	3" x ¾"	Small
98-312	4" x 1¼"	Medium
98-324	4½" x 1⅜"	Large
98-325	4¾" x 1½"	Extra Large



PEDERSON VAGINAL SPECULUMS

98-336	3" x 1½"	Small
98-364	4" x ⅝"	Slender
98-348	4" x ⅞"	Medium
98-362	4¾" x 1"	Large
98-366	6½" x 1½"	Extra Large



SIMPSON UTERINE SOUND

MALLEABLE WITH NODULE JUST BELOW CURVE

98-382	12½"	Centimeter Graduated, Silver Plated
--------	------	-------------------------------------



SIMS UTERINE SOUND

MALLEABLE

98-384	12½"	Non-Graduated, Silver Plated
--------	------	------------------------------

OB/GYN

MERIT™

MERIT

Quality

Service

Price

The most cost-effective line of quality instruments from Pakistan designed for routine surgical use.

- All instruments are forged from surgical grade stainless steel
- Crafted in FDA registered manufacturing facilities by skilled instrument makers
- Designed to meet OR tolerances for surgical quality
- 10 year guarantee
- Same day shipment on orders received by 3:00 PM
- No charge for drop shipments
- Available and invoiced by your distributor
- Backed by Sklar Instruments for value pricing
- Custom pricing available in minutes for volume orders
- Ask about our "Bid Pricing"

MERIT

206 OB/GYN

www.sklarcorp.com



1.800.221.2166

SURGI-OR™


OB/GYN

Needle Holders




DERF NEEDLE HOLDER

95-826 4 3/4" Serrated




WEBSTER NEEDLE HOLDERS

95-809 5 1/4" Smooth
95-813 5 1/4" Serrated




HALSEY NEEDLE HOLDERS

95-810 5 1/4" Smooth
95-812 5 1/4" Serrated



BAUMGARTNER NEEDLE HOLDER

95-832 5 1/4" Serrated




MAYO-HEGAR NEEDLE HOLDERS

SERRATED
95-842 5 1/4"
95-847 6"
95-853 7"
95-869 8"

OB/GYN


SURGI-OR™

Forceps



FOERSTER SPONGE FORCEPS

95-571	7"	Straight	Serrated
95-577	9½"	Straight	Serrated
95-586	9½"	Curved	Serrated




BACKHAUS TOWEL CLAMPS

95-606	3½"
95-611	5¼"




LORNA EDNA TOWEL CLAMPS

95-629	3½"
95-635	5¼"



**PRESBYTERIAN TUBE
OCCLUDING FORCEP**

95-592	7"	Straight, Smooth
--------	----	------------------



MAGILL CATHETER FORCEPS

95-285	8"	Child
95-289	10"	Adult


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SURGI-OR™

OB/GYN


Forceps



ADSON FORCEPS

DRESSING


95-771 4¾" Serrated
95-772 4¾" Serrated, Delicate



ADSON FORCEPS

TISSUE


95-775 4¾" 1x2 Teeth
95-776 4¾" 1x2 Teeth, Delicate



ADSON FORCEPS

ADSON-BROWN


95-778 4¾" Side-Grasping



DRESSING FORCEPS

SERRATED


95-705 4½"
95-710 5"
95-727 5½"
95-731 6"
95-736 8"
95-742 10"



TISSUE FORCEPS

1X2 TEETH

95-753 4½"
95-757 5"
95-762 5½"
95-766 6"
95-770 8"




SPLINTER FORCEPS

95-777 3½"
95-781 4½"

OB/GYN


SURGI-OR™

Forceps




HALSTED MOSQUITO FORCEPS

95-426	5"	Straight
95-431	5"	Curved




HARTMANN MOSQUITO FORCEPS

95-409	3½"	Straight
95-421	3½"	Curved




ROCHESTER-PEAN FORCEPS

95-470	5½"	Straight
95-471	5½"	Curved
95-473	6¼"	Straight
95-477	6¼"	Curved
95-481	7¼"	Straight
95-485	7¼"	Curved
95-492	8"	Straight
95-498	8"	Curved
95-483	9"	Straight
95-487	9"	Curved
95-494	10¼"	Straight



ROCHESTER-OGHSNER FORCEPS

STRAIGHT, 1x2 TEETH		CURVED, 1x2 TEETH	
95-500	5½"	95-501	5½"
95-504	6¼"	95-512	6¼"
95-521	7¼"	95-525	7¼"
95-513	8"	95-514	8"
95-515	9"	95-516	9"
95-517	10"	95-518	10"




ALLIS TISSUE FORCEPS

95-934	6"	5x6 Teeth
95-938	7½"	5x6 Teeth

SURGI-OR™


OB/GYN

Forceps




BOZEMAN UTERINE DRESSING FORCEP

95-374 10¼" S-Curve



SCHROEDER TENACULUM FORCEP


95-386 9½" Straight



BABCOCK FORCEPS

95-952 6¼"

95-441 8"




CRILE FORCEPS

95-447 5½" Straight

95-452 5½" Curved

95-458 6¼" Straight


95-463 6¼" Curved



KELLY FORCEPS

95-436 5½" Straight

95-442 5½" Curved



RANKIN-KELLY FORCEPS

95-465 6¼" Straight

95-469 6¼" Curved

OB/GYN


SURGI-OR™

Scissors



STRABISMUS SCISSORS

95-360 4½" Straight
95-361 4½" Curved




STITCH SCISSORS

95-346 3½" Spencer
95-352 4½" Spencer
95-358 5½" Littauer



LISTER BANDAGE SCISSORS


95-222 3½"
95-226 4½"
95-231 5½"
95-241 7¼"
95-243 8"
95-245 8" With Large Finger Ring



LISTER BANDAGE SCISSORS

WITH POCKET CLIP

95-232 5½"
95-232B 5½" Blue Ring Handle
95-232P 5½" Pink Ring Handle



KNOWLES BANDAGE SCISSOR


95-254 5½"



MULTI-CUT UTILITY SCISSORS

BLACK HANDLE

95-100 6" Serrated
95-105 7½" Serrated



WIRE CUTTING SCISSOR

ANGLED

95-364 4¾" Serrated

SURGI-OR™

OB/GYN

Scissors



OPERATING SCISSORS

STRAIGHT


95-266	4½"	Sharp/Sharp
95-265	4½"	Sharp/Blunt
95-290	5"	Sharp/Sharp
95-275	5½"	Sharp/Sharp
95-272	5½"	Sharp/Blunt
95-278	5½"	Blunt/Blunt
95-284	6½"	Sharp/Sharp
95-281	6½"	Sharp/Blunt



OPERATING SCISSORS


CURVED

95-294	4½"	Curved	Sharp/Sharp
95-295	5½"	Curved	Sharp/Sharp
95-292	5½"	Curved	Sharp/Blunt
95-299	5½"	Curved	Blunt/Blunt




IRIS SCISSORS

95-103	3½"	Straight
95-107	3½"	Curved
95-104	4½"	Straight
95-109	4½"	Curved



MAYO DISSECTING SCISSORS

95-308	5½"	Straight	
95-330	5½"	Curved	
95-321	6¾"	Straight	
95-325	6¾"	Curved	
95-329	6¾"	Straight	Serrated
95-332	6¾"	Curved	Serrated
95-327	9"	Straight	
95-328	9"	Curved	



METZENBAUM DISSECTING SCISSORS

95-335	5¾"	Straight
95-336	5¾"	Curved
95-339	7"	Straight
95-340	7"	Curved



STEVENS DISSECTING SCISSORS

95-160	4⅞"	Straight	Blunt/Blunt
95-161	4⅞"	Curved	Blunt/Blunt

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

SURGI-OR™

Tungsten Carbide




TC HALSEY NEEDLE HOLDERS

95-132 5¼" Smooth
95-133 5¼" Serrated



TC OLSEN-HEGAR NEEDLE HOLDERS

95-146 4¾" Serrated, Delicate
95-147 5½" Serrated
95-148 6½" Serrated
95-149 7½" Serrated



TC BAUMGARTNER NEEDLE HOLDER


95-135 5¼" Serrated



TC MAYO-HEGAR NEEDLE HOLDERS

SERRATED

95-140 5¼"
95-141 6"
95-142 7"
95-143 8"



TC CRILE-WOOD NEEDLE HOLDERS

SERRATED

95-137 6"
95-138 7"
95-139 8"



TC CRILE-MURRAY NEEDLE HOLDER

95-134 6" Serrated

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SURGI-OR™


OB/GYN

Tungsten Carbide




TC MAYO DISSECTING SCISSORS

95-112	6 3/4"	Straight
95-113	6 3/4"	Curved




TC MAYO METZENBAUM SCISSORS

95-119	5 1/4"	Straight	Standard
95-122	5 3/4"	Curved	Standard
95-115	7"	Straight	Standard
95-117	7"	Curved	Standard
95-118	5 3/4"	Straight	Delicate
95-120	5 3/4"	Curved	Delicate
95-114	7"	Straight	Delicate
95-116	7"	Curved	Delicate



TC STRABISMUS SCISSORS

95-126	4 1/2"	Straight	Blunt/Blunt
95-127	4 1/2"	Curved	Blunt/Blunt




TC WIRE CUTTING SCISSOR

95-128	4 3/4"	Angled	Serrated
--------	--------	--------	----------



TC DERF NEEDLE HOLDER

95-136	4 3/4"	Serrated
--------	--------	----------



TC WEBSTER NEEDLE HOLDER

95-130	5 1/4"	Smooth
--------	--------	--------

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

SURGI-OR™

Surgi-OR™ is a full line of high-quality, Mid-Grade OR instruments, hand-crafted and forged from Japanese surgical grade stainless steel. This cost-effective line of instruments are designed to give you maximum performance at affordable cost and offers significant value to the Operating Room as well as other important surgical settings.

Sklar is pleased to offer Surgi-OR™ Tungsten Carbide Instruments that feature this extremely hard alloy resulting in a superior instrument.

Surgi-OR™ Needle Holders feature solid Tungsten Carbide insert jaws for exceptional control and unbeatable strength. These perfectly balanced instruments are readily identified by the dual gold finger rings and afford quality and performance you can trust in a variety of surgical settings.



GRAVES VAGINAL SPECULUMS

95-304	3" x 1/2"	X-small
95-303	3" x 3/4"	Small
95-312	4" x 1 1/4"	Medium
95-324	4 1/2" x 1 3/8"	Large
95-322	4 3/4" x 1 1/2"	X-Large



PEDERSON VAGINAL SPECULUMS

95-334	2 3/4" x 1/2"	X-Small
95-347	3" x 1/2"	Small
95-348	4" x 7/8"	Medium
95-362	4 3/4" x 1"	Large
95-368	4 3/4" x 1 5/8"	X-Large



TC OPERATING SCISSORS

95-110	5 1/2"	Straight	Sharp/Blunt
95-111	5 1/2"	Curved	Sharp/Blunt



TC IRIS SCISSORS

95-124	4 1/2"	Straight	Sharp/Sharp
95-125	4 1/2"	Curved	Sharp/Sharp

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PROCEDURE TRAYS

OB/GYN

Non-Sterile

96-2312 — ECONOMY BABY CARE PACK 50/CASE

DESCRIPTION	QUANTITIES
Aloetouch Wipes, 12/cs	1
Baby Comb	1
Draw String Bag	1
Ivory Soap, 3 oz.	1
Syringe, Bulb Irrigation, 60cc	1
Tape Measure	1



96-2313 — BABY CARE PACK 50/CASE

DESCRIPTION	QUANTITIES
Aloetouch Wipes, 12/cs	1
Baby Comb	1
Digital Thermometer	1
Draw String Bag	1
Ivory Soap, 3 oz.	1
Lubricating Jelly	1
Syringe, Bulb Irrigation, 60cc	1
Tape Measure	1



96-2314 — DELUXE BABY CARE PACK 20/CASE

DESCRIPTION	QUANTITIES
Aloetouch Wipes, 12/cs	1
Baby Comb	1
Diapers	15
Digital Thermometer	1
Draw String Bag	1
Hand/Body Lotion, 2 oz.	1
Ivory Soap, 3 oz.	1
Lubricating Jelly	1
Syringe, Bulb Irrigation, 60cc	1
Tape Measure	1



OB/GYN

PROCEDURE TRAYS

Non-Sterile

96-2309 — ECONOMY MATERNITY PACK 50/CASE

DESCRIPTION	QUANTITIES
Draw String Bag	1
OB Pads	5
Ziplock Bag, 10" x 12"	1



96-2310 — MATERNITY PACK 50/CASE

DESCRIPTION	QUANTITIES
Draw String Bag	1
Peri Cold Packs	3
Stretch Briefs	1



96-2311 — DELUXE MATERNITY PACK 20/CASE

DESCRIPTION	QUANTITIES
Draw String Bag	1
OB Pads	10
OB Pads Maxi	9
Peri Bottle	1
Peri Cold Packs	3
Sani-Bag	12
Stretch Briefs	2
Underpad, 23" x 24"	5



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PROCEDURE TRAYS

OB/GYN

Sterile

96-4306— C-SECTION PACK

1/CASE

DESCRIPTION	QUANTITIES
Baby Receiving Blanket	1
Back Table Cover	1
Bowl, 80 oz.	1
Sponge Bowl, 32 oz.	1
#10 Carbon Scalpel Blade	2
#20 Carbon Scalpel Blade	1
Drape, 24" x 36" Fenestrated	1
Drape, C-Section w/Pouch & Fenestration	1
Gown, XL	1
Gown, XL Reinforced	2
Lap Sponges	15
Mayo Stand Cover	1
Medicine Cup, 2 oz.	1
OB Pads	2
OR Basin, 7 qt.	2
OR Towel, Blue	6
Suction Tubing, 1/4" x 12'	1
Suture Bag	1
Syringe, Ear/Ulcer, 2 oz.	1
Umbilical Clip	1
Yankauer Suction Tube, Non-Vented	1



96-4307— C-SECTION TRAY

2/CASE

DESCRIPTION	QUANTITIES
Baby Cap	1
Baby Receiving Blanket	1
Back Table Cover	1
#10 Carbon Scalpel Blade	2
#20 Carbon Scalpel Blade	1
Drape, C-Section w/Pouch & Fenestration	1
Emesis Basin, 700cc	1
Gauze, 4" x 4"	10
Gauze, 4" x 4" X-Ray	10
Gloves, Pair, Devon Light LF	2
Gown, XL	2
Gown, XL Reinforced	1
Lap Sponges	20
Mayo Stand Cover	1
Medicine Cup, 2 oz.	1
Needle, 18GA x 1 1/2"	1
Needle Counter	1
Non Adherent Dressing, 3" x 4"	1
OB Pads	2
OR Basin, 7 qt.	2
OR Towel, Blue	5
Pitcher, 1200ml	1
Placenta Basin	1
Skin Marker	1
Specimen Cup w/ Lid, 4 oz.	1
Sponge Bowl, 32 oz.	1
Suction Tubing 1/4" x 12'	1
Suture Bag	1
Syringe, 10cc, L/L	1
Syringe, 20cc, L/L	1
Syringe, Ear/Ulcer, 2 oz.	1
Towel Absorbent	2
Tray, Large Foam	1
Umbilical Clip	1
Utility Drape w/Tape	4
Venous Blood Collection Tube, 7ml	1
Venous Blood Collection Tube, 10ml	1
Yankauer Suction Tube, Non-Vented	1



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

PROCEDURE TRAYS

Sterile

96-3171 — LABOR AND DELIVERY TRAY 5/CASE

DESCRIPTION	QUANTITIES
Back Table Cover	1
Bowl, 80 oz.	1
Bowl, Lid	1
Drape, Vaginal Delivery	1
Gauze, 4" x 4" X-Ray	1
Half Drape Sheet	1
Lap Sponges, Pk/5	1
Maternity Pad	2
Needle, 21GA x 1½"	1
OR Towel, Blue	4
Specimen Tube 8ml w/Screw Cap	1
Syringe, Bulb 3 oz.	1
Syringe, 10cc, Control, L/L	1
Umbilical Clip	1



96-4305 — L & D PACK 3/CASE

DESCRIPTION	QUANTITIES
Abdominal Drape	1
Baby Receiving Blanket	1
Back Table Cover	1
Bowl, 16 oz.	1
Bowl, 80 oz.	1
Sponge Bowl, 32 oz.	1
Gauze, 4" x 4" X-Ray	30
Gown, XL	1
Gown, XL Microcool	1
Half Drape Sheet	1
Leggings, 31" x 48"	2
Needle, 22GA x 1½"	1
OB Pads	2
OR Basin, 7 qt.	1
OR Towel, Blue	4
Syringe, 10cc, L/L	1
Syringe, Ear/Ulcer 2 oz.	1
Towel Absorbent	2
Tray, Half Deep	1
Umbilical Clip	1
Under Buttocks Drape w/Pouch	1
Vaginal Packing Gauze, 4" x 36"	1



96-3181 — C-SECTION TRAY 5/CASE

DESCRIPTION	QUANTITIES
Abdominal Pad	1
6" Applicator Prep Stick	3
Back Table Cover	1
Band Bag	1
Drape, Vaginal Delivery	1
Gloves, Pair/Walleted, Latex	1
Gown, Large	1
Lap Sponges	10
Maternity Pad	1
Mayo Stand Cover	1
Needle Counter	1
OR Towel, Blue	4
PVP Prep Bottle, 2 oz.	1
Specimen Tube 8ml w/Screw Cap	1
Suction Tubing, ¼" x 12'	1
Syringe, Bulb, 3 oz.	1
Tray, Half Shallow	1
Tray, Large Foam	1
Tray, Prep Shallow	1



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PROCEDURE TRAYS

OB/GYN

Sterile

96-1916 — LLETZ ELECTROSURGERY TRAY 10/CASE

DESCRIPTION	QUANTITIES
Alcohol Prep Pad	1
Ampule Holder, 5ml	1
Lidocaine 2%, 10ml	1
Lidocaine 2%, 30ml w/Epi	1
Lugol's PVP Solution	1
Monse's Hemostatic Paste	1
Needle, 27GA x 3 1/4"	1
Needle, 30GA x 1"	1
Panty Liner	1
8" Rayon Swab	6
Tray, Half Shallow	1
Underpad, 17" x 24"	1
Ziplock Bag, 4" x 7"	1



96-1917 — PAP SMEAR TRAY

50/CASE

DESCRIPTION	QUANTITIES
Alcohol Swabstick Single	1
Ampule Holder, 5ml	1
Cervix Brush	1
6" Cotton Tip Applicator	2
Glass Slides Frosted	1
Gloves, Pair/Walleted LF	1
Graves Vaginal Speculum, Disposable	1



96-4319 — D & C / GYN PACK

4/CASE

DESCRIPTION	QUANTITIES
6" Applicator Prep Stick	3
Back Table Cover	1
Catheter, Rob-Nel, 14 Fr	1
6" Cotton Tip Applicator	2
Gauze, 4" x 4"	10
Gloves, Pair/Walleted, Powder Free	1
Gown, XL	1
Lap Drape	1
Leggings, Lithotomy, 31" x 43"	2
Lubricating Jelly, Foil Pack 5gram	2
OB Pads	1
OR Towel, Blue	7
PVP Paint Bottle, 4 oz.	1
PVP Scrub Bottle, 4 oz.	1
Specimen Cup w/Lid, 4 oz.	1
Sponge, Wing Small	6
Telfa Pad, 8" x 3"	1
Towel Absorbent	2
Tray, 2 Compartment	1
Tray, 3 Compartment	1
Under Buttocks Drape w/Pouch	1



OB/GYN

PROCEDURE TRAYS

Sterile

96-1801 — ULTRASOUND TRAY

10/CASE

DESCRIPTION	QUANTITIES
Adhesive Bandage	1
6" Applicator Prep Stick	3
Drape, 18" x 26"	3
Drape, 24" x 36" Fenestrated	1
Gauze, 3" x 3"	12
Lidocaine 1%, 5ml Ampule	1
Needle, 20GA x 1½"	1
Needle, 25GA x 1½"	1
6" Paper Ruler	1
Procedure Tray	1
Syringe, 10cc, L/L	3
Tray, Prep Shallow	1



96-1803 — CYST / FNA / CORE TRAY

10/CASE

DESCRIPTION	QUANTITIES
Alcohol Swabsticks Triple	1
Ampule Holder, 5ml	1
Bowl, 16 oz.	1
Drape, 24" x 36" Fenestrated	1
Gauze, 4" x 4"	6
Lidocaine 1%, 5ml Ampule	1
Needle, 18GA x 1½"	1
Needle, 21GA x 1½"	1
Needle, 22GA x 1½"	1
Needle, 25GA x ½"	1
Probe Cover	1
PVP Swabsticks 3's	1
Scalpel, #11, Disposable	1
Syringe, 5cc, L/L	1
Syringe, 10cc, L/L	1
Tray, Half Deep	1
Ziplock Bag, 8" x 10"	1



96-1807 — CYST ASPIRATION TRAY

25/CASE

DESCRIPTION	QUANTITIES
Gauze, 4" x 4"	5
Needle, 21GA x 1½"	1
PVP Swabsticks 3's	1
Syringe, 10cc, L/L	1
Tray, Small Deep Prep	1



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PROCEDURE TRAYS

OB/GYN

Sterile

96-1721 — HYSTEOSALPINGOGRAM BASIC TRAY 10/CASE

DESCRIPTION	QUANTITIES
Drape, 18" x 26" Fenestrated	1
20" Extension Set, L/L	1
Gauze, 4" x 4"	6
Gloves, Pair/Walleted, Powder Free	1
Graves Vaginal Speculum, Disposable	1
Half Drape Sheet	1
Lubricating Jelly	1
Needle, 18GA x 1½"	1
PVP Ointment Packet, 1 gram	1
PVP Solution .75 oz	1
8" Rayon Swab	2
Stopcock, 1-Way LP w/Rot	1
Syringe, 10cc, L/L	1
Syringe, 10cc, Control, L/L	1
Tray, Half Deep	1
Tray, Large Foam	1
Uterine Injector, 2.0	1
Uterine Sound	1
Ziplock Bag, 6" x 6"	1



96-1722 — AMNIOCENTESIS TRAY 10/CASE

DESCRIPTION	QUANTITIES
Adhesive Bandage	1
Drape, 24" x 36" Fenestrated	1
Filter Straw, 5 Micron	1
Gauze, 4" x 4"	8
Lidocaine 1%, 5ml Ampule	1
Needle, 22GA x 1½"	1
Needle, 22GA x 3½" Spinal	1
Needle, 25GA x ½"	1
6" Plastic Ruler	1
PVP Swabsticks 3's	1
Specimen Tube, Clear	3
Syringe, 5cc, L/L	1
Syringe, 10cc, L/L	1
Syringe, 20cc, L/L	1
Towel Absorbent	1
Ultrasound Gel Pak	1



96-1723 — STANDARD CT BIOPSY TRAY 10/CASE

DESCRIPTION	QUANTITIES
Adhesive Bandage	1
Alcohol Prep Pad	1
Ampule Holder 5ml	1
Drape, 24" x 36" Fenestrated	1
Gauze, 4" x 4"	10
Lidocaine 1%, 5ml Ampule	2
Mini Scalpel	1
Needle, 22GA x 1½"	1
Needle, 25GA x ½", Safety Glide	1
6" Paper Ruler	1
PVP Swabsticks 3's	1
Syringe, 10cc, L/L	4
Towel, Absorbent	2
Tray, Half Shallow	1



OB/GYN

DISPOSABLES



SKLAR TRU-PUNCH™ DISPOSABLE VULVAR BIOPSY PUNCH

The SKLAR Tru-punch has a seamless steel tip carefully honed for a razor sharp incision and perfect sample. The ribbed and hollow handle permits tissue to be easily dislodged.

STERILE DISPOSABLE BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION	QUANTITY	PRODUCT NUMBER	DESCRIPTION	QUANTITY
96-0192	1mm	25/bx	96-1092	1mm	50/bx
96-0197	1.5mm	25/bx	96-1098	1.5mm	50/bx
96-1140	2mm	25/bx	96-1100	2mm	50/bx
96-0102	2.5mm	25/bx	96-1102	2.5mm	50/bx
96-1142	3mm	25/bx	96-1105	3mm	50/bx
96-1144	3.5mm	25/bx	96-1110	3.5mm	50/bx
96-1146	4mm	25/bx	96-1115	4mm	50/bx
96-1148	5mm	25/bx	96-1120	5mm	50/bx
96-1150	6mm	25/bx	96-1125	6mm	50/bx
96-1152	8mm	25/bx	96-1130	8mm	50/bx
96-1160	Assorted - 5 of each: 2, 3, 3.5, 4 & 5mm	25/bx	96-0160	Assorted - 5 of each: 2, 3, 3.5, 4 & 5mm	50/bx

DISPOSABLES

OB/GYN

UTERINE MANIPULATOR/ INJECTOR, DISPOSABLE

STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-4120	4.5mm O.D. -33cm Length	12/cs



ENDOMETRIAL BIOPSY CANNULA

STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-4100	For Aspiration Biopsy	25/cs



ENDOMETRIAL SAMPLING SET

STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-4160	3mm	10/bx



NEEDLE EXTENDERS, DISPOSABLE

PLASTIC LUER LOCK ADAPTER

STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-5190	3"	10/bx
96-5192	4"	10/bx
96-5194	5"	10/bx
96-5196	6"	10/bx




NEEDLE EXTENDERS, REUSABLE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-4190	3"	1
96-4192	4"	1
96-4194	5"	1
96-4196	6"	1



OB/GYN

DISPOSABLES



SKLAR CERVIX BRUSH CERVICAL CELL SAMPLER

Sklar's Cervix Brush is made of polyethylene and has a total length of 20 centimeters. The top portion uses a soft flexible brush to obtain cell samples. The brush consists of 57 semi-circular plastic bristles of various lengths. The shape is such that the top edges of the brush follow the contours of the cervix. The longer middle bristles reach deep into the endocervical canal. The shorter bristles touch both the ectocervical area and transformation zone at the same time.

SKLAR CERVIX BRUSH

STERILE

PRODUCT NUMBER	QUANTITIES
96-4140	100/bx

NON-STERILE

PRODUCT NUMBER	QUANTITIES
96-4140NS	100/bx

DISPOSABLES

OB/GYN

PESSARIES

Sklar offers a full line of soft, silicone Pessaries for effective vaginal support and relief of incontinence. Pessaries are uniquely designed for ease of insertion and removal. Properly fitted, Pessaries can be successfully used to treat a variety of gynecologic conditions and abnormalities. Choose from the most popular patterns, available in a variety of shapes and sizes. Latex Free.



RING WITH SUPPORT

For Mild Uterine Prolapse or a Cystocele

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
90-5291	#1 2"	1 Each
90-5292	#2 2¼"	1 Each
90-5293	#3 2½"	1 Each
90-5294	#4 2¾"	1 Each
90-5295	#5 3"	1 Each
90-5296	#6 3¼"	1 Each
90-5297	#7 3½"	1 Each



GELHORN PESSARY WITH DRAIN

For a Mild Cystocele or Moderate Uterine Prolapse

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
90-5282	#2 2"	1 Each
90-5283	#3 2¼"	1 Each
90-5284	#4 2½"	1 Each
90-5285	#5 2¾"	1 Each
90-5286	#6 3"	1 Each



DONUT

For Moderate Uterine Prolapse or a Mild Cystocele or Rectocele

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
90-5271	#1 2¼"	1 Each
90-5272	#2 2½"	1 Each
90-5273	#3 2¾"	1 Each
90-5274	#4 3"	1 Each
90-5275	#5 3¼"	1 Each
90-5276	#6 3½"	1 Each



CUBE WITH DRAIN

For Moderate to Severe Uterine Prolapse or a Mild Cystocele or Rectocele

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
90-5301	#1 29mm	1 Each
90-5302	#2 33mm	1 Each
90-5303	#3 37mm	1 Each
90-5304	#4 41mm	1 Each
90-5305	#5 45mm	1 Each

OB/GYN 185

OB/GYN

DISPOSABLES



SKLAR STERILE DISPOSABLE ENDOCERVICAL CURETTE

SINGLE PEEL POUCHES, SINGLE USE ONLY

PRODUCT NUMBER	QUANTITIES
----------------	------------

96-4143	25/bx
---------	-------

Federal law restricts the sale of this device by or on the order of a physician. Contents sterile unless package has been opened or damaged.




SKLAR OS DILATORS, DISPOSABLE

STERILE

PRODUCT NUMBER	QUANTITIES
----------------	------------

96-4155	50/bx
---------	-------



SKLAR OS DILATOR SET

NON-STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
----------------	-------------	------------

96-4150	Teflon Reusable Set Includes Os Locator, Canal Finder, Fundus Dilator	Set of 3
---------	--	----------

DISPOSABLES

OB/GYN

SKLAR SILICONE VAGINAL DILATORS

Sklar Silicone Vaginal Dilators are used to enlarge the vagina for the purposes of examination, a surgical procedure or for the relief of vaginismus. This is accomplished by supporting or distending the vaginal canal. They are contraindicated in acute genital tract or pelvic infections or with non-compliant patients.


VAGINAL DILATOR SETS

NON-STERILE

PRODUCT NUMBER	DESCRIPTION				QUANTITIES
90-5260	Small Family	Length	Small OD	Large OD	Set of 4
	XL	75.0mm	13.5mm	18.0mm	
	L	65.0mm	12mm	16.0mm	
	M	60.0mm	10.5mm	14.0mm	
	S	50.0mm	9.5mm	13.0mm	
90-5262	Medium Family	Length	Small OD	Large OD	Set of 4
	XL	95.0mm	26.0mm	33.0mm	
	L	86.0mm	22.5mm	28.0mm	
	M	78.0mm	19.5mm	23.0mm	
	S	75.0mm	12.5mm	17.0mm	
90-5264	Large Family	Length	Small OD	Large OD	Set of 4
	XL	150.0mm	30.0mm	38.0mm	
	L	138.0mm	25.5mm	32.0mm	
	M	118.0mm	19.5mm	28.0mm	
	S	120.0mm	19.5mm	24.0mm	



Normal Anatomy



Dilator in Position

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

DISPOSABLES



SKLAR MONSEL'S FERRIC SUBSULFATE PASTE

A ONE-TIME USE pre-mixed paste packaged in convenient ready-to-use bottles. Also includes 12 rayon tip applicators.

Sklar has converted Monsel's solution (ferric subsulfate) into an efficacious and conveniently packaged paste. As a hemostatic agent it is ideal to apply to cervical tissue following colposcopic biopsy, electrosurgical excision (LLETZ or LEEP) and ablative (LASER) surgery. Sklar Monsel's is fabricated to relieve office personnel from the inexact and messy task of mixing their own compound for post-surgical treatment of their female patients.

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-1200	8ml Bottle & Swab	12 Each

SKLAR LUGOL'S STRONG IODINE SOLUTION USP

A ONE-TIME USE pre-mixed iodine solution packaged in convenient ready-to-use bottles. Also includes 12 rayon tip applicators.

This iodine contrast medium is also sealed in airtight bottles which assist in mitigating waste associated with the use of this product. The Lugol's works quickly to assist the Colposcopist to identify potentially diseased cervical tissue (CIN) by changing suspect tissue in the cervical transformation zone into a white "target" for excision or biopsy. Healthy tissue will remain unchanged by the application of Lugol's.

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-1210	8ml Bottle & Swab	12 Each

DISPOSABLES

OB/GYN



PLASTIC VAGINAL SPECULUMS, LIGHT SOURCE STYLE

NON-STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-0827	Small, Blue	100/cs
96-0837	Medium, Green	100/cs
96-0847	Large, White	100/cs



PLASTIC VAGINAL SPECULUMS, CUSCO STYLE STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-2642	Small	50/cs
96-2644	Medium	50/cs
96-2648	Large	50/cs

NON-STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-0557	Small	10/pk, 10pk/cs, 100/cs
96-0567	Medium	10/pk, 10pk/cs, 100/cs
96-0577	Large	10/pk, 10pk/cs, 100/cs

OB/GYN

DISPOSABLES



SKLARSPEC™

DISPOSABLE GRAVES VAGINAL SPECULA

SKLARSPEC plastic specula are clear, disposable, and comfortable. Plastic material allows 360 degree visibility and reduces thermal shock. Provides protection of the doctor, the nurse and the patient from cross contamination. Clean up is eliminated by single-use only. Plastic material avoids unpleasant temperature shock associated with metal. One hand operation with wide range adjustment.

PLASTIC VAGINAL SPECULUMS, GRAVES STYLE

STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-2602	Small	50/cs
96-2604	Medium	50/cs
96-2606	Large	50/cs
96-1857	Small	100/cs
96-1867	Medium	100/cs
96-1877	Large	100/cs

NON-STERILE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-0857	Small	100/cs
96-0867	Medium	100/cs
96-0877	Large	100/cs

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

INSTRUMENT SETS

OB/GYN

98-1401 — RADICAL VULVECTOMY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
50-1472	Adson Forceps	Curved, 7¼"	4
50-3047	Adson Tissue Forceps	1x2 Teeth, 4¾"	2
36-2160	Allis Tissue Forceps	4x5 Teeth, 6"	10
47-2955	Backhaus Towel Clamps	5¼"	2
91-1070	Barrett Tenaculum Forceps	Straight, 7"	2
17-3055	Crile Forceps	Straight, 5½"	4
17-3155	Crile Forceps	Curved, 5½"	18
50-1290	Cushing Nerve/Vein Retractors	9"	2
52-5177	DeBakey Atraumatic Forceps	2mm, 7¾"	2
19-1060	Dressing Forceps	Serrated, 6"	2
19-1080	Dressing Forceps	Serrated, 8"	2
10-1810	E-Z Lock Instrument Stringers	10" x 2½"	2
87-2195	Foerster Sponge Forceps	Straight, Serrated, 9½"	4
55-2880	Gemini Clamps	Delicate, Angled, 8"	4
17-1550	Halsted Mosquito Forceps	Curved, 5"	4
60-1661	Kelly Retractors	9½", 2" x 2½"	2
47-2852	Lorna-Edna Towel Clamps	5¼"	2
60-1676	Richardson Retractors	With SklarGrip™ Handle, Small, 9½", ¾" x 1"	2
60-1681	Richardson Retractors	With SklarGrip™ Handle, Medium, 9½", 1" x 1¼"	2
17-2862	Rochester-Ochsner Forceps	Straight, 1x2 Teeth, 6¼"	8
17-2880	Rochester-Ochsner Forceps	Straight, 1x2 Teeth, 8"	10
17-2262	Rochester-Pean Forceps	Curved, 6¼"	16
06-2903	Scalpel Handles	#3	2
06-2907	Scalpel Handle	#7	1
91-1298	Schroeder Tenaculum Forceps	Rounded Tips, 10"	2
16-1510	Sklar Edge Mayo Dissecting Scissor	Tungsten Carbide, Straight, 6¾"	1
16-1610	Sklar Edge Mayo Dissecting Scissor	Tungsten Carbide, Curved, 6¾"	1
16-1905	Sklar Edge Metzenbaum Lahey Scissor	Tungsten Carbide, Delicate, Curved, 5¾"	1
16-2110	Sklar Edge Metzenbaum Scissor	Tungsten Carbide, Curved, 7"	1
10-4096	SklarLite Full Size Container	23" x 11" x 6" Non-Perforated Bottom	1
10-4290	SklarLite Full Size Wire Basket	21¼" x 10" x 4"	1
91-1465	Stauder-Moore Tenaculum Forceps	6½"	2
21-8031	TC Mayo Hegar Needle Holders	Tungsten Carbide, Serrated, 6"	2
21-8032	TC Mayo Hegar Needle Holders	Tungsten Carbide, Serrated, 7"	2
19-1260	Tissue Forceps	1x2 Teeth, 6"	2
19-1280	Tissue Forceps	1x2 Teeth, 8"	2
60-1085	US Army Retractor	Set of 2, Double End, 8¼"	1
42-1865	Weitlaner Retractors	3x4 Teeth, Blunt, 6½"	2

OB/GYN

INSTRUMENT SETS

90-8557 — COLPOSCOPY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
90-8334	Iris Hook	Quick Release, 4mm, 9"	1
90-6612	Kevorkian-Young Endocervical Curette	With Basket, 12"	1
90-3995	Kogan Endocervical Speculum	Wide Wire Blades, With Gauge	1
90-7520	Tischler Oval Biopsy Forcep	3mm x 7mm	1

98-1400 — GYNECOLOGICAL ABDOMINAL SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
36-2275	Allis Tissue Forceps	5x6 Teeth, 7½"	8
91-1070	Barrett Tenaculum Forceps	Straight, 7"	2
10-1808	E-Z Lock Instrument Stringer	8" x 2½"	1
91-2282	Jacobs Vulsellum Forceps	Straight, 8¼"	2
52-2970	Jorgenson Scissor	Strong Curve, 9"	1
34-1414	Ochsner Trocar	14 French, 10 Fr Catheter	1
90-3025	O'Sullivan O'Connor Abdominal Retractor	With 3 Interchangeable Blades, 2 Small / 1 Large	1
90-2680	Phaneuf Forceps	Angled, 8¼"	4
19-2280	Russian Tissue Forceps	8"	2
91-1298	Schroeder Tenaculum Forceps	Rounded Tips, 10"	2
90-5500	Sims Curettes	#0, Sharp	2
16-1615	Sklar Edge Mayo Dissecting Scissor	Tungsten Carbide, Curved, 9"	1
10-4394	SklarLite Mid-Size Container	18¼" x 11" x 6" Non-Perforated Bottom	1
10-4484	SklarLite Mid-Size Wire Basket	16" x 10" x 4"	1
91-1465	Staude-Moore Tenaculum Forceps	6½"	2

INSTRUMENT SETS

OB/GYN

98-1024 — HYSTERECTOMY SET

(FOR USE WITH MAJOR BASIC SET)

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
36-2295	Allis Tissue Forceps	5x6 Teeth, 9½"	2
60-3090	Deaver Retractor	#0, 1" x 9"	1
60-3112	Deaver Retractor	#1, 1" x 12"	1
60-3212	Deaver Retractor	#3, 1½" x 12"	1
52-5195	DeBakey Tissue Forceps	Straight, 2mm, 9½"	1
87-2195	Foerster Sponge Forceps	Straight, Serrated, 9½"	2
90-2280	Heaney Clamps	Double Tooth, 8"	6
20-4682	Heaney Needle Holders	8¼"	2
90-2480	Heaney-Ballantine Clamps	Straight, Longitudinal Serrations, 8"	2
90-2580	Heaney-Ballantine Clamp	Curved, Longitudinal Serrations, 8"	1
15-2567	Mayo Dissecting Scissors	Curved, 6¾"	1
55-9080	Metz-Nelson Scissors	Curved, 8"	1
55-2894	Mixer Right Angle Forceps	Longitudinal Serrations, 10½"	2
17-2680	Rochester-Ochsner Forceps	Curved, 1x2 Teeth, 8"	4
17-2880	Rochester-Ochsner Forceps	Straight, 1x2 Teeth, 8"	8
19-2280	Russian Tissue Forcep	8"	1
06-2915	Scalpel Handle	#3L	1
74-3075	Schmidt Hemostat Forceps	Half Curved Jaws, 7½"	2

98-1026 — SMALL GYN BIOPSY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
90-8335	Iris Hook	Quick Release, 3mm, 9"	1
90-6612	Kevorkian-Young Endocervical Curette	With Basket, 12"	1
90-3990	Kogan Endocervical Speculum	Narrow Wire Blades, With Gauge	1
90-7525	Sklar Oval Biopsy Forceps	3mm x 5mm	1
91-1290	Schroeder Braun Tenaculum Forcep	10"	1

98-1028 — CESAREAN SECTION SET

(FOR USE WITH MAJOR BASIC SET)

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
92-1012	DeLee OB Forcep	12"	1
87-2195	Foerster Sponge Forcep	Straight, Serrated, 9½"	1
90-3135	Jackson Vaginal Retractor	Medium Blade, 3½" x 1½"	1
20-2910	Masson Needle Holder	Serrated, 10½"	1
90-1590	Mayo Dissecting Scissor	Straight, 9"	1
90-1690	Mayo Dissecting Scissor	Curved, 9"	1
90-2680	Phaneuf Forceps	Angular, 8¼"	2
21-8060	TC Heaney Needle Holder	Curved, Serrated, 8¼"	1

www.sklarcorp.com



1.800.221.2166

OB/GYN 177

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

INSTRUMENT SETS

98-1018 — GYNECOLOGICAL BIOPSY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
90-3526	Auvard Weighted Speculum	Medium, 2½ lbs., 2¾" x 1 ½"	1
80-1410	Bozeman Forcep	Double Curved, 10¼"	1
90-7710	Bozeman Forcep	Double Curved (One Large Ring), 10¼"	1
87-2195	Foerster Sponge Forceps	Straight, Serrated, 9½"	2
87-2395	Foerster Sponge Forceps	Curved, Serrated, 9½"	2
90-7492	Gelhorn Biopsy Forceps	9½"	1
90-3700	Graves Vaginal Speculum	Small, ¾" x 3"	1
90-3701	Graves Vaginal Speculum	Medium, 1¼" x 4"	1
90-3702	Graves Vaginal Speculum	Large, 1¾" x 4½"	1
90-6612	Kevorkian-Young Endocervical Curette	With Basket, 12"	1
90-3995	Kogan Endocervical Speculum	Wide Wire Blades, With Gauge	1
90-6292	Novak Endometrial Curette	4mm, 9½"	1
91-1290	Schroeder Braun Tenaculum Forcep	10"	1
90-5413	Simpson Uterine Sound	CM Graduation	1
90-5423	Sims Curette	#3, Sharp, 11"	1
90-6887	Tischler Cervical Biopsy Forcep	3 x 9.5mm, 8½"	1

98-1020 — OBSTETRICAL SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
36-2380	Allis-Thoms Forcep	6x7 Teeth, 8"	1
47-2955	Backhaus Towel Clamps	5¼"	12
90-1755	Braun Episiotomy Scissor	Angled, 5½"	1
17-3062	Crile Forceps	Straight, 6¼"	3
17-3162	Crile Forceps	Curved, 6¼"	6
19-1060	Dressing Forceps	Serrated, 6"	2
90-3035	Eastman Retractor	Medium Blade, 3½" x 1½"	1
90-3040	Eastman Retractor	Large Blade, 4" x 1½"	1
87-2395	Foerster Sponge Forceps	Curved, Serrated, 9½"	2
91-4310	Iowa Puncturing Forcep	10"	1
92-2912	Laufe Divergent OB Forcep	12½"	1
15-1567	Mayo Dissecting Scissors	Straight, 6¾"	1
92-2215	McLane-Tucker-Luikart OB Forcep	15½"	1
14-1060	Operating Scissors	Straight, Sharp/Blunt, 6"	1
14-2060	Operating Scissors	Curved, Sharp/Blunt, 6"	1
17-2880	Rochester-Ochsner Forcep	Straight, 1x2 Teeth, 8"	1
17-2472	Rochester-Pean Forceps	Straight, 7¼"	3
17-2272	Rochester-Pean Forceps	Curved, 7¼"	3
19-2280	Russian Tissue Forceps	8"	2
92-1114	Simpson OB Forcep	Long, 14¼"	1
92-2314	Simpson-Luikart OB Forcep	14"	1
16-1610	Sklar Edge Mayo Dissecting Scissor	Tungsten Carbide, Curved, 6¾"	1
21-8032	TC Mayo-Hegar Needle Holders	Serrated, 7"	2
21-8033	TC Mayo-Hegar Needle Holder	Serrated, 8"	1
19-1260	Tissue Forceps	1x2 Teeth, 6"	2

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

INSTRUMENT SETS

OB/GYN

98-1016 — GYNECOLOGICAL SURGERY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
36-2275	Allis Tissue Forceps	5x6 Teeth, 7½"	2
36-2295	Allis Tissue Forceps	5x6 Teeth, 9½"	2
90-3526	Auvard Weighted Speculum	Medium, 2½ lbs., 2¾" x 1½"	1
47-2955	Backhaus Towel Clamps	5¼"	8
80-1410	Bozeman Forceps	Double Curved, 10"	2
17-3062	Crile Forceps	Straight, 6¼"	6
17-3162	Crile Forceps	Curved, 6¼"	6
90-2782	DeAlvarez Forceps	Double Curved, 8¼"	2
19-1080	Dressing Forcep	Serrated, 8"	1
90-3030	Eastman Retractor	Small Blade, 3" x 1½"	1
90-3035	Eastman Retractor	Medium Blade, 3½" x 1½"	1
90-3040	Eastman Retractor	Large Blade, 4" x 1½"	1
87-2195	Foerster Sponge Forceps	Straight, Serrated, 9½"	2
90-3700	Graves Vaginal Speculum	Small, ¾" x 3"	1
90-3701	Graves Vaginal Speculum	Medium, 1¼" x 4"	1
90-3702	Graves Vaginal Speculum	Large, 1⅝" x 4½"	1
90-4926	Hank Dilators	Set of 6	1
90-1977	Heaney Forceps	Light, Single Tooth, 7¾"	2
90-2177	Heaney Forceps	Heavy, Single Tooth, 8¼"	2
20-4682	Heaney Needle Holders	8¼"	2
90-1590	Mayo Dissecting Scissor	Straight, 9"	1
90-1690	Mayo Dissecting Scissor	Curved, 9"	1
20-2480	Mayo-Hegar Needle Holder	Serrated, 8"	1
55-9190	Metz-Nelson Scissor	Straight, 9"	1
55-9090	Metz-Nelson Scissor	Curved, 9"	1
14-2065	Operating Scissor	Curved, Sharp/Blunt, 6½"	1
90-2680	Phaneuf Forceps	Angled, 8¼"	2
34-2023	Poole Suction Tubes	23 French, Slip on End	2
17-1055	Pratt T-Shaped Clamps	Straight, 6"	2
17-2880	Rochester-Ochsner Forceps	Straight, 1x2 Teeth, 8"	3
17-2680	Rochester-Ochsner Forceps	Curved, 1x2 Teeth, 8"	3
17-2480	Rochester-Pean Forceps	Straight, 8"	3
17-2280	Rochester-Pean Forceps	Curved, 8"	3
19-2280	Russian Tissue Forcep	8"	1
19-2310	Russian Tissue Forcep	10"	1
06-2903	Scalpel Handle	#3	1
06-2915	Scalpel Handle	#3L	1
06-2904	Scalpel Handle	#4	1
06-2920	Scalpel Handle	#4L	1
91-1290	Schroeder Braun Tenaculum Forceps	10"	2
90-7995	Schweizer Uterine Forcep	Straight, 9½"	1
90-5413	Simpson Uterine Sound	CM Graduation	1
23-1106	Sklarlite Operating Scissor	Straight, Sharp/Blunt, 5½"	1
23-1102	Sklarlite Operating Scissor	Straight, Sharp/Sharp, 5½"	1
97-1560	TC Allis Grasping Forceps	7½"	2
21-8032	TC Mayo-Hegar Needle Holders	Serrated, 7"	2
21-8033	TC Mayo-Hegar Needle Holders	Serrated, 8"	2
51-4382	TC Potts-Smith Dressing Forcep	8¼"	1
19-1280	Tissue Forceps	1x2 Teeth, 8"	2
07-1805	Yankauer Suction Tube	Stainless Steel	1

www.sklarcorp.com



1.800.221.2166

OB/GYN 175

Office address: QUEEN CENTER, OFFICE No.206, 2ND Floor, MT. KHAN Road , Karachi. Page 107 of 372

Prime.instruments@hotmail.com: 03002621673

www.primeinstrumentintl.com

OB/GYN

INSTRUMENT SETS

98-1000 — MAJOR BASIC SURGERY SET (CONTINUED)

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
19-1255	Tissue Forceps	1x2 Teeth, 5½"	3
19-1280	Tissue Forcep	1x2 Teeth, 8"	1
60-1085	US Army Retractor	Set of 2, Double End, 8¼"	2
10-1805	Weinstein Sterilizing Racks	11¼" x 3¼"	3
07-1805	Yankauer Suction Tubes	Stainless Steel	2
60-3598	Zalkind Ribbon Retractors	Malleable, ¾" x 7¼"	2
60-3610	Zalkind Ribbon Retractors	Malleable, 1" x 13"	2
60-3615	Zalkind Ribbon Retractors	Malleable, 1½" x 13"	2

98-1014 — D AND C SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
90-3526	Auward Weighted Speculum	Medium, 2½ lbs., 2¾" x 1½"	1
47-2955	Backhaus Towel Clamps	5¼"	6
80-1410	Bozeman Forcep	Double Curved, 10¼"	1
17-3062	Crile Forceps	Straight, 6¼"	6
17-3162	Crile Forceps	Curved, 6¼"	6
19-1080	Dressing Forcep	Serrated, 8"	1
87-2395	Foerster Sponge Forceps	Curved, Serrated, 9½"	2
90-4713	Goodell Dilator	Large, 13"	1
90-3701	Graves Vaginal Speculum	Medium, 1¼" x 4"	1
90-3702	Graves Vaginal Speculum	Large, 1⅝" x 4½"	1
90-4826	Hegar Dilators	Set of 8	1
90-3130	Jackson Retractor	Small Blade, 3" x 1 ½"	1
90-3135	Jackson Retractor	Medium Blade, 3½" x 1½"	1
92-5095	Javerts Forcep	Serrated, 9½"	1
90-6612	Kevorkian-Younge Endocervical Curette	With Basket, 12"	1
90-3995	Kogan Endocervical Speculum	Wide Wire Blades, With Gauge	1
90-3710	Pederson Vaginal Speculum	Small, ½" x 3"	1
06-2903	Scalpel Handle	#3	1
06-2915	Scalpel Handle	#3L	1
06-2907	Scalpel Handle	#7	1
91-1290	Schroeder Braun Tenaculum Forcep	10"	1
90-5422	Sims Curette	#2, Sharp, 11"	1
90-5423	Sims Curette	#3, Sharp, 11"	1
90-5424	Sims Curette	#4, Sharp, 11"	1
90-5425	Sims Curette	#5, Sharp, 11"	1
90-5312	Sims Sound	Plain, 12½"	1
21-1260	TC Crile-Murray Needle Holder	Serrated, 6"	1
21-8032	TC Mayo-Hegar Needle Holder	Serrated, 7"	1
19-1280	Tissue Forcep	1x2 Teeth, 8"	1

Page 108 of 372

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

INSTRUMENT SETS

OB/GYN

98-1000 — MAJOR BASIC SURGERY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
47-1447	Adson Brown Tissue Forceps	4¾"	2
47-2147	Adson Forcep	Delicate, Serrated, 4¾"	1
36-2160	Allis Tissue Forceps	4x5 Teeth, 6"	3
36-2275	Allis Tissue Forceps	5x6 Teeth, 7½"	6
36-1662	Babcock Intestinal Forceps	6¼"	2
36-1677	Babcock Intestinal Forceps	8"	2
47-2955	Backhaus Towel Clamps	5¼"	12
60-6670	Balfour Abdominal Retractor	7"	1
17-3055	Crile Forceps	Straight, 5½"	12
17-3155	Crile Forceps	Curved, 5½"	12
60-3090	Deaver Retractors	#0, 1" x 9"	2
60-2713	Deaver Retractors	#2, 1" x 13"	2
60-3212	Deaver Retractors	#3, 1½" x 12"	2
60-3312	Deaver Retractor	#4, 2" x 12"	1
19-1055	Dressing Forceps	Serrated, 5½"	2
19-1080	Dressing Forceps	Serrated, 8"	2
87-2195	Foerster Sponge Forceps	Straight, Serrated, 9½"	2
87-2395	Foerster Sponge Forceps	Curved, Serrated, 9½"	2
17-1450	Halsted Mosquito Forceps	Straight, 5"	6
17-1550	Halsted Mosquito Forceps	Curved, 5"	6
17-1235	Hartmann Mosquito Forceps	Straight, Serrated, 3½"	3
17-1335	Hartmann Mosquito Forceps	Curved, Serrated, 3½"	3
15-1555	Mayo Dissecting Scissor	Straight, 5½"	1
15-2567	Mayo Dissecting Scissor	Curved, 6¾"	1
75-5670	Metzenbaum Scissors	Curved, 7"	2
55-9090	Metz-Nelson Scissor	Curved, 9"	1
30-1560	Probes	Double End, 6"	2
30-1960	Probes	With Eye, 6"	2
60-1681	Richardson Retractor	With Sklar-Grip Handle, Medium, 9½", 1" x 1¼"	1
60-1691	Richardson Retractor	With Sklar-Grip Handle, X-Large, 9½", ¾" x 2"	1
60-1695	Richardson-Eastman Retractor	Set of 2	1
17-2862	Rochester-Ochsner Forceps	Straight, 1x2 Teeth, 6¼"	2
17-2662	Rochester-Ochsner Forceps	Curved, 1x2 Teeth, 6¼"	2
17-2880	Rochester-Ochsner Forceps	Straight, 1x2 Teeth, 8"	2
17-2462	Rochester-Pean Forceps	Straight, 6¼"	6
17-2262	Rochester-Pean Forceps	Curved, 6¼"	6
17-2480	Rochester-Pean Forceps	Straight, 8"	3
17-2280	Rochester-Pean Forceps	Curved, 8"	3
06-2860	Ruler	Steel, 6"	1
19-2280	Russian Tissue Forceps	8"	2
06-2903	Scalpel Handles	#3	2
06-2904	Scalpel Handles	#4	2
74-3075	Schmidt Hemostat Forceps	Half Curved Jaws, 7½"	2
60-6164	Senn Retractors	Double End, Blunt, 6¼"	2
60-6162	Senn Retractors	Double End, Sharp, 6¼"	2
23-1106	Sklarite Operating Scissors	Straight, Sharp/Blunt, 5½"	2
23-1130	Sklarite Operating Scissors	Curved, Sharp/Blunt, 5½"	2
23-1157	Sklarite Stevens Tenotomy Scissor	Curved, Blunt/Blunt, 4½"	1
21-8031	TC Mayo-Hegar Needle Holders	Serrated, 6"	2

OB/GYN

SKLAR BLACK™ LASER

Needle Holders



MAYO-HEGAR NEEDLE HOLDERS

90-2287 7"

90-2288 8"



HEANEY NEEDLE HOLDER

90-2330 8 1/4" 2mm Curved



NEWPORT PROBE

45° ANGLE

90-9132 12"



CASTRO CALIPER

90-2264 3 1/4"

SKLAR BLACK™ LASER

OB/GYN

Forceps



TISSUE FORCEPS

1X2 TEETH

90-1020	4½"
90-1021	5"
90-1022	5½"
90-1023	6"



DEBAKEY FORCEPS

90-2310	7¾"	2mm
90-2312	9½"	2mm




POTTS-SMITH FORCEPS

90-2304	7"	Serrated
90-2306	10"	Serrated




GERALD FORCEPS

90-1025	7"	Straight	Serrated
90-1026	7"	Curved	Serrated
90-1027	7"	Straight	1x2 Teeth



TUTTLE/SINGLEY FORCEPS

90-2270	7"	Tuttle
90-2272	9"	Singley




CASTROVIEJO FORCEPS

90-2313	4 ¼"	.5mm Teeth
---------	------	------------

OB/GYN

SKLAR BLACK™ LASER

Forceps



ALLIS TISSUE FORCEPS

90-2276 7½" 5x6 Teeth

90-2279 9½" 5x6 Teeth



PEAN FORCEP

90-2292 7¼" Curved



ADSON FORCEPS

90-2261 4¾" Serrated

90-2262 4¾" 1x2 Teeth

90-1060 4¾" 1x2 Teeth

90-1061 4¾" 2x3 Teeth

90-1062 4¾" 7x7 Teeth



DRESSING FORCEPS

SERRATED

90-1063 4½"

90-1064 5"

90-1065 5½"

90-1066 6"

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SKLAR BLACK™ LASER

Forceps

OB/GYN



SCHWEIZER FORCEP

90-2350 9½" Straight



CLEMETSON FORCEP

NON RATCHETED

90-2348 9½" Straight



LAUFE UTERINE POLYP FORCEP

90-2345 7¾" Curved



**SCHROEDER TENACULUM
FORCEPS**

90-2325 10" Braun Flat Tips
90-2326 10" Pozzi Round Tips



**SCHROEDER VULSELLUM
FORCEPS**

90-1790 9½" Straight
90-1890 9½" Curved Sideways



BABCOCK FORCEPS

90-2278 8"
90-2295 9½"

OB/GYN

SKLAR BLACK™ LASER

Forceps



BOZEMAN FORCEP

90-2342 10¼" S-Curve, One Large Ring



BOZEMAN FORCEP

90-2340 10¼" S-Curve



CAMPION GRASPING FORCEP

90-2298 9¼" Curved



DELEE UTERINE FORCEP

90-7111 11" Heavy S-Curve

SKLAR BLACK™ LASER

OB/GYN


Hooks



GUTHRIE DOUBLE HOOKS

SHARP

- 90-1046 5" 3mm 2 Small Prongs
90-1044 5" 5mm 2 Large Prongs



JOSEPH SKIN HOOKS

SHARP

- 90-1050 6¼" 1 Sharp Prong
90-1051 6¼" 2mm 2 Prong
90-1052 6¼" 5mm 2 Prong
90-1053 6¼" 7mm 2 Prong
90-1054 6¼" 10mm 2 Prong
90-1055 6¼" 12mm 2 Prong



DERMAL RETRACTOR

- 90-1040 5¼" Sharp
90-1042 5¼" Blunt



SENN RETRACTOR


DOUBLE END

- 90-1057 6" Sharp
90-1059 6" Blunt



DESMARRES LID RETRACTORS

- 90-1038 12mm
90-1036 16mm
90-1035 18mm



DEPTH RULER IN MM

- 90-2242 12"

OB/GYN

SKLAR BLACK™ LASER

Hooks



IRIS HOOKS

QUICK RELEASE

90-2126	9"	1mm	90-2128	12"	1mm
90-2122	9"	3mm	90-2124	12"	3mm
90-2118	9"	4mm	90-2120	12"	4mm
90-2114	9"	6mm	90-2116	12"	6mm
90-2160	9"	Set of 4	90-2131	12"	Set of 4

EMMETT TENACULUM HOOKS

90-2171	9"	Style 1
90-2172	9"	Style 2
90-2173	9"	Style 3
90-2174	9"	Style 4
90-2175	9"	Style 5
90-2170	9"	Set of 5



GILLIES SKIN HOOKS

90-1045	7"	3mm	Sharp
90-1047	7"	4mm	Sharp

FRAZIER TISSUE HOOKS

90-1048	5"	Sharp	Extra Delicate
90-1049	5"	Blunt	Extra Delicate

TYRELL IRIS HOOKS

90-1056	5"	3mm	Blunt
90-1058	5"	3mm	Sharp

SKLAR BLACK™ LASER

OB/GYN

Specula



KOGAN ENDOCERVICAL SPECULUM

WITH GAUGE AND SCREWS

90-8855 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH SPRING HANDLE

90-8865 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH RATCHET HANDLE

90-8862 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH GAUGE AND SCREW

90-8858 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH SPRING HANDLE

90-8850 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM


WITH RATCHET HANDLE

90-8860 Wide Wire Blades

OB/GYN

SKLAR BLACK™ LASER

Specula



STEINER-AUVARD VAGINAL SPECULUMS

MEDIUM, EXTRA LONG BLADE

- | | | |
|---------|-----------|---|
| 90-3560 | 5½" x 1¼" | 2½ lbs. without loop |
| 90-3561 | 5½" x 1¼" | 2½ lbs. with loop to hold and direct Disposable Evacuation Tubing |
| 90-3562 | 5½" x 1¼" | 2½ lbs. with Permanent Smoke Evacuation Tube |



COLLIN VAGINAL SPECULUM

- | | | |
|---------|-----------|--|
| 90-3946 | 4¼" x 1¼" | With permanent Smoke Evacuation Tube (not illustrated) |
|---------|-----------|--|



LATERAL VAGINAL RETRACTOR

RATCHET WITH HANDLE

- | | |
|---------|------------------------------------|
| 90-8845 | Blade fits inside Vaginal Speculum |
|---------|------------------------------------|



JACKSON RETRACTORS SE

SINGLE END

- | | | |
|---------|-----------|--|
| 90-9064 | 3" x 1½" | Small with permanent Smoke Evacuation Tube |
| 90-9069 | 3½" x 1½" | Medium without Tube or Loops |
| 90-9065 | 3½" x 1½" | Medium with permanent Smoke Evacuation Tube |
| 90-9066 | 3½" x 1½" | Medium with loop to hold and direct disposable Smoke Evacuation Tubing |
| 90-9067 | 4" x 1½" | Large without Tube or Loops |
| 90-9068 | 4" x 1½" | Large with loop to hold and direct disposable Smoke Evacuation Tubing |
| 90-9070 | 4" x 1½" | Large with permanent Smoke Evacuation Tube |

SKLAR BLACK™ LASER

OB/GYN

Specula



GRAVES VAGINAL SPECULUMS WIDE VIEW PATTERN

WITH LARGE PROXIMAL OPENING FOR
INSTRUMENTATION

- | | | |
|---------|-------------|---|
| 90-9054 | 4" x 1 1/4" | With permanent
Smoke Evacuation Tube |
| 90-9055 | 4" x 1 1/4" | Medium, with loops to
hold and direct Smoke
Evacuation Tubing |



FREEWAY GRAVES VAGINAL SPECULUMS

OPEN SIDE, 90° ANGLE


- | | | |
|---------|-------------|--|
| 90-3754 | 4" x 1 1/4" | With permanent Smoke
Evacuation Tube
(not illustrated) |
| 90-3755 | 4" x 1 1/4" | With loop to hold and
direct Disposable
Smoke Evacuation
Tubing (not illustrated) |



FREEWAY GRAVES VAGINAL SPECULUMS


OPEN SIDE, 45° ANGLE

- | | | |
|---------|-------------|--|
| 90-3814 | 4" x 1 1/4" | With permanent
Smoke Evacuation
Tube (not illustrated) |
| 90-3815 | 4" x 1 1/4" | With loop to hold and
direct Disposable
Smoke Evacuation
Tubing (not illustrated) |




AUVARD VAGINAL SPECULUMS MEDIUM

- | | | |
|---------|-----------------|---|
| 90-9100 | 2 3/4" x 1 1/2" | 2 1/2 lbs. without loop |
| 90-9102 | 2 3/4" x 1 1/2" | 2 1/2 lbs. with loop to
hold and direct
Disposable Evacuation
Tubing |
| 90-9103 | 2 3/4" x 1 1/2" | 2 1/2 lbs. with Permanent
Smoke Evacuation
Tube |



AUVARD VAGINAL SPECULUMS MEDIUM, LONG BLADES

- | | | |
|---------|-------------|---|
| 90-9106 | 4" x 1 3/4" | 2 1/2 lbs. Without loop |
| 90-9108 | 4" x 1 3/4" | 2 1/2 lbs. With loop to
hold and direct
Disposable Evacuation
Tubing |
| 90-9109 | 4" x 1 3/4" | With permanent
Smoke Evacuation
Tube |



DISPOSABLE SMOKE EVACUATION TUBING

- | | | |
|---------|------------------------------|-------|
| 93-2160 | 10' precut 3/16" ID, 1/4" OD | 50/cs |
|---------|------------------------------|-------|

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

SKLAR BLACK™ LASER

Specula



**GRAVES VAGINAL SPECULUMS
WITH PERMANENT SMOKE EVACUATION TUBE**

90-9015	3" x ¾"	Small
90-9016	4" x 1¼"	Medium
90-9018	4½" x 1⅝"	Large
90-9020	4¾" x 1½"	Extra Large



**GRAVES VAGINAL SPECULUMS
WITH LOOPS FOR DISPOSABLE SMOKE EVACUATION
TUBING**

90-9005	3" x ½"	Extra Small
90-9000	3" x ¾"	Small
90-9002	4" x 1¼"	Medium
90-9004	4½" x 1⅝"	Large
90-9006	4¾" x 1½"	Extra Large



**GRAVES VAGINAL SPECULUMS
WITHOUT TUBE OR LOOPS**

90-3812	3" x ½"	Extra Small
90-3804	3" x ¾"	Small
90-3806	4" x 1¼"	Medium
90-3808	4½" x 1⅝"	Large
90-3810	4¾" x 1½"	Extra Large



**PEDERSON VAGINAL SPECULUMS
WITH PERMANENT SMOKE EVACUATION TUBE**

90-9040	3" x ½"	Small
90-9050	4" x ⅝"	Slender
90-9042	4" x ⅞"	Medium
90-9044	4¾" x 1"	Large
90-9048	4¾" x 1⅝"	Extra Large



**PEDERSON VAGINAL SPECULUMS
WITH LOOPS FOR DISPOSABLE SMOKE EVACUATION
TUBING**

90-9035	2¾" x ½"	Extra Small
90-9030	3" x ½"	Small
90-9038	4" x ⅝"	Slender
90-9032	4" x ⅞"	Medium
90-9034	4¾" x 1"	Large



**PEDERSON VAGINAL SPECULUMS
WITHOUT TUBE OR LOOPS**

90-3841	2¾" x ½"	Extra Small
90-3833	3" x ½"	Small
90-3837	4" x ⅝"	Slender
90-3839	4" x ⅞"	Medium
90-3835	4¾" x 1"	Large

ELECTROSURGERY

OB/GYN

Sklar Blue™ Forceps



SCHWEIZER FORCEP

91-5404 9½" Straight



CLEMETSON FORCEP

WITHOUT RATCHET

91-5400 9½" Straight



LAUE UTERINE POLYP FORCEP

91-5408 7¾"



SCHROEDER TENACULUM FORCEPS

91-5450 10" Braun Flat Tips
91-5452 10" Pozzi Round Tips



BABCOCK FORCEPS

91-5268 8"
91-5270 9½"



ALLIS TISSUE FORCEPS

91-5278 7½" 5x6 Teeth
91-5440 9½" 5x6 Teeth

OB/GYN

ELECTROSURGERY

Sklar Blue™ Endocervical Speculums



KOGAN ENDOCERVICAL SPECULUM

WITH GAUGE AND SCREW

91-5230 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH SPRING HANDLE


91-5228 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH RATCHET HANDLE


91-5235 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH GAUGE AND SCREW

91-5232 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH SPRING HANDLE

91-5226 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH RATCHET HANDLE

91-5234 Wide Wire Blades

OB/GYN

SKLARTECH 5000

Laparoscopic Sets

98-1354 - BAND SET

PRODUCT NUMBER	DESCRIPTION	QUANTITY
31-2238-45	Suction/Irrigation System with 45mm Cannula Dual Trumpet Valve	1
31-9012XC-45	Metzenbaum Scissors, Curved - 5mm, 45cm	1
31-9030XC-45	Maryland Dissector - 5mm, 45cm	1
31-9090XC-45	Fenestrated Grasper, 22mm Jaw - 5mm, 45cm	1
31-9095XC-45	Fenestrated Grasper, 40mm Jaw - 5mm, 45cm	1
31-4565	Esophagus Retractor - 5mm	1

98-1355 - CHOLECYSTECTOMY SET

PRODUCT NUMBER	DESCRIPTION	QUANTITY
31-2227	Suction Control Handle for Spatula Electrode	2
31-2235	Spatula Electrode with Suction - 5mm	1
31-2237	L Hook Electrode with Suction - 5mm	1
31-2239	High Volume Suction Irrigation System	1
31-2355	Needle for Injection and Irrigation - 5mm	1
31-9005XC	Metzenbaum Scissors, Straight - 5mm, 33cm	1
31-9010XC	Metzenbaum Scissors, Curved - 5mm, 33cm	1
31-9015XC	Hook Scissors - 5mm, 33cm	1
31-9025XC	Micro Scissors, Curved - 5mm, 33cm	1
31-9030XC	Maryland Dissector - 5mm, 33cm	1
31-9035XC	Mixer Forceps, Diamond Serrations - 5mm, 33cm	1
31-9040XC	Dolphin Forceps with Spoon - 5mm, 33cm	1
31-9050YC	Standard Grasping Forceps - 5cm, 33cm	1
31-9060YC	Wave Grasper - 5mm, 33cm	1
31-9100YC	Maxi (Fundus) Grasper - 5mm, 33cm	1
31-9125XC	Mixer Forceps - 10mm, 33cm	1
31-9130XC	Claw Forceps - 10mm, 33cm	1
31-9135XC	Spoon Forceps - 10mm, 33cm	1

98-1356 - APPENDECTOMY SET

PRODUCT NUMBER	DESCRIPTION	QUANTITY
31-2239	High Volume Suction Irrigation System	1
31-2552	TC Needle Holder, Curved Left - 5mm, 33cm	1
31-9010XC	Metzenbaum Scissors, Curved - 5mm, 33cm	1
31-9025XC	Micro Scissors, Curved - 5mm, 33cm	1
31-9050YC	Standard Grasping Forceps - 5mm, 33cm	1
31-9085YC	Pointed Grasping Forceps with Spoon - 5mm, 33cm	1
31-9125XC	Mixer Forceps - 10mm, 33cm	1
31-9130XC	Claw Forceps - 10mm, 33cm	1

98-1357 - BOWEL RESECTION SET

PRODUCT NUMBER	DESCRIPTION	QUANTITY
31-2552	TC Needle Holder, Curved Left - 5mm, 33cm	1
31-2235	Spatula Electrode with Suction - 5mm	1
31-2227	Suction Control Valve for Spatula Electrode	1
31-9067YC	Allis Forceps TC - 5mm, 33cm	1
31-9080YC	Grasping Forcep with Spoon - 5mm, 33cm	1
31-9051YC	Tapering Grasping Forcep - 5mm, 33cm	1
31-9040XC	Dolphin Forcep with Spoon - 5mm, 33cm	1
31-9140YC	Babcock Forcep - 10mm, 33cm	1
31-9095YC	Fenestrated Grasper, 40mm Jaws - 5mm, 33cm	1
31-9010XC	Metzenbaum Scissors, Curved - 5mm 33cm	1
31-9025XC	Micro Scissors, Curved - 5mm, 33cm	1

SKLARTECH 5000

OB/GYN

Miniature Laparoscopic Accessories

MINIATURE ENDOSCOPE

COMPLETE INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-4301	2.7mm	30cm	Miniature, Autoclavable



MINIATURE TROCAR

COMPLETE INSTRUMENT	REPLACEMENT OBTURATOR	DIAMETER	FEATURES
31-4302	31-4318	3mm	Mini, with Trap Door Cannula, Resusable, Includes 3mm Obturator
31-4320	-	3mm	External Seal, Red For use with Standard and Hasson Cannulas
31-4321	-	3mm	Internal Seal, Gray For use with Trap Door Cannulas only




ELECTROSURGERY

INSULATED ELECTRODES WITHOUT SUCTION

Each Electrode has a permanently fitted adapter for connecting to the Electrosurgical Cable. Insulated Electrodes and Cables sold separately.

INSTRUMENT	DIAMETER	FEATURES
31-4325 (A)	3mm	Mini, Spatula
31-4326 (B)	3mm	Mini, L-Hook
31-4327 (C)	3mm	Mini, J-Hook



ELECTROSURGICAL CABLE

INSTRUMENT	LENGTH	FEATURES
31-2242	277cm	Monopolar



OB/GYN

SKLARTECH 5000

Miniature Laparoscopic Instruments

GRASPING FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-4308YC	31-4308	3mm	33cm	Mini Grasping Forcep, YC Handle
31-4311XC	31-4311	3mm	33cm	Mini Atraumatic, Dissector, Straight, XC Handle
31-4312YC	31-4312	3mm	33cm	Mini Grasping Forcep, Heavy Serrations, YC Handle
31-4313YC	31-4313	3mm	33cm	Mini Cobra, Heavy Serrations, YC Handle
31-4317YC	31-4317	3mm	33cm	Mini Fenestrated, 8mm Double Action Jaws, YC Handle



PUNCH


COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-4305XC	31-4305	3mm	33cm	Mini Punch, Delicate



NEEDLE HOLDER

NON-INSULATED WITH TUNGSTEN CARBIDE JAW INSERTS AND AXIAL HANDLE

COMPLETE INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-9360AX	3mm	33cm	Mini, Straight, Axial Ring Handle
31-9362AX	3mm	33cm	Mini, Curved Right, Axial Ring Handle
31-9364AX	3mm	33cm	Mini, Curved Left, Axial Ring Handle



SKLARTECH 5000

OB/GYN

Miniature Laparoscopic Instruments

METZENBAUM SCISSORS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-4306XC	31-4306	3mm	33cm	Mini, Straight, XC Handle
31-4303XC	31-4303	3mm	33cm	Mini, Curved, XC Handle
31-4307XC	31-4307	3mm	33cm	Mini, Curved, Serrated, XC Handle



HOOK SCISSORS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-4304XC	31-4304	3mm	33cm	Mini, Delicate, Straight, XC Handle



MARYLAND DISSECTOR

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-4310XC	31-4310	3mm	33cm	Mini, Cross Serrations, Curved, XC Handle



FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-4314YC	31-4314	3mm	33cm	Mini Duck Nose, Cross Serrations, Straight, YC Handle
31-4309XC	31-4309	3mm	33cm	Mini Kelly, Double Action Delicate Jaws, XC Handle
31-4315XC	31-4315	3mm	33cm	Mini Spoon, XC Handle
31-4316XC	31-4316	3mm	33cm	Mini Spoon, Serrated, XC Handle



OB/GYN


SKLARTECH 5000

Laparoscopic Electrosurgery

Each coagulating Electrode is able to connect to the suction control valve of your choice (Listed below). Insulated Electrodes and Control Valves sold separately.


INSULATED ELECTRODES, WITH SUCTION

INSTRUMENT	DIAMETER	FEATURES
31-2235 (A)	5mm	Spatula
31-2237 (B)	5mm	L-Hook
31-2241 (C)	5mm	Knife
31-2243 (D)	5mm	J-Hook
31-2244 (E)	5mm	Blunt



HIGH FREQUENCY SUCTION CONTROL VALVES


INSTRUMENT	DIAMETER	FEATURES
31-2227 (A)	5mm	Insulated Suction Control with Single Trumpet Valve
31-2233 (B)	5mm	Insulated Suction Control with Double Trumpet Valve
31-2229 (C)	5mm	Insulated Suction Control with Double Trumpet Valve and Handle



INSULATED ELECTRODES, WITHOUT SUCTION

Each Electrode has a permanently fitted adapter for connecting to the Electrosurgical Cable. Insulated Electrodes and Cables sold separately.

INSTRUMENT	DIAMETER	FEATURES
31-2435 (A)	5mm	Spatula
31-2437 (B)	5mm	L-Hook
31-2443 (C)	5mm	J-Hook
31-2445 (D)	5mm	Button
31-2447 (E)	5mm	Cone



ELECTROSURGICAL CABLE

INSTRUMENT	LENGTH	FEATURES
31-2242	277cm	Monopolar



SKLARTECH 5000

OB/GYN

Laparoscopic Needle Holders



TC NEEDLE HOLDERS

NON-INSULATED WITH TUNGSTEN CARBIDE JAW INSERTS AND AXIAL HANDLE

COMPLETE INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-2550	5mm	33cm	Standard, Straight, Flushable Shaft, 10mm Jaw, Axial Handle
31-2550-45	5mm	45cm	Standard, Straight, Flushable Shaft, 10mm Jaw, Axial Handle
31-2551	5mm	33cm	Standard, Curved Right, Flushable Shaft, 10mm Jaw, Axial Handle
31-2552	5mm	33cm	Standard, Curved Left, Flushable Shaft, 10mm Jaw, Axial Handle
31-2552-45	5mm	45cm	Standard, Curved Left, Flushable Shaft, 10mm Jaw, Axial Handle
31-2553	5mm	33cm	Standard, Straight, Solid Shaft, 8mm Jaw, Axial Handle

STRAIGHT



CURVED RIGHT



CURVED LEFT



OB/GYN

SKLARTECH 5000

Laparoscopic Suction/Irrigation



SUCTION IRRIGATION SYSTEM

COMPLETE INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-2238	5mm	33cm	Dual Trumpet Valve (also fits 10mm Sheath)
31-2238-45	5mm	45cm	Dual Trumpet Valve (also fits 10mm Sheath)
31-2258	5mm	33cm	Single Trumpet Valve (also fits 10mm Sheath)
31-4003	5mm	33cm	Sheath Only
31-4004	10mm	33cm	Sheath Only
31-2205	5mm	33cm	Tube Channel, for Laser-Fiber
31-2238OR	-	-	O-Ring, Replacement Pk/10



HIGH-VOLUME SUCTION IRRIGATION SYSTEM

COMPLETE INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-2239	5mm 10mm	33cm	Dual Trumpet Valve (includes both 5mm and 10mm Sheaths)
31-2347	5mm	33cm	Sheath Only
31-2349	10mm	33cm	Sheath Only

SKLARTECH 5000

OB/GYN

Laparoscopic Retractors




ESOPHAGUS RETRACTOR

INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-4565	5.5mm	33cm	For Gastric Band Application



FAN RETRACTOR

INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-4562	5.5mm	33cm	3 Finger Articulating Fan, Slightly Bent



FAN RETRACTOR

INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-4563	10mm	33cm	5 Finger Articulating Fan, Slightly Bent

OB/GYN

SKLARTECH 5000

Laparoscopic Accessories



VERRES NEEDLE WITH STOPCOCK

INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-2327	2.3mm	100mm	Stopcock
31-2350	2.3mm	120mm	Stopcock
31-2352	2.3mm	150mm	Stopcock



NEEDLE FOR INJECTION AND IRRIGATION

INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-2355	5mm	33cm	Luer Lock, 2mm Hollow Needle
31-2355-45	5mm	45cm	Luer Lock, 2mm Hollow Needle



GRADUATED PROBE


INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-2365	5.5mm	33cm	Graduated

SKLARTECH 5000 Laparoscopic Trocars

OB/GYN

OBTURATORS

INSTRUMENT	LENGTH	FEATURES
31-0004	11mm	Blunt Tip Obturator for 11mm Hasson Cannula
31-1011	11mm	Extremely Sharp Pyramidal Tip Obturator for Standard Trocar
31-1010	5.5mm	Extremely Sharp Pyramidal Tip Obturator for Standard Trocar



THREADED CANNULA SLEEVE

INSTRUMENT	LENGTH	FEATURES
31-1020	5.5mm	Interchangeable with Standard Smooth Cannulas



SNAP ON FLIP REDUCER

INSTRUMENT	LENGTH	FEATURES
31-1042		Reduces from 11mm to 5.5mm



S-RETRACTOR

INSTRUMENT	LENGTH	FEATURES
60-6130	6 1/4"	Set of 2



TROCAR CANNULA PARTS

FOR USE WITH STANDARD & HASSON VALVE CANNULAS				
EXTERNAL SEALS	5.5mm	11mm	11mm (Hasson)	11mm (Hasson)
	RED	BLUE	BLUE	GRAY
	31-0003	31-0001	31-0001	31-0002
FOR USE WITH TRAP DOOR CANNULAS				
INTERNAL SEAL	5.5mm		11mm	
	GRAY		YELLOW	
	31-4323		31-4324	
FOR USE WITH SILICON VALVE CANNULAS				
SILICON VALVE	5.5mm		11mm	
	31-4328		31-4329	

OB/GYN

SKLARTECH 5000

Laparoscopic Trocars

HASSON TROCAR WITH TRUMPET VALVE

COMPLETE INSTRUMENT	REPLACEMENT OBTURATOR	LENGTH	FEATURES
31-1002	31-0004	11mm	Hasson Trumpet Valve Cannula with 11mm Blunt Tip Obturator



STANDARD TROCAR WITH TRUMPET VALVE

COMPLETE INSTRUMENT	REPLACEMENT OBTURATOR	LENGTH	FEATURES
31-1007	31-1010	5.5mm	Trumpet Valve Cannula With 5.5mm Sharp Obturator



STANDARD TROCAR WITH TRAP DOOR

COMPLETE INSTRUMENT	REPLACEMENT OBTURATOR	LENGTH	FEATURES
31-1004	31-1011	11mm	Trap Door Cannula With 11mm Sharp Obturator



STANDARD TROCAR WITH AUTOMATIC VALVE

COMPLETE INSTRUMENT	REPLACEMENT INSERT	LENGTH	FEATURES
31-1015	31-1010	5.5mm	Automatic Valve Cannula with 5.5mm Sharp Obturator, Dismountable
31-1016	31-1011	11mm	Automatic Valve Cannula with 11mm Sharp Obturator, Dismountable



STANDARD TROCAR WITH SILICON VALVE

COMPLETE INSTRUMENT	REPLACEMENT INSERT	LENGTH	FEATURES
31-1017	31-1010	5.5mm	Silicon Valve Cannula with LL Stopcock & 5.5mm Sharp Obturator, Dismountable
31-1018	31-1011	11mm	Silicon Valve Cannula with LL Stopcock & 11mm Sharp Obturator, Dismountable



SKLARTECH 5000

Handles & Cannulas

OB/GYN

SKLAR "XC" CARBON HANDLE

INSTRUMENT

31-9001XC

FEATURES

Standard, non-ratchet, lightweight, carbon fiber reinforced non-conductive safety handle with locking collar handle. Fully rotating with a proximally positioned high-frequency connector.

SKLAR "YC" CARBON HANDLE

INSTRUMENT

31-9002YC

FEATURES

Standard, ratchet, lightweight, carbon fiber reinforced non-conductive safety handle with locking collar handle. Fully rotating with a proximally positioned high-frequency connector.

REPLACEMENT LOCKING COLLAR

INSTRUMENT

31-9160

FEATURES

For Sklar Standard Handles



Instruments may be ordered as complete sets or assembled from components ordered separately. Optional handles for customizing instruments - interchangeable for 3mm, 5mm or 10mm - are available for special order. Contact your Sklar Sales Representative for details.

CANNULAS

INSTRUMENT DIAMETER LENGTH FEATURES

INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-9155	5mm	33cm	Insulated
31-9155-45	5mm	45cm	Insulated

INSTRUMENT	DIAMETER	LENGTH	FEATURES
31-9150	10mm	33mm	Insulated With Permanent Locking Collar Attached
31-9150-45	10mm	45mm	Insulated With Permanent Locking Collar Attached



OB/GYN

SKLARTECH 5000

Laparoscopic Punches

PUNCH

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9046XC	31-9046	5mm	33cm	Probe Excision



BIOPSY PUNCH

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9047XC	31-9047	5mm	33cm	No Tooth, XC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9110YC	31-9110	5mm	33cm	With Tooth, YC Handle
31-9110YC-45	31-9110-45	5mm	45cm	With Tooth, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9115XC	31-9115	5mm	33cm	Spoon with Teeth, XC Handle
31-9115YC	31-9115	5mm	33cm	Spoon with Teeth, YC Handle
31-9115YC-45	31-9115-45	5mm	45cm	Spoon with Teeth, YC Handle



SKLARTECH 5000

OB/GYN

Laparoscopic Grasping Forceps

ATRAUMATIC GRASPERS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9098YC	31-9098	5mm	33cm	Single Action, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9097XC	31-9097	5mm	33cm	Diamond Serrations, XC Handle
31-9097YC	31-9097	5mm	33cm	Diamond Serrations, YC Handle



FENESTRATED GRASPING DISSECTOR

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9042XC	31-9042	5mm	33cm	16mm Double Action Jaw, XC Handle



FENESTRATED GRASPING FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9090YC	31-9090	5mm	33cm	22mm Double Action Jaw, YC Handle
31-9090XC-45	31-9090-45	5mm	45cm	22mm Double Action Jaw, XC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9095XC	31-9095	5mm	33cm	40mm Double Action Jaw, XC Handle
31-9095YC	31-9095	5mm	33mm	40mm Double Action Jaw, YC Handle
31-9095XC-45	31-9095-45	5mm	45mm	40mm Double Action Jaw, XC Handle
31-9095YC-45	31-9095-45	5mm	45mm	40mm Double Action Jaw, YC Handle



OB/GYN

SKLARTECH 5000

Laparoscopic Grasping Forceps

WAVE FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9060YC	31-9060	5mm	33cm	Atraumatic Wavelike Serrations, Single Action, YC Handle
31-9060YC-45	31-9060-45	5mm	45cm	Atraumatic Wavelike Serrations, Single Action, YC Handle
31-9061YC	31-9061	5mm	33cm	Atraumatic Wavelike Serrations, Double Action, YC Handle



DEBAKEY FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9065YC	31-9065	5mm	33cm	40mm Jaws, YC Handle
31-9065XC	31-9065	5mm	33mm	40mm Jaws, XC Handle
31-9065YC-45	31-9065-45	5mm	45mm	40mm Jaws, YC Handle
31-9065XC-45	31-9065-45	5mm	45mm	40mm Jaws, XC Handle



BABCOCK/DEBAKEY FORCEP

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9142XC	31-9142	10mm	33cm	Atraumatic Serrations, XC Handle



COBRA FORCEP

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9075YC	31-9075	5mm	33cm	2x3 Teeth, Double Action, YC Handle



ALLIGATOR FORCEP

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9070YC	31-9070	5mm	33cm	Double Action, YC Handle



SKLARTECH 5000

OB/GYN

Laparoscopic Grasping Forceps

GRASPING FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9050XC	31-9050	5mm	33cm	Standard, Cross Serrations, XC Handle
31-9050YC	31-9050	5mm	33cm	Standard, Cross Serrations, YC Handle
31-9050YC-45	31-9050-45	5mm	45cm	Standard, Cross Serrations, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9051YC	31-9051	5mm	33cm	Tapered, Rounded Jaws, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9055YC	31-9055	5mm	33cm	Endo Clinch, Extremely Atraumatic Serrations, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9080YC	31-9080	5mm	33mm	With Spoon, YC Handle
31-9080XC	31-9080	5mm	33mm	With Spoon, XC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9085YC	31-9085	5mm	33mm	Pointed with Spoon, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9100YC	31-9100	5mm	33mm	Maxi (Fundus), YC Handle



OB/GYN

SKLARTECH 5000

Laparoscopic Forceps

BABCOCK FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9105YC	31-9105	5mm	33cm	YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9140XC	31-9140	10mm	33cm	XC Handle
31-9140YC	31-9140	10mm	33cm	YC Handle
31-9140YC-45	31-9140-45	10mm	45cm	YC Handle



ALLIS FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9067YC	31-9067	5mm	33cm	TC Inserted Jaws, Atraumatic Clamp, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9144YC	31-9144	10mm	33cm	TC Inserted Jaws, Diamond Serrations, YC Handle



STRONG ATRAUMATIC CLAMP

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9145YC	31-9145	10mm	33cm	40mm Jaws, Straight, YC Handle



NEEDLE HOLDER

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9112YC	31-9112	5mm	33cm	10mm Jaws, Pointed, YC Handle



SKLARTECH 5000 Laparoscopic Forceps

OB/GYN

KELLY FORCEP

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9043XC	31-9043	5mm	33cm	Double Action, XC Handle



GLASSMANN FORCEP

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9096YC	31-9096	5mm	33cm	Atraumatic Clamp, YC Handle



TENACULUM FORCEP

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9111YC	31-9111	5mm	33cm	1x1 Sharp Teeth, YC Handle
31-9113YC	31-9113	10mm	33cm	1x1 Sharp Teeth, YC Handle



KOCHER FORCEP

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9128XC	31-9128	10mm	33cm	60mm Jaws, Axial Serration, XC Handle



SPOON FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9135XC	31-9135	10mm	33cm	XC Handle



CLAW FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9045XC	31-9045	5mm	33cm	Traumatic Jaws, XC Handle
31-9045YC	31-9045	5mm	33cm	Traumatic Jaws, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9130XC	31-9130	10mm	33cm	Traumatic Jaws, XC Handle



OB/GYN

SKLARTECH 5000 Laparoscopic Forceps

MARYLAND DISSECTORS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9030XC	31-9030	5mm	33cm	Cross Serration, XC Handle
31-9030XC-45	31-9030-45	5mm	45cm	Cross Serration, XC Handle
COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9031XC	31-9031	5mm	33cm	Diamond Serration, XC Handle
31-9120XC	31-9120	10mm	33cm	Diamond Serration, XC Handle



MIXTER FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9033XC	31-9033	5mm	33cm	90° Angle, Cross Serration, XC Handle
COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9035XC	31-9035	5mm	33cm	90° Angle, Diamond Serration, XC Handle
31-9125XC	31-9125	10mm	33cm	90° Angle, Diamond Serration, XC Handle
31-9125XC-45	31-9125-45	10mm	45cm	90° Angle, Diamond Serration, XC Handle



DISSECTING FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9036XC	31-9036	5mm	33cm	Curved, Serrated, XC Handle



DOLPHIN FORCEPS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9040XC	31-9040	5mm	33cm	Short, XC Handle



PRIME INSTRUMENTS PVT. LTD.

SKLARTECH 5000 Laparoscopic Scissors

OB/GYN

HOOK SCISSORS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9016XC	31-9016	5mm	33cm	Standard, Flush Tips, XC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9015XC	31-9015	5mm	33cm	Sharp Cutting Edge, Overlapping Tips, XC Handle
31-9015XC-45	31-9015-45	5mm	45cm	Sharp Cutting Edge, Overlapping Tips, XC Handle



MICRO SCISSORS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9020XC	31-9020	5mm	33cm	Straight, Extremely Delicate, XC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9025XC	31-9025	5mm	33cm	Curved, Extremely Delicate, XC Handle



OB/GYN

SKLARTECH 5000

Laparoscopic Scissors

PERITONEAL SCISSORS

COMPLETE INSTRUMENT	REPLACEMENT	DIAMETER	LENGTH	FEATURES
31-9026XC	31-9026	5mm	33cm	Straight, Blunt Tips, XC Handle



OPERATING SCISSORS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9027XC	31-9027	10mm	33cm	Straight, XC Handle



METZENBAUM SCISSORS

COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9005XC	31-9005	5mm	33cm	Straight, XC Handle
31-9005YC	31-9005	5mm	33cm	Straight, YC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9010XC	31-9010	5mm	33cm	Curved, XC Handle
31-9010XC-45	31-9010-45	5mm	45cm	Curved, XC Handle



COMPLETE INSTRUMENT	REPLACEMENT INSERT	DIAMETER	LENGTH	FEATURES
31-9012XC	31-9012	5mm	33cm	Curved, Serrated, XC Handle
31-9012XC-45	31-9012-45	5mm	45cm	Curved, Serrated, XC Handle



SKLARTECH 5000

Laparoscopy

OB/GYN

SKLARTECH 5000


The Sklartech 5000 Laparoscopic Instrument line provides a full range of the most popular patterns. Instruments can be purchased as a complete set or individually. When purchased as a complete instrument, Sklartech Laparoscopic Scissors, Dissectors and Forceps include the insert with tip, insulated cannula and carbon fiber handle with locking collar.



OB/GYN


NEEDLE HOLDERS

Econo™ Sterile




DERF NEEDLE HOLDERS

96-2583 4 3/4" Serrated 25/cs
96-2582 4 3/4" Serrated 50/cs




WEBSTER NEEDLE HOLDERS

96-2587 5 1/4" Smooth 25/cs
96-2586 5 1/4" Smooth 50/cs
96-2659 5 1/4" Smooth 100/cs




HALSEY NEEDLE HOLDERS

96-2585 5" Smooth 25/cs
96-2584 5" Smooth 50/cs




BAUMGARTNER NEEDLE HOLDERS

96-2581 5 1/2" Serrated 25/cs
96-2580 5 1/2" Serrated 50/cs




MAYO-HEGAR NEEDLE HOLDERS

96-2595 6" Serrated 25/cs
96-2594 6" Serrated 50/cs



CRILE-WOOD NEEDLE HOLDERS

96-2591 6" Serrated 25/cs
96-2590 6" Serrated 50/cs




WIREFORM NEEDLE HOLDERS

96-8968 5 1/4" Smooth 50/cs
96-8974 5 1/4" Serrated 50/cs

NEEDLE HOLDERS


OB/GYN

SURGI-OR™




DERF NEEDLE HOLDER

95-826 4 3/4" Serrated




WEBSTER NEEDLE HOLDERS

95-809 5 1/4" Smooth
95-813 5 1/4" Serrated




HALSEY NEEDLE HOLDERS

95-810 5 1/4" Smooth
95-812 5 1/4" Serrated



BAUMGARTNER NEEDLE HOLDER

95-832 5 1/4" Serrated




MAYO-HEGAR NEEDLE HOLDERS

SERRATED
95-842 5 1/4"
95-847 6"
95-853 7"
95-869 8"


OB/GYN

NEEDLE HOLDERS

Surgi-OR™ Tungsten Carbide



TC DERF NEEDLE HOLDER
95-136 4¾" Serrated




TC WEBSTER NEEDLE HOLDER
95-130 5¼" Smooth




TC HALSEY NEEDLE HOLDER
95-132 5¼" Smooth
95-133 5¼" Serrated




TC OLSEN-HEGAR NEEDLE HOLDERS
95-146 4¾" Serrated, Delicate
95-147 5½" Serrated
95-148 6½" Serrated, Delicate
95-149 7½" Serrated



TC BAUMGARTNER NEEDLE HOLDER
95-135 5¼" Serrated



TC MAYO-HEGAR NEEDLE HOLDERS
SERRATED
95-140 5¼"
95-141 6"
95-142 7"
95-143 8"




TC CRILE-WOOD NEEDLE HOLDERS
SERRATED
95-137 6" Narrow
95-138 7" DeBaKey
95-139 8" DeBaKey

NEEDLE HOLDERS


OB/GYN

Microsurgical




**MICRO SURGERY
NEEDLE HOLDERS**
ROUND HANDLE WITHOUT LOCK

98-3180	5¼"	Curved
98-3186	7⅞"	Straight
98-3188	7⅞"	Curved




**MICRO SURGERY
NEEDLE HOLDERS**
ROUND HANDLE WITH LOCK

98-3190	7⅞"	Straight
98-3192	7⅞"	Curved




**CASTROVIEJO
NEEDLE HOLDERS**

65-6350	5¼"	Straight w/o Lock
65-6352	5¼"	Straight w/Lock
65-6354	5¼"	Curved w/Lock



**CASTROVIEJO
NEEDLE HOLDERS**
WITH LOCK

65-6357	5½"	Straight
65-6358	5½"	Curved




**CASTROVIEJO KALT
NEEDLE HOLDERS**

65-6055	5½"	Serrated
---------	-----	----------



**BARRAQUER
NEEDLE HOLDERS**
WITHOUT CATCH

65-6452	5¼"	Curved, Delicate
---------	-----	------------------



**BARRAQUER
NEEDLE HOLDERS**
WITH CATCH

65-6448	5¼"	Curved, Delicate
65-6446	5½"	Curved, Delicate

OB/GYN

NEEDLE HOLDERS


General Surgery




NEW ORLEANS NEEDLE HOLDER
20-1870 7" Serrated



CRILE NEEDLE HOLDER
20-1160 6" Serrated




SAROT NEEDLE HOLDER
22-5050 5" Serrated
20-4170 7½" Serrated
20-3910 10¼" Serrated



HEANEY NEEDLE HOLDER
20-4682 8¼" Serrated, Curved



STRATTE NEEDLE HOLDER
20-4795 9" Serrated, S-Curve



FINOCHIETTO NEEDLE HOLDER
20-4510 10¼" Serrated, Angular

NEEDLE HOLDERS

OB/GYN

General Surgery




MATHIEU NEEDLE HOLDER

20-3055 5½" Serrated
20-3075 7½" Serrated



BAUMGARTNER NEEDLE HOLDER

20-1450 5¼" Serrated



MAYO-HEGAR NEEDLE HOLDER

SERRATED


20-2052 5¼"	20-2480 8"
20-2060 6"	20-2490 9½"
20-2070 7"	20-2910 10½" Masson
20-2270 7" Broad Jaw	20-2911 11½" Masson



CRILE-WOOD NEEDLE HOLDER


SERRATED

20-1265 6" Narrow Jaws
20-1270 7"
20-1280 8"
20-1290 9"
20-1295 10½"



CRILE-MURRAY NEEDLE HOLDER

20-1260 6" Serrated




GILLIES NEEDLE HOLDER

20-2662 6" Serrated

OB/GYN

NEEDLE HOLDERS

General Surgery



**COLLIER
NEEDLE HOLDER**

65-6550 5¼" Serrated,
Fenestrated Jaw




DERF NEEDLE HOLDER

20-1547 4¾" Serrated




**WEBSTER
NEEDLE HOLDER**

20-1350 5¼" Smooth
20-1352 5¼" Serrated
20-3090 5¼" Smooth, Heavy




HALSEY NEEDLE HOLDER

20-1750 5¼" Smooth
20-1751 5¼" Serrated



**NEIVERT
NEEDLE HOLDER**

20-1760 5" 1 Offset Ring



**OLSEN-HEGAR
NEEDLE HOLDER**

20-2747 4¾" Serrated, Delicate
20-2748 4¾" Smooth, Delicate
20-2755 5½" Serrated
20-2765 6½" Serrated
20-2775 7" Serrated



**CONVERSE
NEEDLE HOLDER**

20-3350 5" Serrated,
Longitudinal



BROWN NEEDLE HOLDER



HORIZONTAL SERRATIONS

20-4852 5¼" Rounded Jaws
20-4867 6¾" Rounded Jaws

NEEDLE HOLDERS

Premium

OB/GYN



OB/GYN

NEEDLE HOLDERS

Premium



SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™
DERF NEEDLE HOLDER**
23-2124 4¾" Serrated

SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE
EXTRA DELICATE™
DERF NEEDLE HOLDER**
23-2224 4¾" Serrated



SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™
WEBSTER
NEEDLE HOLDER**
23-2125 5¼" Smooth
23-2120 5¼" Serrated

SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE
EXTRA DELICATE™
WEBSTER
NEEDLE HOLDER**
23-2225 5¼" Smooth
23-2220 5¼" Serrated



SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™
HALSEY
NEEDLE HOLDER**
23-2135 5¼" Smooth
23-2130 5¼" Serrated

SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE
EXTRA DELICATE™
HALSEY
NEEDLE HOLDER**
23-2235 5¼" Smooth
23-2230 5¼" Serrated



SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™
OLSEN-HEGAR
NEEDLE HOLDER**
23-2136 5½" Serrated

SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE™
EXTRA DELICATE
OLSEN-HEGAR
NEEDLE HOLDER**
23-2236 5½" Serrated

NEEDLE HOLDERS

Tungsten Carbide

OB/GYN

Sklar 'S' Left-Handed needle Holder S

Now there's a ring-handled needle holder created for the needs of left-handed surgeons.

Sklar's left-handed needle holder features a unique ratchet design that makes it possible for surgeons to easily control the needle holder with their left hand. As a result surgeons experience less fatigue and hand strain than ever before.

Suitable for Left Hand use



TC WEBSTER NEEDLE HOLDER

SMOOTH

21-9012 5" Left Handed



**TC MAYO-HEGAR
NEEDLE HOLDERS**

SERRATED

21-9016 6" Left Handed

21-9018 7" Left Handed



**TC CRILE-WOOD
NEEDLE HOLDERS**

SERRATED

21-9030 6" Left Handed

21-9035 7" Left Handed



TC RYDER NEEDLE HOLDER


SERRATED

21-9045 6" Left Handed

OB/GYN

NEEDLE HOLDERS


Tungsten Carbide



**TC WANGENSTEEN
NEEDLE HOLDER**

SERRATED

21-3510 10½"



**TC FINOCHIETTO
NEEDLE HOLDER**

SERRATED

21-8039 10½"



**TC NEIVERT
NEEDLE HOLDER**

SERRATED

21-8086 5" One Offset Ring



**TC GILLIES
NEEDLE HOLDER**

SERRATED

21-2662 6"




**TC CASTROVIEJO
NEEDLE HOLDERS**

21-8040 5½" Smooth, Straight
21-8046 5½" Serrated, Straight
21-8047 5½" Serrated, Curved



**TC CASTROVIEJO
NEEDLE HOLDERS**

21-8042 5½" Smooth, Straight
21-8043 5½" Smooth, Curved



**TC CASTROVIEJO
NEEDLE HOLDERS**

21-8044 5½" Smooth, Straight
21-8045 5½" Smooth, Curved




**TC JACOBSON
NEEDLE HOLDERS**

21-8055 8¼" Smooth w/Lock
21-8058 8¼" Smooth w/o Lock
21-8059 9¼" Smooth w/Lock

NEEDLE HOLDERS

OB/GYN

Tungsten Carbide



**TC CRILE-MURRAY
NEEDLE HOLDER**

SERRATED


21-8075 6"



**TC RYDER
NEEDLE HOLDERS**

SERRATED

21-8019 5"
21-8020 6"
21-8021 7"
21-8022 8"
21-8023 9"



**TC RYDER
NEEDLE HOLDERS**

MICRO SERRATED

21-8012 5"
21-8013 6"
21-8014 7"
21-8016 9"



**TC NEW ORLEANS
NEEDLE HOLDER**


SERRATED

21-8082 7"




**TC SAROT
NEEDLE HOLDERS**

21-8036 7 1/8" Serrated
21-8037 10 1/2" Serrated




**TC HEANEY
NEEDLE HOLDERS**

21-8060 8 1/4" Serrated, Curved
21-8061 10" Serrated, Curved



**TC STRATTE
NEEDLE HOLDER**

21-8065 9" Serrated,
S-Curve




**TC STERNAL WIRE
TWISTERS**

21-9060 5"
21-9065 7"

OB/GYN


NEEDLE HOLDERS

Tungsten Carbide



**TC DERF
NEEDLE HOLDER**

SERRATED
21-8004 4¾"



**TC WEBSTER
NEEDLE HOLDERS**

21-8002 5¼" Smooth
21-8080 5¼" Serrated



**TC HALSEY
NEEDLE HOLDERS**

21-8000 5¼" Smooth
21-8001 5¼" Serrated



**TC OLSEN-HEGAR
NEEDLE HOLDERS**

21-8025 4¾" Serrated, Delicate
21-8026 4¾" Smooth, Delicate
21-8027 5½" Serrated
21-8028 6½" Serrated
21-8029 7½" Serrated




**TC MATHIEU
NEEDLE HOLDERS**

21-8070 5½"
21-8071 6¾" Narrow Jaw




**TC BAUMGARTNER
NEEDLE HOLDER**

SERRATED
21-8005 5¼"



**TC MAYO-HEGAR
NEEDLE HOLDERS**

SERRATED
21-8030 5½"
21-8031 6"
21-8032 7"
21-8033 8"
21-8051 9½"
21-8034 10½"
21-8035 12"



**TC CRILE-WOOD
NEEDLE HOLDERS**

SERRATED
21-8006 6" Narrow
21-8007 7" DeBakey
21-8008 8" DeBakey
21-8009 9" DeBakey
21-8010 10½" DeBakey
21-8011 12" DeBakey

NEEDLE HOLDERS

OB/GYN

Tungsten Carbide

needle Holders

Tungsten Carbide Needle Holders Page 108

Sklar TC™ needle holders utilize solid tungsten carbide insert jaws with a Sklar pyramid profile to hold suture needles far more securely than any standard instrument. All TC inserts are virtually indestructible with an extremely hard alloy that provides exceptional control, eliminating needle rotation and slippage. These precisely ground needle holders are readily identified by the dual gold finger rings.

Premium Needle Holders Page 112

Sklar offers Premium instruments that have been modified to meet the needs of clinicians for specific procedures and fully guaranteed for a lifetime of dependable use.

Sklarlite™

Premium Needle Holders have been reduced in their overall pattern weight to produce a lighter instrument for improved comfort that minimizes fatigue to clinicians during prolonged procedures.

Sklarlite X™

Extra Delicate Needle Holders provide the most sensitive surgical tool as well as our lightest possible pattern while retaining its needle grasping capability and strength.

General Surgery Instruments Page 114

Our very best operating room, high quality needle holders are forged and finished from optimum grades of German surgical steel with finely detailed craftsmanship and exacting standards of manufacturing, including our exceptional line of Microsurgery instruments.

Surgi-OR™ Tungsten Carbide Instruments Page 118

A Mid-Grade line of high-quality OR instruments, hand-crafted and forged from Japanese surgical grade stainless steel, that features solid tungsten carbide insert jaws for exceptional control and unbeatable strength. Designed to give you maximum performance at affordable cost, Surgi-OR™ TC needle holders offer significant value to the Operating Room and other important surgical settings. Surgi-OR™ TC needle holders are easily identified by the dual gold finger rings.

Surgi-OR™ Instruments Page 119

A Mid-Grade line of high-quality OR instruments, hand-crafted and forged from Japanese surgical grade stainless steel. This cost-effective line of instruments offers significant value to the Operating Room as well as other important surgical settings.


Econo® Sterile Instruments Page 120

Econo™ Sterile instruments are suitable in a number of surgical venues where economical performance is most important. Convenient, ready-to-use sterile instruments provide the very best value with the greatest utility of any surgical instrument.

OB/GYN


FORCEPS

Econo™ Sterile




TISSUE FORCEPS

96-2579	5"	1x2 Teeth	Case/25
96-2578	5"	1x2 Teeth	Case/50
96-2661	5"	1x2 Teeth	Case/100




FINE POINT SPLINTER FORCEPS

96-2412	4½"	Case/25
96-2410	4½"	Case/50




SPONGE FORCEPS

SERRATED			
96-2533	9"	Straight	Case/25
96-2532	9"	Straight	Case/50




BACKHAUS TOWEL CLAMPS

96-2632	3½"	Straight	Case/25
96-2636	5¼"	Curved	Case/25
96-2630	3½"	Straight	Case/50
96-2634	5¼"	Curved	Case/50



LORNA EDNA TOWEL CLAMPS

96-2529	3½"	Case/25
96-2528	3½"	Case/50



PRESBYTERIAN TUBE OCCLUDING FORCEP

SMOOTH			
96-2408	7"	Straight	Case/25
96-2406	7"	Straight	Case/50


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

FORCEPS

OB/GYN

Econo™ Sterile




ROCHESTER-PEAN FORCEPS

96-2545	6 1/4"	Straight	Case/25
96-2547	6 1/4"	Curved	Case/25
96-3006	7 1/4"	Straight	Case/25
96-3008	7 1/4"	Straight	Case/25
96-2548	8 1/2"	Curved	Case/25
96-2544	6 1/4"	Straight	Case/50
96-2546	6 1/4"	Curved	Case/50
96-2549	8 1/2"	Curved	Case/100
96-2664	8 1/2"	Curved	Case/100




ROCHESTER-OSHSNER FORCEPS

96-2566	6 1/4"	Straight	Case/25
96-2567	6 1/4"	Curved	Case/25
96-2707	6 1/4"	Straight	Case/50
96-2708	6 1/4"	Curved	Case/50




ALLIS TISSUE FORCEPS

96-2575	6"	4x5 Teeth	Case/25
96-2576	6"	4x5 Teeth	Case/50
96-2774	7 1/4"	6x7 Teeth	Case/50



ADSON FORCEPS


96-2571	4 3/4"	Dressing	Serrated	Case/25
96-2573	4 3/4"	Tissue	1x2 Teeth	Case/25
96-2570	4 3/4"	Dressing	Serrated	Case/50
96-2572	4 3/4"	Tissue	1x2 Teeth	Case/50
96-2577	4 3/4"	Dressing	Serrated, Delicate	Case/50
96-2775	4 3/4"	Tissue	1x2 Teeth Delicate	Case/50
96-2655	4 3/4"	Dressing	Serrated	Case/100
96-2656	4 3/4"	Tissue	1x2 Teeth	Case/100



ADSON BROWN FORCEPS

SIDE GRASPING

96-2559	4 3/4"	Case/25
96-2558	4 3/4"	Case/50



DRESSING FORCEPS


SERRATED

96-2569	5"	Case/25
96-2921	8"	Case/25
96-2688	4 1/2"	Case/50
96-2568	5"	Case/50
96-2660	5"	Case/100

OB/GYN


FORCEPS

Econo™ Sterile




BOZEMAN FORCEPS

96-2622	10 1/4"	S-Curve	Case/25
96-2620	10 1/4"	S-Curve	Case/50




SCHROEDER TENACULUM FORCEPS

96-2624	9 1/2"	Case/25
96-2626	9 1/2"	Case/50




HALSTED MOSQUITO FORCEPS

96-2537	5"	Straight	Case/25
96-2539	5"	Curved	Case/25
96-2536	5"	Straight	Case/50
96-2538	5"	Curved	Case/50
96-2657	5"	Straight	Case/100
96-2658	5"	Curved	Case/100




CRILE FORCEPS

96-2551	5 1/2"	Straight	Case/25
96-2553	5 1/2"	Curved	Case/25
96-2550	5 1/2"	Straight	Case/50
96-2552	5 1/2"	Curved	Case/50



KELLY FORCEPS

96-2561	5 1/2"	Straight	Case/25
96-2563	5 1/2"	Curved	Case/25
96-2560	5 1/2"	Straight	Case/50
96-2562	5 1/2"	Curved	Case/50
96-2652	5 1/2"	Straight	Case/100
96-2653	5 1/2"	Curved	Case/100




HARTMANN FORCEPS

96-2556	3 1/2"	Straight	Case/25
96-2557	3 1/2"	Curved	Case/25
96-2554	3 1/2"	Straight	Case/50
96-2555	3 1/2"	Curved	Case/50

FORCEPS


Surgi-OR™

OB/GYN




FOERSTER SPONGE FORCEPS

95-571	7"	Straight	Serrated
95-577	9½"	Straight	Serrated
95-586	9½"	Curved	Serrated



BACKHAUS TOWEL CLAMPS

95-606	3½"
95-611	5¼"




LORNA EDNA TOWEL CLAMPS

95-629	3½"
95-635	5¼"



PRESBYTERIAN TUBE OCCLUDING FORCEP

95-592	7"	Straight, Smooth
--------	----	------------------



MAGILL CATHETER FORCEP

95-285	Child
95-289	Adult

OB/GYN

FORCEPS

Surgi-OR™



ADSON DRESSING FORCEPS

95-771 4¾" Serrated

DELICATE

95-772 4¾" Serrated



ADSON TISSUE FORCEPS

95-775 4¾" 1x2 Teeth

DELICATE

95-776 4¾" 1x2 Teeth



ADSON-BROWN TISSUE FORCEP

95-778 4¾" Side-Grasping



DRESSING FORCEPS

95-705 4½" Serrated

95-710 5" Serrated

95-727 5½" Serrated

95-731 6" Serrated

95-736 8" Serrated

95-742 10" Serrated



TISSUE FORCEPS

95-753 4½" 1x2 Teeth

95-757 5" 1x2 Teeth

95-762 5½" 1x2 Teeth

95-766 6" 1x2 Teeth

95-770 8" 1x2 Teeth



SPLINTER FORCEPS


95-777 3½"

95-781 4½"

FORCEPS

Surgi-OR™

OB/GYN




RANKIN-KELLY FORCEPS

95-465	6¾"	Straight
95-469	6¾"	Curved



HARTMANN MOSQUITO FORCEPS

95-409	3½"	Straight
95-421	3½"	Curved




ROCHESTER-PEAN FORCEPS

STRAIGHT		CURVED	
95-470	5½"	95-471	5½"
95-473	6¼"	95-477	6¼"
95-481	7¼"	95-485	7¼"
95-492	8"	95-498	8"
95-483	9"	95-487	9"
95-494	10¼"		



ROCHESTER-OSHSNER FORCEPS

STRAIGHT 1x2 TEETH		CURVED 1x2 TEETH	
95-500	5½"	95-501	5½"
95-504	6¼"	95-512	6¼"
95-521	7¼"	95-525	7¼"
95-513	8"	95-514	8"
95-515	9"	95-516	9"
95-517	10"	95-518	10"



ALLIS TISSUE FORCEPS

95-934	6"	5x6 Teeth
95-938	7½"	5x6 Teeth

OB/GYN


FORCEPS

Surgi-OR™



BOZEMAN UTERINE DRESSING FORCEP

95-374 10 1/4" S-curved




SCHROEDER TENACULUM FORCEP

95-386 9 1/2"




BABCOCK FORCEP

95-952 6 1/4"
95-441 8"




HALSTED MOSQUITO FORCEPS

95-426 5" Straight
95-431 5" Curved



CRILE FORCEPS

95-447 5 1/2" Straight
95-458 6 1/4" Straight
95-452 5 1/2" Curved
95-463 6 1/4" Curved



KELLY FORCEPS

95-436 5 1/2" Straight
95-442 5 1/2" Curved


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

FORCEPS

OB/GYN


General Surgery



DRESSING FORCEPS

SERRATED, DELICATE

22-1161	4½"
22-1162	5"
22-1163	5½"
22-1164	6"




TISSUE FORCEPS

1 X 2 TEETH

19-1245	4½"	19-1280	8"
19-1250	5"	19-1310	10"
19-1255	5½"	19-1312	12"
19-1260	6"		

1 X 2 TEETH, DELICATE

22-1063	5½"
---------	-----



TISSUE FORCEPS

2 X 3 TEETH


19-1450	5"
19-1455	5½"
19-1460	6"
19-1480	8"
19-1510	10"
19-1512	12"

3 X 4 TEETH

19-1450	5"
19-1455	5½"
19-1460	6"


4 X 5 TEETH

19-1855	5½"
19-1910	10"



FINE POINT SPLINTER FORCEPS

19-3035	3½"
19-3040	4"
19-3045	4½"



FOERSTER SPONGE FORCEPS

STRAIGHT

87-2095	9½"	Smooth
87-2195	9½"	Serrated


CURVED

87-2295	9½"	Smooth
87-2395	9½"	Serrated



BACKHAUS TOWEL CLAMPS

47-2935	3½"
47-2955	5¼"



LORNA EDNA TOWEL CLAMP

NON-PERFORATING

47-2837	4"
47-2852	5¼"

OB/GYN

FORCEPS


General Surgery



ALLIS FORCEPS

ATRAUMATIC JAWS

36-2298	6"	Straight
36-2280	8"	Straight
36-2385	10"	Straight



MAGILL CATHETER FORCEPS

07-1777	6½"	Infant
07-1780	8"	Child
07-1797	10"	Adult



ADSON FORCEPS DRESSING

DRESSING

50-3147	4¾"	Serrated
47-2147	4¾"	Serrated, Delicate

TISSUE

50-3047	4¾"	1x2 Teeth
47-2047	4¾"	1x2 Teeth, Delicate
47-2049	4¾"	2x3 Teeth, Delicate

ADSON-BROWN

47-1447	4¾"	Side Grasping
47-2447	4¾"	Side Grasping, Delicate




ADSON FORCEPS

50-3048	4¾"	Serrated, 1x2 Teeth
50-3049	4¾"	1x2 Teeth with Platform



MICRO ADSON FORCEPS

98-2047	4¾"	1x2 Teeth
98-2048	4¾"	Serrated



DRESSING FORCEPS

SERRATED

19-1040	4"	19-1060	6"
19-1045	4½"	19-1080	8"
19-1050	5"	19-1110	10"
19-1055	5½"	19-1112	12"

FORCEPS

OB/GYN

General Surgery



HARTMANN MOSQUITO FORCEPS

SERRATED

17-1235	3½"	Straight
17-1335	3½"	Curved

DELICATE

17-1230	3½"	Straight
17-1232	3½"	Curved




ROCHESTER-PEAN FORCEPS

STRAIGHT

17-2455	5½"
17-2462	6¼"
17-2467	6¾" Heavy
17-2472	7¼"
17-2480	8"
36-2690	9"
36-2710	10¼"

CURVED

17-2255	5½"
17-2262	6¼"
17-2267	6¾" Heavy
17-2272	7¼"
17-2280	8"
36-2790	9"
36-2810	10¼"



ROCHESTER-PEAN FORCEPS

DELICATE

22-3555	5½"	Straight
22-3655	5½"	Curved
17-2362	6¼"	Curved
22-2272	7¼"	Curved
22-2280	8"	Curved



ROCHESTER-OSHSNER FORCEPS

1x2 TEETH

STRAIGHT

17-2855	5½"
17-2862	6¼"
17-2872	7¼"
17-2880	8"

CURVED

17-2655	5½"
17-2662	6¼"
17-2672	7¼"
17-2680	8"



ROCHESTER-OSHSNER FORCEPS


DELICATE, 1x2 TEETH

STRAIGHT

18-1752	5¼"
22-3755	5½"
17-2962	6¼"

CURVED

18-1852	5¼"
22-3855	5½"
17-2762	6¼"



ALLIS TISSUE FORCEPS

36-2060	6"	3x4 Teeth
36-2160	6"	4x5 Teeth
36-2260	6"	5x6 Teeth
36-2275	7½"	5x6 Teeth

DELICATE

18-2140	4"	4x5 Teeth
18-2155	5"	4x5 Teeth
18-2162	6"	4x5 Teeth
18-2164	6"	5x6 Teeth
36-2295	9½"	5x6 Teeth


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

FORCEPS

General Surgery



Z-TYPE HYSTERECTOMY FORCEPS

LONGITUDINAL GROOVES WITH TIP SERRATIONS

90-2785	8 1/4"	Straight	90-2796	12"	Straight
90-2786	8 1/4"	Slight Curve	90-2797	12"	Slight Curve
90-2787	8 1/4"	Full Curve	90-2798	12"	Full Curve
90-2788	8 1/4"	Angled	90-2799	12"	Angled
90-2790	9 1/2"	Straight	90-2801	14"	Straight
90-2792	9 1/2"	Slight Curve	90-2802	14"	Slight Curve
90-2793	9 1/2"	Full Curve	90-2803	14"	Full Curve
90-2794	9 1/2"	Angled	90-2804	14"	Angled




WERTHEIM-CULLEN CLAMP

33-1375	8 1/2"	Pediatric
---------	--------	-----------



SOMER UTERINE ELEVATING FORCEP

90-5297	9"
---------	----



MATHIEU FOREIGN BODY FORCEP

79-3228	7 3/4"
79-3230	11"




HALSTED MOSQUITO HEMOSTATS

17-1450	5"	Straight
17-1650	5"	Straight, 1x2
17-1550	5"	Curved
17-1750	5"	Curved, 1x2

DELICATE

18-1050	5"	Straight
18-1150	5"	Curved



CRILE HEMOSTATIC FORCEPS

17-3055	5 1/2"	Straight
17-3155	5 1/2"	Curved
17-3062	6 1/4"	Straight
17-3162	6 1/4"	Curved
17-3075	7 1/2"	Straight
17-3175	7 1/2"	Curved

DELICATE

18-1952	5 1/2"	Straight
18-2052	5 1/2"	Curved



KELLY HEMOSTATIC FORCEPS

17-2055	5 1/2"	Straight
17-2155	5 1/2"	Curved
74-3270	7"	Straight, Extra Heavy


DELICATE

18-1552	5 1/2"	Straight
18-1652	5 1/2"	Curved

FORCEPS


OB/GYN

General Surgery




LONG FORCEP

90-2690 7½"



LAWTON FORCEP


90-2180 8¼"



HEANEY-BALLENTINE FORCEPS


LONGITUDINAL SERRATIONS

90-2480	8"	Straight
90-2580	8"	Curved




HEANEY-KANTER FORCEP

90-2379



HEANEY REZAK FORCEP

90-1877 8¼"



HEANEY FORCEPS SINGLE TOOTH

90-1977	7¾"
90-2177	8¼"

Light
Heavy




HEANEY FORCEPS HEAVY DOUBLE TOOTH

90-2280	8"
90-2387	8¾"

OB/GYN


FORCEPS

General Surgery




**HILDEBRANDT UTERINE
HEMOSTATIC FORCEP**

90-8065 6½"




DE ALVAREZ FORCEP

90-2782 7¾" S-Curve



**GUSBERG HYSTERECTOMY
FORCEP**

90-2780 8" Curved



GUTGLASS CERVIX FORCEP

90-8090 9½"



MAIER-GROSS UTERINE DRESSING FORCEPS

WITH RATCHET

90-8070	7"	Straight
90-8072	7"	Curved
90-8080	8"	Straight
90-8082	8"	Curved
90-8110	10"	Straight
90-8112	10"	Curved


WITHOUT RATCHET

90-8120	7"	Straight
90-8122	7"	Curved
90-8130	8"	Straight
90-8132	8"	Curved
90-8140	10"	Straight
90-8142	10"	Curved

FORCEPS


OB/GYN

General Surgery




SCHROEDER VULSELLUM FORCEP

91-1590 9½" Straight




SCHROEDER VULSELLUM FORCEPS

91-1690 9" Curved on Flat
91-1890 9½" Curved Sideways




BILLROTH FORCEP

91-3810 10¾" 4x4 Teeth



LAHEY TRACTION FORCEPS

17-3562 6¼"
91-2080 8" (MGH)



GELPI-LOWRIE FORCEP

90-1885 8½"




PHANEUF FORCEPS

90-2682 8¼" Straight
90-2680 8¼" Curved

OB/GYN

FORCEPS


General Surgery




TEALE VULSELLUM FORCEP
CURVED SIDEWAYS
91-2695 9" 3x4 Teeth




HENROTIN VULSELLUM FORCEP
91-2585 8½" 2x3 Teeth




DOYEN VULSELLUM FORCEP
91-2482 8¼"



JACOBS VULSELLUM FORCEPS
SERRATED
91-2282 8¼" Straight
91-2382 8¼" Angled



MUSEUX VULSELLUM FORCEPS
2X2 TEETH
91-2906 9½" 6mm Straight
91-3006 9½" 6mm Curved
91-2908 9½" 8mm Straight
91-3008 9½" 8mm Curved
91-2910 9½" 10mm Straight
91-3010 9½" 10mm Curved



SECOND (COLLIN) FORCEPS
4X4 TEETH
91-3795 9½" Curved on flat
3X3 TEETH
91-3095 9½" Straight
91-3195 9½" Curved sideways


FORCEPS

OB/GYN


General Surgery




DUPLAY TENACULUM FORCEP
91-1111 11" S-Curve




STAUDE-MOORE TENACULUM FORCEPS
91-1455 5½" Straight
91-1465 6½" Straight




HULKA TENACULUM FORCEP
91-1280 11½"



KAHN TENACULUM FORCEP
91-1125 9½"



ADAIR BREAST TENACULUM FORCEP
36-1865 6½"




BABCOCK FORCEPS
22-6155 5½"
36-1662 6¼"
36-1677 8"
36-1692 9½"

OB/GYN


FORCEPS

General Surgery



SCHWEIZER UTERINE FORCEP


90-7995 9½" Straight



CLEMETSON UTERINE FORCEP


NON-RATCHETED

90-7895 9½" Straight




LAUFE UTERINE POLYP FORCEP

90-7977 7¾" Curved




SCHROEDER TENACULUM FORCEPS

91-1290 10" Braun (Flat Tips)
91-1298 10" Pozzi (Round Tips)



BARRETT TENACULUM FORCEPS

18-2257 6½" Delicate
91-1070 7" Straight




JARCHO TENACULUM FORCEP

91-1275 7½" S-Curved

FORCEPS

OB/GYN


General Surgery



BOZEMAN FORCEP


ONE LARGE RING

90-7710 10 3/4" S-Curve




BOZEMAN FORCEPS

90-7610 10 1/2" Straight
80-1410 10 1/4" S-Curve



CAMPION UTERINE FORCEP

90-7898 9 1/4" Curved




DELEE UTERINE FORCEP

90-8111 11" S-Curve



CHERON UTERINE DRESSING FORCEP

90-7910 10 1/2"



WINTER PLACENTA & OVUM FORCEPS

90-5292 11" Straight
90-5293 11" Curved

OB/GYN

FORCEPS

Premium



SKLAR LITE THE LIGHT TOUCH

SKLAR LITEGRIP™ ADSON FORCEPS

50-3148	4¾"	Serrated
50-3050	4¾"	1x2 Teeth

SKLAR LITE XD EXTRA * DELICATE

SKLAR LITEGRIP™ ADSON- BROWN FORCEPS

47-1448	4¾"	7x7 Teeth Side Grasping
---------	-----	----------------------------



SKLAR LITE THE LIGHT TOUCH

SKLAR LITE™ DRESSING FORCEPS

SERRATED

23-2824	4½"
23-2826	5"
23-2828	5½"
23-2830	6"

SKLAR LITE XD EXTRA * DELICATE

SKLAR LITE™ EXTRA DELICATE DRESSING FORCEPS

23-2834	4½"
23-2836	5"
23-2838	5½"



SKLAR LITE THE LIGHT TOUCH

SKLAR LITE™ TISSUE FORCEPS

23-2809	4½"	1x2 Teeth
23-2811	5½"	1x2 Teeth

SKLAR LITE XD EXTRA * DELICATE

SKLAR LITE™ EXTRA DELICATE TISSUE FORCEPS

23-2819	5½"	1x2 Teeth
---------	-----	-----------

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

FORCEPS

Premium

OB/GYN

SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ ROCHESTER-PEAN FORCEPS

23-2532	6 1/4"	Straight
23-2538	6 1/4"	Curved
23-2542	7 1/4"	Curved

SKLARLITE XD
EXTRA * DELICATE

SKLARLITE™ EXTRA DELICATE ROCHESTER-PEAN FORCEPS

23-2552	6 1/4"	Straight
23-2561	6 1/4"	Curved
23-2563	7 1/4"	Curved



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ ROCHESTER-OGHSNER FORCEPS

23-2560	6 1/4"	1x2 Teeth	Straight
23-2568	6 1/4"	1x2 Teeth	Curved

SKLARLITE XD
EXTRA * DELICATE

SKLARLITE™ EXTRA DELCATE ROCHESTER-OGHSNER FORCEPS

23-2577	6 1/4"	1x2 Teeth	Straight
23-2585	6 1/4"	1x2 Teeth	Curved



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ ALLIS TISSUE FORCEPS

23-2851	6"	4x5 Teeth
23-2852	6"	5x6 Teeth

SKLARLITE XD
EXTRA * DELICATE

SKLARLITE™ EXTRA DELICATE ALLIS TISSUE FORCEPS

23-2861	5"	4x5 Teeth
23-2864	6"	4x5 Teeth
23-2865	6"	5x6 Teeth



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ ADSON FORCEPS

DRESSING		
23-2600	4 3/4"	Serrated

TISSUE		
23-2602	4 3/4"	1x2 Teeth

ADSON-BROWN		
23-2603	4 3/4"	Side Grasping

SKLARLITE XD
EXTRA * DELICATE

SKLARLITE™ EXTRA DELICATE ADSON FORCEPS

DRESSING		
23-2604	4 3/4"	Serrated

TISSUE		
23-2606	4 3/4"	1x2 Teeth


ADSON-BROWN		
23-2608	4 3/4"	Side Grasping



OB/GYN

FORCEPS

Premium




PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

FORCEPS


OB/GYN

Tungsten Carbide




TC BABCOCK FORCEPS

97-1565 6¼"
97-1570 8"




TC ALLIS TISSUE FORCEPS

97-1555 6"
97-1560 7½"




TC ADSON DRESSING FORCEPS

51-3447 4¾" Serrated
51-3049 4¾" 1x2 Teeth with platform



TC ADSON-BROWN FORCEPS

51-1447 4¾" Side Grasping



TC LITEGRIP™ ADSON-BROWN FORCEPS

51-1448 4¾" Side Grasping




TC DRESSING FORCEPS

51-3460 6" Serrated
51-3480 8" Serrated

OB/GYN

FORCEPS

Sklar Black™ Laser



GERALD FORCEPS

90-1025	7"	Straight	Serrated
90-1026	7"	Curved	Serrated
90-1027	7"	Straight	1x2 Teeth



TUTTLE/SINGLEY FORCEPS

90-2270	7"	Tuttle
90-2272	9"	Singley



CASTROVIEJO FORCEP

90-2313	4 1/4"	.5mm Teeth
---------	--------	------------

SKLAR DISPOSABLE SCALPELS



ALL SKLAR PRECISION
PRODUCTS ARE LATEX FREE

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3110	#10	10/box
06-3111	#11	10/box
06-3112	#12	10/box
06-3115	#15	10/box
06-3120	#20	10/box
06-3121	#21	10/box
06-3122	#22	10/box
06-3123	#23	10/box

BOXES ARE MARKED WITH LOT NUMBER
AND EXPIRATION DATES FOR TRACEABILITY

FORCEPS

OB/GYN

Sklar Black™ Laser



ADSON TISSUE FORCEPS

90-2262	4½"	1x2 Teeth
90-1060	4¾"	1x2 Teeth
90-1061	4¾"	2x3 Teeth



ADSON-BROWN TISSUE FORCEP

90-1062	4¾"	Side Grasping	7x7 Teeth
---------	-----	---------------	-----------



DRESSING FORCEPS

SERRATED

90-1063	4½"
90-1064	5"
90-1065	5½"
90-1066	6"



TISSUE FORCEPS

90-1020	4½"	1x2 Teeth
90-1021	5"	1x2 Teeth
90-1022	5½"	1x2 Teeth
90-1023	6"	1x2 Teeth



DEBAKEY ATRAUMATIC TISSUE FORCEPS

2.0 MM JAWS

90-2310	7¾"
90-2312	9½"




POTTS-SMITH FORCEPS

90-2304	7"	Serrated
90-2306	10"	Serrated

OB/GYN

FORCEPS

Sklar Black™ Laser



SCHROEDER TENACULUM FORCEPS

90-2325 10" Braun (Flat Tips)
90-2326 10" Pozzi (Round Tips)



SCHROEDER VULSELLUM FORCEPS

90-1790 9½" Straight
90-1890 9½" Curved



BABCOCK FORCEPS

90-2278 8"
90-2295 9½"




ALLIS TISSUE FORCEPS

90-2276 7½" 5x6 Teeth
90-2279 9½" 5x8 Teeth



ROCHESTER-PEAN FORCEP

90-2292 7¼" Curved



ADSON FORCEP

90-2261 4¾" Serrated

FORCEPS

OB/GYN

Sklar Black™ Laser



BOZEMAN FORCEPS

ONE LARGE RING

90-2342 10 1/4" S-Curve



BOZEMAN FORCEPS

90-2340 10 1/4" S-Curve



CAMPION GRASPING FORCEPS

90-2298 9 1/4" Curved



DELEE UTERINE FORCEP

90-7111 11" Heavy



SCHWEIZER FORCEPS

90-2350 9 1/2"



CLEMETSON FORCEPS

NOT RATCHETED

90-2348 9 1/2" Straight



LAUFE UTERINE POLYP FORCEPS

90-2345 7 3/4" Curved

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS

Econo™ Sterile

OB/GYN



UTILITY SCISSORS

96-2701	5½"	Black	25/cs	96-1223	7"	Burgundy	25/cs	96-1228	7"	Lime Green	25/cs
96-2703	7"	Black	25/cs	96-1224	7"	Magenta	25/cs	96-1229	7"	Hunter Green	25/cs
96-1220	7"	Yellow	25/cs	96-1225	7"	Navy Blue	25/cs	96-2702	5½"	Black	50/cs
96-1221	7"	Orange	25/cs	96-1226	7"	Royal	25/cs	96-2704	7"	Black	50/cs
96-1222	7"	Purple	25/cs	96-1227	7"	Teal Blue	25/cs				

STETHOSCOPES

SKLAR OFFERS A COMPLETE
LINE OF HIGH QUALITY
STETHOSCOPES


For a complete listing of Sklar's extensive
Stethoscope product line, contact your
Sklar Account Manager or Distributor
Sales Representative.



OB/GYN

SCISSORS

Econo™ Sterile




STRABISMUS SCISSORS

96-2542	4½" Straight	25/cs
96-2543	4½" Curved	25/cs
96-2540	4½" Straight	50/cs
96-2541	4½" Curved	50/cs




SPENCER STITCH SCISSORS

96-2598	4½"	25/cs
96-2672	3½"	50/cs
96-2596	4½"	50/cs




LITTAUER STITCH SCISSORS

96-2511	3½"	25/cs
96-2513	5½"	25/cs
96-2510	3½"	50/cs
96-2512	5½"	50/cs
96-2663	3½"	100/cs
96-2662	5½"	100/cs




WIRE CUTTING SCISSORS

96-2519	4¾" Angled	25/cs
96-2518	4¾" Angled	50/cs



LISTER BANDAGE SCISSORS

96-2503	5½"	25/cs
96-2501	4½"	50/cs
96-2502	5½"	50/cs
96-2654	5½"	100/cs



KNOWLES BANDAGE SCISSORS

96-2404	5½"	25/cs
96-2402	5½"	50/cs


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS

OB/GYN

Econo™ Sterile



OPERATING SCISSORS

STRAIGHT


96-2597	4½"	Sharp/Blunt	25/cs	96-2514	5½"	Blunt/Blunt	50/cs
96-2592	4½"	Sharp/Sharp	25/cs	96-2520	5½"	Sharp/Blunt	50/cs
96-2521	5½"	Sharp/Blunt	25/cs	96-2522	5½"	Sharp/Sharp	50/cs
96-2523	5½"	Sharp/Sharp	25/cs	96-2726	6½"	Sharp/Blunt	50/cs
96-2728	4½"	Sharp/Blunt	50/cs	96-2649	5½"	Sharp/Blunt	100/cs
96-2593	4½"	Sharp/Sharp	50/cs				



OPERATING SCISSORS


CURVED

96-2526	5½"	Sharp/Blunt	25/cs
96-2527	5½"	Sharp/Sharp	25/cs
96-2524	5½"	Sharp/Blunt	50/cs
96-2525	5½"	Sharp/Sharp	50/cs



MAYO SCISSORS

96-2531	5½"	Straight	25/cs
96-2530	5½"	Straight	50/cs



METZENBAUM SCISSORS

96-2517	5½"	Straight	25/cs
96-2720	5½"	Curved	25/cs
96-2682	7"	Curved	25/cs
96-2516	5½"	Straight	50/cs
96-2721	5½"	Curved	50/cs
96-2680	7"	Curved	50/cs



IRIS SCISSORS

96-2505	4½"	Straight	25/cs
96-2507	4½"	Curved	25/cs
96-2497	3½"	Straight	50/cs
96-2499	3½"	Curved	50/cs
96-2504	4½"	Straight	50/cs
96-2506	4½"	Curved	50/cs
96-2650	4½"	Straight	100/cs
96-2651	4½"	Curved	100/cs


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN


SCISSORS

SURGI-OR™



STRABISMUS SCISSORS

95-360 4½" Straight
95-361 4½" Curved




SPENCER STITCH SCISSORS

95-346 3½"
95-352 4½"



LITTAUER STITCH SCISSOR

95-358 5½"



LISTER BANDAGE SCISSORS


95-222 3½"
95-226 4½"
95-231 5½"
95-241 7¼"
95-243 8"
95-245 8" Large Ring



LISTER BANDAGE SCISSORS


WITH POCKET CLIP

95-232 5½"
95-232B 5½" Blue Ring Handle
95-232P 5½" Pink Ring Handle



KNOWLES BANDAGE SCISSOR


95-254 5½"



MULTI-CUT UTILITY SHEARS

SERRATED

95-100 6" Black Handle
95-105 7½" Black Handle



WIRE CUTTING SCISSOR

SERRATED

95-364 4¾" Angled

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS

OB/GYN


SURGI-OR™



OPERATING SCISSORS

STRAIGHT


95-265	4½"	Sharp/Blunt
95-266	4½"	Sharp/Sharp
95-290	5"	Sharp/Sharp
95-278	5½"	Blunt/Blunt
95-272	5½"	Sharp/Blunt
95-275	5½"	Sharp/Sharp
95-281	6½"	Sharp/Blunt
95-284	6½"	Sharp/Sharp



OPERATING SCISSORS


CURVED

95-294	4½"	Sharp/Sharp
95-292	5½"	Sharp/Blunt
95-295	5½"	Sharp/Sharp
95-299	5½"	Blunt/Blunt



IRIS SCISSORS

95-103	3½"	Straight
95-107	3½"	Curved
95-104	4½"	Straight
95-109	4½"	Curved




MAYO DISSECTING SCISSORS

95-308	5½"	Straight	
95-330	5½"	Curved	
95-321	6¾"	Straight	
95-325	6¾"	Curved	
95-329	6¾"	Straight	Serrated
95-332	6¾"	Curved	Serrated
95-327	9"	Straight	
95-328	9"	Curved	



METZENBAUM DISSECTING SCISSORS

95-335	5¾"	Straight	Lahey
95-336	5¾"	Curved	Lahey
95-339	7"	Straight	
95-340	7"	Curved	



STEVENS TENOTOMY SCISSORS

95-160	4½"	Straight	Blunt/Blunt
95-161	4½"	Curved	Blunt/Blunt

OB/GYN

SCISSORS

SURGI-OR™ Tungsten Carbide



TC OPERATING SCISSORS

95-110	5½"	Straight	Sharp/Blunt
95-111	5½"	Curved	Sharp/Blunt



TC IRIS SCISSORS

95-124	4½"	Straight	
95-125	4½"	Curved	



TC MAYO DISSECTING SCISSORS

95-112	6¾"	Straight	
95-113	6¾"	Curved	



TC MAYO METZENBAUM SCISSORS

95-119	5¾"	Straight	Standard
95-122	5¾"	Curved	Standard
95-115	7"	Straight	Standard
95-117	7"	Curved	Standard
95-118	5¾"	Straight	Delicate
95-120	5¾"	Curved	Delicate
95-114	7"	Straight	Delicate
95-116	7"	Curved	Delicate

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS

OB/GYN


General Surgery



UTILITY SCISSORS

SERRATED

11-1278	6"	Blue
11-1280	7½"	Blue
11-1281	7½"	Black
11-1282	7½"	Green
11-1283	7½"	Red
11-1284	7½"	White
11-1285	7½"	Purple
11-1286	7½"	Orange
11-1287	7½"	Yellow



BRAUN EPISIOTOMY SCISSORS

ANGLED


90-1755	5½"
90-1785	8½"



WALDMANN EPISIOT SCISSOR


SERRATED

90-1770	7"
---------	----




BUSCH UMBILICAL SCISSOR

92-5565	6½"
---------	-----




UMBILICAL SCISSOR

92-5540	4"
---------	----



SCHUMACHER UMBILICAL SCISSOR

92-5560	6¼"
---------	-----



BROPHY SULLIVAN SCISSORS


SERRATED

22-2855	5¾"	Straight
22-2755	5¾"	Curved

OB/GYN


SCISSORS

General Surgery




ESMARCH HEAVY PLASTER SHEARS

11-1380 8"
11-1390 9"




SKLAR HIGH LEVEL BANDAGE SCISSOR

11-1252 5½" Serrated




KNOWLES BANDAGE SCISSOR

11-1255 5½"




MOLESKIN SCISSORS

97-1350 7½" Straight Serrated
97-1351 7½" Curved Serrated



USA GAUZE SCISSOR

11-1180 8"



CROWN SCISSORS

24-2345 4" Straight Smooth
24-2340 4" Straight Serrated
24-2346 4" Curved Smooth
24-2341 4" Curved Serrated

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS

OB/GYN

General Surgery



LISTER BANDAGE SCISSORS


11-1035 3½"
11-1045 4½"
11-1055 5½"
11-1072 7¼"
11-1080 8"
11-1060 8" Large Ring



LISTER BANDAGE SCISSORS


DELICATE

22-1054 5½"
22-1072 7¼"



NURSES BANDAGE SCISSORS

11-1495 5½" Pink
11-1496 5½" Blue



BANTAM BANDAGE SCISSOR

11-1257 5¼" Serrated



COMMANDER UTILITY SHEAR

11-1298 8½" Serrated



HERCULES UTILITY SCISSOR

11-1297 7½" Serrated

OB/GYN


SCISSORS

General Surgery



SUTURE SCISSOR

24-2347 4¾" Straight Serrated Heavy



SUTURE WIRE CUTTING SCISSOR

24-2342 4¾" Angled Serrated



SUTURE WIRE SCISSOR

24-1965 4½" Delicate




SUTURE SCISSORS

22-2935 3½" Shortbent
22-2934 3½" Shortbent Delicate
22-2947 4¾" Northbent



SKLAR STITCH SCISSORS

15-3150 4½" Angled
22-2950 5" Curved Sideways



SKLAR PRECISION SUTURE SCISSORS

22-2820 3½" Straight
22-2822 3½" Curved
22-2824 4½" Straight
22-2825 4½" Curved

SCISSORS

OB/GYN

General Surgery



RIBBON IRIS SCISSORS

47-1250 4" Straight
47-1255 4" Curved



RIBBON STRABISMUS SCISSORS

64-1240 4" Straight
64-1340 4" Curved




RIBBON STEVENS TENOTOMY SCISSORS

64-1241 4" Straight
64-1341 4" Curved




RIBBON STITCH SCISSORS

64-1345 4" Straight
64-1445 4" Curved



UTILITY RIBBON SCISSORS

47-1146 4" Straight Blunt/Blunt
47-1247 4" Curved Blunt/Blunt



SPENCER STITCH SCISSORS

22-2835 3½" Spencer
22-2834 3½" Spencer Delicate
22-2845 4½" Spencer
22-2844 4½" Spencer Delicate
15-1355 5½" Littauer
22-1354 5½" Littauer Delicate

OB/GYN

SCISSORS

General Surgery



IRIS SCISSORS

STRAIGHT

47-1135 3½" Sharp/Sharp
47-1140 3½" Blunt/Blunt
47-1441 4" Sharp/Sharp Heavy
47-1145 4½" Sharp/Sharp
47-1144 4½" Blunt/Blunt

CURVED


47-1235 3½" Sharp/Sharp
47-1240 3½" Blunt/Blunt
47-1442 4" Sharp/Sharp Heavy
47-1245 4½" Sharp/Sharp
47-1244 4½" Blunt/Blunt



IRIS SCISSORS


DELICATE

22-1135 3½" Straight
22-1235 3½" Curved
47-1440 4" Straight
47-1246 4½" Straight
47-1248 4½" Curved




IRIS SCISSORS

64-2042 4¼" Angled Probe Point
64-2142 4½" Angled Sharp/Sharp



STEVENS TENOTOMY SCISSORS


64-1242 4½" Straight Blunt/Blunt
64-1343 4½" Curved Sharp/Sharp
64-1342 4½" Curved Blunt/Blunt



STEVENS TENOTOMY SCISSORS

DELICATE

22-1240 4½" Straight Sharp/Sharp
22-1241 4½" Curved Sharp/Sharp
22-1242 4½" Straight Blunt/Blunt
22-1342 4½" Curved Blunt/Blunt



STRABISMUS SCISSORS

64-1442 4½" Straight
64-1542 4½" Curved


DELICATE

22-1442 4½" Straight
22-1542 4½" Curved

SCISSORS

OB/GYN

General Surgery




KELLY-ADSON SCISSORS

50-4562 6 1/4" Straight Smooth
50-4563 6 1/4" Straight Serrated




KELLY-ADSON SCISSOR

50-4662 6 1/4" Curved Smooth



KELLY-ADSON SCISSOR

50-4462 6 1/4" Angled Smooth




SIMS SCISSORS

90-1280 8" Straight Sharp/Sharp
90-1180 8" Straight Sharp/Blunt
90-1080 8" Straight Blunt/Blunt
90-1580 8" Curved Sharp/Sharp
90-1480 8" Curved Sharp/Blunt
90-1380 8" Curved Blunt/Blunt



PLASTIC SURGERY SCISSOR

47-1345 4 1/2" Curved Sharp/Sharp



DURA SCISSOR

22-2742 4 1/4" Curved Semi-Sharp

OB/GYN

SCISSORS

General Surgery



METZENBAUM SCISSORS

STRAIGHT


75-5570	7"	
55-9180	8"	Nelson
55-9190	9"	Nelson
55-9310	10"	Nelson
55-9311	11"	Nelson



METZENBAUM SCISSORS

CURVED

75-5670	7"	
55-9080	8"	Nelson
55-9090	9"	Nelson
55-9210	10"	Nelson
55-9411	11"	Nelson
55-9612	12"	Nelson
55-9614	14"	Nelson



METZENBAUM SCISSORS


DELICATE

22-1071	7"	Straight	
22-1070	7"	Curved	
55-9184	8"	Straight	Nelson
55-9084	8"	Curved	Nelson
55-9194	9"	Straight	Nelson
55-9094	9"	Curved	Nelson



METZENBAUM (LAHEY) SCISSORS


22-2550	4½"	Straight
22-2545	4½"	Curved
75-5460	5¾"	Straight
47-1360	5¾"	Curved



METZENBAUM (LAHEY) SCISSORS

DELICATE

22-2450	4½"	Straight
22-2650	4½"	Curved
22-2460	5¾"	Straight
22-1057	5¾"	Curved



DOYEN SCISSORS

15-1307	7"	Straight Smooth
15-1308	7"	Straight Serrated
15-2307	7"	Curved Smooth
15-2308	7"	Curved Serrated

SCISSORS

OB/GYN

General Surgery




MAYO DISSECTING SCISSORS

STRAIGHT

15-1555	5½"	
15-1567	6¾"	
90-1590	9"	
55-8711	11"	Harrington

CURVED


15-2555	5½"	
15-2567	6¾"	
90-1690	9"	
55-8611	11"	Harrington



MAYO DISSECTING SCISSORS


DELICATE

22-1556	5½"	Straight
22-2555	5½"	Curved
22-1567	6¾"	Straight
22-2567	6¾"	Curved



MAYO-STILLE DISSECTING SCISSORS

15-1655	5½"	Straight
15-2655	5½"	Curved
15-1667	6¾"	Straight
15-2667	6¾"	Curved




MAYO NOBLE DISSECTING SCISSORS

15-1762	6½"	Straight
15-2762	6½"	Curved



MAYO-NOBLE STILLE DISSECTING SCISSORS

15-1862	6"	Straight
15-2862	6"	Curved



ADA DISSECTING SCISSORS

52-3170	7"	Curved
52-3080	8"	Straight

OB/GYN

SCISSORS

General Surgery



OPERATING SCISSORS

STRAIGHT

13-1045	4½"	Sharp/Sharp
14-1045	4½"	Sharp/Blunt
15-1045	4½"	Blunt/Blunt
13-1050	5"	Sharp/Sharp
14-1050	5"	Sharp/Blunt
15-1050	5"	Blunt/Blunt
13-1055	5½"	Sharp/Sharp
14-1055	5½"	Sharp/Blunt
15-1055	5½"	Blunt/Blunt
13-1060	6"	Sharp/Sharp
14-1060	6"	Sharp/Blunt
15-1060	6"	Blunt/Blunt
13-1065	6½"	Sharp/Sharp
14-1065	6½"	Sharp/Blunt
15-1065	6½"	Blunt/Blunt
13-1072	7¼"	Sharp/Sharp
14-1072	7¼"	Sharp/Blunt

DELICATE

22-1145	4½"	Sharp/Sharp
22-1245	4½"	Sharp/Blunt
22-1345	4½"	Blunt/Blunt
22-1150	5"	Sharp/Sharp
22-1250	5"	Sharp/Blunt
22-1350	5"	Blunt/Blunt
22-1155	5½"	Sharp/Sharp
22-1255	5½"	Sharp/Blunt
22-1355	5½"	Blunt/Blunt



OPERATING SCISSORS

CURVED

13-2045	4½"	Sharp/Sharp
14-2045	4½"	Sharp/Blunt
15-2045	4½"	Blunt/Blunt
13-2050	5"	Sharp/Sharp
14-2050	5"	Sharp/Blunt
15-2050	5"	Blunt/Blunt
13-2055	5½"	Sharp/Sharp
14-2055	5½"	Sharp/Blunt
15-2055	5½"	Blunt/Blunt
13-2060	6"	Sharp/Sharp
14-2060	6"	Sharp/Blunt
15-2060	6"	Blunt/Blunt
13-2065	6½"	Sharp/Sharp
14-2065	6½"	Sharp/Blunt
15-2065	6½"	Blunt/Blunt
13-2072	7¼"	Sharp/Sharp
14-2072	7¼"	Sharp/Blunt

DELICATE

22-1445	4½"	Sharp/Sharp
22-1545	4½"	Sharp/Blunt
22-1645	4½"	Blunt/Blunt
22-1450	5"	Sharp/Sharp
22-1550	5"	Sharp/Blunt
22-1650	5"	Blunt/Blunt
22-1455	5½"	Sharp/Sharp
22-1555	5½"	Sharp/Blunt
22-1655	5½"	Blunt/Blunt

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS

Premium

OB/GYN

SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ SPENCER STITCH SCISSORS

23-1188 3½"
23-1189 4½"
23-1190 5½" Littauer



SKLARLITE XD
EXTRA * DELICATE

SKLARLITE XD™ SPENCER STITCH SCISSORS

23-1193 3½"
23-1194 4½"
23-1195 5½" Littauer

SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ SHORTBENT SUTURE SCISSORS

23-1191 3½"



SKLARLITE XD
EXTRA * DELICATE

SKLARLITE XD™ SHORTBENT SUTURE SCISSORS

23-1196 3½"

SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ PRECISION SUTURE SCISSORS

23-1290 3½" Straight
23-1292 3½" Curved
23-1294 4½" Straight
23-1296 4½" Curved



SKLARLITE XD
EXTRA * DELICATE

SKLARLITE XD™ PRECISION SUTURE SCISSORS

23-1300 3½" Straight
23-1302 3½" Curved
23-1304 4½" Straight
23-1306 4½" Curved

SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ LISTER BANDAGE SCISSORS

23-1199 5½"
23-1200 7¼"



SKLARLITE XD
EXTRA * DELICATE

SKLARLITE XD™ LISTER BANDAGE SCISSORS

23-1202 5½"
23-1203 7¼"

OB/GYN

SCISSORS

Premium

SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™ METZENBAUM
SCISSORS**

23-1180	7"	Straight
23-1182	7"	Curved



SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE XD™ METZENBAUM
SCISSORS**

23-1184	7"	Straight
23-1186	7"	Curved



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ IRIS SCISSORS

23-1148	3½"	Straight
23-1150	3½"	Curved
23-1149	4½"	Straight
23-1151	4½"	Curved



SKLARLITE XD
EXTRA * DELICATE

SKLARLITE™ XD IRIS SCISSORS

23-1152	3½"	Straight
23-1154	3½"	Curved
23-1153	4½"	Straight
23-1155	4½"	Curved



SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™ STEVENS
SCISSORS**

23-1156	4½"	Straight	Blunt/Blunt
23-1157	4½"	Curved	Blunt/Blunt



SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE XD™ STEVENS
SCISSORS**


23-1160	4½"	Straight	Blunt/Blunt
23-1161	4½"	Curved	Blunt/Blunt



SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™ STRABISMUS
SCISSORS**

23-1158	4½"	Straight
23-1159	4½"	Curved



SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE XD™ STRABISMUS
SCISSORS**

23-1162	4½"	Straight
23-1163	4½"	Curved



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS

Premium

OB/GYN

SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ OPERATING SCISSORS

23-1100	4½"	Straight	Sharp/Sharp
23-1104	4½"	Straight	Sharp/Blunt
23-1101	5"	Straight	Sharp/Sharp
23-1105	5"	Straight	Sharp/Blunt
23-1102	5½"	Straight	Sharp/Sharp
23-1106	5½"	Straight	Sharp/Blunt

SKLARLITE XD
EXTRA * DELICATE

SKLARLITE XD™ OPERATING SCISSORS

23-1112	4½"	Straight	Sharp/Sharp
23-1116	4½"	Straight	Sharp/Blunt
23-1113	5"	Straight	Sharp/Sharp
23-1117	5"	Straight	Sharp/Blunt
23-1114	5½"	Straight	Sharp/Sharp
23-1118	5½"	Straight	Sharp/Blunt



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ OPERATING SCISSORS

23-1124	4½"	Curved	Sharp/Sharp
23-1128	4½"	Curved	Sharp/Blunt
23-1126	5½"	Curved	Sharp/Sharp
23-1130	5½"	Curved	Sharp/Blunt

SKLARLITE XD
EXTRA * DELICATE

SKLARLITE XD™ OPERATING SCISSORS

23-1136	4½"	Curved	Sharp/Sharp
23-1140	4½"	Curved	Sharp/Blunt
23-1138	5½"	Curved	Sharp/Sharp
23-1142	5½"	Curved	Sharp/Blunt



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ MAYO DISSECTING SCISSORS

23-1164	5½"	Straight
23-1166	5½"	Curved
23-1165	6¾"	Straight
23-1167	6¾"	Curved

SKLARLITE XD
EXTRA * DELICATE

SKLARLITE XD™ MAYO DISSECTING SCISSORS

23-1168	5½"	Straight
23-1170	5½"	Curved
23-1169	6¾"	Straight
23-1171	6¾"	Curved



SKLARLITE
THE LIGHT TOUCH


SKLARLITE™ METZENBAUM (LAHEY) SCISSORS

23-1172	4½"	Straight
23-1174	4½"	Curved
23-1173	5¾"	Straight
23-1175	5¾"	Curved

SKLARLITE XD
EXTRA * DELICATE

SKLARLITE XD™ METZENBAUM (LAHEY) SCISSORS


23-1176	4½"	Straight
23-1178	4½"	Curved
23-1177	5¾"	Straight
23-1179	5¾"	Curved



OB/GYN

SCISSORS

Premium




SKLARHONE™ OPERATING SCISSORS

STRAIGHT

15-3342	5½"	Smooth	Sharp/Sharp
15-3343	5½"	Smooth	Sharp/Blunt
15-3442	5½"	Serrated	Sharp/Sharp
15-3443	5½"	Serrated	Sharp/Blunt


CURVED

15-3348	5½"	Smooth	Sharp/Sharp
15-3349	5½"	Smooth	Sharp/Blunt
15-3448	5½"	Serrated	Sharp/Sharp
15-3449	5½"	Serrated	Sharp/Blunt




SKLARHONE™ MAYO SCISSORS

15-3325	5½"	Straight	Smooth
15-3330	5½"	Curved	Smooth
15-3430	5½"	Curved	Serrated
15-3326	6¾"	Straight	Smooth
15-3426	6¾"	Straight	Serrated
15-3331	6¾"	Curved	Smooth




SKLARHONE™ IRIS SCISSORS

15-3352	4½"	Straight	Smooth
15-3452	4½"	Straight	Serrated
15-3357	4½"	Curved	Smooth
15-3457	4½"	Curved	Serrated




SKLARHONE™ METZENBAUM (LAHEY) SCISSORS

15-3306	5¾"	Curved	Smooth
15-3406	5¾"	Curved	Serrated
15-3310	7"	Curved	Smooth
15-3410	7"	Curved	Serrated
15-3311	8"	Curved	Smooth
15-3315	9"	Curved	Smooth
15-3412	9"	Curved	Serrated



SKLARHONE™ SKLARLITE™ METZENBAUM (LAHEY) SCISSORS

15-9107	7"	Straight	
15-9117	7"	Curved	
15-9108	8"	Straight	
15-9118	8"	Curved	
15-9109	9"	Straight	
15-9119	9"	Curved	



SKLARHONE™ PLASTIC SURGERY SCISSORS

15-3460	4¾"	Straight	
15-3461	4¾"	Curved	


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS


Premium

OB/GYN



SKLARCUT™ RAGNELL SCISSORS

15-3320 5" Curved




SKLARLITE™ SKLARCUT™ RAGNELL SCISSORS

23-1250 5" Curved




SKLARCUT™ FOMON LATERAL SCISSORS

15-3548 5" Curved




SKLARCUT™ PLASTIC SURGERY SCISSORS

15-3520 4 1/2" Straight
15-3525 4 1/2" Curved
15-3535 4 1/2" Curved
Delicate




SKLARCUT™ GORNEY PLASTIC SCISSORS

15-3577 7 1/4" Straight
15-3578 7 1/4" Curved



SKLARCUT™ GORNEY FACE LIFT SCISSORS

15-3583 7 1/2" Straight
15-3584 7 1/2" Curved



SKLARCUT™ KAYE FACE LIFT SCISSORS

15-3579 6" Curved
15-3580 7 1/2" Curved

SKLARLITE
THE LIGHT TOUCH
SKLARCUT

OB/GYN


SCISSORS

Premium




SKLARLITE™ SKLARCUT™ MAYO SCISSORS

23-1274 6¾" Straight
23-1278 6¾" Curved




SKLARCUT™ METZENBAUM (LAHEY) SCISSORS

15-3569 4½" Straight
15-3571 4½" Curved
15-3570 5¾" Curved
15-3575 7" Curved
15-3574 9" Curved




SKLARLITE™ SKLARCUT™ METZENBAUM (LAHEY) SCISSORS

23-1280 5¾" Curved
23-1267 7" Curved




SKLARCUT™ GOLDMAN-FOX SCISSOR

15-3536 5" Delicate



SKLARCUT™ WAGNER SCISSOR

15-3542 4½" Straight Sharp/Blunt



SKLARLITE™ SKLARCUT™ WAGNER SCISSOR

23-1254 4½" Straight Sharp/Blunt

SCISSORS

Premium

OB/GYN



SKLARCUT™ OPERATING SCISSORS

15-3500	5½"	Straight	Sharp/Sharp
15-3505	5½"	Straight	Sharp/Blunt
15-3510	5½"	Curved	Sharp/Sharp
15-3515	5½"	Curved	Sharp/Blunt



SKLARCUT™ OPERATING SCISSORS

23-1260	5½"	Straight	Sharp/Sharp
23-1262	5½"	Straight	Sharp/Blunt
23-1264	5½"	Curved	Sharp/Sharp
23-1266	5½"	Curved	Sharp/Blunt



SKLARCUT™ IRIS SCISSORS


15-3585	4½"	Straight
15-3590	4½"	Curved



SKLARCUT™ IRIS SCISSORS

23-1268	4½"	Straight
23-1270	4½"	Curved

SKLARLITE
THE LIGHT TOUCH
SKLARCUT™
THE LIGHT TOUCH



SKLARCUT™ STEVENS TENOTOMY SCISSORS

15-3372	4½"	Straight	Blunt/Blunt
15-3370	4½"	Curved	Blunt/Blunt




SKLARCUT™ MAYO SCISSORS

15-3550	5½"	Straight
15-3560	5½"	Curved
15-3555	6¾"	Straight
15-3565	6¾"	Curved
15-3566	9"	Straight
15-3568	9"	Curved

OB/GYN

SCISSORS


Sklar Edge™ Tungsten Carbide



SKLAR EDGE™ TC METZENBAUM SCISSORS

STRAIGHT

16-2010	7"	Standard	
16-2210	9"	Standard	
16-2005	7"	Delicate	
16-2205	9"	Delicate	
16-2015	7"	Delicate	Sharp/Sharp



SKLAR EDGE™ TC METZENBAUM SCISSORS

CURVED

16-2110	7"	Standard	
16-2310	9"	Standard	
16-2311	11"	Standard	
16-2105	7"	Delicate	
16-2305	9"	Delicate	
16-2115	7"	Delicate	Sharp/Sharp




SKLAR EDGE™ TC STRABISMUS SCISSORS

16-2605	4½"	Straight
16-2705	4½"	Curved



SKLAR EDGE™ TC DISSECTING SCISSOR

16-2610	5"	Curved
---------	----	--------



SKLAR EDGE™ TC WIRE CUTTING SCISSOR

16-2800	4¾"	Angled	Serrated
---------	-----	--------	----------



SKLAR EDGE™ TC BLEPHAROPLASTY SCISSOR

47-1148	4½"	Curved
---------	-----	--------

SCISSORS

OB/GYN

Sklar Edge™ Tungsten Carbide



SKLAR EDGE™ TC OPERATING SCISSORS

STRAIGHT

16-1055	5½"	Sharp/Sharp
16-1056	5½"	Sharp/Blunt
16-1057	5½"	Blunt/Blunt



SKLAR EDGE™ TC OPERATING SCISSORS

CURVED

16-1155	5½"	Sharp/Sharp
16-1156	5½"	Sharp/Blunt
16-1157	5½"	Blunt/Blunt



SKLAR EDGE™ TC IRIS SCISSORS

16-2405	4½"	Straight
16-2505	4½"	Curved




SKLAR EDGE™ TC MAYO DISSECTING SCISSORS

16-1505	5½"	Straight
16-1605	5½"	Curved
16-1510	6¾"	Straight
16-1610	6¾"	Curved
16-1515	9"	Straight
16-1615	9"	Curved



SKLAR EDGE™ TC MAYO-STILLE DISSECTING SCISSORS

16-3552	5½"	Straight
16-3562	5½"	Curved
16-3557	6¾"	Straight
16-3567	6¾"	Curved



SKLAR EDGE™ TC METZENBAUM (LAHEY) SCISSORS

16-1805	5¾"	Straight	Delicate
16-1905	5¾"	Curved	Delicate
16-1810	5¾"	Straight	Standard
16-1910	5¾"	Curved	Standard

SKLAR EDGE
TUNGSTEN CARBIDE

OB/GYN

SCISSORS

Scissors

SklarEdge™ Page 53

SklarEdge™ scissors are manufactured using a Sklar Exclusive edging process that permanently bonds Tungsten Carbide to the scissor blades. Tungsten Carbide is an extremely hard alloy that provides an exceptionally strong and sharper edge while maintaining a Sklar 5-year guarantee of durability. These perfectly balanced scissors are readily identified by the dual gold finger rings.

Premium Page 55

Sklar is proud to offer a full line of Premium instruments that have been modified to meet the needs of clinicians for specific procedures. Choose from a wide variety of the most popular patterns.

SklarCut™

SklarCut™ improves the design and function of your most frequently used surgical instruments with the addition of ultra-fine serrations on one blade. This premium feature provides additional surgical control by holding the tissue while cutting more cleanly and with greater precision. SklarCut™ scissors are easily identified by the gold finger ring on the serrated blade handle.

SklarHone™

SklarHone™ scissors combine the cutting edge of a surgical knife with that of a standard scissor. With either smooth or a finely serrated edge, these superior instruments employ the most sophisticated technology to create a smooth, clean cutting action; ideal for delicate surgery or when thick coarse tissue is encountered. SklarHone™ scissors are distinctly identified by their dual black handles.

SklarLite™/SklarLite XD™

SklarLite™ premium scissors provide the clinician a lighter pattern, minimizing fatigue during prolonged procedures. A comfortable instrument that withstands the rigors of frequent use while retaining its structural integrity. SklarLite XD™ extra delicate scissors provide the most sensitive surgical tool as well as our lightest pattern with surprising strength. SklarLite™ and SklarLite XD™ are fully guaranteed for a lifetime of dependable use.

General Surgery Instruments Page 62

Our very best operating room, high quality instruments incorporating optimum grades of German grade surgical steel, finely detailed craftsmanship and exacting standards of manufacturing. Choose from a full line of the most popular patterns.

Surgi-OR™ Instruments Page 72

A Mid-Grade line of high-quality OR instruments, hand-crafted and forged from Japanese surgical grade stainless steel. This cost-effective line of instruments offers significant value to the Operating Room as well as other important surgical settings.

Econo™ Sterile Instruments Page 75

Econo™ Sterile instruments are suitable in a number of surgical venues where economical choice with reliable performance is most important. Convenient, ready-to-use sterile instruments provide the very best value with the greatest utility of any surgical instrument.

COLPOSCOPY SETS

OB/GYN

90-8555 — MINI COLPOSCOPY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
90-3990	Kogan Endocervical Speculum	With Gauge and Screw, Narrow Wire Blades	1
90-8336	Iris Hook	Quick Release, 1mm 9"	1
90-6285	Mini Endocervical Curette		1
90-7530	Mini Tischler Biopsy Forcep	2.3mm x 4.2mm Bite	1



90-8557 — STANDARD COLPOSCOPY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
90-3995	Kogan Endocervical Speculum	With Gauge and Screw, Wide Wire Blades	1
90-8334	Iris Hook	Quick Release, 4mm 9"	1
90-6612	Kevorkian-Younger Endocervical Curette	With Basket, Narrow Tip, 12"	1
90-7520	Tischler Oval Biopsy Forcep	3mm x 7mm Bite	1



OB/GYN

BIOPSY FORCEPS

General Surgery



GAYLOR BIOPSY PUNCH

IMPROVED JAWS

90-7385 9½"



GELHORN BIOPSY PUNCH

90-7492 9½"



SCHUBERT-VAN DOREN BIOPSY PUNCH

90-7511 10" Angular Tip and Shanks

SKLAR DISINFECTANT™ SURFACE WIPES

CONVENIENT DISPENSING CANISTER
READY-TO-USE

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
10-1616	100 Towelettes/Canister	6/CS
10-1617	100 Towelettes/Canister	12/CS

BIOPSY FORCEPS

General Surgery

OB/GYN

SKLAR FULLY ROTATIONAL BIOPSY FORCEPS

FULLY ROTATIONAL, Sklar's line of 360° OR Grade Rotating Biopsy Forceps are designed to obtain a perfect sample without obstructing the biopsy site. The shaft is detachable from the handle allowing for a full assortment of tip-shaft combinations to compliment the most difficult gynecological procedure. Flexible configurations, high-quality and precision craftsmanship are backed by a lifetime guarantee against manufacturing defects.

COMPLETE INSTRUMENT

TISCHLER OVAL BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-7551	3mm x 7mm Bite 7" Shaft Length
90-7571	3mm x 7mm Bite 12" Shaft Length



MINI TISCHLER BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-7553	2.3mm x 4.2mm Bite 7" Shaft Length
90-7573	2.3mm x 4.2mm Bite 12" Shaft Length



OVAL BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-7552	3mm x 5mm Bite 7" Shaft Length
90-7572	3mm x 5mm Bite 12" Shaft Length




KEVORKIAN BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-7554	3mm x 8mm Bite 7" Shaft Length
90-7574	3mm x 8mm Bite 12" Shaft Length



- 360° ROTATING SHAFT
- SIMPLE, TIME-SAVING DETACHABLE SHAFT
- FULL ASSORTMENT OF INTERCHANGEABLE TIP-SHAFT
- WORKABLE LENGTH OF SHAFT MEASURED, FROM HUB TO TIP



CANNULA WITH TIP ONLY

PRODUCT NUMBER	DESCRIPTION
90-7555	7" Tischler Oval
90-7556	7" Oval
90-7557	7" Mini Tischler
90-7558	7" Kevorkian
90-7576	12" Tischler Oval
90-7577	12" Oval
90-7578	12" Mini Tischler
90-7579	12" Kevorkian

HANDLE ONLY

PRODUCT NUMBER	DESCRIPTION
90-7560	Handle for use with 7" Cannula
90-7575	Handle for use with 12" Cannula

OB/GYN

BIOPSY FORCEPS

General Surgery

OVAL BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-7525	3mm x 5mm Bite
90-9565	3mm x 5mm Bite, 60° Opening



KEVORKIAN-YOUNGE BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-6610	3mm x 9.7mm Bite



EPPENDORFER BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-6820	4mm x 8mm Bite



WITTNER UTERINE BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-6600	4mm x 8mm Bite, Straight
90-6602	4mm x 8mm Bite, Angular



BIOPSY FORCEPS

General Surgery

OB/GYN

PRECISION, balance and control are the key features of Sklar's line of OR Grade Biopsy Punches and ensures the highest level of hand-crafted workmanship. A variety of tip-handle combinations provide the flexibility, comfort and control required by the most selective physicians. Select forcep tips are available with a 60° Wide Bite mouth for easy removal of biopsy and safer cleaning and are easily identified by their Dual Gold handles. Guaranteed for a lifetime against manufacturing defects.

TISCHLER KEVORKIAN BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-6887	3mm x 9.5mm Bite
90-9570	3mm x 9.5mm Bite, 60° Opening



MINI TISCHLER BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-7530	2.3mm x 4.2mm Bite
90-9555	2.3mm x 4.2mm, 60° Opening



TISCHLER OVAL BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-7520	3mm x 7mm Bite
90-9560	3mm x 7mm Bite, 60° Opening



KEVORKIAN BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-7540	3mm x 8mm Bite
90-9575	3mm x 8mm Bite, 60° Opening



OB/GYN

BIOPSY FORCEPS

Premium



MICRO-CUT GAYLOR BIOPSY PUNCH

90-8925 9½"



MICRO-CUT GELHORN BIOPSY PUNCH

90-8930 9½"

SKLAR DISINFECTANT™

SPRAY SURFACE DISINFECTANT

2 MINUTE KILL TIME

- Ready to use, non-toxic formula.
- Specially formulated for the control of cross contamination.
- Effective against C. Diff, TB, HBV, HCV, Hydrophilic and Lipophilic viruses including HIV-1 (AIDS Virus), bacteria and fungi on inanimate objects.
- Safe for all areas of the facility including the NICU, operating rooms, isolation rooms, patient care areas and laboratories.
- Satisfies OSHA's Bloodborne Pathogen Standards.
- Non-irritating formula.
- EPA registered.



PRODUCT NUMBER	DESCRIPTION	Quantities
10-1643	24 oz. Spray Bottles	6/cs
10-1653	1 Gallon Bottles	4/cs

BIOPSY FORCEPS

Premium

SKLAR FULLY ROTATIONAL BIOPSY FORCEPS

MICRO-CUT Premium Grade Biopsy Forceps rotate a full 360° and are designed to obtain a perfect sample without obstruction of the biopsy site. Combining fine craftsmanship with the added benefit of a knife honed edge on the upper jaw, Micro-Cut Forceps allow for a cleaner, more precise cut. A full assortment of tip-shaft combinations provide flexible configurations to compliment the most difficult gynecological procedure. Micro-Cut Forceps are easily identifiable by their dual gold handles.

COMPLETE INSTRUMENT

MICRO-CUT TISCHLER OVAL BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8960	3mm x 7mm Bite 7" Shaft Length
90-8981	3mm x 7mm Bite 12" Shaft Length



MICRO-CUT OVAL BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8962	3mm x 5mm Bite 7" Shaft Length
90-8982	3mm x 5mm Bite 12" Shaft Length




MICRO-CUT MINI TISCHLER BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8964	2.3mm x 4.2mm Bite 7" Shaft Length
90-8983	2.3mm x 4.2mm Bite 12" Shaft Length



MICRO-CUT KEVORKIAN BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8966	3mm x 8mm Bite 7" Shaft Length
90-8984	3mm x 8mm Bite 12" Shaft Length



- 360° ROTATING SHAFT
- SIMPLE, TIME-SAVING DETACHABLE SHAFT
- FULL ASSORTMENT OF INTERCHANGEABLE TIP-SHAFT
- WORKABLE LENGTH OF SHAFT MEASURED, FROM HUB TO TIP

CANNULA WITH TIP ONLY

PRODUCT NUMBER	DESCRIPTION
90-8970	7" Tischler Oval
90-8972	7" Oval
90-8974	7" Mini Tischler
90-8976	7" Kevorkian
90-8986	12" Tischler Oval
90-8987	12" Oval
90-8988	12" Mini Tischler
90-8989	12" Kevorkian

HANDLE ONLY

PRODUCT NUMBER	DESCRIPTION
90-8980	Handle for use with 7" Cannula
90-8985	Handle for use with 12" Cannula

OB/GYN

BIOPSY FORCEPS

Premium

MICRO-CUT OVAL BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8948	3mm x 5mm Bite



MICRO-CUT KEVORKIAN-YOUNGE BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8905	3mm x 9.7mm Bite



MICRO-CUT EPPENDORFER BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8010	4mm x 8mm Bite



MICRO-CUT WITTNER UTERINE BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8900	4mm x 8mm Bite, Straight



BIOPSY FORCEPS

Premium

OB/GYN

PRECISION ENGINEERED, hand-crafted and imported from Germany, Sklar Biopsy Punches have a reputation for being the best in the industry. A complete selection of forcep tips combined with three specific handle configurations provide the flexibility, comfort and control required by the most selective physicians.

MICRO-CUT

Micro-Cut Premium Grade Biopsy Punches feature a knife honed edge on the upper jaw for a cleaner, more precise cut. Precisely balanced and superior workmanship delivers the finest value in quality gynecological instruments. Micro-Cut Forceps are easily identified by their dual gold handles.

MICRO-CUT TISCHLER KEVORKIAN BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8920	3mm x 9.5mm Bite



MICRO-CUT MINI-TISCHLER PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8945	2.3mm x 4.2mm Bite



MICRO-CUT TISCHLER OVAL BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8940	3mm x 7mm Bite




MICRO-CUT KEVORKIAN BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION
90-8952	3mm x 8mm Bite




OB/GYN

OBSTETRICAL FORCEPS



ATRAUMATIC OVUM FORCEP

92-5280 10" Straight



GREEN ARMYTAGE FORCEP

90-5420 8 1/4"



COLLYER PELVIMETER

92-4100



MARTIN PELVIMETER

92-4130

SKLAR SURGICAL GOWNS

- Low linting, lightweight and breathable SMS fabric.
- Generous cut provides comfort and freedom of movement during various procedures.
- Superior liquid barrier protection provides confidence against strikethrough.
- Velcro-style neck closure easily adjusts for a perfect fit.
- Meets safety standards for Class I Flame Resistance.
- Meets AAMI Level I requirements.
- Disposable, single use.

AVAILABLE IN NON-STERILE AND BULK, NON-STERILE. SEE PAGE 233 FOR ORDERING INFORMATION.



OBSTETRICAL FORCEPS

OB/GYN



KELLY PLACENTA FORCEP

CURVED


92-5012 12"



JAVERTS PLACENTA FORCEPS


92-5095 9½" Serrated Jaws

92-5195 9½" Smooth Jaws



COREY OVUM FORCEP

92-5295 9½"



BIERER OVUM FORCEPS


WITH RATCHET

92-5298 11" Serrated Jaw

WITHOUT RATCHET


91-3822 11" Serrated Jaw

91-3823 12½" Serrated Jaw



BOND PLACENTA FORCEP

92-5412 12"




BARRETT PLACENTA FORCEP

92-5310 10"


OB/GYN

OBSTETRICAL FORCEPS




**LAUE DIVERGENT
OBSTETRICAL FORCEP**

92-2912 12½"



**MURLESS
HEAD EXTRACTOR**

92-2915 10½"




BILL TRACTION HANDLE
FOR USE ON STANDARD OBSTETRICAL
FORCEPS

92-3350




**LUIKART-BILL TRACTION
HANDLE**

92-3360



BEACHAM AMNIOTOME

91-4210 10¼"



IOWA MEMBRANE FORCEP

S-CURVE

91-4310 10"

OBSTETRICAL FORCEPS

OB/GYN



**MCLANE-TUCKER-LUIKART
OBSTETRICAL FORCEP**

92-2215 15½"



**SIMPSON-LUIKART
OBSTETRICAL FORCEP**

92-2314 14"



**DE LEE OBSTETRICAL
FORCEPS**

92-1012	12"	Short for Caesareans
92-1014	14"	Standard
92-1015	15¼"	Long



**SIMPSON OBSTETRICAL
FORCEPS**

92-1112	12"	Short
92-1114	14¼"	Long



**MCLANE-TUCKER OBSTETRICAL
FORCEP**

92-1415 15¼"



HALE OBSTETRICAL FORCEP

92-1610 10¼"

OB/GYN

OBSTETRICAL FORCEPS



ELLIOT OBSTETRICAL FORCEPS

92-1712 12½"

92-1715 15"



**PIPER OBSTETRICAL FORCEP
FOR BREECH PRESENTATION**

92-2717 17½"



**LAUFE-PIPER OBSTETRICAL
FORCEP**

92-2812 12"



KIELLAND OBSTETRICAL FORCEP

92-1915 15¾"



LUIKART OBSTETRICAL FORCEP

92-2015 15¾"



**KIELLAND-LUIKART
OBSTETRICAL FORCEP**

92-2115 15¾"

HYSTERECTOMY/PARAMETRIUM Z-Type Scissors & Forceps

Sklar's Z-Type Hysterectomy/Parametrium Scissors and Forceps are specifically designed to ideally meet the surgeon's needs during vaginal and abdominal hysterectomy procedures. Manufactured from German grade stainless steel and under the highest quality standards, these specialized instruments provide uncompromised balance, feel and performance.

Z-TYPE SCISSORS

With solid Tungsten Carbide inserts, Z-Type Scissors have extremely sharp, strong and precise cutting edges with lasting durability, and are specifically formed to cut tissue laterally. They are designed to complement the curvature of Z-Type Forceps for optimal access, visibility and control.

PRODUCT NUMBER		DESCRIPTION
16-4000	9"	Z-Type Parametrium Scissors, Slight Curve
16-4001	9"	Z-Type Parametrium Scissors, Curved
16-4002	9"	Z-Type Parametrium Scissors, Full Curve
16-4004	10½"	Z-Type Parametrium Scissors, Slight Curve
16-4005	10½"	Z-Type Parametrium Scissors, Curved
16-4006	10½"	Z-Type Parametrium Scissors, Full Curve
16-4008	12"	Z-Type Parametrium Scissors, Slight Curve
16-4009	12"	Z-Type Parametrium Scissors, Curved
16-4010	12"	Z-Type Parametrium Scissors, Full Curve
16-4012	14"	Z-Type Parametrium Scissors, Slight Curve
16-4013	14"	Z-Type Parametrium Scissors, Curved
16-4014	14"	Z-Type Parametrium Scissors, Full Curve



Z-TYPE FORCEPS


Longitudinal grooves with tip serrations on Z-Type Forceps ensure an atraumatic yet firm grip on tissue, minimizing slippage and tissue trauma. Longer instrument design and off-set handles ensure superior manipulation and exposure especially for deep pelvic surgery or obese patients.

PRODUCT NUMBER		DESCRIPTION
90-2785	8¼"	Z-Type Hysterectomy Forceps, Straight
90-2786	8¼"	Z-Type Hysterectomy Forceps, Slight Curve
90-2787	8¼"	Z-Type Hysterectomy Forceps, Full Curve
90-2788	8¼"	Z-Type Hysterectomy Forceps, Angled
90-2790	9½"	Z-Type Hysterectomy Forceps, Straight
90-2792	9½"	Z-Type Hysterectomy Forceps, Slight Curve
90-2793	9½"	Z-Type Hysterectomy Forceps, Full Curve
90-2794	9½"	Z-Type Hysterectomy Forceps, Angled
90-2796	12"	Z-Type Hysterectomy Forceps, Straight
90-2797	12"	Z-Type Hysterectomy Forceps, Slight Curve
90-2798	12"	Z-Type Hysterectomy Forceps, Full Curve
90-2799	12"	Z-Type Hysterectomy Forceps, Angled
90-2801	14"	Z-Type Hysterectomy Forceps, Straight
90-2802	14"	Z-Type Hysterectomy Forceps, Slight Curve
90-2803	14"	Z-Type Hysterectomy Forceps, Full Curve
90-2804	14"	Z-Type Hysterectomy Forceps, Angled



OB/GYN

CURETTES



RANDALL ENDOMETRIAL CURETTE

90-6295 9½" 4mm



SKLAR MINI ENDOCERVICAL CURETTE

90-6285 Small rounded tip for easier entry into stenosed cervix



KEVORKIAN-YOUNGE ENDOCERVICAL CURETTES

NARROW TIP
90-6611 12" Without Basket
90-6612 12" With Basket



LAUE-NOVAK DIAGNOSTIC CURETTE

WITH SUCTION FINGER CONTROL
90-5834 Small



IUD EXTRACTOR

90-8210 10½"



HOLMAN-MATHIEU CANNULA FOR HYSTEROSALPINGOGRAPHY

90-6582 Small 8¼" (21cm) with Luer-Lock Chrome Plated
90-6590 Large 9" (22.9cm) with Luer-Lock Chrome Plated

CURETTES

OB/GYN



BUMM UTERINE CURETTES

13½", WITH MALLEABLE SHAFTS

SHARP

90-5602 #14

90-5603 #15


90-5604 #16

BLUNT

90-5605 #14

90-5606 #15


90-5607 #16



HOLDEN UTERINE CURETTE

SEMI-SHARP


90-5710 10¾" Extra Small



GREENE UTERINE CURETTE

SERRATED

90-5566 11"



HEANEY UTERINE CURETTES

SERRATED

90-6095 9½" Open Curette

90-6195 9½" With Basket



KELLY-GRAY UTERINE CURETTE


MALLEABLE

90-6081 9½" #1

90-6082 9½" #2


90-6083 9½" #3

90-6084 9½" #4



LOUNSBURY PLACENTA CURETTE

90-5713 13" ¾" Wide



NOVAK ENDOMETRIAL CURETTES

90-6289 9½" 1mm


90-6290 9½" 2mm

90-6291 9½" 3mm

90-6292 9½" 4mm

OB/GYN


CURETTES



THOMAS UTERINE CURETTES

BLUNT, 11"

90-5431 #1	90-5435 #5
90-5432 #2	90-5436 #6
90-5433 #3	90-5438 Set of 6



SIMS UTERINE CURETTES

SHARP, 11"

90-5499 #00	90-5424 #4
90-5500 #0	90-5425 #5
90-5421 #1	90-5426 #6
90-5422 #2	90-5428 Set of 6
90-5423 #3	



RECAMIER UTERINE CURETTES

SHARP, 12½", WITH MALLEABLE SHAFTS

90-5610 #00	90-5616 #5
90-5611 #0	90-5617 #6
90-5612 #1	90-5618 #7
90-5613 #2	90-5619 #8
90-5614 #3	90-5620 #9
90-5615 #4	90-5621 #10

BLUNT, 12½", WITH MALLEABLE SHAFTS

90-5622 #00	90-5628 #5
90-5623 #0	90-5629 #6
90-5624 #1	90-5630 #7
90-5625 #2	90-5631 #8
90-5626 #3	90-5632 #9
90-5627 #4	90-5633 #10

DILATORS/PUNCHES

OB/GYN

General Surgery



WALTHER DILATOR

CHROME

85-2012	12 French	85-2028	28 French
85-2014	14 French	85-2030	30 French
85-2016	16 French	85-2032	32 French
85-2018	18 French	85-2034	34 French
85-2020	20 French	85-2036	36 French
85-2022	22 French	85-2038	38 French
85-2024	24 French	85-2010	Set of 14
85-2026	26 French		



FEMALE CATHETER

CHROME

85-2210	10 French
85-2212	12 French
85-2214	14 French
85-2216	16 French
85-2218	18 French
85-2220	20 French



KEYES VULVAR/DERMAL PUNCH

06-4100	Set of 6, 2mm thru 7mm with Universal Handle and Sterilizable Base
06-4110	KEYES Sharpening Stone



KEYES CUTANEOUS PUNCH


ROUND

06-4142	2mm	06-4146	6mm
06-4143	3mm	06-4147	7mm
06-4144	4mm	06-4148	8mm
06-4145	5mm	06-4110	KEYES Sharpening Stone

OB/GYN


DILATORS

General Surgery



PRATT UTERINE DILATOR


90-5215 13/15 French
90-5219 17/19 French
90-5223 21/23 French
90-5227 25/27 French
90-5231 29/31 French
90-5235 33/35 French



KLEEGMAN DILATOR


NARROW

90-4902 10½"



GOODELL UTERINE DILATOR

90-4713 13"



MINI DILATORS

SET WITH HANDLE


90-6280 Set Includes: 1.0mm
1.5mm
2.0mm
2.5mm
3.0mm

90-6281 Handle for Mini Dilator Set

DILATORS

General Surgery


OB/GYN



HEGAR UTERINE DILATOR

DOUBLE END, 7"


90-4802	1/2 mm
90-4804	3/4 mm
90-4806	5/6 mm
90-4808	7/8 mm
90-4810	9/10 mm
90-4812	11/12 mm
90-4814	13/14 mm
90-4816	15/16 mm
90-4818	17/18 mm
90-4828	27/28 mm
90-4832	31/32 mm
90-4825	Set of 8, 3/4mm thru 17/18mm without Khaki Roll
90-4826	Set of 8, 3/4mm thru 17/18mm with Khaki Roll



HEGAR UTERINE DILATOR

SINGLE END, 7"

90-4784	1.0 mm
90-4785	1.5 mm
90-4786	2.0 mm
90-4787	2.5 mm
90-4788	3.0 mm
90-4789	3.5 mm
90-4790	4.0 mm
90-4791	4.5 mm
90-4792	5.0 mm
90-4793	5.5 mm
90-4794	6.0 mm
90-4795	6.5 mm
90-4796	7.0 mm
90-4797	7.5 mm
90-4798	8.0 mm
90-4799	9.0 mm



HANK UTERINE DILATOR

AMERICAN SCALE

90-4908	7/8
90-4910	9/10
90-4912	11/12
90-4914	13/14
90-4916	15/16
90-4918	17/18
90-4920	19/20
90-4928	27/28
90-4925	Set of 6, 9/10 thru 19/20 without Khaki Roll
90-4926	Set of 6, 9/10 thru 19/20 with Khaki Roll
90-4927	Khaki Roll

OB/GYN

SOUNDS

General Surgery



DITTEL URETHRAL SOUNDS

STRAIGHT, 11 1/2"

85-1308	8 French
85-1310	10 French
85-1312	12 French
85-1314	14 French
85-1316	16 French
85-1318	18 French
85-1320	20 French
85-1322	22 French
85-1324	24 French
85-1326	26 French
85-1328	28 French
85-1330	30 French
85-1332	32 French
85-1334	34 French
85-1336	36 French
85-1338	38 French
85-1340	40 French

SETS IN KHAKI ROLL

85-1342	Set of 6	8, 12, 16, 20, 24, 30FR
85-1345	Set of 17	8-40 French



OTIS DITTEL URETHRAL SOUNDS

CURVED, 11 1/2"

85-1408	8 French
85-1410	10 French
85-1412	12 French
85-1414	14 French
85-1416	16 French
85-1418	18 French
85-1420	20 French
85-1422	22 French
85-1424	24 French
85-1426	26 French
85-1428	28 French
85-1430	30 French
85-1432	32 French
85-1434	34 French
85-1436	36 French
85-1438	38 French
85-1440	40 French

SETS IN KHAKI ROLL

85-1442	Set of 8	16-30 French
85-1446	Set of 17	8-40 French



OTIS BOUGIE A BOULE

STRAIGHT, 9 3/4"

85-1608	8 French
85-1610	10 French
85-1612	12 French
85-1614	14 French
85-1616	16 French
85-1618	18 French
85-1620	20 French
85-1622	22 French
85-1624	24 French
85-1626	26 French
85-1628	28 French
85-1630	30 French
85-1632	32 French

SETS IN KHAKI ROLL

85-1635	Set of 8	16-30 French
85-1636	Set of 13	8-32 French


PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SOUNDS

OB/GYN


General Surgery



SIMS UTERINE SOUNDS

MALLEABLE


90-5312	12½"	Non-Graduated, Silver Plated
90-5314	12½"	Centimeter Graduated, Silver Plated



SIMPSON UTERINE SOUNDS

MALLEABLE WITH NODULE JUST BELOW CURVE

90-5412	12½"	Inch Graduated, Silver Plated
90-5413	12½"	Centimeter Graduated, Silver Plated
90-5415	12½"	Centimeter Graduated, Silver Plated



VAN BUREN URETHRAL SOUNDS

85-1008	8 French	85-1026	26 French
85-1010	10 French	85-1028	28 French
85-1012	12 French	85-1030	30 French
85-1014	14 French	85-1032	32 French
85-1016	16 French	85-1034	34 French
85-1018	18 French	85-1036	36 French
85-1020	20 French	85-1038	38 French
85-1022	22 French	85-1040	40 French
85-1024	24 French		


SETS IN KHAKI ROLL

85-1065	Set of 17	Sizes 8-40 French
85-1070	Set of 6	Sizes 14-24 French
85-1075	Set of 12	Sizes 8-30 French
85-1076	Khaki Roll	Holds up to 12 Sounds

OB/GYN

HOOKS

General Surgery



TYRELL IRIS HOOKS

65-3205	5"	3mm	Blunt
65-3206	5"	3mm	Sharp




GILLIES SKIN HOOKS

47-2570	7"	3mm	Sharp
47-2572	7"	4mm	Sharp



FRAZIER TISSUE HOOKS


47-2555	5"	Sharp	Extra Delicate
47-2556	5"	Blunt	Extra Delicate



GUTHRIE DOUBLE HOOKS

SHARP

47-2651	5"	3mm	2 Small Prongs
47-2650	5"	5mm	2 Large Prongs



JOSEPH SKIN HOOKS


SHARP

47-2500	6 1/4"	-	1 Prong
47-2502	6 1/4"	2mm	2 Prongs
47-2505	6 1/4"	5mm	2 Prongs
47-2507	6 1/4"	7mm	2 Prongs
47-2510	6 1/4"	10mm	2 Prongs
47-2512	6 1/4"	12mm	2 Prongs

HOOKS

OB/GYN


General Surgery



IRIS HOOKS

QUICK RELEASE

90-8336	9"	1mm	90-8348	12"	1mm
90-8335	9"	3mm	90-8346	12"	3mm
90-8334	9"	4mm	90-8344	12"	4mm
90-8333	9"	6mm	90-8342	12"	6mm
90-8332	9"	Set of 4	90-8341	12"	Set of 4



EMMETT HOOKS

90-8301	9"	Style 1
90-8302	9"	Style 2
90-8303	9"	Style 3
90-8304	9"	Style 4
90-8305	9"	Style 5
90-8306	9"	Set of 5

OB/GYN

HOOKS

Sklar Black™ Laser



TYRELL IRIS HOOKS

90-1056	5"	3mm	Blunt
90-1058	5"	3mm	Sharp



GILLIES SKIN HOOKS

90-1045	7"	3mm	Sharp
90-1047	7"	4mm	Sharp



FRAZIER TISSUE HOOKS


90-1048	5"	Sharp	Extra Delicate
90-1049	5"	Blunt	Extra Delicate



GUTHRIE DOUBLE HOOKS

SHARP

90-1046	5"	3mm	2 Small Prongs
90-1044	5"	5mm	2 Large Prongs



JOSEPH SKIN HOOKS


SHARP

90-1050	6 1/4"	-	1 Prong
90-1051	6 1/4"	2mm	2 Prongs
90-1052	6 1/4"	5mm	2 Prongs
90-1053	6 1/4"	7mm	2 Prongs
90-1054	6 1/4"	10mm	2 Prongs
90-1055	6 1/4"	12mm	2 Prongs

HOOKS

OB/GYN

Sklar Black™ Laser



IRIS HOOKS

QUICK RELEASE

90-2126	9"	1mm	90-2128	12"	1mm
90-2122	9"	3mm	90-2124	12"	3mm
90-2118	9"	4mm	90-2120	12"	4mm
90-2114	9"	6mm	90-2116	12"	6mm
90-2160	9"	Set of 4	90-2131	12"	Set of 4




EMMETT TENACULUM HOOKS

90-2171	9"	Style 1
90-2172	9"	Style 2
90-2173	9"	Style 3
90-2174	9"	Style 4
90-2175	9"	Style 5
90-2170	9"	Set of 5

OB/GYN

HOOKS

Sklar Blue™ Electrosurgical



IRIS HOOKS

QUICK RELEASE

91-5244	9"	1mm	91-5248	12"	1mm
91-5243	9"	3mm	91-5247	12"	3mm
91-5242	9"	4mm	91-5246	12"	4mm
91-5241	9"	6mm	91-5245	12"	6mm
91-5236	9"	Set of 4	91-5238	12"	Set of 4

EMMETT HOOKS

91-5251	9"	Style 1
91-5252	9"	Style 2
91-5253	9"	Style 3
91-5254	9"	Style 4
91-5255	9"	Style 5
91-5250	9"	Set of 5

ENDOCERVICAL SPECULA

General Surgery

OB/GYN



NARROW WIRE BLADES

KOGAN ENDOCERVICAL SPECULUM

WITH GAUGE AND SCREW

90-3990 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH SPRING HANDLE

90-3996 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH RATCHET HANDLE

90-3997 Narrow Wire Blades



WIDE WIRE BLADES

KOGAN ENDOCERVICAL SPECULUM

WITH GAUGE AND SCREW

90-3995 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM

WITH SPRING HANDLE

90-3994 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM


WITH RATCHET HANDLE

90-3998 Wide Wire Blades

OB/GYN


ENDOCERVICAL SPECULA

Sklar Black™ Laser



KOGAN ENDOCERVICAL SPECULUM
WITH GAUGE AND SCREW

90-8855 Narrow Wire Blades




KOGAN ENDOCERVICAL SPECULUM
WITH SPRING HANDLE

90-8865 Narrow Wire Blades




KOGAN ENDOCERVICAL SPECULUM
WITH RATCHET HANDLE

90-8862 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM
WITH GAUGE AND SCREW

90-8858 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM
WITH SPRING HANDLE

90-8850 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM
WITH RATCHET HANDLE

90-8860 Wide Wire Blades

ENDOCERVICAL SPECULA


Sklar Blue™ Electrosurgical

OB/GYN



KOGAN ENDOCERVICAL SPECULUM
WITH GAUGE AND SCREW

91-5230 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM
WITH SPRING HANDLE

91-5228 Narrow Wire Blades




KOGAN ENDOCERVICAL SPECULUM
WITH RATCHET HANDLE

91-5235 Narrow Wire Blades



KOGAN ENDOCERVICAL SPECULUM
WITH GAUGE AND SCREW

91-5232 Wide Wire Blades



KOGAN ENDOCERVICAL SPECULUM
WITH SPRING HANDLE

91-5226 Wide Wire Blades




KOGAN ENDOCERVICAL SPECULUM
WITH RATCHET HANDLE

91-5234 Wide Wire Blades

OB/GYN

RETRACTORS

General Surgery



BREISKY RETRACTORS

90-3520	100mm x 20mm
90-3521	100mm x 25mm
90-3522	100mm x 30mm
90-3528	130mm x 40mm
90-3529	156mm x 40mm



DELEE RETRACTOR

90-3452	9 1/2"	2 3/4" x 2"	Universal
---------	--------	-------------	-----------



HEANEY RETRACTORS

60-1616	10 1/2"	Blade 1 3/4" x 1"
60-1617	10 1/2"	Blade 3" x 1"
60-1618	10 1/2"	Blade 3 1/4" x 1 1/8"




GELPI RETRACTORS

22-8850	3 1/2"
22-8853	5 1/2"
22-8855	6 1/4"
60-6570	7"
60-6580	7 1/2" With Stops



RIGBY RETRACTORS

TRU-GRIP RATCHET HANDLE	
60-6370	7"



O'SULLIVAN-O'CONNOR VAGINAL RETRACTORS


SMALL WITH THREE INTERCHANGEABLE BLADES

90-3045	1 1/2" x 3/4"	Small
	3" x 1 1/2"	Medium
	3 1/2" x 1 3/4"	Large

RETRACTORS

General Surgery


OB/GYN



SIMS VAGINAL RETRACTORS


SKLAR-GRIP HANDLE

90-3327	2 3/4" x 1"	Small
90-3330	3" x 1 1/4"	Medium
90-3335	3 1/2" x 1 1/2"	Large




BRANTLEY-TURNER RETRACTOR

90-3470 7" Blade 4" x 1 1/4"




PRYOR-PEAN RETRACTOR

90-3510 10 1/4" Blade 4 1/4" x 1"




DOYEN RETRACTORS

60-2185	9"	Blade 1 3/4" x 2 1/8"	Small
60-2190	9"	Blade 1 3/4" x 3 1/4"	Medium
60-2195	9"	Blade 2 1/4" x 3 1/4"	Large



TAUBER VAGINAL SPATULA

90-3692 10" Blade 3 1/2" x 1 3/8"



SECOND VAGINAL SPATULA

90-3610 10 1/2" Blade 4 3/4" x 1 3/8"

OB/GYN

RETRACTORS

General Surgery



LATERAL VAGINAL RETRACTOR

RATCHET HANDLE


90-3980 Blade fits inside Vaginal Speculum



SIMS VAGINAL RETRACTORS

DOUBLE END

90-3150	7/8" x 2 1/2"	1" x 3"	#1
90-3151	1" x 2 1/2"	1 1/8" x 3"	#2
90-3152	1" x 3"	1 1/4" x 3 1/4"	#3
90-3153	1 1/4" x 3"	1 1/2" x 3 1/4"	#4
90-3154	1 1/4" x 3 1/4"	1 1/2" x 4"	#5
90-3155	1 1/4" x 3 1/4"	1 3/4" x 4"	#6



SIMS VAGINAL RETRACTORS

SINGLE END


90-3225	2 1/2" x 1"	Small
90-3230	3" x 1 1/4"	Medium
90-3235	3 1/2" x 1 1/2"	Large



JACKSON VAGINAL RETRACTORS

SINGLE END

90-3130	3" x 1 1/2"	Small
90-3135	3 1/2" x 1 1/2"	Medium
90-3140	4" x 1 1/2"	Large



EASTMAN RETRACTORS

90-3030	3" x 1 1/2"	Small
90-3035	3 1/2" x 1 1/2"	Medium
90-3040	4" x 1 1/2"	Large



SIMS-KELLY VAGINAL RETRACTORS

SKLAR-GRIP HANDLE

90-3422	2 1/4" x 1"	Small
90-3427	2 3/4" x 1 1/4"	Medium
90-3430	3" x 1 1/2"	Large

RETRACTORS

OB/GYN

Sklar Black™ Laser



LATERAL VAGINAL RETRACTOR

RATCHET HANDLE

90-8845 Blade fits inside Vaginal Speculum



JACKSON RETRACTORS

SINGLE END

90-9064	3" x 1½"	Small with Permanent Smoke Evacuation Tube
90-9065	3½" x 1½"	Medium without Tube or Loops
90-9066	3½" x 1½"	Medium with Loop to hold and direct disposable Smoke Evacuation Tubing
90-9069	3½" x 1½"	Medium with Permanent Smoke Evacuation Tube
90-9067	4" x 1½"	Large without Tube or Loops
90-9068	4" x 1½"	Large with Loop to hold and direct disposable Smoke Evacuation Tubing
90-9070	4" x 1½"	Large with Permanent Smoke Evacuation Tube

SKLAR INSTRUMENT CARE KIT

A TRIO OF SKLAR'S CARE AND CLEANING PRODUCTS
TO KEEP INSTRUMENTS IN TIPTOP SHAPE.

Each Kit Contains:

- 1 Sklar Lube™ 8 oz. Spray Bottle
- 1 Sklar Kleen™ Liquid 8 oz. Bottle
- 1 Sklar Polish™ 8 oz. Jar
- 1 Stainless Steel Bristle Cleaning Brush
- 1 Nylon Bristle Cleaning Brush

PRODUCT NUMBER

10-1625



OB/GYN

RETRACTORS

Sklar Blue™ Electrosurgical



LATERAL VAGINAL RETRACTOR
RATCHET HANDLE

91-5190 Blade fits inside Vaginal Speculum



SIMS VAGINAL RETRACTORS
SINGLE END

91-5200 Medium
91-5202 Medium, with loop to hold and direct Disposable Smoke Evacuation Tubing



SIMS VAGINAL RETRACTORS
DOUBLE END

91-5210 Medium
91-5212 Medium, with loop to hold and direct Disposable Smoke Evacuation Tubing



JACKSON VAGINAL RETRACTORS
SINGLE END

91-5215 3" x 1½" Small
91-5218 3½" x 1½" Medium
91-5220 4" x 1½" Large



TAUBER VAGINAL SPATULA
91-5176 10" Blade 3½" x 1¾"



SECOND VAGINAL SPATULA
91-5178 10½" Blade 4¾" x 1¾"

r retractors

Sklar offers a full line of Retractors in a variety of the most popular styles and sizes to meet the surgeon's need in any situation. Specifically designed to actively separate incisions or safely hold back underlying organs and tissue for optimal access and exposure of the site. Avoid unintended injury to surrounding tissue by selecting blunt retractors or opt for a more secure grip that sharp retractors can afford. Both styles allow for optimum flexibility and precision while working.

Sklar Blue™ Electrosurgical Retractors Page 16

Designed to provide an effective, non-conductive, patient-safe surface, Electrosurgical Retractors feature our exclusive Sklar Blue™ insulated coating that prevents transmission of RF current through the instrument. Manufactured of the highest quality stainless steel, these instruments are fully autoclavable. All Sklar Blue™ electrosurgical coatings are guaranteed for one year against manufacturing defects.

Sklar Black™ Laser Retractors Page 17

Sklar Black™ Laser Retractors use an improved black polymer coating developed specifically for Laser surgery. Testing demonstrates that the Sklar Black™ coating over a sandblasted surface reduced reflective laser energy to less than one-tenth of one percent. The durable coating improves patient safety and provides a reliable instrument you can use with confidence. All Sklar Laser coatings are guaranteed for one year against manufacturing defects.

General Surgery Retractors Page 18

Our very best operating room, high quality instruments incorporating optimum grades of German grade surgical steel, finely detailed craftsmanship and exacting standards of manufacturing. Choose from a full line of the most popular patterns.

OB/GYN


SPECULA

Econo™ Sterile



GRAVES SPECULUMS


96-2603	Small	25/cs
96-2605	Medium	25/cs
96-2607	Large	25/cs



DISPOSABLE GRAVES SPECULUMS


PLASTIC

96-2602	Small	50/cs
96-2604	Medium	50/cs
96-2606	Large	50/cs



PEDERSON SPECULUMS

96-2613	Small	25/cs
96-2615	Medium	25/cs
96-2617	Large	25/cs



DISPOSABLE CUSCO SPECULUMS

PLASTIC


96-2642	Small	50/cs
96-2644	Medium	50/cs
96-2648	Large	50/cs

SPECULA Surgi-OR™

OB/GYN



Hand-crafted and forged from Japanese surgical grade stainless steel, Surgi-OR™ instruments are built to last. Precision engineering and stringent quality control guidelines ensure superior performance. Available in the most popular patterns, this cost effective solution doesn't sacrifice quality for value and convenience in the Operating Room and other important clinical settings.



GRAVES VAGINAL SPECULUMS

95-304	3" x 1½"	X-Small
95-303	3" x ¾"	Small
95-312	4" x 1¼"	Medium
95-324	4½" x 1⅜"	Large
95-322	4¾" x 1½"	X-Large



PEDERSON VAGINAL SPECULUMS

95-334	2¾" x ½"	X-Small
95-347	3" x ½"	Small
95-348	4" x ⅞"	Medium
95-362	4¾" x 1"	Large
95-368	4¾" x 1⅝"	X-Large

SKLAR SOAK™

SURGICAL INSTRUMENT SOAK

Immersing soiled instruments in SKLAR SOAK™ for two minutes substantially reduces the risk of exposure to many infectious organisms.

- Effective Against: H1N1, HIV-1, HPV, VRE, Norovirus, Hep B & C as well as TB, MRSA, A. baumannii, and C. difficile in its vegetative form.
- Non-Glutaraldehyde.
- PH Neutral Formula designed for safe use on all surgical instruments & stainless steelware.
- Non-Caustic Formula allows for use as a manual cleaner or in ultrasonic units.
- EPA Registered.
- Meets all OSHA standards on blood borne pathogen destruction.
- Ready to use.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1603	1 Gallon Bottles	4/case
10-1609	1 Gallon Bottle	1 bottle



OB/GYN

SPECULA

General Surgery



AUVARD WEIGHTED VAGINAL SPECULUMS

90-3525	2¾" x 1½"	2 lbs. Light
90-3526	2¾" x 1½"	2½ lbs. Medium
90-3527	2¾" x 1½"	3 lbs. Heavy
90-3530	4" x 1¾"	2½ lbs. Medium, Long Blade



GARRIGUE WEIGHTED VAGINAL SPECULUMS

90-3535	3" x 1½"	2 lbs. Light
90-3536	3" x 1½"	2½ lbs. Medium
90-3537	3" x 1½"	3 lbs. Heavy



STEINER-AUVARD VAGINAL SPECULUM

90-3559	5½" x 1¼"	2½ lbs. Medium
---------	-----------	----------------



HENROTIN WEIGHTED VAGINAL SPECULUM

90-3555	3¼" x 1¼"	3 lbs. Heavy
---------	-----------	--------------




BERLIND-AUVARD WEIGHTED VAGINAL SPECULUM

90-3558	4¼" x 1⅝"	2½ lbs. Medium
---------	-----------	----------------

SPECULA


OB/GYN

General Surgery



GRAVES OB-VIEW VAGINAL SPECULUMS

90-3740 4" x 1 1/4" Medium
90-3741 4 3/4" x 1 1/2" Large



FREEWAY GRAVES VAGINAL SPECULUM

OPEN SIDE, 90° ANGLE

90-3750 4" x 1 1/4" Medium




FREEWAY GRAVES VAGINAL SPECULUM

OPEN SIDE, 45° ANGLE

90-3752 4" x 1 1/4" Medium



SIDE
SCREW



CENTER
SCREW

CUSCO VAGINAL SPECULUMS

WITH SIDE SCREW

90-3957 3" x 1 1/4" Small
90-3959 3 1/4" x 1 1/4" Medium
90-3960 3 3/4" x 1 3/8" Large


WITH CENTER SCREW

90-3952 3" x 1 1/4" Small
90-3953 3 1/4" x 1 1/4" Medium
90-3954 3 3/4" x 1 3/8" Large



COLLIN VAGINAL SPECULUMS

90-3943 3 1/2" x 1 1/4" Small
90-3945 4 1/4" x 1 1/4" Medium
90-3944 4 1/2" x 1 1/2" Large



DEVILBISS VAGINAL SPECULUM

PARALLEL OPENING, CHROME

90-3865 4" x 1 1/4"

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

OB/GYN

SPECULA

General Surgery



GRAVES VAGINAL SPECULUMS

90-3800	3" x 1/2"	Extra Small
90-3700	3" x 3/4"	Small
90-3701	4" x 1 1/4"	Medium
90-3702	4 1/2" x 1 3/8"	Large
90-3802	4 3/4" x 1 1/2"	Extra Large
90-3803	6" x 1 1/2"	Extra Long



PEDERSON VAGINAL SPECULUMS

90-3849	2 3/4" x 1/2"	Extra Small
90-3710	3" x 1/2"	Small
90-3712	4" x 5/8"	Slender
90-3713	4" x 7/8"	Medium
90-3711	4 3/4" x 1"	Large
90-3851	4 3/4" x 1 5/8"	Extra Large



PEDERSON OPEN-SIDED VAGINAL SPECULUM

OPEN ON RIGHT SIDE

90-3751	4" x 1"
---------	---------



GRAVES VAGINAL SPECULUMS

GOLD SCREWS

90-3720	3/4" x 3"	Small
90-3721	1 1/4" x 4"	Medium
90-3722	1 3/8" x 4 1/4"	Large



PEDERSON VAGINAL SPECULUMS

GOLD SCREWS

90-3730	1/2" x 3"	Small
90-3731	7/8" x 4"	Medium
90-3732	1" x 4 3/4"	Large



GRAVES VAGINAL SPECULUM WIDE VIEW PATTERN


WITH LARGE PROXIMAL OPENING FOR INSTRUMENTATION

90-3705	4" x 1 1/4"	Medium
---------	-------------	--------

SPECULA

OB/GYN


Sklar Black™ Laser



AUVARD VAGINAL SPECULUMS

MEDIUM, LONG BLADES

- 90-9106 4" x 1 3/4" 2 1/2 lbs. without loop
- 90-9108 4" x 1 3/4" 2 1/2 lbs. with loop to hold and direct disposable Smoke Evacuation Tubing
- 90-9109 4" x 1 3/4" 2 1/2 lbs. with permanent Smoke Evacuation Tube



STEINER-AUVARD VAGINAL SPECULUMS

MEDIUM, EXTRA LONG BLADE

- 90-3560 5 1/2" x 1 1/4" 2 1/2 lbs. without loop
- 90-3561 5 1/2" x 1 1/4" 2 1/2 lbs. with loop to hold and direct disposable Smoke Evacuation Tubing
- 90-3562 5 1/2" x 1 1/4" 2 1/2 lbs. with permanent Smoke Evacuation Tube



SKLAR MONSEL'S

FERRIC SUBSULFATE PASTE

A ONE-TIME USE pre-mixed paste packaged in convenient ready-to-use bottles.

Also includes 12 rayon tip applicators.

Product Number

96-1200

SKLAR LUGOL'S

STRONG IODINE SOLUTION USP

A ONE-TIME USE pre-mixed iodine solution packaged in convenient ready-to-use bottles.

Also includes 12 rayon tip applicators.

Product Number

96-1210

OB/GYN

SPECULA

Sklar Black™ Laser



GRAVES VAGINAL SPECULUMS WIDE VIEW PATTERN

WITH LARGE PROXIMAL OPENING FOR
INSTRUMENTATION

- | | | |
|---------|-------------|---|
| 90-9054 | 4" x 1 1/4" | Medium, with permanent Smoke Evacuation Tube |
| 90-9055 | 4" x 1 1/4" | Medium, with loops to hold and direct Smoke Evacuation Tubing |



FREEWAY GRAVES VAGINAL SPECULUMS

OPEN SIDE, 90° ANGLE


- | | | |
|---------|-------------|---|
| 90-3754 | 4" x 1 1/4" | Medium, with permanent Smoke Evacuation Tube (not illustrated) |
| 90-3755 | 4" x 1 1/4" | Medium, with loop to hold and direct disposable Smoke Evacuation Tubing (not illustrated) |



FREEWAY GRAVES VAGINAL SPECULUMS

OPEN SIDE, 45° ANGLE

- | | | |
|---------|-------------|---|
| 90-3814 | 4" x 1 1/4" | With permanent Smoke Evacuation Tube (not illustrated) |
| 90-3815 | 4" x 1 1/4" | With loop to hold and direct disposable Smoke Evacuation Tubing (not illustrated) |




DISPOSABLE SMOKE EVACUATION TUBING

- | | | |
|---------|--|-------|
| 93-2160 | 10" Precut
3/16" Inside Diameter
1/4" Outside Diameter | 50/cs |
|---------|--|-------|



COLLIN VAGINAL SPECULUM

- | | | |
|---------|-----------------|--|
| 90-3946 | 4 1/4" x 1 1/4" | Medium, with permanent Smoke Evacuation Tube (not illustrated) |
|---------|-----------------|--|



AUWARD VAGINAL SPECULUMS

- | | | |
|---------------|-----------------|--|
| MEDIUM | | |
| 90-9100 | 2 3/4" x 1 1/2" | 2 1/2 lbs. without loop |
| 90-9102 | 2 3/4" x 1 1/2" | 2 1/2 lbs. with loop to hold and direct disposable Smoke Evacuation Tubing |
| 90-9103 | 2 3/4" x 1 1/2" | 2 1/2 lbs. with permanent Smoke Evacuation Tube |

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)


SPECULA

OB/GYN

Sklar Black™ Laser



GRAVES VAGINAL SPECULUMS		
WITH LOOPS FOR DISPOSABLE SMOKE EVACUATION TUBING		
90-9005	3" x 1½"	Extra Small
90-9000	3" x ¾"	Small
90-9002	4" x 1¼"	Medium
90-9004	4½" x 1⅝"	Large
90-9006	4¾" x 1½"	Extra Large



GRAVES VAGINAL SPECULUMS		
WITH PERMANENT SMOKE EVACUATION TUBE		
90-9015	3" x ¾"	Small
90-9016	4" x 1¼"	Medium
90-9018	4½" x 1⅝"	Large
90-9020	4¾" x 1½"	Extra Large



GRAVES VAGINAL SPECULUMS		
WITHOUT TUBE OR LOOPS		
90-3812	3" x ½"	Extra Small
90-3804	3" x ¾"	Small
90-3806	4" x 1¼"	Medium
90-3808	4½" x 1⅝"	Large
90-3810	4¾" x 1½"	Extra Large



PEDERSON VAGINAL SPECULUMS		
WITH LOOPS FOR DISPOSABLE SMOKE EVACUATION TUBING		
90-9035	2¾" x ½"	Extra Small
90-9030	3" x ½"	Small
90-9038	4" x ⅝"	Slender
90-9032	4" x ⅞"	Medium
90-9034	4¾" x 1"	Large




PEDERSON VAGINAL SPECULUMS		
WITH PERMANENT SMOKE EVACUATION TUBE		
90-9040	3" x ½"	Small
90-9050	4" x ⅝"	Slender
90-9042	4" x ⅞"	Medium
90-9044	4¾" x 1"	Large
90-9048	4¾" x 1⅝"	Extra Large



PEDERSON VAGINAL SPECULUMS		
WITHOUT TUBE OR LOOPS		
90-3841	2¾" x ½"	Extra Small
90-3833	3" x ½"	Small
90-3837	4" x ⅝"	Slender
90-3839	4" x ⅞"	Medium
90-3835	4¾" x 1"	Large

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



SMOKE EVACUATOR TUBING

ITEM #	DESCRIPTION
H-786T	7/8" x 6" TUBE, NON STERILE, BOX OF 24
H-786TS	7/8" x 6" TUBE, STERILE, BOX OF 24



REDUCER

ITEM #	DESCRIPTION
HNM-36-150722	7/8", TO PRE-FILTER OR CANISTER (FLUTED TUBING) 10/PK



CYTOLOGY BRUSH

ITEM #	DESCRIPTION
HBI-TEAR1011	100/PK



COTTLE TENACULUM HOOK

ITEM #	DESCRIPTION
HNM-1-20-4016	6 1/4", ANG



EMMETT TENACULUM HOOK

ITEM #	DESCRIPTION
HNM-1-20-4201	8 3/4", #1
HNM-1-20-4202	#2
HNM-1-20-4203	8 3/4", #3
HNM-1-20-4204	#4
HNM-1-20-4205	8 3/4", #5

ASPIRATION SET



ASPIRATION SET W/LOCKING 60CC SYRINGE AND CANNULA

ITEM #	DESCRIPTION
HBI-R57-980	3mm, FLEXIBLE CURETTE
HBI-R57-981	4mm, FLEXIBLE CURETTE
HBI-R57-982	5mm, FLEXIBLE CURETTE
HBI-R57-983	6mm, FLEXIBLE CURETTE
HBI-R57-985	7mm, FLEXIBLE CURETTE

GYN SUPPLIES

JACKSON TRACHEAL TENACULUM

ITEM #	DESCRIPTION
HNM-1-20-4014	5 3/4"



UTERINE SOUND



UTERINE SOUND

ITEM #	DESCRIPTION
HBI-BA-100	9 3/4", STR, SINGLE USE, STERILE, BOX OF 25






SIMS UTERINE CURETTE

ITEM #	DESCRIPTION
HNM-1-52-8150	11", #0, MALLEABLE SHAFT, SHARP
HNM-1-52-8151	11", #1, MALLEABLE SHAFT, SHARP
HNM-1-52-8152	11", #2, MALLEABLE SHAFT, SHARP
HNM-1-52-8153	11", #3, MALLEABLE SHAFT, SHARP
HNM-1-52-8154	11", #4, MALLEABLE SHAFT, SHARP
HNM-1-52-8155	11", #5, MALLEABLE SHAFT, SHARP
HNM-1-52-8156	11", #6, MALLEABLE SHAFT, SHARP
HNM-1-52-8150S	11", #0, MALLEABLE SHAFT, SHARP, STND GRADE
HNM-1-52-8152S	11", #2, MALLEABLE SHAFT, SHARP, STND GRADE
HNM-1-52-8156S	11", #6, MALLEABLE SHAFT, SHARP, STND GRADE

SIZE CHART




CURETTES

KEVORKIAN YOUNGE BIOPSY CURETTE

ITEM #	DESCRIPTION
HNM-1-52-7010	W/O BASKET, 11 3/4"
HNM-1-52-7015	W/BASKET, 12 1/8"

NOVAK ENDOMETRIUM BIOPSY SUCTION CURETTE


ITEM #	DESCRIPTION
HNM-1-52-7048	9", DIA 2MM
HNM-1-52-7050	9 7/8", DIA 4MM, W/LUER LOCK HUB AND TEETH



AUVARD VAGINAL SPECULUM

ITEM #	DESCRIPTION
HNM-1-52-2786	2 3/4" x 1 1/2" BLADES, 2LB
HNM-1-52-2787	2 3/4" x 1 1/2" BLADES, 2.5LB
HNM-1-52-2788	2 3/4" x 1 1/2" BLADES, 3LB
HNM-1-52-2793	4" x 1 3/4" BLADES, 2.5LB

ENDOCERVICAL SPECULUM



KOGAN ENDOCERVICAL

ITEM #	DESCRIPTION
HNM-1-52-3955	W/GAUGE AND SET SCREW, 9"
HNM-1-52-3950	W/RACHET, 9"



WEISMAN GRAVES VAGINAL SPECULUM

ITEM #	DESCRIPTION
HNM-1-52-0769	SMALL, 3 1/4" x 1/4" BLADES, RIGHT SIDE OPENING
HNM-1-52-0400	LEFT SIDE OPENING
HNM-36-150759	MEDIUM, RIGHT SIDE OPENING
HNM-36-150760	MEDIUM, 1 1/4" x 4" BLADES, LEFT SIDE OPENING
HNM-1-52-0405	LARGE, LEFT SIDE OPENING
HNM-1-52-0770	LARGE, 4 3/4" x 1 1/2" BLADES, RIGHT SIDE OPENING
HNM-36-150761	EXTRA LARGE, RIGHT SIDE OPENING
HNM-36-150762	EXTRA LARGE, LEFT SIDE OPENING

WEISMAN PEDERSON VAGINAL SPECULUM

ITEM #	DESCRIPTION
HNM-36-150763	MEDIUM, RIGHT SIDE OPENING
HNM-36-150764	MEDIUM, LEFT SIDE OPENING
HNM-36-150766	LARGE, RIGHT SIDE OPENING
HNM-36-150765	LARGE, LEFT SIDE OPENING


VAGINAL SPECULUM

GRAVES VAGINAL SPECULUM


ITEM #	DESCRIPTION
HNM-1-52-0607	SMALL, 3/4" x 3" BLADES
HNM-1-52-0610	MEDIUM, 1 1/4" x 4" BLADES
HNM-1-52-0612	LARGE, 1 3/8" x 4 1/2" BLADES
HNM-36-150755	EXTRA LARGE, 1 1/2" x 6 3/4" BLADES
HNM-36-150756	IMPROVED, MEDIUM, 1 3/8" x 4" BLADES, WIDE ANG

PEDERSON VAGINAL SPECULUM

ITEM #	DESCRIPTION
HNM-1-52-0507	SMALL, 1/2" x 3" BLADES
HNM-1-52-0510	MEDIUM, 7/8" x 4" BLADES
HNM-1-52-0512	LARGE, 1" x 4 3/4" BLADES
HNM-36-151353	VIRGIN SIZE, 5/8" x 2 3/4" BLADES, EXTRA NARROW



GRAVES



PEDERSON




TISCHLER-KEVORKIAN BIOPSY PUNCH

ITEM #	DESCRIPTION
HNM-36-150754	8", 3mm x 9.5mm BITE, SPRING HANDLE




KEVORKIAN-PACIFIC BIOPSY PUNCH

ITEM #	DESCRIPTION
H-233-5365	8", 3mm x 8mm BITE, SPRING HANDLE




BIOPSY PUNCHES



TISCHLER-MORGAN BIOPSY PUNCH

ITEM #	DESCRIPTION
HNM-36-150753	8", 3mm x 7mm BITE, SPRING HANDLE
HNM-1-52-6556	10", 3mm x 7mm BITE, SPRING HANDLE



BABY TISCHLER BIOPSY PUNCH

ITEM #	DESCRIPTION
HNM-36-150750	8", 2mm x 4mm BITE, SPRING HANDLE
HNM-36-150752	10", 2mm x 4mm BITE, SPRING HANDLE



BRAUN TENACULUM FORCEPS

ITEM #	DESCRIPTION
HNM-1-52-4426	10", STR



LAUFE UTERINE POLYP FORCEPS

ITEM #	DESCRIPTION
HNM-1-30-5640	8"




FORCEPS



SCHROEDER TENACULUM FORCEPS

ITEM #	DESCRIPTION
HNM-36-151325	10", TIPS OVERLAP
HNM-36-151324	10", TIPS DO NOT OVERLAP
HNM-1-52-4425	10"
HNM-1-52-4425S	10", STND GRADE





DISPOSABLE SQUARE LOOP ELECTRODES

ITEM #	DESCRIPTION
H-ES14	5mm x 5mm, STERILE, BOX OF 5
H-ES16	10mm x 8mm, STERILE, BOX OF 5
H-OB1010S	10mm x 10mm, STERILE, BOX OF 5




DISPOSABLE BALL ELECTRODES

ITEM #	DESCRIPTION
H-OB3B	3mm, STERILE, BOX OF 5
H-ES07	5mm, STERILE, BOX OF 5




ELECTRODES



DISPOSABLE LOOP ELECTRODES


ITEM #	DESCRIPTION
H-ES08	5mm x 5mm LOOP, STERILE, BOX OF 5
H-ES46	10mm x 5mm, STERILE, BOX OF 5
H-OB-1007	10mm x 7mm, STERILE, BOX OF 5
H-ES52	10mm x 8mm, STERILE, BOX OF 5
H-ES09	10mm x 10mm LOOP, STERILE, BOX OF 5
H-OB-1507	15mm x 7mm, STERILE, BOX OF 5
H-ES11	15mm x 10mm LOOP, STERILE, BOX OF 5
H-ES53	15mm x 12mm, STERILE, BOX OF 5
H-ES43	15mm x 15mm, STERILE, BOX OF 5
H-ES12	20mm x 8mm LOOP, STERILE, BOX OF 5
H-ES49	20mm x 10mm, STERILE, BOX OF 5
H-ES42	20mm x 12mm, STERILE, BOX OF 5
H-ES13	20mm x 15mm LOOP, STERILE, BOX OF 5
H-ES31	20mm x 20mm LOOP, STERILE, BOX OF 5
H-ES47	25mm x 10mm LOOP, STERILE, BOX OF 5

ENDOMETRIAL SAMPLERS



ENDOMETRIAL CURETTE

ITEM #	DESCRIPTION
H-M0014	PREFERRED FIRM DISPOSABLE ENDOMETRIAL CURETTE, 3MM, BOX OF 25
H-M0015	PREFERRED SUREFLEX DISPOSABLE ENDOMETRIAL CURETTE, 3MM, BOX OF 25
H-M0016	SUREFLEX EZ DISPOSABLE ENDOMETRIAL CURETTE, 3.6MM, BOX OF 25



CURETTE ENDOMETRIAL SAMPLER

ITEM #	DESCRIPTION
HBI-RD003	W/TWIST AND LOCK SYRINGE 3mm, BOX OF 25



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



UTERINE INTRODUCING CATHETER

ITEM #	DESCRIPTION
HBI-UIC5F	NON BALLOON SONOHYSTEROGRAPHY CATHETER, BOX OF 25



UTERINE MANIPULATOR/INJECTOR


ITEM #	DESCRIPTION
HUA-1151	4.5MM, BOX OF 12
HNM-36-151532	ADJUSTABLE, BOX OF 12

UTERINE INJECTOR

ITEM #	DESCRIPTION
HNM-36-151530	2mm, BOX OF 12
HNM-36-151531	4mm, BOX OF 12




CATHETERS



WORD BARTHOLIN CATHETER SET

ITEM #	DESCRIPTION
HBI-WC400	W/SYRINGE AND SCALPEL, 6 SETS/PK



HSG CATHETER SETS

ITEM #	DESCRIPTION
HBI-BA05F	W/BALLOON, 5Fr, BOX OF 10
HBI-BA07F	W/BALLOON, 7Fr, BOX OF 10




PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



SILICONE VAGINAL DILATOR SET

ITEM #	DESCRIPTION
HBI-DT-A	LARGE, 4 SIZES, REUSABLE
HBI-DT-B	MEDIUM, 4 SIZES, REUSABLE
HBI-DT-C	SMALL, 4 SIZES, REUSABLE



LAMINARIA NATURAL CERVICAL DILATOR

ITEM #	DESCRIPTION
H-402	2mm, EXTRA SMALL, BOX OF 12
H-403	3mm, SMALL, BOX OF 12
H-404	4mm, MEDIUM, BOX OF 12
H-405	5mm, LARGE, BOX OF 12
H-406	6mm, EXTRA LARGE, BOX OF 12



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)



DENNISTON DILATOR

ITEM #	DESCRIPTION
HNM-36-150523	SOFT PLASTIC, DOUBLE ENDED, SET OF 5



MINI DILATOR SET

ITEM #	DESCRIPTION
HNM-36-151393	(5) PROBES (1mm/1.5mm/2mm/2.5mm/3mm) AND (1) HANDLE




OS FINDER

ITEM #	DESCRIPTION
HBI-OS260BA	OS FINDER COMPLETE SET, SET OF 3 (OS/CANAL/FUNDUS), REUSABLE
HBI-OS265BA-1	DISPOSABLE OS FINDER, BOX OF 25
HBI-OS265BA-2	DISPOSABLE CANAL FINDER, BOX OF 25
HBI-OS265BA-3	DISPOSABLE FUNDUS FINDER, BOX OF 25



DILATORS



HANK UTERINE DILATOR SET

ITEM #	DESCRIPTION
HNM-1-52-3742	11", DOUBLE ENDED, SET OF 6

**Now also available in individual*



PRATT UTERINE DILATOR SET


ITEM #	DESCRIPTION
HNM-1-52-3642	11 1/2", DOUBLE ENDED, SET OF 8

**Now also available in individual*



SMALL CERVICAL DILATOR

ITEM #	DESCRIPTION
HNM-36-150522	7 1/2", SINGLE ENDED, SET OF 6



HEGAR DILATOR SET

ITEM #	DESCRIPTION
HNM-1-52-3542	DOUBLE ENDED, 7 1/2", SET OF 8 (3mm/4mm - 17mm/18mm)

**Now also available in individual*

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

CUBE W/O DRAINS



ITEM #	DESCRIPTION
HNM-36-150682	CUBE #0, W/O DRAINS, 25MM
HNM-36-150683	CUBE #1, W/O DRAINS, 29MM
HNM-36-150684	CUBE #2, W/O DRAINS, 33MM
HNM-36-150685	CUBE #3, W/O DRAINS, 37MM
HNM-36-150686	CUBE #4, W/O DRAINS, 41MM
HNM-36-150687	CUBE #5, W/O DRAINS, 45MM
HNM-36-150688	CUBE #6, W/O DRAINS, 50MM
HNM-36-150689	CUBE #7, W/O DRAINS, 56MM
HNM-36-150690	CUBE #8, W/O DRAINS, 63MM
HNM-36-150691	CUBE #9, W/O DRAINS, 70MM
HNM-36-150692	CUBE #10, W/O DRAINS, 75MM

CUBE W/DRAINS



ITEM #	DESCRIPTION
HNM-36-150671	CUBE #0, W/DRAINS, 25MM
HNM-36-150672	CUBE #1, W/DRAINS, 29MM
HNM-36-150673	CUBE #2, W/DRAINS, 33MM
HNM-36-150674	CUBE #3, W/DRAINS, 37MM
HNM-36-150675	CUBE #4, W/DRAINS, 41MM
HNM-36-150676	CUBE #5, W/DRAINS, 45MM
HNM-36-150677	CUBE #6, W/DRAINS, 50MM
HNM-36-150678	CUBE #7, W/DRAINS, 56MM
HNM-36-150679	CUBE #8, W/DRAINS, 63MM
HNM-36-150680	CUBE #9, W/DRAINS, 70MM
HNM-36-150681	CUBE #10 W/DRAINS, 75MM



INFLATABLE PESSARY

ITEM #	DESCRIPTION
HNM-36-150713	INFLATABLE 2", SMALL
HNM-36-150714	INFLATABLE 2.25", MEDIUM
HNM-36-150715	INFLATABLE 2.5", LARGE
HNM-36-150716	INFLATABLE 2.75", XL

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

GEHRUNG



ITEM #	DESCRIPTION
HNM-36-150661	GEHRUNG #0, 50MM
HNM-36-150662	GEHRUNG #1, 55MM
HNM-36-150663	GEHRUNG #2, 60MM
HNM-36-150664	GEHRUNG #3, 65MM
HNM-36-150665	GEHRUNG #4, 70MM
HNM-36-150666	GEHRUNG #5, 75MM
HNM-36-150667	GEHRUNG #6, 80MM
HNM-36-150668	GEHRUNG #7, 85MM
HNM-36-150669	GEHRUNG #8, 90MM
HNM-36-150670	GEHRUNG #9, 95MM

HODGE W/O SUPPORT



ITEM #	DESCRIPTION
HNM-36-150703	HODGE #0, W/O SUPPORT, 65MM
HNM-36-150704	HODGE #1, W/O SUPPORT, 70MM
HNM-36-150705	HODGE #2, W/O SUPPORT, 75MM
HNM-36-150706	HODGE #3, W/O SUPPORT, 80MM
HNM-36-150707	HODGE #4, W/O SUPPORT, 85MM
HNM-36-150708	HODGE #5, W/O SUPPORT, 90MM
HNM-36-150709	HODGE #6, W/O SUPPORT, 95MM
HNM-36-150710	HODGE #7, W/O SUPPORT, 100MM
HNM-36-150711	HODGE #8, W/O SUPPORT, 105MM
HNM-36-150712	HODGE #9, W/O SUPPORT, 110MM

HODGE W/SUPPORT



ITEM #	DESCRIPTION
HNM-36-150693	HODGE #0, W/SUPPORT, 65MM
HNM-36-150694	HODGE #1, W/SUPPORT, 70MM
HNM-36-150695	HODGE #2, W/SUPPORT, 75MM
HNM-36-150696	HODGE #3, W/SUPPORT, 80MM
HNM-36-150697	HODGE #4, W/SUPPORT, 85MM
HNM-36-150698	HODGE #5, W/SUPPORT, 90MM
HNM-36-150699	HODGE #6, W/SUPPORT, 95MM
HNM-36-150700	HODGE #7, W/SUPPORT, 100MM
HNM-36-150701	HODGE #8, W/SUPPORT, 105MM
HNM-36-150702	HODGE #9, W/SUPPORT, 110MM

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

MARLAND W/SUPPORT



ITEM #	DESCRIPTION
HNM-36-150647	MARLAND #2, W/SUPPORT, 2.25"
HNM-36-150648	MARLAND #3, W/SUPPORT, 2.5"
HNM-36-150649	MARLAND #4, W/SUPPORT, 2.75"
HNM-36-150650	MARLAND #5, W/SUPPORT, 3"
HNM-36-150651	MARLAND #6, W/SUPPORT, 3.25"
HNM-36-150652	MARLAND #7, W/SUPPORT, 3.5"
HNM-36-150653	MARLAND #8, W/SUPPORT, 3.75"

CUP W/O SUPPORT



ITEM #	DESCRIPTION
HNM-36-150639	CUP #0, W/O SUPPORT, 50/18MM
HNM-36-150640	CUP #1, W/O SUPPORT, 55/23MM
HNM-36-150641	CUP #2, W/O SUPPORT, 60/28MM
HNM-36-150642	CUP #3, W/O SUPPORT, 65/33MM
HNM-36-150643	CUP #4, W/O SUPPORT, 70/38MM
HNM-36-150644	CUP #5, W/O SUPPORT, 75/43MM
HNM-36-150645	CUP #6, W/O SUPPORT, 80/48MM
HNM-36-150646	CUP #7, W/O SUPPORT, 85/52MM

CUP W/SUPPORT



ITEM #	DESCRIPTION
HNM-36-150631	CUP #0, W/SUPPORT, 50/18MM
HNM-36-150632	CUP #1, W/SUPPORT, 55/23MM
HNM-36-150633	CUP #2, W/SUPPORT, 60/28MM
HNM-36-150634	CUP #3, W/SUPPORT, 65/33MM
HNM-36-150635	CUP #4, W/SUPPORT, 70/38MM
HNM-36-150636	CUP #5, W/SUPPORT, 75/43MM
HNM-36-150637	CUP #6, W/SUPPORT, 80/48MM
HNM-36-150638	CUP #7, W/SUPPORT, 85/52MM



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

DISH W/O SUPPORT



ITEM #	DESCRIPTION
HNM-36-150613	DISH #0, W/O SUPPORT, 50MM
HNM-36-150614	DISH #1, W/O SUPPORT, 55MM
HNM-36-150615	DISH #2, W/O SUPPORT, 60MM
HNM-36-150616	DISH #3, W/O SUPPORT, 65MM
HNM-36-150617	DISH #4, W/O SUPPORT, 70MM
HNM-36-150618	DISH #5, W/O SUPPORT, 75MM
HNM-36-150619	DISH #6, W/O SUPPORT, 80MM
HNM-36-150620	DISH #7, W/O SUPPORT, 85MM

DISH W/SUPPORT



ITEM #	DESCRIPTION
HNM-36-150621	DISH #0, W/SUPPORT, 50MM
HNM-36-150622	DISH #1, W/SUPPORT, 55MM
HNM-36-150623	DISH #2, W/SUPPORT, 60MM
HNM-36-150624	DISH #3, W/SUPPORT, 65MM
HNM-36-150625	DISH #4, W/SUPPORT, 70MM
HNM-36-150626	DISH #5, W/SUPPORT, 75MM
HNM-36-150627	DISH #6, W/SUPPORT, 80MM
HNM-36-150628	DISH #7, W/SUPPORT, 85MM

MARLAND W/O SUPPORT



ITEM #	DESCRIPTION
HNM-36-150654	MARLAND #2, W/O SUPPORT, 2.25"
HNM-36-150655	MARLAND #3, W/O SUPPORT, 2.5"
HNM-36-150656	MARLAND #4, W/O SUPPORT, 2.75"
HNM-36-150657	MARLAND #5, W/O SUPPORT, 3"
HNM-36-150658	MARLAND #6, W/O SUPPORT, 3.25"
HNM-36-150659	MARLAND #7, W/O SUPPORT, 3.5"
HNM-36-150660	MARLAND #8, W/O SUPPORT, 3.75"



GELLHORN REGULAR STEM (1.5" STEM)



ITEM #	DESCRIPTION
HNM-36-150578	GELLHORN #0, 1.50"
HNM-36-150579	GELLHORN #1, 1.75"
HNM-36-150580	GELLHORN #2, 2"
HNM-36-150581	GELLHORN #3, 2.25"
HNM-36-150582	GELLHORN #4, 2.5"
HNM-36-150583	GELLHORN #5, 2.75"
HNM-36-150584	GELLHORN #6, 3"
HNM-36-150585	GELLHORN #7, 3.25"
HNM-36-150586	GELLHORN #8, 3.5"

GELLHORN SHORT STEM (1.1" STEM)



ITEM #	DESCRIPTION
HNM-36-150587	GELLHORN SHORT STEM #0, 1.50"
HNM-36-150588	GELLHORN SHORT STEM #1, 1.75"
HNM-36-150589	GELLHORN SHORT STEM #2, 2"
HNM-36-150590	GELLHORN SHORT STEM #3, 2.25"
HNM-36-150591	GELLHORN SHORT STEM #4, 2.50"
HNM-36-150592	GELLHORN SHORT STEM #5, 2.75"
HNM-36-150593	GELLHORN SHORT STEM #6, 3"
HNM-36-150594	GELLHORN SHORT STEM #7, 3.25"
HNM-36-150595	GELLHORN SHORT STEM #8, 3.50"

SHAATZ



ITEM #	DESCRIPTION
HNM-36-150596	SHAATZ #0, 1.5"
HNM-36-150597	SHAATZ #1, 1.75"
HNM-36-150598	SHAATZ #2, 2"
HNM-36-150599	SHAATZ #3, 2.25"
HNM-36-150600	SHAATZ #4, 2.5"
HNM-36-150601	SHAATZ #5, 2.75"
HNM-36-150602	SHAATZ #6, 3"
HNM-36-150603	SHAATZ #7, 3.25"
HNM-36-150604	SHAATZ #8, 3.5"

RING W/O SUPPORT



ITEM #	DESCRIPTION
HNM-36-150534	RING #0, W/O SUPPORT, 1.75"
HNM-36-150535	RING #1, W/O SUPPORT, 2"
HNM-36-150536	RING #2, W/O SUPPORT, 2.25"
HNM-36-150537	RING #3, W/O SUPPORT, 2.5"
HNM-36-150538	RING #4, W/O SUPPORT, 2.75"
HNM-36-150539	RING #5, W/O SUPPORT, 3"
HNM-36-150540	RING #6, W/O SUPPORT, 3.25"
HNM-36-150541	RING #7, W/O SUPPORT, 3.5"
HNM-36-150542	RING #8, W/O SUPPORT, 3.75"
HNM-36-150543	RING #9, W/O SUPPORT, 4"

OVAL W/O SUPPORT



ITEM #	DESCRIPTION
HNM-36-150569	OVAL #1, W/O SUPPORT, 2.00"
HNM-36-150570	OVAL #2, W/O SUPPORT, 2.25"
HNM-36-150571	OVAL #3, W/O SUPPORT, 2.5"
HNM-36-150572	OVAL #4, W/O SUPPORT, 2.75"
HNM-36-150573	OVAL #5, W/O SUPPORT, 3"
HNM-36-150574	OVAL #6, W/O SUPPORT, 3.25"
HNM-36-150575	OVAL #7, W/O SUPPORT, 3.5"
HNM-36-150576	OVAL #8, W/O SUPPORT, 3.75"
HNM-36-150577	OVAL #9, W/O SUPPORT, 4"

OVAL W/SUPPORT



ITEM #	DESCRIPTION
HNM-36-150560	OVAL #1, W/SUPPORT, 2"
HNM-36-150561	OVAL #2, W/SUPPORT, 2.25"
HNM-36-150562	OVAL #3, W/SUPPORT, 2.5"
HNM-36-150563	OVAL #4, W/SUPPORT, 2.75"
HNM-36-150564	OVAL #5, W/SUPPORT, 3"
HNM-36-150565	OVAL #6, W/SUPPORT, 3.25"
HNM-36-150566	OVAL #7, W/SUPPORT, 3.5"
HNM-36-150567	OVAL #8, W/SUPPORT, 3.75"
HNM-36-150568	OVAL #9, W/SUPPORT, 4"

RING W/KNOB AND W/SUPPORT



ITEM #	DESCRIPTION
HNM-36-150544	RING #0, W/KNOB AND W/SUPPORT, 1.75"
HNM-36-150545	RING #1, W/KNOB AND W/SUPPORT, 2"
HNM-36-150546	RING #2, W/KNOB AND W/SUPPORT, 2.25"
HNM-36-150547	RING #3, W/KNOB AND W/SUPPORT, 2.5"
HNM-36-150548	RING #4, W/KNOB AND W/SUPPORT, 2.75"
HNM-36-150549	RING #5, W/KNOB AND W/SUPPORT, 3"
HNM-36-150550	RING #6, W/KNOB AND W/SUPPORT, 3.25"
HNM-36-150551	RING #7, W/KNOB AND W/SUPPORT, 3.5"

RING W/KNOB AND W/O SUPPORT



ITEM #	DESCRIPTION
HNM-36-150552	RING #0, W/KNOB AND W/O SUPPORT, 1.75"
HNM-36-150553	RING #1, W/KNOB AND W/O SUPPORT, 2"
HNM-36-150554	RING #2, W/KNOB AND W/O SUPPORT, 2.25"
HNM-36-150555	RING #3, W/KNOB AND W/O SUPPORT, 2.5"
HNM-36-150556	RING #4, W/KNOB AND W/O SUPPORT, 2.75"
HNM-36-150557	RING #5, W/KNOB AND W/O SUPPORT, 3"
HNM-36-150558	RING #6, W/KNOB AND W/O SUPPORT, 3.25"
HNM-36-150559	RING #7, W/KNOB AND W/O SUPPORT, 3.5"

RING W/SUPPORT



ITEM #	DESCRIPTION
HNM-36-150524	RING #0, W/SUPPORT, 1.75"
HNM-36-150525	RING #1, W/SUPPORT, 2"
HNM-36-150526	RING #2, W/SUPPORT, 2.25"
HNM-36-150527	RING #3, W/SUPPORT, 2.5"
HNM-36-150528	RING #4, W/SUPPORT, 2.75"
HNM-36-150529	RING #5, W/SUPPORT, 3"
HNM-36-150530	RING #6, W/SUPPORT, 3.25"
HNM-36-150531	RING #7, W/SUPPORT, 3.5"
HNM-36-150532	RING #8, W/SUPPORT, 3.75"
HNM-36-150533	RING #9, W/SUPPORT, 4"

PESSARY

FITTING SETS

ITEM #	DESCRIPTION
HNM-36-150718	RING WITHOUT SUPPORT FITTING SET (sizes #2-7)
HNM-36-150517	RING WITH SUPPORT FITTING SET (sizes #2-7)
HNM-36-150518	OVAL FITTING SET (sizes #2-7)
HNM-36-150519	CUBE FITTING SET (sizes #2-7)
HNM-36-150520	DISH FITTING SET (sizes #2-7)
HNM-36-150521	GELLHORN FITTING SET (sizes #2-7)



PESSARY REMOVER

ITEM #	DESCRIPTION
HNM-36-150717	PESSARY REMOVER



DONUT



ITEM #	DESCRIPTION
HNM-36-150605	DONUT #0, 2"
HNM-36-150606	DONUT #1, 2.25"
HNM-36-150607	DONUT #2, 2.50"
HNM-36-150608	DONUT #3, 2.75"
HNM-36-150609	DONUT #4, 3"
HNM-36-150610	DONUT #5, 3.25"
HNM-36-150611	DONUT #6, 3.50"
HNM-36-150612	DONUT #7, 3.75"

BUILD YOUR OWN COLPOSCOPY SET

1. Biopsy Punches

Tischler Morgan, Spring Handle, 10"	HNM-1-52-6556
Tischler Kevorikan, Rotating Handle, 10"	HNM-1-52-6541
Baby Tischler, Spring Handle, 8"	HNM-36-150750

2. Endocervical Speculum

Kogan Std Jaws w/Ratchet, 9"	HNM-1-52-3950
Kogan Std Jaws w/Gauge & Set Screw, 9"	HNM-1-52-3955

3. Vaginal Speculum

Graves, Std Opening, Small.....	HNM-1-52-0607
Graves, Std Opening, Medium.....	HNM-1-52-0610
Graves, Std Opening, Large.....	HNM-1-52-0612
Pederson, Std Openinig, Small.....	HNM-1-52-0507
Pederson, Std Openinig, Medium.....	HNM-1-52-0510
Pederson, Std Openinig, Large.....	HNM-1-52-0512

4. Endocervical Curette

Kevorkian Young w/o Basket, 12"	HNM-1-52-7010
Kevorkian Young w/Basket, 12"	HNM-1-52-7015

5. Emmett Tenaculum Hook

Emmett Tenaculum Hook, #1, 8 3/4"	HNM-1-20-4201
Emmett Tenaculum Hook, #2	HNM-1-20-4202
Emmett Tenaculum Hook, #3, 8 3/4"	HNM-1-20-4203
Emmett Tenaculum Hook, #4	HNM-1-20-4204
Emmett Tenaculum Hook, #5, 8 3/4"	HNM-1-20-4205





MedGyn Needle Adapters - Extenders

An effective means of extending ordinary hypodermic needles when longer needles are required. The paracervical and uterosacral areas can easily be reached with an extended needle.

Needle Adapters

030501 Needle Adapter 3" (76mm)

030502 Needle Adapter 4" (102mm)

030503 Needle Adapter 5" (127mm)

030504 Needle Adapter 6" (152mm)



Miscellaneous Products

040004 Emesis Basin Medium 8" (203mm)

040005 Emesis Basin Large 10" (254mm)

040003 Emesis Basin Small 6" (152mm)

040008 Iodine Cup 6oz

040001 Medicine Cup (2 oz.)

040009 Instrument Cleaner 8oz Bottle

040000 Instrument Lubricant Spray - 8oz

040010 Instrument Stain Remover 3oz

040002 Instrument Steel Tray - 12¼" X 7⅞" X 2⅛"
(Graham 3259)



NEEDLE HOLDERS



MedGyn Needle Holders

Webster Needle Holder

040075 Webster Needle Holder 5" (127mm) Serrated Jaws

040076 Webster Needle Holder 5" (127mm) Smooth Jaws



Mayo-Hegar Needle Holder

040093 Mayo-Hegar Needle Holder 5" (127mm) Serrated

040094 Mayo-Hegar Needle Holder 6" (152mm) Serrated

040095 Mayo-Hegar Needle Holder 7" (178mm) Serrated

040096 Mayo-Hegar Needle Holder 8" (203mm) Serrated



Olsen-Hegar Needle Holder/Scissors

040097 Olsen-Hegar Needle Holder/Scissors 4¾" (121mm) Serrated Extra Delicate

040098 Olsen-Hegar Needle Holder/Scissors 4¾" (121mm) Smooth Extra Delicate

040099 Olsen-Hegar Needle Holder/Scissors 5½" (140mm) Serrated

040100 Olsen-Hegar Needle Holder/Scissors 6½" (165mm) Serrated

040101 Olsen-Hegar Needle Holder/Scissors 7¼" (184mm) Serrated



Heaney Needle Holder

040071 Heaney Needle Holder 8½" (216mm) Curved Jaws




Mathieu Needle Holder

040072 Mathieu Needle Holder 5½" (140mm) w/Spring Handles

040073 Mathieu Needle Holder 6¼" (159mm) w/Spring Handles

040074 Mathieu Needle Holder 7½" (191mm) w/Spring Handles



Crile-Wood Needle Holder

040081 Crile-Wood Needle Holder 6" (152mm)

040082 Crile-Wood Needle Holder 7" (178mm)

040083 Crile-Wood Needle Holder 8" (203mm)

040084 Crile-Wood Needle Holder 9" (229mm)

040087 Crile-Wood Needle Holder 10½" (267mm)



Potts-Smith Scissors

- 041015 Potts-Smith Scissors 7" (178mm), 13mm Blades, Delicate - Angled 25°
- 041016 Potts-Smith Scissors 7" (178mm), 13mm Blades, Delicate - Angled 45°
- 041017 Potts-Smith Scissors 7" (178mm), 13mm Blades, Delicate - Angled 60°
- 041018 Potts-Smith Scissors 7½" (191mm), 20mm Blades, Delicate - Angled 25°
- 041019 Potts-Smith Scissors 7½" (191mm), 20mm Blades, Delicate - Angled 45°
- 041020 Potts-Smith Scissors 7½" (191mm), 20mm Blades, Delicate - Angled 60°



Lister Scissors

- 041050 Lister Bandage Scissors 3½" (89mm) Extra Fine
- 041051 Lister Bandage Scissors 4½" (114mm) Extra Fine
- 041052 Lister Bandage Scissors 5½" (140mm) Extra Fine
- 041053 Lister Bandage Scissors 7¼" (184mm) Extra Fine
- 041054 Lister Bandage Scissors 8" (203mm) with One Larger Finger Ring



Littauer Scissors

- 041088 Littauer Stitch Scissors 4½" (114mm) Light
- 041089 Littauer Stitch Scissors 5½" (140mm) Standard
- 042099-ECO Littauer Stitch Scissor 4½" (114mm) Economy



Sweet Scissors

- 041096 Sweet Scissors 8½" (216mm) Angled on Flat



Spencer Scissors

- 041095 Spencer Stitch Scissors 3½" (89mm) Delicate

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SCISSORS



Doyen Scissors

- 041075 Doyen Scissors 7" (178mm) - Curved/Smooth
- 041076 Doyen Scissors 7" (178mm) - Straight/Smooth
- 041077 Doyen Scissors 7" (178mm) - Curved/Serrated
- 041078 Doyen Scissors 7" (178mm) - Straight/Serrated



IUD Scissors

- 041010 IUD Scissors Curved 10 (254mm)



Mayo Scissors

- 041061 Straight Mayo Scissors 5½" (140mm)
- 041063 Straight Mayo Scissors 6¾" (171mm)
- 041070 Straight Mayo Scissor 8" (203mm)
- 041067 Straight Mayo Scissors 9" (229mm)
- 041071 Straight Mayo Scissor 10¼"(260mm)
- 041069 Straight Mayo Scissors 11" (279mm)
- 041060 Curved Mayo Scissors 5½" (140mm)
- 041062 Curved Mayo Scissors 6¾" (171mm)
- 041064 Curved Mayo Scissors 8" (203mm)
- 041066 Curved Mayo Scissors 9" (229mm)
- 041068 Curved Mayo Scissors 11" (279mm)



Iris Scissors

- 041100 Iris Scissors 3½" (89mm), 20mm Blades Delicate Curved
- 041101 Iris Scissors 3½" (89mm), 20mm Blades Delicate Straight
- 041102 Iris Scissors 3¾" (95mm), 14mm Blades Curved Ribbon Type
- 041103 Iris Scissors 3¾" (95mm), 14mm Blades Straight Ribbon Type
- 041104 Iris Scissors 4½" (114mm) Angled on Side
- 041105 Iris Scissors 4½" (114mm) Curved
- 041106 Iris Scissors 4½" (114mm) Curved Sideways
- 041107 Iris Scissors 4½" (114mm) Straight



Kelly-Adson Scissors

- 041085 Kelly-Adson 6¼" (159mm) Smooth/Straight
- 041086 Kelly-Adson 6¼" (159mm) Serrated/Straight
- 041087 Kelly-Adson 6¼" (159mm) Smooth/Curved

MedGyn Scissors



Metzenbaum Scissors - Straight

041026 Metzenbaum 5" (127mm) - Standard Straight
 041028 Metzenbaum 6" (152mm) - Standard Straight
 041030 Metzenbaum 7" (178mm) - Standard Straight
 041032 Metzenbaum 8" (203mm) - Standard Straight
 041034 Metzenbaum 9" (229mm) - Standard Straight
 041036 Metzenbaum 10" (254mm) - Standard Straight
 041038 Metzenbaum 11" (279mm) - Standard Straight
 041040 Metzenbaum 12" (305mm) - Standard Straight
 041026-DEL Metzenbaum 5" (127mm) - Delicate Straight
 041028-DEL Metzenbaum 6" (152mm) - Delicate Straight
 041030-DEL Metzenbaum 7" (178mm) - Delicate Straight
 041032-DEL Metzenbaum 8" (203mm) - Delicate Straight
 041034-DEL Metzenbaum 9" (229mm) - Delicate Straight



Metzenbaum Scissors - Curved

041025 Metzenbaum 5" (127mm) - Standard Curved
 041027 Metzenbaum 6" (152mm) - Standard Curved
 041029 Metzenbaum 7" (178mm) - Standard Curved
 041031 Metzenbaum 8" (203mm) - Standard Curved
 041033 Metzenbaum 9" (229mm) - Standard Curved
 041035 Metzenbaum 10" (254mm) - Standard Curved
 041037 Metzenbaum 11" (279mm) - Standard Curved
 041039 Metzenbaum 12" (305mm) - Standard Curved
 041025-DEL Metzenbaum 5" (127mm) - Delicate Curved
 041027-DEL Metzenbaum 6" (152mm) - Delicate Curved
 041029-DEL Metzenbaum 7" (178mm) - Delicate Curved
 041031-DEL Metzenbaum 8" (203mm) - Delicate Curved
 041033-DEL Metzenbaum 9" (229mm) - Delicate Curved



Operating Scissors - Straight

042000-S Straight Scissors 4½" (114mm) - Blunt/Blunt
 042001-S Straight Scissors 4½" (114mm) - Sharp/Blunt
 042002-S Straight Scissors 4½" (114mm) - Sharp/Sharp
 042004-S Straight Scissors 5" (127mm) - Blunt/Blunt
 042003-S Straight Scissors 5" (127mm) - Sharp/Blunt
 042005-S Straight Scissors 5" (127mm) - Sharp/Sharp
 042006-S Straight Scissors 5½" (140mm) - Blunt/Blunt
 042007-S Straight Scissors 5½" (140mm) - Sharp/Blunt
 042008-S Straight Scissors 5½" (140mm) - Sharp/Sharp
 042009-S Straight Scissors 6" (152mm) - Blunt/Blunt
 042010-S Straight Scissors 6" (152mm) - Sharp/Blunt
 042011-S Straight Scissors 6" (152mm) - Sharp/Sharp
 042012-S Straight Scissors 6½" (165mm) - Blunt/Blunt
 042013-S Straight Scissors 6½" (165mm) - Sharp/Blunt
 042014-S Straight Scissors 6½" (165mm) - Sharp/Sharp
 042015-S Straight Scissors 7½" (191mm) - Sharp/Blunt
 042016-S Straight Scissors 7½" (191mm) - Blunt/Blunt



Operating Scissors - Curved

042000-C Curved Scissors 4½" (114mm) - Blunt/Blunt
 042001-C Curved Scissors 4½" (114mm) - Sharp/Blunt
 042002-C Curved Scissors 4½" (114mm) - Sharp/Sharp
 042004-C Curved Scissors 5" (127mm) - Blunt/Blunt
 042003-C Curved Scissors 5" (127mm) - Sharp/Blunt
 042005-C Curved Scissors 5" (127mm) - Sharp/Sharp
 042006-C Curved Scissors - 5½" (140mm) - Blunt/Blunt
 042007-C Curved Scissors 5½" (140mm) - Sharp/Blunt
 042008-C Curved Scissors 5½" (140mm) - Sharp/Sharp
 042009-C Curved Scissors 6" (152mm) - Blunt/Blunt
 042010-C Curved Scissors 6" (152mm) - Sharp/Blunt
 042011-C Curved Scissors 6" (152mm) - Sharp/Sharp
 042012-C Curved Scissors 6½" (165mm) - Blunt/Blunt
 042013-C Curved Scissors 6½" (165mm) - Sharp/Blunt
 042014-C Curved Scissors 6½" (165mm) - Sharp/Sharp
 042015-C Curved Scissors 7½" (191mm) - Sharp/Blunt

HYSTERECTOMY INSTRUMENTS

MedGyn Hysterectomy Clamps

Hysterectomy ("Z") Clamps



040025 Hysterectomy Clamp 8 1/4" (210mm) - Straight

040026 Hysterectomy Clamp 8 1/4" (210mm) - Slight Curve

040027 Hysterectomy Clamp 8 1/4" (210mm) - Full Curve

040028 Hysterectomy Clamp 8 1/4" (210mm) - Angled

040035 Hysterectomy Clamp 9 1/2" (241mm) - Straight

040036 Hysterectomy Clamp 9 1/2" (241mm) - Slight Curve

040037 Hysterectomy Clamp 9 1/2" (241mm) - Strong Curve

040038 Hysterectomy Clamp 9 1/2" (241mm) - Angled

040045 Hysterectomy Clamp 12" (305mm) - Straight

040046 Hysterectomy Clamp 12" (305mm) - Slight Curve

040047 Hysterectomy Clamp 12" (305mm) - Strong Curve

040048 Hysterectomy Clamp 12" (305mm) - Angled

040055 Hysterectomy Clamp 14" (356mm) - Straight

040056 Hysterectomy Clamp 14" (356mm) - Slight Curve

040057 Hysterectomy Clamp 14" (356mm) - Strong Curve

040058 Hysterectomy Clamp 14" (356mm) - Angled

040015 Hysterectomy Clamp 16" (406mm) - Straight

040016 Hysterectomy Clamp 16" (406mm) - Slight Curve

040017 Hysterectomy Clamp 16" (406mm) - Strong Curve

040018 Hysterectomy Clamp 16" (406mm) - Angled

MedGyn Hysterectomy Forceps



Pean Hysterectomy Forceps - Straight

030999 Pean Hysterectomy Forceps Straight 8" (203mm)

Pean Hysterectomy Forceps - Curved

031000 Pean Hysterectomy Forceps Curved 8" (203mm)

MedGyn Hysterectomy Scissors

041125 Hysterectomy Scissors 9" (229mm) - Slight Curve

041126 Hysterectomy Scissors 9" (229mm) - Curve

041127 Hysterectomy Scissors 9" (229mm) - Full Curve

041128 Hysterectomy Scissors 9" (229mm) - Angled

041135 Hysterectomy Scissors 10 1/2" (267mm) - Slight Curve

041136 Hysterectomy Scissors 10 1/2" (267mm) - Curve

041137 Hysterectomy Scissors 10 1/2" (267mm) - Full Curve

041138 Hysterectomy Scissors 10 1/2" (267mm) - Angled

041145 Hysterectomy Scissors 12" (305mm) - Slight Curve

041146 Hysterectomy Scissors 12" (305mm) - Curve

041147 Hysterectomy Scissors 12" (305mm) - Full Curve

041148 Hysterectomy Scissors 12" (305mm) - Angled

041155 Hysterectomy Scissors 14" (356mm) - Slight Curve

041156 Hysterectomy Scissors 14" (356mm) - Curve

041157 Hysterectomy Scissors 14" (356mm) - Full Curve

041158 Hysterectomy Scissors 14" (356mm) - Angled

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

LONG REACH INSTRUMENTS



Rochester Ochsner Forceps

031184 Rochester Ochsner Forceps 12" (305mm)
Curved

031185 Rochester Ochsner Forceps 14" (356mm)
Curved



Clark Hysterectomy Forceps

031193 Clark Hysterectomy Forceps 12" (305mm)

031194 Clark Hysterectomy Forceps 14" (356mm)



Kocher Forceps

031195 Kocher Forceps Curved 14" (356mm)



Jorgenson Scissors

031177 Jorgenson Scissors 14" (356mm) Curved



Metzenbaum Scissors

040091 Metzenbaum Scissors Curved 14"
(356mm)

040092 Metzenbaum Scissors Straight 14"
(356mm)

LONG REACH INSTRUMENTS



MedGyn Long-Reach Instruments

Allis Tissue Forceps

031099 Allis Forceps 12" (305mm) 5x6 Teeth

031100 Allis Forceps 14" (356mm) 5x6 Teeth



Babcock Tissue Forceps

031188 Babcock Tissue Forceps 12" (305mm)

031189 Babcock Tissue Forceps 14" (356mm)

031190 Babcock Tissue Forceps 18" (457mm)



DeBakey Tissue Forceps

031182 DeBakey Tissue Forceps 12" (305mm)

031183 DeBakey Tissue Forceps 14" (356mm)



Gemini Forceps

031192 Gemini Forceps 14" (356mm) Curved



Crile Wood Needle Holder

040085 Crile Wood Needle Holder 12" (305mm)

040086 Crile Wood Needle Holder 14" (356mm)



Bozeman Uterine Dressing Forceps

031161 Bozeman Uterine Dressing 10½" (267mm) Straight

031162 Bozeman Uterine Dressing 10½" (267mm) Curved


031162-ECO Bozeman Uterine Dressing 10½" (267mm) Curved - Economy



Dressing Forceps

031186 Dressing Forceps 8" (203mm)

031187 Dressing Forceps 10" (254mm)



Fletcher Forceps

031014 Fletcher with Oval Jaws 9½" (241mm) Straight



IUD Hooks


For IUD retrieval; available in single and double notch versions.

031306 IUD Hook 10" (254mm) Single Notch

031307 IUD Hook 10½" (265mm) Double Notch

031308 IUD Hook Universal 12½" (318mm)

FORCEPS



MedGyn Polyp Forceps

Polyp Forceps

031121 Polyp Forceps 10" (254mm) Narrow Jaws



Alligator Forceps

Perfect design for retrieval of IUDs, with hinged serrated jaws.

030112 Alligator Forceps 8" (203mm)

030112-11 Alligator Forceps 11" (279mm)



Foerster Dressing Forceps

031035 Foerster Dressing Forceps 9½" (241mm) Straight - Serrated Jaws

031037 Foerster Dressing Forceps 9½" (241mm) Straight - Smooth Jaws

031036 Foerster Dressing Forceps 9½" (241mm) Curved - Serrated Jaws

031038 Foerster Dressing Forceps 9½" (241mm) Curved - Smooth Jaws



Halsted Mosquito Forceps

031113 Halsted Mosquito Hemostatic Forceps Straight 5" (127mm)

031114 Halsted Mosquito Hemostatic Forceps Straight with Delicate Teeth 5" (127mm)

031115 Halsted Mosquito Hemostatic Forceps Straight with 1x2 Teeth 5" (127mm)

031116 Halsted Mosquito Hemostatic Forceps Curved 5" (127mm)

031118 Halsted Mosquito Hemostatic Forceps Curved with Delicate Teeth 5" (127mm)

031119 Halsted Mosquito Hemostatic Forceps Curved with 1x2 Teeth 5" (127mm)

031116-ECO Halsted Mosquito Hemostatic Forceps Economy 5" (127mm)




Baby Allis Tissue Forceps

- 031111 Baby Allis Forceps 5½" (140mm)
4x5 Teeth-Extra Delicate Jaws
- 031112 Baby Allis Forceps 5" (127mm)
4x5 Teeth-Extra Delicate Jaws



Adson Tissue Forceps

- 031196 Adson Tissue Forceps 4¾" (121mm)
1x2 Teeth



Tissue Forceps

- 031180 Tissue Forceps 8" (203mm) Small 1x2 Teeth
- 031181 Tissue Forceps 10" (254mm) Large 1x2 Teeth



MedGyn Tissue Forceps

Allis Tissue Forceps

- 031107 Allis Tissue Forceps 4" (104mm) 4x5 Teeth
- 031105 Allis Tissue Forceps 5" (127mm) 4x5 Teeth
- 031108 Allis Tissue Forceps 6" (152mm) 4x5 Teeth
- 031109 Allis Tissue Forceps 6" (152mm) 5x6 Teeth
- 031110 Allis Tissue Forceps 7½" (190mm) 5x6 Teeth
- 031103 Allis Tissue Forceps 9½" (241mm) 5x6 Teeth
- 031106 Allis Tissue Forceps 10" (254mm) 5x6 Teeth
- 031104 Allis Tissue Forceps 6" (152mm) 4x5 Teeth
with Angular Jaws
- 031101 Allis Tissue Forceps 6" (152mm) 5x6 Teeth
with Angular Jaws
- 031102 Allis Tissue Forceps 7¼" (184mm) 5x6 Teeth
with Angular Jaws

FORCEPS



MedGyn Placenta Forceps

Hern-Barret Placenta Forceps

Slightly curved placenta forceps with oval jaws.

031151 Hern-Barret Placenta Forceps 13" (330mm)
16mm Serrated Jaws

031152 Hern-Barret Placenta Forceps 13" (330mm)
16mm Smooth Jaws



Kelly Placenta Forceps

Slightly curved placenta forceps with serrated jaws.

031140 Kelly Placenta Forceps 12.5" (318mm)
19mm Serrated Jaws



Hern-Kelly Placenta Forceps

Slightly curved placenta forceps with serrated jaws.

031145 Hern-Kelly 13" (330mm)
19mm Serrated Jaws



MedGyn Myoma Grasping Forceps

Myoma Grasping Forceps

Ideal for transcervical management of polyps and pedunculated myoma.

031210 Myoma Grasping Forceps 13" (330mm) 7mm Jaws

031211 Myoma Grasping Forceps 13" (330mm) 11mm
Jaws

031212 Myoma Grasping Forceps 13" (330mm) 17mm
Jaws



MedGyn Tissue Forceps

Ochsner Tissue Forceps

031169 Ochsner Tissue Forceps



Boston Forceps

Straight forceps with large oval jaws and coarse teeth.

031156 Boston Evacuation Forceps 13½" (343mm)
27mm Jaws - High Fulcrum

031157 Boston Evacuation Forceps 13½" (343mm)
27mm Jaws - Low Fulcrum



Meltzer Forceps

031143 Meltzer Evacuation Forceps 13" (330mm)
19mm Jaws



McClintock Forceps

Slightly curved with fine-tooth jaws.


031122 McClintock Ovum Forceps 9½"
(241mm) 14mm Jaws



Fink Forceps

Slightly curved with jaws and concave tip.

031176 Fink Ovum Forceps 10½" (267mm)
13mm Jaws



Lisa Forceps

031173 Lisa Forceps 10" (254mm)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

FORCEPS



Stubblefield Forceps

Slightly curved, very long forceps with oval jaws.

031144 Stubblefield Forceps 15" (381mm) 14mm Jaws



Agadi Evacuation Forceps

Smaller forceps available in both straight and curved versions with oval jaws.

031170 Agadi Evacuation Forceps 9¾" (248mm) Curved

031171 Agadi Evacuation Forceps 9½" (241mm) Straight



Blumenthal Forceps

Slightly curved forceps with large oval jaws and coarse teeth.

031141 Blumenthal Evacuation Forceps 13" (330mm)
20mm Jaws



Saenger Forceps

Forceps with unique inner-diameter teeth within oval jaws.

031163 Saenger Ovum Forceps 11¼" (286mm) Curved

031164 Saenger Ovum Forceps 10¾" (273mm) Straight
with Ratchet

031165 Saenger Ovum Forceps 10½" (267mm) Curved
with Ratchet


031172 Saenger Ovum Forceps 13" (330mm) Straight
with Ratchet



Jacobson Forceps

Slightly curved forceps, with coarse teeth, oval jaw and enlarged ring handles. With concave tip.

031142 Jacobson Evacuation Forceps 13" (330mm)
Large with 20mm Jaws



Hern-Patterson Forceps

Straight forceps with coarse teeth, oval jaws and ratchet.

031139 Hern Patterson Evacuation Forceps 13" (330mm) 17mm Jaws




Hern Van-Lyth Forceps

Curved forceps with fine teeth.

031166 Hern Van Lyth Evacuation Forceps 9½" (241mm) 9mm Jaws

031167 Hern Van Lyth Evacuation Forceps 10" (254mm) 11mm Jaws



Agadi-Bierer Forceps

031153 Agadi-Bierer Forceps 13" (330mm) - 14mm without Ratchet

031154 Agadi-Bierer Forceps 13" (330mm) - 14mm with Ratchet



Ramplly Ring Forceps

031007 Ramplly with Oval Jaws 10" (254mm) Straight



West-Gluck Forceps

031168 West-Gluck Forceps 11" (279mm)

FORCEPS



MedGyn Ovum Forceps

Sopher Ovum Forceps

Slightly curved forceps with oval, serrated jaws.

- 031123 Sopher Ovum Forceps 11" (279mm) 12mm Jaws
- 031123R Sopher Ovum Forceps 11" (279mm) 12mm Jaws with Ratchet
- 031126 Sopher Ovum Forceps 11" (279mm) 14mm Jaws
- 031149 Sopher Ovum Forceps 11" (279mm) 14mm Jaws with Ratchet
- 031124 Sopher Ovum Forceps 13" (330mm) 12mm Jaws
- 031125 Sopher Ovum Forceps 13" (330mm) 12mm Jaws with H Ratchet
- 031127 Sopher Ovum Forceps 13" (330mm) 14mm Jaws
- 031128 Sopher Ovum Forceps 13" (330mm) 14mm Jaws with Ratchet



Bierer Ovum Forceps

Slightly curved forceps with oval, coarse, serrated jaws.


- 031131 Bierer Ovum Forceps 11" (279mm) 16mm Jaws
- 031132 Bierer Ovum Forceps 11" (279mm) 16mm Jaws with Ratchet
- 031135 Bierer Ovum Forceps 11" (279mm) 19mm Jaws
- 031136 Bierer Ovum Forceps 11" (279mm) 19mm Jaws with Ratchet
- 031133 Bierer Ovum Forceps 13" (330mm) 16mm Jaws
- 031134 Bierer Ovum Forceps 13" (330mm) 16mm Jaws with Ratchet
- 031137 Bierer Ovum Forceps 13" (330mm) 19mm Jaws
- 031138 Bierer Ovum Forceps 13" (330mm) 19mm Jaws with Ratchet



Sopher Original Forceps

Slightly curved forceps with oval jaw. Concave tip design to prevent grasping of fundus.

- 031129 Sopher Original Forceps 10½" (267mm) Small Jaws
- 031130 Sopher Original Forceps 10½" (267mm) Medium Jaws



Hern Evacuation Forceps

Slightly curved forceps with coarse teeth, oval jaws.

- 031146 Hern Evacuation Forceps 13" (330mm) Small 12mm Jaws
- 031147 Hern Evacuation Forceps 13" (330mm) Medium 15mm Jaws
- 031148 Hern Evacuation Forceps 13" (330mm) Large 19mm Jaws



- 030316 Spiral Endocervical Curette Blue Titanium
- 030301 Kevorkian-Young Malleable Curette w/o Basket 12" (305mm)
- 030302 Kevorkian-Young Malleable Curette with Basket 12" (305mm)
- 030303 Kevorkian-Young Rigid Curette w/o Basket 12" (305mm)
- 030304 Kevorkian-Young Rigid Curette with Basket 12" (305mm)
- 030303-H Kevorkian-Young Rigid Curette w/o Basket (Square Tip, Gold Handle, Heavy) 12" (305mm)
- 030304-H Kevorkian-Young Rigid Curette with Basket (Square Tip, Gold Handle, Heavy) 12" (305mm)
- 030315 Kevorkian-Young Rigid Curette w/o Basket (Square Tip, Titanium) 12" (305mm)
- 030311 Kevorkian-Young Rigid Curette with Basket (Square Tip, Titanium) 12" (305mm)
- 030306 Kevorkian-Young Rigid Curette with Basket (Round Tip, Gold Handle, Heavy) 12" (305mm)
- 030310 Kevorkian-Young Rigid Curette with Basket (Round Tip, Titanium) 12" (305mm)

Townsend Curette

Round, tapered tip allows for entry through stenotic os.



- 030305 Townsend Curette 12" (305mm)

Metal Uterine Sounds

MedGyn reusable, graduated metal sounds made of highly malleable metal.



- 030704 Simpson Uterine Sound 3mm Diameter, Graduated in inches, 13" (330mm)
- 030705 Simpson Uterine Sound 3mm Diameter, Graduated in cm, 13" (330mm)

Simpson Uterine Sound



- 030700 Sims Uterine Sound 3mm Diameter, Graduated in inches, 13" (330mm)
- 030701 Sims Uterine Sound 3mm Diameter, Graduated in cm, 13" (330mm)
- 030702 Sims Uterine Sound 2mm Diameter, Graduated in cm, 13" (330mm)
- 030703 Sims Uterine Sound 3mm Plain Probe, Malleable, 13" (330mm)

METAL CURETTES



Sims Curette

Metal curette with sharp blade on malleable or rigid shanks 11" (279mm).

030640 Sims Malleable Curette 8 Piece Set Sizes #00 - #6

030640 12 Sims Malleable Curette 12 Piece Set Sizes #00-10

030641 Sims Malleable Curette #00

030642 Sims Malleable Curette #0

030643 Sims Malleable Curette #1

030644 Sims Malleable Curette #2

030645 Sims Malleable Curette #3

030646 Sims Malleable Curette #4

030647 Sims Malleable Curette #5

030648 Sims Malleable Curette #6

030649 Sims Malleable Curette #7

030637 Sims Malleable Curette #8

030638 Sims Malleable Curette #9

030639 Sims Malleable Curette #10

030650 Sims Rigid Curette 8 Piece Set Sizes #00 - #6

030651 Sims Rigid Curette #00

030652 Sims Rigid Curette #0

030653 Sims Rigid Curette #1

030654 Sims Rigid Curette #2

030655 Sims Rigid Curette #3

030656 Sims Rigid Curette #4

030657 Sims Rigid Curette #5

030658 Sims Rigid Curette #6

030659 Sims Rigid Curette #7

030660 Sims Rigid Curette #8

030661 Sims Rigid Curette #9

030662 Sims Rigid Curette #10

MedGyn Endocervical Curettes



Kevorkian-Youngue Curette

Metal curette available with or without basket.



Kevorkian-Youngue w/o Basket



Kevorkian-Youngue with Basket

MedGyn Endometrial Curettes



Novak Curette

Popular metal curette available in a 3 piece set or in individual sizes.

030600 Novak Curette 3 Piece Set 2mm - 4mm 9³/₄" (248mm)

030602 Novak Curette 2mm 9³/₄" (248mm)

030603 Novak Curette 3mm 9³/₄" (248mm)

030604 Novak Curette 4mm 9³/₄" (248mm)



Evans Curette

030630 Evans Curette 6 Piece Set Sizes #0 - #5 11" (279mm)

030631 Evans Curette #0-4mm 11" (279mm)

030632 Evans Curette #1-6mm 11" (279mm)

030633 Evans Curette #2-8mm 11" (279mm)

030634 Evans Curette #3-10mm 11" (279mm)

030635 Evans Curette #4-12mm 11" (279mm)

030636 Evans Curette #5-14mm 11" (279mm)



Randall Curette

030610 Randall Curette 3 Piece Set #2 - #4 9³/₄" (248mm)

030612 Randall Curette #1 9³/₄" (248mm)

030613 Randall Curette #2 9³/₄" (248mm)

030614 Randall Curette #3 9³/₄" (248mm)



Bumm Curette

Extra-large metal curettes on malleable shanks

030620 4 Piece Set - Sizes #1 - #4 12" (305mm)

030621 Bumm - Size #1 12" (305mm)

030622 Bumm - Size #2 12" (305mm)

030623 Bumm - Size #3 12" (305mm)

030624 Bumm - Size #4 12" (305mm)



Banjo/Hunter Curette

030619 Banjo/Hunter Curette Size 14 - 30mm



Heaney Curette

030680 Heaney Curette Small 9¹/₂" (241mm)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

DILATORS



McMahan Cervico Meter

Angled, used for measuring dilation of cervix.

031300 McMahan Cervico Meter



Mini Dilator Set

6 piece set (5 dilators & handle).

030779 6 Piece Mini Dilator Set:

1mm, 1.5mm, 2mm, 2.5mm, 3mm



030779BT Mini Dilator Set Blue Titanium

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

DILATORS



Hegar Dilator

Stainless steel single-end dilator 7½" (191mm).

030841-1 Hegar Dilator 1mm
030841-1.5 Hegar Dilator 1.5mm
030841-2 Hegar Dilator 2mm
030841-2.5 Hegar Dilator 2.5mm
030841-3 Hegar Dilator 3mm
030841-3.5 Hegar Dilator 3.5mm
030841-4 Hegar Dilator 4mm
030841-4.5 Hegar Dilator 4.5mm
030841-5 Hegar Dilator 5mm
030841-5.5 Hegar Dilator 5.5mm
030842-6 Hegar Dilator 6mm
030842-6.5 Hegar Dilator 6.5mm
030842-7 Hegar Dilator 7mm
030842-7.5 Hegar Dilator 7.5mm
030842-8 Hegar Dilator 8mm
030842-8.5 Hegar Dilator 8.5mm
030842-9 Hegar Dilator 9mm
030842-9.5 Hegar Dilator 9.5mm



030842-10 Hegar Dilator 10mm
030842-10.5 Hegar Dilator 10.5mm
030843-11 Hegar Dilator 11mm
030843-11.5 Hegar Dilator 11.5mm
030843-12 Hegar Dilator 12mm
030843-12.5 Hegar Dilator 12.5mm
030843-13 Hegar Dilator 13mm
030843-13.5 Hegar Dilator 13.5mm
030843-14 Hegar Dilator 14mm
030843-14.5 Hegar Dilator 14.5mm
030843-15 Hegar Dilator 15mm
030843-15.5 Hegar Dilator 15.5mm
030844-16 Hegar Dilator 16mm
030844-16.5 Hegar Dilator 16.5mm
030844-17 Hegar Dilator 17mm
030844-17.5 Hegar Dilator 17.5mm
030844-18 Hegar Dilator 18mm
030844-18.5 Hegar Dilator 18.5mm



Pratt Dilator Plastic

Light-weight, reusable 5 piece plastic dilator set.

030890 Pratt Dilator Plastic 5 Piece Set 5-6mm-
13-14mm



Os Cervical Finder Set

Made of soft, malleable Teflon™. 3 piece set
(Os Finder, Canal Finder & Fundus Finder)

030803 Os Cervical Finder 3 Piece Set

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

DILATORS



Hegar Dilator

Stainless steel double-end dilator 7½" (191mm).

030830 Hegar 8 Piece Dilator Set 3-4mm - 17-18mm

030831 Hegar Dilator 1-2mm

030832 Hegar Dilator 3-4mm

Hegar Dilator 5-6mm

030834 Hegar Dilator 7-8mm

030835 Hegar Dilator 9-10mm

030836 Hegar Dilator 11-12mm

030837 Hegar Dilator 13-14mm

030838 Hegar Dilator 15-16mm

030839 Hegar Dilator 17-18mm



Kleegman Dilator

Double-end dilator available in stainless steel and Teflon™. (Stainless steel pictured) (1.0mm/2.0mm)

030804 Kleegman Stainless Steel Dilator 11" (279mm)

030805 Kleegman Teflon™ Dilator 12" (305mm)

MedGyn Single-End Dilators



Hern Dilator

Teflon single-end dilator with handle 14¾" (375mm).

030807 Hern 21 Piece Dilator Set Sizes 44 Fr - 104 Fr

030808 Hern Dilator 44 Fr

030809 Hern Dilator 47 Fr

030810 Hern Dilator 50 Fr

030811 Hern Dilator 53 Fr

030812 Hern Dilator 56 Fr

030813 Hern Dilator 59 Fr

030814 Hern Dilator 62 Fr

030815 Hern Dilator 65 Fr

030816 Hern Dilator 68 Fr

030817 Hern Dilator 71 Fr

030818 Hern Dilator 74 Fr

030819 Hern Dilator 77 Fr

030820 Hern Dilator 80 Fr

030821 Hern Dilator 83 Fr

030822 Hern Dilator 86 Fr

030823 Hern Dilator 89 Fr

030824 Hern Dilator 92 Fr

030825 Hern Dilator 95 Fr

030826 Hern Dilator 98 Fr

030827 Hern Dilator 101 Fr

030828 Hern Dilator 104 Fr

030829-107 Hern Dilator

107 Fr

030829-110 Hern Dilator

110 Fr

030829-113 Hern Dilator

113 Fr



Hawkin Ambler

Stainless steel single-end, dilator 7¼" (184mm).

030780 16 Piece Dilator Set Sizes 3-6mm - 18-21mm

030781 Hawkin Ambler Dilator 3-6mm

030782 Hawkin Ambler Dilator 4-7mm

030783 Hawkin Ambler Dilator 5-8mm

030784 Hawkin Ambler Dilator 6-9mm

030785 Hawkin Ambler Dilator 7-10mm

030786 Hawkin Ambler Dilator 8-11mm

030787 Hawkin Ambler Dilator 9-12mm

030788 Hawkin Ambler Dilator 10-13mm

030789 Hawkin Ambler Dilator 11-14mm

030790 Hawkin Ambler Dilator 12-15mm

030791 Hawkin Ambler Dilator 13-16mm

030792 Hawkin Ambler Dilator 14-17mm

030793 Hawkin Ambler Dilator 15-18mm

030794 Hawkin Ambler Dilator 16-19mm

030795 Hawkin Ambler Dilator 17-20mm

030796 Hawkin Ambler Dilator 18-21mm

MedGyn Double-End Dilators



Pratt Dilator

Stainless steel double-end dilator 11½" (292mm).

030850 Pratt 8 Piece Dilator Set 13-15 Fr - 41-43 Fr

030850-ECO Pratt 8 Piece Dilator Set 13-15 Fr -41-43 Fr

030859 Pratt Dilator 9-11Fr

030851 Pratt Dilator 13-15 Fr

030852 Pratt Dilator 17-19 Fr

030853 Pratt Dilator 21-23 Fr

030854 Pratt Dilator 25-27 Fr

030855 Pratt Dilator 29-31 Fr

030856 Pratt Dilator 33-35 Fr

030857 Pratt Dilator 37-39 Fr

030858 Pratt Dilator 41-43 Fr



Hank Dilator

Designed with guide rings for ease-of-use 11" (279mm)

030801 Hank 6 Piece Dilator Set 11-12 Fr - 21-22 Fr

030801-7 Hank 7 Piece Dilator Set 9-10 Fr - 21-22 Fr

030801-7/8 Hank Dilator 7-8 Fr

030801-9/10 Hank Dilator 9-10 Fr

030801-11/12 Hank Dilator 11-12 Fr

030801-13/14 Hank Dilator 13-14 Fr

030801-15/16 Hank Dilator 15-16 Fr

030801-17/18 Hank Dilator 17-18 Fr

030801-19/20 Hank Dilator 19-20 Fr

030801-21/22 Hank Dilator 21-22 Fr



Pratt Extra Large Dilator

Pratt style in larger sizes 11½" (292mm).

030870 Pratt Extra Large 9 Piece Dilator Set 45-47 Fr - 77-79 Fr

030870-ECO Pratt Extra Large 9 Piece Dilator Set 45-47 Fr - 77-79 Fr

030871 Pratt Dilator Extra Large 45-47 Fr

030872 Pratt Dilator Extra Large 49-51 Fr

030873 Pratt Dilator Extra Large 53-55 Fr

030874 Pratt Dilator Extra Large 57-59 Fr

030875 Pratt Dilator Extra Large 61-63 Fr

030876 Pratt Dilator Extra Large 65-67 Fr

030877 Pratt Dilator Extra Large 69-71 Fr

030878 Pratt Dilator Extra Large 73-75 Fr

030879 Pratt Dilator Extra Large 77-79 Fr

030880 Pratt Extra Large, Light Weight 9 Piece Dilator Set 45-47 Fr - 77-79 Fr

030881 Pratt Extra Large, Light Weight 45-47 Fr

030882 Pratt Extra Large, Light Weight 49-51 Fr

030883 Pratt Extra Large, Light Weight 53-55 Fr

030884 Pratt Extra Large, Light Weight 57-59 Fr

030885 Pratt Extra Large, Light Weight 61-63 Fr

030886 Pratt Extra Large, Light Weight 65-67 Fr

030887 Pratt Extra Large, Light Weight 69-71 Fr

030888 Pratt Extra Large, Light Weight 73-75 Fr

030889 Pratt Extra Large, Light Weight 77-79 Fr

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

TENACULUMS



Atraumatic Tenaculum

- 031001 Atraumatic Straight 1/4" (6mm) Jaws, 9" (229mm)
- 031002 Atraumatic Straight 1/8" (3mm) Jaws, 9" (229mm)
- 031003 Atraumatic Angled 1/4" (6mm) Jaws, 9 1/2" (241mm)



Shah Tenaculum

- 031016 Shah Tenaculum 9 1/2" (241mm)



Iris Hook

- 030406 Iris Hook 10" (254mm)



Emmett Tenaculum Hooks

- 031401 Emmett Curved Hook 9 1/2" (241mm)
- 031402 Emmett Half-Curved Hook 9 3/4" (248mm)
- 031403 Emmett Hook L-Shaped 9" (229mm)
- 031404 Emmett Hook Acute Angle 9" (229mm)
- 031405 Emmett Hook Double Angled 9" (229mm)





Schroeder Double-Tooth Tenaculum

- 031011 Schroeder Double-Tooth - 10" Straight (254mm)
- 031012 Schroeder Double-Tooth - 10" Curved (254mm)
- 031013 Schroeder Double-Tooth - 10" Side Curve (254mm)



Jacobs Double-Tooth Tenaculum

- 031197 Jacobs Double-Tooth 8½" Straight - 2x2 Teeth (215mm)
- 031198 Jacobs Double-Tooth 8½" Side Curve - 2x2 Teeth (215mm)



Bierer Tenaculum

- 031010 Bierer Angled 5x6 Teeth, 9" (229mm)
- 031033 Bierer Angled 5x6 Teeth, Economy 9" (229mm)



Teale Tenaculum

- 031005 Teale Curved, 3x4 Teeth, 5mm Wide, 10" (254mm)



Kelly Tenaculum

- 031017 Kelly Tenaculum 4X5 Teeth, 9" (229mm)

TENACULUMS



MedGyn Tenaculums

Schroeder Braun Single-Tooth Tenaculum

031004 Schroeder Braun Single-Tooth - 9¾" (248mm)

031019 Schroeder Braun Tenaculum Light-Touch Straight - 9¾" (248mm)

031034 Schroeder Braun Tenaculum Economy - 9¾" (248 mm)

031018 Braun Tenaculum Light-Touch Round - 7¼" (184mm)



Skene Single-Tooth Tenaculum

031006 Skene Single-Tooth-Side-Curve - 10" (254mm)



Duplay Single-Tooth Tenaculum

031008 Duplay Single-Tooth Double-Curve - 11" (279mm)



Barrett Tenaculum

031009 Barrett Straight - 7¼" (184mm)



Jacobson Double-Tooth Tenaculum

031199 Jacobson Double-Tooth 10" Straight - 2x2 Teeth with Larger Rings (254mm)

031200 Jacobson Double-Tooth 10" Side Curve - 2x2 Teeth with Larger Rings (254mm)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

BIOPSY PUNCHES



Van Doren Biopsy Forceps

- 030155 Van Doren Biopsy Forceps with Ratchet 10" (254mm)
- 030166 Van Doren Biopsy Forceps without Ratchet 10" (254mm)



Schubert Biopsy Forceps

- 030157 Schubert Biopsy Forceps 10" (254mm) Curved



Wittner Biopsy Punch

- 030117 Wittner Biopsy Punch 8½" (216mm)



Rotating Biopsy Punch Options

- 030778 Vulvar Keyes Biopsy Punch Set (Handle + 6 Tips)
- 030104 Rotating Tischler Biopsy Punch
- 030145 Rotating Tischler Biopsy Punch - Angled Down
- 030146 Rotating Tischler Biopsy Punch - Angled Up
- 030105 Rotating Baby Tischler Biopsy Punch
- 030105-DOWN Rotating Baby Tischler Biopsy Punch - Angled Down
- 030105-UP Rotating Baby Tischler Biopsy Punch - Angled Up
- 030106 Rotating Kevorkian Biopsy Punch
- 030107 Rotating Burke Biopsy Punch
- 030111-BB Rotating Shaft - Baby Tischler Tip Only
- 030111-BRK Rotating Shaft - Burke Tip Only
- 030111-KEV Rotating Shaft - Kevorkian Tip Only
- 030111-TIS Rotating Shaft - Tischler Tip Only
- 030147 Rotating Shaft - Tischler Tip - Up
- 030148 Rotating Shaft - Tischler Tip - Down
- 030108 Rotating Biopsy Punch Handle Only
- 030100 Rotating Shaft (with any tip)



BIOPSY PUNCHES



Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma.

Bite Size: 4.2mm x 2.3mm

- 030102-10 Baby Tischler Biopsy Punch 10 (254mm)
- 029811 Baby Tischler Biopsy Punch 10 (254mm) Angled Up
- 029812 Baby Tischler Biopsy Punch 10 (254mm) Angled Down
- 029813 Baby Tischler Biopsy Punch 10 (254mm) Angled Left
- 029814 Baby Tischler Biopsy Punch 10 (254mm) Angled Right
- 030102 Baby Tischler Biopsy Punch 8½ (216mm)
- 030140 Baby Tischler Biopsy Punch 8½ (216mm) Angled Up
- 030141 Baby Tischler Biopsy Punch 8½ (216mm) Angled Down
- 030142 Baby Tischler Biopsy Punch 8½ (216mm) Angled Left
- 030143 Baby Tischler Biopsy Punch 8½ (216mm) Angled Right



Kevorkian Biopsy Punch

(Kevorkian-Youngs Biopsy Forceps displayed)

- 030103-10 Kevorkian Biopsy Punch 10 (254mm)
- 030103 Kevorkian Biopsy Punch 8½ (216mm)
- 030169 Kevorkian-Youngs Biopsy Forceps without Spring 8" (203mm)



Tischler-Kevorkian Cervical Biopsy Punch

- 030116 Tischler-Kevorkian Biopsy Punch 10 (254mm)
- 030115 Tischler-Kevorkian Biopsy Punch 8½ (216mm)



Burke Biopsy Punch

Open-handle design for optimal control with small cutting jaws and sharp teeth.

- 030109 Burke Biopsy Punch 8½



Eppendorfer Biopsy Forceps

- 030110 Eppendorfer Biopsy Forceps 8" (203mm)

BIOPSY PUNCHES



Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

- 030173 Baby Tischler Biopsy Punch Titanium 10 (254mm)
Heavy-Grip Handle
- 030121 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- 030134 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- Angled Up
- 030135 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- Angled Down
- 030161 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- Angled Left
- 030160 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- Angled Right
- 030119 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- 030132 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- Angled Up
- 030133 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
-Angled Down
- 030159 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- Left
- 030158 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- Right



MedGyn Stainless Steel

Tischler Biopsy Punch

Classic Tischler bite designed for drawing in maximum tissue near Os with minimal trauma.

Bite Size 7mm x 3mm

- 030101-10 Tischler Biopsy Punch 10 (254mm)
- 029801 Tischler Biopsy Punch 10 (254mm) Angled Up
- 029802 Tischler Biopsy Punch 10 (254mm) Angled Down
- 029803 Tischler Biopsy Punch 10 (254mm) Angled Left
- 029804 Tischler Biopsy Punch 10 (254mm) Angled Right
- 030101 Tischler Biopsy Punch 8½ (216mm)
- 030113 Tischler Biopsy Punch 8½ (216mm) Angled Up
- 030114 Tischler Biopsy Punch 8½ (216mm) Angled Down
- 030136 Tischler Biopsy Punch 8½ (216mm) Angled Left
- 030137 Tischler Biopsy Punch 8½ (216mm) Angled Right

BIOPSY PUNCHES



MedGyn Titanium

Tischler Biopsy Punch

Classic Tischler bite, with durable coating.

Bite Size: 7mm x 3mm

030172 Tischler Biopsy Punch Titanium 10 (254mm) Heavy-Grip Handle

030120 Tischler Biopsy Punch Titanium 10 (254mm)

030124 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Up

030125 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Down

030162 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Left

030163 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Right

030118 Tischler Biopsy Punch Titanium 8½ (216mm)

030122 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Up

030123 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Down

030138 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Right

030139 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Left



Baby Tischler Biopsy Punch (Forceps Style)

Classic Baby Tischler bite with durable coating. Bite Size:

4.2mm x 2.3mm

030170 Baby Tischler Biopsy Punch Titanium Gold Handle (Forceps Style) 10" (254mm)



Oval Tischler Biopsy Punch

030171 Oval Tischler Biopsy Punch Titanium Heavy-Grip Handle 10 (254mm)



Kevorkian Biopsy Punch

Popular design with strong, narrow jaw with durable coating. Bite Size: 8mm x 3mm.

030174 Kevorkian Biopsy Punch Titanium Heavy-Grip Handle 10" (254mm)

030127 Kevorkian Biopsy Punch Titanium 10 (254mm)

030126 Kevorkian Biopsy Punch Titanium 8½ (216mm)



MedGyn Blue Titanium

Classic Tischler bite with durable coating.
Bite Size: 7mm x 3mm

030176 Tischler Biopsy Punch Blue Titanium
Heavy-Grip Handle 10 (254mm)

030181 Tischler Biopsy Punch Blue Titanium 10
(254mm)



Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

030177 Baby Tischler Biopsy Punch Blue Titanium
Heavy-Grip Handle 10 (254mm)

030129 Baby Tischler Biopsy Punch Blue Titanium
8½ (216mm)

030180 Baby Tischler Biopsy Punch Blue Titanium
10 (254mm)



Baby Tischler Biopsy Punch (Forceps Style)

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

030179 Baby Tischler Biopsy Punch Blue Titanium
Gold Handle (Forceps Style) 10" (254mm)



Kevorkian Biopsy Punch

Popular design with strong, narrow jaw with durable coating and a heavy-grip handle. Bite Size: 8mm x 3mm

030178 Kevorkian Biopsy Punch Blue Titanium
Heavy-Grip Handle 10 (254mm)



Oval Tischler Biopsy Punch

Oval Tischler bite with durable coating and a heavy-grip handle.

030175 Oval Tischler Biopsy Punch Blue Titanium
Heavy-Grip Handle 10 (254mm)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SPECULUMS

Weighted Speculum



- 030921 Auward 2.0lbs, 75° angle,
Blade 1½" x 3¼" (38mm x 83mm)
- 030922 Auward 2.5lbs, 75° angle,
Blade 1½" x 3¼" (38mm x 83mm)
- 030923 Auward 3.0lbs, 75° angle,
Blade 1½" x 3¼" (38mm x 83mm)
- 030924 Auward 2.5lbs, 45° angle,
Blade 1½" x 4" (35mm x 102mm)
- 030926 Auward 2.5lbs, 90° angle,
Blade 1" x 5½" (25mm x 140mm)

Collins Speculum

The unique hands-free design offers the advantages of an open-sided speculum while permitting multiple positioning for optimal exposure with blades that open laterally instead of vertically.



- 030946 1¼" x 3¾" (32mm x 95mm)
- 030947 1½" x 4" (35mm x 102mm)
- 030948 1½" x 4½" (41mm x 114mm)


Lateral Vaginal Wall Retractor

Used for the retraction of vaginal wall and made to fit inside any standard speculum to expand past blades of speculum. Allows for greater exposure of cervix for safer, more accurate procedures.




- 030970 Lateral Vaginal Wall Retractor Blue Titanium
- 030966 Open ½" (13mm), 10" (254mm),
Blade Size 1" x 2½" (25mm x 64mm)
- 030967 Open 1" (25mm), 10" (254mm),
Blade Size 1" x 2½" (25mm x 64mm)
- 030968 Wide Open, 10" (254mm),
Blade Size 1" x 2½" (25mm x 64mm)
- 030969 Closed, 10" (254mm),
Blade Size 1" x 2½" (25mm x 64mm)

Collins Speculum



Endocervical Speculum



- 030205 Endocervical Speculum Blue Titanium
- 030915 Kogan with Screw 2-4mm
- 030916 Kogan with Screw 4-8mm
- 030201 Kogan with Gauge & Ratchet 2-4mm
- 030202 Kogan with Gauge & Ratchet 4-8mm
- 030203 Kogan with Gauge & Ratchet 8-12mm
- 030204 Townsend Endocervical Speculum 8-12mm

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SPECULUMS

Klopper Speculum

Wider opening (4cm) and short blade



030931 1 3/8" x 3 1/2" (35mm x 89mm)

Grundman Speculum

Wide open design with screw adjustment/lever



030949 Small 1 1/8" x 3 1/8" (34mm x 85mm)

030950 Medium 1 3/8" x 3 3/4" (35mm x 95mm)

030951 Large 1 1/2" x 4 3/4" (38mm x 121mm)

Cusco Speculum

Unique handle-design speculum with center screw or side screw options.



030960 Small w/Center Screw 1 1/8" x 3 1/2" (34mm x 89mm)

030961 Medium w/Center Screw 1 3/8" x 4" (35mm x 102mm)


030962 Large w/Center Screw 1 3/8" x 4 1/2" (35mm x 114mm)

030963 Small w/Side Screw 1 1/8" x 3 1/2" (34mm x 89mm)

030964 Medium w/Side Screw 1 3/8" x 4" (35mm x 102mm)

030965 Large w/Side Screw 1 3/8" x 4 1/2" (35mm x 114mm)

Klopper, Grundman, Cusco, For-Way Speculum



Four-Way Speculum

The design of this speculum combines the characteristics of both a standard speculum and a vaginal wall retractor.



038095 4-Way Lateral Expander Speculum - Medium
7/8" x 4 1/4" (22mm x 108mm)

038096 4-Way Lateral Expander Speculum - Large
7/8" x 5 1/4" (22mm x 133mm)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SPECULUMS

Pederson Vu-Max Speculum

Largest view area (7cm vs. standard 3cm) with durable coating.



030983 Medium 4½" (114mm)

030984 Large 5" (127mm)

Pederson Speculum

Narrow blade design



030941 XNarrow ⅝" x 4" (17mm x 102mm)

030913 XSmall ½" x 2¾" (13mm x 70mm)

030906 Small ½" x 3" (13mm x 76mm)


030907 Medium ⅞" x 4" (22mm x 102mm)

030908 Large 1" x 4¾" (25mm x 121mm)

030914 XLarge 1" x 6½" (25mm x 165mm)

030917 Extra Long Snowman Clinton 1" x 6" (25mm x 152mm)

Pederson Speculum



Pederson Vu-More Speculum

Larger view area (4cm vs. standard 3cm) with durable coating.



030981 Medium 4½" (114mm)

030982 Large 5" (127mm)

Pederson - Open Side Speculum

Narrow blade design with open side



030972 45° Medium - Right ⅞" x 4" (22mm x 102mm)

030973 45° Large - Right 1" x 4¾" (25mm x 121mm)

030942 XNarrow - Right ⅝" x 4" (17mm x 102mm)

030943 XNarrow - Left ⅝" x 4" (17mm x 102mm)

030939 Medium - Right ⅞" x 4" (22mm x 102mm)

030939-LT Medium - Left ⅞" x 4" (22mm x 102mm)

030908R Large - Right 1" x 4¾" (25mm x 121mm)

030908-LT Large - Left 1" x 4¾" (25mm x 121mm)

030914-L-XL XLarge - Left 1" x 6½" (25mm x 165mm)

030914-R-XL XLarge - Right 1" x 6½" (25mm x 165mm)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SPECULUMS

Moore Graves Speculum
Short blade design



030920 1 1/4" x 3 1/2" (32mm x 89mm)

Graves Speculum
Wider blade design




030905 XSmall 3/4" x 2 3/4" (19mm x 70mm)
030901 Small 3/4" x 3" (19mm x 76mm)
030902 Medium 1 3/8" x 4" (35mm x 102mm)
030903 Large 1 3/8" x 4 3/4" (35mm x 121mm)
030904 XLarge 1 1/2" x 6 1/2" (38mm x 165mm)
030900 Extra Long Snowman Clinton 1 1/2" x 6" (38mm x 152mm)

Graves Wide Angled Speculum

Wider Graves blade design with 45° open angle.



030909 Medium 1 3/8" x 4" (35mm x 102mm)
030910 Right Side Open Medium 1 3/8" x 4" (35mm x 102mm)
030971 Right Side Open Large 1 3/8" x 4 3/4" (35mm x 121mm)



PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

SPECULUMS

Graves Vu-More Speculum

Wider blade design with larger view area (4cm vs. standard 3cm) with durable coating.



030985 Medium 4½" (114mm)
030986 Large 5" (127mm)


Graves Vu-More Speculum

Wider blade design with larger view area (4cm vs. standard 3cm).



030952 Small 3¾" (95mm)
030953 Medium 4½" (114mm)
030954 Large 5" (127mm)
030955 XLong 6" (152mm)
030956 XLong 6½" (161mm)
030957 XLong 6¾" (171mm)
030958 XLong 7" (178mm)

Graves Speculum



Graves Vu-Max Speculum

Wider blade design with largest view area (7cm vs. standard 3cm) with durable coating.



030987 Medium 4½" (114mm)
030988 Large 5" (127mm)

Graves Vu-Max Speculum

Wider blade design with largest view area (7cm vs. standard 3cm).



030975 Medium 4½" (114mm)
030976 Large 5" (127mm)
030977 XLong 6" (152mm)
030978 XLong 6½" (161mm)
030979 XLong 6¾" (171mm)
030980 XLong 7" (178mm)

Weisman Graves Speculum

Open-side speculum



030911-SM Small - Right ¾" x 3"
(10mm x 76mm)
030911 Medium - Right 1 ¾" x 4"
(35mm x 102mm)
030911-LG Large - Right 1 ¾" x 4 ¾"
(35mm x 121mm)
030911-XL XLarge - Right 1 ½" x 6 ½"
(38mm x 165mm)
030912 Medium - Left 1 ¾" x 4"
(35mm x 102mm)
030912-LG Large - Left 1 ¾" x 4 ¾"
(35mm x 121mm)
030912-XL XLarge - Left 1 ½" x 6 ½"
(38mm x 165mm)



MedGyn Oval with Support

- 050111 Oval w/ Support #1 (51mm, 2in)
- 050112 Oval w/ Support #2 (57mm, 2.25in)
- 050113 Oval w/ Support #3 (64mm, 2.5in)
- 050114 Oval w/ Support #4 (69mm, 2.7in)
- 050115 Oval w/ Support #5 (76mm, 3in)
- 050116 Oval w/ Support #6 (83mm, 3.25in)
- 050117 Oval w/ Support #7 (89mm, 3.5 in)
- 050118 Oval w/ Support #8 (95mm, 3.75in)
- 050119 Oval w/ Support #9 (102mm, 4in)

MedGyn Oval w/o Support

- 050101 Oval w/o Support #1 (51mm, 2in)
- 050102 Oval w/o Support #2 (57mm, 2.25in)
- 050103 Oval w/o Support #3 (64mm, 2.5in)
- 050104 Oval w/o Support #4 (69mm, 2.7in)
- 050105 Oval w/o Support #5 (76mm, 3in)
- 050106 Oval w/o Support #6 (83mm, 3.25in)
- 050107 Oval w/o Support #7 (89mm, 3.5in)
- 050108 Oval w/o Support #8 (95mm, 3.75 in)
- 050109 Oval w/o Support #9 (102mm, 4in)



MedGyn Dish with Support

- 050137 Dish w/ Support #0 (50mm, 2in)
- 050138 Dish w/ Support #1 (55mm, 2.2in)
- 050139 Dish w/ Support #2 (60mm, 2.4in)
- 050140 Dish w/ Support #3 (63mm, 2.5in)
- 050141 Dish w/ Support #4 (68mm, 2.7in)
- 050142 Dish w/ Support #5 (73mm, 2.9in)
- 050143 Dish w/ Support #6 (80mm, 3.1in)
- 050144 Dish w/ Support #7 (85mm, 3.3in)

MedGyn Dish w/o Support

- 050129 Dish w/o Support #0 (50mm, 2in)
- 050130 Dish w/o Support #1 (55mm, 2.2in)
- 050131 Dish w/o Support #2 (60mm, 2.4in)
- 050132 Dish w/o Support #3 (63mm, 2.5in)
- 050133 Dish w/o Support #4 (68mm, 2.7in)
- 050134 Dish w/o Support #5 (73mm, 2.9in)
- 050135 Dish w/o Support #6 (80mm, 3.1in)
- 050136 Dish w/o Support #7 (85mm, 3.3in)



MedGyn Pessary Remover

Designed to aid in the removal of many pessaries, the firm grip and flexibility allows patients to safely and easily remove their own pessaries. The high-quality plastic is not cold, unlike metal removers, greatly increasing patient comfort.

031315 Pessary Remover

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PESSARIES



MedGyn Cube with Drain

050066 Cube w/ Drain #0 (25mm, 0.98in)
050067 Cube w/ Drain #1 (29mm, 1.1in)
050068 Cube w/ Drain #2 (32mm, 1.2in)
050069 Cube w/ Drain #3 (36mm, 1.4in)
050070 Cube w/ Drain #4 (40mm, 1.6in)
050071 Cube w/ Drain #5 (45mm, 1.8in)
050072 Cube w/ Drain #6 (49mm, 1.9in)
050073 Cube w/ Drain #7 (55mm, 2.2in)
050074 Cube w/ Drain #8 (63mm, 2.5in)
050075 Cube w/ Drain #9 (69mm, 2.7in)
050076 Cube w/ Drain #10 (74mm, 2.9in)



MedGyn Cube w/o Drain

050055 Cube w/o Drain #0 (25mm, 0.98in)
050056 Cube w/o Drain #1 (29mm, 1.1in)
050057 Cube w/o Drain #2 (32mm, 1.2in)
050058 Cube w/o Drain #3 (36mm, 1.4in)
050059 Cube w/o Drain #4 (40mm, 1.6in)
050060 Cube w/o Drain #5 (45mm, 1.8in)
050061 Cube w/o Drain #6 (49mm, 1.9in)
050062 Cube w/o Drain #7 (55mm, 2.2in)
050063 Cube w/o Drain #8 (63mm, 2.5in)
050064 Cube w/o Drain #9 (69mm, 2.7in)
050065 Cube w/o Drain #10 (74mm, 2.9in)



MedGyn Hodge w/o Support

050035 Hodge w/o Support #0 (65mm, 2.6in)
050036 Hodge w/o Support #1 (70mm, 2.8in)
050037 Hodge w/o Support #2 (75mm, 3in)
050038 Hodge w/o Support #3 (80mm, 3.2in)
050039 Hodge w/o Support #4 (85mm, 3.4in)
050040 Hodge w/o Support #5 (90mm, 3.6in)
050041 Hodge w/o Support #6 (95mm, 3.8in)
050042 Hodge w/o Support #7 (100mm, 4in)
050043 Hodge w/o Support #8 (105mm, 4.2in)
050044 Hodge w/o Support #9 (110mm, 4.4in)



MedGyn Gehrung with Support

050145 Gehrung w/ Support #0 (50mm, 1.9in)
050146 Gehrung w/ Support #1 (55mm, 2.2in)
050147 Gehrung w/ Support #2 (60mm, 2.4in)
050148 Gehrung w/ Support #3 (65mm, 2.6in)
050149 Gehrung w/ Support #4 (70mm, 2.8in)
050150 Gehrung w/ Support #5 (75mm, 3in)
050151 Gehrung w/ Support #6 (80mm, 3.2in)
050152 Gehrung w/ Support #7 (85mm, 3.4in)
050153 Gehrung w/ Support #8 (90mm, 3.6in)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PESSARIES



MedGyn Gellhorn Short Stem

Stem Height 28 mm, 1.1 in

- 050155S Gellhorn Short Stem #1 (44mm, 1.7in)
- 050156S Gellhorn Short Stem #2 (51mm, 2in)
- 050157S Gellhorn Short Stem #3 (57mm, 2.2in)
- 050158S Gellhorn Short Stem #4 (63mm, 2.5in)
- 050159S Gellhorn Short Stem #5 (70mm, 2.8in)
- 050160S Gellhorn Short Stem #6 (76mm, 3in)
- 050161S Gellhorn Short Stem #7 (82mm, 3.2in)
- 050162S Gellhorn Short Stem #8 (89mm, 3.5in)



MedGyn Shaatz with Drain

- 050120 Shaatz w/ Drain #0 (38mm, 1.5in)
- 050121 Shaatz w/ Drain #1 (44mm, 1.75in)
- 050122 Shaatz w/ Drain #2 (51mm, 2in)
- 050123 Shaatz w/ Drain #3 (57mm, 2.2in)
- 050124 Shaatz w/ Drain #4 (64mm, 2.5in)
- 050125 Shaatz w/ Drain #5 (70mm, 2.8in)
- 050126 Shaatz w/ Drain #6 (76mm, 3in)
- 050127 Shaatz w/ Drain #7 (82mm, 3.2in)
- 050128 Shaatz w/ Drain #8 (88mm, 3.5in)



MedGyn Cup with Support

- 050171 Cup w/ Support #0 (48mm, 1.8in)
- 050172 Cup w/ Support #1 (52mm, 2in)
- 050173 Cup w/ Support #2 (57mm, 2.25in)
- 050174 Cup w/ Support #3 (63mm, 2.5in)
- 050175 Cup w/ Support #4 (68mm, 2.7in)
- 050176 Cup w/ Support #5 (73mm, 2.9in)
- 050177 Cup w/ Support #6 (77mm, 3in)
- 050178 Cup w/ Support #7 (82mm, 3.2in)

MedGyn Cup w/o Support

- 050163 Cup w/o Support #0 (48mm, 1.8in)
- 050164 Cup w/o Support #1 (52mm, 2in)
- 050165 Cup w/o Support #2 (57mm, 2.25in)
- 050166 Cup w/o Support #3 (63mm, 2.5in)
- 050167 Cup w/o Support #4 (68mm, 2.7in)
- 050168 Cup w/o Support #5 (73mm, 2.9in)
- 050169 Cup w/o Support #6 (77mm, 3in)
- 050170 Cup w/o Support #7 (82mm, 3.2in)




MedGyn Hodge with Support

- 050045 Hodge w/ Support #0 (65mm, 2.6in)
- 050046 Hodge w/ Support #1 (70mm, 2.8in)
- 050047 Hodge w/ Support #2 (75mm, 3in)
- 050048 Hodge w/ Support #3 (80mm, 3.2in)
- 050049 Hodge w/ Support #4 (85mm, 3.4in)
- 050050 Hodge w/ Support #5 (90mm, 3.6in)
- 050051 Hodge w/ Support #6 (95mm, 3.8in)
- 050052 Hodge w/ Support #7 (100mm, 4in)
- 050053 Hodge w/ Support #8 (105mm, 4.2in)
- 050054 Hodge w/ Support #9 (110mm, 4.4in)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PESSARIES



MedGyn Ring with Knob w/Support

- 050026K Ring w/ Knob - w/ Support #1 (51mm, 2in)
- 050027K Ring w/ Knob - w/ Support #2 (57mm, 2.2in)
- 050028K Ring w/ Knob - w/ Support #3 (63mm, 2.5in)
- 050029K Ring w/ Knob - w/ Support #4 (70mm, 2.8in)
- 050030K Ring w/ Knob - w/ Support #5 (76mm, 3in)
- 050031K Ring w/ Knob - w/ Support #6 (83mm, 3.3in)
- 050032K Ring w/ Knob - w/ Support #7 (89mm, 3.5in)

MedGyn Ring w/ Knob w/o Support

- 050017K Ring w/ Knob - w/o Support #1 (51mm, 2in)
- 050018K Ring w/ Knob - w/o Support #2 (57mm, 2.2in)
- 050019K Ring w/ Knob - w/o Support #3 (63mm, 2.5in)
- 050020K Ring w/ Knob - w/o Support #4 (70mm, 2.8in)
- 050021K Ring w/ Knob - w/o Support #5 (76mm, 3in)
- 050022K Ring w/ Knob - w/o Support #6 (83mm, 3.3in)
- 050023K Ring w/ Knob - w/o Support #7 (89mm, 3.5in)



MedGyn Gellhorn with Drain

- 050155 Gellhorn w/ Drain #1 (44mm, 1.7in)
- 050156 Gellhorn w/ Drain #2 (51mm, 2in)
- 050157 Gellhorn w/ Drain #3 (57mm, 2.2in)
- 050158 Gellhorn w/ Drain #4 (63mm, 2.5in)
- 050159 Gellhorn w/ Drain #5 (70mm, 2.8in)
- 050160 Gellhorn w/ Drain #6 (76mm, 3in)
- 050161 Gellhorn w/ Drain #7 (82mm, 3.2in)
- 050162 Gellhorn w/ Drain #8 (89mm, 3.5in)




MedGyn Marland with Support

- 050084 Marland w/ Support #2 (57mm, 2.25in)
- 050085 Marland w/ Support #3 (63mm, 2.5in)
- 050086 Marland w/ Support #4 (70mm, 2.75in)
- 050087 Marland w/ Support #5 (76mm, 3in)
- 050088 Marland w/ Support #6 (83mm, 3.25in)
- 050089 Marland w/ Support #7 (89mm, 3.5in)
- 050090 Marland w/ Support #8 (95mm, 3.75in)

Medgyn Marland w/o Support

- 050077 Marland w/o Support #2 (57mm, 2.25in)
- 050078 Marland w/o Support #3 (63mm, 2.5in)
- 050079 Marland w/o Support #4 (70mm, 2.75in)
- 050080 Marland w/o Support #5 (76mm, 3in)
- 050081 Marland w/o Support #6 (83mm, 3.25in)
- 050082 Marland w/o Support #7 (89mm, 3.5in)
- 050083 Marland w/o Support #8 (95mm, 3.75in)



MedGyn Donut

- 050010 Donut #0 (51mm, 2in)
- 050011 Donut #1 (57mm, 2.2in)
- 050012 Donut #2 (64mm, 2.5in)
- 050013 Donut #3 (69mm, 2.7in)
- 050014 Donut #4 (76mm, 3in)
- 050015 Donut #5 (82mm, 3.2in)
- 050016 Donut #6 (88mm, 3.5in)
- 050199 Donut #7 (95mm, 3.7in)

MedGyn Silicone Pessaries

Pessaries are an effective tool for the non-surgical treatment of a variety of gynecologic problems. Commonly used for management of pelvic support defects, pessaries are also used to treat stress urinary incontinence. MedGyn offers an extensive line of quality medical-grade silicone pessaries and easy-to-use pessary fitting sets.



051001 MedGyn Pessary Case
(28 different pessaries including Ring Fitting Set)



051000 Ring w/o Support
Pessary Fitting Set of 6



051005 Ring w/ Support
Pessary Fitting Set of 6



051006 Shaatz Pessary
Fitting Set of 6



MedGyn Ring w/o Support

- 050017 Ring w/o Support #1 (51mm, 2in)
- 050018 Ring w/o Support #2 (57mm, 2.2in)
- 050019 Ring w/o Support #3 (62mm, 2.4in)
- 050020 Ring w/o Support #4 (70mm, 2.75in)
- 050021 Ring w/o Support #5 (74mm, 2.9in)
- 050022 Ring w/o Support #6 (82mm, 3.2in)
- 050023 Ring w/o Support #7 (89mm, 3.5in)
- 050024 Ring w/o Support #8 (96mm, 3.8in)
- 050025 Ring w/o Support #9 (102mm, 4in)

MedGyn Ring with Support

- 050026 Ring w/ Support #1 (51mm, 2in)
- 050027 Ring w/ Support #2 (57mm, 2.2in)
- 050028 Ring w/ Support #3 (62mm, 2.4in)
- 050029 Ring w/ Support #4 (70mm, 2.75in)
- 050030 Ring w/ Support #5 (74mm, 2.9in)
- 050031 Ring w/ Support #6 (82mm, 3.2in)
- 050032 Ring w/ Support #7 (89mm, 3.5in)
- 050033 Ring w/ Support #8 (96mm, 3.8in)
- 050034 Ring w/ Support #9 (102mm, 4in)



MedGyn Modified Cup w/ Drain

- 050210M Modified Cup #1 w/ Drain (65mm, 2.6in)
- 050211M Modified Cup #2 w/ Drain (70mm, 2.8in)
- 050212M Modified Cup #3 w/ Drain (65mm, 2.6in)
- 050213M Modified Cup #4 w/ Drain (70mm, 2.8in)
- 050214M Modified Cup #5 w/ Drain (65mm, 2.6in)
- 050215M Modified Cup #6 w/ Drain (70mm, 2.8in)

MedGyn Modified Cup w/o Drain

- 050191M Modified Cup #1 w/o Drain (65mm, 2.6in)
- 050192M Modified Cup #2 w/o Drain (70mm, 2.8in)
- 050193M Modified Cup #3 w/o Drain (65mm, 2.6in)
- 050194M Modified Cup #4 w/o Drain (70mm, 2.8in)
- 050195M Modified Cup #5 w/o Drain (65mm, 2.6in)
- 050196M Modified Cup #6 w/o Drain (70mm, 2.8in)

DISPOSABLES KITS

MedGyn Colposcopy Set

Everything in one sterile package for an efficient and safe procedure.

030090 Disposable Colposcopy Set

030091 Disposable Colposcopy Set with Rotating Biopsy Punch.



MedGyn IUD Kit

This sterile, single use kit contains everything needed to insert or remove an IUD (intrauterine device).

We have included the following products:

022365 MedGyn IUD Kit consists of Cherron Forceps, Pozzy Forceps, Long Littauer Scissors, Uterine Sound (25/Box)



022365-M IUD Kit consisting of Metal Cherron Forceps, Pozzy Forceps, Long Littauer Scissors, Uterine Sound (25/box)





MedGyn Monsel's Solution

Ferric Subsulfate Solution 20% is an effective hemostatic solution for use after cervical biopsies.

022610 Monsel's Solution - 8mL bottles
(12/Box)

022612 Monsel's Solution - 500 mL bottle



MedGyn's Disposable Derma Curette

022797-3 Disposable Derma- Curette 3.0mm- STERILE

022797-4 Disposable Derma- Curette 4.0mm- STERILE

022797-7 Disposable Derma- Curette 7.0mm- STERILE



MedGyn's Disposable Derma (Keyes Biopsy) Punch

MedGyn's latest products include a disposable derma (keyes biopsy) punch made of quality stainless steel and a ribbed sturdy handle for the ejection of lodged specimen.

022727-1 Disposable Derma (Keyes Biopsy) Punch 1mm- STERILE

022727-2 Disposable Derma (Keyes Biopsy) Punch 2mm- STERILE

022727-3 Disposable Derma (Keyes Biopsy) Punch 3mm- STERILE

022727-4 Disposable Derma (Keyes Biopsy) Punch 4mm- STERILE

022727-5 Disposable Derma (Keyes Biopsy) Punch 5mm- STERILE



MedGyn's Disposable Scalpel

MedGyn's latest products include disposable scalpels made of quality stainless steel blades with sturdy, non-slip plastic handles.

022798-11 Disposable Scalpel #11- STERILE

022798-15 Disposable Scalpel #15- STERILE



MedGyn Word Bartholin's Catheter

For use following an incision and drainage of a Bartholin's cyst or abscess. Catheter kit includes balloon catheter, syringe for balloon inflation and scalpel.

022719 MedGyn Word Bartholin's Catheter (6/Box)



Other

040065 Video Print Paper - 110mm x 20mm rolls
(10/Box)

DISPOSABLES



MedGyn Single-Use Forceps

029701-D Tischler Biopsy Single-Use Forceps, 9 3/4" (10/Pack)



029702-D Baby Tischler Biopsy Single-Use Forceps, 9 3/4" (10/Pack)



029703-D Kevorkian-Young Biopsy Single-Use Forceps, 9 3/4" (10/Pack)



030112-D Alligator Single-Use Forceps 8" (10/Pack)



MedGyn Rotating Biopsy Forceps

MedGyn's Rotating Biopsy Punch is our latest single-use device with a plastic handle and stainless steel shaft to sample cervical tissue.

030100-D Disposable Rotating Biopsy Forceps



MedGyn Needle Extender

An effective means of extending ordinary hypodermic needles when longer needles are required.

022600 MedGyn Disposable Needle Extender



MedGyn Aspiration Kits

MedGyn Aspiration Kit and Aspiration Kit IV are effective devices for simultaneous endometrial curettage and aspiration.

022514 MedGyn Aspiration Kit with 4mm Cannula-Curette (10/Pack)

022515 MedGyn Aspiration Kit with 5mm Cannula-Curette (10/Pack)

022516 MedGyn Aspiration Kit with 6mm Cannula-Curette (10/Pack)

022517 MedGyn Aspiration Kit with 7mm Cannula-Curette (10/Pack)

022511 MedGyn Aspiration Kit IV with 4mm, 5mm, 6mm and 7mm Cannula-Curette (10/Pack)



MedGyn Aspiration Kit Double Valve

MedGyn Aspiration Kit Double Valve is used for uterine aspiration/uterine evacuation in obstetrics and gynecologic procedures. The device fits the MedGyn Flexible Karman Curettes and a wide variety of MedGyn Rigid and Straight Curettes.

022526 MedGyn Aspiration Kit Double Valve



MedGyn s IUD Retriever

MedGyn s IUD Retriever is designed to assist in the removal of the intrauterine device from the uterus.

022366 Disposable IUD Retriever (25/Box)



MedGyn Os Finders

Ideal for cervical dilation, uterine sounding or other intrauterine procedures. Made of semi-flexible plastic and packaged individually sterile.

022731 MedGyn Disposable Os Finder (50/Pack)

DISPOSABLES



022310



022450

022352

022343

022340

022341

MedGyn Collection Sets, Collection Jars and Bags

MedGyn Disposable Collection Sets include 6' tubing (either 3/8" or 1/2" I.D.) attached to a control handle with slip-ring for finger-tip vacuum control. The Collection Sets are individually packaged and sterile.

Collection Sets & Bottles

- 022310 3/8" Collection Set with Handle (50/Pack)
- 022310NS 3/8" Collection Set w/Handle & Adapter (Non-Sterile)
- 022310-WO 3/8" Collection Set with Handle (Without Adapter, Sterile) (50/Pack)
- 022314 3/8" Collection Set x 8 Ft with Handle & Adaptor
- 022312 1/2" Collection Set with Handle & Adapter (25/Pack)
- 022312NS 1/2" Collection Set w/Handle & Adapter (Non-Sterile)
- 022331 3/8" Disp. Collection Jar - with Cap (10/Box)
- 022332 1/2" Disp. Collection Jar - with Cap (10/Box)
- 022323 Disp. Tops Large Non-Conductive (10/Box)
- 022324 Disp. Tops Small Non-Conductive (10/Box)

Metal Swivel Handles

Stainless steel, autoclavable, with slip-ring control

- 022450 3/8" Handle for Curettes up to 14mm with small base
- 022451 1/2" Handle for 14mm, 15mm & 16mm Curettes
- 022340 Stainless Steel 3/8" Female Fitting
- 022341 Stainless Steel 3/8" Male Fitting
- 022343 Stainless Steel 1/2" Male Fitting
- 022352 Stainless Steel 3/8" Handle Extender
- 022353 Stainless Steel 1/2" Handle Extender


Accessories

- 022316-3/8-6 Tubing 3/8" Precut to 6'
- 022311 Tubing 3/8", Custom Cut Per Foot
- 022316-1/2-6 Tubing 1/2" Precut to 6'
- 022313 Tubing 1/2", Custom Cut Per Foot
- 022317 Tubing 16" Connector
- 022318 Tubing 20" Connector
- 022317-ADAP Tube Fitting - Male 3/8" Plastic
- 022317-INT Tube - Interconnecting with Adapter (10/Box)
- 022319-ADAP Tube Fitting - Male 1/2" Plastic
- 022350 Plastic Adapter for 1/2" Base Curettes (Sterile)
- 022351 Plastic Adapter for 1/2" Base Curettes (Non-Sterile)
- O-RING-FAD 3/8" O Ring for Female Adapter (10/Pack)
- 022320 Disposable Filter Assembly (10/Box)
- ADAP-REDC 1/2" to 3/8" Adapter for Aspiration Handle (For Small Curette) (10Box)
- 022354 3/8" to 1/2" Adapter for Aspiration Handle (For Large Curette)

Permeable Gauze Sacks

Gauze collection sack easily attaches to collection canister for simple and efficient separation of tissue from aspirated fluid.

- 022460 Gauze Sacks Tissue Trap (10/Box)
- O-RING3/8 O Ring for 3/8" Sack (10/Box)
- O-RING1/2 O Ring for 1/2" Sack (10/Pack)



MedGyn Vacuum Cannulas-Curettes

MedGyn offers a complete line of high-quality, disposable vacuum cannulas-curettes in the industry. Available in Rigid Curved, Rigid Straight, Flexible and Flexible Karman options, our cannulas-curettes are individually packaged and sterilized.

Rigid Curved Cannula-Curette

022106 6mm (50/Pack)	022113 13mm - 3/8 Base (25/Pack)
022107 7mm (50/Pack)	022114 14mm - 1/2 Base (25/Pack)
022108 8mm (50/Pack)	022144 14mm - 1/2 Base w/Adapter (25/Pack)
022109 9mm (50/Pack)	022145 14mm - 3/8 Base (25/Pack)
022110 10mm (50/Pack)	022115 15mm - 1/2 Base (25/Pack)
022111 11mm (50/Pack)	022116 16mm - 1/2 Base (25/Pack)
022112 12mm (50/Pack)	



Rigid Straight Cannula-Curette

022206 6mm (50/Pack)	022213 13mm - 3/8 Base (25/Pack)
022207 7mm (50/Pack)	022214 14mm - 1/2 Base (25/Pack)
022208 8mm (50/Pack)	022244 14mm - 1/2 Base w/Adapter (25/Pack)
022209 9mm (50/Pack)	022245 14mm - 3/8 Base (25/Pack)
022210 10mm (50/Pack)	022215 15mm - 1/2 Base (25/Pack)
022211 11mm (50/Pack)	022216 16mm - 1/2 Base (25/Pack)
022212 12mm (50/Pack)	



Flexible/Karman Cannula-Curette

022004K 4mm (10/Pack)
022005K 5mm (10/Pack)
022006K 6mm (10/Pack)
022007K 7mm (10/Pack)
022008K 8mm (10/Pack)
022009K 9mm (10/Pack)
022010K 10mm (10/Pack)
022011K 11mm (10/Pack)
022012K 12mm (10/Pack)

Flexible Cannula-Curette with Dual Port Design

022004 4mm (50/Pack)
022005 5mm (50/Pack)
022006 6mm (50/Pack)
022007 7mm (50/Pack)
022008 8mm (50/Pack)
022009 9mm (50/Pack)
022010 10mm (50/Pack)
022011 11mm (50/Pack)
022012 12mm (50/Pack)



DISPOSABLES



MedGyn Laminaria

Laminaria (type of hygroscopic dilator) is a 100% natural sterile cervical dilator made from a sea-grown plant containing no synthetic materials. Safely utilized for over a century, Laminaria provides gentle, gradual dilation. After harvesting, the plant is dehydrated. Once inserted into the cervical canal, Laminaria absorbs liquid, swelling to several times its initial dried diameter, gently dilating the cervix.

Standard, approximate length, 60mm

- 021002 2mm Extra Small/Extra Thin (25/Box)
- 021003 3mm Small (25/Box)
- 021004 4mm Medium (25/Box)
- 021005 5mm Large/Thick (25/Box)
- 021006 6mm Extra Large/Extra Thick (25/Box)
- 021008 8mm Jumbo (25/Box)
- 021010 10mm Extra Jumbo (25/Box)

Extra Long, approximate length, 70mm

- 021014 4mm Long 70mm (25/Box)
- 021015 5mm Long 70mm (25/Box)
- 021016 6mm Long 70mm (25/Box)
- 021018 8mm Long 70mm (25/Box)
- 021003L Laminaria #3 Long-String
- 021004L Laminaria #4 Long-String
- 021005L Laminaria #5 Long-String
- 021006L Laminaria #6 Long-String
- 021008L Laminaria #8 Long-String

With buttons to prevent over-insertion

- 021022 2mm Extra Small with Buttons (12/Box)
- 021023 3mm Small with Buttons (12/Box)
- 021024 4mm Medium with Buttons (12/Box)
- 021025 5mm Large with Buttons (12/Box)
- 021026 6mm Extra Large withButtons (12/Box)

Laminaria Inserter

- 031201 MedGyn Laminaria Inserter
- 031202 McMahon Laminaria Inserter



MedGyn Pregnancy Test Kits

MedGyn One-Step Pregnancy Test Kit is an easy-to-use clinical pregnancy test that delivers fast and accurate results. Our tests have built-in test controls and highly visible color.

022704 MedGyn Pregnancy HCG Cassette Urine Test (25/Box)

022705 MedGyn Pregnancy HCG Combo Urine/Serum Test (25/Box)

022706 MedGyn Pregnancy HCG Dipstick Urine Test (25/Box)



MedGyn Amnio Hook

MedGyn Amnio Hook is designed to perforate the amniotic membrane with a rigid talon tip.

020301 MedGyn Disposable Amnio Hook (50/Pack)



MedGyn Umbilical Cord

MedGyn Umbilical Cord Clamp has grooved clamping and lock to secure the umbilical cord prior to cutting.

020302 MedGyn Disposable Umbilical Cord Clamp (100/Pack)



MedGyn Uterine Manipulator

This device is used for laparoscopic tubal occlusion, diagnostic laparoscopy, mini laparotomy, salpingoplasty and fertility exams. Pear-shaped balloon provides optimal hold and a perfect seal at the internal OS. Smaller cannula size allows for easier insertion and comfort.

022701 MedGyn Uterine Manipulator(12/Box)

DISPOSABLES



MedGyn LED Speculum

MedGyn LED Speculum has a cordless LED light to improve visualization of the exam area.

020227 MedGyn Graves/Pederson LED Large (80/Box)

020226 MedGyn Graves/Pederson w/ LED Medium (100/Box)

020225 MedGyn Graves/Pederson w/ LED Small (100/Box)

020225STRL MedGyn Sterile Graves/Pederson w/ LED Small (100/Box)

020226STRL MedGyn Sterile Graves/Pederson w/ LED Medium (100/Box)

020227STRL MedGyn Sterile Graves/Pederson w/ LED Large (80/Box)



Large



Medium



Small



MedGyn Vaginal Speculum

MedGyn Disposable Speculum is available in both Graves and Pederson styles and easily locks in a variety of positions to reduce examination time.

020205 MedGyn Disposable Graves Large (100/Box)

020204 MedGyn Disposable Graves Medium (100/Box)

020201 MedGyn Disposable Pederson Small (100/Box)



Large



Medium



Small



MedGyn Balloon Catheter

MedGyn Postpartum Balloon Catheter is a disposable, multiple lumen catheter attached to an inflatable balloon system designed to provide tamponade for controlling hemorrhage from the uterus and vagina.

022703 MedGyn Balloon Catheter



MedGyn HSG Catheters

MedGyn HSG Catheter is used to inject contrast media or saline into the uterus and fallopian tubes for hysterosalpingography and sonohysterography procedures. It has a latex-free balloon catheter with a rigid insertion sheath.

022726 MedGyn HSG Catheters - 5Fr (10/Box)

022728 MedGyn HSG Catheters - 7Fr (10/Box)



MedGyn HVP-100 Hand-Held Vacuum Delivery Pump and Vacuum Assisted Delivery Cups

MedGyn Hand-Held Vacuum Delivery Pump may be used to safely assist medical personnel during vaginal deliveries. MedGyn Vacuum Delivery Cups offer optimal adhesion to ensure proper vacuum levels.

020322 MedGyn Vacuum Delivery Cups 50mm with Release Valve

020323 MedGyn Vacuum Delivery Cups 60mm with Release Valve

020316 MedGyn Vacuum Delivery Cups 70mm with Release Valve

020320 Vacuum Delivery Cup 50mm

020321 Vacuum Delivery Cup 60mm

016160 MedGyn HVP-100 Hand-Held Vacuum Delivery Pump

DISPOSABLES



MedGyn Endocervical Curette

The MedGyn Endocervical Curette provides an efficient sampling of the cervical canal to test for abnormal, precancerous conditions, or cervical cancer.

The Kevorkian-Young endocervical curette is a sterile, single-use device with a basket and a square tip.



022732 MedGyn Disposable Endocervical Curette (25/Box)




MedGyn IUI Catheters

MedGyn Insemination Catheter is a sterile, flexible catheter designed for use in intrauterine insemination. Tip of catheter is rounded and beveled for patient comfort.

022723 MedGyn Intrauterine Insemination Catheter - Regular/Straight (25/Box)

022724 MedGyn Intrauterine Insemination Catheter - Mini (25/Box)

022725 MedGyn Intrauterine Insemination Catheter - Curved (25/Box)



Regular

Mini

Curved



MedGyn Ezinject Catheter

A simple and effective way to allow visibility of the uterine cavity.

Balloonless catheter allows for a more complete view

Long tip at distal end improves placement in the cervical canal

022797 MedGyn Ezinject Catheter (10/box)





MedGyn Uterine Sounds

MedGyn Disposable Uterine Sounds are made of malleable plastic with centimeter graduations for accurate measurement of uterine cavity.

020250 MedGyn Disposable Uterine Sounds O.D. 3mm 10" Length with Marking (25/Box)



MedGyn Cell Sweep®

The MedGyn Cell Sweep® makes pap smears easier than ever by simultaneously sampling both the endocervix and ectocervix, eliminating the need for both a cytobrush and spatula.

- Optimal cell harvesting and transfer
- Separation of cells in one-slide technique
- Minimizes false negatives, red blood cells, and blood artifacts

- 022362 MedGyn Cell Sweep® Slide Holder (100/pack)
- 022361 MedGyn Cell Sweep® Non-Sterile (100/pack)
- 022361K MedGyn Cell Sweep® Non-Sterile with Kit (Slide and Slide Holder) (100/pack)
- 022363 MedGyn Cell Sweep® Sterile (100/pack)
- 022363K MedGyn Cell Sweep® Sterile with Kit (Slide and Slide Holder) (100/pack)



MedGyn Cytobrush and Pap Cell

- 022710 MedGyn Cytobrush (100/Box)
- 022364 MedGyn Pap Cell (50/Box)

DISPOSABLES



MedGyn Endosampler™

The unique design of the MedGyn Endosampler™ allows for superior accuracy and a more reliable sample. Available with one or two holes for the right amount of suction for each patient.

- 3mm curette and 10cc syringe
- Highly-visible graduated markings
- Semi-rigid curette
- Curved for improved manipulation
- Prevent back-flow with syringe lock-spring
- Sharp slot in curette tip
- Closed system prevents sample contamination



Endosampler™ EndosamplerII™

022720 MedGyn Endosampler™ (25/Box)

022722 MedGyn EndosamplerII™ (25/Box)



MedGyn Pipette™

The MedGyn Pipette™, Pipette II™, and Pipette IV™ are unique, sterile, disposable cannulas for sampling the endometrium.

- Semi-rigid and tapered for easy insertion and safety
- Flexible enough collect tissue from all areas

022721 MedGyn Pipette 3.1mm O.D. (25/box)

022734 MedGyn Pipette II 3.1mm O.D. (25/box)

022730 MedGyn Pipette IV 3.1mm O.D. (25/box)



Pipette™



Pipette II™



Pipette IV™



MedGyn Mucous Sampler

The MedGyn Mucous Sampler is a single-use, sterile, endocervical aspirator.

- 23 cm long, clear, flexible polypropylene tube
- 2.2mm diameter circular opening
- Distal end rounded for safety



022729 MedGyn Disposable Mucous Sampler (25/pack)



MedGyn LED-100 Examination Lamp

MedGyn LED-100 Examination Lamp is used to assist in performing clinical and gynecological exams. The Lamp has a flexible, chrome-plated, top base with four casters for excellent maneuverability and stability. The LED-100 provides three times better light output than standard lamps with exceptional longevity and brightness

- Highly flexible LED light is easy to maneuver
- True tissue color rendition and pinpoint accuracy
- Power dual voltage 110V or 240V, 60/50Hz
- Illuminance of 20,000 lux (Distance of 50 cm)
- Minimum expected life of LED is 35,000 hours
- Net weight: 5 kg (10.2 lb)
- Variable height of 110cm and 150 cm (95.5lb x 130lb)
- Electrical Standards: NBR-IEC 60.601-1 and NBR-IEC 60.601-1-2



010551 MedGyn LED-100 Examination Lamp
(Includes Rolling base with 4 Casters)



MedGyn Mini-Microscope

MedGyn Mini-Microscope is a high-magnification, illuminated microscope designed for cytological review of cervical or endometrial cells.

- Compact design
- Up to 60x magnification
- Linear Field > 4mm
- Working distance 15 18mm

017020 MedGyn Mini-Microscope

VACUUM PUMPS



Senator Vacuum Pump

Senator Vacuum Pump is made of high-quality components and with Swiss mechanical precision to guarantee excellent functionality. The Senator Vacuum Pump is ideal for removal of secretions, vacuum extractions, uterine suction and for pleural drainage.

- Suction capacity of 30 L/min
- Suction strength can be adjusted continuously and precisely
- Compact and user-friendly pump with large, integrated recessed grip
- Reliable 3-guard protective system with float levers integrated in bottle lids and hydrophobic bacterial filter



016100 Senator Vacuum Pump (Accessories and Trolley Sold Separately)

016150 Senator Vacuum Pump System (Includes Footswitch, Trolley, and Bottle Holder)
110 Volt Only



MedGyn MVP-100 Vacuum Pump

MedGyn MVP-100 Vacuum Pump is designed for safety and convenience. It features a push-switch linked directly to the pump and a trap bottle to control overflow.

- Applied oil-less piston vacuum pump
- Two collection bottles of 3,000 cc each
- Built-in anti-overflow device
- Push-switch (red button) linked directly to pump for ease of checking pump status
- Control low-voltage for user safety (under DC 12v)

016155 MedGyn MVP-100: Vacuum Pump
(Includes Foot Switch, Collection Bottle)



MedGyn Portable Vacuum Pump

- Oil-less piston requiring no lubrication and no maintenance
- 1000 mL autoclavable collection jar
- 12V rechargeable battery for two hours of operation

016156 MedGyn Portable Vacuum Pump
022316 Suction Tube 2m for MedGyn Portable Pump
022321 Suction Filter with Side Tube for MedGyn Portable Pump
022322 Suction Bottle with Bottle Cap for MedGyn Portable Pump




Eco 3 Ultrasound

12" LED Monitor, 2 Probe Connectors,
3.5MHz Convex Probe
Ergonomic Shape
USB Connection
THI (Tissue-Harmonic Imaging)
SRA (Speckle-Reduction Algorithm)

015090 Eco 3 Ultrasound System

Eco3 Ultrasound Accessories

- 015091 3.5MHz Convex Probe C3-A for Eco-3
- 015092 5.3 - 10.0 MHz Linear Probe L7M-A for Eco-3
- 015093 5.3 - 11.0 MHz Linear Probe L7S-A for Eco-3
- 015094 3.0MHz Phased Array Probe P3-A for Eco-3
- 015095 6.0MHz Trans-vaginal Probe V6-A for Eco-3
- 015096 7.5MHz Trans-rectal Probe R7-A for Eco-3
- 015097 6.0MHz Micro-convex Probe (Pediatrics) MC6-A for Eco-3
- 015098 3.0MHz Micro-convex Probe MC3-A for Eco-3
- 015099 5.0MHz Micro-convex Probe MC5V-A for Eco-3
- 015084 4D Complete Package (Includes: 4.0 MHz 4D Volume Probe, 4D Motor Board , 4D Software)
- 015083 DICOM 3.0
- 015100 DICOM for Eco-3



Q5 Color Ultrasounds

Q5: 15" LCD Monitor, 2 Probe Connectors,
Convex and Linear Probes
Double-Phase Digital Beam Forming (DPDBF)
Tissue-Harmonic Imaging (THI) on all probes
Speckle-Reduction Algorithm (SRA)
Memory Storage via 2 USB Ports
Built-in Software with Measurement and Reporting Functions

Q5 Ultrasound Accessories

- 015076 3.5 MHz Convex Probe
- 015077 7.5 MHz Linear Probe
- 015078 6.0 MHz Trans-vaginal Probe
- 015079 7.0 MHz Trans-vaginal Probe (180°)
- 015080 3.0 MHz Micro-convex Probe (adult cardiac)
- 015081 5.0 MHz Micro-convex Probe
- 015082 6.0 MHz Pediatric Micro-convex Probe

DOPPLERS



MedGyn Ultrasonic Pocket Doppler

MedGyn Pocket Doppler combines reliable and accurate detection with easy-to-use, light-weight design. Available with rechargeable batteries and probe frequency options.

- Easy-to-read, back-lit LCD display
- Built-in recorder
- Accurate FHR detection with sound
- Automatic and manual FHR counting
- Interchangeable waterproof probes



- 017100 MedGyn Ultrasonic Pocket Doppler (2MHz Probe)
- 017101 Ultrasonic Pocket Doppler with 2MHz Probe
- 017126 Ultrasonic Pocket Doppler with 3MHz Probe
- 017102 Ultrasonic Pocket Doppler with 2MHz Probe Rechargeable
- 017127 Ultrasonic Pocket Doppler with 3MHz Probe Rechargeable

Ultrasonic Pocket Doppler Accessories

- 017128 2MHz Probe
- 017129 3MHz Probe
- 017130 4 MHz Probe
- 017132 5 MHz Probe
- 017131 8 MHz Probe



SD5 Table-Top Doppler (Wired)

MedGyn Ultrasonic Table-Top Doppler

MedGyn Table-Top Dopplers are innovative additions to our Doppler category in a portable, light-weight package. Available in both wireless/infrared and wired models.

- Easy-to-read, back-lit LCD display
- Built-in recorder
- Accurate FHR calculation
- Portable lightweight design



SD6 Table-Top Doppler (Wireless)

- 017150 SD5 Table-Top Doppler with 2MHz probe, wired
- 017151 SD6 Table-Top Doppler with 2MHz probe, wireless

Ultrasonic Table-Top Doppler Accessories

- 017155 2MHz Probe, Wired for SD5
- 017157 3MHz Probe, Wired for SD5
- 017156 2MHz Probe, Wireless for SD6
- 017158 3MHz Probe, Wireless for SD6
- 017159 Carry Bag for SD5/SD6
- 017160 Rolling Stand for SD5/SD6
- 017199 Doppler Gel 8oz



F6 Fetal Monitor



F9 Fetal Monitor

MedGyn Fetal Monitors "F" Series

MedGyn's F6 & F9 Fetal and Maternal Monitor series combine simultaneous mother and baby monitoring. These monitors offer powerful patient information management and centralized monitoring system compatibility with a color display.

- 10.2" foldable wide color screen
- 24-hour CTG storage
- Built-in thermal array printer
- Data management software for real-time data transmission to PC
- Intrapartum monitoring (DECG/IUP) (Optional)
- The F6 Express offers the same superior features as the F6 System but with maternal monitoring included

017118 MedGyn F6 Fetal Monitor (FM10)

017119 MedGyn F6 Fetal Monitor (FM10) Twin FHR

017120 MedGyn F6 Express Fetal Monitor (FM10)

The F9 Fetal Monitor offers the same superior features as the F6 Fetal Monitor plus a foldable 12" color screen

017122 MedGyn F9 Fetal Monitor

017123 MedGyn F9 Express Fetal Monitor

Fetal Monitor Accessories

017105 Strap for Fetal Monitor F6

017106 Battery for Fetal Monitor F6

017191 Recording Paper for Fetal Monitor F6


017192 Recording Paper for Fetal Monitor F9

017125 Fetal Monitor Stand (Adjustable with Rolling Base)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)


ELECTROSURGERY



Lletz Kogan Endocervical Speculum

038057 Speculum 2-4mm, 11" (279mm)

038058 Speculum 4-8mm, 11" (279mm)



Lletz Tenaculum


038070 Single-Tooth Tenaculum 9 3/4" (248mm)



Lletz Tenaculum Hook

038080 Emmet Tenaculum 9" (229mm)



038081 Iris Tenaculum Hook 9" (229mm)



Lletz Dressing Forceps

038082 Lletz Dressing Forceps Small 8" (203mm)


038092 Lletz Dressing Forceps Large 10" (254mm)



Lletz Tissue Forceps

038090 Lletz Tissue Forceps Small 8" (203mm)

038091 Lletz Tissue Forceps Large 10" (254mm)



Lletz Foerster Forceps

038071 Lletz Foerster Forceps - Straight 9 1/2" (241mm)

038072 Lletz Foerster Forceps - Curved 9 1/2" (241mm)





Pederson Vu-More Speculum, (40mm) Opening

- 038100 Medium 4½" (114mm)
- 038101 Large 4¾" (124mm)
- 038102 Extra Large 5" (149mm)



Pederson Vu-Max Speculum, (70mm) Opening

- 038110 Medium 4½" (114mm)
- 038111 Large 4¾" (124mm)
- 038112 Extra Large 5" (149mm)



Weisman Open Side Speculum Closed Tube

- 038050 Right Side Open Medium 1⅝" x 4" (35mm x 102mm)
- 038050-ANG Right Side Open Medium Angled 1⅝" x 4" (35mm x 102mm)
- 038052 Right Side Open Large 1⅝" x 4¾" (35mm x 121mm)
- 038052-ANG Right Side Open Large Angled 1⅝" x 4¾" (35mm x 121mm)
- 038051 Left Side Open Medium 1⅝" x 4" (35mm x 102mm)
- 038051-ANG Left Side Open Medium Angled 1⅝" x 4" (35mm x 102mm)
- 038053 Left Side Open Large 1⅝" x 4¾" (35mm x 121mm)
- 038053-ANG Left Side Open Large Angled 1⅝" x 4¾" (35mm x 121mm)



Collins Speculum Closed Tube

- 038047 Small 1¼" x 3¾" (32mm x 95mm)
- 038048 Medium 1⅝" x 4" (35mm x 102mm)
- 038049 Large 1⅝" x 4½" (41mm x 114mm)



Four-Way Lateral Speculum

- 038093 Medium ⅞" x 4¼" (22mm x 108mm)
- 038094 Large ⅞" x 5¼" (22mm x 133mm)



Lateral Vaginal Wall Retractor

- 038060 Closed, 10" (254mm), Blade Size 1" x 2½" (25mm x 64mm)
- 038061 Open ½" (13mm) 10" (254mm), Blade Size 1" x 2½" (25mm x 64mm)
- 038062 Open 1" (25mm) 10" (254mm) Blade Size 1" x 2½" (25mm x 64mm)
- 038063 Wide Open, 10" (254mm) Blade Size 1" x 2½" (25mm x 64mm)



Klopfer Speculum

- 038054 1⅝" x 3½" (35mm x 89mm)

ELECTROSURGERY

MedGyn Lletz Instruments - Blue Coated

All MedGyn Lletz instruments have a durable coating to insulate against transmission of electric current, ensuring the highest level of protection for patients and physicians.



Graves Speculum Closed

- 038030 Small, $\frac{3}{4}$ " x 3" (19mm x 76mm)
- 038031 Medium, $1\frac{3}{8}$ " x 4" (35mm x 102mm)
- 038032 Large, $1\frac{3}{8}$ " x $4\frac{3}{4}$ " (35mm x 121mm)
- 038033 Extra Large, $1\frac{1}{2}$ " x $6\frac{1}{2}$ " (38mm x 165mm)
- 030900-LTZ Extra Long Snowman, Clinton $1\frac{1}{2}$ " x 6" (38mm x 152mm)



Graves Vu-More Speculum Closed

- 038055 Medium $1\frac{1}{4}$ " x 4" (32mm x 102mm)
- 038056 Large $1\frac{3}{8}$ " x $4\frac{1}{2}$ " (35mm x 114mm)



Graves Vu-Max Speculum Closed Tube

- 038046 Large $1\frac{3}{8}$ " x $4\frac{1}{2}$ " (35mm x 114mm)



Graves Vu-More Speculum, (40mm) Opening

- 038115 Small $3\frac{3}{8}$ " (98mm)
- 038116 Medium $4\frac{1}{2}$ " (114mm)
- 038117 Large $4\frac{3}{8}$ " (124mm)
- 038118 Extra Large $5\frac{7}{8}$ " (149mm)
- 038119 Extra Large $6\frac{1}{8}$ " (161mm)
- 038120 Extra Large $6\frac{3}{8}$ " (169mm)
- 038121 Extra Large 7" (178mm)




Graves Vu-Max Speculum, (70mm) Opening

- 038125 Small $3\frac{3}{8}$ " (98mm)
- 038126 Medium $4\frac{1}{2}$ " (114mm)
- 038127 Large $4\frac{3}{8}$ " (124mm)
- 038128 Extra Large $5\frac{7}{8}$ " (149mm)
- 038129 Extra Large $6\frac{1}{8}$ " (161mm)
- 038130 Extra Large $6\frac{3}{8}$ " (169mm)
- 038131 Extra Large 7" (178mm)

MedGyn Disposable Electrodes

Made of high-quality Tungsten wire, for superior shape and integrity, MedGyn offers a full range of shapes and sizes to suit most electrosurgical procedures.



- 028067 Disposable Loop Electrode 5mm x 5mm (5/Pack)
- 028069 Disposable Loop Electrode 10mm x 5mm (5/Pack)
- 028070 Disposable Loop Electrode 10mm x 7mm (5/Pack)
- 028071 Disposable Loop Electrode 10mm x 8mm (5/Pack)
- 028055 Disposable Loop Electrode 10mm x 10mm (5/Pack)
- 028056 Disposable Loop Electrode 10mm x 13mm (5/Pack)
- 028057 Disposable Loop Electrode 15mm x 5mm (5/Pack)
- 028087 Disposable Loop Electrode 15mm X 7mm (5/Pack)
- 028058 Disposable Loop Electrode 15mm x 8mm (5/Pack)
- 028068 Disposable Loop Electrode 15mm x 10mm (5/Pack)
- 028076 Disposable Loop Electrode 15mm x 13mm (5/Pack)
- 028059 Disposable Loop Electrode 15mm x 15mm (5/Pack)
- 028060 Disposable Loop Electrode 20mm x 8mm (5/Pack)
- 028063 Disposable Loop Electrode 20mm x 10mm (5/Pack)
- 028061 Disposable Loop Electrode 20mm x 12mm (5/Pack)
- 028062 Disposable Loop Electrode 20mm x 13mm (5/Pack)
- 028064 Disposable Loop Electrode 20mm x 15mm (5/Pack)
- 028065 Disposable Loop Electrode 20mm x 20mm (5/Pack)
- 028066 Disposable Loop Electrode 25mm x 10mm (5/Pack)

- 028041 Disposable Ball Electrode 3mm (5/Pack)
- 028042 Disposable Ball Electrode 5mm (5/Pack)
- 028081 Disposable Blade Electrode (5/Pack)
- 028053 Disposable Needle 5mm (5/Pack)

- 028075 Disposable Square Loop Electrode 5mm x 5mm (5/Pack)
- 028072 Disposable Square Loop Electrode 10mm x 4mm (5/Pack)
- 028073 Disposable Square Loop Electrode 10mm x 8mm (5/Pack)
- 028074 Disposable Square Loop Electrode 10mm x 10mm (5/Pack)



ELECTROSURGERY



MedGyn Electrosurgical Generator System

MedGyn Electrosurgical Generator Complete System includes the ES Generator, Smoke Evacuator, and trolley.

018011 MedGyn Electrosurgical Generator Complete System includes the ES Generator, Smoke Evacuator, and Trolley

018102 MedGyn Electrosurgical Generator Trolley



MedGyn Smoke Evacuators

Smoke filtration system minimizes airborne microbes and viruses while removing smoke from site for improved visualization for use during all electrosurgical procedures.

Smoke Evacuator Filter
Smoke Evacuator Tube
Reducer Fitting

018100 MedGyn Smoke Evacuator System (110V)

018109M MedGyn Smoke Evacuator System SA-01 (220V)

018109 MedGyn Smoke Evacuator (220V)



Smoke Evacuator Accessories

018105 Smoke Filter Tube - For MedGyn Smoke Evacuator (110V or 220V) (5/Box)

018106 Disposable Smoke Evacuator (DSE) Male tubing 8 (Non-Sterile) (50/Box)

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

ELECTROSURGERY



028030



028032



028033

- FS02-01R Double Foot Switch - 4Pin Plug Cable 4m
- FS01-02R Single Foot Switch - 3Pin Plug Cable 4m
- HD02-01D Disposable Twin Button Handle - 3Pin Plug Cable 3m
- HD01-01R Mono Polar Handle & Cable - Silicon Cable 4ø 3m For ES-200/300/400
- HD01-02R Mono Polar Handle & Cable - Silicon Cable 4ø 3m For ES-100
- PL02-04R Silicon Patient Plate - 240mm x 150mm Valley Lab Type Plug Cable 3m
- PL03-02D Disposable Patient Plate (2ea) - Adult Dual (120cm²) For ES-200/300/400
- PL03-03D Disposable Patient Plate (2ea) - Adult Dual (68cm²) For ES-100
- CA02-02R Patient Return Plate Cable - 6.3 Pin Plug Silicon Cable 3m For ES-100
- CA02-03R Patient Return Plate Cable - Valley Lab Type Plug Silicon Cable 3m For ES-200/300/400
- CA03-01R Ground Cable - Cable 4ø, 3m
- CA01-01R Bipolar Forceps Cable - 2pin Plug Cable 3m EN Type Connector
- BF01-01R Bipolar Forceps - Bayonet 17.8cm 1.0 Tip
- EL01-02D Electrode - Knife 2.4 * 70mm Black
- EL02-02D Electrode - Ball 5mm Black
- EL03-02D Electrode - Loop 6mm Black
- EL04-02D Electrode - Needle 2.4 * 70mm Black
- EL04-04D Electrode - Needle Angled 45° 2.4 * 70mm, Black
- ST03-01R Latex Strap Perforated - 600mm x 1250mm, 60 hole
- PL03-03D Disposable Patient Plate - Pediatric Dual (66cm²)
- BF01-02R Bipolar Forceps - Bayonet 15.2cm 1.0 Tip
- BF02-01R Bipolar Forceps - Straight 9.0cm 0.5 Tip
- BF02-02R Bipolar Forceps - Straight 12.1cm 1.0 Tip
- BF02-03R Bipolar Forceps - Straight 17.8cm 1.0 Tip
- BF03-01R Bipolar Forceps - Angled 9.0cm 0.5 Tip
- BF03-02R Bipolar Forceps - Angled 17.8cm 1.0 Tip
- CT01-01R Cart - W390xD480xH850(mm)

ELECTROSURGERY



ES-200, ES-300, ES-400



ES-100

MedGyn Electrosurgical Generators (ES Series)

The MedGyn ES-100, ES-200, ES-300, ES-400 Electrosurgical Generator Series meet all demands for a LEEP/ LLETZ procedure. MedGyn's generators feature both monopolar and bipolar functions to satisfy all surgical demands with safety, flexibility, reliability, and convenience.

- Pure cut and blend cut functions
- Contact coagulation and bipolar coagulation
- Displays output power digitally
- Allows operator to check output power easily during the procedure
- Power source: AC120V or AC230V, 50/60 Hz (manual switch)

- 018015 MedGyn ES-100: Electrosurgical Generator (110-230V)
(Includes Single Foot Switch, Disposable Twin Button Handle, Silicone Patient Plate)
- 018120 MedGyn ES-200: Electrosurgical Generator (220V)
(Includes Single and Double Foot Switch, Disposable Twin Button Handle, Silicone Patient Plate)
- 018130 MedGyn ES-300: Electrosurgical Generator (220V)
(Includes Single and Double Foot Switch, Disposable Twin Button Handle, Silicone Patient Plate)
- 018135 MedGyn ES-400: Electrosurgical Generator (220V)
(Includes Single and Double Foot Switch, Disposable Twin Button Handle, Silicone Patient Plate)
- 028030 Disposable Handswitch Pencil Single Use (10/Pack)
- 028032 Handswitch Pencil Reusable
- 028033 Disposable Grounding Pad with Cable
- 028031 Disposable Grounding Pads without Cable
- 028034 Reusable Cable only for Disposable Grounding pad



MedGyn MTA-100 Thermal Ablation System

In 2019, WHO issued new guidelines for the use of thermal ablation for cervical pre-cancer lesions. While they continue to recommend cryotherapy, the WHO has offered thermal ablation as an alternative method of treatment for specific environments. MedGyn is pleased to announce its latest device: MedGyn's MTA-100 Thermal Ablation System. The device is a portable, lightweight device for the treatment of pre-cancerous lesions optimized for 20, 40, or 60 second intervals of treatment.

- Optimized to offer 20, 40, or 60 second intervals for procedure
- Bright LED lights for clear viewing of cervical area.
- Sound and light activation for beginning and end of procedure.
- Digital display to detailing temperature at all stages of procedure
- Durable and ergonomic handle for ease of use
- 25mm tip available offering larger treatment area

018003 MedGyn MTA-100 Thermal Ablation System

* Patented

	018003/2500 Exocervical, 25mm, Flat
	018003/2507 Endo/Exocervical, Large
	018003/1910 Endo/Exocervical, Small
	018003/1905 Exocervical, Convex, 19mm x 5mm
	018003/1900 Exocervical, 19mm, Flat

019022 Plastic Probe Protection Bag

MedGyn MTA-100 Tips

Safe, interchangeable, autoclaveable tips are made from a conductive gold alloy, to allow for various pre-cancerous lesion treatment sizes.

CRYOTHERAPY

	019002/2507 Cryotherapy Tip GT-2507 Endo/Exocervical, Large
	019002/2500 Cryotherapy Tip GT-2500 Exocervical, 25mm Flat
	019002/1910 Cryotherapy Tip GT-1910 Endo/Exocervical, Small
	019002/1905 Cryotherapy Tip GT-1905 Exocervical, Convex (19mm X 5mm)
	019002/1900 Cryotherapy Tip GT-1900 Exocervical, 19mm Flat
	019002/0811 Cryotherapy Tip GT-0811 Skin Lesion, 8mm Diam, 45° angle
	019002/0800 Cryotherapy Tip GT-0800 Endocervical, Round
	019002/0519 Cryotherapy Tip GT-0519 Endocervical, Nulliparous
	019002/0507 Cryotherapy Tip GT-0507 Skin Lesion, 5mm Diam, 45° angle
	019002/0219 Cryotherapy Tip GT-0219 Micro, 2mm Diam
	019002/0832 Cryotherapy Tip GT-0832 Anorectal

MedGyn Cryotherapy Tips

Safe, interchangeable, autoclavable tips are made from a conductive gold alloy, resulting in very rapid freeze and defrost capabilities. Designed to allow for change-out without the need to shut off gas flow, reducing risk of blown O ring.



Carrying Case



Cryotherapy System

MedGyn MGC-200 Cryotherapy System

The MedGyn Cryotherapy System is designed by physicians and engineered by cryogenic experts to ensure safety, performance, and convenience. The 3-position single trigger delivers optimal cryosurgical temperatures with Nitrous Oxide (N₂O) or Carbon Dioxide (CO₂) gas. The system has standard-specific adaptors and a wide variety of autoclavable tips available.

- Single-hand control with 3 settings: freeze, off, and defrost
- Instant defrost
- Change tips without shutting off the gas
- Built-in regulators control pressure at tip for added safety and gas economy

019000 MedGyn Cryotherapy System MGC-200 (tips not included)


Cryotherapy Equipment & Accessories

- 019019 Carrying Case for Cryotherapy System
- 019020 Protective Cap for Cryotherapy System
- 019022 Plastic Probe Protection Bag
- 019021 Rubber Plug 6448K92 for MedGyn Cryo Tip
- 019030 450 Purifier Housing with 1/4 npt Male Thread Adapter
- 019031 541 Filter Cartridge for Purifier
- 019001 Cart for Cylinder, 6lb or 20lb
- 019003 "O" Ring for Cryo System (3/pack)
- 019008 French Adapter for CO₂
- 019016 French Adapter For N₂O

- 019009 German Adapter for CO₂
- 019012 German Adapter for N₂O
- 019010 British Adapter for CO₂
- 019011 British Adapter for N₂O
- 019015 Chinese Adapter for CO₂
- 019014 US Adapter for CO₂ (20lb)
- 019013 US Adapter for N₂O (20 lb)
- 019017 US adapter for N₂O (6lb)
- 019018 US adapter for CO₂ (6lb)



COLPOSCOPES



MedGyn AL-102 and AL-101 Colposcopes

MedGyn's classic stereoscopic colposcopes AL-102 and AL-101 are the industry standard. These models combine high-quality lens systems with easy maneuverability at an attractive price.

010200 MedGyn AL-102 Colposcope

- Interchangeable eyepieces, 10x, 13.5x, 18x
- 3-position illumination settings
- Video and camera compatible
- Green filter
- Locking caster wheels

010199 MedGyn AL-101 Colposcope

- 1-step magnification: 16x
- Crystal clear image
- 12.5x eyepieces, 1 fixed and 1 moveable
- Green filter
- Adjustable stand

AL-102 Accessories

010204 AL-102 10x Eyepiece (set of 2)

010205 AL-102 18x Eyepiece (set of 2)

010206 AL-102 13.5x Eyepiece (set of 2)



MedGyn AL-105 and AL-105Z Colposcopes

Traditional optical colposcopes compatible the the MedGyn Image & Data Management Software.

Five-step magnification 3x, 4x, 7x, 11x, 17x
LED light source with continuous adjustment
High-contrast images with excellent stereoscopic effect
Great depth of field
Green filters

010505 MedGyn AL-105Z Colposcope with Swing Arm

010500 MedGyn AL-105 Colposcope

AL-105 Accessories


090025 Power Cord #209-000 - UK

090026 Power Cord #199-000 - Europe

090027 Power Cord #364-010 - Italy

090028 Power Cord #456-010 - South Africa

COLOSCOPES



MedGyn Image & Data Management Software

Designed for use with our AL-106 series of colposcopes, MedGyn's Image & Data Management Software captures and manages colposcopy images in a self-contained patient record. The user-friendly software includes numerous templates and is easily customizable to your clinic's requirements. The software automatically creates a variety of reports which can include colposcopy images.

- Image annotation and measurement
- Full video capture
- Statistical analysis
- Complete system security, internal password protected
- Includes capture card, security key, and foot switch

010477 MedGyn Image and Data Management Software

010476 MedGyn Image and Data Management Software
with Laptop (monitor not included)

010475 MedGyn Image and Data Management Software
with Laptop (monitor Included)

MedGyn Colposcopy Instrument Set

All the instruments you need for a thorough colposcopic examination.

- Tischler Biopsy Punch
- Endocervical Speculum
- Graves Speculum
- Endocervical Curette
- Emmet Hook



030350 MedGyn Colposcopy Instrument Set



MedGyn Digital Video Colposcopes

The MedGyn AL-106B and AL-106 BWL are more affordable options to the AL-106 HD. The AL-106 BWL is a battery operated colposcope useful in remote areas with limited access to electricity.

- Battery operated and wireless
- High-resolution of 440K pixels
- LED light
- 1-54x magnification (36x for AL-106 BWL)
- Ultra-fast auto focus and manual focus option
- Video and image freeze
- Green light filters

- 010471 MedGyn AL-106 BWL Digital Video Colposcope, Battery Operated
- 010469 MedGyn AL-106B Digital Video Colposcope
- *Monitor Sold Separately

COLPOSCOPES



MedGyn AL-106 HD™ Digital Video Colposcopes

Different from the traditional optical colposcopes, the AL-106HD™ is a digital colposcope with an LED light source for ultra-bright visualization and a long life cycle. Its unique optical precision and multiple green light filters ensure the highest quality images and improved visualization of blood vessels.

- High-resolution: 1.43M px for AL-106HD, 440K for AL-106
- LED light
- 1-40x magnification
- Ultra-fast auto focus and manual focus option
- Video and image freeze
- Green light filters

- 010466-HD MedGyn AL-106 High Definition Digital Video Colposcope*
- 010466SA-HD MedGyn AL-106 High Definition Digital Video Colposcope with Swing Arm*
- 010010 Colposcope Monitor (Used for AL-101, AL-102, AL-401, AL-105, AL-106)*
- 010490 Colposcope Monitor Stand (Adjustable with Rolling Base)

Neurological Instruments


Catalog #	Description
MH1-226	Wartenberg Neurological Pin Wheel, 7" (178 mm)




MH1-226

Obstetrics and Gynecology Instruments - Forceps and Scissors

Catalog #	Description
MH30-1700	Heaney Hysterectomy Forceps, 8½" (213 mm), heavy pattern, single tooth, curved
MH30-1710	Heaney Hysterectomy Forceps, 8½" (215 mm), heavy pattern, double tooth, curved
MH30-1716	Heaney-Ballentine Hysterectomy Forceps, 8½" (216 mm), longitudinal serrations, single tooth, curved
MH30-1718	Ballentine Hysterectomy Forceps, 8¾" (220 mm), longitudinal serrations, single tooth, straight
MH30-2190	Braun Episiotomy Scissors, 5½" (140 mm), angled to side, guarded lower blade
MH30-2195	Braun Episiotomy Scissors, 8½" (216 mm), angled to side, guarded lower blade



MH30-1700, MH30-1710



MH30-1716, MH30-1718



MH30-2190, MH30-2195

Obstetrics and Gynecology Instruments - Biopsy Forceps

Catalog #	Description
MH30-1482	Kevorkian-Young Uterine Biopsy Forceps, 9 1/2" (241 mm) shaft, 3.5 x 8 mm bite, teeth on lower jaw
MH30-1485	Eppendorfer Uterine Biopsy Forceps, 8 1/2" (216 mm) shaft
MH30-1405	Thoms-Gaylor Uterine Biopsy Forceps, 8 1/2" (216 mm), angled shanks, curved jaws, 5 mm bite, interfitted cups
MH301442WL	Tischler Biopsy Forceps, 8" (203 mm), with lock
MH301443WL	Burke Baby Tischler Biopsy Forceps, 7 3/4" (197 mm), with lock



MH30-1482



MH30-1485



MH30-1405

Obstetrics and Gynecology Instruments - Biopsy Curettes

Catalog #	Description
MH30-1382	Kevorkian-Young Endocervical Biopsy Curette, 12 1/4" (310 mm), loop size 3 x 12 mm, with basket
MH30-1383	Kevorkian-Young Endocervical Biopsy Curette, 12 1/4" (310 mm), loop size 3 x 12 mm, without basket
MH30-1385	Heaney Uterine Biopsy Curette, 9" (229 mm), serrated loop



MH30-1382



MH30-1383



MH30-1385

Obstetrics and Gynecology Instruments

Catalog #	Description
MH30-997	Teale Uterine Vulsellum Forceps, 9" (230 mm), curved sideways, 3 x 4 teeth, 5 mm wide
MH30-1000	Schroeder Uterine Vulsellum Forceps, 10" (254 mm), straight
MH30-1025	Jacobs Uterine Vulsellum Forceps, 8 1/4" (210 mm), straight, serrated jaws, 2 x 2 teeth
MH30-1205-0	Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 0
MH30-1205-1	Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 1
MH30-1205-2	Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 2
MH30-1205-3	Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 3
MH30-1205-4	Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 4
MH30-1205-5	Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 5
MH30-1205-6	Sims Uterine Curette, 11" (279 mm), sharp blades on malleable shank, size 6




MH30-997



MH30-1000



MH30-1025



0

1

2

3

4

5


6




MH30-1205-0 to MH30-1205-6

Obstetrics and Gynecology Instruments - Tenaculum Forceps

Catalog #	Description
MH30-965	Schroeder-Braun Uterine Tenaculum Forceps, 9 3/4" (248 mm), straight
MH30-965ATR	Braun Tenaculum Forceps, 10" (254 mm), straight/squared jaw with non-overlapping atraumatic tips
MH30-966	Schroeder Uterine Tenaculum Forceps, 10" (254 mm), round jaw
MH30-966ATR	Schroeder Tenaculum Forceps, 10" (254 mm), rounded jaw with non-overlapping atraumatic tips
MH30-967	Hulka Tenaculum Forceps, 11" (279 mm), with 3.5 mm probe
MH30-985	Barrett Uterine Tenaculum Forceps, 7 1/4" (184 mm), straight




MH30-965



MH30-966



MH30-967



MH30-985

Obstetrics and Gynecology Instruments - Dilators

Catalog #	Description
MH30-561	Pratt Uterine Dilators, 11½" (292 mm), double end, set of 8 Dilators in Khaki Roll, sizes 13-15 Fr. (4.3 to 4.6 mm) to 41-43 Fr. (13.5-14.2 mm)
MH30-560-1315	Pratt Uterine Dilator, 11½" (292 mm), double end, 13-15 Fr. (4.3-4.6 mm)
MH30-560-1719	Pratt Uterine Dilator, 11½" (292 mm), double end, 17-19 Fr. (5.6-6.3 mm)
MH30-560-2123	Pratt Uterine Dilator, 11½" (292 mm), double end, 21-23 Fr. (7-7.6 mm)
MH30-560-2527	Pratt Uterine Dilator, 11½" (292 mm), double end, 25-27 Fr. (8.3-9 mm)
MH30-560-2931	Pratt Uterine Dilator, 11½" (292 mm), double end, 29-31 Fr. (9.6-10.2 mm)
MH30-560-3335	Pratt Uterine Dilator, 11½" (292 mm), double end, 33-35 Fr. (10.9-11.6 mm)
MH30-560-3739	Pratt Uterine Dilator, 11½" (292 mm), double end, 37-39 Fr. (12.2-12.9 mm)
MH30-560-4143	Pratt Uterine Dilator, 11½" (292 mm), double end, 41-43 Fr. (13.5-14.2 mm)



MH30-561



MH30-560-1315 to MH30-560-4143

Obstetrics and Gynecology Instruments - Dilators

Catalog #	Description
MH30-550	Hegar Uterine Dilators, 7 $\frac{3}{4}$ " (197 mm), double end, set of 8 Dilators in Khaki Roll, sizes 3-4 mm to 17-18 mm
MH30-555-34	Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 3-4 mm
MH30-555-56	Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 5-6 mm
MH30-555-78	Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 7-8 mm
MH30-555-910	Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 9-10 mm
MH30-555-1112	Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 11-12 mm
MH30-555-1314	Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 13-14 mm
MH30-555-1516	Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 15-16 mm
MH30-555-1718	Hegar Uterine Dilator, 7 $\frac{3}{4}$ " (197 mm), double end, 17-18 mm



MH30-550



MH30-555-34 to MH30-555-1718

Obstetrics and Gynecology Instruments - Dilators

Catalog #	Description
MH30-500	Hank Uterine Dilators, 10½" (267 mm), double end, set of 6 Dilators in Khaki Roll, Sizes 11-12, (3.6-4 mm), 13-14 (4.3-4.6 mm), 15-16 (5-5.3 mm), 17-18 (5.6-6 mm), 19-20 (6.3-6.6 mm), and 21-22 (7-7.3 mm) French
MH30-505-910	Hank Uterine Dilator, 10½" (267 mm), double end, 9-10 Fr. (3-3.3 mm)
MH30-505-1112	Hank Uterine Dilator, 10½" (267 mm), double end, 11-12 Fr. (3.6-4 mm)
MH30-505-1314	Hank Uterine Dilator, 10½" (267 mm), double end, 13-14 Fr. (4.3-4.6 mm)
MH30-505-1516	Hank Uterine Dilator, 10½" (267 mm), double end, 15-16 Fr. (5-5.3 mm)
MH30-505-1718	Hank Uterine Dilator, 10½" (267 mm), double end, 17-18 Fr. (5.6-6 mm)
MH30-505-1920	Hank Uterine Dilator, 10½" (267 mm), double end, 19-20 Fr. (6.3-6.6 mm)



MH30-500



MH30-505-910 to MH30-505-1920

Obstetrics and Gynecology Instruments - Specula

Catalog #	Description
MH30-10	Graves Vaginal Speculum, original model, small size, ¾" (19 mm) x 3" (76 mm)
MH30-15	Graves Vaginal Speculum, original model, medium size, 1⅜" (35 mm) x 4" (102 mm)
MH30-20	Graves Vaginal Speculum, original model, large size, 1½" (38 mm) x 4½" (114 mm)
MH30-30	Graves Open Sided Vaginal Speculum, medium size, 1⅜" (35 mm) x 4" (102 mm), wide angle blades
MH30-50	Pederson Vaginal Speculum, narrow blades, small size, ½" (13 mm) x 3" (76 mm)
MH30-55	Pederson Vaginal Speculum, narrow blades, medium size, ⅞" (22 mm) x 4" (102 mm)
MH30-60	Pederson Vaginal Speculum, narrow blades, large size, 1" (25 mm) x 4¾" (121 mm)
MH30-62	Pederson Vaginal Speculum, extra narrow blades, ⅝" (16 mm) x 4" (102 mm)



MH30-10 to MH30-20



MH30-50 to MH30-62

Suction Tubes

Catalog #	Description
MH2-104SS	Yankauer Suction Tube, 11½" (292 mm), with removable tip and tubing connection
MH2-109	Yankauer Pediatric Suction Tube, 8½" (216 mm), with removable tip, delicate
MH2-110SS	Andrews-Pynchon Suction Tube 10" (254 mm), delicate
MH10-310	Poole Suction Tube, 24 French, (8 mm), curved
MH10-312	Poole Suction Tube, 30 French, (10 mm), straight



MH2-104SS



MH2-110SS



MH10-310



MH10-312

OB/GYN

CARE AND CLEANING

TROUBLE SHOOTING

CORROSION

- Problem:** Debris and residue in box locks, serrations and ratchets.
- Solution:** Serrations, box locks and ratchets should be carefully cleaned and rinsed. A soft bristle brush is effective in these difficult areas.
- Problem:** Moisture remaining on the instruments within the sterile wrapped packs; excessive moisture in autoclave during sterilization and termination of cycle.
- Solution:** Preheat autoclave, drying time cannot be accelerated, open autoclave door slightly. Routinely inspect valves. Malfunctioning valves can precipitate corrosion.
- Problem:** Some local water supplies contain excessive quantities of alkali. This can manifest itself by depositing corrosive matter on the instruments.
- Solution:** It is generally accepted that cold distilled or demineralized water be used during sterilization. An alternative is to use the local water supply, while instituting a program of preventive maintenance. Carefully clean surfaces with SKLAR KLEEN™, using SKLAR POLISH™ to remove difficult stains and SKLAR LUBE™ for proper lubrication.

PITTING

- Problem:** Instruments are inevitably exposed to chemical compounds that initiate pitting. The most common are saline solutions and blood.
- Solution:** It is essential that instruments be rinsed with distilled or demineralized water and then cleaned routinely.
- Problem:** Some Non-Sklar detergents contain a chloride or have an acid base. If rinsing is not thorough, the residual detergent can precipitate pitting. Hydrochloric acid can be formed by the combination of chloride in the cleaner and steam, causing the instrument to pit.
- Solution:** Evaluate detergents without a chloride base. SKLAR KLEEN™ is recommended.
- Problem:** Instruments may be pitted when the detergent used has an excessively high pH. If instruments are not rinsed thoroughly and then autoclaved, the heat increases the effect of the solution on the instruments.
- Solution:** Use only detergents specifically developed for instrument cleaning. SKLAR KLEEN™ products are pH balanced.

BROWNISH STAIN

- Problem:** This stain is more obvious on dull finish instruments. Most likely, it is a deposit of chromium oxide. The thin hard layer protects stainless steel from corrosion precipitated by the atmosphere. The brownish film (sometimes dull blue) occurs naturally when stainless is heated. Occasionally this stain can be traced to the sterilizer. Improper cleaning compounds could contain poly phosphates that cause the copper components in the sterilizer to be deposited by electrolytic action on the stainless steel instruments.
- Solution:** Water supply to autoclave should be cold and demineralized. Use the correct amount of SKLAR KLEEN™ as recommended for best results.

BLUISH GRAY STAINS

- Problem:** Cold sterilizing solutions.
- Solution:** In the preparation of the solution make certain the manufacturer's directions are followed accurately. To minimize discoloration use distilled water and add a rust inhibitor. Using the solution beyond recommended time limitations will cause the solution to become corrosive.

RUST DEPOSIT

- Problem:** Sterilizing instruments of different metals in the same cycle will cause rusting. An electrolytic action carries the carbon particles from the exposed metal and deposits them on the stainless steel instruments. The rapid oxidation of molecules results in a surface film of rust.
- Solution:** The ideal technique would be to sterilize the two types of instruments separately. It is recommended that chrome plated instruments that indicate signs of peeling or flaking be replaced with stainless steel instruments.
- Problem:** The rust film on the surface of stainless steel instruments can be caused by chemicals in the detergent. Excess amounts of iron or other minerals from the local water supply may cause the rust deposits.
- Solution:** It is recommended that instruments be rinsed with distilled or demineralized water. Surface rust can be removed by light buffing or an application of SKLAR POLISH™.

LIGHT AND DARK COLORED SPOTS

- Problem:** Spotting caused by slow evaporation of water drops identified with the mineral content of the water.
- Solution:** Eliminate water droplets by adhering to autoclave manufacturer's operating instructions. Autoclave door should not be opened until steam has been completely vented.

PRIME Surgical-GYNECOLOGY & OBSTETRICS

(Surgical Instruments, Medical Equipments, Hospital Supplies)

PRIME SURGICAL INTERNATIONAL