



PRIME SURGICAL INTERNATIONAL

UK OFFICE

24 King Street
Manchester M2 4LG United Kingdom
Phone: +44 20 7946 0123
www.primesurgicalinternational.com

SINGAPORE OFFICE

20 Cecil Street, #05-03 PLUS
Singapore 049705
+65 9076 3821
Zia@primesurgicalinternational.com
info@primesurgicalinternational.com

PAKISTAN OFFICE

OFFICE # 206, 2ND FLOOR,
QUEEN CENTER, M.T. KHAN ROAD,
KARACHI, PAKISTAN
+92 300 2621673
ziarehman790@gmail.com

Damage Control Surgery - Advanced Instrument List

Code: PSI-SET-A1434

Category: Damage Control Surgery

Set Level: Advanced

SR	Specialty	Instrument Name	Size / Specification	Qty	Code
1	Emergency Trauma	PRIME Scalpel Handle	No. 4	2	
2	Emergency Trauma	PRIME Trauma Shears	19 cm	2	
3	Emergency Trauma	PRIME Mayo Scissors	17 cm straight	2	
4	Emergency Trauma	PRIME Metzenbaum Scissors	18 cm curved	2	
5	Emergency Trauma	PRIME Artery Forceps	Kelly 14 cm curved	8	
6	Emergency Trauma	PRIME Kocher Forceps	16 cm straight	6	
7	Emergency Trauma	PRIME Allis Forceps	15 cm	4	
8	Emergency Trauma	PRIME Needle Holder	Mayo-Hegar 18 cm	3	
9	Emergency Trauma	PRIME Towel Clamp	Backhaus 11 cm	6	
10	Emergency Trauma	PRIME Suction Tip	Yankauer	2	
11	Emergency Trauma	PRIME Retractor	Army-Navy	2	
12	Emergency Trauma	PRIME Skin Staple Remover/Forceps	standard	1	
13	Emergency Trauma	PRIME Rochester-Carmalt Forceps	20 cm curved	4	
14	Emergency Trauma	PRIME Rochester-Pean Forceps	18 cm	4	
15	Emergency Trauma	PRIME Lahey Forceps	18 cm	2	
16	Emergency Trauma	PRIME Bone Holding Forceps	Lane 20 cm	2	
17	Emergency Trauma	PRIME Periosteal Elevator	Mayo 17 cm	2	
18	Emergency Trauma	PRIME Volkman Retractor	4 prong	2	
19	Emergency Trauma	PRIME Self-Retaining Retractor	Weittlaner 13 cm	2	
20	Emergency Trauma	PRIME Suction Tube	Poole	2	
21	Emergency Trauma	PRIME Debridement Forceps	Adson 12 cm	2	



PRIME SURGICAL INTERNATIONAL

UK OFFICE

24 King Street
 Manchester M2 4LG United Kingdom
 Phone: +44 20 7946 0123
 www.primesurgicalinternational.com

SINGAPORE OFFICE

20 Cecil Street, #05-03 PLUS
 Singapore 049705
 +65 9076 3821
 Zia@primesurgicalinternational.com
 info@primesurgicalinternational.com

PAKISTAN OFFICE

OFFICE # 206, 2ND FLOOR,
 QUEEN CENTER, M.T. KHAN ROAD,
 KARACHI, PAKISTAN
 +92 300 2621673
 ziarehman790@gmail.com

Damage Control Surgery - Advanced Instrument List

SR	Specialty	Instrument Name	Size / Specification	Qty	Code
22	Emergency Trauma	PRIME Bowel Clamp	Doyen 23 cm	2	
23	Emergency Trauma	PRIME Tourniquet Clamp/Forceps	standard	2	
24	Emergency Trauma	PRIME Instrument Tray	trauma medium	1	
25	Emergency Trauma	PRIME Rib Spreader	Finochietto adult	1	
26	Emergency Trauma	PRIME Vascular Clamp	Satinsky 24 cm	2	
27	Emergency Trauma	PRIME Vascular Clamp	DeBakey assorted	4	
28	Emergency Trauma	PRIME Amputation Saw	Gigli saw with handles	1	
29	Emergency Trauma	PRIME Orthopedic Drill Guide	basic set	1	
30	Emergency Trauma	PRIME External Fixator Instrument Set	basic	1	
31	Emergency Trauma	PRIME Large Bone Retractor	Hohmann assorted	4	
32	Emergency Trauma	PRIME Hemostatic Clip Applicator	medium-large	1	
33	Emergency Trauma	PRIME Damage Control Laparotomy Pack	instrument support set	1	
34	Emergency Trauma	PRIME Thoracic Suction Tube	large bore	2	