



PRIME SURGICAL INTERNATIONAL

UK OFFICE

24 King Street
Manchester M2 4LG United Kingdom
Phone: +44 20 7946 0123
www.primesurgicalinternational.com

SINGAPORE OFFICE

20 Cecil Street, #05-03 PLUS
Singapore 049705
+65 9076 3821
Zia@primesurgicalinternational.com
info@primesurgicalinternational.com

PAKISTAN OFFICE

OFFICE # 206, 2ND FLOOR,
QUEEN CENTER, M.T. KHAN ROAD,
KARACHI, PAKISTAN
+92 300 2621673
ziarehman790@gmail.com

Endovascular Aneurysm Repair (EVAR) - Basic Instrument List

Code: PSI-SET-A1180

Category: Endovascular Aneurysm Repair (EVAR)

Set Level: Basic

| SR | Specialty | Instrument Name | Size / Specification | Qty | Code |
|----|-------------------------|---------------------------|--------------------------|-----|------|
| 1 | Cardiovascular/Vascular | PRIME Scalpel Handle | No. 3 & 4 | 2 | |
| 2 | Cardiovascular/Vascular | PRIME Metzenbaum Scissors | 18 cm curved | 2 | |
| 3 | Cardiovascular/Vascular | PRIME Mayo Scissors | 17 cm straight | 1 | |
| 4 | Cardiovascular/Vascular | PRIME Adson Forceps | 12 cm fine toothed | 2 | |
| 5 | Cardiovascular/Vascular | PRIME DeBakey Forceps | 15 cm atraumatic | 2 | |
| 6 | Cardiovascular/Vascular | PRIME Mosquito Forceps | 12.5 cm curved | 6 | |
| 7 | Cardiovascular/Vascular | PRIME Needle Holder | Crile-Wood 15 cm | 2 | |
| 8 | Cardiovascular/Vascular | PRIME Bulldog Clamp | straight/curved assorted | 4 | |
| 9 | Cardiovascular/Vascular | PRIME Vessel Loop Passer | small | 1 | |
| 10 | Cardiovascular/Vascular | PRIME Suction Tip | Frazier 10 Fr | 1 | |
| 11 | Cardiovascular/Vascular | PRIME Retractor | Senn-Miller | 2 | |
| 12 | Cardiovascular/Vascular | PRIME Towel Clamp | Backhaus 9 cm | 4 | |