



# PRIME SURGICAL INTERNATIONAL

## UK OFFICE

24 King Street  
Manchester M2 4LG United Kingdom  
Phone: +44 20 7946 0123  
www.primesurgicalinternational.com

## SINGAPORE OFFICE

20 Cecil Street, #05-03 PLUS  
Singapore 049705  
+65 9076 3821  
Zia@primesurgicalinternational.com  
info@primesurgicalinternational.com

## PAKISTAN OFFICE

OFFICE # 206, 2ND FLOOR,  
QUEEN CENTER, M.T. KHAN ROAD,  
KARACHI, PAKISTAN  
+92 300 2621673  
ziarehman790@gmail.com

## Living Donor Surgery - Advanced Instrument List

**Code: PSI-SET-A609**

**Category: Living Donor Surgery**

**Set Level: Advanced**

| SR | Specialty  | Instrument Name                | Size / Specification | Qty | Code |
|----|------------|--------------------------------|----------------------|-----|------|
| 1  | Transplant | PRIME Scalpel Handle           | No. 3                | 2   |      |
| 2  | Transplant | PRIME Scalpel Handle           | No. 4                | 1   |      |
| 3  | Transplant | PRIME Mayo Scissors            | 14 cm straight       | 2   |      |
| 4  | Transplant | PRIME Metzenbaum Scissors      | 18 cm curved         | 2   |      |
| 5  | Transplant | PRIME Dissecting Forceps       | 14 cm toothed 1x2    | 2   |      |
| 6  | Transplant | PRIME Dissecting Forceps       | 14 cm plain          | 2   |      |
| 7  | Transplant | PRIME Mosquito Forceps         | 12.5 cm curved       | 6   |      |
| 8  | Transplant | PRIME Kelly Hemostatic Forceps | 14 cm curved         | 6   |      |
| 9  | Transplant | PRIME Allis Tissue Forceps     | 15 cm 4x5 teeth      | 4   |      |
| 10 | Transplant | PRIME Babcock Forceps          | 16 cm                | 2   |      |
| 11 | Transplant | PRIME Needle Holder            | Mayo-Hegar 16 cm     | 2   |      |
| 12 | Transplant | PRIME Towel Clamp              | Backhaus 9 cm        | 4   |      |
| 13 | Transplant | PRIME Retractor                | Langenbeck small     | 2   |      |
| 14 | Transplant | PRIME Sponge Holding Forceps   | 24 cm                | 1   |      |
| 15 | Transplant | PRIME Balfour Retractor        | complete             | 1   |      |
| 16 | Transplant | PRIME Kocher Forceps           | 18 cm                | 6   |      |
| 17 | Transplant | PRIME Right Angle Clamp        | Mixer 20 cm          | 2   |      |
| 18 | Transplant | PRIME Gall Duct Forceps        | Desjardins           | 1   |      |
| 19 | Transplant | PRIME Common Duct Dilator      | Bakes set            | 1   |      |
| 20 | Transplant | PRIME Intestinal Clamp         | non-crushing         | 4   |      |
| 21 | Transplant | PRIME Long Tissue Forceps      | 20-25 cm             | 2   |      |



# PRIME SURGICAL INTERNATIONAL

## UK OFFICE

24 King Street  
 Manchester M2 4LG United Kingdom  
 Phone: +44 20 7946 0123  
 www.primesurgicalinternational.com

## SINGAPORE OFFICE

20 Cecil Street, #05-03 PLUS  
 Singapore 049705  
 +65 9076 3821  
 Zia@primesurgicalinternational.com  
 info@primesurgicalinternational.com

## PAKISTAN OFFICE

OFFICE # 206, 2ND FLOOR,  
 QUEEN CENTER, M.T. KHAN ROAD,  
 KARACHI, PAKISTAN  
 +92 300 2621673  
 ziarehman790@gmail.com

## Living Donor Surgery - Advanced Instrument List

| SR | Specialty  | Instrument Name                             | Size / Specification             | Qty | Code |
|----|------------|---|----------------------------------|-----|------|
| 22 | Transplant | PRIME Malleable Retractor                   | 30 cm                            | 2   |      |
| 23 | Transplant | PRIME DeBakey Forceps                       | 20 cm atraumatic                 | 2   |      |
| 24 | Transplant | PRIME Potts-Smith Scissors                  | 18 cm angled                     | 2   |      |
| 25 | Transplant | PRIME Satinsky Vascular Clamp               | 24 cm                            | 2   |      |
| 26 | Transplant | PRIME Cooley Vascular Clamp                 | 20-24 cm                         | 2   |      |
| 27 | Transplant | PRIME Organ Retrieval Instrument Set        | complete                         | 1   |      |
| 28 | Transplant | PRIME Vascular Clamp Set                    | DeBakey/Cooley/Satinsky assorted | 1   |      |
| 29 | Transplant | PRIME Microvascular Instrument Set          | transplant                       | 1   |      |
| 30 | Transplant | PRIME Long DeBakey Forceps                  | 25-30 cm                         | 4   |      |
| 31 | Transplant | PRIME Vessel Dilator Set                    | Garrett/micro                    | 1   |      |
| 32 | Transplant | PRIME Perfusion Cannula Support Instruments | assorted                         | 1   |      |
| 33 | Transplant | PRIME Deep Abdominal Retractor System       | Bookwalter/Balfour               | 1   |      |
| 34 | Transplant | PRIME Bile Duct/ureteric Fine Instruments   | transplant                       | 1   |      |
| 35 | Transplant | PRIME Donor Organ Bench Set                 | advanced                         | 1   |      |
| 36 | Transplant | PRIME Atraumatic Bulldog Clamp Set          | assorted                         | 1   |      |